

Maestro Command Reference Manual

Version 9.3, September 2012

For inquiries about Maestro TM contact:

Schrödinger
101 SW Main Street, Suite 1300
Portland, OR 97204-3212
503-299-1150
503-299-4532 fax
Email: help@schrodinger.com

Copyright © 2012 Schrödinger, LLC. All rights reserved.

Maestro, CombiGlide, Epik, Glide, Impact, Jaguar, Liaison, Lig-Prep, Phase, Prime, PrimeX, QikProp, QikFit, QikSim, QSite, and Strike are trademarks of Schrödinger, LLC. MacroModel and Schrödinger are registered trademarks of Schrödinger, LLC.

To the maximum extent permitted by applicable law, this document is provided “as is” without warranty of any kind. This document may contain trademarks of other companies.

Please note that any third party programs (“Third Party Programs”) or third party Web sites (“Linked Sites”) referred to in this document may be subject to third party license agreements and fees. Schrödinger, LLC and its affiliates have no responsibility or liability, directly or indirectly, for the Third Party Programs or for the Linked Sites or for any damage or loss alleged to be caused by or in connection with use of or reliance thereon. Any warranties that we make regarding our own products and services do not apply to the Third Party Programs or Linked Sites, or to the interaction between, or interoperability of, our products and services and the Third Party Programs. Referrals and links to Third Party Programs and Linked Sites do not constitute an endorsement of such Third Party Programs or Linked Sites.

1 The Maestro Command Language

This manual contains a listing of all supported Maestro commands, their functions, and their usage options. Use Maestro commands by entering them in the Command Input Area of the Main Application Window.

The Maestro command language syntax is:

keyword [*options*] [operand]

For example:

entryimport *wsreplace=true all=true format=maestro ligands.mae*

In the example given, “entryimport” is the keyword. As this example illustrates, the first item entered must be a keyword, and this keyword must match a known keyword completely, i.e., to all characters. However, Maestro does have an alias facility, which allows users to map keywords to shorter names, for example, “entryimport” to simply “import”. See the “alias” command for details.

Command options are used to set values that generally correspond to states within the program. As the name “options” suggests, inclusion of these items in a Maestro command is not mandatory. However, any number of option *name*=⟨value⟩ pairs can be specified. Options can be of type: Boolean (yes/no, y/n, t/f, and true/false are valid values), string, integer, or real.

Once an option is set, its value will persist for all subsequent commands until it is specified again. A command consisting only of a keyword and options will set the states corresponding to those option values, but will not perform any other action. For example, the command:

entryimport *format=maestro all=true*

will not import structures from any file. It will only set the format for subsequent file reading and specify that all structures contained in the subsequently specified file be imported. If the following command is then issued:

entryimport *ligands.mae*

all structures from the *ligands.mae* Maestro file will be imported.

The operand is the remaining element of a Maestro command. The type of operand required is dependent on the specific command. For “entryimport”, the operand must be a file name. For other commands, it may be one or more atom numbers, or a string describing arbitrary atom sets, specified using Maestro’s Atom Specification Language (ASL). For more information about ASL, see [Chapter 3 \[ASL\]](#), page 5.

2 Conventions Used in this Manual

This manual uses a number of typographical conventions to describe the command language.

The names of keywords are printed in boldface:

entryimport

The names of options are printed in italic, and literal choices for that option are in normal roman:

wsreplace=true

Values which the user is to replace with appropriate values are represented in angle brackets:

entryimport \langle file_name \rangle

Where there is a choice of a number of values, these are shown separated by vertical bars:

entryimport *format*=maestro|mmod|sd|pdb|mol2

If a value is optional then it is enclosed in square brackets:

print [\langle file_name \rangle]

The following conventions are used to describe acceptable values for options:

\langle n \rangle means an integer

\langle x \rangle means a real number

yes|no means a boolean value (could also be true|false, y|n, on|off)

\langle text \rangle represents a string value. If this contains spaces, it must be enclosed in double quotation marks

The following are some conventions used for operands:

\langle atom_number \rangle
represents an atom number

\langle atom1 \rangle represents the first atom number in the operations, \langle atom2 \rangle the second, and so on

\langle ASL \rangle represents a valid string in the Atom Specification Language

\langle ESL \rangle represents a valid string in the Entry Specification Language

Chapter 2: Conventions Used in this Manual

3 The Atom Specification Language

This chapter documents the syntax of Maestro’s Atom Specification Language (ASL).

3.1 Why an Atom Specification Language?

- To provide a flexible way to define sets of atoms in complex macromolecular systems
- To serve as the basis of a “Sets” facility in Maestro
- To allow atom specification from textual input which some users find faster than picking atoms from the main structure window.

3.2 The ASL Hierarchy

`entry > molecule > chain > residue > atom`

There are five classes which make up the atom specification language. Each is listed below.

In this section the minimum acceptable abbreviation is shown outside the square brackets. So for example `m[olecule]` means that the minimal acceptable abbreviation for `molecule` is `m`.

`e[ntry]` This is the top level class in the language. An entry is all atoms in the workspace associated with a single entry in the currently open project.

`m[olecule]` The term molecule is used in the normal chemical sense meaning all atoms which are connected by a single covalent path.

`c[hain]` This corresponds to a chain as specified in the PDB file format. Note that this chain may be a subset of a molecule, e.g., when chains are linked by disulphide linkages.

`r[esidue]` An arbitrary collection of one or more covalently bound atoms within a molecule, such as the monomer units in a polymer.

`a[tom]` A single atom.

Each class is optional. If absent all entities of that type are matched.

3.3 Atom Specification

A complete specification is a class name and some property specified by a property name and property list. The syntax is:

class.property ⟨propertylist⟩

Items in a property list may be separated by comma, whitespace or both. Ranges (lower-upper) may be used where appropriate. Unterminated ranges are taken to include all available numbers. For example, if there are four molecules in the system then the specifications:

- `mol. 2, 3, 4`
- `mol. >=2`
- `mol. 2-4`
- `mol. >1`

are equivalent.

In a similar manner,

- `mol. 1, 2, 3`
- `mol. <=3`
- `mol. 1-3`
- `mol. <4`

are equivalent.

All names of properties and characters in property lists are treated in a case insensitive manner.

Wildcards are supported for atom and set names. A ‘*’ will match zero or more characters and a ‘?’ will match any single character. You can include comments in a specification by placing a ‘#’ character before the text you wish to hide.

In this section the minimum acceptable abbreviation is shown outside the square brackets. So, for example, `a[tom].pt[type]` means that the minimal acceptable abbreviation for `atom.ptype` is `a.pt`.

`e[ntry]`.

[name]

Because entry names are the only entry properties specifiable using ASL expressions, the word **name** can be completely omitted, i.e., `entry.name ename` and `entry. ename` are equivalent. Note, however, that the ‘.’ is still required even if the property name is not included. A valid property list for the ‘entry name’ property is a list of entry names. Wildcard characters are permitted. For example:

```
entry. e1
entry.name recep, lig*
```


- `m[olecule]`.
 `[number]` Because numbers are the only molecular properties specifiable using ASL expressions, the word **number** can be completely omitted, i.e., `mol.number 1` and `mol. 1` are equivalent. Note, however, that the `'.'` is still required even if the property name is not included. A valid property list for the 'molecule number' property is a set of numbers or a range. For example:
 `mol. 1-4`
 `mol. 1,2,3,4`
- `m[odulo]` a property which can be used to select every *n*th molecule. For example: `mol.mod 10 1` will select molecule 1, 11, 21, etc.
- `e[ntrynum]`
 a property which can be used to select molecules based on their entry-relative numbering. For example: `mol.entrynum 1` will select the first molecule in each entry.
- `a[toms]` a property which can be used to select molecules based on the number of atoms they contain. For example: `mol.atoms 200` will select molecules that contain exactly 300 atoms. Other examples are:
 `mol.atoms 200-300`
 `mol.a > 200`
- `w[eight]` a property which can be used to select molecules based on their molecular weight. For example: `mol.weight 200.12` will select the atoms that have a molecular weight of exactly 200.12. Other examples are:
 `mol.weight <=300.0`
 `mol.weight > 200.0`
- `c[hain]` This class designation allows you to specify atoms using chain attributes. Combine with the 'name' property.
- `[name]` Because names are the only chain properties that are specifiable by ASL, the word **name** can be omitted. For example, the specifications `chain.name A` and `chain. A` are equivalent. Note, however, that the `'.'` is required even if the property name is not included. A valid property list for the 'chain name' property is a single character representing a PDB

chain name. Some examples of equivalent acceptable 'chain name' expressions are:

```
chain.name A
chain. A
c. A
```

r[esidue]

This ASL class designation allows you to specify atoms based on residue properties. Combine with one of the following property specifications.

[name] or **[number]**

Either residue names or numbers (but not both) can be used in property lists which do not specify a property name. For example, the following are valid ASL expressions:

```
res. ala val leu
res. 1 2 3
```

and will return atoms which are either in alanine, valine, or leucine residues, or atoms in residues 1, 2 or 3, respectively. Residue numbers can be specified with a range, e.g., **residue. 1-4** and may include negative values or zero.

pt[type]

The three-letter PDB code for the residue. This is the default for non-numeric characters in the property list, so the expression **res. arg** and **res.ptype arg** are equivalent. A valid property list for 'ptype' is comprised of three-character tokens. For example:

```
res.ptype gly, val, ala
res. gly val ala
```

m[type]

The one letter residue codes as used in Maestro. A valid property list for 'mtype' is comprised of one-character tokens. For example:

```
res.mtype g, v, a
res.m g, v, a
```

po[larity]

The polarity of the residue. The property list must consist only of the following descriptor types:

h[ydrophobic]

returns atoms in hydrophobic residues

pol[ar]

returns atoms in polar residues

`pos[itive]`
returns atoms in residues with positive formal charges

`n[egative]`
returns atoms in residues with negative formal charges

For example,

```
residue.polarity hydrophobic
residue.pol pos,neg
res.pol h pos neg
```

`sec[ondary_structure]`

The secondary structure of the residue. The property list must consist only of the following descriptor types:

`h[helix]` returns atoms in helical regions

`s[trand]` returns atoms in strand regions

`l[oop]` returns atoms in loop regions

For example,

```
residue.sec helix
residue.sec hel, str
res.sec l, s
```

`pos[ition]`

The fractional position of the residue. The property list must include two real numbers representing a fractional range of residue numbers. For example, if there are 100 residues numbered from 1-100, the specification:

```
residue.pos 0.0 0.1
```

will return residues 1 to 10.

`i[nscore]`

The insertion code of the residue. A property list should include one-character tokens representing insertion codes. For example:

```
residue.inscode a
```

will get all residues with insertion code 'a', while:

```
res. 25 and res.inscode b
```

will get residue 25b.

a[atom] Atom ptypes and numbers may be mixed in property lists where no explicit property is specified. For example, the following is valid:

```
atom. 1,2,3,CA
```

and returns atoms 1,2 and 3 and any alpha carbons.

pt[type] The PDB atom names. A valid property list for 'ptype' consists of acceptable PDB names. This is the default property for non-numeric components of property lists and, as such, the word 'ptype' may be omitted. Note that the '.' is required. The following specifications are equivalent and return the "backbone" atoms in a structure:

```
atom.ptype N,CA,C,O
```

```
atom. N,CA,C,O
```

```
a. n,ca,c,o
```

Note: see below for a discussion of how PDB atom names are specified and matched.

Wildcards as described above can be applied to ptypes.

na[me] The atom names. The property list must contain valid atom names. A valid atom name could be a string of any length which:

1. must have at least one non-digit character
2. must not contain any control characters (ASCII value < 0x20)
3. must not contain spaces or equal signs, unless the name is quoted

Wildcards as described above can be applied to atom names. Examples:

```
atom.name the_36th_carbon
```

```
atom.na C15, O:66, H-77
```

```
atom.na C* (returns atoms with name  
starting with C)
```

```
atom.nam ??0* (returns atoms whose  
name's 3rd character is '0')
```

n[umber] The atom numbering. The property list must be a list or a range of numbers. This is the default property for numeric components of property lists. The following expressions are equivalent and return the atoms numbered 1, 2, 3, and 4.

```
atom.num 1,2,3,4
a. 1 2 3 4
atom. 1-4
```

mo[lnum] The atom numbering in the “by molecule” scheme. The property list must be a list or a range of numbers. For example, the expression:

```
atom.molnum 1
```

returns the first atom in each molecule, while the specification:

```
atom.molnum 1-10
```

returns the first 10 atoms in each molecule.

en[trynum] The atom numbering in the “by entry” scheme. The property list must be a list or a range of numbers. For example, the expression:

```
atom.entrynum 1
```

returns the first atom in each entry, while the specification:

```
atom.entrynum 1-10
```

returns the first 10 atoms in each entry.

m[type] The Maestro atom type. A valid property list for the ‘mtype’ property consists of Maestro atom types. The following expression is valid and specifies sp2 carbons and oxygens.

```
atom.mtype C2,O2
```

e[lement] The element symbol for the atom. A valid property list for the element is a list of standard periodic table symbols. To define all carbons and oxygens:

```
atom.ele C,O
```

att[achments] The number of bonds the atom has to it. The property list must be a number in the range 0-6, but greater than ‘>’, less than ‘<’, and equals ‘=’ signs may also be used. The expression:

```
atom.att 1
```

returns all terminal atoms. The specification:

```
atom.att <=2
```

returns all terminal atoms and other atoms with 2 or fewer bonds. The specification:

`atom.att 0`

returns all isolated atoms.

`ato[micnumber]`

The atomic number. The property list must contain numbers only. Ranges of integers and greater than '>', less than '<', and equals '=' signs may also be used. The expression:

`atom.atomicnum 1`

returns all hydrogen atoms. The specification:

`atom.ato 1-6`

returns all atoms in the range H to C.

`c[harge]` The partial charge on the atom. A valid property list contains a value or range of floating point values. The expression:

`atom.charge 0.400`

returns atoms with partial charges of 0.400. The expression:

`atom.charge -0.6--0.4`

returns atoms with partial charges -0.6 to -0.4,

`atom.charge <0.0`

returns atoms with negative partial charges, and

`atom.charge >=0.5`

returns atoms with charges of 0.5 or greater.

`f[ormalcharge]`

The formal charge on the atom. A valid property list contains a value or range of integer values. The expression:

`atom.formal 0`

returns atoms with formal charges of 0. The expression:

`atom.formal -2 -- 1`

returns atoms with formal charges -2 to -1,

`atom.formal <0`

returns atoms with negative formal charges, and

`atom.formal >=1`

returns atoms with formal charges of 1 or greater.

d[isplayed]

Whether or not the atom is currently displayed in Maestro. No property list is used. For example, the expressions:

`atom.displayed` returns the set of all displayed atoms
`not atom.displayed` returns all atoms not currently displayed
`sidechain and fillres atom.displayed` returns the side chains of residues where at least one atom is displayed.

s[elect]

Whether or not the atom is currently selected in Maestro. No property list is used. For example, the expression:

`atom.selected` returns the set of all atoms selected in the Workspace

Using generalized atom properties. Some structures may have additional properties available. These are referenced directly by their data names appended onto the "atom." class. These properties are either of integer, real, boolean or string type and the datanames are encoded as beginning with `i_`, `r_`, `b_` or `s_` respectively. It is possible to use these atom properties in conjunction with any other ASL expression. Any atoms which don't have these properties associated with them will never match. Some examples of using the ASL to address these properties are:

`atom.i_my_integer_prop 1-4` `atom.b_my_boolean_prop`
`atom.r_my_real_prop < 4.0` `atom.s_my_string_prop LIG_`

3.4 Operators

A number of operators are supported:

3.4.1 The Boolean and Operator

Boolean AND (set intersection). The syntax for this operation is:

`spec1 and spec2`

where `spec1` and `spec2` are valid atom specifications. This operation will return the set of atoms which meets the specifications `spec1` and `spec2`. For example, the expression:

`mol. 1 and atom. CA`

will return the set of all the alpha carbons of molecule 1. The specification:

```
res.num 1-100 and res. ala
```

returns all alanines in residues with numbers in the range 1-100.

3.4.2 The Boolean or Operator

The Boolean OR operator. The syntax for this operations is:

```
spec1 or spec2
```

where **spec1** and **spec2** are valid atom specifications. This returns the set of atoms which meet either specification **spec1** or **spec2**. For example, the expression:

```
mol. 1 or atom.ptype CA
```

returns the set of all atoms that are in molecule number 1, or are alpha carbons. The specification:

```
res.num 1-100 or res.ptype ala
```

returns all residues either with numbers in the range 1-100, or any alanines.

3.4.3 The Boolean not Operator

The Boolean NOT operator. The syntax for this operation is:

```
not spec1
```

where **spec1** is a valid atom specification. This returns the set of atoms that are not part of those defined by **spec1**. For example, the expression:

```
not atom. CA,C,N,O
```

will return a set containing all side chain atoms.

3.4.4 The fillres and fillmol Operators

Two special operations, **fillres** and **fillmol**, can be used to “fill out” the atoms defined by a atom specification to complete residue or molecule boundaries. For example:

```
fillres atom.num 1,100,40
```

will return all the atoms in residues of which atoms 1,100 and 40 are members. In a similar way:

```
fillmol atom.num 1,100,40
```

will return all the atoms in molecules of which atoms 1,100 and 40 are members.

3.4.5 The `within` and `beyond` Operators

The operators `within` and `beyond` can be used to define sets of atoms based on their distance from a set of atoms defined by an atom specification. The syntax for these operators is:

`within distance spec1`

which returns the atoms within, i.e., less than or equal to the distance in Angstroms of the set defined by `spec1`, and

`beyond distance spec2`

which returns the atoms which are further than the distance in Angstroms from the set defined by `spec2`.

For example, the expression:

`within 5.0 mol. 1`

returns the set of all atoms that are within 5 Å of molecule 1. The expression:

`beyond 5.0 mol. 2`

returns all atoms that are farther than 5 Å from molecule 2.

The combination of `fillres` and `within` or `beyond` is especially powerful.

`fillres within 5.0 mol. 1`

will produce a set containing the atoms of all complete residues that have atoms within 5 Å of molecule 1. Note that the `within` operator will also return the reference set of atoms:

`within 5.0 mol. 1`

returns the reference set of all atoms that are within 5 Å of molecule 1 and those that are part of molecule 1.

The `and` operator, when used with `within` and `beyond`, can be used to allow more specificity:

`mol. 2 and within 5.0 mol. 1`

returns the set of all atoms of molecule 2 that are within 5 Å of molecule 1.

3.4.6 The `withinbonds` Operator

A special type of `within`, one that finds all atoms within a certain number of bonds of the reference set:

`withinbonds <num_bonds> <spec>`

For example:

`withinbonds 4 atom. 1`

will find all the atoms that are within four bonds of atom 1.

3.4.7 The `beyondbonds` Operator

A special type of `within`, one that finds all atoms beyond a certain number of bonds of the reference set:

```
beyondbonds <num_bonds> <spec>
```

For example:

```
beyondbonds 4 atom. 1
```

will find all the atoms that are in the same molecule as atom 1 but beyond four bonds of atom 1.

3.5 Operator Priority

The order of priority of operators is (in decreasing order):

- `not/fillres/fillmol`
- `and/or`
- `within/beyond`

At equal levels of priority the expression will be evaluated left to right.

Examples:

```
within 5.0 mol. 1 or mol. 2
```

returns the set of all atoms that are within 5.0 Å of either molecule 1 or molecule 2 (or has higher priority).

```
not atom.ptype CA,C,O,N and mol. 1
```

returns the side chain atoms of molecule 1.

```
atom.ptype CA or mol. 1 and not res.pol polar
```

returns all alpha carbons and atoms in hydrophobic residues of molecule 1.

Parentheses can be used to override the order of evaluation:

```
not (atom.ptype CA,C,O,N or mol. 1)
```

produces all atoms either not in the backbone or not in molecule 1.

3.6 Implicit Operators

When no operator is specified, the following operations are assumed:

- Specifications within a property list have a Boolean `or` relationship. So the following specification:

```
atom.ptype CA,CB
```

matches any atom which has a PDB name of either CA or CB.

- When the classes are in order of decreasing priority in the class hierarchy, the `and` operator is assumed:

```
mol. 1 chain. A atom.ptype CA
```

is equivalent to:

```
mol. 1 and chain. A and atom.ptype CA
```

and both return the alpha carbons in chain A of molecule 1.

3.7 Creating New Sets from Existing Ones

Maestro supports a mechanism for defining and naming atom sets via its “Sets” panel. The names of existing sets may be used in expressions if they are prefixed with the word **set**. For example, if there are two sets defined as:

```
set S1 mol. 1 (molecule number 1)
set S2 atom.ptype C,O,N,CA (all backbone atoms)
```

the following are then valid atom specifications:

- **set S1 and set S2**, all backbone atoms in molecule 1
- **set S1 or set S2**, all backbone atoms or atoms in molecule 1
- **within 5.0 set S1**, all atoms within 5.0 Å of molecule 1.

3.8 Special Specifications

- The specification **all** matches everything. You can use this in any Maestro command which expects an ASL operand in order to apply the command to all the atoms. For example to color all atoms green:

```
coloratom color=green all
```

Note that the syntax for this is just **all**, not **atom. all** or **molecule. all**.

- Anything enclosed in **/ /** is treated as a string in the linear substructure notation - see the BatchMin Reference Manual for a complete description of this. Some examples are:

```
/C3(-H1)(-H1)(-H1)/ specifies all methyl groups
/C2(=O2)-N2-H2/ specifies all amide groups
/C0-N0/ specifies all pairs bound C-N atoms, and
/O0-S0/ specifies any atom bound to a sulfur.
```

Specifications made in this way can be treated in the same way as any other specification in the ASL, e.g. they can have operators applied to them.

3.9 Matching PDB Atom Names

Because PDB atom names are four characters wide we need to employ the following strategy in order to conveniently specify PDB atom names:

- Unquoted names which begin with a non-numeric character have a blank character inserted in front of them and are padded with blanks from the right to make up four characters before matching. Examples:
 - Property List: CA,C
Actually Matched: " CA "," C "
 - Property List: CG1,CG2
Actually Matched: " CG1 "," CG2"
- Unquoted names which begin with a number do not have an initial blank character inserted but are right to padded four characters. Example:
 - Property List: CA,1HB,2HB
Actually Matched: " CA ","1HB ","2HB "
- Quoted names (double or single quotes are acceptable) are treated as literals but will be right padded with blanks to make up four characters. Examples:
 - Property List: CA,"CA"
Actually Matched: " CA ","CA " (alpha carbons or calciums)
 - Property List: " N A"
Actually Matched: " N A" (heme atoms)

3.10 Miscellaneous

3.10.1 Atoms not yet present in a structure

If a structure in the Workspace has only 100 atoms when the ASL definition: `atom.num 1,8,44,101,103` is issued, Maestro will simply match the atoms numbered 1,8 and 44. If additional atoms are subsequently added to the structure, the atoms bearing the numbers 101 and 103 will be added to the previously defined set.

3.10.2 Aliasing

Maestro allows you to define your own aliases, using either the Command Input Area or the Command Aliases panel. Maestro converts all aliases into their corresponding commands before performing operations involving the aliased commands. Users must ensure that aliases produce sensible results. Some aliases are supplied with the distribution. They are:

Operator: **and**
Aliases: intersection, INTERSECTION, &

Operator: **or**
Aliases: UNION, union, |

Operator: **not**
Aliases: !

Class Designator: **mol.**
Aliases: MOL, mol

Class Designator: **atom.**
Aliases: ATOM, atom

Class Designator: **res.**
Aliases: RES, res

Class Designator: **chain.**
Aliases: CHAIN, chain

ASL Definition: **atom. ca,c,n,h,o**
Aliases: BACKBONE, backbone

ASL Definition: **not (atom.pt ca,c,n,h,o)**
Aliases: SIDECHAIN, sidechain

ASL Definition: **"/H2-O3-H2/ or atom.mtype OW"**
Aliases: WATER, water

3.11 Useful Hints when using ASL with the Project Facility

Maestro 4.1 introduced its Project Facility. Entries can be included into and excluded from the Workspace. The order in which this is performed affects the molecule numbers. For example, if you have two entries in your Project Table called "A" and "B" and you include into an empty Workspace first A and then B, the molecule numbers will be 1 for A and 2 for B; however, if you first include B and then A, the molecule numbers will be 1 for B and 2 for A.

This means that an expression such as **mol. 1** will match different atoms in each of the above cases. In earlier versions of Maestro, it often made sense to use the **mol.** notation because there was no Project Facility. But with the Project Facility, Maestro is entry-centric and in most cases it makes more sense to use entry names.

For example, if you have an inhibitor and a receptor that are in different entries and wish to have a ribbon appear on only the receptor, use the entry

name in the ASL expression, not the molecule number. This will ensure that when the receptor is included that it, and only it, will be used to generate the ribbon(s). A different inclusion order of entries in the Workspace will then result in the same matching atoms. So for ribbons with a receptor called "receptor" it is more useful to use `entry.name receptor` as the ASL definition.

3.12 ASL Examples

This section gives some examples of the use of the ASL in real-life situations. Note that while these examples all use lower-case, the ASL expressions themselves are not case sensitive.

i) Defining a set to refer to a ligand and/or receptor.

The exact command will depend on the nature of your system. If the ligand and the receptor are separate entries then it will suffice to use

```
set ligand entry.name <ligand_name>
```

where <ligand_name> is the name of the entry that contains the ligand. Similarly

```
set receptor entry.name <receptor_name>
```

for the receptor with entry called <receptor_name>.

In order to define sets that will work with multiple ligands it's also possible to define the ligand as everything that is not part of the receptor. A definition of:

```
set ligand not set receptor
```

will identify the ligand as anything that's not part of the receptor.

If the ligand and the receptor are part of the same entry then molecule numbers are the best way to define the ligand and the receptor. Assuming the receptor is molecule 1 and the ligand molecule 2:

```
set ligand mol.num 2
```

```
set receptor mol.num 1
```

Note however that the use of molecule numbers in set definitions should be avoided where possible as these depend on the order in which the project entries are included into the Workspace. If it is possible to use entry names, then these should be used.

The subsequent examples assume that sets for the receptor and the ligand have been defined using one of the methods defined above.

ii) The set of atoms within a given distance of the ligand.

One common task is to do something with the set of atoms within a given distance of the ligand. For example to only display those atoms or to include them in a substructure region for a MacroModel calculation. These examples

will use the "displayonlyatom" command but the ASL which follows can be used with any other command that uses ASL.

To only display atoms within 5.0 Angstroms of the ligand:

```
displayonlyatom within 5.0 set ligand
```

A common variation is to display complete residues which have any of their atoms within a given distance of the ligand:

```
displayonlyatom fillres within 5.0 set ligand
```

It's also possible to restrict the expression so that it only applies to receptor atoms within a given distance of the ligand. Here the Boolean "and" operator is used to restrict the displayed atoms to the receptor only:

```
displayonlyatom set receptor and fillres within 5.0 set  
ligand
```

Because this is a lengthy expression it's often convenient to make this into a set itself:

```
set active_site set receptor and fillres within 5.0 set  
ligand
```

An equivalent form of this is:

```
set active_site (! set ligand) & fillres within 5.0 set  
ligand
```

Note that "!" is a standard alias for "not" and similarly "&" for "and".

iii) Sidechain and backbone.

The ASL has standard aliases for the definition of sidechain and backbone atoms in proteins. For example to only display the atoms of the backbone:

```
displayonlyatom backbone
```

These aliases can be used with operators to build up more complicated expressions. For example to only display the sidechain of the receptor:

```
displayonlyatom sidechain and set receptor
```

To display only the sidechains of the atoms within 5.0 Angstroms of the ligand:

```
displayonlyatom sidechain and set receptor and fillres  
within 5.0 set ligand
```

iv) Atoms of a given type.

There are a variety of ways to specify atoms of a given type. For example to specify all carbons, nitrogens and oxygens the following is used:

```
atom.ele C,N,O
```

To specify non-hydrogen atoms:

```
not atom.ele H
```

To specify the alpha carbons in a protein:

```
atom.ptype CA
```

Chapter 3: The Atom Specification Language

To specify all sp² carbons there are two choices. The first relies on knowing that the MacroModel atom type for such an atom is "C2" and using:

```
atom.mtype C2
```

The other (assuming no formally charged or radical carbons are present) uses the number of attachments to the atom:

```
atom.ele C and atom.att 3
```

To specify polar hydrogens:

```
atom.ele H and not /C0-H0/
```

or

```
atom.ele H and not atom.mtype H1
```

v) Water molecules.

The ASL has a standard alias "water". For example to delete all water molecules the Maestro command is:

```
delete atom water
```

vii) Restricting an operation to the atoms that are currently displayed in the workspace.

Often a user will be working with only a subset of the atoms in the workspace displayed. If an operation is to be performed only on the atoms that are displayed then the "atom.displayed" property can be used. For example to change the color to green of all the atoms currently displayed in the Workspace and to leave alone the undisplayed Workspace atoms:

```
coloratom color=green atom.displayed
```

To only do it for the atoms that are displayed and in the receptor:

```
coloratom color=green atom.displayed and receptor
```

viii) Specifying molecules.

All molecules with between 30 and 100 atoms:

```
mol.atoms 30-100
```

All molecules with over 100 atoms:

```
mol.atoms >100
```

All molecules with a molecular weight over 300:

```
mol.weight > 300.0
```

All molecules which contain a halogen:

```
fillmol atom.ele F,Cl,Br,I
```

ix) Specifying atoms based on a linear-substructure notation.

The ASL supports the use of a SMILES-like linear substructure notation to specify atoms with a particular bonding arrangement. The atoms are referred to by MacroModel atom types, but there are wildcard types that can be used to allow the expression to apply to any atoms of a given element type.

Some examples:

Any five-membered ring:

```
/00-00-00-00-00-1/
```

Aromatic six-membered carbon rings(C2 is sp2 carbon)

```
/C2-C2*C2-C2*C2-C2*1/
```

Amide groups:

```
/C2(*O2)-N2/
```

Methyl groups:

```
/C3(-H1)(-H1)(-H1)/
```

Water:

```
/H2-O3-H2/
```

Guanadinium group:

```
/N2(-H3)-C2(*N4(-H4)(-H4))-N2(-H3)(-H3)/
```

x) Using wildcard characters.

Most string-type property values can use wildcard characters. Some examples:

All PDB atom names beginning with C

```
atom.ptype C*
```

All forms of the histidine residue:

```
res.ptype HI*
```

All entries that begin with "lig":

```
entry.name lig*
```

xi) Addressing residues with negative residue numbers.

From version 6.1 of Maestro it is possible to use residue numbers that are negative or zero.

```
res. -8
```

will match residues with residue number -8.

xii) Use SMARTS expressions.

From version 7.0 of Maestro it is possible to use SMARTS expressions as part of an ASL expression.

Some examples:

```
smarts. CCC
```

will match all three-carbon subsequences

```
smarts. [R] and atom.ele N
```

will match all ring nitrogens

```
smarts. C1CCCCC1
```

will match all six-membered carbon rings

4 The Entry Specification Language

This chapter documents the syntax of Maestro's Entry Specification Language (ESL).

4.1 Why an Entry Selection Language?

The Entry Selection Language (ESL) is used for selecting entries in a Maestro project based on the properties of those entries. ESL expressions can be of arbitrary complexity and can use parentheses and logical operators. While the ESL is similar in appearance to the Atom Specification Language (ASL), they are used for quite different purposes - the ASL for specifying sets of atoms in the workspace, the ESL for selecting entries in the project. A typical application for the ESL would be for filtering. An ESL expression can be used to select only those entries that meet a particular criteria, usually based on entry properties. Once selected, those entries can be displayed in the workspace or exported to an external file.

4.2 Entry Properties

There are a number of properties that can be used in ESL expressions:

entry The name of the entry is the most useful property at present. The syntax is:

```
entry <name>
entry_re <name-expression>
```

Entry names that contain spaces need to be enclosed in double quotes. The following wildcard characters can be used for the **entry_re** expression:

- * matches any number of any characters
- ? matches any single character
- # matches any single digit

Examples:

```
entry anentry
entry crambin-1
entry_re cramb*
entry_re ligand##
```

Note: **name** can be used as an alias for **entry** and **name_re** can be used as an alias for **entry_re**.

selected This allows the selection of entries based on whether or not they are already selected in the project. Examples:

	<code>entry_re lig* and not selected</code> <code>not selected</code>
<code>included</code>	This allows the selection of entries based on whether or not they are already included for display in the workspace. <code>entry_re lig* and not included</code> <code>not included</code>
<code>all</code>	This is a quick way to select all entries.

4.3 Logical Operators

The following logical operators are supported (in order of decreasing priority)

- `not`
- `and` `or`

At equal levels of priority expression will be evaluated from left to right. Parentheses can be used to override the default order of priority.

4.4 Entry Property Comparisons

ESL expressions may also contain comparisons involving entry properties. Properties may be defined by their data names or their user names. The data names are the names that appear in the Maestro files. An example is “`i_qp_n_stars`”. User names are the names that appear in the first row of Maestro’s Project Table. An example of how this would appear in Maestro’s Project Table Header is `n_stars`. User names are specified in commands as:

```
user(<propertyname>)
```

For example:

```
user(volume)
```

and so,

```
user(n_stars) and  
i_qp_n_stars
```

are equivalent names for properties.

For integer and real type properties the following comparison operators are supported: `==`, `!=`, `<`, `>`, `<=`, `>=`.

For boolean (logical) properties, the name of the property can just be specified and will match any entries where that property is true. It is also possible to use `==` and `!=` in explicit comparisons with the values `true` and `false`.

Comparisons involving string properties have the form:

`<property name> <string literal>`

For example:

```
s_m_string_prop left
s_m_string_prop left*
```

String literals may include the *, ? and # wildcard characters available for entry names as described above.

4.5 Examples

Some examples of the use of the ESL are:

```
entry anentry
entry_re ligand*
entry_re lig* and selected
(entry_re lig# and selected ) and not included
user(volume) < 100.0 and selected
included and user(n stars) == 2 and user(n amine) == 0
user(log s) > 0.5 and user(log S) < 1.0
```


5 Commands

1ddataset

Controls the appearance of the a single dataset on a plot of a dihedral drive of a single dihedral.

Syntax:

```
1ddataset ccolor=black | red | green | blue | purple | orange |
          blue_green | light_green | red_purple | yellow | cyan
          curve=solid | dashed | noline cwidth=<n> scolor=black | red
          | green | blue | purple | orange | blue_green | light_green |
          red_purple | yellow | cyan ssize=<n> symbol=filled_rectangle
          | hollow_rectangle | filled_circle | hollow_circle | cross | point
          | filled_diamond | hollow_diamond | no_symbol
          <data_set_name> [<grd_file_name>]
```

Options:

<i>ccolor</i>	The color of the curve used on the graph for this dataset.
Valid values:	black red green blue purple orange blue_green light_green red_purple yellow cyan
Default value:	black
<i>curve</i>	The type of curve which will join the points on the graph.
Valid values:	solid dashed noline
Default value:	solid
<i>cwidth</i>	The width of the line used on the graph for this dataset.
Valid values:	integers
Default value:	1

Chapter 5: Commands

	Minimum:	1
	Maximum:	5
<i>scolor</i>	The color of the symbols used on the graph for this dataset.	
	Valid values:	black red green blue purple orange blue_green light_green red_purple yellow cyan
	Default value:	black
<i>ssize</i>	The size of the symbol used on this plot for this dataset.	
	Valid values:	integers
	Default value:	3
	Minimum:	2
	Maximum:	5
<i>symbol</i>	The type of symbol used for this data set on the graph.	
	Valid values:	filled_rectangle hollow_rectangle filled_circle hollow_circle cross point filled_diamond hollow_diamond no_symbol
	Default value:	filled_rectangle

Operands:

$\langle \text{data_set_name} \rangle$ [$\langle \text{grd_file_name} \rangle$]

The first operand is the name of the dataset which is to be created or to have its properties changed. If a second operand is present then that is treated as the name of a file from which the data is to be read. The full name of the file, including the suffix, should be included.

1dplot

Controls the appearance of the plot the energy as a function of a drive of a single dihedral angle.

Syntax:

```
1dplot amax=<x> amin=<x> emax=<x> emin=<x>
      escale=absolute | relative orientation=landscape | portrait
      papersize=letter | a4 showlegend=yes | no
      squareplot=yes | no title=<text> units=kj | kcal
      xlabel=<text> xdecimal=<n> ylabel=<text>
      ydecimal=<n> [<postscript_file_name>]
```

Options:

<i>amax</i>	The maximum angle value displayed on the graph Valid values: reals Default value: 0
<i>amin</i>	The minimum angle value displayed on the graph Valid values: reals Default value: 0
<i>emax</i>	The maximum energy value displayed on the graph Valid values: reals Default value: 0
<i>emin</i>	The minimum energy value displayed on the graph Valid values: reals Default value: 0
<i>escale</i>	The energy scale (absolute/relative) for the graph Valid values: absolute relative Default value: absolute
<i>orientation</i>	The paper orientation for the postscript output Valid values: landscape portrait Default value: landscape
<i>papersize</i>	The paper size for the postscript output Valid values: letter a4 Default value: letter

Chapter 5: Commands

showlegend

Whether the legend is displayed for the 1D plot

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

squareplot A boolean which controls whether the plot will be constrained to be square.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

title The plot title

Valid values: text strings

Default value:

units The energy units for the graph

Valid values: kj

 kcal

Default value: **kj**

xaxislabel The label for the X-axis.

Valid values: text strings

Default value: **Coordinate**

xdecimal The number of decimal places for the X axis.

Valid values: integers

Default value: **0**

Minimum: 0

Maximum: 5

yaxislabel The label for the Y-axis.

Valid values: text strings

Default value: **Energy**

ydecimal The number of decimal places for the Y axis.

Valid values: integers

Default value: **1**

Minimum: 0

Maximum: 5

Operands:

[<postscript_file_name>]

If an operand is given then this will be treated as the name of a file to which a postscript representation of the plot will be written. The full name of the file, including any suffix, must be included.

1drescale

Rescales the 1D plot axes back to be able to view the entire data range.

Syntax:

1drescale

1dtable

Write out a table of the data from the 1D plot to an external file.

Syntax:

1dtable

2ddataset

Controls the appearance of the dataset plotted on the 2D grid contour plot.

Syntax:

2ddataset *colormap*= $\langle \text{text} \rangle$ *cwidth*= $\langle \text{n} \rangle$ *emax*= $\langle \text{x} \rangle$ *emin*= $\langle \text{x} \rangle$
fillcontours=yes | no *labelcontours*=yes | no
negativedashed=yes | no *numcontours*= $\langle \text{n} \rangle$
showlegend=yes | no *title*= $\langle \text{text} \rangle$ *xaxislabel*= $\langle \text{text} \rangle$
yaxislabel= $\langle \text{text} \rangle$

Options:

<i>colormap</i>	The colormap to be used for displaying the contour plot Valid values: text strings Default value: jet
<i>cwidth</i>	The width of the contours drawn on the plot Valid values: integers Default value: 1 Minimum: 1 Maximum: 5
<i>emax</i>	The maximum energy value displayed on the graph Valid values: reals Default value: 0

Chapter 5: Commands

<i>emin</i>	The minimum energy value displayed on the graph Valid values: reals Default value: 0
<i>fillcontours</i>	Whether contours are filled or drawn as lines Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>labelcontours</i>	Whether contours are labelled directly on the plot. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>negativedashed</i>	A boolean which controls whether contours corresponding to negative values will be drawn with a dashed line. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>numcontours</i>	The number of contours used on the graph Valid values: integers Default value: 10 Minimum: 2 Maximum: 20
<i>showlegend</i>	Whether the legend is displayed for the contour plot Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>title</i>	The plot title Valid values: text strings Default value:
<i>xaxislabel</i>	The label for the X-axis. Valid values: text strings Default value: Coordinate 1
<i>yaxislabel</i>	The label for the Y-axis. Valid values: text strings Default value: Coordinate 2

2dplot

Controls the appearance of a plot of the energy as a function of two dihedral angles.

Syntax:

```
2dplot contdecimal=<n> escale=absolute | relative
      orientation=landscape | portrait papersize=letter | a4
      squareplot=yes | no units=kj | kcal xdecimal=<n>
      ydecimal=<n> [(postscript_file_name)]
```

Options:

contdecimal

The number of decimal places for the contour labels

Valid values: integers

Default value: **0**

Minimum: 0

Maximum: 5

escale

The energy scale (absolute/relative) for the graph

Valid values: absolute
relative

Default value: **absolute**

orientation

The paper orientation for the postscript output

Valid values: landscape
portrait

Default value: **landscape**

papersize

The paper size for the postscript output

Valid values: letter
a4

Default value: **letter**

squareplot

A boolean which controls whether the plot will be constrained to be square.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

units

The energy units for the graph

Valid values: kj
kcal

Default value: **kj**

Chapter 5: Commands

xdecimal The number of decimal places for the X axis.

Valid values: integers

Default value: **0**

Minimum: 0

Maximum: 5

ydecimal The number of decimal places for the Y axis.

Valid values: integers

Default value: **0**

Minimum: 0

Maximum: 5

Operands:

[*(postscript_file_name)*]

If an operand is given, this is treated as the name of a file to which a postscript representation of the file is to be written.

2drescale

Rescales the plot axes back to be able to view the entire data range.

Syntax:

2drescale

addatomproperty

Creates a new atom-level property for the atoms which are selected in the Workspace.

Syntax:

addatomproperty *value*=*(text)* *(property)*

Options:

value An optional initial value for the property.

Valid values: text strings

Default value:

Operands:

⟨property⟩

The M2IO data name of the new property.

addfromhold

This command retrieves the structure from the hold set and adds it to the current on-screen structure. The added hold structure is placed at the periphery of the on-screen structure and then a scale-to-screen is done.

Syntax:

addfromhold ⟨hold_name⟩

Operands:

⟨hold_name⟩

The name of the hold. This must be the name which was specified when the hold was created using the “hold” command.

adjustangle

Adjust the bond angle specified by the 3 atoms to the given value.

Syntax:

adjustangle *angle*=⟨x⟩ *move*=attached | terminal | single
 ⟨atom1⟩ ⟨atom2⟩ ⟨atom3⟩

Options:

angle Value to which to set the angle

Valid values: reals

Default value: **0**

Minimum: 0.0

Maximum: 180.0

move This option sets how to move other atoms that attached to the selected moving atom. There are three options: (1) move all attached atoms, (2) move all attached terminal atoms, and (3) move only single atom. Default option is (1).

Valid values:	attached terminal single
Default value:	attached

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle \langle \text{atom3} \rangle$

Three atoms which are used to adjust an angle. Note that direction does matter. So a-b-c is different from c-b-a. The third atom 'points' to the part of the structure that will be adjusted.

adjustchirality

Adjustl chiral center(s) in the specified set of atoms.

Syntax:

adjustchirality $\langle \text{ASL} \rangle$

Operands:

$\langle \text{ASL} \rangle$

An ASL specification of the atoms which are to have their chiral centers inverted. Because of the way this works, this set of atom should involve at least whole molecules.

adjustdihedral

Adjust the dihedral specified by 4 atoms to the given value.

Syntax:

adjustdihedral *dihedral*= $\langle x \rangle$ *dihedralterm*= $\langle x \rangle$ *move*=attached
| terminal | single $\langle \text{atom1} \rangle \langle \text{atom2} \rangle \langle \text{atom3} \rangle \langle \text{atom4} \rangle$

Options:

<i>dihedral</i>	Value to which to set the torsion
Valid values:	reals
Default value:	0
Minimum:	-180.1
Maximum:	180.1

dihedralterm

Value to which to terminally attached to the dihedral

Valid values: reals
 Default value: **0**
 Minimum: -180.1
 Maximum: 180.1

move

This option sets how to move other atoms that attached to the selected moving atom. There are three options: (1) move all attached atoms, (2) move all attached terminal atoms, and (3) move only single atom. Default option is (1).

Valid values: attached
 terminal
 single
 Default value: **attached**

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle \langle \text{atom3} \rangle \langle \text{atom4} \rangle$

Four atoms which are used to adjust a torsion. Note that direction does matter. So a-b-c-d is different from d-c-b-a. The fourth atom 'points' to the part of the structure that will be adjusted.

adjustdistance

Adjust the distance specified by 2 atoms to the given value.

Syntax:

adjustdistance *distance*= $\langle x \rangle$ *move*=attached | terminal | single
 $\langle \text{atom1} \rangle \langle \text{atom2} \rangle$

Options:

distance Value to which to set the torsion

Valid values: reals
 Default value: **0**
 Minimum: 0.0

move

This option sets how to move other atoms that attached to the selected moving atom. There are three options: (1) move all attached atoms, (2) move all attached terminal atoms, and (3) move only single atom. Default option is (1).

Valid values:	attached terminal single
Default value:	attached

Operands:

`<atom1> <atom2>`

Two atoms which are used to adjust a distance. Note that direction does matter. So a-b is different from b-a. The second atom 'points' to the part of the structure that will be adjusted.

alias

Define an alias for a command. This allows an abbreviated symbol to be defined for any command.

Syntax:

alias `<alias_name> <definition>`

Operands:

`<alias_name> <definition>`

The first operand is the name of the alias. If this contains embedded spaces then it must be enclosed in double quotes. The remaining operands are the definition of the alias. When command processing takes place, all occurrences of the alias name (the first operand) will be replaced by the definition of the alias.

alignaxis

Specifies a pair of atoms of a molecule to align the molecule to one of X, Y, or Z axis.

Syntax:

alignaxis *axis*=x | y | z `<atom1> <atom2>`

Options:

axis This option determines to which axis will the molecule be aligned: x=1, y=2, z=3.

Valid values:	x y z
Default value:	x

Operands:

⟨atom1⟩ ⟨atom2⟩

The first atom will be put at the origine and the second atom will be at whatever the distance is on the chosen axis.

alignplane

Define a plane of all atoms in the set described by the ASL operand and align the plane to the XY plane.

Syntax:

alignplane *plane*=xy | yz | zx ⟨ASL⟩

Options:

plane This option determines to which plane will the molecule be aligned: xy=1, yz=2, zx=3.

Valid values:	xy yz zx
Default value:	xy

Operands:

⟨ASL⟩

A string in the atom specification language which describes the set of atoms which are used to define a plane to be aligned.

angle

Specifies a triplet of atoms to have their bond angle measured and displayed.

Syntax:

angle *xoffset*=⟨x⟩ *yoffset*=⟨x⟩ ⟨atom1⟩ ⟨atom2⟩ ⟨atom3⟩

Options:

xoffset Specifies the X offset in Angstroms. Any angle created after this value is set will use this new value. Can be applied to an existing measurement by respecifying the measurement. Does not affect any already created angles.

Valid values: reals

Default value: **0**

yoffset Specifies the Y offset in Angstroms. Any angle created after this value is set will use this new value. Can be applied to an existing measurement by respecifying the measurement. Does not affect any already created angles.

Valid values: reals

Default value: **0**

Operands:

⟨atom1⟩ ⟨atom2⟩ ⟨atom3⟩

The three atoms between which the angle is to be measured. Note that the specifying a-b-c is the same as specifying c-b-a.

appendribbons

Generate a ribbon for the parts that are not currently in ribbon.

Syntax:

appendribbons

application

Determines which backend application is currently active and which application-specific menu is displayed. For example the command application Impact will result in the Impact menu being displayed in the main menu bar.

Syntax:

application \langle application_name \rangle

Operands:

\langle application_name \rangle

The name of the backend application which Maestro is currently interacting with. Current options are “none”, “macromodel” or “impact”.

assigncharges

This keyword is used to assign partial charges to all atoms in the workspace. The forcefield used is OPLS2005.

Syntax:

assigncharges

atom

Create a new atom in space with the current type and at the position specified by the operands.

Syntax:

atom *by*=element | type *element*= \langle text \rangle *type*= \langle n \rangle \langle x \rangle \langle y \rangle \langle z \rangle

Options:

<i>by</i>	This options determines whether atoms are to be described by element symbols (the <i>element</i> = option) or by atom types (type= option).
	Valid values: element type
	Default value: element
<i>element</i>	This option sets the current element of all atoms to be created (or retyped) if the <i>by</i> option == element
	Valid values: text strings
	Default value: C
<i>type</i>	This option sets the current type of all atoms to be created (or retyped) if the <i>by</i> option == type
	Valid values: integers
	Default value: 3

Operands:

$\langle x \rangle \langle y \rangle \langle z \rangle$

The operands are three real numbers which specify the x, y and z coordinates for the new atom.

atomname

Set the PDB atom name for all atoms which match the ASL specification.

Syntax:

atomname $\langle \text{PDBNAME} \rangle \langle \text{ASL} \rangle$

Operands:

$\langle \text{PDBNAME} \rangle \langle \text{ASL} \rangle$

The first operand is the PDB atom name which will be used for all atoms which match the specification. The second operand is a valid ASL string which defines the set of atoms which are to have their atom names changed.

atomproplabel

Turns the display of additional (non-standard) atom properties on or off. These are the properties read from the structure file at run-time as opposed to the standard properties controlled by the 'labelatom' command.

Syntax:

atomproplabel *display*=yes | no $\langle \text{property display name} \rangle$

Options:

display A boolean option which determines if the atom property is to be displayed or not in the current label definition.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

Operands:

$\langle \text{property display name} \rangle$

The name of the atom property. It's not an error if this property does not exist. Note this can be the 'Display Name' as it appears in the atom labels panel or the 'Data Name' as it appears in the Maestro format file.

attach

Attach the currently selected fragment to the specified atom. The bond to this atom (and the atom itself) will be replaced by a bond from the incoming fragment.

Syntax:

```
attach <atom_num>
```

Operands:

```
<atom_num>
```

The atom number defines the bond to which the fragment be attached. The atom specified must be a terminally attached atom.

attachmentmarkerdump

Print out the current option values of the attachment marker command.

Syntax:

```
attachmentmarkerdump
```

attachmentmarkersettings

Set graphical data of attachment markers.

Syntax:

```
attachmentmarkersettings ambient=<x> blue=<x>
                           cornradius=<x> cylinderheight=<x> cylinderradius=<x>
                           diffuse=<x> drawstyle=solid | line emission=<x> green=<x>
                           linewidth=<n> reagentradius=<x> red=<x> selectblue=<x>
                           selectgreen=<x> selectred=<x> shininess=<x> sliceline=<n>
                           slicesolid=<n> specular=<x> stackline=<n> stacksolid=<n>
                           transparency=<x>
```

Options:

ambient Set material property - ambient, to its red, green, and blue components, for front face.

Valid values: reals

Chapter 5: Commands

	Default value: 0.5
	Minimum: 0.0
	Maximum: 1.0
<i>blue</i>	The blue component of attachment markers.
	Valid values: reals
	Default value: 0
	Minimum: 0.0
	Maximum: 1.0
<i>cornradius</i>	The radius of corn of attachment markers.
	Valid values: reals
	Default value: 0.55
	Minimum: 0.0
<i>cylinderheight</i>	The cylinder height ratio of attachment markers.
	Valid values: reals
	Default value: 0.6
	Minimum: 0.0
	Maximum: 1.0
<i>cylinderradius</i>	The radius of cylinder of attachment markers.
	Valid values: reals
	Default value: 0.26
	Minimum: 0.0
<i>diffuse</i>	Set material property - diffuse, to its red, green, and blue components, for front face.
	Valid values: reals
	Default value: 0.4
	Minimum: 0.0
	Maximum: 1.0
<i>drawstyle</i>	The styles of rendering attachment markers, they are: 1 - solid, and 2 - lines. Default is solid.
	Valid values: solid line
	Default value: solid
<i>emission</i>	Set material property - emission, to its red, green, and blue components, for front face.
	Valid values: reals
	Default value: 0.05
	Minimum: 0.0
	Maximum: 1.0

<i>green</i>	The green component of attachment markers. Valid values: reals Default value: 0.6 Minimum: 0.0 Maximum: 1.0
<i>linewidth</i>	Set the width of lines in drawing attachment. Valid values: integers Default value: 2 Minimum: 1
<i>reagentradius</i>	The radius of sphere for attachment markers having associated reagents. Valid values: reals Default value: 0.8 Minimum: 0.0
<i>red</i>	The red component of attachment markers. Valid values: reals Default value: 1 Minimum: 0.0 Maximum: 1.0
<i>selectblue</i>	The blue component of selected attachment markers. Valid values: reals Default value: 1 Minimum: 0.0 Maximum: 1.0
<i>selectgreen</i>	The green component of selected attachment markers. Valid values: reals Default value: 0.9 Minimum: 0.0 Maximum: 1.0
<i>selectred</i>	The red component of selected attachment markers. Valid values: reals Default value: 0.2 Minimum: 0.0 Maximum: 1.0
<i>shininess</i>	Set material property - shininess, for front face. Valid values: reals Default value: 80 Minimum: 0.0 Maximum: 128.0

Chapter 5: Commands

<i>sliceline</i>	Set the slices of drawing line attachment. Valid values: integers Default value: 10 Minimum: 2
<i>slicesolid</i>	Set the slices of drawing solid attachment. Valid values: integers Default value: 36 Minimum: 2
<i>specular</i>	Set material property - specular, to its red, green, and blue components, for front face. Valid values: reals Default value: 0.2 Minimum: 0.0 Maximum: 1.0
<i>stackline</i>	Set the stacks of drawing line attachment. Valid values: integers Default value: 8 Minimum: 2
<i>stacksolid</i>	Set the stacks of drawing solid attachment. Valid values: integers Default value: 18 Minimum: 2
<i>transparency</i>	The transparency of rendering attachment markers. Valid values: reals Default value: 20 Minimum: 0.0 Maximum: 100.0

autosetup

Performs the automatic setup for the conformational search. The setup is only performed on those atoms which match the ASL definition given in the operand.

Syntax:

```

autosetup chiralatoms=yes | no compatoms=yes | no
             rings=yes | no torsionchecks=yes | no <ASL>

```

Options:

chiralatoms

A boolean option which determines whether the the automatic setup will include finding chiral atoms.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

compatoms

A boolean option which determines whether the the automatic setup will include finding comparison atoms.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

rings

A boolean option which determines whether the the automatic setup will include finding ring closures.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

torsionchecks

A boolean option which determines whether the the automatic setup will include finding torsion angles to be checked.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

Operands:

<ASL>

A string in the atom specification language. The automatic setup will be performed for only that atoms which match this ASL specification.

beginundoblock

Begin a new undoable command block.

Syntax:

```
beginundoblock
```

bmincomfile

Write a MacroModel command file with the current energy settings.

Syntax:

bmincomfile

bond

Create a new bond between two atoms or increments the bond order between two existing atoms.

Syntax:

bond $\langle \text{atom1} \rangle \langle \text{atom2} \rangle$

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle$

The operands are the numbers of the two atoms which are to be connected by a bond or have the bond order of an existing bond incremented.

bondorder

Increment, decrement or set the bond order for the specified bond.

Syntax:

bondorder increment|decrement| $\langle \text{bond_order} \rangle \langle \text{atom1} \rangle \langle \text{atom2} \rangle$

Operands:

increment|decrement| $\langle \text{bond_order} \rangle \langle \text{atom1} \rangle \langle \text{atom2} \rangle$

The first operand is either the word “increment”, the word “decrement”, or an integer value representing the bond order to be applied to the specified bond. Note that “0” is a valid bond order. The second and third operands are the atoms which define the bond which is to have its bond order changed.

bondoutline

Control features related to bond outlines

Syntax:

```
bondoutline autoscale=yes | no fixedwidth=<x> largewidth=<x>
           mediumwidth=<x> smallwidth=<x> technique=2pass | 3pass
           usewhenmoving=never | automatic | always
           visibleatomcutoff=<n> whenmovingcutoff=<n>
```

Options:

<i>autoscale</i>	The width of the outline border is automatically scaled when this is enabled. If this is off, then the value of width is used at all scales. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>fixedwidth</i>	Fixed width bond outline. Ignore if repall bondoutline is off or when bondoutline is on but autoscale is on. Valid values: reals Default value: 2 Minimum: 0.01 Maximum: 20.0
<i>largewidth</i>	The border width to use when the scale is large, meaning you are zoomed in and things appear large. Valid values: reals Default value: 4 Minimum: .5 Maximum: 20.0
<i>mediumwidth</i>	The border width to use when the scale is medium, meaning you are partially zoomed out. Valid values: reals Default value: 3 Minimum: .5 Maximum: 20.0
<i>smallwidth</i>	The border width to use when the scale is small, meaning you are zoomed out and things appear small. Valid values: reals Default value: 1 Minimum: .5 Maximum: 20.0
<i>technique</i>	Use a 2-pass + stenciling technique. The results are slightly different. Currently, only the 2-pass is used. 3-pass is not currently used and should do the same as using 2-pass. When bond

outlines are enabled you may not notice any performance degradation depending on your graphics system and what you're viewing. It's possible, however, that you may see a significant performance hit (roughly 5x).

Valid values: 2pass
 3pass
Default value: **2pass**

usewhenmoving

When set to always bond outlines will be drawn when the structure is being rotated. If never, then no bond outlines will be drawn when the structure is rotated. If set to automatic, then bond outlines are drawn but only when the number of atoms is less than the value of 'whenmovingcutoff'.

Valid values: never
 automatic
 always
Default value: **never**

visibleatomcutoff

If the number of visible atoms is larger than this amount, then do not display bond borders.

Valid values: integers
Default value: **1500**
Minimum: 1

whenmovingcutoff

If usewhenmoving is set to automatic, then see if the number of visible atoms is less than this amount. If so, then display bond borders. If it is greater than this amount, then do not display bond borders when the user is rotating or translating the structure.

Valid values: integers
Default value: **1500**

bondtonew

Create a new atom and bond it with a single bond to an existing atom.

Syntax:

bondtnew $\langle \text{atom_num} \rangle \langle x \rangle \langle y \rangle \langle z \rangle$

Operands:

$\langle \text{atom_num} \rangle \langle x \rangle \langle y \rangle \langle z \rangle$

The first operand is the number of an existing atom from which the bond is to be drawn. The remaining three operands are the x, y and z coordinates at which the new atom is to be placed.

buildoptions

Sets options associated with structure building. These are persistent and will be retained between Maestro sessions.

Syntax:

buildoptions *adjustbondlengths*=yes | no
 adjustnumhydrogens=yes | no *useunitedatoms*=yes | no

Options:

adjustbondlengths

When an atom type is changed, bond lengths around that atom will be set to ideal values.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

adjustnumhydrogens

During bond order, drawing or formal charge change operations, rectify the number of hydrogens to be consistent with the new bonding or formal charge.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

useunitedatoms

Will allow the generation of united atom types while structure building.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

calcenergy

Used to set options associated with a single point energy calculation.

Syntax:

calcenergy *listing*=none | minimal | complete

Options:

listing This option determines the extent of the listing of energy components.

Valid values: none
minimal
complete

Default value: **none**

Aliases:

ecalc (see [\[ecalc\]](#), page 127)

canonicalname

Sets atom names to one generated from atom properties (such as element or PDB names) for all the atoms matched by the ASL operand.

Syntax:

canonicalname <separator> <option_flag> <ASL>

Operands:

<separator> <option_flag> <ASL>

The first operand is the separator symbol to be used in generating the canonical names of all atoms. The second operand is the option_flag, which indicates which two atom properties will be used to construct the canonical names. Possible values are “3” - from element and atom number and “12” PDB chain, residue and atom name. The final operand is the ASL expression to indicate which atoms are to have a canonical name set.

caption

Create a Workspace caption and select it

Syntax:

`caption`

captionbackground

Caption backgroud settings

Syntax:

`captionbackground type=automatic | auto | none | user`

Options:

<i>type</i>	Background caption type
Valid values:	automatic auto none user
Default value:	automatic

captioncolor

Color of the caption text. Default color is yellow

Syntax:

`captioncolor alpha=<x> blue=<x> green=<x> red=<x>`

Options:

<i>alpha</i>	Alpha component. 0 is invisible. 1 is opaque.
Valid values:	reals
Default value:	0.67
Minimum:	0.0
Maximum:	1.0
<i>blue</i>	Blue component.
Valid values:	reals
Default value:	0
Minimum:	0.0
Maximum:	1.0
<i>green</i>	Green component.
Valid values:	reals

	Default value:	1
	Minimum:	0.0
	Maximum:	1.0
<i>red</i>	Red component.	
	Valid values:	reals
	Default value:	1
	Minimum:	0.0
	Maximum:	1.0

captiondelete

Delete the selected captions.

Syntax:

captiondelete

captionfont

Font settings for the caption

Syntax:

captionfont *name*=⟨text⟩ *size*=⟨x⟩ *style*=regular | bold | italic
| italicbold

Options:

<i>name</i>	Font name for the caption (helvetica, etc.)
	Valid values: text strings
	Default value: helvetica
<i>size</i>	Font size for the caption
	Valid values: reals
	Default value: 14
	Minimum: 3
	Maximum: 96
<i>style</i>	Font style for the caption (regular, bold, etc.)
	Valid values: regular
	bold
	italic
	italicbold

Default value: **regular**

captionhide

Hide the caption

Syntax:

captionhide

captionposition

Position of the caption.

Syntax:

captionposition $x=\langle x \rangle$ $y=\langle x \rangle$

Options:

x	Left-side X position of the caption. A fraction from zero to one. X origin is left-hand side of the window.
	Valid values: reals
	Default value: 0.4
	Minimum: 0.0
	Maximum: 0.98
y	Bottom left Y position of the caption. A fraction from zero to one. Y origin is bottom of window.
	Valid values: reals
	Default value: 0.1
	Minimum: 0.0
	Maximum: 0.98

captionselect

Select only the named caption

Syntax:

captionselect

captionselectadd

Add the named caption to the selected captions

Syntax:

captionselectadd

captionshow

Show the caption

Syntax:

captionshow

captionsoff

Hide all captions but leave their visible state alone.

Syntax:

captionsoff

captionson

Display all captions which are visible

Syntax:

captionson

captiontext

Text for the caption to be displayed in the Workspace

Syntax:

captiontext

captionunselect

Unselect the specified caption

Syntax:

captionunselect

captionunselectall

Unselect all captions

Syntax:

captionunselectall

cascadepanels

Cascade visible panels

Syntax:

cascadepanels

cd

This is a standard alias for **changedirectory** (see [\[changedirectory\]](#), page 68).

cellsmarkerdump

Print out the current option values of the Glide hydrophobic Cells markers command.

Syntax:

cellsmarkerdump

cellsmarkersettings

Set graphical data of Glide hydrophobic cells markers.

Syntax:

```
cellsmarkersettings activeblue=⟨x⟩ activegreen=⟨x⟩  
activered=⟨x⟩ ambient=⟨x⟩ blue=⟨x⟩ diffuse=⟨x⟩  
emission=⟨x⟩ green=⟨x⟩ highlightblue=⟨x⟩  
highlightgreen=⟨x⟩ highlightred=⟨x⟩ labelregions=yes | no  
red=⟨x⟩ regionblue=⟨x⟩ regiongreen=⟨x⟩ regionred=⟨x⟩  
roundingeffect=⟨x⟩ shininess=⟨x⟩ specular=⟨x⟩ step=⟨n⟩  
transparency=⟨x⟩
```

Options:

activeblue The blue color component of cell active frame and region label.

Valid values: reals
Default value: **0**
Minimum: 0.0
Maximum: 1.0

activegreen

The green color component of cell active frame and region label.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

activered The red color component of cell active frame and region label.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

ambient

Set material property - ambient, to its red, green, and blue components, for front face.

Valid values: reals
Default value: **0.5**
Minimum: 0.0
Maximum: 1.0

<i>blue</i>	<p>The blue color component of cell markers.</p> <p>Valid values: reals</p> <p>Default value: 0.75</p> <p>Minimum: 0.0</p> <p>Maximum: 1.0</p>
<i>diffuse</i>	<p>Set material property - diffuse, to its red, green, and blue components, for front face.</p> <p>Valid values: reals</p> <p>Default value: 0.5</p> <p>Minimum: 0.0</p> <p>Maximum: 1.0</p>
<i>emission</i>	<p>Set material property - emission, to its red, green, and blue components, for front face.</p> <p>Valid values: reals</p> <p>Default value: 0.1</p> <p>Minimum: 0.0</p> <p>Maximum: 1.0</p>
<i>green</i>	<p>The green color component of cell markers.</p> <p>Valid values: reals</p> <p>Default value: 0.75</p> <p>Minimum: 0.0</p> <p>Maximum: 1.0</p>
<i>highlightblue</i>	<p>The blue color component of cell markers if the cell is of a highlight region.</p> <p>Valid values: reals</p> <p>Default value: 0.3</p> <p>Minimum: 0.0</p> <p>Maximum: 1.0</p>
<i>highlightgreen</i>	<p>The green color component of cell markers if the cell is of a highlight region.</p> <p>Valid values: reals</p> <p>Default value: 0.3</p> <p>Minimum: 0.0</p> <p>Maximum: 1.0</p>
<i>highlightred</i>	<p>The red color component of cell markers if the cell is of a highlight region.</p> <p>Valid values: reals</p>

Chapter 5: Commands

	Default value: 1 Minimum: 0.0 Maximum: 1.0
<i>labelregions</i>	This option determines whether a region should be labeled (true) or not (false). Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>red</i>	The red color component of cell markers. Valid values: reals Default value: 0.75 Minimum: 0.0 Maximum: 1.0
<i>regionblue</i>	The blue color component of cell markers if the cell is of a region. Valid values: reals Default value: 0 Minimum: 0.0 Maximum: 1.0
<i>regiongreen</i>	The green color component of cell markers if the cell is of a region. Valid values: reals Default value: 0 Minimum: 0.0 Maximum: 1.0
<i>regionred</i>	The red color component of cell markers if the cell is of a region. Valid values: reals Default value: 1 Minimum: 0.0 Maximum: 1.0
<i>roundingeffect</i>	This determines the rounding effect of edges of a cell, for front face. Valid values: reals Default value: 8 Minimum: 0.0
<i>shininess</i>	Set material property - shininess, for front face. Valid values: reals Default value: 80 Minimum: 0.0 Maximum: 128.0

specular Set material property - specular, to its red, green, and blue components, for front face.

Valid values: reals
 Default value: **0.1**
 Minimum: 0.0
 Maximum: 1.0

step The step domain tolerance of cells.

Valid values: integers
 Default value: **3**
 Minimum: 1

transparency

The transparency of phobic cell markers.

Valid values: reals
 Default value: **30**
 Minimum: 0.0
 Maximum: 100.0

centeratom

Set global center of rotation to the given atom

Syntax:

centeratom *atom*=⟨n⟩

Options:

atom Atom to center transformations on

Valid values: integers
 Default value: **1**
 Minimum: 1

centerbond

Set global center of rotation to the given bond

Syntax:

centerbond *at1*=⟨n⟩ *at2*=⟨n⟩

Options:

<i>at1</i>	Atom1 of bond to be center of transformation
	Valid values: integers
	Default value: 1
	Minimum: 1
<i>at2</i>	Atom 2 of bond to be center of transformation
	Valid values: integers
	Default value: 1
	Minimum: 1

centercoordinates

Set global center of rotation to the given coordinate

Syntax:

centercoordinates *x*=⟨x⟩ *y*=⟨x⟩ *z*=⟨x⟩

Options:

<i>x</i>	X center
	Valid values: reals
	Default value: 0
<i>y</i>	Y center
	Valid values: reals
	Default value: 0
<i>z</i>	Z center
	Valid values: reals
	Default value: 0

centroid

This command takes the previously defined centroid atoms, averages the x, y, z co-ordinates and then creates a new atom of type 61 at the average position.

Syntax:

centroid

centroidatom

Define a single atom which will be used to define the centroid. Associated with this atom is a single atom marker.

Syntax:

centroidatom \langle atom_num \rangle

Operands:

\langle atom_num \rangle

The number of the atom to be used to define the centroid.

centroidatomset

Specifies a set of atoms to be added to the current centroid definition.

Syntax:

centroidatomset \langle ASL \rangle

Operands:

\langle ASL \rangle

A string in the atom specification language which describes the set of atoms which are to be added to the current centroid definition.

cglidedockconstraintposition

Specifies a constraint position in the receptor for a Glide calculation.

Syntax:

Chapter 5: Commands

```
cglidedockconstraintposition feature=⟨n⟩ index=⟨n⟩  
                                radius=⟨x⟩ radiusnoe=⟨x⟩ type=⟨n⟩ use1=yes | no  
                                use2=yes | no use3=yes | no use4=yes | no  
                                usenoe=yes | no x=⟨x⟩ y=⟨x⟩ z=⟨x⟩
```

Options:

<i>feature</i>	The constraint feature of position in docking. Valid values: integers Default value: -1
<i>index</i>	The index of position in docking. Valid values: integers Default value: 0
<i>radius</i>	The radius of a Glide constraint position. Valid values: reals Default value: 1 Minimum: 0.0001
<i>radiusnoe</i>	The minimum NOE distance of a Glide constraint position. Valid values: reals Default value: 0 Minimum: 0.0
<i>type</i>	The constraint type of position in docking. Valid values: integers Default value: 0
<i>use1</i>	The flag indicates if this position will be used in docking for group 1. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use2</i>	The flag indicates if this position will be used in docking for group 2. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use3</i>	The flag indicates if this position will be used in docking for group 3. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use4</i>	The flag indicates if this position will be used in docking for group 4. Valid values: boolean (true false; yes no; y n; on off) Default value: false

<i>usenoe</i>	Use the NOE constraint or not. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>x</i>	The X coordinate of a Glide constraint position. Valid values: reals Default value: 0
<i>y</i>	The Y coordinate of a Glide constraint position. Valid values: reals Default value: 0
<i>z</i>	The Z coordinate of a Glide constraint position. Valid values: reals Default value: 0

cglidedockconstraintregion

Specifies a constraint region in the receptor for a Glide calculation.

Syntax:

```
cglidedockconstraintregion atoms=<n> feature=<n>
index=<n> type=<n> use1=yes | no use2=yes | no
use3=yes | no use4=yes | no <region_name>
```

Options:

<i>atoms</i>	The number of required ligand atoms of constraint region in docking. Valid values: integers Default value: 1
<i>feature</i>	The constraint feature of constraint region in docking. Valid values: integers Default value: -1
<i>index</i>	The index of constraint region in docking. Valid values: integers Default value: 0
<i>type</i>	The constraint type of constraint region in docking. Valid values: integers Default value: 0

Chapter 5: Commands

<i>use1</i>	The flag indicates if this constraint region will be used in docking for group 1. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use2</i>	The flag indicates if this constraint region will be used in docking for group 2. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use3</i>	The flag indicates if this constraint region will be used in docking for group 3. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use4</i>	The flag indicates if this constraint region will be used in docking for group 4. Valid values: boolean (true false; yes no; y n; on off) Default value: false

Operands:

⟨region_name⟩

The name which will be applied to the region. If the name contains embedded spaces then it must be enclosed in double quotation marks.

chainname

Set the chain name for all atoms which match the ASL specification.

Syntax:

chainname ⟨CHAINNAME⟩ ⟨ASL⟩

Operands:

⟨CHAINNAME⟩ ⟨ASL⟩

The first operand is the PDB chain name (a single character) that will be used for all atoms which match the specification. The second operand is a valid ASL string which defines the set of atoms which are to have their residue names changed.

changedirectory

Change the current directory to that given by the operand of this command.

Syntax:

changedirectory \langle new_directory \rangle

Operands:

\langle new_directory \rangle

The name of the directory to be made the current directory.

Aliases:

cd (see [\[cd\]](#), page 59), **chdir** (see [\[chdir\]](#), page 69)

chdir

This is a standard alias for **changedirectory** (see [\[changedirectory\]](#), page 68).

chiralatom

Specifies an atom which is to be marked as “chiral” in the conformational search

Syntax:

chiralatom \langle atom-number \rangle

Operands:

\langle atom-number \rangle

The number of an atom which is to be added to the list of atoms which are to be considered as chiral in the conformational search.

clip

Adjust the clipping planes

Syntax:

Chapter 5: Commands

```
clip back=⟨x⟩ backselect=⟨x⟩ backsurface=⟨x⟩ boxoffset=⟨x⟩  
      front=⟨x⟩ frontselect=⟨x⟩ frontsurface=⟨x⟩  
      leftslopesurface=⟨x⟩ leftsurface=⟨x⟩ objects=all | surfaces |  
      surfacesselected | selected rightslopesurface=⟨x⟩  
      rightsurface=⟨x⟩
```

Options:

<i>back</i>	Position at which to set the back clipping plane Valid values: reals Default value: 0
<i>backselect</i>	Position at which to set the back clipping plane for selected atoms Valid values: reals Default value: 0
<i>backsurface</i>	Position at which to set the back clipping plane for surfaces Valid values: reals Default value: 0
<i>boxoffset</i>	The offset of surface clipping box for clipping surfaces to selected atoms. Valid values: reals Default value: 0
<i>front</i>	Position at which to set the front clipping plane Valid values: reals Default value: 0
<i>frontselect</i>	Position at which to set the front clipping plane for selected atoms Valid values: reals Default value: 0
<i>frontsurface</i>	Position at which to set the front clipping plane for surfaces Valid values: reals Default value: 0
<i>leftslopesurface</i>	Position at which to set the left slope clipping plane for surfaces Valid values: reals Default value: 0
<i>leftsurface</i>	Position at which to set the left clipping plane for surfaces Valid values: reals Default value: 0

objects What type of objects is to be clipped.

Valid values: all
 surfaces
 surfacesselected
 selected

Default value: **all**

rightslopesurface

Position at which to set the right slope clipping plane for surfaces

Valid values: reals

Default value: **0**

rightsurface

Position at which to set the right clipping plane for surfaces

Valid values: reals

Default value: **0**

clipatomset

Syntax:

clipatomset \langle ASL \rangle

Operands:

\langle ASL \rangle

A string in the atom specification language. This set defines clipping atoms for clipping atoms.

clipsurfaceset

Syntax:

clipsurfaceset \langle ASL \rangle

Operands:

\langle ASL \rangle

A string in the atom specification language. This set defines clipping atoms for clipping surfaces.

cliptoallatoms

Move clipping planes to all atoms.

Syntax:

cliptoallatoms

cliptoligand

Move clipping planes to ligand atoms.

Syntax:

cliptoligand

cliptoselectedatoms

Move clipping planes to selected atoms.

Syntax:

cliptoselectedatoms

clusteratom

Defines a single comparison atom to be used for XCluster jobs where clustering is based on Atomic RMS differences.

Syntax:

clusteratom \langle atom \rangle

Operands:

\langle atom \rangle

The number of an atom which is to added to the list of comparison atoms for a clustering calculation.

clusterheavyatoms

Specifies the atoms to be used in comparisons of conformers in XCluster jobs where clustering is based on Atomic RMS differences. The atoms specified are all of the heavy atoms and optionally hydrogens on oxygen and sulfur.

Syntax:

```
clusterheavyatoms osh=yes | no
```

Options:

<i>osh</i>	A boolean which controls whether hydrogens on oxygen and sulfur are included with heavy atoms.
Valid values:	boolean (true false; yes no; y n; on off)
Default value:	false

clusterset

Specifies a set of atoms to be used in comparisons of conformers in XCluster jobs where clustering is based on Atomic RMS differences.

Syntax:

```
clusterset <ASL>
```

Operands:

<ASL>

A string in the atom specification language. All atoms which match this description will be added to the list of XCluster comparison atoms.

clustertorsion

Specifies four atoms which define a torsion to be compared in XCluster jobs where clustering is based on Torsional RMS differences.

Syntax:

```
clustertorsion <atom1> <atom2> <atom3> <atom4>
```

Operands:

<atom1> <atom2> <atom3> <atom4>

The numbers of four atoms which define a torsion angle to be compared in XCluster jobs. Note that specifying a-b-c-d is the same as specifying d-c-b-a.

color

Specifies a color or set color options.

Syntax:

color *alpha*=⟨n⟩ *blue*=⟨n⟩ *green*=⟨n⟩ *red*=⟨n⟩

Options:

<i>alpha</i>	Alpha color component. Valid values: integers Default value: 128 Minimum: 0 Maximum: 255
<i>blue</i>	Blue color component. Valid values: integers Default value: 128 Minimum: 0 Maximum: 255
<i>green</i>	Green color component. Valid values: integers Default value: 128 Minimum: 0 Maximum: 255
<i>red</i>	Red color component. Valid values: integers Default value: 128 Minimum: 0 Maximum: 255

coloratom

Set the color for a group of atoms defined by the ASL operand.

Syntax:

coloratom *index*=⟨n⟩ *color*=⟨text⟩ ⟨ASL⟩

Options:

index An integer which indicates color index which is to be used for the atoms

Valid values: integers

Default value: **2**

Minimum: 1

Maximum: 256

color A string which is the color name for atom coloring. Valid color names are described in the file \$SCHRODINGER/maestro-vX.X/data/res/colors.res

Valid values: text strings

Default value:

Operands:

⟨ASL⟩

A string in the atom specification language. All atoms which match this specification will have their color changed to the current color.

coloratombonds

Set the color of all bonds, specified in ASL Change the color used to draw an on-screen bond.

Syntax:

coloratombonds ⟨ASL⟩

Operands:

⟨ASL⟩

A string in the atom specification language.

coloratomrgb

Set the color for a group of atoms defined by the ASL operand.

Chapter 5: Commands

Syntax:

coloratomrgb *blue*=⟨n⟩ *green*=⟨n⟩ *red*=⟨n⟩ ⟨ASL⟩

Options:

<i>blue</i>	Blue color component. Valid values: integers Default value: 128 Minimum: 0 Maximum: 255
<i>green</i>	Green color component. Valid values: integers Default value: 128 Minimum: 0 Maximum: 255
<i>red</i>	Red color component. Valid values: integers Default value: 128 Minimum: 0 Maximum: 255

Operands:

⟨ASL⟩

A string in the atom specification language. All atoms which match this specification will have their color changed to the current color.

colorbond

Set the color of a bond, specified by two atoms Change the color used to draw an on-screen bond.

Syntax:

colorbond *cindex*=⟨n⟩ *color*=⟨text⟩ *reset*=yes | no ⟨atom1⟩
⟨atom2⟩

Options:

<i>cindex</i>	An integer which indicates color index which is to be used for bond coloring. Valid values: integers
---------------	---

	Default value: 2
	Minimum: 1
	Maximum: 256
<i>color</i>	A string which is the color name for bond coloring. Valid color names are described in the file \$SCHRODINGER/maestro-vX.X/data/res/colors.res
	Valid values: text strings
	Default value:
<i>reset</i>	Resets the bond color to use atom colors.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false

Operands:

⟨ atom1 ⟩ ⟨ atom2 ⟩

The operands represent the numbers of the two atoms which define the bond to have its color changed.

colorscheme

Apply a predefined color scheme to the group of atoms defined by the ASL spec.

Syntax:

colorscheme *scheme*=⟨ text ⟩ ⟨ ASL ⟩

Options:

<i>scheme</i>	The name of the current color scheme.
	Valid values: text strings
	Default value: atype

Operands:

⟨ ASL ⟩

A string in the atom specification language. All atoms which match this specification will have their color changed to match the currently selected color scheme.

combilibenum

Start the job of Combinatorial Library Enumeration.

Syntax:

combilibenum

combilibenumaddattachment

Adds an attachment to the core molecule using the given atoms.

Syntax:

combilibenumaddattachment *atom1*=⟨n⟩ *atom2*=⟨n⟩
 ⟨attachment name⟩

Options:

atom1 The atom number of the atom in the original core to set as an attachment point. This is the atom which will be kept.

Valid values: integers

Default value: **1**

Minimum: 1

atom2 The atom number of the atom in the original core to set as an attachment point. This is the atom which will be removed.

Valid values: integers

Default value: **1**

Minimum: 1

Operands:

⟨attachment name⟩

The name of the attachment.

combilibenumclearreagentfile

Clears the reagent file for the selected rows.

Syntax:

combilibenumclearreagentfile

combilibenumdeleteattachment

Deletes all of the selected attachments.

Syntax:

combilibenumdeleteattachment

combilibenumexportdefinition

Stores the current core molecule and attachments in a file.

Syntax:

combilibenumexportdefinition <file name>

Operands:

<file name>

The name of the file to store the core definition in.

combilibenumimportdefinition

Reads a core molecule and attachments from the given file.

Syntax:

combilibenumimportdefinition <file name>

Operands:

<file name>

The name of the file to read the core definition from.

combilibenumoptions

This command holds general options for Combinatorial Library Enumeration.

Syntax:

combilibenumoptions *untangle*=yes | no

Options:

untangle An option which allows post-combgen minimization (for library enumeration).

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

combilibenumrefreshstructure

This function refreshes the structure in the Workspace from the current core structure in CombiGlide.

Syntax:

combilibenumrefreshstructure *viewcappedcore*=yes | no

Options:

viewcappedcore

An option which allows viewing of the minimally capped core, rather than the original core, in the Workspace.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

combilibenumrenameattachment

Renames the attachment in CombiLibEnum to the new name.

Syntax:

combilibenumrenameattachment *row*=⟨n⟩ ⟨new name⟩

Options:

row The row to rename.

Valid values: integers

Default value: **1**

Minimum: 1

Operands:

⟨ new name ⟩

The new name for the attachment.

combilibenumselectextendtablerow

Extends the selection to this row in the attachments table in the CombiGlide Library Enumeration.

Syntax:

```
combilibenumselectextendtablerow ⟨ row ⟩
```

Operands:

⟨ row ⟩

The row number to extend the select to.

combilibenumselectonlytablerow

Selects only this row in the attachments table in the CombiGlide Library Enumeration.

Syntax:

```
combilibenumselectonlytablerow ⟨ row ⟩
```

Operands:

⟨ row ⟩

The row number to select only in the table row.

combilibenumselecttablerow

Selects the given row in the attachments table in the CombiGlide Library Enumeration.

Syntax:

combilibenumselecttablerow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to select in the table.

combilibenumsetmolecule

Sets the core molecule for the CombiLibEnum to the molecule containing the given atom.

Syntax:

combilibenumsetmolecule *title*= $\langle \text{text} \rangle$ $\langle \text{atom number} \rangle$

Options:

title This option sets the title for the core molecule.

Valid values: text strings

Default value: **core**

Operands:

$\langle \text{atom number} \rangle$

The atom number of the molecule.

combilibenumsetreagentfile

Sets the reagent file for the selected rows.

Syntax:

combilibenumsetreagentfile $\langle \text{reagent name} \rangle$

Operands:

$\langle \text{reagent name} \rangle$

The name of the reagent file.

combilibenumunselecttablerow

Unselects the given row in the attachments table in the CombiGlide Library Enumeration.

Syntax:

```
combilibenumunselecttablerow <row>
```

Operands:

<row>

The row number to unselect in the table.

compareatom

Defines a single comparison atom to be used during a multiple minimization or conformational search. Comparison atoms are used to make comparisons in the process of determining if conformers are unique.

Syntax:

```
compareatom <atom>
```

Operands:

<atom>

The number of an atom which is to added to the list of comparison atoms during a conformational search.

compareset

Specifies a set of atoms to be used in comparisons of conformers during a multiple minimization or conformational search.

Syntax:

```
compareset <ASL>
```

Operands:

<ASL>

A string in the atom specification language. All atoms which match this description will be added to the list of comparison atoms.

confelim

This keyword is used to set various options associated with starting Redundant Conformer Elimination jobs from Maestro.

Syntax:

```
confelim compare_in_place=yes | no energy_source=none | jaguar  
| mm2* | mm3* | amber* | opls* | amber94 | mmff | mmffs |  
oplsaa | opls2005 | opls2008 incorporate=append | replace |  
ignore | appendungrouped | workspace
```

Options:

compare_in_place

A boolean which controls whether ConfElim will compare structures without first doing a superposition.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

energy_source

Source of energy for the comparisons

Valid values: none
 jaguar
 mm2*
 mm3*
 amber*
 opls*
 amber94
 mmff
 mmffs
 oplsaa
 opls2005
 opls2008

Default value: **opls2005**

incorporate

How the results are to be incorporated into the project. This can be done with replacement of the existing entries or by appending as new entries to the project or by ignoring the final results.

Valid values: append
 replace
 ignore
 appendungrouped
 workspace

Default value: **append**

confelimstart

Start a Redundant Conformer Elimination job with the current settings.

Syntax:

```
confelimstart
```

confelimwrite

Write a Redundant Conformer Elimination input file with the current settings.

Syntax:

```
confelimwrite
```

confgenltsearch

Defines settings for ConfGen conformational searching in MacroModel.

Syntax:

```
confgenltsearch amidebonds=vary | retain | trans
               compareatoms=none | heavy | heavy_polar_h
               distinguishenantiomers=yes | no eliminate=atom_deviation |
rmsd limitsave=yes | no max_rmsd=<x> maxdist=<x>
               maxringconf=<n> maxtorsdiff=<x> numsave=<n>
               numsteps=<n> samplerings=yes | no searchmode=standard |
               rapid | complete | thorough usenumsteps=<n>
               usesearchmoves=yes | no window=<x>
```

Options:

amidebonds

This determines whether to vary amide bond conformation, retain original amide bond conformation, or set amide bond conformation to trans.

Valid values: vary
 retain
 trans

Default value: **vary**

Chapter 5: Commands

compareatoms

This determines which kinds of comparison atoms are automatically identified for judging redundant conformers.

Valid values: none
 heavy
 heavy_polar_h

Default value: **heavy_polar_h**

distinguishenantiomers

A boolean for whether to retain different enantiomers. If true, an additional opcode, NANT, is written to the .com file.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

eliminate

The method to use for eliminating redundant conformers: maximum atom deviation or RMSD.

Valid values: atom_deviation
 rmsd

Default value: **rmsd**

limitsave

A boolean which controls whether number of structures saved is limited by numsave option (if true) or by numsteps (if false).

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

max_rmsd

Maximum RMSD for considering two structures equal.

Valid values: reals

Default value: **1**

Minimum: 0.0

maxdist

Maximum distance between atoms in equal structures.

Valid values: reals

Default value: **1.5**

Minimum: 0.0

maxringconf

The maximum number of ring conformations generated by the search.

Valid values: integers

Default value: **16**

Minimum: 1

maxtorsdiff

Maximum torsional angle difference between polar hydrogens in equal structures.

Valid values: reals

Default value: **60**

	Minimum: 0.0
	Maximum: 180.0
<i>numsave</i>	The number of structures that will be saved at the end of each search.
	Valid values: integers
	Default value: 1000
	Minimum: 0
<i>numsteps</i>	An option which sets the number of steps which will be performed during the ConfGen conformational search.
	Valid values: integers
	Default value: 1000
	Minimum: 0
<i>samplerings</i>	A boolean which controls whether to sample rings.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>searchmode</i>	This determines whether rapid (standard) or thorough (complete) search will be used.
	Valid values: standard
	rapid
	complete
	thorough
	Default value: standard
<i>usenumsteps</i>	An option which sets the number of steps which will be performed during the ConfGen conformational search.
	Valid values: integers
	Default value: 5
	Minimum: 1
<i>usesearchmoves</i>	A boolean which controls whether number of search moves is limited by maxsearch moves option (if true).
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>window</i>	The energy window (in kJ/mol) within which structures will be saved.
	Valid values: reals
	Default value: 100
	Minimum: 0.0

confgenmini

Used to set values associated with a MacroModel energy minimization for Ligand Torsional Search (ConfGen)

Syntax:

```
confgenmini converge=nothing | energy | gradient | movement  
             method=sd | prcg | osvm | fmnr | tncg | lbfgs | optimal  
             postmaxiter=⟨ n ⟩ premaxiter=⟨ n ⟩ threshold=⟨ x ⟩
```

Options:

converge This option determines which convergence criterion will be used during an energy minimization.

Valid values: nothing
 energy
 gradient
 movement

Default value: **gradient**

method This option determines which minimization method will be used.

Valid values: sd
 prcg
 osvm
 fmnr
 tncg
 lbfgs
 optimal

Default value: **tncg**

postmaxiter

This option determines the maximum number of iterations for post-minimization of generated structures.

Valid values: integers
Default value: **0**
Minimum: 0
Maximum: 9999999

premaxiter

This option determines the maximum number of iterations for pre-minimization of input structures.

Valid values: integers
Default value: **100**
Minimum: 0
Maximum: 9999999

threshold This option determines what the convergence threshold will be.

Valid values:	reals
Default value:	0.05
Minimum:	0.0

confgenpotential

Set various options associated with the definition of the potential energy to be used in a MacroModel job.

Syntax:

```
confgenpotential cele=⟨x⟩ charges=force_field | structure_file
                   chnd=⟨x⟩ cutoff=normal | extended | user_defined | none
                   cvdw=⟨x⟩ dielectric=⟨x⟩ electrostatics=field_field | constant
                   | distance_dependant field=mm2* | mm3* | amber* | op1s* |
                   amber94 | mmff | mmffs | op1saa | op1s2005 | op1s2008
                   solvent=none | water | chl3 | octanol
                   suppressbond=yes | no
```

Options:

cele This option determines what cutoff will be used for the electrostatic part of the energy calculation.

Valid values:	reals
Default value:	12
Minimum:	0.0
Maximum:	99999.0

charges This option determines where the charges to be used in the energy calculation will come from.

Valid values:	force_field structure_file
Default value:	force_field

chnd This option determines what cutoff will be used for the hydrogen bond part of the energy calculation.

Valid values:	reals
Default value:	4
Minimum:	0.0
Maximum:	99999.0

cutoff This option determines what type of non-bonded cutoff will be used in the energy calculation.

Chapter 5: Commands

	Valid values:	normal extended user_defined none
	Default value:	normal
<i>cvdw</i>	This option determines what VDW cutoff will be used in the energy calculation.	
	Valid values:	reals
	Default value:	7
	Minimum:	0.0
	Maximum:	99999.0
<i>dielectric</i>	The dielectric constant to be used in the electrostatic part of the energy calculation.	
	Valid values:	reals
	Default value:	4
	Minimum:	0.9999999999
<i>electrostatics</i>	The electrostatic treatment to be used in the energy calculation.	
	Valid values:	field_field constant distance_dependant
	Default value:	distance_dependant
<i>field</i>	The force field to be used for the energy calculation.	
	Valid values:	mm2* mm3* amber* opls* amber94 mmff mmffs oplsaa opls2005 opls2008
	Default value:	opls2005
<i>solvent</i>	The solvent model to be used for the energy calculation	
	Valid values:	none water chcl3 octanol
	Default value:	none

suppresshbond

A boolean which controls whether to suppress hydrogen bonding electrostatics.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

confgenreadpotential

Read potential settings from a command file.

Syntax:

confgenreadpotential <file_name>

Operands:

<file_name>

The name of the file from which the potential settings will be read. If no name is specified, the default settings will be used.

confgenstart

Start a Conformer Generation (Ligand Torsion Search) job with the current settings.

Syntax:

confgenstart

confgenstd

Launches a Standard ConfGen job.

Syntax:

confgenstd

confgenstdsettings

Used to set values associated with a Standard ConfGen job.

Syntax:

```
confgenstdsettings minimize_input=yes | no
                    minimize_output=yes | no strategy=very_fast | fast |
                    intermediate | comprehensive | fast_cf | phase_fast |
                    phase_quality structure_file=<text>
                    structure_source=selected_entries | workspace | file
```

Options:

minimize_input

Indicates whether or not the input structures will be minimized before the search.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

minimize_output

Indicates whether or not the output conformers will be minimized.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

strategy

This option determines which search strategy will be used.

Valid values: very_fast
 fast
 intermediate
 comprehensive
 fast_cf
 phase_fast
 phase_quality

Default value: **fast**

structure_file

If the source of structures is set to file, then this option contains the file to use.

Valid values: text strings

Default value:

structure_source

Whether to use the selected entries in the current project, or a specified file with multiple structures as structure input for the job.

Valid values:	selected_entries workspace file
Default value:	selected_entries

confgenstdwrite

Write out the Standard ConfGen job files.

Syntax:

```
confgenstdwrite
```

confgenwrite

Write a Conformer Generation (Ligand Torsion Search) input file with the current settings.

Syntax:

```
confgenwrite
```

confsearch

Defines settings for conformational searching in MacroModel.

Syntax:

```

confsearch distinguishenantiomers=yes | no
              eliminate=atom_deviation | rmsd enableauto=yes | no
              max_distance_lmcs=⟨x⟩ max_rmsd=⟨x⟩ maxdist=⟨x⟩
              method=mcm | serial_mcm | summ | mixed_lmcs |
              serial_mcm_lmcs | pure_lmcs | serial_lmcs | large_lmcs |
              mixed_large_lmcs min_distance_lmcs=⟨x⟩
              multiligand=yes | no numsteps=⟨n⟩ numstructures=⟨n⟩
              probability_tors_lmcs=⟨x⟩ searchvariable=ringclose | torsrot |
              moltrans | compatom | chiralatom | distcheck | torscheck |
              ligbond torsionsampling=restricted | intermediate | enhanced
              | extended usenumsteps=⟨n⟩ usesteps=yes | no window=⟨x⟩

```

Options:

distinguishenantiomers

This option determines to consider conformers more than once

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

eliminate

The method to use for eliminating redundant conformers: maximum atom deviation or RMSD.

Valid values: atom_deviation
 rmsd

Default value: **atom_deviation**

enableauto

Toggle whether or not to use AUTO setup for processing multiple structures.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

max_distance_lmcs

The maximum possible distance for the fastest moving atom in each LMCS move.

Valid values: reals

Default value: **6**

Minimum: 0.0

max_rmsd

Maximum RMSD for considering two structures equal.

Valid values: reals

Default value: **0.5**

Minimum: 0.0

maxdist

Maximum distance between atoms in equal structures.

Valid values: reals

Default value: **0.5**

Minimum: 0.0

- method* This determines which method will be used to perform the conformational search.
- Valid values: mmmm
serial_mmmm
summ
mixed_lmcs
serial_mmmm_lmcs
pure_lmcs
serial_lmcs
large_lmcs
mixed_large_lmcs
- Default value: **serial_mmmm_lmcs**
- min_distance_lmcs* The minimum possible distance for the fastest moving atom in each LMCS move.
- Valid values: reals
Default value: **3**
Minimum: 0.0
- multiligand* An option determines to use multiple ligands or not.
- Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**
- numsteps* An option which sets the number of steps which will be performed during the conformational search.
- Valid values: integers
Default value: **1000**
Minimum: 0
- numstructures* The number of structures that will be saved at the end of each search.
- Valid values: integers
Default value: **0**
Minimum: 0
- probability_tors_lmcs* The probability that a TORS/MOLS move will be made during an LMCS search.
- Valid values: reals
Default value: **0.5**
Minimum: 0.0
Maximum: 1.0

Chapter 5: Commands

searchvariable

Controls which search variable will be used.

Valid values: ringclose
 torsrot
 moltrans
 compatom
 chiralatom
 distcheck
 torscheck
 ligbond

Default value: **ringclose**

torsionsampling

Controls how sampling of torsions is done during the automatic setup for torsional sampling methods.

Valid values: restricted
 intermediate
 enhanced
 extended

Default value: **intermediate**

usenumsteps

An option which sets the number of steps which will be performed during the conformational search.

Valid values: integers
Default value: **100**
Minimum: 1

usesteps

An option which sets the number of steps which will be performed during the conformational search.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

window

The energy window (in kJ/mol) within which structures will be saved.

Valid values: reals
Default value: **21**
Minimum: 0.0

connect

Connect atom pairs.

Syntax:

connect

connectfuseatom

Define an atom pair for which will be connected in a subsequent connect operation.

Syntax:

connectfuseatom $\langle \text{atom1} \rangle \langle \text{atom2} \rangle$

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle$

Two atom numbers which represent an atom pair to be connected by a subsequent connect or fuse command. The two atoms must be from different molecules and all connectfuseatom commands must specify molecules in the same order.

constrainedangle

Specifies a angle between three atoms which is to be constrained by a harmonic constraint during a MacroModel calculation.

Syntax:

constrainedangle *angle*= $\langle x \rangle$ *constant*= $\langle x \rangle$ *select*=yes | no
width= $\langle x \rangle$ $\langle \text{atom1} \rangle \langle \text{atom2} \rangle \langle \text{atom3} \rangle$

Options:

<i>angle</i>	The angle at which the atoms are to be constrained. Valid values: reals Default value: -1
<i>constant</i>	The force constant for the harmonic constraint to be applied to the angle. Valid values: reals Default value: 100 Minimum: 0.0
<i>select</i>	Selection state of the model.

	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	false
<i>width</i>	The half-width for a flat-bottomed constraint.	
	Valid values:	reals
	Default value:	0
	Minimum:	0.0

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle \langle \text{atom3} \rangle$

Three atoms which are to have the angle between them constrained by a harmonic constraint. The three atoms are not necessarily connected. Note that specifying a-b-c will be treated by the program as the same as specifying c-b-a.

constrainedatom

Specifies a single atom to be constrained at its current positions during a MacroModel energy calculation. The atom will be constrained by the use of a harmonic constraint.

Syntax:

```
constrainedatom constant= $\langle x \rangle$  frozen=yes | no select=yes | no
                  width= $\langle x \rangle$   $\langle \text{atom\_number} \rangle$ 
```

Options:

<i>constant</i>	The force constant (in kJ/mol/Angs) for the harmonic constraint.	
	Valid values:	reals
	Default value:	0
<i>frozen</i>	The constraint atom is frozen or not	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	false
<i>select</i>	Selection state of the model.	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	false
<i>width</i>	The half-width of a flat bottomed restraint.	
	Valid values:	reals
	Default value:	0
	Minimum:	0

Operands:

$\langle \text{atom_number} \rangle$

The number of the atom to which a constraint is to be applied.

constraineddist

Specifies a distance between two atoms which is to be constrained by a harmonic constraint during a MacroModel calculation.

Syntax:

```
constraineddist constant= $\langle x \rangle$  distance= $\langle x \rangle$ 
                remove_nb=yes | no select=yes | no width= $\langle x \rangle$   $\langle \text{atom1} \rangle$ 
                 $\langle \text{atom2} \rangle$ 
```

Options:

<i>constant</i>	The force constant for the harmonic constraint to be applied to the distance. Valid values: reals Default value: 100 Minimum: 0.0
<i>distance</i>	The distance at which the atoms are to be constrained. Valid values: reals Default value: -1
<i>remove_nb</i>	A boolean option which controls whether the non-bonded interaction (if any) between the two specified atoms will be removed. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>select</i>	Selection state of the model. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>width</i>	The half-width for a flat-bottomed constraint. Valid values: reals Default value: 0 Minimum: 0.0

Operands:

$\langle \text{atom1} \rangle$ $\langle \text{atom2} \rangle$

Two atoms which are to have the distance between them constrained by a harmonic constraint. Note that specifying a-b will be treated by the program as the same as specifying b-a.

constrainedset

Specifies a set of atoms to be constrained at their current positions. The atoms will be constrained by the use of a harmonic constraint.

Syntax:

constrainedset *constant*=⟨x⟩ *width*=⟨x⟩ ⟨ASL⟩

Options:

constant The force constant (in kJ/mol/Angs) for the harmonic constraint.

Valid values:	reals
Default value:	100
Minimum:	0.0

width The half-width of a flat bottomed restraint

Valid values:	reals
Default value:	0
Minimum:	0

Operands:

⟨ASL⟩

The operand must be a valid string in the atom specification language. Any atoms which match will have constraints generated for them.

constrainedtorsion

Specifies a torsion between four atoms which is to be constrained by a harmonic constraint during a MacroModel calculation.

Syntax:

constrainedtorsion *constant*=⟨x⟩ *select*=yes | no *torsion*=⟨x⟩
width=⟨x⟩ ⟨atom1⟩ ⟨atom2⟩ ⟨atom3⟩ ⟨atom4⟩

Options:

constant The force constant for the harmonic constraint to be applied to the torsion.

Valid values: reals
 Default value: **100**
 Minimum: 0.0

select Selection state of the model.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

torsion The torsion at which the atoms are to be constrained.

Valid values: reals
 Default value: **500**

width The half-width for a flat-bottomed constraint.

Valid values: reals
 Default value: **0**
 Minimum: 0.0

Operands:

⟨atom1⟩ ⟨atom2⟩ ⟨atom3⟩ ⟨atom4⟩

Four atoms which are to have the angle between them constrained by a harmonic constraint. Note that specifying a-b-c-d is treated by the program the same as specifying d-c-b-a.

contactcriteria

Specify the criteria for calculating good, bad, ugly contacts. There is no upper-bound to C. $C = \{\text{dist between the two atoms}\} \text{ divided by } \{\text{vdW radius first atom} + \text{vdW radius of second atom}\}$.

Chapter 5: Commands

Syntax:

```
contactcriteria bad=<x> display=yes | no displaybad=yes | no
displaygood=yes | no displayugly=yes | no
exclude14interactions=yes | no excludehbond=yes | no
good=<x> ugly=<x>
```

Options:

bad Criteria for bad contact

Valid values: reals

Default value: **0.89**

Minimum: 0.0

display This option determines whether contact markers will be displayed.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

displaybad If contact markers are shown, this option determines whether the bad contact markers will be displayed.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

displaygood

If contact markers are shown, this option determines whether the good contact markers will be displayed.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

displayugly

If contact markers are shown, this option determines whether the ugly contact markers will be displayed.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

exclude14interactions

If this option is true, atoms that have 1,4 interactions should not be considered to have bad or ugly contacts.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

excludehbond

If this option is true, atoms that are H-bonded should not be considered to have bad or ugly contacts.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

<i>good</i>	Criteria for good contact
	Valid values: reals
	Default value: 1.3
	Minimum: 0.0
<i>ugly</i>	Criteria for ugly contact
	Valid values: reals
	Default value: 0.75
	Minimum: 0.0

contactset1

Specify the first set of atoms used in determining good, bad, and ugly contacts.

Syntax:

contactset1 <ASL>

Operands:

<ASL>

A string in the atom specification language. Typical usage is to define contactset1 and contactset2. This set, contactset1, defines the “from” atoms. The contactset2 atoms define the “to” atoms. Contacts are calculated between these two sets. That is, the contacts are inter-set contacts. No intra-set contacts are calculated. If contactset2’s ASL string is empty, then contacts are calculated for all atoms in contactset1.

contactset2

Specify the second set of atoms used in determining good, bad, and ugly contacts.

Syntax:

contactset2 <ASL>

Operands:

<ASL>

A string in the atom specification language. Typical usage is to define contactset1 and contactset2. This set, contactset2, defines the “to” atoms. The

contactset1 atoms define the “from” atoms. Contacts are calculated between these two sets. That is, the contacts are inter-set contacts. No intra-set contacts are calculated. If contactset2’s ASL string is empty, then contacts are calculated for all atoms in contactset1.

coordinatescan

Used to set a distance or angle or dihedral coordinate to be scanned. The four operands are the atom numbers defining the dihedral to be driven. A maximum of two angles can be driven. The four operands are the atom numbers defining the coordinate. Depending on the type of scan coordinate set in MM_SCAN_MODE: For distance the third and fourth atom numbers will be 0 and for angle the fourth one will be 0. In case of dihedral, all the four atom numbers are used. A maximum of two coordinates can be driven.

Syntax:

```
coordinatescan finishangle= $\langle x \rangle$  finishdihedral= $\langle x \rangle$ 
                finishdist= $\langle x \rangle$  incrementangle= $\langle x \rangle$  incrementdihedral= $\langle x \rangle$ 
                incrementdist= $\langle x \rangle$  startangle= $\langle x \rangle$  startdihedral= $\langle x \rangle$ 
                startdist= $\langle x \rangle$   $\langle \text{atom1} \rangle$   $\langle \text{atom2} \rangle$   $\langle \text{atom3} \rangle$   $\langle \text{atom4} \rangle$ 
```

Options:

finishangle

Specifies the finishing angle for the angle scan.

Valid values: reals

Default value: **360**

finishdihedral

Specifies the finishing angle for the dihedral scan.

Valid values: reals

Default value: **360**

finishdist

Specifies the finishing distance for the distance scan.

Valid values: reals

Default value: **10**

incrementangle

Specifies the angle increment for the angle scan.

Valid values: reals

Default value: **30**

incrementdihedral

Specifies the angle increment for the dihedral scan.

Valid values: reals
 Default value: **30**

incrementdist

Specifies the distance increment for the distance scan.

Valid values: reals
 Default value: **1**

startangle Specifies the starting angle for the angle scan.

Valid values: reals
 Default value: **0**

startdihedral

Specifies the starting angle for the dihedral scan.

Valid values: reals
 Default value: **0**

startdist Specifies the starting distance for the distance scan.

Valid values: reals
 Default value: **1**

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle \langle \text{atom3} \rangle \langle \text{atom4} \rangle$

The four operands are treated as atom numbers which define the distance/angle/dihedral to be driven. The MM_SCAN_MODE option will determine the type of scan.

coupling

Specifies a pair of atoms to have their coupling measured and displayed.

Syntax:

coupling *xoffset*= $\langle x \rangle$ *yoffset*= $\langle x \rangle$ $\langle \text{atom1} \rangle \langle \text{atom2} \rangle$

Options:

xoffset Specifies the X offset in Angstroms. Any coupling created after this value is set will use this new value. Can be applied to an existing measurement by respecifying the measurement. Does not affect any already created couplings.

Valid values: reals
 Default value: **0**

yoffset Specifies the Y offset in Angstroms. Any coupling created after this value is set will use this new value. Can be applied to an existing measurement by respecifying the measurement. Does not affect any already created couplings.

Valid values: reals

Default value: **0**

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle$

The two hydrogen atoms between which the coupling is to be measured. Note that the specifying a-b is the same as specifying b-a.

createcustompropertygroup

This command creates a custom property group for the given space separated property names.

Syntax:

```
createcustompropertygroup  $\langle \text{propertygroupname} \rangle$   
                           $\langle \text{propertynames} \rangle$ 
```

Operands:

$\langle \text{propertygroupname} \rangle \langle \text{propertynames} \rangle$

The name of the custom property group. The space separated property names use to create the custom property group.

createlibrary

Start the creation of the library in the Analyse Library step of CombiGlide.

Syntax:

```
createlibrary
```

createpropsubset

This command creates a property subset for the given table using the given property names.

Syntax:

createpropsubset \langle table \rangle \langle propertynames \rangle

Operands:

\langle table \rangle \langle propertynames \rangle

The number of the table to create the subset for. If the table operand is missing, no resize will be done. The names of the properties to use to create the subset.

createsubset

Creates a subset in the project table using the currently selected entries. This function also switches to subset view.

Syntax:

createsubset

defaultfc

Requires a single operand which is the default force constant to be used for the constrainedatom command. The default value is used only when no value has been specified by the constant= option of the “constrainedatom” command.

Syntax:

defaultfc \langle default_force_constant_value \rangle

Operands:

\langle default_force_constant_value \rangle

The value of the default force constant in kJ/mol/Angstrom

delete

Delete a named object. The object type is the same as the command which is used to create that type of object. For example to delete a set named “set1” use: delete set set1. All instances of any type of object can be deleted with

the “all” name. There are three special objects: “atom” and “bond” for deleting atoms and bonds respectively

For example: delete set set1

delete constrainedatom 1

delete set all delete res. 1-5 delete bond 1 2

Syntax:

```
delete includeterminal=yes | no <object_type>  
      <object_name> | all | <ASL> | <bond>
```

Options:

includeterminal

This option is used to control the behavior when deleting atoms and bonds. If this is “true” then all terminally attached atoms will also be deleted when deleting bonds and atoms.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

Operands:

```
<object_type> <object_name> | all | <ASL> | <bond>
```

The first operand determines what type of object is to be deleted. This can be the name of any of the named objects in Maestro (sets, filters, holds etc.) or can be “atom” or “bond”. The subject operands depend on the first operand. If the first operand names an object then the second operand must be the name which was specified when that object was created. If the first operand is “atom” then the second operand must be a valid string in the atom specification language and all atoms which match that string will be deleted. Finally if the first operand is “bond” then the second and third operands are expected to be the numbers of the two atoms which define that bond.

deletecustompropertygroup

This command deletes a custom property group.

Syntax:

```
deletecustompropertygroup <propertygroupname>
```

Operands:

```
<propertygroupname>
```

The name of the custom property group.

deleteemptyentrygroups

Deletes all the empty entry groups.

Syntax:

deleteemptyentrygroups

deleteproperty

This is a standard alias for **propertydelete** (see [\[propertydelete\]](#), page 604).

deleteselectedproperties

Deletes the selected properties from the atom-property table and the Workspace.

Syntax:

deleteselectedproperties

deletezeroorderbonds

Delete zero order bonds connected to atoms in the specified set of atoms. The zero order bond will be deleted if it connects any atom within the set.

Syntax:

deletezeroorderbonds \langle ASL \rangle

Operands:

\langle ASL \rangle

The set of atoms which will have zero order bonds deleted.

densityblobupdate

Update electron density blobs for current density map.

Syntax:

densityblobupdate *minvol*=⟨x⟩

Options:

<i>minvol</i>	Threshold for displaying density blobs. Only blobs whose volumes are greater than or equal to this number will be displayed.
Valid values:	reals
Default value:	0.05
Minimum:	0.0

dialogoptions

Sets options for dialog boxes.

Syntax:

dialogoptions *show*=yes | no

Options:

<i>show</i>	If this option is set to false, then dialog boxes will not be displayed.
Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

dihedral

Specifies a quartet of atoms to have their dihedral angle measured and displayed.

Syntax:

dihedral *xoffset*=⟨x⟩ *yoffset*=⟨x⟩ ⟨atom1⟩ ⟨atom2⟩ ⟨atom3⟩
 ⟨atom4⟩

Options:

xoffset Specifies the X offset in Angstroms. Any dihedral created after this value is set will use this new value. Can be applied to an existing measurement by respecifying the measurement. Does not affect any already created dihedrals.

Valid values: reals

Default value: **0**

yoffset Specifies the Y offset in Angstroms. Any dihedral created after this value is set will use this new value. Can be applied to an existing measurement by respecifying the measurement. Does not affect any already created dihedrals.

Valid values: reals

Default value: **0**

Operands:

⟨atom1⟩ ⟨atom2⟩ ⟨atom3⟩ ⟨atom4⟩

The four atoms between which the dihedral angle is to be measured. Note that the specifying a-b-c-d is the same as specifying d-c-b-a.

dihedraldrive

Used to set a dihedral angle to be driven. The four operands are the atom numbers defining the dihedral to be driven. A maximum of two angles can be driven.

Syntax:

dihedraldrive *finish*=⟨x⟩ *increment*=⟨x⟩ *start*=⟨x⟩ ⟨atom1⟩
 ⟨atom2⟩ ⟨atom3⟩ ⟨atom4⟩

Options:

finish Specifies the finishing angle for the dihedral drive.

Valid values: reals

Default value: **360**

increment Specifies the angle increment for the dihedral drive.

Valid values: reals

Default value: **30**

Chapter 5: Commands

start Specifies the starting angle for the dihedral drive.
Valid values: reals
Default value: **0**

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle \langle \text{atom3} \rangle \langle \text{atom4} \rangle$

The four operands are treated as atom numbers which define the dihedral to be driven.

dipolemoment

Set dipole moment options.

Syntax:

dipolemoment *ambient*= $\langle x \rangle$ *arrowwidth*= $\langle x \rangle$ *chargeblue*= $\langle x \rangle$
chargegreen= $\langle x \rangle$ *chargered*= $\langle x \rangle$ *diffuse*= $\langle x \rangle$ *diskblue*= $\langle x \rangle$
diskgreen= $\langle x \rangle$ *diskoffset*= $\langle x \rangle$ *diskred*= $\langle x \rangle$ *diskwidth*= $\langle x \rangle$
display=yes | no *emission*= $\langle x \rangle$ *frompositive*=yes | no
jagusrblue= $\langle x \rangle$ *jagusrgreen*= $\langle x \rangle$ *jagusrred*= $\langle x \rangle$
negativeblue= $\langle x \rangle$ *negativegreen*= $\langle x \rangle$ *negativered*= $\langle x \rangle$
positiveblue= $\langle x \rangle$ *positivegreen*= $\langle x \rangle$ *positivered*= $\langle x \rangle$
ringblue= $\langle x \rangle$ *ringgreen*= $\langle x \rangle$ *ringred*= $\langle x \rangle$ *ringustep*= $\langle n \rangle$
ringustep= $\langle n \rangle$ *ringwidth*= $\langle x \rangle$ *scale*= $\langle x \rangle$ *shaftwidth*= $\langle x \rangle$
shininess= $\langle x \rangle$ *showdisk*=yes | no *showpositive*=yes | no
specular= $\langle x \rangle$ *transparency*= $\langle x \rangle$ *ustep*= $\langle n \rangle$ *vstep*= $\langle n \rangle$

Options:

ambient Set material property - ambient, to its red, green, and blue components, for front face.
Valid values: reals
Default value: **0.5**
Minimum: 0.0
Maximum: 1.0

arrowwidth The width of dipole moment arrow.
Valid values: reals
Default value: **0.2**
Minimum: 0.0001

chargeblue The blue color component for dipole moments defined by partial atomic charges.

	Valid values: reals
	Default value: 0.5
	Minimum: 0.0
	Maximum: 1.0
<i>chargegreen</i>	The green color component for dipole moments defined by partial atomic charges.
	Valid values: reals
	Default value: 1
	Minimum: 0.0
	Maximum: 1.0
<i>chargedred</i>	The red color component for dipole moments defined by partial atomic charges.
	Valid values: reals
	Default value: 0.8
	Minimum: 0.0
	Maximum: 1.0
<i>diffuse</i>	Set material property - diffuse, to its red, green, and blue components, for front face.
	Valid values: reals
	Default value: 0.4
	Minimum: 0.0
	Maximum: 1.0
<i>diskblue</i>	The blue color component of disk.
	Valid values: reals
	Default value: 0
	Minimum: 0.0
	Maximum: 1.0
<i>diskgreen</i>	The green color component of disk.
	Valid values: reals
	Default value: 1
	Minimum: 0.0
	Maximum: 1.0
<i>diskoffset</i>	The offset (percentage) of distance between the dipole moment disk and the positive end.
	Valid values: reals
	Default value: 0.05
	Minimum: 0.0
	Maximum: 1.0
<i>diskred</i>	The red color component of disk.

Chapter 5: Commands

	Valid values: reals Default value: 1 Minimum: 0.0 Maximum: 1.0
<i>diskwidth</i>	The width of dipole moment disk. Valid values: reals Default value: 0.3 Minimum: 0.0001
<i>display</i>	Whether or not display dipole moment in workspace. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>emission</i>	Set material property - emission, to its red, green, and blue components, for front face. Valid values: reals Default value: 0.15 Minimum: 0.0 Maximum: 1.0
<i>frompositive</i>	Set to true if the arrow is drawn from positive to negative. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>jagusrblue</i>	The blue color component for dipole moments defined by Jaguar properties. Valid values: reals Default value: 1 Minimum: 0.0 Maximum: 1.0
<i>jagusrgreen</i>	The green color component for dipole moments defined by Jaguar properties. Valid values: reals Default value: 0.7 Minimum: 0.0 Maximum: 1.0
<i>jagusrred</i>	The red color component for dipole moments defined by Jaguar properties. Valid values: reals Default value: 0.8 Minimum: 0.0 Maximum: 1.0

negativeblue

The blue color component of dipole moment markers for negative region.

Valid values: reals
 Default value: **0**
 Minimum: 0.0
 Maximum: 1.0

negativegreen

The green color component of dipole moment markers for negative region.

Valid values: reals
 Default value: **0**
 Minimum: 0.0
 Maximum: 1.0

negativered

The red color component of dipole moment markers for negative region.

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 1.0

positiveblue

The blue color component of dipole moment markers for positive region.

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 1.0

positivegreen

The green color component of dipole moment markers for positive region.

Valid values: reals
 Default value: **0**
 Minimum: 0.0
 Maximum: 1.0

positivered

The red color component of dipole moment markers for positive region.

Valid values: reals
 Default value: **0**
 Minimum: 0.0
 Maximum: 1.0

Chapter 5: Commands

<i>ringblue</i>	The blue color component of the ring. Valid values: reals Default value: 1 Minimum: 0.0 Maximum: 1.0
<i>ringgreen</i>	The green color component of the ring. Valid values: reals Default value: 1 Minimum: 0.0 Maximum: 1.0
<i>ringred</i>	The red color component of the ring. Valid values: reals Default value: 0 Minimum: 0.0 Maximum: 1.0
<i>ringustep</i>	The U step domain tolerance for the ring. Valid values: integers Default value: 5 Minimum: 2
<i>ringvstep</i>	The V step domain tolerance for the ring. Valid values: integers Default value: 3 Minimum: 2
<i>ringwidth</i>	The width of the ring around the disk. Valid values: reals Default value: 0.05 Minimum: 0.0001
<i>scale</i>	The scale to convert unit from Debye to Angstrong for the length of dipole moment arrow. Valid values: reals Default value: 1 Minimum: 0.0001
<i>shaftwidth</i>	The width of dipole moment shaft. Valid values: reals Default value: 0.1 Minimum: 0.0001
<i>shininess</i>	Set material property - shininess, for front face. Valid values: reals Default value: 80

	Minimum:	0.0
	Maximum:	128.0
<i>showdisk</i>	Set to true if the disk for positive sign is shown.	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	true
<i>showpositive</i>	Set to true if positive and negative colors are used.	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	true
<i>specular</i>	Set material property - specular, to its red, green, and blue components, for front face.	
	Valid values:	reals
	Default value:	0.7
	Minimum:	0.0
	Maximum:	1.0
<i>transparency</i>	The transparency of dipole moment markers.	
	Valid values:	reals
	Default value:	0.75
	Minimum:	0.0
	Maximum:	1.0
<i>ustep</i>	The U step domain tolerance for the cylinder.	
	Valid values:	integers
	Default value:	3
	Minimum:	2
<i>vstep</i>	The V step domain tolerance for the cylinder.	
	Valid values:	integers
	Default value:	5
	Minimum:	2

displayatom

Display all atoms in the set described by the ASL operand.

Syntax:

displayatom \langle ASL \rangle

Operands:

\langle ASL \rangle

A string in the atom specification language which describes the set of atoms which are to be displayed.

displayonlyatom

Display only the set described by the ASL operand.

Syntax:

displayonlyatom \langle ASL \rangle

Operands:

\langle ASL \rangle

A string in the atom specification language which describes the set of atoms which are to be displayed.

displayopt

Set properties for the display, such as the background color.

Syntax:

displayopt *adjustclip*=yes | no *aligndepth*= \langle x \rangle
angledependenttransparency=yes | no *bgcindex*= \langle n \rangle
bgcolor= \langle text \rangle *eye_sepf*= \langle x \rangle *fog*=auto | on | true | yes | off
| false | no *fogatomlabels*=yes | no *fogcutoff*= \langle n \rangle
fogdensity= \langle x \rangle *fogendz*= \langle x \rangle *fogstartz*= \langle x \rangle *fogtype*=linear |
exp | exp2 | automatic *fsaa*=yes | no *invertstereo*=yes | no
maximize=yes | no *opengl_level*= \langle text \rangle *perspective*=yes | no
perspectivescale= \langle x \rangle *rubber_banding*=xor | main *sepf*= \langle x \rangle
sizef= \langle x \rangle *start_str*= \langle text \rangle *stereo*=yes | no
stereomethod=hardware | crosseyed | walleyed | fullscreen |
interlaced *stop_str*= \langle text \rangle

Options:

adjustclip A boolean option which determines whether to adjust clipping to preserve stereo depth.

- Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**
- aligndepth* Relative z-depth (1 for back clipping plane, 0 for front clipping plane) at which stereo images appear to converge. Moves region inside clipping planes in or out of the screen.
- Valid values: reals
 Default value: **0.5**
 Minimum: 0.0
 Maximum: 1.0
- angledependenttransparency*
 A bool which determines whether transparent objects should use angle dependent transparency or not
- Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**
- bgcindex* An integer which indicates color index that may be used to set the background color. Since the background color can now be an arbitrary RGB color (via MM_OPTION_BACKGROUND_COLOR), and this MM_OPTION_BACKGROUND_CINDEX value is only the nearest indexed color, not the actual color, this option is deprecated.
- Valid values: integers
 Default value: **1**
 Minimum: 1
 Maximum: 256
- bghcolor* A string which is the color name or 6-hexadecimal-digit RGB string for background color. Valid color names are described in the file \$SCHRODINGER/mmshare-vX.X/data/colors.res
- Valid values: text strings
 Default value: **000000**
- eye_sepf* Eye separation for stereo as fraction of normal. Actual separation may be made smaller to limit maximum stereo disparity to about 1 inch at clipping planes.
- Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 2.0
- fog* Set fog to be displayed, to not be displayed, or to toggle automatically based on the number of atoms in the workspace.

	Valid values:	auto on true yes off false no
	Default value:	auto
<i>fogatomlabels</i>	A Boolean option which determines whether to display fog on atom labels (only displayed if fog is already on in general).	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	true
<i>fogcutoff</i>	Specifies the minimum number of atoms in the workspace for which fog will be displayed, when fogging is set to auto.	
	Valid values:	integers
	Default value:	40
	Minimum:	0
<i>fogdensity</i>	Factor controlling optical density of fog.	
	Valid values:	reals
	Default value:	1.5
	Minimum:	0.0
	Maximum:	4.0
<i>fogendz</i>	Factor controlling end of linear fog ramp.	
	Valid values:	reals
	Default value:	1.25
	Minimum:	1.0
	Maximum:	2.0
<i>fogstartz</i>	Factor controlling start of linear fog ramp.	
	Valid values:	reals
	Default value:	0.05
	Minimum:	0.0
	Maximum:	0.999
<i>fogtype</i>	Specifies method used for calculating fog.	
	Valid values:	linear exp exp2 automatic
	Default value:	automatic
<i>fsaa</i>	This option is used to toggle full scene anti aliasing on the systems that support it.	

	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>invertstereo</i>	A boolean option which determines whether to swap stereo images
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>maximize</i>	A boolean option which determines whether to maximize the display window (on) or restore to the original window size (off)
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>opengl_level</i>	Sets an OpenGL level.
	Valid values: text strings
	Default value:
<i>perspective</i>	A Boolean option which when on enables perspective and when off uses orthogonal projection in the Workspace.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>perspectivescale</i>	Ratio between the actual size and the projected size of an object placed 1 foot behind display screen. This has a corresponding viewing distance, $\text{view_dist} = (1 \text{ ft}) / (\text{perspectivescale} - 1)$. For a known viewing distance, $\text{perspectivescale} = 1 + (1 \text{ ft} / \text{view_dist})$. Therefore, for a normal viewing distance of 2 feet, $\text{perspectivescale} = 1 + (1 \text{ ft} / 2 \text{ ft}) = 1.5$.
	Valid values: reals
	Default value: 1.35
	Minimum: 1.0
	Maximum: 4.0
<i>rubber_banding</i>	Specifies method used for drawing temporary graphics to track cursor motion (such as lasso rectangle or rubber band line). Valid values are “xor” or “main”. These draw to the front buffer in XOR mode, or draw to the back buffer at the end of main window drawing (before the background and foreground GL buffers are swapped).
	Valid values: xor main
	Default value: xor

Chapter 5: Commands

<i>sepf</i>	Image separation as fraction of maximum horizontal separation (window width minus image width) Valid values: reals Default value: 0 Minimum: 0.0 Maximum: 1.0
<i>sizef</i>	Image size as fraction of available width Valid values: reals Default value: 1 Minimum: 0.0 Maximum: 1.0
<i>start_str</i>	The system command to start hardware stereo. Valid values: text strings Default value:
<i>stereo</i>	A boolean option which determines whether to display in stereo Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>stereomethod</i>	Specifies method used for stereo viewing Valid values: hardware crosseyed walleyed fullscreen interlaced Default value: hardware
<i>stop_str</i>	The system command to stop hardware stereo. Valid values: text strings Default value:

distance

Specifies a pair of atoms to have their distance measured and displayed.

Syntax:

distance *surfacename1*= $\langle \text{text} \rangle$ *surfacename2*= $\langle \text{text} \rangle$
surfacex1= $\langle x \rangle$ *surfacex2*= $\langle x \rangle$ *surfacey1*= $\langle x \rangle$ *surfacey2*= $\langle x \rangle$
surfacez1= $\langle x \rangle$ *surfacez2*= $\langle x \rangle$ *xoffset*= $\langle x \rangle$ *yoffset*= $\langle x \rangle$
 $\langle \text{atom1} \rangle$ $\langle \text{atom2} \rangle$

Options:

surfacename1

Specifies the name of surface on which atom1 sits.

Valid values: text strings

Default value:

surfacename2

Specifies the name of surface on which atom2 sits.

Valid values: text strings

Default value:

surfacex1

Specifies the X coordinate for atom1 on surface.

Valid values: reals

Default value: **0**

surfacex2

Specifies the X coordinate for atom2 on surface.

Valid values: reals

Default value: **0**

surfacey1

Specifies the Y coordinate for atom1 on surface.

Valid values: reals

Default value: **0**

surfacey2

Specifies the Y coordinate for atom2 on surface.

Valid values: reals

Default value: **0**

surfacez1

Specifies the Z coordinate for atom1 on surface.

Valid values: reals

Default value: **0**

surfacez2

Specifies the Z coordinate for atom2 on surface.

Valid values: reals

Default value: **0**

xoffset

Specifies the X offset in Angstroms. Any distance created after this value is set will use this new value. Can be applied to an existing measurement by respecifying the measurement. Does not affect any already created distances.

Valid values: reals

Default value: **0**

yoffset Specifies the Y offset in Angstroms. Any distance created after this value is set will use this new value. Can be applied to an existing measurement by respecifying the measurement. Does not affect any already created distances.

Valid values: reals
Default value: **0**

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle$

The two atoms between which the distance is to be measured. Note that the specifying a-b is the same as specifying b-a.

distancecheck

Specifies two atoms which define a distance to be checked during a conformational search.

Syntax:

distancecheck *allowed*= $\langle x \rangle$ *maximum*= $\langle x \rangle$ *minimum*= $\langle x \rangle$
value= $\langle x \rangle \langle \text{atom1} \rangle \langle \text{atom2} \rangle$

Options:

allowed The maximum allowed variation from the “target” value during the distance check.

Valid values: reals
Default value: **0**
Minimum: 0.0

maximum The maximum distance for the distance check (Angstroms)

Valid values: reals
Default value: **0**
Minimum: 0.0

minimum The minimum distance for the distance check (Angstroms)

Valid values: reals
Default value: **0**
Minimum: 0.0

value The target value for the distance check (Angstroms)

Valid values: reals
Default value: **0**
Minimum: 0.0

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle$

The numbers of two atoms which are to be considered are to have their distance checked during a conformational search. Note that specifying a-b is the same as specifying b-a.

dockgridbox

Settings for the grid boxes for docking.

Syntax:

```
dockgridbox innersizex= $\langle x \rangle$  innersizey= $\langle x \rangle$  innersizez= $\langle x \rangle$ 
outersizex= $\langle x \rangle$  outersizey= $\langle x \rangle$  outersizez= $\langle x \rangle$ 
show=yes | no xcent= $\langle x \rangle$  ycent= $\langle x \rangle$  zcent= $\langle x \rangle$ 
```

Options:

*innersize*x The size of the inner box in X-coordinate.

Valid values: reals
 Default value: **5**
 Minimum: 0.0

*innersize*y The size of the inner box in Y-coordinate.

Valid values: reals
 Default value: **5**
 Minimum: 0.0

*innersize*z The size of the inner box in Z-coordinate.

Valid values: reals
 Default value: **5**
 Minimum: 0.0

*outersize*x The size of the outer box in X-coordinate.

Valid values: reals
 Default value: **10**
 Minimum: 0.0

*outersize*y The size of the outer box in Y-coordinate.

Valid values: reals
 Default value: **10**
 Minimum: 0.0

*outersize*z The size of the outer box in Z-coordinate.

Valid values: reals

	Default value: 10
	Minimum: 0.0
<i>show</i>	An option which determines whether to show grid boxes or not.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>xcent</i>	The X-coordinate of the box center.
	Valid values: reals
	Default value: 0
<i>ycent</i>	The Y-coordinate of the box center.
	Valid values: reals
	Default value: 0
<i>zcent</i>	The Z-coordinate of the box center.
	Valid values: reals
	Default value: 0

doplanealignment

This command actually does the plane alignment.

Syntax:

doplanealignment

dynamics

Used to set options associated with MacroModel molecular or stochastic dynamics procedures.

Syntax:

dynamics *equilibration*=⟨x⟩ *method*=molecular | stochastic
minimize_sampled=yes | no *sample*=⟨n⟩ *shake*=nothing |
hydrogens | all *simulation*=⟨x⟩ *temperature*=⟨x⟩
timestep=⟨x⟩

Options:

equilibration

The time used for equilibration for the equilibration part of the simulation in ps.

	Valid values: reals
	Default value: 1
	Minimum: 0.0
<i>method</i>	This option determines which dynamics method will be used.
	Valid values: molecular stochastic
	Default value: stochastic
<i>minimize_sampled</i>	Minimize the sampled structure.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>sample</i>	The number of structures to be sampled during the simulation.
	Valid values: integers
	Default value: 0
	Minimum: 0
<i>shake</i>	This option determines how the SHAKE algorithm will be applied during the simulation.
	Valid values: nothing hydrogens all
	Default value: nothing
<i>simulation</i>	The total time to be used in the simulation.
	Valid values: reals
	Default value: 10
	Minimum: 0.0
<i>temperature</i>	The temperature at which the simulation is to be run in Kelvin.
	Valid values: reals
	Default value: 300
	Minimum: 0.0
<i>timestep</i>	The timestep which is to be used for the simulation in fs.
	Valid values: reals
	Default value: 1.5
	Minimum: 0.0

ecalc

This is a standard alias for **calcenergy** (see [\[calcenergy\]](#), page 53).

editatomproperty

Sets the given property for the given atom in the atom-property table to the new value.

Syntax:

editatomproperty *column*= $\langle n \rangle$ *row*= $\langle n \rangle$ \langle new value \rangle

Options:

column The column number in the atom-property table.

Valid values: integers

Default value: **2**

Minimum: 2

row The row number in the atom-property table. If set to -1, it means all selected rows.

Valid values: integers

Default value: **1**

Minimum: -1

Operands:

\langle new value \rangle

The operand is the new value for the atom property.

ejob

This keyword is used to set various options associated with running Macro-Model jobs.

Syntax:

```
ejob host=⟨text⟩ incorporate=append | replace | ignore |
      appendungrouped | workspace input_file=⟨text⟩ job=⟨text⟩
      login=⟨text⟩ rce_structure_source=selected_entries | workspace
      | file structure_source=selected_entries | workspace | file
```

Options:

host The name of the host for the MacroModel job.

Valid values: text strings

Default value:

incorporate

How the results are to be incorporated into the project. This can be done with replacement of the existing entries or by appending as new entries to the project or by ignoring the final results.

Valid values: append
 replace
 ignore
 appendungrouped
 workspace

Default value: **append**

input_file The name of the structure input file.

Valid values: text strings

Default value:

job The name for the MacroModel job.

Valid values: text strings

Default value: **mmodtmp**

login The login name under which a MacroModel will be run.

Valid values: text strings

Default value:

rce_structure_source

Whether to use the selected entries in the current project, or a specified file with multiple structures as structure input for the job.

Valid values: selected_entries
 workspace
 file

Default value: **selected_entries**

structure_source

Whether to use the selected entries in the current project or what is in the workspace as input for the job.

Valid values:	selected_entries workspace file
Default value:	workspace

elementlabels

Deprecated and has no effect. It used to set the element label options. Now it is always off.

Syntax:

elementlabels *labelscheme*=heterohydro | hetero | off

Options:

labelscheme

Deprecated and always set to off. We used to supported three options for displaying element labels: 1) on all non-carbon atoms, including hydrogen. 2) on all non-carbon atoms, except hydrogen (i.e., not on hydrogen) 3) off (no element labels displayed).

Valid values:	heterohydro hetero off
---------------	------------------------------

Default value:	off
----------------	------------

embrace

Used to control the eBMrAcE feature of MacroModel.

Syntax:

embrace *energymode*=interaction | energydiff *inputfile*=⟨text⟩
ligandsmode=inputfile | selectedentries *ouput*=complexes |
ligands | all *receptor*=⟨text⟩ *receptormode*=firstitem | entry

Options:

energymode

This option determines which energy mode will be used by eBMrAcE.

	Valid values:	interaction energydiff
	Default value:	energydiff
<i>inputfile</i>	The name of the Embrace ligands input file.	
	Valid values:	text strings
	Default value:	
<i>ligandsmode</i>	This option determines which source of ligands will be used by eBMrAcE.	
	Valid values:	inputfile selectedentries
	Default value:	inputfile
<i>ouput</i>	This option determines what structural output will be produced.	
	Valid values:	complexes ligands all
	Default value:	complexes
<i>receptor</i>	This option specifies which entry in the project is to be treated as the receptor.	
	Valid values:	text strings
	Default value:	
<i>receptormode</i>	This option determines which kind of receptor will be used by eBMrAcE.	
	Valid values:	firstitem entry
	Default value:	firstitem

embracesearch

Defines settings for Embrace conformational searching in MacroModel.

Syntax:

```
embracecsearch eliminate=atom_deviation | rmsd
               max_distance_lmcs=<x> max_rmsd=<x> maxdist=<x>
               method=embrace_mcmm | embrace_lmcs | embrace_mixed_lmcs
               min_distance_lmcs=<x> numconfs=<n> numsteps=<n>
               numstructures=<n> probability_tors_lmcs=<x> window=<x>
```

Options:

- eliminate* The method to use for eliminating redundant conformers: maximum atom deviation or RMSD.
 Valid values: atom_deviation
 rmsd
 Default value: **atom_deviation**
- max_distance_lmcs* The maximum possible distance for the fastest moving atom in each LMCS move.
 Valid values: reals
 Default value: **6**
 Minimum: 0.0
- max_rmsd* Maximum RMSD for considering two structures equal.
 Valid values: reals
 Default value: **0.5**
 Minimum: 0.0
- maxdist* Maximum distance between atoms in equal structures.
 Valid values: reals
 Default value: **0.25**
 Minimum: 0.0
- method* This determines which method will be used to perform the Embrace conformational search.
 Valid values: embrace_mcmm
 embrace_lmcs
 embrace_mixed_lmcs
 Default value: **embrace_mcmm**
- min_distance_lmcs* The minimum possible distance for the fastest moving atom in each LMCS move.
 Valid values: reals
 Default value: **3**
 Minimum: 0.0
- numconfs* The number of input conformations available to seed each search.
 Valid values: integers

	Default value: 1
	Minimum: 0
<i>numsteps</i>	An option which sets the number of steps which will be performed during the Embrace conformational search.
	Valid values: integers
	Default value: 100
	Minimum: 0
<i>numstructures</i>	The number of structures that will be saved at the end of each search.
	Valid values: integers
	Default value: 1
	Minimum: 0
<i>probability_tors_lmcs</i>	The probability that a TORS/MOLS move will be made during an LMCS search.
	Valid values: reals
	Default value: 0.5
	Minimum: 0.0
	Maximum: 1.0
<i>window</i>	The energy window (in kJ/mol) within which structures will be saved.
	Valid values: reals
	Default value: 500
	Minimum: 0.0

emptyviewreset

Resets the viewing transform to empty workspace.

Syntax:

```
emptyviewreset
```

endundoblock

End the current undoable command block.

Syntax:

endundoblock

energykill

Specifies a running job to be killed. This will force the program running the job to be stopped.

Syntax:

energykill <job_id> | <job_name> | <job_name> <project_name>

Operands:

<job_id> | <job_name> | <job_name> <project_name>

The ID or name of the job which is to be put to killed. This must be the ID or name of a job which is currently running or sleeping.

Aliases:

kill (see [\[kill\]](#), page 310)

energymonitor

Specifies a job to monitor.

Syntax:

energymonitor <job_id> | <job_name> | <job_name>
 <project_name>

Operands:

<job_id> | <job_name> | <job_name> <project_name>

The ID or name of the job which is to be monitored. If no operand is specified then the job which was most recently monitored will be monitored.

Aliases:

monitor (see [\[monitor\]](#), page 346)

energysleep

Specifies a job to be put to sleep

Syntax:

energysleep $\langle \text{job_id} \rangle$ | $\langle \text{job_name} \rangle$ | $\langle \text{job_name} \rangle \langle \text{project_name} \rangle$

Operands:

$\langle \text{job_id} \rangle$ | $\langle \text{job_name} \rangle$ | $\langle \text{job_name} \rangle \langle \text{project_name} \rangle$

The ID or name of the job which is to be put to sleep. This must be the ID or name of a job which is currently running.

Aliases:

sleep (see [\[sleep\]](#), page 707)

energystart

Start a MacroModel job with the current energy settings.

Syntax:

energystart *onscreen*=yes | no $\langle \text{job_name} \rangle$

Options:

onscreen Determines whether the on-screen structure will be written and used as the .dat file for the MacroModel job.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

$\langle \text{job_name} \rangle$

The name of the job which is to be started. This operand is not the name of a file and should not contain any suffix. From this name the files “job_name”.com and “job_name”.mae will be created.

energystop

Specifies a running job to be stopped. The program running the job will be requested to stop but may take some time to respond to the command.

Syntax:

energystop <job_id> | <job_name> | <job_name> <project_name>

Operands:

<job_id> | <job_name> | <job_name> <project_name>

The ID or name of the job which is to be stopped. This must be the ID or name of a job which is currently running or sleeping.

Aliases:

stop (see [\[stop\]](#), page 709)

energytask

Determines which energy task is currently being set up.

Syntax:

energytask ecalc | mini | moldyn | drive | mult | csearch | mcsd | mintal
embrace | embracecsearch

Operands:

ecalc | mini | moldyn | drive | mult | csearch | mcsd | mintal | embrace | embracecsearch

The type of energy task which is to be set up. The operand cannot be abbreviated and must be given in full.

energyupdate

Specifies a running job to be updated.

Syntax:

energyupdate <job_id> | <job_name> | <job_name> <project_name>

Operands:

<job_id> | <job_name> | <job_name> <project_name>

The ID or name of the job which is to be put to updated. This is not the name of a file and no suffix is required. At present only MacroModel conformational searches will pay any attention to a request to update, however it is no error for an update to be issued for any type of job.

Aliases:

update (see [\[update\]](#), page 775)

energyupdatejobstatus

update the job status of all the jobs in the monitor panel.

Syntax:

energyupdatejobstatus

energywake

Specifies a sleeping job to be woken up.

Syntax:

energywake $\langle \text{job_id} \rangle$ | $\langle \text{job_name} \rangle$ | $\langle \text{job_name} \rangle \langle \text{project_name} \rangle$

Operands:

$\langle \text{job_id} \rangle$ | $\langle \text{job_name} \rangle$ | $\langle \text{job_name} \rangle \langle \text{project_name} \rangle$

The ID of the job which is to be put to woken up. This must be the ID of a job which is currently sleeping.

Aliases:

wake (see [\[wake\]](#), page 818)

enhance3d

Fit clipping planes to displayed structure. If display stereo is on, this may help prevent conflicts between the clipping range and the stereo angle that lead to clipping of the displayed structure or flatter appearance of the structure in stereo, especially when zooming in. If display perspective is on, this may help prevent conflicts between the front clipping plane and the viewing distance (which corresponds to the perspective scale) that lead to clipping of the displayed structure or reduced perspective scaling, especially when zooming in. If display fog is on, this helps maximize the shading contrast between near and far atoms and helps make the front part of the displayed structure more visible.

Syntax:

enhance3d \langle ASL \rangle

Operands:

\langle ASL \rangle

If present, we do enhance3d on ASL, if not present we do enhance3d on all atoms.

enterposeviewermode

This command enters into the pose viewer mode.

Syntax:

enterposeviewermode

entryadddisplay

Adds the given type of data to be displayed.

Syntax:

entryadddisplay \langle display properties \rangle

Operands:

\langle display properties \rangle

A bitwise combination of flags for data to be displayed.

entrycolorbyproperty

Colors the entry property values in Project Table based on the options.

Syntax:

entrycolorbyproperty *abovevalue*=⟨x⟩ *belowvalue*=⟨x⟩
colorbetween=solid | interpolate *exactvalue*=⟨x⟩
interpolatesteps=⟨n⟩ *useabove*=yes | no *usebelow*=yes | no
usebetween=yes | no *useexact*=yes | no ⟨property_name⟩

Options:

abovevalue

Valid values: reals
 Default value: **10**

belowvalue

Valid values: reals
 Default value: **1**

colorbetween

Valid values: solid
 interpolate
 Default value: **interpolate**

exactvalue

Valid values: reals
 Default value: **10**

interpolatesteps

Valid values: integers
 Default value: **10**

useabove

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

usebelow

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

usebetween

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

useexact

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

Operands:

⟨property_name⟩

Indicates the property name for which color applied to the property value in Project Table.

entrycombine

Combine selected project entries into a single entry.

Syntax:

```
entrycombine order=entryid | row replace=yes | no <new_name>  
                <ESL>
```

Options:

<i>order</i>	This determines whether the entries should be combined based on entryid order or table row order in the project table. Valid values: entryid row Default value: entryid
<i>replace</i>	This determines whether the combined entry is allowed to replace an existing entry having the name <new_name>. If the name of the combined entry matches the name of an entry in the workspace, the entry in the workspace is also replaced. If the replace option is off, the user will be asked to accept a unique name based upon <new_name>. Valid values: boolean (true false; yes no; y n; on off) Default value: false

Operands:

<new_name> <ESL>

The name to be given to the single combined project entry, followed by a valid ESL expression. If no name is specified, no rename will be done. The ESL expression specifies which entries are to be combined.

entrycopyprop

Set property values for selected project entries from another property. The to property value will be copied from the from property for the entries selected by the ESL expression. If the to property does not already exist, it will be created. When the data type for the two properties is not the same, a reasonable attempt is made to convert between them. If a selected entry has no value for the from property, or the conversion from a string value fails, the value is generally cleared for the to value for that entry. Entry names are never cleared.

Syntax:

entrycopyprop *from*= \langle text \rangle *to*= \langle text \rangle \langle ESL \rangle

Options:

<i>from</i>	The name of the property which provides the values to be copied. This can be either the user name or the m2io data name for the property. Valid values: text strings Default value:
<i>to</i>	The name of the property to be modified. This can be either the user name or the m2io data name for the property, if the to property exists. If the to property does not exist, it will be created using the data type of the from property, with user as the author. If the property is the entry name (e.g. Entry Name or s_m_entry_name), an entryrename with replace=no will be done. If the property is the 'included in Workspace' property (In or b_m_entry_in_workspace), then entrywsinclude or entrywsexclude will be done, if needed. Valid values: text strings Default value:

Operands:

\langle ESL \rangle

The ESL expression specifies for which entries the property values are to be copied.

entrycreatesmoothed

Create new smoothed entries using the entries specified in the ESL. The entries from the ESL are duplicated and the intermediates are placed in between these duplicates. All these new entries are selected.

Syntax:

entrycreatesmoothed *intermediates*= \langle n \rangle \langle ESL \rangle

Options:

intermediates

Specifies the number of intermediate entries to create between each entry pair. The first half of the intermediates are duplicated

from the current key entry and the second half from the next key entry. If this number is odd, then the second half gets the additional entry.

Valid values: integers
Default value: **48**
Minimum: 1

Operands:

⟨ESL⟩

A valid ESL specification.

entrydelete

Delete each selected project entry.

Syntax:

entrydelete ⟨ESL⟩

Operands:

⟨ESL⟩

A valid ESL specification. Delete those entries which match the ESL description.

entrydisassociate

For each selected project entry, create an entry for each molecule except water, and one for all waters in the entry. Each of these new entries will be given a unique name, based upon the name of its originating entry. Any of the originating entries that are in the workspace are replaced by their disassociated molecule entries.

Syntax:

entrydisassociate ⟨ESL⟩

Operands:

⟨ESL⟩

A valid ESL specification. Disassociate those entries which match the ESL description.

entrydisassociatebychain

Create a new entry for each chain in each selected project entry. Each of these new entries will be given a unique name, based upon the name of its originating entry. Any of the originating entries that are in the workspace are replaced by their disassociated chain entries.

Syntax:

entrydisassociatebychain \langle ESL \rangle

Operands:

\langle ESL \rangle

A valid ESL specification. Disassociate those entries which match the ESL description.

entrydisassociatebyligand

Create a new entry for each chain in each selected project entry. Each of these new entries will be given a unique name, based upon the name of its originating entry. Any of the originating entries that are in the workspace are replaced by their disassociated chain entries. For each selected project entry, create an entry for each ligand-sized molecule, an entry for all waters in the selected entry, and an entry for all other atoms in the selected entry. Each of these new entries will be given a unique name, based upon the name of its originating entry. Any of the originating entries that are in the workspace are replaced by their disassociated molecule entries.

Syntax:

entrydisassociatebyligand \langle ESL \rangle

Operands:

\langle ESL \rangle

A valid ESL specification. Disassociate those entries which match the ESL description.

entrydisplaysurfaces

Display surfaces for the entry.

Syntax:

entrydisplaysurfaces $\langle \text{entry_index} \rangle$

Operands:

$\langle \text{entry_index} \rangle$

The index of entry.

entrydragselection

Move current entry selection to the specified table row.

Syntax:

entrydragselection $\langle \text{row} \rangle \langle \text{table} \rangle$

Operands:

$\langle \text{row} \rangle \langle \text{table} \rangle$

The destination row number. In place of the row number, “top” may be used to specify the first row of the table (row 1) and “bottom” may be used to specify placement after the last unselected row ($\text{last_row} + 1$). The selected entries will be moved as a block, placing the first selected entry below the first unselected entry which is above the destination row. If there is no unselected entry above the destination row, the selected entries will be moved to the top of the table. The name of the table to use as the source for the row number. If the table operand is missing, the current or default table (1) will be used.

entryduplicate

Create a new duplicate entry for each selected project entry. Each of these new entries will be given a unique name, based upon the name of its originating entry.

Syntax:

entryduplicate $\langle \text{ESL} \rangle$

Operands:

$\langle \text{ESL} \rangle$

A valid ESL specification. Duplicate those entries which match the ESL description.

entryexport

Export selected entries from the current project, or workspace entries, to a file or a number of files. If the export command is issued with options, but without a filename, then nothing is written but the options are updated. An export command with both options and a file name specified will result in the export being performed with the new options to the specified file(s). When exporting each entry to an individual file, the given file name will be the base file name which is used to name those individual files by adding their associated postfixes to the base file name.

Syntax:

```
entryexport append=yes | no displayedonly=yes | no files=single
| individual format=mmod | pdb | mol2 | maestro | sd |
jaguarinput | gaussianinput | biograp | xyz | jaguaroutput |
gaussian92 | gaussian94 | gamess | mopaccartesian |
mopacinternal | mopacoutput | babelpdb | mdl | babelmol |
babelmol2 | spartan | spartansemi | spartanmm | gamessinput
| gausscartesian | gaussianz | jaguarzmatrix | jaguarcartesian
| any | reagentprep | mcproinput | cms | maegz |
sdcompressed | glideposeviewer | smiles | smilesgz | common
graphical=yes | no names=withentry | withnumber | entry |
withentryid | title properties=all | subset
reorderbyresnum=yes | no sdfv3000=yes | no
seqheaders=yes | no source=workspace | selected
writedupconectrecs=yes | no [{filename}]
```

Options:

append This option determines whether to append to the file which is going to be written. It will be grayed out when PDB format is selected.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

displayedonly

If this option is set to true, then only the displayed atoms in the Workspace will be exported.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

Chapter 5: Commands

files This option sets the number of export files. The two options are: export entries to a single file (0 default) or export each entry to an individual file.

Valid values: single
 individual
Default value: **single**

format This option sets the format of the file to be written. Valid values are “mmod”, “pdb”, “mol2”, “maestro” or “sd”.

Valid values: mmod
 pdb
 mol2
 maestro
 sd
 jaguarinput
 gaussinput
 biograf
 xyz
 jaguaroutput
 gaussian92
 gaussian94
 gamess
 mopaccartesian
 mopacinternal
 mopacoutput
 babelpdb
 mdl
 babelmol
 babelmol2
 spartan
 spartansemi
 spartanmm
 gamessinput
 gausscartesian
 gaussianz
 jaguarzmatrix
 jaguarcartesian
 any
 reagentprep
 mcproinput
 cms
 maegz
 sdcompressed
 glideposeviewer
 smiles
 smilesgz
 common

Default value: **maestro**

graphical This option determines whether to export graphical information when an entry is being saved to disk.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

Chapter 5: Commands

<i>names</i>	<p>This option sets output file names when a number of files are exported. There are three options: file name + entry name, file name + number, just entry name. This option will be grayed out if only one file is exported.</p> <p>Valid values: withentry withnumber entry withentryid title</p> <p>Default value: withentry</p>
<i>properties</i>	<p>This option sets the properties to be exported. The value can either be all or subset .</p> <p>Valid values: all subset</p> <p>Default value: all</p>
<i>reorderbyresnum</i>	<p>This option determines whether to reorder the records by residue number while writing pdb files to disk</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: true</p>
<i>sdfv3000</i>	<p>If set, then SDF files will be written in v3000 format.</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>
<i>seqheaders</i>	<p>This option determines whether to export secondary structure information and seqres headers when an entry is being written in pdb format to disk</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>
<i>source</i>	<p>Whether to export selected entries in the current project table or to export what is in the workspace, including scratch entry.</p> <p>Valid values: workspace selected</p> <p>Default value: selected</p>
<i>writedupconectrecs</i>	<p>This option determines whether to write duplicate CONECT record information while writing pdb files to disk.</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: true</p>

Operands:

[*<filename>*]

The name of the file to which entries will be written. If no name is specified, then no export will be done.

entryexportcanvas

Export selected entries from the current project to a canvas project.

Syntax:

entryexportcanvas *canvasimport*=*<text>* *<canvas_project>*

Options:

canvasimport

This sets import option for canvas, if the canvas project already exists. The option can be append or replace . If blank, this defaults to append .

Valid values: text strings

Default value:

Operands:

<canvas_project>

The path of the canvas project to which selected entries will be exported. If no path is specified, then no export will be done.

entryexportspreadsheet

Export selected entries from the current project to a file for use in a spreadsheet. This can be comma-separated value (.csv) format or tab-delimited format.

Syntax:

entryexportspreadsheet *columns*=all | subset *delimiter*=*<text>*
rows=all | subset | selected | included *<filename>*

Options:

columns This option sets the columns to be output. The columns can either be all or subset .

	Valid values:	all subset
	Default value:	all
<i>delimiter</i>	This option sets the delimiter to use to separate columns.	
	Valid values:	text strings
	Default value:	,
<i>rows</i>	This option sets the rows to be written out. The rows can be all , subset , selected , or included .	
	Valid values:	all subset selected included
	Default value:	all

Operands:

⟨ filename ⟩

The name of the file to which entries will be written. If no name is specified, then no export will be done.

entryextendselection

Extend current entry selection to encompass the specified entry. This is slower than using entryextendselectrow, because of the need to locate the entry row, given the entry name.

Syntax:

entryextendselection ⟨ entry_name ⟩

Operands:

⟨ entry_name ⟩

The name of an entry in the current project. A range selection will be done between currently selected entries and the specified entry.

entryextendselectrow

Extend current entry selection to encompass the specified entry table row.

Syntax:

entryextendselectrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

A row number between 1 and the total number of rows in the project. A range selection will be done for the entries between currently selected entry rows and the specified entry row number.

entryextendwsinclude

Extends the included entries to encompass the selected project entry.

Syntax:

entryextendwsinclude $\langle \text{ESL} \rangle$

Operands:

$\langle \text{ESL} \rangle$

A valid ESL specification. Extends included entries to encompass the entries which match the ESL description into the workspace.

entrygroupbyproperty

Create new groups, if required; and re-group the entries according to their property values.

Syntax:

entrygroupbyproperty *alignon*= $\langle x \rangle$ *basename*= $\langle \text{text} \rangle$
maximum= $\langle x \rangle$ *minimum*= $\langle x \rangle$ *stepsize*= $\langle x \rangle$
usealign=yes | no *usemaximum*=yes | no
useminimum=yes | no $\langle \text{property_name} \rangle$ $\langle \text{ESL} \rangle$

Options:

alignon

Valid values:	reals
Default value:	0
Minimum:	0.0

Chapter 5: Commands

basename

Valid values: text strings
Default value: **group**

maximum

Valid values: reals
Default value: **0**

minimum

Valid values: reals
Default value: **0**

stepsize

Valid values: reals
Default value: **1**
Minimum: 0.0

usealign

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

usemaximum

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

useminimum

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

Operands:

⟨property_name⟩ ⟨ESL⟩

entrygroupcollapse

Collapses the given group.

Syntax:

entrygroupcollapse *isexp*=yes | no ⟨group_name⟩

Options:

isexp Specifies whether the <group_name> is actual group name or it is expression such as 'selected'.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

Operands:

`<group_name>`

Group name. Collapses the given group and hides all the entries in that group.

entrygroupcountcontacts

Record contact counts for entries in a group, using the first group entry as the reference structure (receptor).

Syntax:

entrygroupcountcontacts `<group_name>`

Operands:

`<group_name>`

An existing group name. The first entry in the group is treated as the reference, for measuring contacts. Its contact count property values are deleted, and the counts for all other entries in the group are set.

entrygroupcounthbonds

Record H-bond counts for entries in a group, using the first group entry as the reference structure (receptor).

Syntax:

entrygroupcounthbonds `<group_name>`

Operands:

`<group_name>`

An existing group name. The first entry in the group is treated as the reference, for measuring H-bonds. Its H-bond count property value is deleted, and the count for all other entries in the group is set.

entrygroupcreate

Creates a group for given entries.

Syntax:

entrygroupcreate <group_name> <ESL>

Operands:

<group_name> <ESL>

The new name to be given to the group of entries that match given ESL expression. The name should be unique. The new group name should be followed by a valid ESL expression to specify which entries are to be grouped. If no entry matches the ESL expression, no new group will be created. No empty group can exists.

entrygroupdelete

Deletes the given group and all the entries in it.

Syntax:

entrygroupdelete *isexp*=yes | no <group_name>

Options:

isexp Specifies whether the <group_name> is actual group name or it is expression such as 'selected'.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

<group_name>

Name of the group to be deleted. Deletes the given group and all the entries in that group.

entrygroupdragselection

Moves the selected entry groups to the specified table row.

Syntax:

entrygroupdragselection $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

the destination row number. In place of the row number, “top” may be used to specify the first row of the table (row 1) and “bottom” may be used to specify placement after the last row (`last_row + 1`). If the target row corresponds to an ungrouped entry, the group will be placed just after that last ungrouped entry section. If the target row corresponds to a grouped entry, the group will be placed just before the group to which the target entry belongs to. All the selected groups will be moved as a block.

entrygroupduplicate

Creates a new duplicate entry for each entry of the source group. Each of these new entries will be given a unique name, based upon the name of its originating entry. All these newly created entries will be grouped under a new group; and that new group will be given a unique name, based upon the name of source group.

Syntax:

entrygroupduplicate *isexp*=yes | no $\langle \text{group_name} \rangle$

Options:

isexp Specifies whether the $\langle \text{group_name} \rangle$ is actual group name or it is expression such as 'selected'.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

$\langle \text{group_name} \rangle$

The name of the group to be duplicated. Duplicates the given group along with all its entries.

entrygroupexpand

Expands the given group.

Syntax:

entrygroupexpand *isexp*=yes | no <group_name>

Options:

isexp Specifies whether the <group_name> is actual group name or it is expression such as 'selected'.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

<group_name>

Group name. Expands the given group and displays all the entries in that group in PT.

entrygroupexpandonly

Expands the given group and collapses all other groups

Syntax:

entrygroupexpandonly *isexp*=yes | no <group_name>

Options:

isexp Specifies whether the <group_name> is actual group name or it is expression such as 'selected'.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

<group_name>

Group name. Expands the given group, display all the entries in that group in PT and collapses all the other groups.

entrygroupextendselection

Extend current entry selection to encompass all the entries in the specified group.

Syntax:

entrygroupextendselection <group_name>

Operands:

<group_name>

The name of a entry group in the current project. A range selection will be done between currently selected entries and the specified group; all the entries in that group will also be selected.

entrygroupinvertselection

Invert the selection of entries of an existing group.

Syntax:

entrygroupinvertselection *isexp*=yes | no <group_name>

Options:

isexp Specifies whether the <group_name> is actual group name or it is expression such as 'selected'.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

<group_name>

An existing group name. Inverts the selection state of entries of the given group.

entrygroupmove

Moves the given entry group to the specified table row.

Syntax:

entrygroupmove <group_name> <row>

Operands:

<group_name> <row>

Name of the group to be moved. Group name and the destination row number. In place of the row number, “top” may be used to specify the first row of the table (row 1) and “bottom” may be used to specify placement after the last row (last_row + 1). If the target row corresponds to an ungrouped entry, the group will be placed just after that last ungrouped entry section. If the target row corresponds to a grouped entry, the group will be placed just before the group to which the target entry belongs to. All the entries of the group will be moved as a block.

entrygrouprename

Rename the given group.

Syntax:

```
entrygrouprename <old_name> <new_name>
```

Operands:

<old_name> <new_name>

The name of the group which we want to rename followed by the new name for the group.

entrygroupselect

Select all the entries of given group.

Syntax:

```
entrygroupselect isexp=yes | no <group_name>
```

Options:

isexp Specifies whether the <group_name> is actual group name or it is expression such as 'selected'.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

<group_name>

Group name. Selects all the entries of given group does not change the selection state of the other entries.

entrygroupselectonly

Selects the all entries of the given group and unselects all other entries.

Syntax:

```
entrygroupselectonly <group_name>
```

Operands:

<group_name>

Group name. Selects all the entries of given group and unselects all other entries.

entrygroupsettitle

Set the title for the given group.

Syntax:

```
entrygroupsettitle <group_name> <title>
```

Operands:

<group_name> <title>

The name of the group whose title we want to set followed by the title.

entrygroupungroup

Ungroups the entries in the given group.

Syntax:

```
entrygroupungroup isexp=yes | no <group_name>
```

Options:

isexp Specifies whether the <group_name> is actual group name or it is expression such as 'selected'.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

`<group_name>`

`<group_name>` Remove the group and move all the entries in that group to the end of ungrouped section. In project table, all the ungrouped entries (if any) will be present at the top i.e. before all the groups.

entrygroupunselect

Unselects all the entries of given group.

Syntax:

entrygroupunselect *isexp*=yes | no `<group_name>`

Options:

isexp Specifies whether the `<group_name>` is actual group name or it is expression such as 'selected'.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

`<group_name>`

Group name. Unselects all the entries of given group and does not change the selection state of the other entries.

entrygroupviewposes

Prepare to view group poses (using eplayer), using the first group entry as the reference structure (receptor).

Syntax:

entrygroupviewposes *contactdisplay*=yes | no
hbonddisplay=yes | no `<group_name>`

Options:

contactdisplay

This option determines whether Contact markers will be displayed for entry group view poses.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

hbonddisplay

This option determines whether H-bond markers will be displayed for entry group view poses.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

Operands:

⟨group_name⟩

An existing group name. The first entry in the group is treated as the reference, for viewing H-bonds and contacts, and is fixed in the Workspace so that it remains while stepping through other group entries with the eplayer.

entrygroupsexclude

Exclude all the entries of given group from the workspace.

Syntax:

entrygroupsexclude *isexp*=yes | no ⟨group_name⟩

Options:

isexp Specifies whether the <group_name> is actual group name or it is expression such as 'selected'.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

Operands:

⟨group_name⟩

A valid group name. Excludes those entries which belong to the given group from the workspace.

entrygroupsinclude

Include all entries of given group into the workspace.

Chapter 5: Commands

Syntax:

entrygroupwsinclude *isexp*=yes | no <group_name>

Options:

isexp Specifies whether the <group_name> is actual group name or it is expression such as 'selected'.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

<group_name>

A valid group name. Include all the entries which belong to the given group.

entrygroupwsincludeonly

Include only the entries of given group into the workspace. Exclude all others.

Syntax:

entrygroupwsincludeonly *isexp*=yes | no <group_name>

Options:

isexp Specifies whether the <group_name> is actual group name or it is expression such as 'selected'.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

<group_name>

A valid group name. Include only those entries which belong to the given group and exclude all other entries.

entryimport

Import structures into the current project.

Syntax:

```
entryimport all=yes | no allfiles=yes | no creategroups=multiple
| all | none end=yes | no extradata=yes | no
firstonly=yes | no fittoscreen=yes | no format=mmod | pdb
| mol2 | maestro | sd | jaguarinput | gaussianinput | biograf |
xyz | jaguaroutput | gaussian92 | gaussian94 | gamess |
mopaccartesian | mopacinternal | mopacoutput | babelpdb |
mdl | babelmol | babelmol2 | spartan | spartansemi |
spartanmm | gamessinput | gausscartesian | gaussianz |
jaguarzmatrix | jaguarcartesian | any | reagentprep |
mcproinput | cms | maegz | sdcompressed | glideposeviewer |
smiles | smilesgz | common glideposeviewer=yes | no
graphical=yes | no hidewarnings=yes | no
openprepwizard=yes | no readalternate=yes | no
sdttitle=moleculename | propertyname
sdttitlepropertyname=<text> start=<n> total=<n>
wsinclude=none | first | all wsreplace=yes | no [<filename>]
```

Options:

- | | |
|---------------------|--|
| <i>all</i> | This determines whether all structures will be imported, or just a specified range.
Valid values: boolean (true false; yes no; y n; on off)
Default value: true |
| <i>allfiles</i> | This determines whether all files in the directory given in the operand will be imported.
Valid values: boolean (true false; yes no; y n; on off)
Default value: false |
| <i>creategroups</i> | This option determines which of the imported structures are to be grouped. Valid values are “multiple”, “all” or “none”.
Valid values: multiple
all
none
Default value: multiple |
| <i>end</i> | This determines if all structures in the file are to be imported starting from the structure specified by start.
Valid values: boolean (true false; yes no; y n; on off)
Default value: false |
| <i>extradata</i> | This option determines whether to import data files associated with imported structures. The additional data file associations |

Chapter 5: Commands

are expected to be listed in .smap files having the same base name as the structure files.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

firstonly This option imports only the first structure incase the file format is of Desmond (.cms) type.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

fittoscreen This determines whether the newly imported and included structures should be fit to Workspace or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

format This option sets the format of the file to be read. Valid values are “maestro”, “mmod”, “pdb”, “mol2”, or “sd”.

Valid values:	mmod pdb mol2 maestro sd jaguarinput gaussinput biograf xyz jaguaroutput gaussian92 gaussian94 gamess mopaccartesian mopacinternal mopacoutput babelpdb mdl babelmol babelmol2 spartan spartansemi spartanmm gamessinput gausscartesian gaussianz jaguarzmatrix jaguarcartesian any reagentprep mcproinput cms maegz sdcompressed glideposeviewer smiles smilesgz common
Default value:	common

glideposeviewer

This option determines whether to setup pose viewer files for viewing poses with the ePlayer.

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

Chapter 5: Commands

- graphical* This determines whether the file will be read with just the geometric information or if the visual information in the file will be read and used for displaying
- Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**
- hidewarnings* This flag turns off the Maestro warnings about invalid structures on import.
- Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**
- openprepwizard* This option determines whether to open protein prepwizard panel after importing PDB files.
- Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**
- readalternate* This option determines whether to read alternate locations when importing PDB files.
- Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**
- sdttitle* The source of titles from a SD format file. Valid values are “moleculename” or “propertyname” .
- Valid values: moleculename
propertyname
Default value: **moleculename**
- sdttitlepropertyname* The property to be used to construct titles from a SD format file, if sdttitle option is “propertyname” .
- Valid values: text strings
Default value:
- start* This option sets the number of the first structure to be imported, if not importing all.
- Valid values: integers
Default value: **1**
Minimum: 1
- total* The total number of structures to be imported from the file, if not importing all structures.
- Valid values: integers
Default value: **1**
Minimum: 1

- wsinclude* This option determines which of the imported structures are to be included in the workspace. Valid values are “none”, “first”, or “all”.
- Valid values: none
 first
 all
- Default value: **first**
- wsreplace* This determines whether the structures currently in the workspace will be replaced by the included imported structures, or whether they will be included in the workspace also.
- Valid values: boolean (true|false; yes|no; y|n; on|off)
- Default value: **true**

Operands:

[<filename>]

The name of the file from which structures will be imported. If no name is specified, then no import will be done.

entryimportspreadsheet

Import entries from a file into the project table.

Syntax:

```
entryimportspreadsheet applyto=selected | all delimiter=comma
| tab | userdefined userdelimiter=<text> <filename>
<import_key> <proj_key>
```

Options:

- applyto* This option sets whether the data needs to be imported to all or selected entries.
- Valid values: selected
 all
- Default value: **selected**
- delimiter* This option sets the delimiter to use to separate columns.
- Valid values: comma
 tab
 userdefined
- Default value: **comma**

userdelimiter

This option sets the delimiter defined by the user to separate columns.

Valid values: text strings

Default value:

Operands:

⟨ filename ⟩ ⟨ import_key ⟩ ⟨ proj_key ⟩

The name of the file from which entries will be read. The key to be considered as the reference from file. Property in project that maps to import_key.

entryimportvibration

Import vibration data for the given entry

Syntax:

entryimportvibration *entry*=⟨ text ⟩ ⟨ filename ⟩

Options:

entry This is the entry which the vibration data will be attached to.

Valid values: text strings

Default value:

Operands:

⟨ filename ⟩

The name of the file from which vibration data will be imported.

entryinvertselection

Invert the selection state of all project entries.

Syntax:

entryinvertselection

entrymergeprop

Merge properties from a destination, appending to a source property. The from property values will be appended to the to property for the entries selected by the ESL expression. If the to property does not already exist, it will be created. The to property must be a string type. The separator value will be placed between the to and each from property during the merge operation.

Syntax:

entrymergeprop *from*=⟨text⟩ *separator*=⟨text⟩ *to*=⟨text⟩ ⟨ESL⟩

Options:

<i>from</i>	<p>The name of the property(s) which provides the values to be merged. This can be either the user names or the m2io data names for the property.</p> <p>Valid values: text strings</p> <p>Default value:</p>
<i>separator</i>	<p>The string that will be placed between each value during the merge operation.</p> <p>Valid values: text strings</p> <p>Default value: :</p>
<i>to</i>	<p>The name of the property to be modified. This can be either the user name or the m2io data name for the property, if the to property exists. If the to property does not exist, it will be created using the data type of the from property, with user as the author. The to property must be a string type.</p> <p>Valid values: text strings</p> <p>Default value:</p>

Operands:

⟨ESL⟩

The ESL expression specifies for which entries the property values are to be merged.

entrymovetogroup

Move the given entries to an existing group.

Syntax:

entrymovetogroup <group_name> <ESL>

Operands:

<group_name> <ESL>

An existing group name and a valid ESL expression. Moves the entries that match the ESL description to the given group. The entries are removed from the previous groups.

entryremovedisplay

Removes the given type of data to be displayed.

Syntax:

entryremovedisplay <display properties>

Operands:

<display properties>

A bitwise combination of flags for data to be displayed.

entryrename

Rename selected project entries.

Syntax:

entryrename *replace*=yes | no <new_name> <ESL>

Options:

replace This determines whether a single selected entry is allowed to replace an existing entry having the name <new_name>. Otherwise, the user will be asked to accept a unique name based upon <new_name>.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

<new_name> <ESL>

The new name to be given to a single selected project entry, or the basename for multiple selected project entries from which unique names will be derived for each. The new name should be followed by a valid ESL expression to specify which entries are to be renamed. If no name is specified, no rename will be done.

entryresetcolor

Reset background color of selected project entries to default.

Syntax:

```
entryresetcolor <ESL>
```

Operands:

<ESL>

A valid ESL specification. Reset the background color of those entries which match the ESL description to default.

entryselect

Select specified entries in current project.

Syntax:

```
entryselect <ESL>
```

Operands:

<ESL>

A valid ESL specification. Selects those entries which match the ESL description, adding to the currently selected entries.

entryselectall

Select all project entries.

Syntax:

entryselectall

entryselectonly

Select only the specified entries in current project.

Syntax:

entryselectonly \langle ESL \rangle

Operands:

\langle ESL \rangle

A valid ESL specification. Selects only those entries which match the ESL description, unselecting all other entries.

entryselectonlyrow

Select only the entry corresponding to the given row in the current project.

Syntax:

entryselectonlyrow \langle row \rangle

Operands:

\langle row \rangle

A row number between 1 and the total number of rows. Selects only the entry corresponding to the given row, unselecting all other entries.

entryselectrandom

Selects a random set of entries.

Syntax:

entryselectrandom *entries*=selected | all *numentries*= \langle n \rangle
percentage= \langle x \rangle *selectoption*=nentries | npercentofentries None

Options:

<i>entries</i>	This determines whether the random entries will be chosen from the selected entries or from all entries. Valid values: selected all Default value: selected
<i>numentries</i>	This is the number of entries to select. Valid values: integers Default value: 1 Minimum: 0
<i>percentage</i>	This is the percentage of entries to select. Valid values: reals Default value: 50 Minimum: 0.0 Maximum: 100.0
<i>selectoption</i>	This determines whether the percentage of entries has to be selected or the number of entries has to be selected. Valid values: nentries npercentofentries Default value: npercentofentries

Operands:

None

None

entryselectrow

Select the entry corresponding to the given row in the current project.

Syntax:

entryselectrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

A row number between 1 and the total number of rows. Selects the entry corresponding to the given row.

entrysetcolor

Set background color for selected project entries.

Syntax:

```
entrysetcolor blue=<n> green=<n> red=<n> <ESL>
```

Options:

<i>blue</i>	Blue value.	
	Valid values:	integers
	Default value:	0
	Minimum:	0
	Maximum:	255
<i>green</i>	Green value.	
	Valid values:	integers
	Default value:	0
	Minimum:	0
	Maximum:	255
<i>red</i>	Red value.	
	Valid values:	integers
	Default value:	0
	Minimum:	0
	Maximum:	255

Operands:

<ESL>

A valid ESL specification. Set the background color of those entries which match the ESL description with given RGB value.

entrysetdeletable

Allow deletion of selected project entries.

Syntax:

```
entrysetdeletable <ESL>
```

Operands:

<ESL>

A valid ESL specification. Allow deletion of entries that match the ESL description.

entrysetnondeletable

Disallow deletion of selected project entries.

Syntax:

entrysetnondeletable \langle ESL \rangle

Operands:

\langle ESL \rangle

A valid ESL specification. Disallow deletion of entries that match the ESL description.

entrysetprop

Set property value for selected project entries.

Syntax:

entrysetprop *property*= \langle text \rangle *value*= \langle text \rangle \langle ESL \rangle

Options:

property The name of the property to be modified. This can be either the user name or the m2io data name for the property. If the property is the entry name (e.g. Entry Name or s.m.entry_name), an entryrename with replace=no will be done. If the property is the 'included in Workspace' property (In or b.m.entry_in_workspace), then entrywsinclude or entrywsexclude will be done, if needed.

Valid values: text strings

Default value:

value The value to be set for the selected entries. For Boolean properties, legal values are yes, no, true, false, on, off, and 1, 0.

Valid values: text strings

Default value:

Operands:

⟨ESL⟩

The specified property value will be assigned to the entries selected by the ESL expression.

entrysetreadonly

Disallow modification of selected project entries.

Syntax:

entrysetreadonly ⟨ESL⟩

Operands:

⟨ESL⟩

A valid ESL specification. Disallow modification to entries that match the ESL description.

entrysettitle

Set entry title.

Syntax:

entrysettitle ⟨title⟩ ⟨ESL⟩

Operands:

⟨title⟩ ⟨ESL⟩

New entry title and a valid ESL specification. Set title for those entries which match the ESL description from the workspace.

entrysetwritable

Allow modification of selected project entries.

Syntax:

entrysetwritable $\langle \text{ESL} \rangle$

Operands:

$\langle \text{ESL} \rangle$

A valid ESL specification. Allow modification to entries that match the ESL description.

entryshowall

Creates a entry subset in the project table consisting of all the entries in project.

Syntax:

entryshowall

entrystars

Sets the number of stars for the given entries

Syntax:

entrystars *stars*= $\langle n \rangle$ $\langle \text{ESL} \rangle$

Options:

stars The number of stars to set for the given entries.

Valid values: integers

Default value: **1**

Minimum: 0

Operands:

$\langle \text{ESL} \rangle$

Which entries to set stars for Sets the number of stars for the given entries.

entrytable

This keyword is used to set various options associated with the project entry table (Project Table).

Syntax:

entrytable *swapnametitle*=yes | no

Options:

swapnametitle

This determines whether the table is shown and otherwise treated as if the entry name and title columns are swapped with respect to their positions in the project table.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

entryundisplaysurfaces

Undisplay surfaces for the entry.

Syntax:

entryundisplaysurfaces \langle entry_index \rangle

Operands:

\langle entry_index \rangle

The index of entry.

entryunselect

Unselect specified entries in current project.

Syntax:

entryunselect \langle ESL \rangle

Operands:

\langle ESL \rangle

A valid ESL specification. Unselect those entries which match the ESL description.

entryunselectall

Unselect all project entries.

Syntax:

```
entryunselectall
```

entryunselectrow

Unselect the entry corresponding to the given row.

Syntax:

```
entryunselectrow <row>
```

Operands:

<row>

A row number between 1 and the total number of rows in the project. Unselects the entry corresponding to the given row number.

entrywatermapexaminerresults

Attempt to view WaterMap results for the specified entry, if it has WaterMap files.

Syntax:

```
entrywatermapexaminerresults <entry_name>
```

Operands:

<entry_name>

The name (entry ID) of the entry, whose WaterMap data is to be viewed.

entrywscreate

Create project entry from atoms in the workspace.

Syntax:

entrywsc*create* *replace*=yes | no <new_name> <ASL>

Options:

replace This option is ignored starting with version 60105 of Maestro.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

Operands:

<new_name> <ASL>

The name to be given to the created project entry, followed by a valid ASL expression to specify which atoms constitute the new project entry. If no name is specified, no operation will be done.

entrywsexclude

Exclude selected project entries from the workspace.

Syntax:

entrywsexclude <ESL>

Operands:

<ESL>

A valid ESL specification. Excludes those entries which match the ESL description from the workspace.

entrywsexcludeunfixed

Exclude all unfixed entries from the workspace.

Syntax:

entrywsexcludeunfixed

entrywsinclude

Include selected project entries into the workspace.

Syntax:

entrywsinclude \langle ESL \rangle

Operands:

\langle ESL \rangle

A valid ESL specification. Include those entries which match the ESL description into the workspace.

entrywsincludelock

Locks the given entries into the workspace.

Syntax:

entrywsincludelock \langle ESL \rangle

Operands:

\langle ESL \rangle

A valid ESL specification. Locks the given entries into the workspace.

entrywsincludeonly

Include only selected project entries in the workspace. Exclude all others.

Syntax:

entrywsincludeonly \langle ESL \rangle

Operands:

\langle ESL \rangle

A valid ESL specification. Include only those entries which match the ESL description in the workspace. Exclude all others.

entrywsincludeunlock

Unlocks the given entries in the Workspace.

Syntax:

entrywsincludeunlock \langle ESL \rangle

Operands:

\langle ESL \rangle

A valid ESL specification. Unlocks the given entries in the Workspace.

epik

This keyword is used to set various options associated with running epik jobs.

Syntax:

epik *analysis_mode*=query | predict | scan *gen_tautomers*=yes | no
input_file= \langle text \rangle *max_output_struct*= \langle n \rangle
metal_binding=yes | no *metal_params*= \langle text \rangle
original_ion_state=yes | no *original_tautomer*=yes | no
ph= \langle x \rangle *ph_tolerance*= \langle x \rangle *pka_file*= \langle text \rangle *solvent*=h2o |
dmso *structure_source*=selected_entries | workspace | file
tautomer_file= \langle text \rangle

Options:

analysis_mode

In predict mode it generates the new structure and calculates the pka values for the ionizable atoms in these new structures. In Query mode it just calculates pka values for certain ionizable atoms in the structure.

Valid values: query
 predict
 scan

Default value: **predict**

gen_tautomers

Generate tautomers

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

input_file

The name of the structure input file.

Valid values: text strings
Default value:

max_output_struct

Maximum output structure to be generated.

Valid values: integers

Default value: **16**

Minimum: 0

metal_binding

Set this to true to turn on generation of metal-binding states.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

metal_params

This is the file name of parameters to use for generating metal binding states.

Valid values: text strings

Default value:

original_ion_state

Include original ionization state

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

original_tautomer

Include original Tautomer.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

ph

Ionization pH

Valid values: reals

Default value: **5**

Minimum: 0.0

Maximum: 14.0

ph_tolerance

Ionization pH tolerance

Valid values: reals

Default value: **2**

Minimum: 0.0

Maximum: 7.0

pka_file

The name of the custom pKa parameter file.

Valid values: text strings

Default value:

solvent

Specify the solvent used for calculating pka values for the structure.

Valid values: h2o
dms0

Default value: **h2o**

structure_source

Whether to use the selected entries in the current project, or what is in the workspace, or a specified file with multiple structures as structure input for the job.

Valid values: selected_entries
 workspace
 file

Default value: **workspace**

tautomer_file

The name of the custom tautomer parameter file.

Valid values: text strings

Default value:

epikreadinput

Read the given Epik options file and set the Epik panel options.

Syntax:

epikreadinput <filename>

Operands:

<filename>

The name of the epik options file to read.

epikstart

Start a epik job with the current settings.

Syntax:

epikstart

epikwriteinput

Write the Epik options to the given file.

Syntax:

epikwriteinput `<filename>`

Operands:

`<filename>`

The name of the epik file to write epik options.

eplayergoto

Go to the specified entry in the ordered sequence of selected project entries and pause there. The table specified in the `eplayersettings` command determines the sequence order.

Syntax:

eplayergoto `<entry_name>`

Operands:

`<entry_name>`

The name of the entry, within the ordered sequence of selected project entries, which is to be included in the workspace. This frame is recorded as an option to the `eplayersettings` command.

eplayergotofirst

Go to the first entry in the ordered sequence of selected project entries. The table specified in the `eplayersettings` command determines the sequence order.

Syntax:

eplayergotofirst

eplayergotolast

Go to the last entry in the ordered sequence of selected project entries. The table specified in the `eplayersettings` command determines the sequence order.

Syntax:

`eplayergotolast`

eplayernextgroup

Go to the first selected entry in the next group that has any selected entries, if there is one, after the group for the current `eplayersettings` frame entry. If the current frame is an ungrouped entry, go to the first selected entry that is in an entry group. The table specified in the `eplayersettings` command determines the sequence order.

Syntax:

`eplayernextgroup`

eplayerpreviousgroup

Go to the first selected entry in the previous group that has any selected entries, if there is one, before the group for the current `eplayersettings` frame entry. This will never go to the ungrouped entries at the top of the table. The table specified in the `eplayersettings` command determines the sequence order.

Syntax:

`eplayerpreviousgroup`

eplayersettings

Set eplayer state variables.

Syntax:

```

eplayersettings frame=<text> frameduration=<x>
intermediates=<n> playmode=loop | reverse | once
playsync=yes | no referentry=<text> script=noaction |
current | file scriptfile=<text> superimpose=none | previous |
reference table=<text> title=<text>
useintermediates=yes | no

```

Options:

frame This option sets the name of the entry being displayed.

Valid values: text strings

Default value:

frameduration

This option determines the minimum duration, in seconds, of each displayed (entry) frame during continuous play. The actual frame duration may be longer than the specified value, due to time required for drawing and screen update.

Valid values: reals

Default value: **0**

Minimum: 0.0

Maximum: 5.0

intermediates

Set the number of interpolated intermediates to generate between each real entry. This does not create entries. Recording movies has a similar option to generate smoother movies.

Valid values: integers

Default value: **20**

playmode This option sets the mode for continuous play. Valid values are “loop”, “reverse”, or “once”. These cause play to continue, change direction, or stop, respectively, when reaching either end of the current entry selection.

Valid values: loop
reverse
once

Default value: **once**

playsync If true, save Workspace changes during continuous play.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

referentry This option sets the name of the reference entry used when the superimpose option is set to “reference”. The reference entry must exist in the current project in order to use this option.

Valid values: text strings

	Default value:
<i>script</i>	This option sets the action for each step in the ePlayer. Valid values are “noaction”, “current”, or “file”. At each step, ePlayer may take no action, run current command script, or run command script from a script file. Valid values: noaction current file Default value: noaction
<i>scriptfile</i>	This option sets the name of the script file used when the script option is set to “Execute Command Script From File”. Valid values: text strings Default value:
<i>superimpose</i>	This option sets the type of superposition done for entries displayed by the eplayer. Valid values are “none”, “previous”, or “reference”. These cause the incoming entry to be displayed with its current coordinates, with its atoms superimposed upon those of the outgoing entry, or with its atoms superimposed upon those of the reference structure for registration. Valid values: none previous reference Default value: none
<i>table</i>	This option identifies the project table which is used to define to sequence order of the selected entries in the project. Valid values: text strings Default value:
<i>title</i>	This option sets the title of the entry being displayed. Valid values: text strings Default value:
<i>useintermediates</i>	If true, temporarily generate and visualize ‘intermediate’ number of additional frames between each real entry when playing forward. Between any two real entries are n intermediates. The first half are duplicates of the first entry, the second half are duplicates of the second entry. All have their coordinates interpolated between the 2 entries. Note that this only works if the entries have the same number of atoms though they should really be conformers. If false, then we don’t generate any intermediates and simply eplay forward the actual entries. Note

that eplaying backwards does not display intermediates even if this option is on. These are temporary intermediates that are for display only in the Workspace.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

eplayerstepahead

Go to the next entry in the ordered sequence of selected project entries, if there is one after the frame specified in the eplayersettings command. The table specified in the eplayersettings command determines the sequence order.

Syntax:

eplayerstepahead

eplayerstepback

Go to the previous entry in the ordered sequence of selected project entries, if there is one before the frame specified in the eplayersettings command. The table specified in the eplayersettings command determines the sequence order.

Syntax:

eplayerstepback

errorcheck

Specifies optional error checking to be performed.

Syntax:

errorcheck *ct*=yes | no

Options:

ct An option which determines whether mmct error checking is enabled during certain (primarily graphical) operations.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

excludedvolumesmarkersettings

Set graphical data of Phase excluded volume markers.

Syntax:

```
excludedvolumesmarkersettings ambient=<x> diffuse=<x>
    emission=<x> radius=<x> shininess=<x> slices=<n>
    specular=<x> stacks=<n> transparency=<x>
```

Options:

<i>ambient</i>	Set material property - ambient, to its red, green, and blue components, for front face. Valid values: reals Default value: 0.4 Minimum: 0.0 Maximum: 1.0
<i>diffuse</i>	Set material property - diffuse, to its red, green, and blue components, for front face. Valid values: reals Default value: 0.4 Minimum: 0.0 Maximum: 1.0
<i>emission</i>	Set material property - emission, to its red, green, and blue components, for front face. Valid values: reals Default value: 0.1 Minimum: 0.0 Maximum: 1.0
<i>radius</i>	The radius of QSAR markers. Valid values: reals Default value: 1 Minimum: 0.1
<i>shininess</i>	Set material property - shininess, for front face. Valid values: reals Default value: 80

	Minimum:	0.0
	Maximum:	128.0
<i>slices</i>	Set the slices of drawing sphere.	
	Valid values:	integers
	Default value:	18
	Minimum:	2
<i>specular</i>	Set material property - specular, to its red, green, and blue components, for front face.	
	Valid values:	reals
	Default value:	0.1
	Minimum:	0.0
	Maximum:	1.0
<i>stacks</i>	Set the stacks of drawing sphere.	
	Valid values:	integers
	Default value:	9
	Minimum:	2
<i>transparency</i>	The transparency of QSAR markers.	
	Valid values:	reals
	Default value:	20
	Minimum:	0.0
	Maximum:	100.0

exportrotatablegroups

Exports the receptor plus optional poses, setting the position of the receptor's movable atoms from the current pose.

Syntax:

```
exportrotatablegroups poses=none | current | matching
    <file name>
```

Options:

<i>poses</i>	Determine what to export.	
	Valid values:	none current matching
	Default value:	none

Operands:

⟨file name⟩

The file name to export the receptor and poses to.

extendselecteditatom

Extends the selection to the given row in the edit atom-property table.

Syntax:

extendselecteditatom ⟨row number⟩

Operands:

⟨row number⟩

The operand is the row number to use to extend the selection.

fieldqsarmarkersettings

Set graphical data of Phase field-based QSAR markers.

Syntax:

fieldqsarmarkersettings *ambient*=⟨x⟩ *diffuse*=⟨x⟩
emission=⟨x⟩ *field*=⟨text⟩ *negative_color*=⟨text⟩
negativecoefficient=⟨x⟩ *negativesaturation*=⟨x⟩
numberpls=⟨n⟩ *positive_color*=⟨text⟩ *positivecoefficient*=⟨x⟩
positivesaturation=⟨x⟩ *radius*=⟨x⟩ *shininess*=⟨x⟩ *slice*=⟨n⟩
specular=⟨x⟩ *stack*=⟨n⟩ *transparency*=⟨x⟩

Options:

<i>ambient</i>	Set material property - ambient, to its red, green, and blue components, for front face. Valid values: reals Default value: 0.4 Minimum: 0.0 Maximum: 1.0
<i>diffuse</i>	Set material property - diffuse, to its red, green, and blue components, for front face. Valid values: reals

	Default value: 0.4
	Minimum: 0.0
	Maximum: 1.0
<i>emission</i>	Set material property - emission, to its red, green, and blue components, for front face.
	Valid values: reals
	Default value: 0.1
	Minimum: 0.0
	Maximum: 1.0
<i>field</i>	The field to set coefficients and saturation for.
	Valid values: text strings
	Default value: ff_e
<i>negative_color</i>	Negative coefficient color
	Valid values: text strings
	Default value: red
<i>negativecoefficient</i>	Set the QSAR visualization option of negative coefficient threshold.
	Valid values: reals
	Default value: -0.044
	Maximum: 0.0
<i>negativesaturation</i>	Set the QSAR visualization option of negative saturation threshold.
	Valid values: reals
	Default value: -0.044
	Maximum: 0.0
<i>numberpls</i>	Set the QSAR visualization option of number of PLS factors.
	Valid values: integers
	Default value: 1
	Minimum: 1
<i>positive_color</i>	Positive coefficient color
	Valid values: text strings
	Default value: blue
<i>positivecoefficient</i>	Set the QSAR visualization option of positive coefficient threshold.

Chapter 5: Commands

	Valid values: reals Default value: 0.044 Minimum: 0.0
<i>positivesaturation</i>	Set the QSAR visualization option of positive saturation threshold. Valid values: reals Default value: 0.044 Minimum: 0.0
<i>radius</i>	Radius in angstroms of the spheres used for visualizing intensities. Valid values: reals Default value: 0.1 Minimum: 0.1
<i>shininess</i>	Set material property - shininess, for front face. Valid values: reals Default value: 80 Minimum: 0.0 Maximum: 128.0
<i>slice</i>	Number of slices to use for visualizing intensities. Valid values: integers Default value: 20 Minimum: 1
<i>specular</i>	Set material property - specular, to its red, green, and blue components, for front face. Valid values: reals Default value: 0.1 Minimum: 0.0 Maximum: 1.0
<i>stack</i>	Number of stacks to use for visualizing intensities. Valid values: integers Default value: 20 Minimum: 1
<i>transparency</i>	The transparency of QSAR markers. Valid values: reals Default value: 50 Minimum: 0.0 Maximum: 100.0

fileread

Read a structure file. If a filename is given then structures will be read from that file. If an explicit “start=” option is not included then starting structure is incremented at each read.

For example:

```
fileread start=1 fileread mmodtmp.dat fileread start=1 test.dat
```

Syntax:

```
fileread delete=yes | no format=mmod | pdb | mol2 | maestro |
sd | jaguarinput | gaussianinput | biograp | xyz | jaguaroutput |
gaussian92 | gaussian94 | gamess | mopaccartesian |
mopacinternal | mopacoutput | babelpdb | mdl | babelmol |
babelmol2 | spartan | spartansemi | spartanmm | gamessinput |
gausscartesian | gaussianz | jaguarzmatrix | jaguarcartesian |
any | reagentprep | mcproinput | cms | maegz | sdcompressed |
glideposeviewer | smiles | smilesgz | common ginfo=yes | no
start=⟨n⟩ tile=yes | no total=⟨n⟩ [⟨file_name⟩]
```

Options:

delete This determines whether the currently displayed structure will be deleted before the new structure is read in.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

format This option sets the format of the file to be read. Valid values are “maestro”, “mmod”, “pdb”, “mol2”, or “sd”.

Valid values:	mmod pdb mol2 maestro sd jaguarinput gaussinput biograf xyz jaguaroutput gaussian92 gaussian94 gamess mopaccartesian mopacinternal mopacoutput babelpdb mdl babelmol babelmol2 spartan spartansemi spartanmm gamessinput gausscartesian gaussianz jaguarzmatrix jaguarcartesian any reagentprep mcproinput cms maegz sdcompressed glideposeviewer smiles smilesgz common
Default value:	maestro

ginfo

This determines whether the file will be read with just the geometric information or if the visual information in the file will be read and used for displaying the structures.

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	false

<i>start</i>	This option sets the number of the first structure to be read
	Valid values: integers
	Default value: 1
	Minimum: 1
<i>tile</i>	This determines whether the structures displayed on screen will be “tiled”.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>total</i>	The total number of structures to be read from the file
	Valid values: integers
	Default value: 1
	Minimum: 1

Operands:

[<file_name>]

The name of the file from which the structure will be read. If no name is specified, then no file read will be done.

Aliases:

read (see [\[read\]](#), page 666)

filewrite

If the write command is issued with options, but without a filename then nothing is written but the options are updated. If the write command alone is issued then writing is performed to the currently open file with the current options. A write command with both options and a file name specified will result in the write being performed with the new options to the specified file.

Syntax:

```
filewrite append=yes | no displayed_atoms=yes | no
        format=mmod | pdb | mol2 | maestro | sd | jaguarinput |
        gaussianinput | biograf | xyz | jaguaroutput | gaussian92 |
        gaussian94 | gamess | mopaccartesian | mopacinternal |
        mopacoutput | babelpdb | mdl | babelmol | babelmol2 |
        spartan | spartansemi | spartanmm | gamessinput |
        gausscartesian | gaussianz | jaguarzmatrix | jaguarcartesian |
        any | reagentprep | mcproinput | cms | maegz | sdcompressed
        | glideposeviewer | smiles | smilesgz | common
graphical=yes | no separate=yes | no title=<text>
<file_name>
```

Options:

append This option determines whether to append to the file which is going to be written

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

displayed_atoms

This option determines whether the displayed atoms will be saved. This option is not currently supported.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

format This option sets the format of the file to be read. Valid values are “mmod”, “pdb” or “mol2”.

Valid values: mmod
 pdb
 mol2
 maestro
 sd
 jaguarinput
 gaussinput
 biograf
 xyz
 jaguaroutput
 gaussian92
 gaussian94
 gamess
 mopaccartesian
 mopacinternal
 mopacoutput
 babelpdb
 mdl
 babelmol
 babelmol2
 spartan
 spartansemi
 spartanmm
 gamessinput
 gausscartesian
 gaussianz
 jaguarzmatrix
 jaguarcartesian
 any
 reagentprep
 mcproinput
 cms
 maegz
 sdcompressed
 glideposeviewer
 smiles
 smilesgz
 common

Default value: **maestro**

graphical This option determines whether to write out graphical information when a structure is being saved to disk.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

separate This option determines whether to write out the CT as one structure or separate structures. This option is not currently supported.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

title The title for the output structure.

Valid values: text strings
Default value:

Operands:

⟨ file_name ⟩

The name of the file to which the structure will be written.

Aliases:

write (see [\[write\]](#), page 824)

filter

Creates a new entry filter. The filter name must be a single token (or “quoted” if multiple tokens). A filter can be redefined by specifying a new definition.

Syntax:

filter ⟨ filter_name ⟩ ⟨ ESL-definition ⟩

Operands:

⟨ filter_name ⟩ ⟨ ESL-definition ⟩

The name which will be applied to the filter. If the name contains embedded spaces then it must be enclosed in double quotation marks.

find

Find atom(s), residue(s), chain, or molecule described by the ASL operand.

Syntax:

```

find atombyname=⟨text⟩ atombynum=⟨n⟩
      atomnumberstring=⟨text⟩ byasl=⟨text⟩ center=yes | no
      chainname=⟨text⟩ elementstring=⟨text⟩
      findmethod=bynumber | byname findtype=atom | residue |
      chain | molecule | asl fit=yes | no inscode=⟨text⟩
      label=yes | no markall=yes | no matchesnum=⟨text⟩
      mode=resnum | atomnum | smarts | restype | secstruct |
      element | pdbatom | userdef molnum=⟨n⟩
      pdbatomnamestring=⟨text⟩ residuenumberstring=⟨text⟩
      residuetypeststring=⟨text⟩ resnum=⟨n⟩
      secondarystructure=helix | loop | strand
      showmarkers=yes | no smartsstring=⟨text⟩
      userdefinedstring=⟨text⟩ zoom=yes | no ⟨ASL⟩

```

Options:

atombyname

This option determines name of the atom to be found.

Valid values: text strings

Default value:

atombynum

This option determines value of the atom number to be found

Valid values: integers

Default value: **1**

Minimum: 1

atomnumberstring

The string of Atom number mode.

Valid values: text strings

Default value:

byasl

This option determines the ASL specification fo the objects to be found (atom, residue, chain ot molecule).

Valid values: text strings

Default value:

center

This option determines whether the structure will center on the found atom or on the centroid of the atoms [of residue, chain or molecule].

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

chainname

This option determines the chain name to be found, or the chain name in which to find a specified residue.

Valid values: text strings

Chapter 5: Commands

	Default value:
<i>elementstring</i>	The string of Element mode. Valid values: text strings Default value:
<i>findmethod</i>	This option determines the method used to specify an atom to be found, including By Number and By Name . Valid values: bynumber byname Default value: bynumber
<i>findtype</i>	This option determines the find type, including five options of atom, residue, chain, molecule and asl. Valid values: atom residue chain molecule asl Default value: atom
<i>fit</i>	This option determines whether to zoom and center the found item. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>inscode</i>	This option determines value of the insertion code of the residue to be found. Valid values: text strings Default value:
<i>label</i>	This option determines whether the found atoms will be labeled. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>markall</i>	This option determines whether all the found atom(s) or atom(s) of residue, chain or molecule are marked (TRUE), or only the principal atom(s) are marked (FALSE). Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>matchesnum</i>	This option shows the number of matched atoms over the total number of atoms in the structure. Valid values: text strings Default value: /

<i>mode</i>	<p>This option determines the find mode.</p> <p>Valid values: resnum atomnum smarts restype secstruct element pdbatom userdef</p> <p>Default value: resnum</p>
<i>molnum</i>	<p>This option determines the number of the molecule to be found</p> <p>Valid values: integers</p> <p>Default value: 1</p> <p>Minimum: 1</p>
<i>pdbatomnamestring</i>	<p>The string of PDB atom name mode.</p> <p>Valid values: text strings</p> <p>Default value:</p>
<i>residuenumberstring</i>	<p>The string of Residue number mode.</p> <p>Valid values: text strings</p> <p>Default value:</p>
<i>residuetypestring</i>	<p>The string of residue type.</p> <p>Valid values: text strings</p> <p>Default value:</p>
<i>resnum</i>	<p>This option determines value of the residue number to be found.</p> <p>Valid values: integers</p> <p>Default value: 1</p>
<i>secondarystructure</i>	<p>This option determines the secondary structure type.</p> <p>Valid values: helix loop strand</p> <p>Default value: helix</p>
<i>showmarkers</i>	<p>This option determines whether the markers for the principal atoms should be shown (TRUE) or not (FALSE).</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: true</p>

Chapter 5: Commands

smartsstring

The string of SMARTS mode.

Valid values: text strings

Default value:

userdefinedstring

The string of User defined mode.

Valid values: text strings

Default value:

zoom

This option determines whether or not the structure will be zoomed to fit the found atoms.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

Operands:

⟨ ASL ⟩

A string in the atom specification language which describes the set of atoms which are to be found.

findnext

Select the next item of found items.

Syntax:

findnext

findprevious

Select the previous item of found items.

Syntax:

findprevious

fit

Fit structures to the Workspace.

Syntax:

fit \langle ASL \rangle

Operands:

\langle ASL \rangle

If present, we do fit on the ASL, if not present we do fit on all atoms.

fitbox

Fit the box to the Workspace.

Syntax:

fitbox \langle xmin xmax ymin ymax zmin zmax \rangle

Operands:

\langle xmin xmax ymin ymax zmin zmax \rangle

Fit the box to the Workspace.

fitligand

Fit a ligand to the Workspace.

Syntax:

fitligand

fitnextligand

Fit next ligand to the Workspace.

Syntax:

fitnextligand

fitoptions

Options for the fit command

Syntax:

```
fitoptions move_clip_planes=yes | no zoom_topview=yes | no
```

Options:

move_clip_planes

If set to true, then the clipping planes will be moved when we do a fit to Workspace. If it is set to false, we won't.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

zoom_topview

If set to true, then the clipping planes view will be zoomed when we fit to Workspace.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

fitpreviousligand

Fit previous ligand to the Workspace.

Syntax:

```
fitpreviousligand
```

flexibleligandalignment

Align ligands in a flexible manner against a specified template. This always operates on selected entries in the project but allows the template to be specified either as the first selected entry, the included entry or a specific entry.

Syntax:


```
flexibleligandalignment template=first_selected | included |
entry [ templat_entry_id ]
```

Options:

template This option determines how the template for alignment is specified. It's either the first selected entry, the included entry or a specified entry.

Valid values: first_selected
 included
 entry

Default value: **first_selected**

Operands:

```
[ templat_entry_id ]
```

If the template is being specified by a *entry_id*, the operand is used as the id of the template entry.

forcefield

Used to display the force field in the force field viewer.

Syntax:

```
forcefield
```

forcefieldbend

A keyword which controls the display of the bend interactions in the FF viewer.

Syntax:

```
forcefieldbend select=⟨n⟩ show=all | high | medium | low
                 sortbyenergy=yes | no
```

Options:

select Set the interaction which is currently selected

Valid values: integers

Default value: **0**

Minimum: 0

show Sets a filter which determines which interactions are shown depending on the quality of their parameters.

Valid values: all
high
medium
low

Default value: **all**

sortbyenergy

Determines if the interactions are sorted by energy

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

forcefielddele

A keyword which controls the display of the electrostatic interactions in the FF viewer.

Syntax:

```
forcefielddele select=<n> show=all | high | medium | low  
sortbyenergy=yes | no
```

Options:

select Set the interaction which is currently selected

Valid values: integers

Default value: **0**

Minimum: 0

show Sets a filter which determines which interactions are shown depending on the quality of their parameters.

Valid values: all
high
medium
low

Default value: **all**

sortbyenergy

Determines if the interactions are sorted by energy

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

forcefieldgbsolv

A keyword which controls the display of the GBSOLV interactions in the FF viewer.

Syntax:

```
forcefieldgbsolv select=<n> show=all | high | medium | low
                  sortbyenergy=yes | no
```

Options:

<i>select</i>	Set the interaction which is currently selected
	Valid values: integers
	Default value: 0
	Minimum: 0
<i>show</i>	Sets a filter which determines which interactions are shown depending on the quality of their parameters.
	Valid values: all high medium low
	Default value: all
<i>sortbyenergy</i>	Determines if the interactions are sorted by energy
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false

forcefieldimproper

A keyword which controls the display of the improper torsion interactions in the FF viewer.

Syntax:

```
forcefieldimproper select=<n> show=all | high | medium | low
                  sortbyenergy=yes | no
```

Options:

<i>select</i>	Set the interaction which is currently selected
	Valid values: integers
	Default value: 0

	Minimum:	0
<i>show</i>	Sets a filter which determines which interactions are shown depending on the quality of their parameters.	
	Valid values:	all high medium low
	Default value:	all
<i>sortbyenergy</i>	Determines if the interactions are sorted by energy	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	false

forcefieldsasolv

A keyword which controls the display of the SASOLV interactions in the FF viewer.

Syntax:

```
forcefieldsasolv select=<n> show=all | high | medium | low
                  sortbyenergy=yes | no
```

Options:

<i>select</i>	Set the interaction which is currently selected	
	Valid values:	integers
	Default value:	0
	Minimum:	0
<i>show</i>	Sets a filter which determines which interactions are shown depending on the quality of their parameters.	
	Valid values:	all high medium low
	Default value:	all
<i>sortbyenergy</i>	Determines if the interactions are sorted by energy	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	false

forcefieldstretch

A keyword which controls the display of the stretch interactions in the FF viewer.

Syntax:

```
forcefieldstretch select=⟨n⟩ show=all | high | medium | low
                  sortbyenergy=yes | no
```

Options:

<i>select</i>	Set the interaction which is currently selected Valid values: integers Default value: 0 Minimum: 0
<i>show</i>	Sets a filter which determines which interactions are shown depending on the quality of their parameters. Valid values: all high medium low Default value: all
<i>sortbyenergy</i>	Determines if the interactions are sorted by energy Valid values: boolean (true false; yes no; y n; on off) Default value: false

forcefieldtorsion

A keyword which controls the display of the torsion interactions in the FF viewer.

Syntax:

```
forcefieldtorsion select=⟨n⟩ show=all | high | medium | low
                  sortbyenergy=yes | no
```

Options:

<i>select</i>	Set the interaction which is currently selected Valid values: integers Default value: 0
---------------	--

	Minimum:	0
<i>show</i>	Sets a filter which determines which interactions are shown depending on the quality of their parameters.	
	Valid values:	all high medium low
	Default value:	all
<i>sortbyenergy</i>	Determines if the interactions are sorted by energy	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	false

forcefieldvdw

A keyword which controls the display of the VDW interactions in the FF viewer.

Syntax:

```
forcefieldvdw select=<n> show=all | high | medium | low
               sortbyenergy=yes | no
```

Options:

<i>select</i>	Set the interaction which is currently selected	
	Valid values:	integers
	Default value:	0
	Minimum:	0
<i>show</i>	Sets a filter which determines which interactions are shown depending on the quality of their parameters.	
	Valid values:	all high medium low
	Default value:	all
<i>sortbyenergy</i>	Determines if the interactions are sorted by energy	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	false

forcefieldview

Used to specify a .mmo file which contains the interactions to be displayed in the FF viewer.

Syntax:

```
forcefieldview <mmo_file_name>
```

Operands:

```
<mmo_file_name>
```

The name of the .mmo file which is to be opened and displayed. The full filename, including the .mmo suffix must be specified.

forcefieldwilson

A keyword which controls the display of the wilson interactions in the FF viewer.

Syntax:

```
forcefieldwilson select=<n> show=all | high | medium | low  
                  sortbyenergy=yes | no
```

Options:

select Set the interaction which is currently selected

Valid values: integers

Default value: **0**

Minimum: 0

show Sets a filter which determines which interactions are shown depending on the quality of their parameters.

Valid values: all
 high
 medium
 low

Default value: **all**

sortbyenergy

Determines if the interactions are sorted by energy

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

formalcharge

Increment, decrement or set the formal charge for atoms which match the ASL specification.

Syntax:

formalcharge increment|decrement|⟨formal_charge⟩⟨ASL⟩

Operands:

increment|decrement|⟨formal_charge⟩⟨ASL⟩

The first operand is either the word “increment”, the word “decrement”, or an integer value represening the formal charge to be used for all atoms which match the specification. “increment” and “decrement” change the atom’s current formal charge by +1 or -1 respectively. The second operand is a valid ASL string which defines the set of atoms which are to have their atom names changed. An error will be issued for the first atom in the set which cannot have its formal charge set (usually because no suitable charged type exists) and no further atoms will have their charge changed.

fragment

This command sets the fragment mode and the current fragment used for growing, placing and replacing.

Syntax:

fragment ⟨fragment_mode⟩ [⟨current_fragment⟩]

Operands:

⟨fragment_mode⟩ [⟨current_fragment⟩]

The first operand is the name of the fragment mode, this must be a valid name from the file \$SCHRODINGER/mmshare-vX.X/data/res/mmfrag.ini. The optional second argument is the name of a fragment within that mode which is to be made the current fragment.

fragmentfuse

Fuse the currently selected fragment on to the atoms given in the operand. If the atoms are not bound then it’s an error. It is also an error if the currently

selected fragment is not one for which fusion is defined (mostly just small rings in Organic mode).

Syntax:

fragmentfuse at1 at2

Operands:

at1 at2

The two atoms which define the bond where fusion of the currently selected fragment should take place.

fragmenttype

This command sets the fragment type for the current fragment mode.

Syntax:

fragmenttype <fragment_type>

Operands:

<fragment_type>

The first operand is the a qualifier for the fragment which some modes (such as furanose and pyranose) have. For example alpha or beta for pyranose or furanose sugars. The fragment_type must be a valid type for the current fragment mode.

frozenatom

Specifies a single atom to be frozen at its current position in a MacroModel calculation. Other atoms in the molecule will continue to interact with this atom. The frozen atom however will not “feel” any other atoms in the molecule.

Syntax:

frozenatom <atom_number>

Operands:

<atom_number>

The number of an atom which is to be treated as frozen during a MacroModel calculation.

frozenset

Specifies a set of atoms to be frozen at their current positions in a MacroModel calculation. Other atoms in the molecule will continue to interact with these atoms. The frozen atoms however will not “feel” any other atoms in the molecule.

Syntax:

frozenset <ASL>

Operands:

<ASL>

A string in the atom specification language which describes the set of atoms which are to be treated as frozen in a MacroModel calculation.

fuse

Fuse atom list pairs

Syntax:

fuse

geometrycleanup

Cleans up the geometry

Syntax:

geometrycleanup

glideactivesiteres

Defines a residue to be included in the active site for a Glide grid calculation.

Syntax:

```
glideactivesiteres <chain>:<molnum>:<resnum>:<insertioncode>
```

Operands:

```
<chain>:<molnum>:<resnum>:<insertioncode>
```

A colon separated list giving the PDB chain name, the residue number and insertion code for a residue to be included in the list of active site residues.

glideactivesiteset

Defines a set from atoms, from which all residues will be added as part of the active site for Glide calculations.

Syntax:

```
glideactivesiteset <ASL>
```

Operands:

```
<ASL>
```

A valid ASL specification. All residues which have any atoms in this set will be added to the list of residues that define the active site.

glidecalcboxfromligand

Uses the information specified in the `glideligandgrid` command to calculate the grid box for a Glide calculation.

Syntax:

```
glidecalcboxfromligand
```

glideconstraintatomlabel

Specifies a constraint atom label in the receptor for a Glide calculation.

Syntax:

glideconstraintatomlabel \langle atom-number \rangle \langle label \rangle

Operands:

\langle atom-number \rangle \langle label \rangle

The name which will be applied to the atom.

glideconstraintposition

Specifies a constraint position in the receptor for a Glide calculation.

Syntax:

glideconstraintposition *radius*= \langle x \rangle *radiusnoe*= \langle x \rangle
usenoe=yes | no *x*= \langle x \rangle *y*= \langle x \rangle *z*= \langle x \rangle

Options:

<i>radius</i>	The radius of a Glide constraint position. Valid values: reals Default value: 1 Minimum: 0.0001
<i>radiusnoe</i>	The minimum NOE distance of a Glide constraint position. Valid values: reals Default value: 0 Minimum: 0.0
<i>usenoe</i>	Use the NOE constraint or not. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>x</i>	The X coordinate of a Glide constraint position. Valid values: reals Default value: 0
<i>y</i>	The Y coordinate of a Glide constraint position. Valid values: reals Default value: 0
<i>z</i>	The Z coordinate of a Glide constraint position. Valid values: reals Default value: 0

glideconstraintregion

Specifies a constraint region in the receptor for a Glide calculation.

Syntax:

```
glideconstraintregion <region_name> <cell>
```

Operands:

<region_name> <cell>

The name which will be applied to the region. If the name contains embedded spaces then it must be enclosed in double quotation marks.

glideconstraintregionactivecell

Set an active cell for the region.

Syntax:

```
glideconstraintregionactivecell <region_name> <active_cell>
```

Operands:

<region_name> <active_cell>

The name which will be applied to the region. If the name contains embedded spaces then it must be enclosed in double quotation marks.

glideconstraintregiongrow

Grow cells for a region.

Syntax:

```
glideconstraintregiongrow <region_name>
```

Operands:

<region_name>

The name which will be applied to the region. If the name contains embedded spaces then it must be enclosed in double quotation marks.

glideconstraintregioninvisible

Set a constraint region as invisible in workspace.

Syntax:

```
glideconstraintregioninvisible <region_name>
```

Operands:

<region_name>

The name which will be applied to the region. If the name contains embedded spaces then it must be enclosed in double quotation marks.

glideconstraintregionshrink

Shrink cells for a region.

Syntax:

```
glideconstraintregionshrink <region_name>
```

Operands:

<region_name>

The name which will be applied to the region. If the name contains embedded spaces then it must be enclosed in double quotation marks.

glideconstraintregionvisible

Set a constraint region as visible in workspace.

Syntax:

```
glideconstraintregionvisible <region_name>
```

Operands:

<region_name>

The name which will be applied to the region. If the name contains embedded spaces then it must be enclosed in double quotation marks.

glidecoremolecule

Defines a core molecule to be treated as the reference ligand for a Glide docking.

Syntax:

```
glidecoremolecule <molecule_number>
```

Operands:

<molecule_number>

The number of a molecule to be included as the ligand.

glidecorermssubset

Sets the subset for RMSD calculation in Glide docking.

Syntax:

```
glidecorermssubset <atom number>
```

Operands:

<atom number>

The atom number of the RMSD subset.

glidedisplayreceptor

Displays the receptor in the Workspace.

Syntax:

```
glidedisplayreceptor
```

glidedockconstraintposition

Specifies a constraint position in the receptor for a Glide calculation.

Syntax:

```
glidedockconstraintposition feature=<n> index=<n>
radius=<x> radiusnoe=<x> type=<n> use1=yes | no
use2=yes | no use3=yes | no use4=yes | no
usenoe=yes | no x=<x> y=<x> z=<x>
```

Options:

<i>feature</i>	The constraint feature of position in docking. Valid values: integers Default value: -1
<i>index</i>	The index of position in docking. Valid values: integers Default value: 0
<i>radius</i>	The radius of a Glide constraint position. Valid values: reals Default value: 1 Minimum: 0.0001
<i>radiusnoe</i>	The minimum NOE distance of a Glide constraint position. Valid values: reals Default value: 0 Minimum: 0.0
<i>type</i>	The constraint type of position in docking. Valid values: integers Default value: 0
<i>use1</i>	The flag indicates if this position will be used in docking for group 1. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use2</i>	The flag indicates if this position will be used in docking for group 2. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use3</i>	The flag indicates if this position will be used in docking for group 3. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use4</i>	The flag indicates if this position will be used in docking for group 4. Valid values: boolean (true false; yes no; y n; on off) Default value: false

<i>usenoe</i>	Use the NOE constraint or not. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>x</i>	The X coordinate of a Glide constraint position. Valid values: reals Default value: 0
<i>y</i>	The Y coordinate of a Glide constraint position. Valid values: reals Default value: 0
<i>z</i>	The Z coordinate of a Glide constraint position. Valid values: reals Default value: 0

glidedockconstraintregion

Specifies a constraint region in the receptor for a Glide calculation.

Syntax:

```
glidedockconstraintregion atoms=<n> feature=<n>
                        index=<n> type=<n> use1=yes | no use2=yes | no
                        use3=yes | no use4=yes | no <region_name>
```

Options:

<i>atoms</i>	The number of required ligand atoms of constraint region in docking. Valid values: integers Default value: 1
<i>feature</i>	The constraint feature of constraint region in docking. Valid values: integers Default value: -1
<i>index</i>	The index of constraint region in docking. Valid values: integers Default value: 0
<i>type</i>	The constraint type of constraint region in docking. Valid values: integers Default value: 0

Chapter 5: Commands

<i>use1</i>	The flag indicates if this constraint region will be used in docking for group 1. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use2</i>	The flag indicates if this constraint region will be used in docking for group 2. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use3</i>	The flag indicates if this constraint region will be used in docking for group 3. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use4</i>	The flag indicates if this constraint region will be used in docking for group 4. Valid values: boolean (true false; yes no; y n; on off) Default value: false

Operands:

`<region_name>`

The name which will be applied to the region. If the name contains embedded spaces then it must be enclosed in double quotation marks.

glidedockexcludedvolume

Specifies an excluded volume in the receptor for a Glide docking calculation.

Syntax:

glidedockexcludedvolume *radius*=`<x>` *x*=`<x>` *y*=`<x>` *z*=`<x>`

Options:

<i>radius</i>	The radius of a Glide excluded volume. Valid values: reals Default value: 3 Minimum: 0.0001
<i>x</i>	The X coordinate of a Glide excluded volume. Valid values: reals Default value: 0

y The Y coordinate of a Glide excluded volume.

Valid values: reals

Default value: **0**

z The Z coordinate of a Glide excluded volume.

Valid values: reals

Default value: **0**

glidedockingreset

Resets all the options related to glide ligand docking to there defaults.

Syntax:

glidedockingreset

glidedockregionnumatoms

Specifies the required ligand atoms for this constraint region.

Syntax:

glidedockregionnumatoms \langle region_name \rangle \langle num_atoms \rangle

Operands:

\langle region_name \rangle \langle num_atoms \rangle

The name of the region and the minimum number of atoms permitted for this region. If the name contains embedded spaces then it must be enclosed in double quotation marks.

glidedockregionselect

Selects a hydrophobic constraint region to be used as a constraint in the docking job.

Syntax:

glidedockregionselect *group*=⟨*n*⟩ ⟨*region_name*⟩

Options:

<i>group</i>	The constraint group for region.
Valid values:	integers
Default value:	1
Minimum:	0
Maximum:	3

Operands:

⟨*region_name*⟩

The name of the region to be selected as a constraint in the docking job. If the name contains embedded spaces then it must be enclosed in double quotation marks.

glidedockregionunselect

Unselects a hydrophobic constraint region so it will not be used as a constraint in the docking job.

Syntax:

glidedockregionunselect *group*=⟨*n*⟩ ⟨*region_name*⟩

Options:

<i>group</i>	The constraint group for region.
Valid values:	integers
Default value:	1
Minimum:	0
Maximum:	3

Operands:

⟨*region_name*⟩

The name of the region to be unselected as a constraint in the docking job. If the name contains embedded spaces then it must be enclosed in double quotation marks.

glideexcludedvolume

Specifies an excluded volume in the receptor for a Glide calculation.

Syntax:

```
glideexcludedvolume radius=<x> select=yes | no x=<x> y=<x>
                        z=<x>
```

Options:

<i>radius</i>	The radius of a Glide excluded volume. Valid values: reals Default value: 3 Minimum: 0.0001
<i>select</i>	Selection state of the excluded volume. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>x</i>	The X coordinate of a Glide excluded volume. Valid values: reals Default value: 0
<i>y</i>	The Y coordinate of a Glide excluded volume. Valid values: reals Default value: 0
<i>z</i>	The Z coordinate of a Glide excluded volume. Valid values: reals Default value: 0

glidegridhydrophobic

Used to set all the options associated with the hydrophobic constraints in Glide grid generation.

Syntax:

```
glidegridhydrophobic num_vertices=<n> threshold=<x>
```

Options:

<i>num_vertices</i>	At least this number of cell vertices must meet the threshold in order for the cell to be displayed.
---------------------	--

Valid values:	integers
Default value:	5
Minimum:	1
Maximum:	8

threshold This is cutoff for which grid cells are considered hydrophobic.

Valid values:	reals
Default value:	-0.5
Minimum:	-0.65
Maximum:	-0.01

glidegridreset

Resets all the options related to glide grid generation to their defaults.

Syntax:

glidegridreset

glideligand

Specifies settings about the ligand(s) to be used in a Glide job.

Syntax:

glideligand *amidebondrotations*=yes | no *bondrotation*=⟨text⟩
bondrotationtype=all | amides | none
conformations=usesupplied | generate | inplace | mininplace
definereference=yes | no *dockdisplayed*=yes | no
dockfromfile=yes | no *dockrange*=all | specified *econfcut*=⟨x⟩
endlig=⟨n⟩ *format*=maestro | sd | mol2 | pdb
inputring=yes | no *lig_ccut*=⟨x⟩ *lig_vscale*=⟨x⟩
ligandsfile=⟨text⟩ *ligandsource*=neither | entries | extfile
ninvert=yes | no *reflig*=yes | no *ringconf*=yes | no
startlig=⟨n⟩ *useinputcharges*=yes | no

Options:

amidebondrotations

An option which allows amide bond rotations.

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

bondrotation

The option of amide bond rotation.

Valid values: text strings

Default value: **penal**

bondrotationtype

The type of bond rotation.

Valid values: all
amides
none

Default value: **amides**

conformations

Whether to use only supplied conformations, to generate them with confgen or to treat them “in-place”.

Valid values: usesupplied
generate
inplace
mininplace

Default value: **generate**

definereference

An option which determines if a reference ligand is to be defined.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

dockdisplayed

An option which determines if the currently displayed ligand will be docked.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

dockfromfile

[NOTE: This option is no longer used.] An option which determines ligands from external files are to be docked:

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

dockrange

The range of structures in the file to be docked - whether to dock all from the file or a specified range.

Valid values: all
specified

Default value: **all**

econfcut

The energy window for keeping conformations.

Valid values: reals

Chapter 5: Commands

	Default value: 12 Minimum: 0.0
<i>endlig</i>	The final structure from the file to be docked. Valid values: integers Default value: 1000 Minimum: 0
<i>format</i>	The format of the ligands. Valid values: maestro sd mol2 pdb Default value: maestro
<i>inputring</i>	An option which allows input ring conformation Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>lig_ccut</i>	The partial atomic charge below which ligand atoms are considered to be non-polar and will have their VDW radii scaled. Valid values: reals Default value: 0.15 Minimum: 0.00000001
<i>lig_vscale</i>	The scaling factor for the VDW radii of non-polar ligand atoms. Valid values: reals Default value: 0.8 Minimum: 0.00000001
<i>ligandsfile</i>	The file containing one or more ligands. Valid values: text strings Default value:
<i>ligandsource</i>	The source of the ligands to be docked Valid values: neither entries extfile Default value: extfile
<i>ninvert</i>	An option which allows nitrogen inversions Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>reflig</i>	An option which determines if a reference ligand is to be used in this Glide job. Valid values: boolean (true false; yes no; y n; on off) Default value: false

<i>ringconf</i>	An option which allows ring flips
Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true
<i>startlig</i>	The first structure from the file to be docked.
Valid values:	integers
Default value:	1
Minimum:	1
<i>useinputcharges</i>	An option which determines if Impact uses Maestro charges (true) or force field derived ones (false).
Valid values:	boolean (true false; yes no; y n; on off)
Default value:	false

glideligandfile

Specifies a file containing one or more ligands to be used for a Glide docking run.

Syntax:

glideligandfile `<ligand_file_name>`

Operands:

`<ligand_file_name>`

The name of a file containing one or more ligands.

glideligandgrid

Defines the ligand to be used to calculate the grid for the glide calculation.

Syntax:

glideligandgrid *ligandentry*=`<text>` *ligandformat*=maestro | sd |
mol2 | pdb *ligandindex*=`<n>` *ligandsource*=displayed |
entryname | extfile | none *usedisplayed*=yes | no
`<ligand file name>`

Options:

Chapter 5: Commands

ligandentry

The name of the entry that will be used to define the grid for the Glide calculation.

Valid values: text strings

Default value:

ligandformat

The format of the file that will be used to calculate the grid box.

Valid values: maestro
sd
mol2
pdb

Default value: **maestro**

ligandindex

The index of the structure that is to be read as a ligand.

Valid values: integers

Default value: **1**

Minimum: 1

ligandsource

The source of the ligand to be used to define the Grid for the Glide calculation.

Valid values: displayed
entryname
extfile
none

Default value: **displayed**

usedisplayed

[NOTE: This option is no longer used.] An option which determines if the currently displayed ligand will be used to calculate the grid for the Glide calculation

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

⟨ ligand file name ⟩

The name of the file containing the ligand to be used to define the grid.

glidelocatehydrophobiccells

Runs hppmap and displays hydrophobic cells.

Syntax:

`glidelocatehydrophobiccells`

glideoutput

Controls a number of options for how the poses are output during a Glide calculation.

Syntax:

```
glideoutput applystrain=yes | no delpose=<x>
           interactionradius=<x> maxperlig=<n> nreport=<n>
           numposes=<n> outputdisp=poseviewer | ligandsonly
           outputformat=maestro | sd posedist=<x> postdock=yes | no
           residueinteraction=range | picking rmsdtoinputgeom=yes | no
           strainscale=<x> strainthreshold=<x> tolerance=<x>
           writeinteraction=yes | no writereport=yes | no
```

Options:

applystrain

An option which determines if the Glide job will apply strain correction terms.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

delpose

The maximum displacement for discarding duplicate poses.

Valid values: reals

Default value: **1.3**

Minimum: 0.0

interactionradius

The radius of grid center for per-residue interaction scores.

Valid values: reals

Default value: **12**

Minimum: 0.0

maxperlig

The maximum number of poses per ligand to be written.

Valid values: integers

Default value: **1**

Minimum: 0

nreport

The number of coordinate sets to write to disk.

Valid values: integers

Chapter 5: Commands

	Default value: 10000 Minimum: 1
<i>numposes</i>	The number of poses per ligand to be included. Valid values: integers Default value: 5 Minimum: 0
<i>outputdisp</i>	Whether the output structure file will be intended to be viewed with the Maestro pose viewer (and include the receptor structure) or is just to be the docked ligands. Valid values: poseviewer ligandsonly Default value: poseviewer
<i>outputformat</i>	Whether the output structure file format will be the Maestro format or SD format. Valid values: maestro sd Default value: maestro
<i>posedist</i>	The distance criterion for distinct poses. Valid values: reals Default value: 0.5 Minimum: 0.0
<i>postdock</i>	An option which determines if the Glide job will do post-docking minimization. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>residueinteraction</i>	Whether writing er-residue interaction scores for 1) residues within a distance of grid center or 2) picking residues to include. Valid values: range picking Default value: range
<i>rmsdtoinputgeom</i>	An option which determines if the Glide job will calculate the RMSD against the user supplied input geometry. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>strainscale</i>	The amount to scale ligand strain by. Valid values: reals

Default value: **0.25**
 Minimum: 0.0
 Maximum: 1.0

strainthreshold

Threshold for correcting ligand strain.

Valid values: reals
 Default value: **4**
 Minimum: 0.0

tolerance

For XP docking, a tolerance value that controls the window size where a post-docking minimized structure is kept, even though its GlideScore is less favorable than the original pose.

Valid values: reals
 Default value: **0.5**
 Minimum: 0.0
 Maximum: 10000.0

writeinteraction

An option which determines if the Glide job will write per-residue interaction scores for residues.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

writereport

Turn this on to generate a .rept or .scor file

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

glidereceptorligand

Defines a on-screen molecule or entry to be treated as the ligand for a Glide calculation.

Syntax:

```
glidereceptorligand ligand=molecule | entry  

  { molecule_number | entry_id }
```

Options:

ligand Determine the ligand is 1) a molecule (operand is molecule number), or 2) an entry (operand is entry ID).
 Valid values: molecule
 entry

Default value: **molecule**

Operands:

`< molecule_number | entry_id >`

The molecule number or entry ID to be included as the ligand.

glidereferenceligand

Specifies a file containing the reference ligand for a glide docking run.

Syntax:

```
glidereferenceligand index=<n> ligandentry=<text>  
    ligandfile=<text> ligandformat=maestro | sd | mol2 | pdb  
    ligandsource=displayed | entryname | extfile | none  
    use_displayed=yes | no
```

Options:

index The index of the reference ligand in the reference ligand file.
Valid values: integers
Default value: **1**
Minimum: 1

ligandentry The name of the entry that will be used as the reference ligand.
Valid values: text strings
Default value:

ligandfile The name of the file containing the ligand structure that will be
used as the reference ligand.
Valid values: text strings
Default value:

ligandformat The format of the file that will be used to calculate the grid box.
Valid values: maestro
 sd
 mol2
 pdb
Default value: **maestro**

ligandsource The source of the ligand to be used as the reference ligand.

Valid values: displayed
 entryname
 extfile
 none
 Default value: **none**

use_displayed

[NOTE: This option is no longer used.] An option which determines if the reference ligand is that which is already displayed.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

glideringsampling

Controls a number of options for how rings are sampled.

Syntax:

glideringsampling *energywindow*=⟨x⟩

Options:

energywindow

Set the cutoff, in kcal/mol, for ring-conformation energies in Glide.

Valid values: reals
 Default value: **2.5**
 Minimum: 0.0

gliderotatablegroups

Settings for rotatable groups in a grid generation job.

Syntax:

gliderotatablegroups *allowrotation*=yes | no

Options:

allowrotation

An option which controls whether or not we allow rotation of receptor hydroxyl groups.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

gliderotatablegroupsextenduserow

Extend the used groups to include the given row

Syntax:

gliderotatablegroupsextenduserow \langle row \rangle

Operands:

\langle row \rangle

The row to extend to in the rotatable groups table.

gliderotatablegroupstoggleatom

Toggles the row containing the given atom.

Syntax:

gliderotatablegroupstoggleatom \langle row \rangle

Operands:

\langle row \rangle

The atom to toggle in the rotatable groups table.

gliderotatablegroupstogglerow

Toggles the given row

Syntax:

gliderotatablegroupstogglerow \langle row \rangle

Operands:

\langle row \rangle

The row to toggle in the rotatable groups table.

gliderotatablegroupsuseexisting

Select only the groups which have existing definitions in the Workspace.

Syntax:

gliderotatablegroupsuseexisting

glidescoring

Controls a number of options for how the poses are scored.

Syntax:

glidescoring *cvcutoff*=⟨x⟩ *dielectric*=⟨x⟩
expandedsampling=yes | no *hbcutoff*=⟨x⟩ *itmax*=⟨n⟩
maxkeep=⟨n⟩ *maxref*=⟨n⟩ *mlcutoff*=⟨x⟩ *posedist*=⟨x⟩
scorecut=⟨x⟩ *short_distance*=anneal | softcore

Options:

<i>cvcutoff</i>	Reject poses if the Coulomb/VDW energy exceeds this value Valid values: reals Default value: 0
<i>dielectric</i>	The dielectric constant used in the grid energy calculation. Valid values: reals Default value: 2 Minimum: 0.9999999999
<i>expandedsampling</i>	Whether or not to use expanded sampling. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>hbcutoff</i>	Reject poses if the H-bond energy exceeds this value Valid values: reals Default value: 0 Maximum: 0.00000001
<i>itmax</i>	The maximum number of conjugate gradient steps Valid values: integers Default value: 100 Minimum: 0
<i>maxkeep</i>	The number of ligand poses on the coarse grid.

Chapter 5: Commands

	Valid values: integers Default value: 5000 Minimum: 1
<i>maxref</i>	The maximum number of poses to keep after rough score refinement. Valid values: integers Default value: 400 Minimum: 1
<i>mlcutoff</i>	Reject poses if the metal-ligand interaction score exceeds this value. Valid values: reals Default value: 10
<i>posedist</i>	The distance criterion for distinct poses. Valid values: reals Default value: 0.5 Minimum: 0.0
<i>scorecut</i>	The window of rough scores for keeping poses. Valid values: reals Default value: 100 Minimum: 0.0
<i>short_distance</i>	The short distance behavior of the potentials. Valid values: anneal softcore Default value: anneal

glidesettings

Allows the settings of some values that determine how the overall Glide job runs.

Syntax:

```
glidesettings enhancedplanarity=yes | no
epikstatepenalties=yes | no function=run | setup
griddirectory=<text> gridfilename=<text> gridfiletype=<text>
maxatom=<n> maxreceptoratom=<n> maxrotbonds=<n>
mode=normal | throughput | accurate numjobs=<n>
numprocs=<n> numreggroup1=<n> numreggroup2=<n>
numreggroup3=<n> numreggroup4=<n> penaltylevel=large |
medium | small receptor=alone | withligand | ligandonly
recyclesteps=<n> reqmodegroup1=all | at least
reqmodegroup2=all | at least reqmodegroup3=all | at least
reqmodegroup4=all | at least showvolumes=yes | no
source=structure | file splitligands=yes | no
testsatisfaction=yes | no volumepenalty=yes | no
writegrid=yes | no writekeyvalue=yes | no
writeexpviz=yes | no
```

Options:

enhancedplanarity

An option which determines if Glide will use the enhanced planarity force field parameters.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

epikstatepenalties

An option which determines if Glide will use the state penalties calculated for the input ligands by Epik as part of the Glide docking score.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

function

Determines whether the Glide job is actually to dock some ligands or just calculate a grid for a receptor.

Valid values: run
setup

Default value: **setup**

griddirectory

The directory where the grid files are located.

Valid values: text strings

Default value:

gridfilename

The base name for the file which the receptor grid is to be written to or read from.

Chapter 5: Commands

	Valid values: text strings Default value:
<i>gridfiletype</i>	The type for the file which the receptor grid is to be written to or read from. The type is either zip or grd . Valid values: text strings Default value: zip
<i>maxatom</i>	Any ligands in the input with more than this number of atoms will be skipped. Valid values: integers Default value: 300 Minimum: 1 Maximum: 300
<i>maxreceptoratom</i>	Receptor in the input with more than this number of atoms will be stopped. Valid values: integers Default value: 80000 Minimum: 1 Maximum: 80000
<i>maxrotbonds</i>	Any ligands in the input with more than this number of rotatable bonds will be skipped. Valid values: integers Default value: 50 Minimum: 1 Maximum: 50
<i>mode</i>	Determines how a number of other values are set to ensure the job runs either faster than normal or more accurate than normal. Valid values: normal throughput accurate Default value: normal
<i>numjobs</i>	Number of jobs to split this docking job into Valid values: integers Default value: 1 Minimum: 1
<i>numprocs</i>	Number of processors to use for each job. Valid values: integers Default value: 1 Minimum: 1

numreqgroup1

Number of constraints to be required for group 1 in docking.

Valid values: integers

Default value: **1**

Minimum: 0

Maximum: 4

numreqgroup2

Number of constraints to be required for group 2 in docking.

Valid values: integers

Default value: **1**

Minimum: 0

Maximum: 4

numreqgroup3

Number of constraints to be required for group 3 in docking.

Valid values: integers

Default value: **1**

Minimum: 0

Maximum: 4

numreqgroup4

Number of constraints to be required for group 4 in docking.

Valid values: integers

Default value: **1**

Minimum: 0

Maximum: 4

penaltylevel

Determine the penalty specification for Glide excluded volumes: large, medium, or small.

Valid values: large
medium
small

Default value: **large**

receptor

Determines whether the structure displayed is just the receptor or a receptor + ligand.

Valid values: alone
withligand
ligandonly

Default value: **alone**

recyclesteps

Specifies the number of ligand recycling iterations in high accuracy docking jobs (glidesettings mode=accurate).

Chapter 5: Commands

Valid values: integers
Default value: **5**
Minimum: 1

reqmodegroup1

The mode determines how to set the number of required constraints for group 1 in docking.

Valid values: all
atleast
Default value: **atleast**

reqmodegroup2

The mode determines how to set the number of required constraints for group 2 in docking.

Valid values: all
atleast
Default value: **atleast**

reqmodegroup3

The mode determines how to set the number of required constraints for group 3 in docking.

Valid values: all
atleast
Default value: **atleast**

reqmodegroup4

The mode determines how to set the number of required constraints for group 4 in docking.

Valid values: all
atleast
Default value: **atleast**

showvolumes

An option which determines if showing excluded volumes or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

source

Where the grid files come from: calculated from the structure or from previously calculated files.

Valid values: structure
file
Default value: **structure**

splitligands

Whether to split the ligand structures into separate files per job

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

testsatisfaction

An option which determines if Glide tests constraint satisfaction after docking or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

volumepenalty

An option which determines if applying excluded volumes penalties or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

writegrid

An option which determines if the grid will be written following the job.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

writekeyvalue

An option which determines whether to write the input file in traditional format or key value format.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

writexpviz

An option which determines if the Glide job will write out a special data file for Glide XP visualizer.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

glidesite

Settings for the site in a glide job.

Syntax:

```
glidesite boxcenter=ligand | residues | supplied boxsize=<n>
          boxsizex=<n> boxsizey=<n> boxsizez=<n>
          display_center_box=yes | no display_enclosing_box=yes | no
          enclosingbox=fitlargest | fitdisplayed | fitlength
          maxliglength=<x> recep_ccut=<x> recep_cscale_per_atom=<x>
          recep_scale_per_atom=none | file | asl recep_vscale=<x>
          recep_vscale_per_atom=<x> useinputcharges=yes | no
          xcent=<x> ycent=<x> zcent=<x>
```

Options:

boxcenter Specifies how the box center is defined - by specifying residues in the active site, by specifying a ligand or by using the values as input by the user.

Valid values: ligand
 residues
 supplied

Default value: **ligand**

boxsize The size of the enclosing box.

Valid values: integers

Default value: **10**

Minimum: 6

Maximum: 40

boxsizex The length in X of the inner box.

Valid values: integers

Default value: **10**

Minimum: 6

Maximum: 40

boxsizey The length in Y of the inner box.

Valid values: integers

Default value: **10**

Minimum: 6

Maximum: 40

boxsizez The length in Z of the inner box.

Valid values: integers

Default value: **10**

Minimum: 6

Maximum: 40

display_center_box

An option which determines if the box which contains the ligand center of motion will be displayed.

Valid values: boolean (true|false; yes|no; y|n; on|off)

- Default value: **true**
- display_enclosing_box* An option which determines if the box which encloses the ligand center of motion box will be displayed.
- Valid values: boolean (true|false; yes|no; y|n; on|off)
- Default value: **true**
- enclosing_box* Specifies how the enclosing box is sized - based on the largest ligand, the displayed ligand or ligand length
- Valid values: fitlargest
fitdisplayed
fitlength
- Default value: **fitdisplayed**
- maxliglength* The maximum ligand length expected.
- Valid values: reals
- Default value: **20**
- Minimum: 0.0000000001
- recep_ccut* The partial atomic charge below which receptor atoms are considered to be non-polar and will have their VDW radii scaled.
- Valid values: reals
- Default value: **0.25**
- Minimum: 0.00000001
- recep_cscale_per_atom* The scaling factor for the Charge of receptor per atoms.
- Valid values: reals
- Default value: **1**
- Minimum: 0.0
- recep_scale_per_atom* An option which specifies the source from which backend will pick the VdW Radius scale factor and Charge scaling factor.
- Valid values: none
file
asl
- Default value: **none**
- recep_vscale* The scaling factor for the VDW radii of non-polar receptor atoms.
- Valid values: reals
- Default value: **1**
- Minimum: 0.00000001

Chapter 5: Commands

recep_vscale_per_atom

The scaling factor for the VDW radii of receptor per atoms.

Valid values: reals

Default value: **1**

Minimum: 0.0

useinputcharges

An option which determines if Impact uses Maestro charges (true) or force field derived ones (false).

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

xcent

The X-coordinate of the enclosing box.

Valid values: reals

Default value: **1e-07**

ycent

The Y-coordinate of the enclosing box.

Valid values: reals

Default value: **1e-07**

zcent

The Z-coordinate of the enclosing box.

Valid values: reals

Default value: **1e-07**

glidetorsionalconstraint

Specifies overall options for torsional constraints used in Glide ligand docking.

Syntax:

```
glidetorsionalconstraint apply=yes | no  
                        displaystructure=yes | no
```

Options:

apply

The flag indicates whether torsional constraints should be applied for Glide ligand docking.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

displaystructure

The flag indicates whether structure for selected pattern should be displayed in Workspace, when pattern has an associated structure and torsional constraints are applied.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

glidetorsionalconstraintallbonds

Specify whether all bonds in the selected SMARTS pattern should be constrained to their current torsion angle values.

Syntax:

```
glidetorsionalconstraintallbonds true|false
```

Operands:

true|false

true to constrain all bonds in the selected SMARTS pattern, false to constrain only the torsions defined in the torsions table for the selected SMARTS pattern.

glidetorsionalconstraintpatternadd

Add a new SMARTS pattern, if pattern is valid. This adds a new row to the patterns table, with 0 torsions defined, and selects it. No check is made for duplication of an existing SMARTS pattern.

Syntax:

```
glidetorsionalconstraintpatternadd <pattern>
```

Operands:

<pattern>

The new SMARTS pattern.

glidetorsionalconstraintpatterndelete

Delete the currently selected pattern row (including its table of torsions).

Syntax:

glidetorsionalconstraintpatterndelete

glidetorsionalconstraintpatterndeleteall

Delete all torsion constraint pattern rows (and all torsion tables associated with them).

Syntax:

glidetorsionalconstraintpatterndeleteall

glidetorsionalconstraintpatterndrag

Move the selected pattern row to the specified row.

Syntax:

glidetorsionalconstraintpatterndrag $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The new row number for the selected SMARTS pattern.

glidetorsionalconstraintpatternreplace

Replace existing pattern for currently selected pattern row, if new SMARTS pattern is valid.

Syntax:

glidetorsionalconstraintpatternreplace $\langle \text{pattern} \rangle$

Operands:

$\langle \text{pattern} \rangle$

The new SMARTS pattern.

glidetorsionalconstraintpatternselect

Select a SMARTS pattern, specified by row index.

Syntax:

```
glidetorsionalconstraintpatternselect <index>
```

Operands:

<index>

The 1-based table row index for the pattern.

glidetorsionalconstrainttorsionadd

Add and select a new torsion row, in the torsion table for the currently selected SMARTS pattern, if the four specified atoms define a proper torsion that does not conflict with any previously constrained torsion.

Syntax:

```
glidetorsionalconstrainttorsionadd <atom1> <atom2>  
                                     <atom3> <atom4>
```

Operands:

<atom1> <atom2> <atom3> <atom4>

Atom numbers, relative to the SMARTS pattern, for four atoms to define the torsion to be constrained.

glidetorsionalconstrainttorsiondelete

Delete the currently selected torsion row of the torsion table for the currently selected pattern.

Syntax:

```
glidetorsionalconstrainttorsiondelete
```

glidetorsionalconstrainttorsiondeleteall

Delete all rows from the torsion table for the currently selected pattern.

Syntax:

```
glidetorsionalconstrainttorsiondeleteall
```

glidetorsionalconstrainttorsiondeleteinvalid

Delete all rows from the torsion table having Status that is not OK, but keep the lowest-numbered torsion for any group of torsions that conflict with each other.

Syntax:

```
glidetorsionalconstrainttorsiondeleteinvalid
```

glidetorsionalconstrainttorsionreplace

Replace atoms for currently selected torsion row, if new atoms define a proper torsion.

Syntax:

```
glidetorsionalconstrainttorsionreplace <atom1> <atom2>  
                                     <atom3> <atom4>
```

Operands:

```
<atom1> <atom2> <atom3> <atom4>
```

Atom numbers, relative to the SMARTS pattern, for four atoms to define the torsion to be constrained.

glidetorsionalconstrainttorsionselect

Select a torsion, specified by row index.

Syntax:

```
glidetorsionalconstrainttorsionselect <index>
```

Operands:

```
<index>
```

The 1-based table row index for the torsion.

glidetorsionalconstrainttorsionsetangle

Specify whether the currently selected torsion should be set to a specified angle, or constrained to its current value.

Syntax:

```
glidetorsionalconstrainttorsionsetangle true|false
    [<angle>]
```

Operands:

```
true|false [<angle>]
```

true to set the angle to a specified value, false to constrain the torsion angle to its current value. This argument can be followed by an angle, in degrees, to specify the value to which the torsion angle should be set.

glideundisplayreceptor

Undisplays the receptor in the Workspace.

Syntax:

```
glideundisplayreceptor
```

googleprofiling

Turn on or off google cpu profiling. This uses the google perftools mechanism and Maestro needs to have been compiled with this feature to work. Otherwise, enabling profiling does nothing. You can specify an output file as the operand. If no operand is specified and you enable profiling, then the file name is set to maestro.perf . If profiling is enabled and you change the profile file name, then the current experiment is terminated, the output flushed to the current file and a new experiment started with the output going to the new file name. The environment variable CPUPROFILE is not used.

Syntax:

```
googleprofiling enable=yes | no
```

Options:

enable Enable or disable google profiling.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

graphicsdebugdisable

Do not allow graphics debug messages to be printed for specified item.

Syntax:

graphicsdebugdisable <item>

Operands:

<item>

If present, we disable debug messages for the specified item. If not present, or item is all , disable debug messages for all items. If all debug messages are enabled, then disabling individual items will have no effect until all is disabled.

graphicsdebugenable

Allow graphics debug messages to be printed for specified item.

Syntax:

graphicsdebugenable <item>

Operands:

<item>

If present, we enable debug messages for the specified item (without disabling previously enabled items). If not present, or item is all , enable debug messages for all items.

grow

Grows the currently selected fragment onto the currently selected grow bond.

Syntax:

grow

growbond

Specifies a grow bond for growing. A subsequent grow command will grow the currently selected fragment onto the this grow bond.

Syntax:

growbond $\langle \text{from_atom} \rangle \langle \text{to_atom} \rangle$

Operands:

$\langle \text{from_atom} \rangle \langle \text{to_atom} \rangle$

The two atoms which define the grow bond. In the growing procedure the “to” atom will be replaced by a suitable atom from the incoming fragment.

growbond2

Specifies a second grow bond for growing. A subsequent grow command will grow the currently selected fragment onto the this grow bond.

Syntax:

growbond2 $\langle \text{from_atom} \rangle \langle \text{to_atom} \rangle$

Operands:

$\langle \text{from_atom} \rangle \langle \text{to_atom} \rangle$

The two atoms which define the grow bond. In the growing procedure the “to” atom will be replaced by a suitable atom from the incoming fragment.

growdirection

Sets the grow direction in the current fragment mode. The direction name must be one of the valid directions for the current fragment mode.

Syntax:

growdirection \langle direction_name \rangle

Operands:

\langle direction_name \rangle

The name of the grow direction which is to be made current. This must be a valid grow direction name within the current fragment.

growname

Set the growname to that specified for all atoms which match the ASL specification.

Syntax:

growname \langle grow_name \rangle \langle ASL \rangle

Operands:

\langle grow_name \rangle \langle ASL \rangle

The first operand is the grow name which is to be applied to the atom. Only the first four characters of the growname will be used. The second operand is the ASL specification for all the atoms which are to have the growname applied.

happly

This is a standard alias for **hydrogenapply** (see [\[hydrogenapply\]](#), page 283).

hbondcriteria

Specify the criteria for calculating H-bonds in the workspace. The acceptor atom (A), H-bonded to the donor Hydrogen (H), is Oxygen, Nitrogen, Sulfur, or Fluorine.

Syntax:

hbondcriteria *acceptorangle*= $\langle x \rangle$ *display*=yes | no *distance*= $\langle x \rangle$
donorangle= $\langle x \rangle$

Options:

acceptorangle

Minimum H — A-R angle, in degrees.

Valid values: reals

Default value: **90**

Minimum: 0.0

Maximum: 180.0

display

This option determines whether H-bond markers will be displayed.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

distance

Maximum separation between donor hydrogen atom and acceptor atom (H — A distance, in Angstroms).

Valid values: reals

Default value: **2.5**

Minimum: 0.0

donorangle

Minimum D-H — A angle, in degrees.

Valid values: reals

Default value: **120**

Minimum: 0.0

Maximum: 180.0

hbondset1

Specify the first set of atoms used in finding H-bonds in the 3D workspace.

Syntax:

hbondset1 $\langle \text{ASL} \rangle$

Operands:

$\langle \text{ASL} \rangle$

A string in the atom specification language. Typical usage is to define **hbondset1** and **hbondset2**. This set, **hbondset1**, defines the “from” atoms. The **hbondset2** atoms define the “to” atoms. H-bonds are calculated between these two sets. That is, the H-bonds are inter-set H-bonds. No intra-set

H-bonds are calculated. If `hbondset2`'s ASL string is empty, then H-bonds are calculated for all atoms in `hbondset1`.

hbondset2

Specify the second set of atoms used in finding H-bonds in the 3D workspace.

Syntax:

hbondset2 \langle ASL \rangle

Operands:

\langle ASL \rangle

A string in the atom specification language. Typical usage is to define `hbondset1` and `hbondset2`. This set, `hbondset2`, defines the “to” atoms. The `hbondset1` atoms define the “from” atoms. H-bonds are calculated between these two sets. That is, the H-bonds are inter-set H-bonds. No intra-set H-bonds are calculated. If `hbondset2`'s ASL string is empty, then H-bonds are calculated for all atoms in `hbondset1`.

help

This is a standard alias for **helpsearch** (see [\[helpsearch\]](#), page 259).

helpauto

Turns on or off the automatic popup help.

Syntax:

helpauto *delay*= \langle n \rangle on|off

Options:

delay The delay (in seconds) before the automatic help is popped up after the mouse enters the widget.

Valid values: integers

Default value: **1**

Operands:

on|off

If the operand is “on” then the popup help will be turned on. If off then it will be turned off.

helpcategory

Changes the current help category

Syntax:

helpcategory <category_name>

Operands:

<category_name>

The name of the category for which the help topics are to be displayed.

helpsearch

Searches the help file for the specified search string. The search can be performed on all text or just the titles.

Syntax:

helpsearch *casesensitive*=yes | no *regexexpression*=yes | no
titleonly=yes | no <search_string>

Options:

casesensitive

A boolean option which determines if the search is performed in a case sensitive manner

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

regexexpression

A boolean option which determines if the search string is to be treated as a regular expression

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

titleonly A boolean option which determines if the search applies to the title of the help topics or the entire text

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

⟨search_string⟩

The text which the help file is to be searched for.

Aliases:

help (see [\[help\]](#), page 258)

helptopic

Loads the topic given by the operand

Syntax:

helptopic ⟨topic_name⟩

Operands:

⟨topic_name⟩

The name of the help topic to be displayed in the help panel. This must match the name of a topic exactly.

hidedockablepanels

Hide all visible dockable panels

Syntax:

hidedockablepanels

hideentries

Hide the selected entried in the project table and create a subset of unselected entries.This function also switches to the subset view.

Syntax:

hideentries

hidemarkers

This command hides the given marker.

Syntax:

hidemarkers ⟨marker_name⟩

Operands:

⟨marker_name⟩

The name of the marker.

hidepanel

Hide the panel whose name is given by the operands.

Syntax:

hidepanel ⟨panel_name⟩

Operands:

⟨panel_name⟩

The name of the panel which is to be hidden. The names available for use in the “hidepanel” command are displayed in parentheses after each item in the main menu bar.

hidepanels

Hide all visible panels

Syntax:

hidepanels

hideproperty

This command hides the given property by creating a subset of all remaining properties in show state. This function also switch to the property subset view.

Syntax:

hideproperty `<propertyname>`

Operands:

`<propertyname>`

The name of the property to hide.

hidetoolbar

Hide the toolbar of given id displayed under given panel.

Syntax:

hidetoolbar `<panel_name:toolbar_id>`

Operands:

`<panel_name:toolbar_id>`

The first part is the name of the panel under which toolbar would be marked as hidden. The name must match to all characters and it is a case sensitive. The second part is the name of a toolbar id which needs to be marked as hidden. Id must match to all characters and it is a case sensitive. Any toolbar can not be displayed under any panel. There is a fix set of toolbars which can be displayed under given panel.

highlight

Create a Workspace highlight

Syntax:

highlight

highlightarrowatoms

Specify the ASL for the named highlight arrow. The ASL specifies the atoms in the Workspace to which the arrow is associated. This can change as the ASL is re-applied when composition of the Workspace changes.

Syntax:

highlightarrowatoms

highlightarrowcolor

Color of the arrow highlight. Default color is yellow

Syntax:

highlightarrowcolor *blue*=⟨x⟩ *green*=⟨x⟩ *red*=⟨x⟩

Options:

<i>blue</i>	Blue component.
	Valid values: reals
	Default value: 0
	Minimum: 0.0
	Maximum: 1.0
<i>green</i>	Green component.
	Valid values: reals
	Default value: 1
	Minimum: 0.0
	Maximum: 1.0
<i>red</i>	Red component.
	Valid values: reals
	Default value: 1
	Minimum: 0.0
	Maximum: 1.0

highlightarrowsettings

Global settings that affect all highlight arrows

Syntax:

highlightarrowsettings *headdirection*=pointing | flat

Options:

headdirection

How the arrow head should be oriented. Currently this is either 'pointing' which means it points in the direction the tail points or 'flat' which means the arrow head is flattened into the screen so it is maximally visible no matter the orientation.

Valid values: pointing
 flat

Default value: **flat**

highlightatoms

Specify the ASL for the named highlight. *selected* can be used to specify the highlight currently selected in the highlights table. The ASL specifies the atoms in the Workspace to which the highlight is associated. This can change as the ASL is re-applied when composition of the Workspace changes

Syntax:

highlightatoms

highlightbackground

Highlight background settings

Syntax:

highlightbackground *type*=automatic | auto | none

Options:

type Background highlight type

Valid values: automatic
 auto
 none

Default value: **automatic**

highlightcolor

Color of the silhouette highlight.

Syntax:

```
highlightcolor alpha=⟨x⟩ blue=⟨x⟩ green=⟨x⟩ red=⟨x⟩
```

Options:

alpha Alpha component. 0 is invisible. 1 is opaque.

Valid values: reals

Default value: **1**

Minimum: 0.0

Maximum: 1.0

blue Blue component.

Valid values: reals

Default value: **0**

Minimum: 0.0

Maximum: 1.0

green Green component.

Valid values: reals

Default value: **1**

Minimum: 0.0

Maximum: 1.0

red Red component.

Valid values: reals

Default value: **1**

Minimum: 0.0

Maximum: 1.0

highlightdelete

Delete the named highlight

Syntax:

```
highlightdelete
```

highlightdirection

Direction in which the highlight points. Only applies to those highlights that have direction.

Syntax:

highlightdirection $x=\langle x \rangle$ $y=\langle x \rangle$ $z=\langle x \rangle$

Options:

x	X component of highlight direction vector Valid values: reals Default value: 0
y	Y component of highlight direction vector Valid values: reals Default value: 0
z	Z direction of the highlight. Valid values: reals Default value: 0

highlightduplicate

Duplicate the named highlight.

Syntax:

highlightduplicate

highlighthide

Hide the highlight

Syntax:

highlighthide

highlightmethod

Create or modify a Workspace highlight.

Syntax:

highlightmethod *type*=silhouette | contrast | saturation

Options:

<i>type</i>	Sets the type of highlight
Valid values:	silhouette contrast saturation
Default value:	silhouette

highlightpreference

Modify various highlight preferences

Syntax:

highlightpreference *contrast*=⟨*x*⟩ *saturation*=⟨*x*⟩

Options:

<i>contrast</i>	Increases or decreases the contrast of the atoms and bonds that are not in any contrast highlight. Atoms and bonds in a contrast highlight definition are left alone.
Valid values:	reals
Default value:	0.55
Minimum:	0.0
Maximum:	1.0
<i>saturation</i>	Increases or decreases the saturation of the atoms and bonds that are not in any saturation highlight. Atoms and bonds in a saturation highlight definition are left alone.
Valid values:	reals
Default value:	0.3
Minimum:	0.0
Maximum:	1.0

highlightramp

Control over a highlight's intensity ramp. Only applies to highlights that use a ramp.

Chapter 5: Commands

Syntax:

highlightramp *ramp*=constant | linear | exponential *step*=⟨x⟩

Options:

ramp Ramp type to use for those highlights that change over time. Those that don't ignore this option (effectively setting it to constant).

Valid values: constant
 linear
 exponential

Default value: **linear**

step Step size for linear ramp. 0 causes the highlight to not cycle and to use its initial color values.

Valid values: reals
Default value: **0.025**
Minimum: 0.0

highlightrename

Rename the specified highlight

Syntax:

highlightrename

highlightselect

Select only the named highlight

Syntax:

highlightselect

highlightshow

Show the highlight

Syntax:

highlightshow

highlightsoff

Hide all ighlights

Syntax:

highlightsoff

highlightson

Display all highlights which are visible

Syntax:

highlightson

highlighttext

Optional text for the highlight to be displayed in the Workspace

Syntax:

highlighttext

highlighttextcolor

Color of the highlight text text. Default color is yellow

Syntax:

highlighttextcolor *blue*=⟨x⟩ *green*=⟨x⟩ *red*=⟨x⟩

Options:

blue Blue component.
Valid values: reals

	Default value: 0
	Minimum: 0.0
	Maximum: 1.0
<i>green</i>	Green component.
	Valid values: reals
	Default value: 1
	Minimum: 0.0
	Maximum: 1.0
<i>red</i>	Red component.
	Valid values: reals
	Default value: 1
	Minimum: 0.0
	Maximum: 1.0

highlighttextfont

Font settings for the highlight text

Syntax:

```
highlighttextfont name=<text> size=<x> style=regular | bold  
                | italic | italicbold
```

Options:

<i>name</i>	Font name for the highlight text (helvetica, etc.)
	Valid values: text strings
	Default value: helvetica
<i>size</i>	Font size for the highlight text
	Valid values: reals
	Default value: 14
	Minimum: 3
	Maximum: 96
<i>style</i>	Font style for the highlight text (regular, bold, etc.)
	Valid values: regular bold italic italicbold
	Default value: regular

highlighttextposition

Position of the optional text for a highlight. If there is text an associated, an arrow is rendered. Its tail is always attached to the text. The head position is controlled by the headx, heady, headz option values. These are automatically set and updated based on the ASL definition of the highlight. If the text is empty, no arrow is drawn.

Syntax:

```
highlighttextposition headx= $\langle x \rangle$  heady= $\langle x \rangle$  headz= $\langle x \rangle$ 
                    x= $\langle x \rangle$  y= $\langle x \rangle$ 
```

Options:

<i>headx</i>	3D X position of the arrow head Valid values: reals Default value: 0
<i>heady</i>	3D Y position of the arrow head Valid values: reals Default value: 0
<i>headz</i>	3D Z position of the arrow head Valid values: reals Default value: 0
<i>x</i>	Fractional X window position of the optional highlight text. Valid values: reals Default value: 0.75 Minimum: 0.0 Maximum: 1.0
<i>y</i>	Fractional Y window position of the optional highlight text. Valid values: reals Default value: 0.5 Minimum: 0.0 Maximum: 1.0

highlighttextselect

Select only the named highlight text

Syntax:

highlighttextselect

highlighttextselectadd

Add the named highlight text to the selected highlight texts

Syntax:

highlighttextselectadd

highlighttextsettings

Modify settings for the named highlight

Syntax:

highlighttextsettings *whenhighlightedatoms*=yes | no

Options:

whenhighlightedatoms

When false always display any text. When true only display text if there are highlighted atoms for the named highlight. This means that there have to be atoms that matched the highlight ASL.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

highlighttextunselect

Unselect the specified highlight text

Syntax:

highlighttextunselect

highlighttextunselectall

Unselect all highlight texts

Syntax:

```
highlighttextunselectall
```

highlightwidth

Width of the border around the bond's line or tube.

Syntax:

```
highlightwidth line=<x> sphere=<x> tube=<x>
```

Options:

<i>line</i>	Border width around the bond when it is in line representation
	Valid values: reals
	Default value: 0.15
	Minimum: 0.01
	Maximum: 0.5
<i>sphere</i>	Border width around an atom sphere (when it is in CPK or Ball representation)
	Valid values: reals
	Default value: 0.15
	Minimum: 0.01
	Maximum: 0.5
<i>tube</i>	Border width around the bond when it is in tube representation
	Valid values: reals
	Default value: 0.15
	Minimum: 0.01
	Maximum: 0.5

historyvisible

Determines which types of commands will be visible in the command history list. Issuing a “historyvisible” alone will cause the list to be rebuilt with the current settings. A historyvisible all will make all commands in the list visible. A historyvisible readwrite=n will make all the read,write,sread and swrite commands hidden. Note that this only affects the visibility of the commands in the command history list of the script editor, not which commands are logged to a file.

Chapter 5: Commands

Syntax:

```
historyvisible 1dplot=yes | no 2dplot=yes | no adjust=yes | no
alias=yes | no application=yes | no atomprop=yes | no
beginundoblock=yes | no bondprop=yes | no build=yes | no
clip=yes | no coloratom=yes | no debug=yes | no
delete=yes | no displayatom=yes | no displayopt=yes | no
entryexport=yes | no entryimport=yes | no
entrywscree=yes | no errorcheck=yes | no ffview=yes | no
find=yes | no fit=yes | no glide=yes | no helpauto=yes | no
helpcategory=yes | no helpsearch=yes | no helptopic=yes | no
historyvisible=yes | no hold=yes | no htreat=yes | no
impact=yes | no labelatom=yes | no liaison=yes | no
ligrep=yes | no macromodel=yes | no measurements=yes | no
monitor=yes | no pause=yes | no pausecommands=yes | no
prefer=yes | no print=yes | no project=yes | no
qikprop=yes | no qsite=yes | no quit=yes | no
rename=yes | no rep=yes | no ribbons=yes | no
rotate=yes | no saveimage=yes | no savelayout=yes | no
saverestoreview=yes | no script=yes | no sets=yes | no
showhide=yes | no spotcenter=yes | no superimpose=yes | no
system=yes | no table=yes | no tile=yes | no
transformation=yes | no translate=yes | no
undoredo=yes | no zoom=yes | no [all]
```

Options:

1dplot

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

2dplot

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

adjust

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

alias

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

application

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

<i>atomprop</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>beginundoblock</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>bondprop</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>build</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>clip</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>coloratom</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>debug</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>delete</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>displayatom</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>displayopt</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>entryexport</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>entryimport</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true

Chapter 5: Commands

entrywscree

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

errorcheck

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

ffview

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

find

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

fit

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

glide

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

helpauto

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

helpcategory

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

helpsearch

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

helptopic

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

historyvisible

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

hold

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

<i>htreat</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>impact</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>labelatom</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>liaison</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>ligrep</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>macromodel</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>measurements</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>monitor</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>pause</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>pausecommands</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>prefer</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>print</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true

Chapter 5: Commands

<i>project</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>qikprop</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>qsite</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>quit</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>rename</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>rep</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>ribbons</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>rotate</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) false
<i>saveimage</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>savelayout</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>saverestoreview</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>script</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true

<i>sets</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>showhide</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>spotcenter</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>superimpose</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>system</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>table</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>tile</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>transformation</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) false
<i>translate</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) false
<i>undoredo</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) true
<i>zoom</i>	Valid values: Default value:	boolean (true false; yes no; y n; on off) false

Operands:

[all]

If “all” is specified then all commands will be visible in the command history list.

hold

The hold command permits holding the current structures off screen.

Syntax:

hold <hold_name>

Operands:

<hold_name>

The name which will be applied to the hold set. This name will reference the hold in subsequent “addfromhold” and “replacefromhold” commands. If the name contains embedded spaces then it must be enclosed in quotation marks.

hppmap

Sets the parameters for hppmap

Syntax:

hppmap *boxmargin*=<x> *boxtype*=box | ligand *cutoff*=<x>
cutoffscheme=atom | neutral *gridspacing*=standard | high
incorporate=append | replace | ignore | appendungrouped |
workspace method=h2o | dipole *source*=selected_entries |
workspace | file

Options:

boxmargin This value is used to increase the size of the box.

Valid values: reals
Default value: **6**
Minimum: 0.

boxtype This option indicates whether to treat the box ASL as defining a box or as defining a ligand.

Valid values: box
ligand
Default value: **box**

cutoff This specifies the cutoff value to use for the hppmap calculation, in angstroms.

Valid values: reals
Default value: **20**
Minimum: 0.

cutoffscheme

This option controls the cutoff scheme.

Valid values: atom
 neutral

Default value: **atom**

gridspacing

This option indicates what type of grid spacing (standard or high) to use for the hppmap.

Valid values: standard
 high

Default value: **standard**

incorporate

This option controls the incorporation of the results.

Valid values: append
 replace
 ignore
 appendungrouped
 workspace

Default value: **replace**

method

This option controls the method.

Valid values: h2o
 dipole

Default value: **dipole**

source

This is the source for the structure.

Valid values: selected_entries
 workspace
 file

Default value: **workspace**

hppmapbox

Sets the box parameters for hppmap

Syntax:

hppmapbox <ASL>

Operands:

<ASL>

The ASL defines the box to use for the hppmap calculation.

hppmapset

This command is used to specify the atoms for which the hppmap calculation will be performed.

Syntax:

hppmapset \langle ASL \rangle

Operands:

\langle ASL \rangle

The ASL defines the set of atoms for which the hppmap calculation is to be performed.

hppmapstart

Starts a hppmap job with the current parameters.

Syntax:

hppmapstart

hppmapwrite

Writes the hppmap input file with the current parameters.

Syntax:

hppmapwrite \langle file name \rangle

Operands:

\langle file name \rangle

Name of the file to write to. If the operand is blank, then use the default file name.

htreat

This is a standard alias for **hydrogentreat** (see [\[hydrogentreat\]](#), page 283).

hydrogenapply

Add or remove hydrogens and lone pairs to the atoms defined by the ASL operand according to the currently set treatment.

Syntax:

hydrogenapply \langle ASL \rangle

Operands:

\langle ASL \rangle

A string in the atom specification language. All atoms which match this specification will have hydrogens added or removed to become consistent with the current treatment.

Aliases:

happly (see [\[happly\]](#), page 256)

hydrogentreat

Choose a treatment for hydrogen atom addition/deletion.

Syntax:

hydrogentreat \langle treatment_name \rangle

Operands:

\langle treatment_name \rangle

The operand is the name of a treatment which is to be made current. Note: treatment names are defined in `$SCHRODINGER/mmshare-vX.X/data/mmhtreat.ini`. There are both long and short names for each treatments and either can be used.

Aliases:

htreat (see [\[htreat\]](#), page 282)

impactbufferedatom

Specifies an atom to be buffered in an Impact calculation.

Syntax:

```
impactbufferedatom <atom_number>
```

Operands:

<atom_number>

The number of an atom which is to be treated as buffered during an Impact calculation.

impactbufferedset

Specifies a set of atoms to be buffered in an Impact calculation.

Syntax:

```
impactbufferedset <ASL>
```

Operands:

<ASL>

The operand must be a valid string in the atom specification language.

impactconstraints

Used to set all the options associated with constraints in impact

Syntax:

```
impactconstraints bonds=yes | no buffer_force=<x>  
                  hmc_bonds=yes | no hmc_solvent=yes | no shake_tolerance=<x>  
                  solvent=yes | no
```

Options:

<i>bonds</i>	If true, then all bonds will be constrained during the Impact calculation.
Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

<i>buffer_force</i>	The buffered atom force constant Valid values: reals Default value: 25
<i>hmc_bonds</i>	If true, then all bonds will be constrained during the HMC calculation. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>hmcsolvent</i>	If true, then all solvent molecules will be held rigid during the HMC calculation. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>shake_tolerance</i>	The tolerance for the SHAKE (or RATTLE) algorithms. Valid values: reals Default value: 1e-07 Minimum: 0.0
<i>solvent</i>	[NOTE: This option is no longer used.] If true, then all solvent molecules will be held rigid during the Impact calculation. Valid values: boolean (true false; yes no; y n; on off) Default value: true

impactcontinuumsolvent

Used to set all the options associated with the continuum solvent in Impact.

Syntax:

```
impactcontinuumsolvent pbf_cutoff=<x> qsitepbf_cutoff=<x>
  qsite_resolution=low | medium | high qsite_sgb_cutoff=<x>
  qsite_type=sgb | pbf | agbnp resolution=low | medium | high
  sgb_cutoff=<x> type=sgb | pbf | agbnp
```

Options:

<i>pbf_cutoff</i>	The type displacement threshold for the PBF calculation. Valid values: reals Default value: 0.1 Minimum: 0.0
-------------------	--

Chapter 5: Commands

qsitepbfcutoff

The type displacement threshold for the PBF calculation.
shake_tolerance=d

Valid values: reals
Default value: **0.1**
Minimum: 0.0

qsiteresolution

The resolution for the PBF solvation calculation.

Valid values: low
medium
high
Default value: **low**

qsitesgbcutoff

The type displacement threshold for the SGB calculation.

Valid values: reals
Default value: **0.1**
Minimum: 0.0

qsitetype

The type of continuum solvent which is to be used in the calculation.

Valid values: sgb
pbf
agbnp
Default value: **sgb**

resolution

The resolution for the PBF solvation calculation.

Valid values: low
medium
high
Default value: **low**

sgbcutoff

The type displacement threshold for the SGB calculation.

Valid values: reals
Default value: **0.1**
Minimum: 0.0

type

The type of continuum solvent which is to be used in the calculation.

Valid values: sgb
pbf
agbnp
Default value: **sgb**

impactdynamics

Settings associated with molecular dynamics simulations in Impact.

Syntax:

```
impactdynamics effectivedensity=⟨x⟩ ensemble=nvt | nve | npt
               inittempgauss=⟨x⟩ initvelo=yes | no isothercomp=⟨x⟩
               numbermdsteps=⟨n⟩ targetpress=⟨x⟩ targettemp=⟨x⟩
               tautemp=⟨x⟩ tauvol=⟨x⟩ timestep=⟨x⟩
               volumescaling=centerofmass | atom
```

Options:

effectivedensity

The effective density used during an NPT Impact dynamics simulation.

Valid values: reals

Default value: **1**

Minimum: 0.00000000000001

ensemble

The target ensemble to be achieved during an Impact dynamics simulation.

Valid values: nvt

nve

npt

Default value: **nvt**

inittempgauss

The temperature used to initialize the velocities from a Gaussian distribution during an NPT Impact dynamics simulation.

Valid values: reals

Default value: **298.15**

Minimum: 0.00000000000001

initvelo

An option which determines if the velocities are to be initialized from a Gaussian distribution during an Impact dynamics simulation.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

isothercomp

The solvent isothermal compressibility to be used during an NPT Impact dynamics simulation.

Valid values: reals

Default value: **4.96e-05**

Minimum: 0.00000000000001

Chapter 5: Commands

numbermdsteps

The total number of time steps to be performed during an Impact dynamics simulation.

Valid values: integers

Default value: **100**

Minimum: 1

targetpress

The target pressure to be used during an NPT Impact dynamics simulation.

Valid values: reals

Default value: **1**

targettemp

The target temperature to be used during a NVT or NPT Impact dynamics simulation.

Valid values: reals

Default value: **298.15**

Minimum: 0.0000000001

tautemp

The temperature relaxation time to be used during a NVT or NPT Impact dynamics simulation.

Valid values: reals

Default value: **0.01**

Minimum: 0.0000000001

tauvol

The volume relaxation time to be used during an NPT Impact dynamics simulation.

Valid values: reals

Default value: **0.01**

Minimum: 0.0000000001

timestep

The time step to be used (in ps) during an Impact dynamics simulation.

Valid values: reals

Default value: **0.001**

Minimum: 0.0000001

volumescaling

The method of volume scaling to be used during an NPT Impact Dynamics simulation.

Valid values: centerofmass

atom

Default value: **centerofmass**

impactfastmultipole

Used to set all the options associated with the fast multipole method in Impact.

Syntax:

```
impactfastmultipole level=<n> maximum=<n>
                    smoothing=yes | no
```

Options:

level The level for the fast multipole method.

Valid values: integers

Default value: **2**

Minimum: 1

maximum The maximum for the fast multipole method.

Valid values: integers

Default value: **7**

Minimum: 4

Maximum: 20

smoothing Whether or not to use smoothing with the fast multipole method.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

impactfrozenatom

Specifies a single atom to be frozen at its current position in an Impact calculation.

Syntax:

```
impactfrozenatom <atom_number>
```

Operands:

<atom_number>

The number of an atom which is to be treated as frozen during an Impact calculation.

impactfrozenet

Specifies a set of atoms to be frozen at their current positions in an Impact calculation.

Syntax:

```
impactfrozenet <ASL>
```

Operands:

<ASL>

A string in the atom specification language which describes the set of atoms which are to be treated as frozen in an Impact calculation.

impacthybridmc

Settings associated with hybrid Monte Carlo simulations in Impact.

Syntax:

```
impacthybridmc ncycles=<n> nmdmc=<n> timestep=<x>
```

Options:

<i>ncycles</i>	The total number of HMC cycles to be performed during an Impact hybrid Monte Carlo simulation. Valid values: integers Default value: 100 Minimum: 1
<i>nmdmc</i>	The number of MD steps per HMC cycle to be performed during an Impact hybrid Monte Carlo simulation. Valid values: integers Default value: 4 Minimum: 1
<i>timestep</i>	The time step to be used (in ps) during an Impact hybrid Monte Carlo simulation Valid values: reals Default value: 0.001 Minimum: 0.0000001

impactjob

This keyword is used to set various options associated with running Impact jobs.

Syntax:

```
impactjob host=<text> incorporate=append | replace | ignore |
appendungrouped | workspace job=<text> login=<text>
structure_source=selected_entries | workspace | file
```

Options:

host The name of the host for the Impact job.

Valid values: text strings

Default value:

incorporate

How the results are to be incorporated into the project. This can be done with replacement of the existing entries, by appending as new entries to the project or by ignoring the final results.

Valid values: append
 replace
 ignore
 appendungrouped
 workspace

Default value: **append**

job The name for the Impact job.

Valid values: text strings

Default value: **impacttmp**

login The login name under which a Impact will be run.

Valid values: text strings

Default value:

structure_source

Whether to use the selected entries in the current project or what is in the workspace as input for the job.

Valid values: selected_entries
 workspace
 file

Default value: **workspace**

impactmdparams

Used to set parameters associated with molecular dynamics simulations in Impact.

Syntax:

```
impactmdparams every=<n> fastfreq=<n> hmcintegrator=verlet |  
                rrespa integrator=verlet | rrespa mdstatistics=yes | no  
                mediumfreq=<n> nprint=<n> slowfreq=<n> stoprot=yes | no  
                traj=yes | no trajfile=<text> trajvelocities=yes | no
```

Options:

every An option which determines how many MD steps between saving a frame of the trajectory for an Impact dynamics simulation.

Valid values: integers

Default value: **5**

Minimum: 1

fastfreq RRESPA frequency for the fast forces.

Valid values: integers

Default value: **4**

Minimum: 1

hmcintegrator

Which integration technique is to be used during the Hybrid Monte Carlo simulation

Valid values: verlet

 rrespa

Default value: **rrespa**

integrator Which integration technique is to be used during the dynamics simulation

Valid values: verlet

 rrespa

Default value: **verlet**

mdstatistics

An option which determines if statistics are to be gathered during an Impact dynamics simulation.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

mediumfreq

RRESPA frequency for the medium forces.

Valid values: integers

	Default value: 2 Minimum: 1
<i>nprint</i>	The number of steps at which output will be written during an Impact dynamics simulation. Valid values: integers Default value: 5 Minimum: 0
<i>slowfreq</i>	RRESPA frequency for the slow forces. Valid values: integers Default value: 1 Minimum: 1
<i>stoprot</i>	An option which determines if overall rotation and translation is stopped during the dynamics simulation. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>traj</i>	An option which determines if the trajectory is to be saved during an Impact dynamics simulation. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>trajfile</i>	The name of the trajectory file. Valid values: text strings Default value: trajectory.trj
<i>trajvelocities</i>	An option which determines if the velocities are to be saved during an Impact dynamics simulation. Valid values: boolean (true false; yes no; y n; on off) Default value: true

impactminimization

Used to set up a minimization in Impact.

Syntax:

Chapter 5: Commands

```
impactminimization algorithm=newton | conjugate | steepest  
                    convergence=eadng | energy | gradient cutoff_1=⟨x⟩  
                    cutoff_2=⟨x⟩ energy_change=⟨x⟩ gradient=⟨x⟩  
                    initial_step_size=⟨x⟩ maximum_step_size=⟨x⟩ ncycles=⟨n⟩  
                    nfull_1=⟨n⟩ nfull_2=⟨n⟩ qsite_energychange=⟨x⟩  
                    qsite_gradient=⟨x⟩ qsite_ncycles=⟨n⟩  
                    qsiteinitialalgorithm=newton | conjugate | steepest
```

Options:

<i>algorithm</i>	Which algorithm is to be used for the minimization
	Valid values: newton conjugate steepest
	Default value: newton
<i>convergence</i>	How convergence is to be established during the minimization.
	Valid values: eadng energy gradient
	Default value: eadng
<i>cutoff_1</i>	Minimization parameter for Truncated Newton. Long range force cutoff.
	Valid values: reals
	Default value: 10
	Minimum: 0.0
<i>cutoff_2</i>	Minimization parameter for Truncated Newton. Long range force cutoff.
	Valid values: reals
	Default value: 10
	Minimum: 0.0
<i>energy_change</i>	The energy change convergence criterion.
	Valid values: reals
	Default value: 1e-07
	Minimum: 0.0
<i>gradient</i>	The gradient convergence criterion.
	Valid values: reals
	Default value: 0.01
	Minimum: 0.0
<i>initial_step_size</i>	The initial steps size for the minimization

	Valid values: reals
	Default value: 0.05
	Minimum: 0.0
<i>maximum_step_size</i>	The maximum step size for the minimization.
	Valid values: reals
	Default value: 1
	Minimum: 0.0
<i>ncycles</i>	The maximum number of minimization cycles to be performed in non-QSite jobs.
	Valid values: integers
	Default value: 100
	Minimum: 0
<i>nfull_1</i>	Minimization parameter for Truncated Newton. Update long range forces every X steps.
	Valid values: integers
	Default value: 10
	Minimum: 0
<i>nfull_2</i>	Minimization parameter for Truncated Newton. Update long range forces every X steps.
	Valid values: integers
	Default value: 10
	Minimum: 0
<i>qsite_energychange</i>	The energy change convergence criterion for QSite.
	Valid values: reals
	Default value: 0.1
	Minimum: 0.0
<i>qsite_gradient</i>	The gradient convergence criterion for QSite.
	Valid values: reals
	Default value: 0.01
	Minimum: 0.0
<i>qsite_ncycles</i>	The maximum number of minimization cycles to be performed for QSite jobs.
	Valid values: integers
	Default value: 1000
<i>qsite_minialgorithm</i>	Which algorithm is to be used for the minimization for qsite jobs.

Valid values:	newton conjugate steepest
Default value:	newton

impactperiodicboundary

Used to set all the options associated with the periodic boundary conditions in Impact.

Syntax:

```
impactperiodicboundary alpha=<x> ewald=yes | no
                        kvectormax=<n> xsize=<x> ysize=<x> zsize=<x>
```

Options:

<i>alpha</i>	The alpha factor for the Ewald long-range correction. Valid values: reals Default value: 0.25 Minimum: 0.0
<i>ewald</i>	Whether or not to use the Ewald long-range correction with the fast multipole method. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>kvectormax</i>	The maximum length of the k-space vectors in the Ewald long-range correction. Valid values: integers Default value: 5 Minimum: 1
<i>xsize</i>	The size of the box in the X-dimension. Valid values: reals Default value: 18.65 Minimum: 18.65
<i>ysize</i>	The size of the box in the Y-dimension. Valid values: reals Default value: 18.65 Minimum: 18.65
<i>zsize</i>	The size of the box in the Z-dimension.

Valid values:	reals
Default value:	18.65
Minimum:	18.65

impactpotential

Set various options associated with the definition of the potential energy to be used in a Impact job.

Syntax:

```
impactpotential continuum_solvent=yes | no dielectric=⟨x⟩
               electrostatics=constant | distance_dependant
               fast_multipole=yes | no field=oplsaa | opsl1999 | opsl2001 |
               opsl2005 | opsl2008 force_field_checks=yes | no
               paramfile=⟨text⟩ partial_charges=yes | no
               periodic_boundary=yes | no qsitecontinuum_solvent=yes | no
               qsitedielectric=⟨x⟩ qsiteelectrostatics=constant |
               distance_dependant qsitefield=oplsaa | opsl1999 | opsl2001 |
               opsl2005 | opsl2008 qsitetruncate=yes | no truncate=yes | no
```

Options:

continuum_solvent

Whether or not to use the continuum solvent

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	false

dielectric The dielectric constant to be used in the electrostatic part of the energy calculation.

Valid values:	reals
Default value:	1
Minimum:	0.9999999999

electrostatics

The electrostatic treatment to be used in the Impact calculation.

Valid values:	constant distance_dependant
Default value:	constant

fast_multipole

Whether or not to use the fast multipole method for electrostatic interactions.

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	false

Chapter 5: Commands

<i>field</i>	The force field to be used for the Impact calculation. Valid values: oplsaa opls1999 opls2001 opls2005 opls2008 Default value: opls2005
<i>force_field_checks</i>	Whether or not to skip force field checks. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>paramfile</i>	The name of the parameter file. Valid values: text strings Default value: paramstd.dat
<i>partial_charges</i>	Whether or not to use atomic partial charges in structure file. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>periodic_boundary</i>	Whether or not to use the periodic boundary conditions Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>qsitecontinuum_solvent</i>	Whether or not to use the continuum solvent Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>qsitedielectric</i>	The dielectric constant to be used in the electrostatic part of the energy calculation. Valid values: reals Default value: 1 Minimum: 0.9999999999
<i>qsiteelectrostatics</i>	The electrostatic treatment to be used in the QSite calculation. Valid values: constant distance_dependant Default value: constant
<i>qsitefield</i>	The force field to be used for the QSite calculation.

Valid values: oplsaa
 opls1999
 opls2001
 opls2005
 opls2008
 Default value: **opls2005**

qsitettruncate

Whether or not to truncate the non-bonded interactions.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

truncate

Whether or not to truncate the non-bonded interactions.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

impactread

Read an Impact input file and set corresponding panels.

Syntax:

impactread <filename>

Operands:

<filename>

The name of an impact input file.

impactrepexch

Settings associated with Replica Exchange in Impact

Syntax:

impactrepexch *dorxmd*=yes | no *nexch*=<n> *nrepl*=<n>
 reftempid=<n> *restart*=yes | no *restartfile*=<text>
 targettemp=<text>

Options:

dorxmd Set the mode to Replica Exchange

Valid values: boolean (true|false; yes|no; y|n; on|off)

Chapter 5: Commands

	Default value: false
<i>nexch</i>	The number of MD steps between Exchanges Valid values: integers Default value: 250 Minimum: 1
<i>nrepl</i>	The number of Replicas Valid values: integers Default value: 2 Minimum: 2 Maximum: 100
<i>reftempid</i>	The index of temperature of interest Valid values: integers Default value: 1 Minimum: 1
<i>restart</i>	Set the mode to restart a previous Replica Exchange job. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>restartfile</i>	Restart input file name. Valid values: text strings Default value:
<i>targettemp</i>	List of target temperatures separated by commas Valid values: text strings Default value: 300.00, 350.00

impactselectextendtablerow

Extends the selection to this row in the given orbital table in the Surfaces Property for Impact.

Syntax:

```
impactselectextendtablerow table=alpha | beta <row>
```

Options:

<i>table</i>	Indicates whether to select a row in the alpha or beta orbital table. Valid values: alpha beta
--------------	--

Default value: **alpha**

Operands:

⟨row⟩

The row number to extend the select to.

impactselectonlytablerow

Selects only this row in the given orbital table in the Surface Property for Impact.

Syntax:

`impactselectonlytablerow table=alpha | beta ⟨row⟩`

Options:

<i>table</i>	Indicates whether to select a row in the alpha or beta orbital table.
Valid values:	alpha beta
Default value:	alpha

Operands:

⟨row⟩

The row number to select only in the table row.

impactstart

Start an Impact input file with the current settings.

Syntax:

`impactstart`

impacttask

Determines which impact task is currently being set up.

Syntax:

impacttask soak | mini | dynamics | hmc | qsite | glide | liaison

Operands:

soak | mini | dynamics | hmc | qsite | glide | liaison

The type of impact task which is to be set up. The operand cannot be abbreviated and must be given in full.

impacttruncation

Used to set all the options associated with the truncation of non-bonded interactions in Impact.

Syntax:

impacttruncation *distance*=⟨x⟩ *qsitedistance*=⟨x⟩
qsiteupdatefrequency=⟨n⟩ *updatefrequency*=⟨n⟩

Options:

distance The truncation distance for residue-based cutoffs.

Valid values: reals
Default value: **12**
Minimum: 0.0

qsitedistance

The truncation distance for residue-based cutoffs.

Valid values: reals
Default value: **12**
Minimum: 0.0

qsiteupdatefrequency

The number of steps between update of the neighbor list.

Valid values: integers
Default value: **10**
Minimum: 1

updatefrequency

The number of steps between update of the neighbor list.

Valid values: integers
Default value: **10**
Minimum: 1

impactwrite

Write an Impact input file with the current settings.

Syntax:

```
impactwrite
```

invert

Inverts the chirality around a chiral atom.

Syntax:

```
invert <chiral_atom> <non_moving_atom1> <non_moving_atom2>
```

Operands:

```
<chiral_atom> <non_moving_atom1> <non_moving_atom2>
```

Three atom numbers. The first is the atom around which the chirality is to be inverted. The second and third are atoms which are attached to the chiral atom but are not to be moved in the inversion process.

invertset

Inverts all chiral centers in the specified set of atoms.

Syntax:

```
invertset <ASL>
```

Operands:

```
<ASL>
```

An ASL specification of the atoms which are to have their chiral centers inverted. Because of the way this works, this set of atom should involve at least whole molecules.

jaguarassignatomnames

Sets the atom names of the matching atoms as per the jaguar standards i.e, element name + atom number. Incase of any duplicates, integer part of the name is incremented till it is unique in the Workspace

Syntax:

```
jaguarassignatomnames <ASL>
```

Operands:

<ASL>

The ASL expression describing the set of atoms which are to have their atom names made unique.

jaguarimportgeometry

Import structures into the current project from a Jaguar input file.

Syntax:

```
jaguarimportgeometry all=yes | no end=yes | no start=<n>  
total=<n> wsinclude=none | first | all wsreplace=yes | no  
[<filename>]
```

Options:

<i>all</i>	This determines whether all structures will be imported, or just a specified range. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>end</i>	This determines if all structures in the file are to be imported starting from the structure specified by start. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>start</i>	This option sets the number of the first structure to be imported, if not importing all. Valid values: integers Default value: 1 Minimum: 1
<i>total</i>	The total number of structures to be imported from the file, if not importing all structures.

	Valid values: integers
	Default value: 1
	Minimum: 1
<i>wsinclude</i>	This option determines which of the imported structures are to be included in the workspace. Valid values are “none”, “first”, or “all”.
	Valid values: none first all
	Default value: none
<i>wsreplace</i>	This determines whether the structures currently in the workspace will be replaced by the included imported structures, or whether they will be included in the workspace also.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true

Operands:

[<filename>]

The name of the file from which structures will be imported. If no name is specified, then no import will be done.

jaguarinputfilesjob

This keyword is used to set various options associated with running Jaguar input files jobs.

Syntax:

jaguarinputfilesjob *input_files*=<text>

Options:

<i>input_files</i>	Specifies existing structure file to be used as input files for running the job
	Valid values: text strings
	Default value:

jaguarjob

This keyword is used to set various options associated with running Jaguar jobs.

Syntax:

```
jaguarjob structure_files=<text> structure_source=selected_entries  
          | workspace | file
```

Options:

structure_files

If structure_source is Selected structure files , specifies existing structure file to be used as input files for running the job.

Valid values: text strings

Default value:

structure_source

Whether to use the selected entries in the current project or what is in the workspace as input for the job.

Valid values: selected_entries
workspace
file

Default value: **workspace**

jaguarselectextendtablerow

Extends the selection to this row in the given orbital table in the Surfaces Property for Jaguar.

Syntax:

```
jaguarselectextendtablerow table=alpha | beta <row>
```

Options:

table Indicates whether to select a row in the alpha or beta orbital table.

Valid values: alpha
beta

Default value: **alpha**

Operands:

⟨row⟩

The row number to extend the select to.

jaguarselectonlytablerow

Selects only this row in the given orbital table in the Surface Property for Jaguar.

Syntax:

```
jaguarselectonlytablerow table=alpha | beta ⟨row⟩
```

Options:

table Indicates whether to select a row in the alpha or beta orbital table.

Valid values: alpha
 beta

Default value: **alpha**

Operands:

⟨row⟩

The row number to select only in the table row.

jobcleanup

Specifies a job to cleanup.

Syntax:

```
jobcleanup files=jobandmonitor | all ⟨job_id⟩
```

Options:

files This option determines which files are removed when the job is cleaned up. Either just the job record and monitoring files or all files associated with the job.

Valid values: jobandmonitor
 all

Default value: **jobandmonitor**

Operands:

`<job_id>`

The ID of the job which is to be cleaned up.

jobsettings

This keyword is used to set various options associated with running a backend job.

Syntax:

```
jobsettings compress=yes | no directory=<text>
               distributesubjobs=maxprocessors | specifiedprocessors
               gpuindices=<text> host=<text> hostlist=<text>
               hostsubjobs=<n> incorporate=appendentries | replaceentries |
               ignoreentries | ingoreentries | appendentriesungrouped |
               workspace jobname=<text> login=<text> maxstructures=<n>
               numcpus=<n> numsubjobs=<n> subjobsprocessors=<n>
               title=<text> tmpdir=<text> writedirectory=<text>
               <model_name>
```

Options:

- | | |
|--------------------------|---|
| <i>compress</i> | An option which determines if output files are generated in an mmzip archive.
Valid values: boolean (true false; yes no; y n; on off)
Default value: true |
| <i>directory</i> | The directory to write the output to.
Valid values: text strings
Default value: |
| <i>distributesubjobs</i> | How the subjobs are to be distributed over the available processors
Valid values: maxprocessors
specifiedprocessors
Default value: maxprocessors |
| <i>gpuindices</i> | The GPU indices for running jobs.
Valid values: text strings
Default value: |
| <i>host</i> | The name of the host for the backend job. |

	Valid values: text strings Default value: host
<i>hostlist</i>	The string of hosts and number of CPUs which will be passed to the job. Valid values: text strings Default value:
<i>hostsubjobs</i>	The number of subjobs this current job is to be run as. Valid values: integers Default value: 1 Minimum: 1
<i>incorporate</i>	How the results are to be incorporated into the project. This can be done by appending as new entries to the project, replacing entries, or by ignoring the final results. Valid values: appendentries replaceentries ignoreentries ingoreentries appendentriesungrouped workspace Default value: appendentries
<i>jobname</i>	The name for the backend job. Valid values: text strings Default value:
<i>login</i>	The login name under which the job will be run. Valid values: text strings Default value: user
<i>maxstructures</i>	The maximum number of structures to be incorporated into the project. Valid values: integers Default value: 100 Minimum: 0
<i>numcpus</i>	The number of CPUs to be used for jobs that support multiple CPUs. Valid values: integers Default value: 1 Minimum: 1

Chapter 5: Commands

numsubjobs

The number of subjobs this current job is to be run as.

Valid values: integers

Default value: **1**

Minimum: 1

subjobsprocessors

The number of processors to distribute the subjobs over

Valid values: integers

Default value: **1**

Minimum: 1

title

Set the input entry's title.

Valid values: text strings

Default value:

tmpdir

The path name of the scratch directory.

Valid values: text strings

Default value:

writedirectory

The directory to write job input files to.

Valid values: text strings

Default value:

Operands:

`<model_name>`

The name of a model which has job settings like host, login, and so on.

kill

This is a standard alias for **energykill** (see [\[energykill\]](#), page 134).

labelatom

Displays labels for the atoms specified by the ASL operand.

Syntax:


```

labelatom 1charge=yes | no 1chargeformatstring=<text>
          2charge=yes | no 2chargeformatstring=<text>
          acolor=yes | no anum=yes | no anumformatstring=<text>
          atomicnumber=yes | no atomname=yes | no
          atomnameformatstring=<text> atype=yes | no
          atypeformatstring=<text> chain=yes | no
          chainformatstring=<text> chirality=yes | no
          chiralityformatstring=<text> cindex=<n> color=<text>
          compositionfields=<text> dmsopka=yes | no
          dmsopkaformatstring=<text> element=yes | no
          elementformatstring=<text> entryname=yes | no
          entrynameformatstring=<text>
          entrypropertyformatstring=<text> entrypropertynames=<text>
          font=<text> font_size=<n> font_style=normal | italic | bold |
          bolditalic formalcharge=yes | no
          formalchargeformatstring=<text> grownname=yes | no
          grownnameformatstring=<text> h2opka=yes | no
          h2opkaformatstring=<text> headings=yes | no
          inscode=yes | no keeplabels=on_top | with_atom
          mode=replace | append | clear molnum=yes | no
          molnumentry=yes | no molnumentryformatstring=<text>
          molnumformatstring=<text> numentry=yes | no
          numentryformatstring=<text> nummol=yes | no
          nummolformatstring=<text> occupancy=yes | no
          occupancyformatstring=<text> oneletter=yes | no
          pdbbfactor=yes | no pdbbfactorformatstring=<text>
          pdbname=yes | no pdbnameformatstring=<text>
          reapply=yes | no reapplylabels=yes | no resname=yes | no
          resnameformatstring=<text> resnum=yes | no
          resnumformatstring=<text> separator=<text>
          showlabel=yes | no stereochemistry=yes | no
          stereochemistryformatstring=<text> title=yes | no
          titleformatstring=<text> user=yes | no
          useratompropertyformatstring=<text> utext=<text>
          xoffset=<x> xyz=yes | no xyzformatstring=<text>
          yoffset=<x> <ASL>

```

Options:

1charge A boolean option which determines if the charge 1 value will be included in the atom label

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

Chapter 5: Commands

1chargeformatstring

A string which determines which format is to be used for the charge 1 value labels

Valid values: text strings

Default value: **%C1.3**

2charge

A boolean option which determines if the charge 2 value will be included in the atom label

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

2chargeformatstring

A string which determines which format is to be used for the charge 2 value labels

Valid values: text strings

Default value: **%C2.3**

acolor

A boolean option which determines if labels will be colored the same as the atoms they label

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

anum

A boolean option which determines if the atom numbers will be included in the label

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

anumformatstring

A string which determines which format is to be used for the atom number labels

Valid values: text strings

Default value: **%NU**

atomicnumber

A boolean option which determines if the atomic number will be included in the atom label

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

atomname

A boolean option which determines if the atom name will be included in the atom label

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

atomnameformatstring

A string which determines which format is to be used for the atom name labels

	Valid values: text strings Default value: %AT
<i>atype</i>	A boolean option which determines if the MacroModel atom types will be included in the label Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>atypeformatstring</i>	A string which determines which format is to be used for the MacroModel atom type labels Valid values: text strings Default value: %TY
<i>chain</i>	A boolean option which determines if the chain name will be included in the label Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>chainformatstring</i>	A string which determines which format is to be used for chain name labels Valid values: text strings Default value: %CH
<i>chirality</i>	A boolean option which determines if the atom chirality (R or S) will be included in the label. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>chiralityformatstring</i>	A string which determines which format is to be used for the atom chirality labels Valid values: text strings Default value: %CY
<i>cindex</i>	An integer which indicates color index which is to be used for the labels. This will be ignored unless the <i>acolor</i> option is off Valid values: integers Default value: 2 Minimum: 1 Maximum: 256
<i>color</i>	A string which is the color name for atom labels. This will be ignored unless the <i>acolor</i> option is off Valid values: text strings Default value:

Chapter 5: Commands

compositionfields

This option sets the names of the properties (semicolon separated list) to be used to compose the labels for selected atoms.

Valid values: text strings

Default value:

dmsopka

A boolean option which determines if the pka value in DMSO solvent will be included in the atom label

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

dmsopkaformatstring

A string which determines which format is to be used for the labels of pka value in DMSO solvent

Valid values: text strings

Default value: **%DP**

element

A boolean option which determines if the element symbol will be included in the atom label

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

elementformatstring

A string which determines which format is to be used for the element symbol labels

Valid values: text strings

Default value: **%EL**

entryname

A boolean option which determines if the entry name will be included in the atom label

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

entrynameformatstring

A string which determines which format is to be used for the entry name labels

Valid values: text strings

Default value: **%EN**

entrypropertyformatstring

A string which determines which format is to be used for the ct property labels

Valid values: text strings

Default value: **%EP**

entrypropertynames

This option sets the names of the entry properties (semicolon separated list) to be used to compose the labels for selected atoms.

Valid values: text strings

Default value:

font

A string which determines which font is to be used for the labels

Valid values: text strings

Default value:

font_size

An integer determines which font size to be used for the labels

Valid values: integers

Default value: **14**

Minimum: 0

font_style

An integer determines which font style to be used for the labels

Valid values: normal
italic
bold
bolditalic

Default value: **bold**

formalcharge

A boolean option which determines if the formal charge will be included in the atom label

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

formalchargeformatstring

A string which determines which format is to be used for the formal charge labels

Valid values: text strings

Default value: **%FC**

growname

A boolean option which determines if the grow name will be included in the atom label

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

grownameformatstring

A string which determines which format is to be used for the grow name labels

Valid values: text strings

Default value: **%GN**

h2opka

A boolean option which determines if the pka value in H2O solvent will be included in the atom label

Chapter 5: Commands

	Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>h2opkaformatstring</i>	A string which determines which format is to be used for the labels of pka value in H2O solvent Valid values: text strings Default value: %HP
<i>headings</i>	A boolean option which determines if the label fields will have headers in the labels (anum=1, atype=C2 and so on). Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>inscode</i>	A boolean option which determines if the pdb residue insertion code will be included in the label string Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>keeplabels</i>	This option controls whether labels are kept on top of atoms or whether they are kept at the atom level. Keeping the labels at atom level means labels will be obscured by other atoms and objects in the Workspace. Valid values: on_top with_atom Default value: on_top
<i>mode</i>	Specifies one of two modes: 1) replace: replace the existing label string with the new one. 2) append: append to the existing label string Valid values: replace append clear Default value: replace
<i>molnum</i>	A boolean option which determines if the molecule number will be included in the label Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>molnumentry</i>	A boolean option which determines if the molecule number by entry will be included in the label Valid values: boolean (true false; yes no; y n; on off) Default value: false

molnumentryformatstring

A string which determines which format is to be used for molecule number by entry labels

Valid values: text strings

Default value: **%ME**

molnumformatstring

A string which determines which format is to be used for the molecule number labels

Valid values: text strings

Default value: **%MN**

numentry

A boolean option which determines if the atom number by entry will be included in the label

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

numentryformatstring

A string which determines which format is to be used for the atom number by entry labels

Valid values: text strings

Default value: **%NE**

nummol

A boolean option which determines if the atom number by molecule will be included in the label

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

nummolformatstring

A string which determines which format is to be used for the atom number by molecule labels

Valid values: text strings

Default value: **%NM**

occupancy

A boolean option which determines if the partial occupancy data will be included in the label string

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

occupancyformatstring

A string which determines which format is to be used for the partial occupancy data labels

Valid values: text strings

Default value: **%OC**

oneletter

A boolean option which determines whether single letter pdb residue name is displayed or three letters.

Chapter 5: Commands

	Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>pdbbfactor</i>	A boolean option which determines if the pdb bfactor will be included in the label string Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>pdbbfactorformatstring</i>	A string which determines which format is to be used for the pdb bfactor labels Valid values: text strings Default value: %BF
<i>pdbname</i>	A boolean option which determines if the pdb atom name will be included in the label string Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>pdbnameformatstring</i>	A string which determines which format is to be used for pdb name labels Valid values: text strings Default value: %PA
<i>reapply</i>	A boolean option which determines if reapplying current atom labels when Workspace changes, by using last used option in the main toolbar Label atoms button menu Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>reapplylabels</i>	A boolean option which determines if reapplying current atom labels when Workspace changes, by using the options specified in the Atom Labels panel Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>resname</i>	A boolean option which determines if the pdb residue name will be included in the label string Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>resnameformatstring</i>	A string which determines which format is to be used for pdb residue name labels Valid values: text strings Default value: %RT

<i>resnum</i>	<p>A boolean option which determines if the residue number will be included in the label string</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>
<i>resnumformatstring</i>	<p>A string which determines which format is to be used for the residue number labels</p> <p>Valid values: text strings</p> <p>Default value: %RN</p>
<i>separator</i>	<p>A string option which is the inserted between each field in the label string.</p> <p>Valid values: text strings</p> <p>Default value:</p>
<i>showlabel</i>	<p>A boolean option which determines if an atom label is displayed when the atom itself is undisplayed.</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>
<i>stereochemistry</i>	<p>A boolean option which determines if the double bond stereochemistry (E or Z) will be included in the label.</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>
<i>stereochemistryformatstring</i>	<p>A string which determines which format is to be used for the stereochemistry labels</p> <p>Valid values: text strings</p> <p>Default value: %ST</p>
<i>title</i>	<p>A boolean option which determines if the entry title will be included in the atom label</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>
<i>titleformatstring</i>	<p>A string which determines which format is to be used for the entry title labels</p> <p>Valid values: text strings</p> <p>Default value: %ET</p>
<i>user</i>	<p>A boolean option which determines if the user-defined text will be included in the label</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>

Chapter 5: Commands

useratompropertyformatstring

A string which determines which format is to be used for the user specified atom property labels

Valid values: text strings

Default value: **%UA**

utext

This option sets the user defined text for atom labels

Valid values: text strings

Default value:

xoffset

A double option which sets the x-offset from the atom center at which the labels are to appear.

Valid values: reals

Default value: **0.05**

xyz

A boolean option which determines if the xyz location of the atom will be included in the label.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

xyzformatstring

A string which determines which format is to be used for the xyz location labels

Valid values: text strings

Default value: **%XY**

yoffset

A double option which sets the y-offset from the atom center at which the labels are to appear.

Valid values: reals

Default value: **-0.15**

Operands:

⟨ ASL ⟩

A string in the atom specification language. All atoms which match this specification will be labeled with the current label settings.

labelatomfieldsnone

Allows user to turn off all the atom label fields in one simple command instead of having to issue a labelatom command with all of the options set to false: anum=false etc.

Syntax:

labelatomfieldsnone

labelatomoffset

Set the offset for the specified atoms. If not specified, then Maestro uses the offset specified in preferences (i.e. last value specified to the labelatom command).

Syntax:

labelatomoffset $x=\langle x \rangle$ $y=\langle x \rangle$ $\langle \text{ASL} \rangle$

Options:

x	Specifies x offset
	Valid values: reals
	Default value: 0
y	Specifies y offset
	Valid values: reals
	Default value: 0

Operands:

$\langle \text{ASL} \rangle$

A string in the atom specification language. All atoms which match this specification will have their offset modified.

labelatomoffsetreset

Remove any customized label offsets for the specified atoms. Note: any entry included after this command that has customized label offsets will display with those customized offsets. This command only removes offsets for those labels (structures) in the Workspace.

Syntax:

labelatomoffsetreset

labelclear

Clears labels from the atoms specified by the ASL operand.

Syntax:

labelclear \langle ASL \rangle

Operands:

\langle ASL \rangle

A string in the atom specification language. All atoms which match this specification will have their labels removed.

labeldisplayusertext

Displays user text labels for the atoms specified by the ASL operand.

Syntax:

labeldisplayusertext \langle ASL \rangle

Operands:

\langle ASL \rangle

A string in the atom specification language. All atoms which match this specification will have their previously defined user text label shown in the workspace.

labelformat

Sets the label format string which will be used in subsequent labeling operations.

Syntax:

labelformat \langle format_string \rangle

Operands:

\langle format_string \rangle

A string which defines the current labeling format. See on-line help for details.

labelupdate

Apply the current labelatom format to all the atoms that have a non-empty atom label format. This does not apply to atoms with a heteroatom element label format that is generated automatically from the elementlabels scheme (i.e., to Hetatom Labels). The update is based on the atom label format, regardless of whether any non-blank label is displayed for the atom in the Workspace.

Syntax:

labelupdate [NON_EL]

Operands:

[NON_EL]

An optional string, to indicate that only atom labels, not elementlabels, are relabelled. This argument is now ignored.

labelusertextatomupdate

Update the user text for the specified atoms

Syntax:

labelusertextatomupdate *utext*=⟨text⟩ ⟨ASL⟩

Options:

utext This option sets the user defined text for atom relabeling
 Valid values: text strings
 Default value:

Operands:

⟨ASL⟩

ASL of atoms to relabel

liaisonanalysis

Used to do the analysis

Chapter 5: Commands

Syntax:

```
liaisonanalysis liaisonalpha=⟨x⟩ liaisonanalysistype=fit |  
predict liaisonbefile=⟨text⟩ liaisonbeta=⟨x⟩  
liaisonenerymodeltype=equation | glidescore  
liaisonfitname=⟨text⟩ liaisonfixalpha=yes | no  
liaisonfixbeta=yes | no liaisonfixgamma=yes | no  
liaisonfixintercept=yes | no liaisonfixslope=yes | no  
liaisongamma=⟨x⟩ liaisonintercept=⟨x⟩ liaisonlistfile=⟨text⟩  
liaisonslope=⟨x⟩ liaisonspecligs=readligs | enterligs
```

Options:

liaisonalpha

This option sets the van der Waals (Alpha) coefficient.

Valid values: reals

Default value: **0**

*liaisonanalysis*type

This option determines the analysis type, including two options of Predict (predict) and Fit (fit).

Valid values: fit
predict

Default value: **fit**

liaisonbefile

This option determines the name of the ligand binding energy file

Valid values: text strings

Default value:

liaisonbeta

This option sets the electrostatic (Beta) coefficient.

Valid values: reals

Default value: **0**

liaisonenerymodeltype

This option determines which binding energy model is used, LIA equation or Glidescore

Valid values: equation
glidescore

Default value: **equation**

liaisonfitname

This option determines the names of ligands to be predicted

Valid values: text strings

Default value:

liaisonfixalpha

Allows the alpha value to be used as a constraint.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

liaisonfixbeta

Allows the beta value to be used as a constraint.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

liaisonfixgamma

Allows the gamma value to be used as a constraint

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

liaisonfixintercept

Allows the Glidescore intercept to be used as a constraint

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

liaisonfixslope

Allows the glidescore slope to be used as a constraint.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

liaisongamma

This option sets cavity (Gamma) coefficient.

Valid values: reals

Default value: **0**

liaisonintercept

This option sets the Glidescore intercept

Valid values: reals

Default value: **0**

liaisonlistfile

This option determines the name of file that contains names of ligands to be predicted

Valid values: text strings

Default value:

liaisonslope

This option sets the Glidescore slope

Valid values: reals

Default value: **0**

liaisonspecligs

This option determines the methods used to specify the ligands either by reading ligand names from a text file (readligs), or entering a comma separated list of ligand names (enterligs).

Valid values: readligs
 enterligs

Default value: **readligs**

liaisonparameters

Used to specify simulation parameters in a Liaison (Linear Interaction Approximation).

Syntax:

```
liaisonparameters algorithm=newton | conjugate | steepest
liaisonboundgcut=<x> liaisonboundheat=<x>
liaisonboundint=<n> liaisonboundmini=<n>
liaisonboundprod=<x> liaisonboundpremini=<n>
liaisonfreeint=<n> liaisonliggcut=<x> liaisonligheat=<x>
liaisonligmini=<n> liaisonligpremini=<n> liaisonligprod=<x>
liaisonmethod=mini | hmc | dyn liaisonreltime=<x>
liaisonrescut=<x> liaisontemp=<x> minibound=yes | no
minilig=yes | no
```

Options:

algorithm Which algorithm is to be used for the minimization

Valid values: newton
 conjugate
 steepest

Default value: **newton**

liaisonboundgcut

This option determines the RMS gradient for convergence (in kcal/mol/Å)

Valid values: reals

Default value: **0.05**

liaisonboundheat

This option determines the time for heating/ equilibration (ps)

Valid values: reals

Default value: **5**

Minimum: 0.0

liaisonboundint

This option determines the number of steps for collecting Liaison statistics

Valid values: integers
 Default value: **10**
 Minimum: 0

liaisonboundmini

This option determines the maximum number of minimization steps

Valid values: integers
 Default value: **500**
 Minimum: 0

liaisonboundpremini

This option determines the maximum number of minimization steps

Valid values: integers
 Default value: **500**
 Minimum: 0

liaisonboundprod

This option determines the time for data collection (in ps)

Valid values: reals
 Default value: **5**
 Minimum: 0.0

liaisonfreeint

This option determines the number of steps for collecting Liaison statistics

Valid values: integers
 Default value: **10**
 Minimum: 0

liaisonliggcut

This option determines the RMS gradient for convergence

Valid values: reals
 Default value: **0.01**
 Minimum: 0.0

liaisonligheat

This option determines the time for heating/ equilibration (ps)

Valid values: reals
 Default value: **5**
 Minimum: 0.0

Chapter 5: Commands

liaisonligmini

This option determines the maximum number of minimization steps

Valid values: integers
Default value: **500**
Minimum: 0

liaisonligpremini

This option determines maximum number of minimization steps

Valid values: integers
Default value: **1000**
Minimum: 0

liaisonligprod

This option determines the time for data collection (in ps)

Valid values: reals
Default value: **5**
Minimum: 0.0

liaisonmethod

This option determines which sampling method will be used during a Liaison analysis. The 3 options are minimization (mini), Hybrid Monte Carlo (hmc) and Molecular dynamics (dyn).

Valid values: mini
hmc
dyn
Default value: **mini**

liaisonreltime

This option determines the temperature relaxation time

Valid values: reals
Default value: **0.01**
Minimum: 0.0

liaisonrescut

This option determines the residue-based cutoff for nonbonded interactions

Valid values: reals
Default value: **15**
Minimum: 0.0

liaisontemp

This option determines the simulation temperature

Valid values: reals
Default value: **300**
Minimum: 0.0

minibound

This option determines whether the ligand- receptor complex is minimized before running simulation.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

minilig

This option determines whether the ligand is minimized before running the simulation

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

liaisonselectlig

Defines a on-screen molecule to be treated as the ligand for a Liaison calculation

Syntax:

liaisonselectlig <molecule_number>

Operands:

<molecule_number>

The number of a molecule to be included as the ligand.

liaisonsettings

Used to set values associated with the Linear Interaction Approximation (Liaison)

Syntax:

liaisonsettings *liaisonsubdir*=<text> *liaison*type=sim | analyze
numproc=<n>

Options:

liaisonsubdir

This option determines the subdirectory to use for Liaison jobs.

Valid values: text strings

Default value: **liaison**

liaisontype

This option determines the job type for Liaison analysis. There are two types: (1) Simulate (sim); (2) Analyze results of earlier simulations (analyze)

Valid values: sim
 analyze

Default value: **sim**

numproc

This option determines the number of processors the job is to be run on.

Valid values: integers

Default value: **1**

Minimum: 1

liaisonsystem

Used to specify receptor and ligand(s) to be simulated in the Linear Interaction Approximation (Liaison)

Syntax:

```
liaisonsystem haslialig=yes | no liaisonligfmt=m2io | sd | mol2 |  
              pdb liaisonligfname=<text> liaisonsimtype=mult | single  
              liaisonsname=<text> ligftype=ligstruct | liglist  
              usepartialcharges=yes | no
```

Options:

haslialig

This option determines whether the displayed structure includes the ligands

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

liaisonligfmt

This option for selecting the format of the ligand file. The four possible formats are MAESTRO (m2io), MDL SD (sd), MOL2 (mol2) and PDB (pdb)

Valid values: m2io
 sd
 mol2
 pdb

Default value: **m2io**

liaisonligfname

This option determines the name of the ligand file

Valid values: text strings
 Default value:

liaisonsimtype

This option determines the type of ligand(s)-receptor structures to be simulated, having two options: (1) Multiple ligands, single receptor; (2) Single ligand, single receptor.

Valid values: mult
 single
 Default value: **mult**

liaisonsname

This option determines the name to be used for this ligand

Valid values: text strings
 Default value: **Lig**

ligftype

This option determines the source of ligand structure(s), having two options: (1) File containing a single ligand structure (`ligstruct`); and (2) File containing a list of ligand names and structure files (`liglist`).

Valid values: `ligstruct`
`liglist`
 Default value: **ligstruct**

usepartialcharges

This option determines whether to use partial charges from the Maestro file

Valid values: boolean (`true|false`; `yes|no`; `y|n`; `on|off`)
 Default value: **false**

ligandbond

Specifies two atoms which define a bond to be treated as a ligand bond during a conformational search.

Syntax:

ligandbond `<atom1>` `<atom2>`

Operands:

`<atom1>` `<atom2>`

The numbers of two atoms which are to be treated as a ligand bond during a conformational search. Note that specifying a-b is the same as specifying b-a.

light

Specifies a light or set light options.

Syntax:

```
light ambient= $\langle n \rangle$  diffuse= $\langle n \rangle$  specular= $\langle n \rangle$  use=yes | no  
x= $\langle x \rangle$  y= $\langle x \rangle$  z= $\langle x \rangle$ 
```

Options:

<i>ambient</i>	Ambient lighting component. Valid values: integers Default value: 20 Minimum: 0 Maximum: 100
<i>diffuse</i>	Diffuse lighting component. Valid values: integers Default value: 70 Minimum: 0 Maximum: 100
<i>specular</i>	Specular lighting component. Valid values: integers Default value: 30 Minimum: 0 Maximum: 100
<i>use</i>	Use the light or not (switch it ON or OFF). Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>x</i>	The X coordinate of light position. Valid values: reals Default value: -1
<i>y</i>	The Y coordinate of light position. Valid values: reals Default value: 1
<i>z</i>	The Z coordinate of light position. Valid values: reals Default value: 3

lightambient

Sets the ambient intensity of the light.

Syntax:

```
lightambient intensity=<x>
```

Options:

<i>intensity</i>	Ambient lighting component.
Valid values:	reals
Default value:	0.2
Minimum:	0.0
Maximum:	1.0

lightdiffuse

Sets the diffuse intensity of the light.

Syntax:

```
lightdiffuse intensity=<x>
```

Options:

<i>intensity</i>	Diffuse lighting component.
Valid values:	reals
Default value:	0.7
Minimum:	0.0
Maximum:	1.0

lightposition

Sets the position of the light.

Syntax:

```
lightposition <x> <y> <z>
```

Operands:

```
<x> <y> <z>
```

The position of the light.

lightspecular

Sets the specular intensity of the light.

Syntax:

```
lightspecular intensity=<x>
```

Options:

<i>intensity</i>	Specular lighting component.
Valid values:	reals
Default value:	0.3
Minimum:	0.0
Maximum:	1.0

ligprep

This keyword is used to set various options associated with running LigPrep jobs.

Syntax:

```
ligprep desalt=yes | no epikmetalbinding=yes | no  
filter_file=<text> forcefield=mmffs | opls2005 | opls2008  
gen_conform=<n> gen_stereo=<n> gen_tautomers=yes | no  
includeoriginalstate=yes | no input_file=<text>  
ionization=generate | neutralize | retain  
ionizationmethod=ionizer | epik output_format=maestro | sdf  
ph=<x> ph_tolerance=<x> stereoisomers=retain | determine |  
generate structure_source=selected_entries | workspace | file
```

Options:

<i>desalt</i>	Desalt
Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true
<i>epikmetalbinding</i>	Add metal binding states

	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>filter_file</i>	The name of the filter criteria file.
	Valid values: text strings
	Default value:
<i>forcefield</i>	Force field to use for minimization and filtering
	Valid values: mmffs opls2005 opls2008
	Default value: opls2005
<i>gen_conform</i>	Generate low energy ring conformations
	Valid values: integers
	Default value: 1
	Minimum: 1
<i>gen_stereo</i>	Generate stereoisomers (maximum)
	Valid values: integers
	Default value: 32
	Minimum: 1
<i>gen_tautomers</i>	Generate tautomers
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>includeoriginalstate</i>	Includes original state.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>input_file</i>	The name of the structure input file.
	Valid values: text strings
	Default value:
<i>ionization</i>	How to generate the ionization states
	Valid values: generate neutralize retain
	Default value: generate
<i>ionizationmethod</i>	Which ionization computation method to use (default or Epik)
	Valid values: ionizer epik
	Default value: ionizer

Chapter 5: Commands

output_format

Whether to output the file in Maestro or SDF format.

Valid values: maestro
 sdf

Default value: **maestro**

ph

Ionization pH

Valid values: reals

Default value: **7**

Minimum: 0.0

Maximum: 14.0

ph_tolerance

Ionization pH tolerance

Valid values: reals

Default value: **2**

Minimum: 0.0

Maximum: 7.0

stereoisomers

How to generate the stereoisomers

Valid values: retain
 determine
 generate

Default value: **retain**

structure_source

Whether to use the selected entries in the current project, or what is in the workspace, or a specified file with multiple structures as structure input for the job.

Valid values: selected_entries
 workspace
 file

Default value: **file**

ligprepread

Read the given Ligprep options file and set the Ligprep panel options.

Syntax:

ligprepread \langle filename \rangle

Operands:

\langle filename \rangle

The name of the ligprep options file to read.

ligprepstart

Start a LigPrep job with the current settings.

Syntax:

ligprepstart

ligprepwrite

Write a LigPrep input file with the current settings.

Syntax:

ligprepwrite

localatoms

Specify the group of local (transformation) atoms

Syntax:

localatoms \langle ASL \rangle

Operands:

\langle ASL \rangle

For specifying which atoms are to be locally transformed. A list of atoms is created at the time the command is issued. Use the transform command to specify local scope, then use the rotate and/or translate commands to perform local transformations. To get back to global scope use the transform command to set the scope back to global.

localbitset

Specify whether local bitsets (atoms and center) are updated whenever contents of Workspace change. Unless the lock is enabled, these bitsets are re-evaluated each time atoms are changed, including changes in atom coordinates due to local adjustment.

Syntax:

```
localbitset lock=yes | no
```

Options:

lock Enable/disable lock on local bitsets for as long as ASL expressions and number of atoms in Workspace remain unchanged.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

localcenter

Set local transformation center to given centroid

Syntax:

```
localcenter
```

logfile

This is a standard alias for **scriptlogfile** (see [\[scriptlogfile\]](#), page 698).

logp

Support calculation of LogP of MacroModel.

Syntax:

```
logp enable=yes | no secsolvent=none | water | chcl3 | octanol
```

Options:

enable Determines whether logp will be used with this multiple minimization.

	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	false
<i>secsolvent</i>	This option determines the second solvent for logp.	
	Valid values:	none water chcl3 octanol
	Default value:	water

macrodefine

Define a macro. Definition is a semicolon separated string of maestro commands.

Syntax:

macrodefine \langle macro_name \rangle \langle definition \rangle

Operands:

\langle macro_name \rangle \langle definition \rangle

The first operand is the name of the macro. If this contains embedded spaces then it must be enclosed in double quotes. The second operand is the definition. If this contains embedded spaces then it must be enclosed in double quotes. To run the macro, execute `macrorun <macro_name>`

macrorun

Invokes an already defined macro

Syntax:

macrorun \langle macro_name \rangle

Operands:

\langle macro_name \rangle

The operand is the name of the macro. If this contains embedded spaces then it must be enclosed in double quotes. As a shortcut macros can also be run without preceding them with `macrorun`. Using `macrorun` has the advantage that if the macro doesn't exist the program will tell you the macro doesn't exist. If you omit `macrorun` and Maestro cannot find the macro, then

Maestro will report that the command does not exist. This will be misleading and possibly confusing because it's really the macro that will not exist.

makedirectory

Create new directory.

Syntax:

makedirectory \langle directory_name \rangle

Operands:

\langle directory_name \rangle

The name of the directory to created.

markers

The markers command creates a new marker.

Syntax:

markers *blue*= \langle x \rangle *green*= \langle x \rangle *red*= \langle x \rangle \langle marker_name \rangle

Options:

blue The blue component of the color for the new marker.

Valid values: reals

Default value: **1**

Minimum: 0.0

Maximum: 1.0

green The green component of the color for the new marker.

Valid values: reals

Default value: **1**

Minimum: 0.0

Maximum: 1.0

red The red component of the color for the new marker.

Valid values: reals

Default value: **1**

Minimum: 0.0

Maximum: 1.0

Operands:

`<marker_name>`

The name is the name to use for subsequent `showmarkername` and `hide-markername` commands.

material

Specifies a material or set material options.

Syntax:

```
material ambient=<n> changed=yes | no diffuse=<n>
          emission=<n> minimize=yes | no shininess=<n>
          specular=<n>
```

Options:

<i>ambient</i>	Ambient material component. Valid values: integers Default value: 14 Minimum: 0 Maximum: 100
<i>changed</i>	If any material option changed or not. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>diffuse</i>	Diffuse material component. Valid values: integers Default value: 100 Minimum: 0 Maximum: 100
<i>emission</i>	Emission material component. Valid values: integers Default value: 40 Minimum: 0 Maximum: 100
<i>minimize</i>	Use <code>glColorMaterial</code> to minimize performance cost or not. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>shininess</i>	Shininess material component. Valid values: integers

	Default value:	55
	Minimum:	0
	Maximum:	100
<i>specular</i>	Specular material component.	
	Valid values:	integers
	Default value:	50
	Minimum:	0
	Maximum:	100

materialdefault

A command which sets default material settings.

Syntax:

materialdefault

materialmolecular

Syntax:

materialmolecular *shininess*=⟨x⟩ *specular*=⟨x⟩

Options:

<i>shininess</i>	Shininess component for material of molecules	
	Valid values:	reals
	Default value:	0.5
	Minimum:	0.0
	Maximum:	1.0
<i>specular</i>	Specular component for material of molecules	
	Valid values:	reals
	Default value:	0.7
	Minimum:	0.0
	Maximum:	1.0

mcsd

A command which defines settings for the MC/SD (Monte Carlo/Stochastic Dynamics) simulation.

Syntax:

mcsd *maxdof*=⟨n⟩ *mindof*=⟨n⟩ *ratio*=⟨n⟩ *temperature*=⟨x⟩

Options:

<i>maxdof</i>	The maximum number of torsions which will be varied at each MC trial. Valid values: integers Default value: 1 Minimum: 1
<i>mindof</i>	The minimum number of torsions which will be varied at each MC trial. Valid values: integers Default value: 1 Minimum: 1
<i>ratio</i>	The ratio of Monte Carlo to Stochastic Dynamics steps in the simulation. Valid values: integers Default value: 1 Minimum: 1
<i>temperature</i>	The temperature at which the MCSD simulation is to be performed. [NOTE: This option is no longer used.] Valid values: reals Default value: 300 Minimum: 0.0

mini

This is a standard alias for **minienergy** (see [\[minienergy\]](#), page 343).

minienergy

Used to set values associated with a MacroModel energy minimization

Syntax:

```
minienergy converge=nothing | energy | gradient | movement  
             maxiter=⟨n⟩ method=sd | prcg | osvm | fmnr | tncg | lbfgs  
             | optimal threshold=⟨x⟩
```

Options:

<i>converge</i>	This option determines which convergence criterion will be used during an energy minimization. Valid values: nothing energy gradient movement Default value: gradient
<i>maxiter</i>	This option determines the maximum number of iterations which will be performed during an energy minimization. Valid values: integers Default value: 500 Minimum: 0 Maximum: 9999999
<i>method</i>	This option determines which minimization method will be used. Valid values: sd prcg osvm fmnr tncg lbfgs optimal Default value: prcg
<i>threshold</i>	This option determines what the convergence threshold will be. Valid values: reals Default value: 0.05 Minimum: 0.0

Aliases:

mini (see [\[mini\]](#), page 343)

minimize

Minimizes the structure in the Workspace, or just the selected atoms if there's a selection.

Syntax:

minimize

minta

Used to set up options associated with BatchMin molecular energy calculation.

Syntax:

```
minta hardlimit=⟨x⟩ numenergy=⟨n⟩ numint=⟨n⟩ softlimit=⟨n⟩
        temperature=⟨x⟩ ⟨input_file_name⟩
```

Options:

hardlimit The hard limit for sampling along normal modes.

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 3.0

numenergy

This option determines the number of energy evaluations per MINTA integration.

Valid values: integers
 Default value: **2000**
 Minimum: 1

numint

This option determines the number of MINTA iteration.

Valid values: integers
 Default value: **5**
 Minimum: 1

softlimit

The soft limit for sampling along normal modes.

Valid values: integers
 Default value: **3**
 Minimum: 1
 Maximum: 3

temperature

The temperature for MINTA calculation.

Valid values: reals
 Default value: **300**
 Minimum: 0.0

Operands:

$\langle \text{input_file_name} \rangle$

The name of the input file. This name must be given in full, including any suffix.

monitor

This is a standard alias for **energymonitor** (see [\[energymonitor\]](#), page 134).

monitorangle

Specifies a triplet of atoms to have their bond angle monitored during a MacroModel dynamics simulation.

Syntax:

monitorangle $\langle \text{atom1} \rangle \langle \text{atom2} \rangle \langle \text{atom3} \rangle$

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle \langle \text{atom3} \rangle$

The atom numbers of three atoms which are to have the angle between them monitored during a dynamics simulation. Note that specifying a-b-c is the same as c-b-a

monitordistance

Specifies a pair of atoms to have their distance monitored during a MacroModel dynamics simulation.

Syntax:

monitordistance $\langle \text{atom1} \rangle \langle \text{atom2} \rangle$

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle$

The atom numbers of two atoms which are to have the distance between them monitored during a dynamics simulation. Note that specifying a-b is the same as b-a

monitorhbond

Specifies a quartet of atoms which define an H-bond to be monitored during a dynamics simulation. The atoms are specified as X-H...Y-Z e.g. N2-H3...O2=C2 for a amide-amide bond.

Syntax:

```
monitorhbond distance=<x> hyzangle=<x> xhyangle=<x>
             <xatom> <hatom> <yatom> <zatom>
```

Options:

distance Specifies the maximum H...Y distance for an acceptable hydrogen bond.

Valid values: reals

Default value: **2.5**

Minimum: 0.0

hyzangle Specifies the minimum H...Y-Z angle for an acceptable hydrogen bond.

Valid values: reals

Default value: **90**

Minimum: 0.0

Maximum: 180.0

xhyangle Specifies the minimum X-H...Y angle for an acceptable hydrogen bond.

Valid values: reals

Default value: **120**

Minimum: 0.0

Maximum: 180.0

Operands:

```
<xatom> <hatom> <yatom> <zatom>
```

The first operand is the atom number of the heavy atom (X) to which the donor hydrogen is attached. The second operand is the atom number of the donor hydrogen itself (H). The third operand is the atom number of acceptor atom (Y) and the fourth operand is the atom number of the heavy atom (Z) attached to the acceptor.

monitorsetsurf

Specifies a set of atoms to be have their surface areas monitored during a MacroModel molecular dynamics simulation.

Syntax:

```
monitorsetsurf <ASL>
```

Operands:

<ASL>

A string in the atom specification language which describes the set of atom which are to have their surface areas monitored during a dynamics simulation.

monitorsettings

Controls settings for how the monitor panel in Maestro functions. At present it's just for whether jobs for the current project are displayed or all jobs.

Syntax:

```
monitorsettings monitorinterval=<x> monitorlines=<n>  
                showjobs=current | active | all
```

Options:

monitorinterval

This option determines the minimum interval, in seconds, between monitoring updates while monitoring a job. The actual interval may be longer than the specified value, due to time spent in panel and Workspace updates.

Valid values:	reals
Default value:	1
Minimum:	1.0
Maximum:	300.0

monitorlines

This option determines the maximum number of lines to be displayed in the monitor panel File tab in each monitoring interval while monitoring a job.

Valid values:	integers
Default value:	100
Minimum:	100

showjobs This determines whether the jobs listed in the monitor panel are limited to those for the current project or all known jobs for the current user.

Valid values: current
 active
 all

Default value: **current**

monitorsurf

Specifies a single atom to have its surface area monitored during a MacroModel molecular dynamics simulation.

Syntax:

monitorsurf $\langle \text{atom_number} \rangle$

Operands:

$\langle \text{atom_number} \rangle$

The number of an atom which is to have its surface areas monitored during a dynamics simulation.

monitortorsion

Specifies a quartet of atoms to have their torsion angle monitored during a MacroModel dynamics simulation.

Syntax:

monitortorsion $\langle \text{atom1} \rangle \langle \text{atom2} \rangle \langle \text{atom3} \rangle \langle \text{atom4} \rangle$

Operands:

$\langle \text{atom1} \rangle \langle \text{atom2} \rangle \langle \text{atom3} \rangle \langle \text{atom4} \rangle$

The atom numbers of four atoms which are to have the torsion between them monitored during a dynamics simulation. Note that specifying a-b-c-d is the same as d-c-b-a

mouse

Specify how mouse movements are to be interpreted

Syntax:

mouse *rotate*=xy | x | y | z | adjust *translate*=xy | x | y | z

There are no operands for this command.

Options:

rotate How mouse movements are interpreted for rotation. Default is rotate in xy.

Valid values: xy
 x
 y
 z
 adjust

Default value: **xy**

translate How mouse movements are interpreted for translation. Default is translate in xy.

Valid values: xy
 x
 y
 z

Default value: **xy**

Operands:

There are no operands for this command.

move

Moves the specified atoms by the offsets specified in x, y and z values

Syntax:

move <atom_num> <xinc> <yinc> <zinc>

Operands:

<atom_num> <xinc> <yinc> <zinc>

The first operand is the number of the atom which is to be moved. The following three real numbers represent the offsets to be applied to the x, y and z coordinates in Angstroms.

multiplemini

Used to set up a multiple minimization. The operand is the name of a file which contains a number of structures.

Syntax:

```
multiplemini distinguishenantiomers=yes | no
              eliminate=atom_deviation | rmsd max_rmsd=<x>
              maxdist=<x> savemaximum=<n> setsavemaximum=yes | no
              window=<x> <input_file_name>
```

Options:

distinguishenantiomers

A boolean for whether to retain different enantiomers. If true, an additional opcode, NANT, is written to the .com file.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

eliminate

The method to use for eliminating redundant conformers: maximum atom deviation or RMSD.

Valid values: atom_deviation
rmsd

Default value: **atom_deviation**

max_rmsd

Maximum RMSD for considering two structures equal.

Valid values: reals

Default value: **0.5**

Minimum: 0.0

maxdist

Maximum distance between atoms in equal structures.

Valid values: reals

Default value: **0.5**

Minimum: 0.0

savemaximum

Maximum number of structures to save

Valid values: integers

Default value: **100**

Minimum: 1

setsavemaximum

A boolean for whether to set a maximum number of structures to save.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Chapter 5: Commands

window Specifies the energy window (in kJ/mol) for saving structures during a multiple minimization.

Valid values: reals

Default value: **21**

Minimum: 0.0

Operands:

`<input_file_name>`

The name of the input file. This name will contain all the structures which are to be minimized. This name must be given in full, including any suffix.

mutate

Replace the sidechain of any residue in set `<ASL>` with the sidechain of the currently selected fragment

Syntax:

mutate `<ASL>`

Operands:

`<ASL>`

The set of atoms which are to be mutated. Each residue represented in this set will be mutated.

nextresidue

If only a single residue is selected in the Workspace, then find the next one.

Syntax:

nextresidue

optimizefog

Fit clipping planes to displayed structure to maximize the shading contrast between near and far atoms.

Syntax:

optimizefog

partialcharge

Set the partial atomic charges for all atoms which match the ASL specification.

Syntax:

partialcharge $\langle \text{charge1} \rangle \langle \text{charge2} \rangle \langle \text{ASL} \rangle$

Operands:

$\langle \text{charge1} \rangle \langle \text{charge2} \rangle \langle \text{ASL} \rangle$

The first operand must be a valid real number which will be used for the charge1 field of all matching atoms. The second operand must be a valid real number which will be used for the charge2 field of all matching atoms. The final operand is a valid ASL string which specifies which atoms are to have their charges changed.

pause

Pauses immediately for the length of time specified (in seconds, unless other units are specified). Setting a pause of zero (e.g. “0” or “0.00”) will mean there is no pause. If no operand is specified, or a negative time is given, then Maestro will remain paused indefinitely, until it is interrupted. All pauses can be interrupted by pressing any key on the keyboard or any mouse button.

Syntax:

pause $\langle \text{pause_time} \rangle [\text{sec} | \text{min} | \text{hour} | \text{day}]$

Operands:

$\langle \text{pause_time} \rangle [\text{sec} | \text{min} | \text{hour} | \text{day}]$

The length of time to pause - floating point value in seconds, accurate to within about 0.05 second.

pausecommands

Sets the time to pause after execution of every command in script (in seconds, unless other units are specified). Setting a pause of zero (e.g. “0” or “0.00”) will mean there is no pause after each command. If no operand is specified, or a negative time is given, then the pause after each command will be indefinite, pausing until interrupted. All pauses can be interrupted by pressing any key on the keyboard or any mouse button.

Syntax:

pausecommands \langle pause_time \rangle [sec | min | hour | day]

Operands:

\langle pause_time \rangle [sec | min | hour | day]

The time to pause after each command - floating point value in seconds, accurate to within about 0.05 second.

peptidecistrans

Perform a cis/trans interconversion of the peptide in the residue which contains the specified atom.

Syntax:

peptidecistrans \langle atom \rangle

Operands:

\langle atom \rangle

The atom which defines the residue to have the cis/trans inter-conversion done.

peptiderotate

Rotate the plane of the peptide in the residue which contains the specified atom

Syntax:

peptiderotate *carbonyl_only*=yes | no *increment*= $\langle x \rangle$ $\langle atom \rangle$

Options:

carbonyl_only

When this is true then only the carbonyl of the peptide is rotated.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

increment Value by which to rotate the plane of the peptide from its current value.

Valid values: reals

Default value: **2**

Minimum: -180.1

Maximum: 180.1

Operands:

$\langle atom \rangle$

The atom which defines the residue to have a peptide rotated.

perresiduedisplay

Specify whether or not display per-residue interactions.

Syntax:

perresiduedisplay *display*=yes | no

Options:

display If this option is true, display per-residue interactions.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

perresiduesubset

Set subset of per-residue interactions by the ASL operand.

Syntax:

perresiduesubset $\langle \text{ASL} \rangle$

Operands:

$\langle \text{ASL} \rangle$

A string in the atom specification language which describes the subset of per-residue interactions.

perresiduetype

Specify the type of displaying per-residue interactions.

Syntax:

perresiduetype *type*=eint | vdw | coulomb | hbond | distance

Options:

type Specify the type of displaying per-residue interactions, Eint, vdw, Coulomb, H-Bond energy, or Minimum distance.

Valid values: eint
 vdw
 coulomb
 hbond
 distance

Default value: **eint**

phaseaddconstraint

Adds the given constraint.

Syntax:

phaseaddconstraint *value*= $\langle x \rangle$ $\langle \text{constraint} \rangle$

Options:

value The value to set for the given constraint.

Valid values: reals

Default value: **1**

Operands:

⟨constraint⟩

The constraint to add. Should be in the form: <site 1>-<site 2> and so on, up to four sites.

phaseaddcustomfeature

Adds the given custom feature.

Syntax:

phaseaddcustomfeature *code*=⟨text⟩ ⟨name⟩

Options:

code The name of the custom feature to add.

Valid values: text strings

Default value: **A**

Operands:

⟨name⟩

The name for the feature.

phaseaddhypotheses

Copies the selected hypotheses from Develop Common Pharmacophore panel (in either the Score Hypothesis or Build QSAR model steps) to the Project Table, and opens the Project Table.

Syntax:

phaseaddhypotheses

phaseaddligands

Adds the given ligands to the table.

Syntax:

```
phaseaddligands activity=⟨text⟩ convert_activity=yes | no  
                convert_scale=⟨x⟩ ⟨ESL⟩
```

Options:

activity This determines which property (if any) to use as the activity property for the ligands.

Valid values: text strings

Default value:

convert_activity

Set to true if the activity values should be converted from concentration to -log[concentration]

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

convert_scale

A scale factor for conversion.

Valid values: reals

Default value: **1**

Operands:

⟨ESL⟩

The entries to add as ligands.

phaseaddligandsfromfile

Adds the ligands in the given files to the current Phase run.

Syntax:

```
phaseaddligandsfromfile activity=⟨text⟩  
                convert_activity=yes | no convert_scale=⟨x⟩ ⟨file 1⟩ ⟨file 2⟩
```

Options:

activity Which property (if any) to use as the activity property for the ligands.

Valid values: text strings

Default value:

convert_activity

Set to true if the activity values should be converted from concentration to $-\log[\text{concentration}]$

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

convert_scale

A scale factor for conversion

Valid values: reals

Default value: **1**

Operands:

$\langle \text{file 1} \rangle \langle \text{file 2} \rangle$

The files to add ligands from.

phaseaddligandsfromrun

Adds the given ligands to the table.

Syntax:

phaseaddligandsfromrun *runname*= $\langle \text{text} \rangle \langle \text{ligand_names} \rangle$

Options:

runname This determines which run to use as the reference to copy the corresponding properties of the ligands.

Valid values: text strings

Default value:

Operands:

$\langle \text{ligand_names} \rangle$

The ligand names to be added separated by semi-colon.

phaseaddsite

Adds a site to the current freestyle hypothesis

Syntax:

phaseaddsite $x=\langle x \rangle$ $y=\langle x \rangle$ $z=\langle x \rangle$ $\langle \text{site type} \rangle$

Options:

x	X-coordinate of the new site. Valid values: reals Default value: 0
y	X-coordinate of the new site. Valid values: reals Default value: 0
z	X-coordinate of the new site. Valid values: reals Default value: 0

Operands:

$\langle \text{site type} \rangle$

The single-letter feature type.

phasealignhypotheses

Sets whether or not to view non-model ligands and whether or not to have a site mask.

Syntax:

phasealignhypotheses $\text{alignment_type}=\text{sites} \mid \text{rmsd}$
 $\text{must_match}=\langle n \rangle$ $\text{tolerance}=\langle x \rangle$ $\langle \text{hypothesis} \rangle$

Options:

alignment_type

Indicates how to choose the best alignment

Valid values: sites
rmsd

Default value: **sites**

must_match

The minimum number of equivalent sites that have to match to do the alignment.

Valid values: integers
Default value: **3**
Minimum: 3

tolerance The tolerance in angstroms.

Valid values: reals

Default value: **2**

Minimum: 0.0

Operands:

⟨hypothesis⟩

The name of the hypothesis to set the options for.

phasebuildqsar

Sets the QSAR options and launches a Build QSAR job.

Syntax:

```
phasebuildqsar color_ligands=yes | no eliminatetvalue=yes | no
gridspacing=⟨x⟩ keep_training_set=yes | no
maxplsfactor=⟨n⟩ modeltype=atom_based | phase_based
sample_uniformly=yes | no tolerance_a=⟨x⟩ tolerance_d=⟨x⟩
tolerance_h=⟨x⟩ tolerance_n=⟨x⟩ tolerance_p=⟨x⟩
tolerance_r=⟨x⟩ tolerance_x=⟨x⟩ tolerance_y=⟨x⟩
tolerance_z=⟨x⟩ tvalue=⟨x⟩
```

Options:

color_ligands

Indicates whether or not to color ligands to match the coloring in the ligands table.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

eliminatetvalue

Indicates whether to use t-value filter before launching QSAR job.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

gridspacing

The grid spacing to be for the QSAR model grid. The valid range is 0.5 to 2.0 angstrom.

Valid values: reals

Default value: **1**

Minimum: 0.5

Maximum: 2.0

Chapter 5: Commands

keep_training_set

Indicates whether or not to keep the actives and inactives in the training set.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

maxplsfactor

This option value range depends on the ligands in the training set. The valid range is 1 to N/5, where N is the number of ligands in the training set. If N is less than 5, the maximum number of factors is 1, but no QSAR model can be constructed.

Valid values: integers

Default value: **1**

Minimum: 1

modeltype The value of the model type in the Build QSAR Model step options dialog.

Valid values: atom_based
 phase_based

Default value: **atom_based**

sample_uniformly

Indicates whether or not to sample uniformly over the activity coordinate when assigning to test and training sets.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

tolerance_a

The feature matching tolerance for the hydrogen bond acceptor feature. The valid range is 0.0 to 100.0

Valid values: reals

Default value: **1**

Minimum: 0.0

Maximum: 100.0

tolerance_d

The feature matching tolerance for the hydrogen bond donor feature. The valid range is 0.0 to 100.0

Valid values: reals

Default value: **1**

Minimum: 0.0

Maximum: 100.0

tolerance_h

The feature matching tolerance for the hydrophobic feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1.5**
 Minimum: 0.0
 Maximum: 100.0

tolerance_n

The feature matching tolerance for the negative feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **0.75**
 Minimum: 0.0
 Maximum: 100.0

tolerance_p

The feature matching tolerance for the positive feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **0.75**
 Minimum: 0.0
 Maximum: 100.0

tolerance_r

The feature matching tolerance for the aromatic ring feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1.5**
 Minimum: 0.0
 Maximum: 100.0

tolerance_x

The feature matching tolerance for the custom(X) feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 100.0

tolerance_y

The feature matching tolerance for the custom(Y) feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 100.0

tolerance_z

The feature matching tolerance for the custom(Z) feature. The valid range is 0.0 to 100.0

Valid values:	reals
Default value:	1
Minimum:	0.0
Maximum:	100.0

tvalue This option value indicates the values to be used to filter t-value.

Valid values:	reals
Default value:	2
Minimum:	0.01

phasechangesitetype

Changes the type of the selected sites

Syntax:

phasechangesitetype \langle site type \rangle

Operands:

\langle site type \rangle

The new site type.

phasechooseactiveset

Sets the Active flag based on whether the activity value is larger than the cutoff value.

Syntax:

phasechooseactiveset \langle cutoff \rangle

Operands:

\langle cutoff \rangle

The cutoff value for deciding between active and inactive.

phasecleanupstructures

Launches a Cleanup Structures job for Phase.

Syntax:

phasecleanupstructures

phaseclearcentroidatoms

Removes all atoms from the centroid list for excluded volumes.

Syntax:

phaseclearcentroidatoms

phaseclusterhypotheses

Sets the Cluster Hypotheses options and launches the job.

Syntax:

phaseclusterhypotheses *linkage*=complete | average | single

Options:

<i>linkage</i>	The linkage type.	
	Valid values:	complete average single
	Default value:	complete

phaseconfgen

Defines settings for Phase Generate Conformers job.

Syntax:

```
phaseconfgen amidebonds=vary | retain | trans
              eliminate=atom_deviation | rmsd field=mmffs | opls2005
              incorporate=append | replace | ignore | appendungrouped |
              workspace max_rmsd=<x> maxdist=<x> method=default |
              mixed minimizationsteps=<n> numrotatablesteps=<n>
              numsteps=<n> postmaxiter=<n> postprocessing=yes | no
              postprocessingmethod=mini | filter | rce premaxiter=<n>
              preprocessing=yes | no sampling=standard | rapid | complete
              | thorough solvation=gbsa | distance_dependent window=<x>
```

Options:

amidebonds

This determines whether to vary amide bond conformation, retain original amide bond conformation, or set amide bond conformation to trans.

Valid values: vary
 retain
 trans

Default value: **vary**

eliminate

The method to use for eliminating redundant conformers: maximum atom deviation or RMSD.

Valid values: atom_deviation
 rmsd

Default value: **rmsd**

field

This determines which force field mmffs|mmff|opls2001 is used. Currently we always use mmffs, so it will have only one option value.

Valid values: mmffs
 opls2005

Default value: **opls2005**

incorporate

This option controls the incorporation of the results (replace or append).

Valid values: append
 replace
 ignore
 appendungrouped
 workspace

Default value: **replace**

max_rmsd

Maximum RMSD for considering two structures equal.

Valid values: reals

	Default value: 1
	Minimum: 0.0
<i>maxdist</i>	Maximum distance between atoms in equal structures.
	Valid values: reals
	Default value: 2
	Minimum: 0.0
<i>method</i>	This determines whether MacroModel uses the ligand torsion search method (default) or the mixed MCMM/LMOD search method (mixed) to generate conformers.
	Valid values: default mixed
	Default value: default
<i>minimizationsteps</i>	This option determines the maximum number of minimization steps for Mixed MCMM/LMOD generation.
	Valid values: integers
	Default value: 100
	Minimum: 1
<i>numrotatablesteps</i>	An option which sets the number of steps which will be performed during the ConfGen conformational search.
	Valid values: integers
	Default value: 100
	Minimum: 1
<i>numsteps</i>	An option which sets the number of steps which will be performed during the conformational search. This also limits number of conformations generated.
	Valid values: integers
	Default value: 1000
	Minimum: 0
<i>postmaxiter</i>	This option determines the maximum number of iterations for post-minimization of generated structures.
	Valid values: integers
	Default value: 50
	Minimum: 0
	Maximum: 9999999
<i>postprocessing</i>	Indicates whether or not to perform MacroModel postprocessing.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

postprocessingmethod

This determines which type of postprocessing method to use (minimization, filtering and redundant conformer elimination, or redundant conformer elimination only).

Valid values: mini
 filter
 rce
 Default value: **mini**

premaxiter

This option determines the maximum number of iterations for pre-minimization of input structures.

Valid values: integers
 Default value: **100**
 Minimum: 0
 Maximum: 9999999

preprocessing

Indicates whether or not to perform MacroModel preprocessing.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

sampling

This determines whether rapid (standard) or thorough (complete) sampling will be used.

Valid values: standard
 rapid
 complete
 thorough
 Default value: **standard**

solvation

This determines whether GB/SA Water (gbsa) or Distance Dependent Dielectric (distance_dependent) solvation treatment is used.

Valid values: gbsa
 distance_dependent
 Default value: **distance_dependent**

window

The energy window (in kcal/mol) within which structures will be saved.

Valid values: reals
 Default value: **10**
 Minimum: 0.0

phasecopyhypothesisfrombuilder

Copy the given hypothesis from the hypotheses table in the Edit Hypothesis panel.

Syntax:

phasecopyhypothesisfrombuilder \langle row number \rangle

Operands:

\langle row number \rangle

The row number of the hypothesis to copy.

phasecreateevactives

Creates excluded volumes from actives and inactives.

Syntax:

phasecreateevactives

phasecreateevactivessettings

Settings for creating excluded volumes from actives and inactives

Syntax:

phasecreateevactivessettings *actives_file*= \langle text \rangle
distance= \langle x \rangle *inactives_file*= \langle text \rangle *min_inactives*= \langle n \rangle
radii= \langle x \rangle

Options:

actives_file

The path to the actives file.

Valid values: text strings

Default value:

distance

The minimum distance between the active surface and the excluded volumes.

Valid values: reals

Default value: **1**

Minimum: 0.0

Chapter 5: Commands

inactives_file

The path to the inactives file.

Valid values: text strings

Default value:

min_inactives

The minimum number of inactives that must experience a clash.

Valid values: integers

Default value: **1**

Minimum: 1

radii

The radii for the excluded volumes.

Valid values: reals

Default value: **1**

Minimum: 0.0

phasecreateevreceptor

Creates excluded volumes from receptor atoms.

Syntax:

phasecreateevreceptor

phasecreateevreceptoratoms

Specify the atoms to use to define an excluded volume shell.

Syntax:

phasecreateevreceptoratoms \langle ASL \rangle

Operands:

\langle ASL \rangle

A string in the atom specification language. Defines a set of atoms to use for building an excluded volume shell.

phasecreateevreceptorsettings

Settings for creating excluded volumes from receptor atoms.

Syntax:

```
phasecreateevreceptorsettings fixed_radius=<x>
    fixed_scaling=<x> ignore_atoms=yes | no ignore_distance=<x>
    limit_thickness=yes | no radii_scaling=fixed | property
    radii_sizes=vdw | fixed | property radius_property=<text>
    scaling_property=<text> thickness=<x> <ASL>
```

Options:

fixed_radius

The fixed radius to use for excluded volumes.

Valid values: reals

Default value: **1**

Minimum: 0.0

fixed_scaling

The fixed scaling to apply to excluded volume radii.

Valid values: reals

Default value: **1**

Minimum: 0.0

ignore_atoms

Ignore receptor atoms which are too close to the ligand surface.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

ignore_distance

The distance to consider as too close to the ligand atoms.

Valid values: reals

Default value: **1**

Minimum: 0.0

limit_thickness

Limit excluded volume shell thickness to a specified number.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

radii_scaling

Indicates how to scale the radii of the excluded volumes.

Valid values: fixed
property

Default value: **fixed**

radii_sizes

Indicates how to determine the radii of the excluded volumes.

Valid values: vdw
fixed
property

Default value: **vdw**

radius_property

The atom property to use as the radius for excluded volumes.

Valid values: text strings

Default value:

scaling_property

The atom property to use to scale the excluded volume radii by.

Valid values: text strings

Default value:

thickness The thickness to limit the excluded volume shell to.

Valid values: reals

Default value: **5**

Minimum: 0.0

Operands:

⟨ ASL ⟩

An ASL to use as the receptor atoms.

phasecreateevreference

Creates excluded volumes from reference structures.

Syntax:

phasecreateevreference

phasecreateevreferenceatoms

Specify the atoms to use to define an excluded volume shell.

Syntax:

phasecreateevreferenceatoms ⟨ ASL ⟩

Operands:

⟨ ASL ⟩

A string in the atom specification language. Defines a set of atoms to use for building an excluded volume shell.

phasecreateevreferencesettings

Settings for creating excluded volumes from reference structure atoms.

Syntax:

```
phasecreateevreferencesettings distance=<x> radii=<x>
                             source=ligand | atoms <ASL>
```

Options:

<i>distance</i>	The minimum distance between the ligand surface and the excluded volume shell.
	Valid values: reals
	Default value: 1
	Minimum: 0.0
<i>radii</i>	The excluded volume sphere radii.
	Valid values: reals
	Default value: 1
	Minimum: 0.0
<i>source</i>	Indicates the source of the reference structures.
	Valid values: ligand atoms
	Default value: ligand

Operands:

<ASL>

An ASL to use as the reference atoms.

phasecreateexcludedvolume

Creates an excluded volume from the current centroid atoms.

Syntax:

```
phasecreateexcludedvolume
```

phasecreatefreestylehypothesis

Sets the given entry to have the current freestyle hypothesis.

Syntax:

phasecreatefreestylehypothesis <entry id>

Operands:

<entry id>

The entry to add the freestyle hypothesis to.

phasedbaddconfset

Adds the selected confset to the Phase DB

Syntax:

phasedbaddconfset

phasedbaddligands

Adds ligands from the specified file to the database. The ligands will go through a cleanup process.

Syntax:

phasedbaddligands *alreadycleaned*=yes | no
 generatetautomers=yes | no *ionization*=retain | neutralize |
 generate *maxisomers*=<n> *maxtautomers*=<n>
 stereoisomers=retain | determine | all *target-ph*=<x>
 <file name>

Options:

alreadycleaned

Indicates whether or not the structures have already been cleaned.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

generatetautomers

Indicates whether to generate tautomers or not

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

ionization Sets the handling of ionization.

Valid values: retain
 neutralize
 generate
 Default value: **retain**

maxisomers

The maximum number of stereoisomers to return for any ligand.

Valid values: integers
 Default value: **10**
 Minimum: 1

maxtautomers

The maximum number of stereoisomers to return for any ligand.

Valid values: integers
 Default value: **8**
 Minimum: 1

stereoisomers

Sets the handling of chiralities.

Valid values: retain
 determine
 all
 Default value: **all**

target_ph

The target pH for generating ionization states.

Valid values: reals
 Default value: **7**

Operands:

⟨ file name ⟩

The name of the file of structure to be added. These will go through the cleanup process with the options specified below. The ligands can be in Maestro or SD format.

phasedbaddligandsfromdb

Adds ligands from the specified database

Syntax:

phasedbaddligandsfromdb

phasedbconfgen

Defines settings for PhaseDB Generate Conformers job.

Syntax:

```
phasedbconfgen amidebonds=vary | retain | trans
                eliminate=atom_deviation | rmsd field=mmffs | opls2005
                incorporate=append | replace | ignore | appendungrouped |
workspace max_rmsd=<x> maxdist=<x> method=default |
mixed minimizationsteps=<n> numrotatablesteps=<n>
numsteps=<n> postmaxiter=<n> postprocessing=yes | no
postprocessingmethod=mini | filter | rce premaxiter=<n>
preprocessing=yes | no sampling=standard | rapid | complete
| thorough solvation=gbsa | distance_dependent window=<x>
```

Options:

amidebonds

This determines whether to vary amide bond conformation, retain original amide bond conformation, or set amide bond conformation to trans.

Valid values: vary
 retain
 trans

Default value: **vary**

eliminate

The method to use for eliminating redundant conformers: maximum atom deviation or RMSD.

Valid values: atom_deviation
 rmsd

Default value: **rmsd**

field

This determines which force field mmffs|mmff|opls2001 is used. Currently we always use mmffs, so it will have only one option value.

Valid values: mmffs
 opls2005

Default value: **opls2005**

incorporate

This option controls the incorporation of the results (replace or append).

	Valid values:	append replace ignore appendungrouped workspace
	Default value:	replace
<i>max_rmsd</i>	Maximum RMSD for considering two structures equal.	
	Valid values:	reals
	Default value:	1
	Minimum:	0.0
<i>maxdist</i>	Maximum distance between atoms in equal structures.	
	Valid values:	reals
	Default value:	2
	Minimum:	0.0
<i>method</i>	This determines whether MacroModel uses the ligand torsion search method (default) or the mixed MCMM/LMOD search method (mixed) to generate conformers. Currently database creation always uses the default method, so it will have only one option value.	
	Valid values:	default mixed
	Default value:	default
<i>minimizationsteps</i>	This option determines the maximum number of iterations for minimization of generated structures in Mixed MCMM/LMOD.	
	Valid values:	integers
	Default value:	100
	Minimum:	1
<i>numrotatablesteps</i>	An option which sets the number of steps which will be performed during the ConfGen conformational search.	
	Valid values:	integers
	Default value:	10
	Minimum:	1
<i>numsteps</i>	An option which sets the number of steps which will be performed during the conformational search. This also limits number of conformations generated.	
	Valid values:	integers
	Default value:	100
	Minimum:	0

Chapter 5: Commands

postmaxiter

This option determines the maximum number of iterations for post-minimization of generated structures.

Valid values: integers

Default value: **50**

Minimum: 0

Maximum: 9999999

postprocessing

Indicates whether or not to perform MacroModel postprocessing.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

postprocessingmethod

This determines which type of postprocessing method to use (minimization, filtering and redundant conformer elimination, or redundant conformer elimination only).

Valid values: mini

filter

rce

Default value: **rce**

premaxiter

This option determines the maximum number of iterations for pre-minimization of input structures.

Valid values: integers

Default value: **100**

Minimum: 0

Maximum: 9999999

preprocessing

Indicates whether or not to perform MacroModel preprocessing.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

sampling

This determines whether rapid (standard) or thorough (complete) sampling will be used.

Valid values: standard

rapid

complete

thorough

Default value: **standard**

solvation

This determines whether GB/SA Water (gbsa) or Distance Dependent Dielectric (distance_dependent) solvation treatment is used.

	Valid values:	gbsa distance_dependent
	Default value:	distance_dependent
<i>window</i>	The energy window (in kcal/mol) within which structures will be saved.	
	Valid values:	reals
	Default value:	10
	Minimum:	0.0

phasedbcreatesites

Launches a Create Sites job for a Phase database.

Syntax:

phasedbcreatesites *withligands*=all | withsites | selected

Options:

withligands

Controls which ligands have sites regenerated.

Valid values: all
withsites
selected

Default value: **withsites**

phasedbcreatesubsetfromhitfile

Create a new subset based on a Maestro structure file which is the result of a phase DB search.

Syntax:

phasedbcreatesubsetfromhitfile <subsetname> <hit file>

Operands:

<subsetname> <hit file>

The name of the new subset and a Maestro structure file which is the result of a Phase database search.

phasedbcreatesubsetfromselected

Create a new subset from those ligands which are currently selected in the ligand table.

Syntax:

phasedbcreatesubsetfromselected `<subsetname>`

Operands:

`<subsetname>`

The name of the new subset.

phasedbcreatesubsetfromtextfile

Create a new subset based on a file containing ligand titles - one per line.

Syntax:

phasedbcreatesubsetfromtextfile `<subset name>` `<text file>`

Operands:

`<subset name>` `<text file>`

The name of the subset and a file which contains a list of ligand names - one per line.

phasedbdeleteligands

Launches a Delete Ligands job for selected ligands in a Phase DB.

Syntax:

phasedbdeleteligands *whichstructures*=selected | all

Options:

whichstructures

Determines which structures this operation will be applied to.

Valid values: selected
 all

Default value: **selected**

phasedbdeletesubset

Delete the named subset from the DB

Syntax:

phasedbdeletesubset \langle subset name \rangle

Operands:

\langle subset name \rangle

The name of the subset to be deleted.

phasedbdisplayrange

Control the range of ligands which are displayed in the table.

Syntax:

phasedbdisplayrange *end*= \langle n \rangle *start*= \langle n \rangle

Options:

<i>end</i>	Sets the end point or the display range
	Valid values: integers
	Default value: 1000
	Minimum: 1
<i>start</i>	Sets the starting point or the display range
	Valid values: integers
	Default value: 1
	Minimum: 1

phasedbexportligands

Export the ligands corresponding to selected ligands to the project table.

Syntax:

phasedbexportligands *confs*=firstconf | allconfs

Options:

<i>confs</i>	Controls what is exported - all conformations or just the first of each ligand.
Valid values:	firstconf allconfs
Default value:	firstconf

phasedbfilter

Replace currently filtered structures with ones which match the expression

Syntax:

phasedbfilter \langle filter_expr \rangle

Operands:

\langle filter_expr \rangle

The filter string which is to be applied to the ligand names to determine which ligands are to be displayed in the database table.

phasedbfilteradd

Add to the currently filtered structures with ones which match the expression.

Syntax:

phasedbfilteradd \langle filter_expr \rangle

Operands:

\langle filter_expr \rangle

The filter string which is to be applied to the ligand names to determine which ligands are to be added to the database table.

phasedbgenerateconformers

Launches a Generate Conformers job for a Phase database.

Syntax:

phasedbgenerateconformers *whichstructures*=selected | all

Options:

whichstructures

Determines which structures this operation will be applied to.

Valid values: selected
 all

Default value: **selected**

phasedbnewfromfile

Create a new Phase 3D database and open it in Maestro.

Syntax:

phasedbnewfromfile <db_dir_path> <name>

Operands:

<db_dir_path> <name>

The path (location) of the directory to be created as the new Phase 3D database directory and the name of the DB

phasedbopen

Open an existing Phase 3D database into Maestro.

Syntax:

phasedbopen <db_dir_path> <db_name>

Operands:

<db_dir_path> <db_name>

The path (location) and name of the Phase 3D database directory to be opened.

phasedbremoveconformers

Launches a Remove Conformers job for selected ligands in a Phase DB.

Syntax:

phasedbremoveconformers *whichstructures*=selected | all

Options:

whichstructures

Determines which structures this operation will be applied to.

Valid values: selected
 all

Default value: **selected**

phasedbselectextendablerow

Extend the selection in the Phase Database table from the selected table row to joining up with an existing selection.

Syntax:

phasedbselectextendablerow <row_number>

Operands:

<row_number>

The row number in the table from which the selection is to begin.

phasedbselectonlyablerow

Select a row from the phase database table.

Syntax:

phasedbselectonlyablerow <row_number>

Operands:

<row_number>

The row number in the table which is to be selected.

phasedbselectsuubsetrow

Select the row in the subset table.

Syntax:

phasedbselectsuubsetrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number of the subset to be selected.

phasedbselecttablerow

Selects a row from the phase database table.

Syntax:

phasedbselecttablerow $\langle \text{row_number} \rangle$

Operands:

$\langle \text{row_number} \rangle$

The row number in the table which is to be selected.

phasedbunselecttablerow

Unselects a row from the Phase Database table.

Syntax:

phasedbunselecttablerow $\langle \text{row_number} \rangle$

Operands:

$\langle \text{row_number} \rangle$

The row number in the table which is to be unselected.

phasedeletcustomfeature

Deletes the given custom feature.

Syntax:

phasedeletecustomfeature $\langle \text{code} \rangle$

Operands:

$\langle \text{code} \rangle$

The code for the feature to delete.

phasedeleteexcludedvolumes

Deletes the selected excluded volumes.

Syntax:

phasedeleteexcludedvolumes

phasedeletehypothesis

Deletes the given hypothesis from the hypotheses table in the Score Hypotheses or Build QSAR Model step.

Syntax:

phasedeletehypothesis $\langle \text{ID} \rangle$

Operands:

$\langle \text{ID} \rangle$

The hypothesis to delete.

phasedeletehypothesisfrombuilder

Deletes the given hypothesis from the hypotheses table in the Edit Hypothesis and Find Matches panels.

Syntax:

phasedeletehypothesisfrombuilder \langle row number \rangle

Operands:

\langle row number \rangle

The row number of the hypothesis to delete.

phasedeleteselectedconstraints

Deletes all of the selected constraints.

Syntax:

phasedeleteselectedconstraints

phasedeleteselectedligands

Deletes the selected ligands in the current Phase step.

Syntax:

phasedeleteselectedligands

phasedeletesites

Deletes the selected sites from the current freestyle hypothesis.

Syntax:

phasedeletesites

phasedisplayproperty

Displays the given property from the table.

Syntax:

phasedisplayproperty \langle property name \rangle

Operands:

\langle property name \rangle

The property to display.

phaseevoptions

Options for Phase excluded volumes

Syntax:

phaseevoptions *radius*= \langle x \rangle

Options:

radius Radius for creating new excluded volumes.

Valid values: reals

Default value: **1**

Minimum: 0.0

phaseexcludetablerow

Excludes the given row in the first table in the step from the Workspace.

Syntax:

phaseexcludetablerow \langle row \rangle

Operands:

\langle row \rangle

The row number to exclude in the Workspace.

phaseexportalignmentstofile

Exports the selected alignments from the ligands table in the Score Hypotheses or Build QSAR steps.

Syntax:

phaseexportalignmentstofile <file name>

Operands:

<file name>

The file name to export the selected alignments to.

phaseexportconformerstofile

Exports the selected conformers from the ligands table in the Prepare Ligands or Create Sites steps.

Syntax:

phaseexportconformerstofile <file name>

Operands:

<file name>

The file name to export the selected conformers to.

phaseexportfeature

Exports the features to the given file name.

Syntax:

phaseexportfeature <file name>

Operands:

<file name>

A path to a write the features file to.

phaseexporthypothesis

Exports the included hypothesis from the hypotheses table in the Score Hypotheses or Build QSAR Model step.

Syntax:

phaseexporthypothesis `<file name>`

Operands:

`<file name>`

The file name to export the included hypothesis to.

phaseexporthypothesisfrombuilder

Exports the selected hypothesis from the hypotheses table in the Edit Hypotheses panel.

Syntax:

phaseexporthypothesisfrombuilder `<file name>`

Operands:

`<file name>`

The file name to export the selected hypothesis to.

phaseexportselectedalignments

Exports the selected alignments from the current Phase step to the project table.

Syntax:

phaseexportselectedalignments

phaseexportselectedhypotheses

Exports all selected hypotheses from the hypotheses table in the Score Hypotheses or Build QSAR Model step.

Syntax:

phaseexportselectedhypotheses <directory name>

Operands:

<directory name>

The directory name to export all selected hypothesis to.

phaseexportselectedhypothesesfrombuilder

Exports the selected hypotheses from the hypotheses table in the Edit Hypotheses panel.

Syntax:

phaseexportselectedhypothesesfrombuilder <directory name>

Operands:

<directory name>

The directory to export the selected hypotheses to.

phaseexportselectedligands

Exports the selected ligands from the current Phase step to the project table.

Syntax:

phaseexportselectedligands

phaseextendsite

Extends the selection to the given site in the freestyle hypothesis

Syntax:

phaseextendsite <site index>

Operands:

<site index>

The index of the site to select.

phasefieldqsar

Command for controlling Phase Field-Based QSAR Models

Syntax:

```
phasefieldqsar contourstyle=solid | mesh | dots  
               crossvalidation=⟨ n ⟩ electrostaticcutoff=⟨ x ⟩  
               eliminatetvalue=yes | no extendgrid=⟨ x ⟩  
               featureddefinitions=⟨ text ⟩ fieldstyle=ff | gauss | custom  
               grid_spacing=⟨ x ⟩ ignoredistance=⟨ x ⟩ ignorefields=yes | no  
               max_factors=⟨ n ⟩ random=⟨ n ⟩ random_seed=⟨ n ⟩  
               stddevcutoff=⟨ x ⟩ stericcutoff=⟨ x ⟩  
               truncateelectrostatic=yes | no truncatesteric=yes | no  
               tvaluecutoff=⟨ x ⟩ useinputpartialcharges=yes | no
```

Options:

contourstyle

What style to use for contours.

Valid values: solid
 mesh
 dots

Default value: **solid**

crossvalidation

Leave-n-out cross validation.

Valid values: integers
Default value: **1**
Minimum: 1

electrostaticcutoff

Value for truncating electrostatic fields in kcal/mol.

Valid values: reals
Default value: **30**
Minimum: 0.1

eliminatetvalue

Controls whether or not to eliminate variables based on t-value

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

extendgrid

Indicates how much (in angstroms) to extend the grid beyond the limits of the training set ligands.

Valid values: reals
Default value: **3**
Minimum: 0.0

featuredefinitions

The file name of a feature definitions file to use with a custom field style.

Valid values: text strings

Default value:

fieldstyle

What type of field to use for generating the model.

Valid values: ff
gauss
custom

Default value: **gauss**

grid_spacing

The grid spacing to be for the QSAR model grid. The valid range is 0.5 to 4.0 angstrom.

Valid values: reals

Default value: **1**

Minimum: 0.5

Maximum: 4.0

ignoredistance

The distance with which to ignore fields.

Valid values: reals

Default value: **2**

Minimum: 0.0

ignorefields

Indicates whether or not Phase should ignore Force field fields.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

max_factors

This option value range depends on the ligands in the training set. The valid range is 1 to N/5, where N is the number of ligands in the training set. If N is less than 5, the maximum number of factors is 1, but no QSAR model can be constructed.

Valid values: integers

Default value: **1**

Minimum: 1

random

The percentage of the ligands to set as the training set.

Valid values: integers

Default value: **50**

Minimum: 0

Maximum: 100

Chapter 5: Commands

random_seed

The seed used to compute the random training / test set in the Individual QSAR Model panel. Zero means to use a completely random seed—any other value is used explicitly.

Valid values: integers
Default value: **0**
Minimum: 0

stddevcutoff

Threshold for standard deviation

Valid values: reals
Default value: **0.01**
Minimum: 0.001

stericcutoff

Value for truncating steric fields in kcal/mol.

Valid values: reals
Default value: **30**
Minimum: 0.1

truncateelectrostatic

Whether or not to truncate the Force field electrostatic fields.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

truncatesteric

Controls whether or not Phase truncates steric fields.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

tvaluecutoff

The cutoff value for eliminating variables based on t-value.

Valid values: reals
Default value: **2**

useinputpartialcharges

Whether to use the input structures' partial charges.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

phasefieldqsaraddligands

Adds the given ligands to the Individual QSAR Model panel.

Syntax:

```
phasefieldqsaraddligands activity=<text>
                        convert_activity=yes | no convert_scale=<x>
                        set_property=<text> test_value=<text> training_value=<text>
                        <ESL>
```

Options:

activity This determines which property (if any) to use as the activity property for the ligands.

Valid values: text strings

Default value:

convert_activity

Set to true if the activity values should be converted from concentration to $-\log[\text{concentration}]$

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

convert_scale

A scale factor for conversion.

Valid values: reals

Default value: **1**

set_property

This property is used to read the training / test set information from.

Valid values: text strings

Default value:

test_value The value which is considered part of the test set.

Valid values: text strings

Default value:

training_value

The value which is considered part of the training set.

Valid values: text strings

Default value:

Operands:

<ESL>

The entries to add as ligands.

phasefieldqsaraddligandsfromfile

Adds the ligands in the given files to the Phase Field-Based QSAR Model panel.

Syntax:

```
phasefieldqsaraddligandsfromfile <file 1> <file 2>
```

Operands:

<file 1> <file 2>

The files to add ligands from.

phasefieldqsaraddtoproject

Adds the selected ligands from the Field-Based QSAR Model panel to the Project Table.

Syntax:

```
phasefieldqsaraddtoproject
```

phasefieldqsarbuild

Build a QSAR model from the current ligands.

Syntax:

```
phasefieldqsarbuild
```

phasefieldqsardelete

Delete the selected ligands.

Syntax:

```
phasefieldqsardelete
```

phasefieldqsardeleteall

Deletes all of the ligands from the field-based QSAR panel.

Syntax:

```
phasefieldqsardeleteall
```

phasefieldqsarexport

Exports the QSAR model.

Syntax:

```
phasefieldqsarexport <file name>
```

Operands:

<file name>

The name of the file to save the QSAR model to.

phasefieldqsarexportligandstofile

Exports the selected ligands from the Field-Based QSAR model.

Syntax:

```
phasefieldqsarexportligandstofile <file name>
```

Operands:

<file name>

The file name to export the selected conformers to.

phasefieldqsarextendinclude

Extends the inclusion to the given ligand in the Field-Based QSAR Model panel.

Syntax:

phasefieldqsarextendinclude $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to include.

phasefieldqsarextendselectligandrow

Extends the selection to the given row in the ligands table in the Field-Based QSAR Model panel.

Syntax:

phasefieldqsarextendselectligandrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to extend-select.

phasefieldqsarimport

Imports a QSAR model from the given file.

Syntax:

phasefieldqsarimport $\langle \text{file} \rangle$

Operands:

$\langle \text{file} \rangle$

The QSAR model file to import.

phasefieldqsarinclude

Includes just the given ligand in the Field-Based QSAR Model panel in the Workspace.

Syntax:

phasefieldqsarinclude $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to include.

phasefieldqsarinvertselection

Invert the selection in the ligands table.

Syntax:

phasefieldqsarinvertselection

phasefieldqsarplotactivity

Creates a scatter plot of predicted versus actual activities for the Field-Based QSAR Model.

Syntax:

phasefieldqsarplotactivity *draw_line*=yes | no
pls_factor= $\langle n \rangle$ $\langle \text{all} | \text{selected} | \text{training} | \text{test} \rangle$

Options:

draw_line If set, then the plot will include a line at 45 degrees.
 Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

pls_factor This is which set of PLS factors to plot.
 Valid values: integers
 Default value: **1**
 Minimum: 1

Operands:

$\langle \text{all} | \text{selected} | \text{training} | \text{test} \rangle$

Which ligands to plot: all, selected, the training set, or the test set.

phasefieldqsarpredict

Predicts activities for the given entries using the Field-Based QSAR Model.

Syntax:

```
phasefieldqsarpredict <ESL>
```

Operands:

<ESL>

The entries to predict activities for.

phasefieldqsarrandomtraining

Sets the random training set for the Field-Based QSAR Model panel.

Syntax:

```
phasefieldqsarrandomtraining
```

phasefieldqsarselectall

Selects all of the ligands.

Syntax:

```
phasefieldqsarselectall
```

phasefieldqsarselectcustom

Selects just the given field in the Custom Field Style dialog

Syntax:

```
phasefieldqsarselectcustom <field>
```

Operands:

<field>

The custom field to select

phasefieldqsarselectligandrow

Selects the given row in the ligands table on the Field-Based QSAR Model panel.

Syntax:

```
phasefieldqsarselectligandrow <row>
```

Operands:

<row>

The row number to select.

phasefieldqsarsetactivity

Sets the activity for the given row in the Field-Based QSAR Model panel.

Syntax:

```
phasefieldqsarsetactivity activity=<x> <row>
```

Options:

activity The activity value for the given ligand.

Valid values: reals

Default value: **1**

Operands:

<row>

The row number of the ligand to set the activity for.

phasefieldqsarsettrainingrows

Toggles the Training Set property on or off for the given row in the ligands table on the Field-Based QSAR Model panel.

Syntax:

```
phasefieldqsarsettrainingrows value=training | test | none  
    <rows>
```

Options:

value The value of the training/test column in the ligands table in the Field-Based QSAR Model panel.

Valid values: training

test

none

Default value: **training**

Operands:

<rows>

The row numbers to toggle, two row numbers should be separated by , or if the user wants to specify a range then it can be given like 1-5 separated by a - . eg 1,2,5 or 1,4-7 .

phasefieldqsarsort

Sort the Field-Based QSAR Model ligands table based on the data in the specified column

Syntax:

```
phasefieldqsarsort <column_name>
```

Operands:

<column_name>

The name of the column to be sorted.

phasefieldqsartest

Tests a QSAR model on the current ligand test set.

Syntax:

```
phasefieldqsartest
```

phasefieldqsartoggleinclude

Toggles the given ligand in the Field-Based QSAR Model panel into or out of the Workspace.

Syntax:

```
phasefieldqsartoggleinclude <row>
```

Operands:

<row>

The row number to toggle.

phasefieldqsartoggleligandrow

Toggles the selection for the given row in the ligands table in the Field-Based QSAR Model panel.

Syntax:

```
phasefieldqsartoggleligandrow <row>
```

Operands:

<row>

The row number to toggle-select.

phasefieldqsartogglerow

Toggles the Training Set property on or off for the given row in the Field-Based QSAR Model panel.

Syntax:

```
phasefieldqsartogglerow <row>
```

Operands:

<row>

The row number to toggle.

phasefieldqsartoggleselectcustom

Toggles the given field in the Custom Field Style dialog.

Syntax:

phasefieldqsartoggleselectcustom $\langle \text{field} \rangle$

Operands:

$\langle \text{field} \rangle$

The custom field to toggle select

phasefieldqsarvisselectcontour

Selects just the given field contour to visualize.

Syntax:

phasefieldqsarvisselectcontour $\langle \text{field} \rangle$

Operands:

$\langle \text{field} \rangle$

The field contour to select for visualization

phasefieldqsarvisselectintensity

Selects just the given field intensity to visualize.

Syntax:

phasefieldqsarvisselectintensity $\langle \text{field} \rangle$

Operands:

$\langle \text{field} \rangle$

The field intensity to select for visualization

phasefieldqsarvistoggleselectcontour

Toggle-selects the given field contour to visualize.

Syntax:

phasefieldqsarvistoggleselectcontour \langle field \rangle

Operands:

\langle field \rangle

The field contour to toggle-select for visualization

phasefindmatches

Launches a Find Matches job for start operand. Writes Find Matches job input files for write operand.

Syntax:

phasefindmatches *align_cutoff*= \langle x \rangle *align_weight*= \langle x \rangle
apply_constraints=yes | no *apply_excluded_volumes*=yes | no
apply_must_match=yes | no *apply_permitted*=yes | no
apply_tolerances=yes | no *computeincludedvolume*=yes | no
conformers=existing | generate *cpu_limit*= \langle x \rangle
database= \langle text \rangle *distance_tolerance*= \langle x \rangle *entry_name*= \langle text \rangle
extfile= \langle text \rangle *generate_sites*=yes | no *hits_molecule*= \langle n \rangle
hits_total= \langle n \rangle *included_volume_cutoff*= \langle x \rangle
included_volume_weight= \langle x \rangle *jobinput*=extfile | selectedentries
| 3ddatabase *matches*=new | existing *minsites*= \langle n \rangle
prefer_big_matches=yes | no *refine*=yes | no
remote_database=yes | no *saveallmatches*=yes | no
score_in_place=yes | no *sorthits*=yes | no *subset*= \langle text \rangle
use_align_cutoff=yes | no *use_cpu_limit*=yes | no
use_distance_tolerance=yes | no
use_included_volume_cutoff=yes | no
use_vector_cutoff=yes | no *use_volume_cutoff*=yes | no
use_volume_groups=yes | no *usedbkeys*=yes | no
useqsar=yes | no *vector_cutoff*= \langle x \rangle *vector_weight*= \langle x \rangle
verbose=yes | no *volume_cutoff*= \langle x \rangle *volume_weight*= \langle x \rangle
 \langle write|start \rangle

Options:

align_cutoff

The cutoff for the alignment score.

Valid values: reals

Default value: **1.2**

Minimum: 0.0001

Chapter 5: Commands

align_weight

The weight for the alignment score

Valid values: reals

Default value: **1**

Minimum: 0.0

apply_constraints

Whether or not to apply constraints.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

apply_excluded_volumes

Indicates whether or not to use the excluded volumes information.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

apply_must_match

Whether or not to apply site must-match values.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

apply_permitted

Whether or not to apply site-based permitted and prohibited matching criteria.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

apply_tolerances

Whether or not to apply site tolerances.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

computeincludedvolume

Indicates whether or not to compute the included volume scores.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

conformers

Whether to use existing conformers or generate them during the search

Valid values: existing
generate

Default value: **generate**

cpu_limit

The CPU limit for each molecule, in seconds.

Valid values: reals

	Default value: 0.1
	Minimum: 0.1
<i>database</i>	The source database file for Phase Find Matches.
	Valid values: text strings
	Default value:
<i>distance_tolerance</i>	
	Valid values: reals
	Default value: 2
<i>entry_name</i>	
	The entry name to use as the hypothesis to search against.
	Valid values: text strings
	Default value:
<i>extfile</i>	The source structure .sd or .mae file for Phase Find Matches.
	Valid values: text strings
	Default value:
<i>generate_sites</i>	
	Use hypothesis feature definitions to generate sites in memory for each database structure searched.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>hits_molecule</i>	
	The maximum number of hits per molecule to return.
	Valid values: integers
	Default value: 1
	Minimum: 1
<i>hits_total</i>	The maximum number of hits to return.
	Valid values: integers
	Default value: 1000
	Minimum: 1
<i>included_volume_cutoff</i>	
	The cutoff for the included volume score.
	Valid values: reals
	Default value: 0
	Minimum: 0.0
	Maximum: 1.0
<i>included_volume_weight</i>	
	The weight for included volume scores.
	Valid values: reals
	Default value: 1
	Minimum: 0.0

Chapter 5: Commands

<i>jobinput</i>	<p>The source of job input of Phase Find Matches.</p> <p>Valid values: extfile selectedentries 3ddatabase</p> <p>Default value: 3ddatabase</p>
<i>matches</i>	<p>Whether to search for new or existing matches.</p> <p>Valid values: new existing</p> <p>Default value: new</p>
<i>minsites</i>	<p>The minimum number of sites to match.</p> <p>Valid values: integers</p> <p>Default value: 5</p> <p>Minimum: 1</p>
<i>prefer_big_matches</i>	<p>If this option is true, there is a preference for partial matches that involve a larger number of sites.</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: true</p>
<i>refine</i>	<p>Whether or not to refine matches</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>
<i>remote_database</i>	<p>If set, then this means that we have a remote database which isn't necessarily accessible to the local machine.</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>
<i>saveallmatches</i>	<p>If this option is true, all the matches will be written to disk, to the file matches.out.</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: true</p>
<i>score_in_place</i>	<p>Indicates whether or not the structures get aligned to the pharmacophore before being scored.</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>
<i>sorthits</i>	<p>Sort hits by decreasing fitness</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: true</p>

- subset* The name of the subset for the Phase Find Matches.
Valid values: text strings
Default value:
- use_align_cutoff*
The flag for filtering hits based on align score. If this option is true, reject hits with align score greater than align_cutoff.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**
- use_cpu_limit*
If set, then use a per-molecule CPU time limit.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**
- use_distance_tolerance*
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**
- use_included_volume_cutoff*
The flag for filtering hits based on included volume score. If this option is true, reject hits with an included volume score less than the cutoff.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**
- use_vector_cutoff*
The flag for filtering hits based on vector score. If this option is true, reject hits with vector score less than vector_cutoff.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**
- use_volume_cutoff*
The flag for filtering hits based on volume score. If this option is true, reject hits with volume score less than volume_cutoff.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**
- use_volume_groups*
If this option is true, volume scores will be higher for alignments that overlap atoms of the same MacroModel type.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**
- usedbkeys* If this option is true, phase_dbsearch uses database keys to perform a single-CPU pre-screen of the database to eliminate molecules.

Chapter 5: Commands

	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>useqsar</i>	The flag of using QSAR model or not.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>vector_cutoff</i>	The cutoff for the vector score.
	Valid values: reals
	Default value: -1
	Minimum: -1.0
	Maximum: 1.0
<i>vector_weight</i>	The weight for the vector score.
	Valid values: reals
	Default value: 1
	Minimum: 0.0
<i>verbose</i>	If set, the Find Matches job will produce more verbose output.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>volume_cutoff</i>	The cutoff for the volume score.
	Valid values: reals
	Default value: 0
	Minimum: 0.0
	Maximum: 1.0
<i>volume_weight</i>	The weight for the volume score.
	Valid values: reals
	Default value: 1
	Minimum: 0.0

Operands:

⟨ write|start ⟩

This should be set to Find Matches .

phasefindpharmacophores

Sets the Find Common Pharmacophores options and launches the job.

Syntax:

phasefindpharmacophores *finalboxsize*=⟨ x ⟩ *maxtreedepth*=⟨ n ⟩
minintersitedistance=⟨ x ⟩

Options:

finalboxsize

The final box size in angstrong.

Valid values: reals

Default value: **1**

maxtreedepth

The maximum depth of the tree.

Valid values: integers

Default value: **5**

minintersitedistance

The minimum distance between the sites in angstrong.

Valid values: reals

Default value: **2**

phasesgenerateconformers

Launches a Generate Conformers job for Phase.

Syntax:

phasesgenerateconformers

phasesgroupligands

Group together ligands

Syntax:

phasesgroupligands

phasesgroupligandsbytitle

Group together ligands by title

Syntax:

`phasesgroupligandsbytitle`

phasehookimport

Sets up the Import panel and dialog for Phase to import ligands and conformers. This command is deprecated as of Maestro 8.0. Use `phaseaddligandsfromfile`

Syntax:

`phasehookimport`

phasehypothesiscreateexcludedvolume

Creates an excluded volume from the current centroid atoms for the currently selected hypothesis.

Syntax:

`phasehypothesiscreateexcludedvolume`

phasehypothesisdeleteexcludedvolumes

Deletes the selected excluded volumes for the currently selected hypothesis.

Syntax:

`phasehypothesisdeleteexcludedvolumes`

phasehypothesisselectevrow

Selects the given row in the excluded volumes table for the currently selected hypothesis.

Syntax:

phasehypothesisselectevrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to select in the table.

phasehypothesisselectextendevrow

Extends the selection to this row in the excluded volumes table for the currently selected hypothesis.

Syntax:

phasehypothesisselectextendevrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to extend the select to.

phasehypothesisselectonlyevrow

Selects only this row in the excluded volumes table for the currently selected hypothesis.

Syntax:

phasehypothesisselectonlyevrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to select only in the table row.

phasehypothesisselectrow

Toggles the given hypothesis into or out of the Workspace based on the currently selected hypothesis.

Syntax:

phasehypothesisselectrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to toggle inclusion state for.

phasehypothesissetexcludedvolumes

Sets the value for the given cell for the currently selected hypothesis.

Syntax:

phasehypothesissetexcludedvolumes *column*= $\langle n \rangle$ *row*= $\langle n \rangle$
 $\langle \text{value} \rangle$

Options:

column The column number of the cell to change.

Valid values: integers

Default value: **1**

row The row number of the cell to change.

Valid values: integers

Default value: **1**

Operands:

$\langle \text{value} \rangle$

The value for the given cell.

phasehypothesissetid

Sets the value for the hypothesis ID for the currently selected hypothesis.

Syntax:

phasehypothesissetid *row*= $\langle n \rangle$ $\langle \text{value} \rangle$

Options:

row The row number of the hypothesis ID to change.
 Valid values: integers
 Default value: **1**

Operands:

⟨ value ⟩

The value for the hypothesis ID.

phasehypothesissetsitemask

Sets the itemask for the currently selected hypothesis

Syntax:

phasehypothesissetsitemask ⟨ value ⟩

Operands:

⟨ value ⟩

The value for the given cell.

phasehypothesissettolerances

Sets the feature matching tolerances for the currently selected hypothesis

Syntax:

phasehypothesissettolerances *tolerance_a*=⟨ x ⟩
 tolerance_d=⟨ x ⟩ *tolerance_h*=⟨ x ⟩ *tolerance_n*=⟨ x ⟩
 tolerance_p=⟨ x ⟩ *tolerance_r*=⟨ x ⟩ *tolerance_x*=⟨ x ⟩
 tolerance_y=⟨ x ⟩ *tolerance_z*=⟨ x ⟩ *usetolerances*=yes | no

Options:

tolerance_a

The feature matching tolerance for the hydrogen bond acceptor feature. The valid range is 0.0 to 100.0

Valid values: reals

Default value: **1**

Minimum: 0.0

Maximum: 100.0

tolerance_d

The feature matching tolerance for the hydrogen bond donor feature. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 100.0

tolerance_h

The feature matching tolerance for the hydrophobic feature. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **1.5**
Minimum: 0.0
Maximum: 100.0

tolerance_n

The feature matching tolerance for the negative feature. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **0.75**
Minimum: 0.0
Maximum: 100.0

tolerance_p

The feature matching tolerance for the positive feature. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **0.75**
Minimum: 0.0
Maximum: 100.0

tolerance_r

The feature matching tolerance for the aromatic ring feature. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **1.5**
Minimum: 0.0
Maximum: 100.0

tolerance_x

The feature matching tolerance for the custom(X) feature. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 100.0

tolerance_y

The feature matching tolerance for the custom(Y) feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 100.0

tolerance_z

The feature matching tolerance for the custom(Z) feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 100.0

usetolerances

When using feature matching tolerances, this option is true otherwise false.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

phasehypothesistoggleinclude

Toggles the given hypothesis into or out of the Workspace.

Syntax:

phasehypothesistoggleinclude <row>

Operands:

<row>

The row number to toggle inclusion state for.

phasehypothesisunselectevrow

Unselects the given row in the excluded volumes table for the currently selected hypothesis.

Syntax:

phasehypothesisunselectevrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to unselect in the table.

phaseimportfeature

Imports the features from the given file name. This replaces the existing features.

Syntax:

phaseimportfeature $\langle \text{file name} \rangle$

Operands:

$\langle \text{file name} \rangle$

A path to a valid features file.

phaseimporthypothesis

Imports a hypothesis from the selected file into the Find Matches and Edit Hypotheses panels.

Syntax:

phaseimporthypothesis *isphase*=yes | no $\langle \text{file name} \rangle$

Options:

isphase Indicates whether or not the hypothesis is a Phase hypothesis or an external hypothesis.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

$\langle \text{file name} \rangle$

The file name to import a hypothesis from.

phaseincludeextendtablerow

Extends the rows included in the workspace to include this one.

Syntax:

```
phaseincludeextendtablerow <row>
```

Operands:

<row>

The row number to include in the Workspace.

phaseincludeonlytablerow

Includes only the given row in the first table in the step into the Workspace.

Syntax:

```
phaseincludeonlytablerow <row>
```

Operands:

<row>

The row number to include in the Workspace.

phaseincludetablerow

Includes the given row in the first table in the step into the Workspace.

Syntax:

```
phaseincludetablerow <row>
```

Operands:

<row>

The row number to include in the Workspace.

phaseinverttableselection

Inverts the row selection in the first table in the step.

Syntax:

`phaseinverttableselection`

phaseligprep

Defines settings for Phase LigPrep job.

Syntax:

`phaseligprep` *gen_stereo*=⟨n⟩ *ionization*=generate | neutralize |
retain *ph*=⟨x⟩ *stereoisomers*=retain | determine | generate

Options:

<i>gen_stereo</i>	Generate stereoisomers (maximum)
Valid values:	integers
Default value:	32
Minimum:	1
<i>ionization</i>	How to generate the ionization states
Valid values:	generate neutralize retain
Default value:	retain
<i>ph</i>	The ionization pH
Valid values:	reals
Default value:	7
Minimum:	0.0
<i>stereoisomers</i>	How to generate the stereoisomers
Valid values:	retain determine generate
Default value:	retain

phaselocatesites

Launches a Locate Sites job.

Syntax:

~~phaselocatesites~~
phasemarkerdump

Print out the current option values of the phase marker command.

Syntax:

~~phasemarkerdump~~
phasemarkersettings

Set graphical data of feature markers.

Syntax:

```

phasesettings acceptorarrowblue=<x>
    acceptorarrowgreen=<x> acceptorarrowred=<x>
    acceptorcornradius=<x> acceptorcylinderheight=<x>
    acceptorcylinderradius=<x> acceptorsiteblue=<x>
    acceptorsitegreen=<x> acceptorsiteradius=<x>
    acceptorsitered=<x> ambient=<x> aromaticringradius=<x>
    aromaticringradius5=<x> aromaticsiteblue=<x>
    aromaticsitegreen=<x> aromaticsitered=<x>
    aromatictuberadius=<x> aromaticustep=<n>
    aromaticvstep=<n> diffuse=<x> donorarrowblue=<x>
    donorarrowgreen=<x> donorarrowred=<x>
    donorcornradius=<x> donorcylinderheight=<x>
    donorcylinderradius=<x> donorsiteblue=<x>
    donorsitegreen=<x> donorsiteradius=<x> donorsitered=<x>
    drawstyle=solid | line emission=<x> featurexsiteblue=<x>
    featurexsitegreen=<x> featurexsiteradius=<x>
    featurexsitered=<x> featureysiteblue=<x>
    featureysitegreen=<x> featureysiteradius=<x>
    featureysitered=<x> featurezsiteblue=<x>
    featurezsitegreen=<x> featurezsiteradius=<x>
    featurezsitered=<x> hydrophobicsiteblue=<x>
    hydrophobicsitegreen=<x> hydrophobicsiteradius=<x>
    hydrophobicsitered=<x> labelblue=<x> labelgreen=<x>
    labelred=<x> linewidth=<n> negativesiteblue=<x>
    negativesitegreen=<x> negativesiteradius=<x>
    negativesitered=<x> positivesiteblue=<x>
    positivesitegreen=<x> positivesiteradius=<x>
    positivesitered=<x> projectedsiteblue=<x>
    projectedsitegreen=<x> projectedsiteradius=<x>
    projectedsitered=<x> shininess=<x> sitedimpercentage=<x>
    sliceline=<n> slicesolid=<n> specular=<x> stackline=<n>
    stacksolid=<n> transparency=<x>

```

Options:

acceptorarrowblue

The arrow color blue component of acceptor markers.

Valid values: reals
 Default value: **0.5**
 Minimum: 0.0
 Maximum: 1.0

acceptorarrowgreen

The arrow color green component of acceptor markers.

Valid values: reals
 Default value: **0.5**

Minimum: 0.0
Maximum: 1.0

acceptorarrowred

The arrow color red component of acceptor markers.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

acceptorcornradius

The radius of arrow corn of acceptor markers.

Valid values: reals
Default value: **0.35**
Minimum: 0.0

acceptorcylinderheight

The height of arrow cylinder of acceptor markers.

Valid values: reals
Default value: **0.6**
Minimum: 0.0
Maximum: 1.0

acceptorcylinderradius

The radius of arrow cylinder of acceptor markers.

Valid values: reals
Default value: **0.15**
Minimum: 0.0

acceptorsiteblue

The site color blue component of acceptor markers.

Valid values: reals
Default value: **0.5**
Minimum: 0.0
Maximum: 1.0

acceptorsitegreen

The site color green component of acceptor markers.

Valid values: reals
Default value: **0.5**
Minimum: 0.0
Maximum: 1.0

acceptorsiteradius

The radius of acceptor site markers.

Valid values: reals
Default value: **0.8**
Minimum: 0.0

Chapter 5: Commands

acceptorsitered

The site color red component of acceptor markers.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

ambient

Set material property - ambient, to its red, green, and blue components, for front face.

Valid values: reals
Default value: **0.5**
Minimum: 0.0
Maximum: 1.0

aromaticringradius

The radius of the ring of the ring torus, for the ring constructed by 6 or more atoms.

Valid values: reals
Default value: **0.5**
Minimum: 0.0

aromaticringradius5

The radius of the ring of the ring torus, for the ring constructed by 5 or less atoms.

Valid values: reals
Default value: **0.4**
Minimum: 0.0

aromaticsiteblue

The site color blue component of aromatic markers.

Valid values: reals
Default value: **0.2**
Minimum: 0.0
Maximum: 1.0

aromaticsitegreen

The site color green component of aromatic markers.

Valid values: reals
Default value: **0.51**
Minimum: 0.0
Maximum: 1.0

aromaticsitered

The site color red component of aromatic markers.

Valid values: reals
Default value: **0.92**

Minimum: 0.0
Maximum: 1.0

aromatictuberadius

The radius of the tube of the ring torus.

Valid values: reals
Default value: **0.15**
Minimum: 0.0

aromaticustep

The U step domain tolerance.

Valid values: integers
Default value: **5**
Minimum: 2

aromaticvstep

The V step domain tolerance.

Valid values: integers
Default value: **3**
Minimum: 2

diffuse

Set material property - diffuse, to its red, green, and blue components, for front face.

Valid values: reals
Default value: **0.4**
Minimum: 0.0
Maximum: 1.0

donorarrowblue

The arrow color blue component of donor markers.

Valid values: reals
Default value: **0.9**
Minimum: 0.0
Maximum: 1.0

donorarrowgreen

The arrow color green component of donor markers.

Valid values: reals
Default value: **0.8**
Minimum: 0.0
Maximum: 1.0

donorarrowred

The arrow color red component of donor markers.

Valid values: reals
Default value: **0.5**
Minimum: 0.0
Maximum: 1.0

Chapter 5: Commands

donorcornradius

The radius of arrow corn of donor markers.

Valid values: reals
Default value: **0.35**
Minimum: 0.0

donorcylinderheight

The height of arrow cylinder of donor markers.

Valid values: reals
Default value: **0.6**
Minimum: 0.0
Maximum: 1.0

donorcylinderradius

The radius of arrow cylinder of donor markers.

Valid values: reals
Default value: **0.15**
Minimum: 0.0

donorsiteblue

The site color blue component of donor markers.

Valid values: reals
Default value: **0.9**
Minimum: 0.0
Maximum: 1.0

donorsitegreen

The site color green component of donor markers.

Valid values: reals
Default value: **0.8**
Minimum: 0.0
Maximum: 1.0

donorsiteradius

The radius of donor site markers.

Valid values: reals
Default value: **0.8**
Minimum: 0.0

donorsitered

The site color red component of donor markers.

Valid values: reals
Default value: **0.5**
Minimum: 0.0
Maximum: 1.0

- drawstyle* The styles of rendering feature markers, they are: 1 - solid, and 2 - lines. Default is solid.
- Valid values: solid
line
- Default value: **solid**
- emission* Set material property - emission, to its red, green, and blue components, for front face.
- Valid values: reals
- Default value: **0.05**
- Minimum: 0.0
- Maximum: 1.0
- featurexsiteblue*
- The site color blue component of featurex markers.
- Valid values: reals
- Default value: **1**
- Minimum: 0.0
- Maximum: 1.0
- featurexsitegreen*
- The site color green component of featurex markers.
- Valid values: reals
- Default value: **1**
- Minimum: 0.0
- Maximum: 1.0
- featurexsiteradius*
- The radius of featurex site markers.
- Valid values: reals
- Default value: **0.8**
- Minimum: 0.0
- featurexsitered*
- The site color red component of featurex markers.
- Valid values: reals
- Default value: **0**
- Minimum: 0.0
- Maximum: 1.0
- featureysiteblue*
- The site color blue component of featurey markers.
- Valid values: reals
- Default value: **1**
- Minimum: 0.0
- Maximum: 1.0

Chapter 5: Commands

featureysitegreen

The site color green component of featurey markers.

Valid values: reals
Default value: **0**
Minimum: 0.0
Maximum: 1.0

featureysiteradius

The radius of featurey site markers.

Valid values: reals
Default value: **0.8**
Minimum: 0.0

featureysitered

The site color red component of featurey markers.

Valid values: reals
Default value: **0.3**
Minimum: 0.0
Maximum: 1.0

featurezsiteblue

The site color blue component of featurez markers.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

featurezsitegreen

The site color green component of featurez markers.

Valid values: reals
Default value: **0**
Minimum: 0.0
Maximum: 1.0

featurezsiteradius

The radius of featurez site markers.

Valid values: reals
Default value: **0.8**
Minimum: 0.0

featurezsitered

The site color red component of featurez markers.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

hydrophobicsiteblue

The site color blue component of hydrophobic markers.

Valid values: reals
 Default value: **0.3**
 Minimum: 0.0
 Maximum: 1.0

hydrophobicsitegreen

The site color green component of hydrophobic markers.

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 1.0

hydrophobicsiteradius

The radius of hydrophobic site markers.

Valid values: reals
 Default value: **0.8**
 Minimum: 0.0

hydrophobicsitered

The site color red component of hydrophobic markers.

Valid values: reals
 Default value: **0.3**
 Minimum: 0.0
 Maximum: 1.0

labelblue

The blue component of feature labels.

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 1.0

labelgreen

The green component of feature labels.

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 1.0

labelred

The red component of feature labels.

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 1.0

linewidth

Set the width of lines in drawing sphere.

Valid values: integers

Chapter 5: Commands

Default value: **2**
Minimum: 1

negativesiteblue

The site color blue component of negative markers.

Valid values: reals
Default value: **0.3**
Minimum: 0.0
Maximum: 1.0

negativesitegreen

The site color green component of negative markers.

Valid values: reals
Default value: **0.3**
Minimum: 0.0
Maximum: 1.0

negativesiteradius

The radius of negative site markers.

Valid values: reals
Default value: **0.8**
Minimum: 0.0

negativesitered

The site color red component of negative markers.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

positivesiteblue

The site color blue component of positive markers.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

positivesitegreen

The site color green component of positive markers.

Valid values: reals
Default value: **0.6**
Minimum: 0.0
Maximum: 1.0

positivesiteradius

The radius of positive site markers.

Valid values: reals

Default value: **0.8**
 Minimum: 0.0

positivesitered

The site color red component of positive markers.

Valid values: reals
 Default value: **0.5**
 Minimum: 0.0
 Maximum: 1.0

projectedsiteblue

The site color blue component of projected markers.

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 1.0

projectedsitegreen

The site color green component of projected markers.

Valid values: reals
 Default value: **0**
 Minimum: 0.0
 Maximum: 1.0

projectedsiteradius

The radius of projected site markers.

Valid values: reals
 Default value: **0.8**
 Minimum: 0.0

projectedsitered

The site color red component of projected markers.

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 1.0

shininess Set material property - shininess, for front face.

Valid values: reals
 Default value: **80**
 Minimum: 0.0
 Maximum: 128.0

sitedimpercentage

The percentage to dim Phase sites by.

Valid values: reals
 Default value: **0.5**

Chapter 5: Commands

	Minimum:	0.0
	Maximum:	1.0
<i>sliceline</i>	Set the slices of drawing line sphere.	
	Valid values:	integers
	Default value:	10
	Minimum:	2
<i>slicesolid</i>	Set the slices of drawing solid sphere.	
	Valid values:	integers
	Default value:	36
	Minimum:	2
<i>specular</i>	Set material property - specular, to its red, green, and blue components, for front face.	
	Valid values:	reals
	Default value:	0
	Minimum:	0.0
	Maximum:	1.0
<i>stackline</i>	Set the stacks of drawing line sphere.	
	Valid values:	integers
	Default value:	8
	Minimum:	2
<i>stacksolid</i>	Set the stacks of drawing solid sphere.	
	Valid values:	integers
	Default value:	18
	Minimum:	2
<i>transparency</i>	The transparency of rendering feature markers.	
	Valid values:	reals
	Default value:	20
	Minimum:	0.0
	Maximum:	100.0

phasemarkfeature

Marks the given feature in the Workspace.

Syntax:

phasemarkfeature $\langle \text{feature} \rangle$

Operands:

$\langle \text{feature} \rangle$

A single letter (A-Z) indicating the feature to mark.

phasemergestereoisomers

Merges stereoisomers for the selected conformers.

Syntax:

phasemergestereoisomers

phasemergetitles

Merges the selected conformers by title.

Syntax:

phasemergetitles

phaseoptions

Sets the options for Phase.

Syntax:

```

phaseoptions convert_hypothesis=yes | no
                include_reference_ligand=yes | no keep_ligands=yes | no
                prompt_convert=yes | no random_seed=⟨ n ⟩
                separate_stereoisomers=yes | no show_unmatched=yes | no
                showhypothesisangles=yes | no
                showhypothesisconstraints=yes | no
                showhypothesisdistanes=yes | no
                showhypothesisexcludedvolumes=yes | no
                showhypothesislabels=yes | no
                showhypothesistolerance=yes | no split_by_title=yes | no

```

Options:

convert_hypothesis

If the *prompt_convert* option is turned off, then this option controls whether or not hypotheses are converted to the current format as needed.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

include_reference_ligand

If true then includes the reference ligand when including hypothesis from DPM panel.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

keep_ligands

Controls whether or not the same ligands are kept in the Workspace when switching hypotheses in the Score Hypotheses and Build QSAR steps.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

prompt_convert

If set, then Find Matches will prompt for whether or not to convert hypotheses as needed.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

random_seed

The seed used to compute the random training / test set in the Build QSAR step. Zero means to use a completely random seed—any other value is used explicitly.

Valid values: integers
 Default value: **0**
 Minimum: 0

separate_stereoisomers

Whether to separate stereoisomers or combine them into a single confset. Deprecated.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

show_unmatched

Controls whether or not unmatched ligands are displayed in the alignments table in Score Hypotheses.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

showhypothesisangles

Controls whether or not angles are displayed when hypotheses are included in the Workspace.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

showhypothesisconstraints

Controls whether or not hypothesis constraints are displayed in the Workspace.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

showhypothesisdistances

Controls whether or not inter-site distances are displayed when hypotheses are included in the Workspace.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

showhypothesisexcludedvolumes

Controls whether or not excluded volumes are displayed when hypotheses are included in the Workspace.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

showhypothesislabels

Controls whether or not hypothesis labels are displayed in the Workspace.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

showhypothesistolerance

Controls whether or not hypotheses are shown as tolerance markers.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

split_by_title

Controls whether or not conformers are split by title.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

phasepatterndelete

Deletes the given pattern.

Syntax:

phasepatterndelete *pattern*=⟨text⟩ ⟨feature⟩

Options:

pattern The SMARTS pattern to remove.

Valid values: text strings

Default value:

Operands:

⟨feature⟩

The feature to remove the pattern from.

phasepatternedit

Updates the given pattern with the new data.

Syntax:

phasepatternedit *geometry*=group | point | vector *group*=⟨text⟩
group_all=yes | no *pattern*=⟨text⟩ *point*=⟨n⟩ *vector*=⟨n⟩
⟨feature⟩

Options:

geometry The type of geometry for this pattern.

Valid values: group
 point
 vector

Default value: **group**

<i>group</i>	The atom numbers for the pattern if it is group geometry. Valid values: text strings Default value:
<i>group_all</i>	Set to true if the group should include all atom numbers. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>pattern</i>	The SMARTS pattern to update. Valid values: text strings Default value:
<i>point</i>	The atom number for the pattern if it is point geometry. Valid values: integers Default value: 1 Minimum: 1
<i>vector</i>	The atom number for the pattern if it is vector geometry. Valid values: integers Default value: 1 Minimum: 1

Operands:

⟨ feature ⟩

The feature to which the pattern is associated with.

phasepatternmovedown

Moves the given pattern down lower in the list for the feature.

Syntax:

phasepatternmovedown *pattern*=⟨ text ⟩ ⟨ feature ⟩

Options:

<i>pattern</i>	The SMARTS pattern to move.
Valid values:	text strings
Default value:	

Operands:

⟨ feature ⟩

The feature which contains the pattern.

phasepatternmoveup

Moves the given pattern up higher in the list for the feature.

Syntax:

phasepatternmoveup *pattern*=⟨text⟩ ⟨feature⟩

Options:

pattern The SMARTS pattern to move.
Valid values: text strings
Default value:

Operands:

⟨feature⟩

The feature which contains the pattern.

phasepatternnew

Adds the given pattern.

Syntax:

phasepatternnew *geometry*=group | point | vector *group*=⟨text⟩
group_all=yes | no *pattern*=⟨text⟩ *point*=⟨n⟩ *vector*=⟨n⟩
⟨feature⟩

Options:

geometry The type of geometry for this pattern.
Valid values: group
 point
 vector
Default value: **group**

group The atom numbers for the pattern if it is group geometry.
Valid values: text strings
Default value:

group_all Set to true if the group should include all atom numbers.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

<i>pattern</i>	The SMARTS pattern to add. Valid values: text strings Default value:
<i>point</i>	The atom number for the pattern if it is point geometry. Valid values: integers Default value: 1 Minimum: 1
<i>vector</i>	The atom number for the pattern if it is vector geometry. Valid values: integers Default value: 1 Minimum: 1

Operands:

⟨ feature ⟩

The feature to add the pattern to.

phasepatternsetoptions

Sets the options for the given pattern.

Syntax:

phasepatternsetoptions *exclude*=yes | no *ignore*=yes | no
mark=yes | no *pattern*=⟨ text ⟩ ⟨ feature ⟩

Options:

<i>exclude</i>	Indicates whether or not to exclude the pattern from the feature. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>ignore</i>	Indicates whether or not to ignore the pattern in the feature. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>mark</i>	Indicates whether or not to mark the pattern in the Workspace. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>pattern</i>	The SMARTS pattern to update. Valid values: text strings Default value:

Operands:

⟨ feature ⟩

The feature which contains the pattern.

phaseplotactivity

Creates a scatter plot of predicted versus actual activities.

Syntax:

```
phaseplotactivity draw_line=yes | no pls_factor=⟨ n ⟩  
                  ⟨ all | selected | training | test ⟩
```

Options:

<i>draw_line</i>	If set, then the plot will include a line at 45 degrees. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>pls_factor</i>	This is which set of PLS factors to plot. Valid values: integers Default value: 1 Minimum: 1

Operands:

⟨ all | selected | training | test ⟩

Which ligands to plot: all, selected, the training set, or the test set.

phasepreview

Launches a Preview job.

Syntax:

```
phasepreview
```

phaseqsar

Command for controlling Phase Individual QSAR Models

Syntax:

```
phasesar crossvalidation=<n> eliminatetvalue=yes | no
        grid_spacing=<x> max_factors=<n> random=<n>
        random_seed=<n> tvalue=<x>
```

Options:

crossvalidation

Leave-n-out cross validation.

Valid values: integers

Default value: **1**

Minimum: 1

eliminatetvalue

Indicates whether to use t-value filter before launching QSAR job.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

grid_spacing

The grid spacing to be for the QSAR model grid. The valid range is 0.5 to 2.0 angstrom.

Valid values: reals

Default value: **1**

Minimum: 0.5

Maximum: 2.0

max_factors

This option value range depends on the ligands in the training set. The valid range is 1 to N/5, where N is the number of ligands in the training set. If N is less than 5, the maximum number of factors is 1, but no QSAR model can be constructed.

Valid values: integers

Default value: **1**

Minimum: 1

random

The percentage of the ligands to set as the training set.

Valid values: integers

Default value: **50**

Minimum: 0

Maximum: 100

random_seed

The seed used to compute the random training / test set in the Individual QSAR Model panel. Zero means to use a completely random seed—any other value is used explicitly.

	Valid values:	integers
	Default value:	0
	Minimum:	0
<i>tvalue</i>	This option value indicates the values to be used to filter t-value.	
	Valid values:	reals
	Default value:	2
	Minimum:	0.01

phaseqsaraddligands

Adds the given ligands to the Individual QSAR Model panel.

Syntax:

```
phaseqsaraddligands activity=<text> convert_activity=yes | no
                    convert_scale=<x> set_property=<text> test_value=<text>
                    training_value=<text> <ESL>
```

Options:

<i>activity</i>	This determines which property (if any) to use as the activity property for the ligands.
	Valid values: text strings
	Default value:
<i>convert_activity</i>	Set to true if the activity values should be converted from concentration to -log[concentration]
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>convert_scale</i>	A scale factor for conversion.
	Valid values: reals
	Default value: 1
<i>set_property</i>	This property is used to read the training / test set information from.
	Valid values: text strings
	Default value:
<i>test_value</i>	The value which is considered part of the test set.
	Valid values: text strings
	Default value:

training_value

The value which is considered part of the training set.

Valid values: text strings

Default value:

Operands:

⟨ESL⟩

The entries to add as ligands.

phaseqsaraddligandsfromfile

Adds the ligands in the given files to the Phase Individual QSAR Model panel.

Syntax:

phaseqsaraddligandsfromfile ⟨file 1⟩ ⟨file 2⟩

Operands:

⟨file 1⟩ ⟨file 2⟩

The files to add ligands from.

phaseqsaraddtohypothesis

Adds the current Individual QSAR Model to the given entries, which must be hypotheses.

Syntax:

phaseqsaraddtohypothesis ⟨ESL⟩

Operands:

⟨ESL⟩

The entries to add the current QSAR model to.

phaseqsaraddtoproject

Adds the selected ligands from the Individual QSAR Model panel to the Project Table.

Syntax:

```
phaseqsaraddtoproject
```

phaseqsarbuild

Build a QSAR model from the current ligands.

Syntax:

```
phaseqsarbuild
```

phaseqsardelete

Delete the selected ligands.

Syntax:

```
phaseqsardelete
```

phaseqsardeleteall

Deletes all of the ligands from the atom-based QSAR panel.

Syntax:

```
phaseqsardeleteall
```

phaseqsarexport

Exports the QSAR model.

Syntax:

phaseqsarexport \langle file name \rangle

Operands:

\langle file name \rangle

The name of the file to save the QSAR model to.

phaseqsarexportligandstofile

Exports the selected conformers from the ligands table in the Prepare Ligands or Create Sites steps.

Syntax:

phaseqsarexportligandstofile \langle file name \rangle

Operands:

\langle file name \rangle

The file name to export the selected conformers to.

phaseqsarextendinclude

Extends the inclusion to the given ligand in the Individual QSAR Model panel.

Syntax:

phaseqsarextendinclude \langle row \rangle

Operands:

\langle row \rangle

The row number to include.

phaseqsarextendselectligandrow

Extends the selection to the given row in the ligands table on the Individual QSAR Model panel.

Syntax:

```
phaseqsarextendselectligandrow <row>
```

Operands:

<row>

The row number to extend-select.

phaseqsarimport

Imports a QSAR model from the given file.

Syntax:

```
phaseqsarimport <file>
```

Operands:

<file>

The QSAR model file to import.

phaseqsarinclude

Includes just the given ligand in the Individual QSAR Model panel in the Workspace.

Syntax:

```
phaseqsarinclude <row>
```

Operands:

<row>

The row number to include.

phaseqsarinvertselection

Invert the selection in the ligands table.

Syntax:

phaseqsarinvertselection

phaseqsarplotactivity

Creates a scatter plot of predicted versus actual activities for the Individual QSAR Model.

Syntax:

phaseqsarplotactivity *draw_line*=yes | no *pls_factor*=⟨n⟩
 ⟨all|selected|training|test⟩

Options:

draw_line If set, then the plot will include a line at 45 degrees.
 Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

pls_factor This is which set of PLS factors to plot.
 Valid values: integers
 Default value: **1**
 Minimum: 1

Operands:

⟨all|selected|training|test⟩

Which ligands to plot: all, selected, the training set, or the test set.

phaseqsarpredict

Predicts activities for the given entries using the Individual QSAR Model.

Syntax:

phaseqsarpredict ⟨ESL⟩

Operands:

⟨ESL⟩

The entries to predict activities for.

phaseqsarrandomtraining

Sets the random training set for the Individual QSAR Model panel.

Syntax:

phaseqsarrandomtraining

phaseqsarsearch

Transfers the current hypothesis from the Build QSAR step to the Find Matches panel. Also opens the Find Matches panel. Deprecated - use phasesearchformatches instead.

Syntax:

phaseqsarsearch

phaseqsarselectall

Selects all of the ligands.

Syntax:

phaseqsarselectall

phaseqsarselecthypothesis

Selects the given hypothesis.

Syntax:

phaseqsarselecthypothesis \langle ID \rangle

Operands:

\langle ID \rangle

The ID of the hypothesis to select in the Build QSAR step. This will populate the Alignment table. This function is single select. If the operand is blank, the alignments table will be filled in for the currently selected hypothesis, if any.

phaseqsarselectligandrow

Selects the given row in the ligands table on the Individual QSAR Model panel.

Syntax:

```
phaseqsarselectligandrow <row>
```

Operands:

<row>

The row number to select.

phaseqsarsetactivity

Sets the activity for the given row in the Individual QSAR Model panel.

Syntax:

```
phaseqsarsetactivity activity=<x> <row>
```

Options:

activity The activity value for the given ligand.

Valid values: reals

Default value: **1**

Operands:

<row>

The row number of the ligand to set the activity for.

phaseqsarsettrainingrows

Toggles the Training Set property on or off for the given row in the ligands table on the Individual QSAR Model panel.

Syntax:

phaseqsarsettrainingrows *value*=training | test | none <rows>

Options:

value The value of the training/test column in the ligands table in the Individual QSAR Model panel.

Valid values: training
 test
 none

Default value: **training**

Operands:

<rows>

The row numbers to toggle, two row numbers should be separated by , or if the user wants to specify a range then it can be given like 1-5 separated by a - . eg 1,2,5 or 1,4-7 .

phaseqsarsort

Sort the Individual QSAR Model ligands table based on the data in the specified column

Syntax:

phaseqsarsort <column_name>

Operands:

<column_name>

The name of the column to be sorted.

phaseqsartest

Tests a QSAR model on the current ligand test set.

Syntax:

phaseqsartest

phaseqsartoggleinclude

Toggles the given ligand in the Individual QSAR Model panel into or out of the Workspace.

Syntax:

```
phaseqsartoggleinclude <row>
```

Operands:

<row>

The row number to toggle.

phaseqsartoggleligandrow

Toggles the selection for the given row in the ligands table on the Individual QSAR Model panel.

Syntax:

```
phaseqsartoggleligandrow <row>
```

Operands:

<row>

The row number to toggle-select.

phaseqsartogglerow

Toggles the Training Set property on or off for the given row in the Individual QSAR Model panel.

Syntax:

```
phaseqsartogglerow <row>
```

Operands:

<row>

The row number to toggle.

phaserandomtraining

Sets the random training set.

Syntax:

phaserandomtraining

phaserandomtrainingoptions

Holds the options for the random training set.

Syntax:

phaserandomtrainingoptions *percentage*=⟨n⟩

Options:

percentage The percentage of ligands to use as a training set.

Valid values:	integers
Default value:	50
Minimum:	1
Maximum:	100

phaserefreshfrequencytable

Restores the default set of min/max values to the feature frequencies table in the Find Common Pharamcophores step.

Syntax:

phaserefreshfrequencytable

phaseremoveligands

Removes the given ligands from the table.

Syntax:

phaseremoveligands $\langle \text{ESL} \rangle$

Operands:

$\langle \text{ESL} \rangle$

The entries to remove as ligands.

phaserenameligand

Changes the name for the given ligand.

Syntax:

phaserenameligand *row*= $\langle n \rangle$ $\langle \text{new name} \rangle$

Options:

<i>row</i>	The row number of the ligand to rename.
Valid values:	integers
Default value:	1
Minimum:	1

Operands:

$\langle \text{new name} \rangle$

The new name for the given ligand.

phaserescorehypotheses

Rescores all the hypotheses in the hypotheses table in the Score Hypotheses step.

Syntax:

phaserescorehypotheses

phaserescoreweighting

Sets the weighting factors for rescoring hypotheses in the Score Hypotheses step.

Chapter 5: Commands

Syntax:

```
phaserescoreweighting activity=<x> energy=<x> inactive=<x>  
match=<x> selectivity=<x> site=<x> vector=<x>  
volume=<x>
```

Options:

<i>activity</i>	The reference ligand activity scoring factor. The valid range is 0.0 to 100.0 Valid values: reals Default value: 0 Minimum: 0.0 Maximum: 100.0
<i>energy</i>	The reference ligand relative conformational energy scoring factor. The valid range is 0.0 to 100.0 Valid values: reals Default value: 0 Minimum: 0.0 Maximum: 100.0
<i>inactive</i>	The inactive match scoring factor. The valid range is 0.0 to 100.0 Valid values: reals Default value: 1 Minimum: 0.0 Maximum: 100.0
<i>match</i>	The number of matches scoring factor. The valid range is 1.0 to infinity, but values only slightly above 1.0 might generate huge values for a large number of matched ligands. Valid values: reals Default value: 1 Minimum: 1.0
<i>selectivity</i>	The selectivity scoring factor. The valid range is 0.0 to 100.0 Valid values: reals Default value: 0 Minimum: 0.0 Maximum: 100.0
<i>site</i>	The aligned sites scoring factor. The valid range is 0.0 to 100.0 Valid values: reals Default value: 1 Minimum: 0.0 Maximum: 100.0

<i>vector</i>	The vector scoring factor. The valid range is 0.0 to 100.0
	Valid values: reals
	Default value: 1
	Minimum: 0.0
	Maximum: 100.0
<i>volume</i>	The volume scoring factor. The valid range is 0.0 to 100.0
	Valid values: reals
	Default value: 1
	Minimum: 0.0
	Maximum: 100.0

phaseresetfeatures

Resets the features to the installation defaults.

Syntax:

```
phaseresetfeatures
```

phaseruncreate

Creates the run with the given name.

Syntax:

```
phaseruncreate <run name>
```

Operands:

```
<run name>
```

The name of the new run to create.

phaserundelete

Deletes the current run from the project.

Syntax:

phaserundelete

phaserunopen

Opens the run with the given name.

Syntax:

phaserunopen ⟨run name⟩

Operands:

⟨run name⟩

The name of the run to open.

phaserunrename

Changes the current run's name to the given name.

Syntax:

phaserunrename ⟨run name⟩

Operands:

⟨run name⟩

The name to change the current run's name to.

phaserunsaveas

Saves a copy of the current run under the given name.

Syntax:

phaserunsaveas ⟨run name⟩

Operands:

⟨run name⟩

The name of the run to save as.

phaserunsetseed

Sets the random seed of the current run based on the current Phase options.

Syntax:

```
phaserunsetseed
```

phaserunsetstereo

Sets the stereoisomer behavior of the current run based on the current Phase options.

Syntax:

```
phaserunsetstereo
```

phasesavematching

Saves the advanced matching options for the given hypothesis.

Syntax:

```
phasesavematching <entry ID>
```

Operands:

<entry ID>

The ID of the hypothesis to save the advanced matching options for.

phasescorehypotheses

Launches a Score Hypotheses job.

Syntax:

```
phasescorehypotheses activityweight=<x> aligncutoff=<x>
energyweight=<x> matchreward=<x> maxboxes=<n>
minboxes=<n> selectivitybasedalignment=yes | no
selectivityweight=<x> sitetoppercentage=<n> siteweight=<x>
tolerance_a=<x> tolerance_d=<x> tolerance_h=<x>
tolerance_n=<x> tolerance_p=<x> tolerance_r=<x>
tolerance_x=<x> tolerance_y=<x> tolerance_z=<x>
usetolerances=yes | no vectorlowercutoff=<x>
vectorweight=<x> volumebasedalignment=yes | no
volumeweight=<x>
```

Options:

activityweight

For calculating the Survival score, the weighting factor of the reference ligand activity. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **0**
 Minimum: 0.0
 Maximum: 100.0

aligncutoff

In case of score by site-based alignment, keep the hypotheses with RMSD values below the value specified by this option. The valid range is 0.0001 to infinity

Valid values: reals
 Default value: **1.2**
 Minimum: 0.0001

energyweight

For calculating the Survival score, the weighting factor of the reference ligand relative conformational energy. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **0**
 Minimum: 0.0
 Maximum: 100.0

matchreward

For calculating the Survival score, increase the score by this value, raised to the power of the number of matched ligands. The valid range is 1.0 to infinity, but values only slightly above 1.0 might generate huge values for a large number of matched ligands.

Valid values: reals
 Default value: **1**

	Minimum: 1.0
<i>maxboxes</i>	In case of score by site-based alignment, keep at most the number specified by this option. The valid range is 1 to infinity Valid values: integers Default value: 50 Minimum: 1
<i>minboxes</i>	In case of score by site-based alignment, keep at least the number specified by this option. The valid range is 1 to infinity Valid values: integers Default value: 10 Minimum: 1
<i>selectivitybasedalignment</i>	For score by selectivity-based alignment, this option is true otherwise false. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>selectivityweight</i>	For calculating the Survival score, the weighting factor of the selectivity score. The valid range is 0.0 to 100.0 Valid values: reals Default value: 0 Minimum: 0.0 Maximum: 100.0
<i>sitetoppercentage</i>	In case of score by site-based alignment, keep the hypotheses in the percentage specified by this option. The valid range is 0 to 100 Valid values: integers Default value: 10 Minimum: 0 Maximum: 100
<i>siteweight</i>	For calculating the Survival score, the weighting factor given for the site score. The valid range is 0.0 to 100.0 Valid values: reals Default value: 1 Minimum: 0.0 Maximum: 100.0
<i>tolerance_a</i>	The feature matching tolerance for the hydrogen bond acceptor feature. The valid range is 0.0 to 100.0

Chapter 5: Commands

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 100.0

tolerance_d

The feature matching tolerance for the hydrogen bond donor feature. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 100.0

tolerance_h

The feature matching tolerance for the hydrophobic feature. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **1.5**
Minimum: 0.0
Maximum: 100.0

tolerance_n

The feature matching tolerance for the negative feature. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **0.75**
Minimum: 0.0
Maximum: 100.0

tolerance_p

The feature matching tolerance for the positive feature. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **0.75**
Minimum: 0.0
Maximum: 100.0

tolerance_r

The feature matching tolerance for the aromatic ring feature. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **1.5**
Minimum: 0.0
Maximum: 100.0

tolerance_x

The feature matching tolerance for the custom(X) feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 100.0

tolerance_y

The feature matching tolerance for the custom(Y) feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 100.0

tolerance_z

The feature matching tolerance for the custom(Z) feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 100.0

usetolerances

When using feature matching tolerances, this option is true otherwise false.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

vectorlowercutoff

In case of score by vector-based alignment, keep the hypotheses that score above the value specified by this option. The valid range is 0.0 to 1.0

Valid values: reals
 Default value: **0.5**
 Minimum: 0.0
 Maximum: 1.0

vectorweight

For calculating the Survival score, the weighting factor given for the vector score. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 100.0

volumebasedalignment

For score by volume-based alignment, this option is true otherwise false.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

volumeweight

For calculating the Survival score, the weighting factor given for the volume score. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 100.0

phasescoreinactives

Launches a Score Inactives job.

Syntax:

phasescoreinactives *matchscore*= $\langle x \rangle$

Options:

matchscore

For calculating the adjusted survival score, the weighting factor given for the inactive match score. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 100.0

phasescoreselecthypothesis

Selects the given hypothesis.

Syntax:

phasescoreselecthypothesis *alignall*=yes | no $\langle ID \rangle$

Options:

alignall This option is used to control the behavior when filling in hypothesis alignments table. If this is “true” then alignments are

generated for all ligands, including inactives, that match the hypothesis on at least 3 sites. If this is “false” then the table is filled in with only the ligands from the active set that match all of the hypothesis sites.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

`< ID >`

The ID of the hypothesis to select. This will populate the Alignment table. This function is single select. If the operand is blank, the alignments table will be filled in for the currently selected hypothesis, if any.

phasesearchformatches

Copies the selected hypotheses from Develop Common Pharmacophore panel (in either the Score Hypothesis or Build QSAR model steps) to the Project Table, opens the Find Matches panel, and sets the first new entry as the chosen hypothesis entry in that panel.

Syntax:

phasesearchformatches

phaseselectalltablerows

Selects all rows in the first table in the step.

Syntax:

phaseselectalltablerows

phaseselectconstraintrow

Selects only this row in the constraints table.

Syntax:

phaseselectconstraintrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to select in the constraint table.

phaseselectevrow

Selects the given row in the excluded volumes table.

Syntax:

phaseselectevrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to select in the table.

phaseselectexcludedfeature

Excludes the given feature from Find Common Pharmacophores.

Syntax:

phaseselectexcludedfeature $\langle \text{feature} \rangle$

Operands:

$\langle \text{feature} \rangle$

The feature to select for excluding.

phaseselectextendconstraintrow

Extends the selection to this row in the constraint table.

Syntax:

phaseselectextendconstraintrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to extend the select to.

phaseselecttextendevrow

Extends the selection to this row in the excluded volumes table.

Syntax:

phaseselecttextendevrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to extend the select to.

phaseselecttextendhypothesisrow

Extends the selection to this row in the hypothesis table.

Syntax:

phaseselecttextendhypothesisrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to extend the select to.

phaseselecttextendtablerow

Extends the selection to this row in the table.

Syntax:

phaseselectextendtablerow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to extend the select to.

phaseselecthypothesisrow

Selects the given row in the hypothesis table in the step.

Syntax:

phaseselecthypothesisrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to select in the hypothesis table.

phaseselectonlyevrow

Selects only this row in the excluded volumes table.

Syntax:

phaseselectonlyevrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to select only in the table row.

phaseselectonlyhypothesisrow

Selects only this row in the hypothesis table.

Syntax:

phaseselectonlyhypothesisrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to select only in the hypothesis table.

phaseselectonlytablerow

Selects only this row in the table.

Syntax:

phaseselectonlytablerow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to select only in the table row.

phaseselectsite

Selects the given site in the freestyle hypothesis

Syntax:

phaseselectsite $\langle \text{site index} \rangle$

Operands:

$\langle \text{site index} \rangle$

The index of the site to select.

phaseselecttablerow

Selects the given row in the first table in the step.

Syntax:

phaseselecttablerow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to select in the table.

phaseseparatestereoisomers

Separate stereoisomers for the selected conformers.

Syntax:

phaseseparatestereoisomers

phaseseparatetitles

Separate the selected conformers by title.

Syntax:

phaseseparatetitles

phasesetactiverows

Toggles the Pharm Set property on or off for the given row in the ligands table in Prepare Ligands or Create Sites.

Syntax:

phasesetactiverows *value*=active | inactive | none $\langle \text{rows} \rangle$

Options:

value The value of the Pharm set column in the ligands table in Prepare Ligands or Create Sites.

Valid values: active
 inactive
 none

Default value: **active**

Operands:

$\langle \text{rows} \rangle$

The row numbers to toggle, two row numbers should be separated by , or if the user wants to specify a range then it can be given like 1-5 separated by a - . eg 1,2,5 or 1,4-7 .

phasesetactivity

Sets the activity for the given row in the current step.

Syntax:

phasesetactivity *activity*= $\langle x \rangle$ $\langle \text{row} \rangle$

Options:

activity The activity value for the given ligand.

Valid values: reals

Default value: **1**

Operands:

$\langle \text{row} \rangle$

The row number of the ligand to set the activity for.

phasesetactivityproperty

Sets the activity property.

Syntax:

phasesetactivityproperty $\langle \text{property} \rangle$

Operands:

$\langle \text{property} \rangle$

The property to get the activity values from.

phasesetactivitythresholds

Sets the activity thresholds of the current run based on the current Phase options.

Syntax:

```
phasesetactivitythresholds active_threshold=<text>
                           inactive_threshold=<text>
```

Options:

active_threshold

The value used as the cutoff for assigning a ligand to the active Pharm Set. Can be left blank.

Valid values: text strings

Default value:

inactive_threshold

The value used as the cutoff for assigning a ligand to the inactive Pharm Set. Can be left blank.

Valid values: text strings

Default value:

phasesetalignmentoptions

Sets whether or not to view non-model ligands and whether or not to have a site mask.

Syntax:

```
phasesetalignmentoptions alignnonmodel=yes | no
                           sitemask=<text> <hypothesis>
```

Options:

alignnonmodel

Indicates whether or not to align non-model ligands.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

sitemask

The site mask is a set of 0s and 1s indicating which sites to include in the mask. A 1 in a position means to include that site, and a 0 means to ignore it.

Valid values: text strings

Default value:

Operands:

\langle hypothesis \rangle

The name of the hypothesis to set the options for.

phasesetconstrainttolerance

Sets the tolerance for this row in the constraints table.

Syntax:

phasesetconstrainttolerance *tolerance*= \langle x \rangle \langle row \rangle

Options:

tolerance This is the tolerance to set for the given row.

Valid values: reals

Default value: **1**

Minimum: 0.0

Operands:

\langle row \rangle

The row number to set the tolerance for in the constraint table.

phasesetexcludedvolumes

Sets the value for the given cell.

Syntax:

phasesetexcludedvolumes *column*= \langle n \rangle *row*= \langle n \rangle \langle value \rangle

Options:

column The column number of the cell to change.

Valid values: integers

Default value: **1**

row The row number of the cell to change.

Valid values: integers

Default value: **1**

Operands:

⟨ value ⟩

The value for the given cell.

phasesetfrequency

Sets the minimum or maximum frequency for a given feature.

Syntax:

phasesetfrequency *column*=⟨ n ⟩ *value*=⟨ n ⟩ ⟨ row ⟩

Options:

column The table column to set.

Valid values: integers

Default value: **2**

Minimum: 1

value The minimum or maximum frequency.

Valid values: integers

Default value: **0**

Minimum: 0

Operands:

⟨ row ⟩

The row number of the min / max value to set.

phasesettrainingrows

Toggles the Training Set property on or off for the given row in the Alignments table under Build QSAR step.

Syntax:

phasesettrainingrows *value*=training | test | none ⟨ rows ⟩

Options:

value The value of the training/test column in the alignment table under Build QSAR Model step.

Valid values:	training test none
Default value:	training

Operands:

$\langle \text{rows} \rangle$

The row numbers to toggle, two row numbers should be separated by , or if the user wants to specify a range then it can be given like 1-5 separated by a - . eg 1,2,5 or 1,4-7 .

phaseshowclusters

Shows hypothesis clusters in the Score Hypotheses step.

Syntax:

phaseshowclusters *similarity*= $\langle x \rangle$

Options:

similarity How similar hypotheses need to be to be clustered.

Valid values:	reals
Default value:	0.9

phasesiteoptions

Sets the site options for a Find Pharmacophores job.

Syntax:

phasesiteoptions *match*=all | minimum *minimum*= $\langle n \rangle$
minsites= $\langle n \rangle$ *numsites*= $\langle n \rangle$

Options:

match Indicates whether to match against all or a minimum set of active ligands.

Valid values:	all minimum
Default value:	all

<i>minimum</i>	The minimum number of active ligands which a pattern has to match against. Valid values: integers Default value: 2 Minimum: 2
<i>minsites</i>	The minimum number of sites to match. Valid values: integers Default value: 5 Minimum: 3 Maximum: 7
<i>numsites</i>	The maximum number of sites to match. Valid values: integers Default value: 5 Minimum: 3 Maximum: 7

phasesorttable

Resort the given Phase table based on the data in the specified column

Syntax:

phasesorttable *table*=⟨n⟩ ⟨column_name⟩

Options:

<i>table</i>	The table to set. Valid values: integers Default value: 10
--------------	---

Operands:

⟨column_name⟩

The name of the column to be sorted.

phasestepforward

Moves forward to the next Phase step. Deletes any steps after the current step, then creates the next step, using the data from previous steps.

Syntax:

phasetstepforward

phasetstepgoto

Moves to an existing step in the current project.

Syntax:

phasetstepgoto \langle step name \rangle

Operands:

\langle step name \rangle

The name of the step to switch to.

phasettoggleactivetablerow

Toggles the Active property on or off for the given row in the first table in the step.

Syntax:

phasettoggleactivetablerow \langle row \rangle

Operands:

\langle row \rangle

The row number to toggle.

phasettogglecentroidatom

Adds or removes the given atom from the current centroid list for excluded volumes.

Syntax:

phasetogglecentroidatom \langle atom index \rangle

Operands:

\langle atom index \rangle

The index of the atom to toggle.

phasetoggleconstraintmark

Toggles the mark for this row in the constraints table.

Syntax:

phasetoggleconstraintmark \langle row \rangle

Operands:

\langle row \rangle

The row number to toggle the mark for in the constraint table.

phasetoggleexcludedfeature

Toggles the given feature for including / excluding from Find Common Pharmacophores.

Syntax:

phasetoggleexcludedfeature \langle feature \rangle

Operands:

\langle feature \rangle

The feature to toggle for excluding.

phasetoglerequiredmatch

Toggles the Required Match property for the given row.

Syntax:

phasetoglerequiredmatch $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number of the ligand to toggle the Required Match property for.

phasetoggleselectconstraintrow

Toggles the selection of the given row in the constraints table.

Syntax:

phasetoggleselectconstraintrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to toggle in the constraints table.

phasetoggleselectsite

Toggles the selection of the given site in the freestyle hypothesis

Syntax:

phasetoggleselectsite $\langle \text{site index} \rangle$

Operands:

$\langle \text{site index} \rangle$

The index of the site to toggle.

phasetoggletrainingrow

Toggles the Training Set property on or off for the given row in the Build QSAR step.

Syntax:

phasetoggletrainingrow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to toggle.

phasetranslateexcludedvolumes

Translate in Angstroms the selected excluded volumes.

Syntax:

phasetranslateexcludedvolumes $x=\langle x \rangle$ $y=\langle x \rangle$ $z=\langle x \rangle$

Options:

x Amount in Angstroms to translate in X

Valid values: reals

Default value: **0**

y Amount in Angstroms to translate in Y

Valid values: reals

Default value: **0**

z Amount in Angstroms to translate in Z

Valid values: reals

Default value: **0**

phasetranslatesite

Translate in Angstroms the selected Phase freestyle hypothesis site.

Syntax:

phasetranslatesite $x=\langle x \rangle$ $y=\langle x \rangle$ $z=\langle x \rangle$

Options:

x Amount in Angstroms to translate in X

Valid values: reals

Default value: **0**

<i>y</i>	Amount in Angstroms to translate in Y
	Valid values: reals
	Default value: 0
<i>z</i>	Amount in Angstroms to translate in Z
	Valid values: reals
	Default value: 0

phaseundisplayproperty

Undisplays the given property from the table.

Syntax:

phaseundisplayproperty \langle property name \rangle

Operands:

\langle property name \rangle

The property to undisplay.

phaseungroupactives

Ungroups a ligand group

Syntax:

phaseungroupactives

phaseunmarkfeature

Unmarks the given feature in the Workspace.

Syntax:

phaseunmarkfeature \langle feature \rangle

Operands:

\langle feature \rangle

A single letter (A-Z) indicating the feature to unmark.

phaseunselectevrow

Unselects the given row in the excluded volumes table.

Syntax:

```
phaseunselectevrow <row>
```

Operands:

<row>

The row number to unselect in the table.

phaseunselecthypothesisrow

Unselects the given row in the hypothesis table in the step.

Syntax:

```
phaseunselecthypothesisrow <row>
```

Operands:

<row>

The row number to unselect in the hypothesis table.

phaseunselecttablerow

Unselects the given row in the first table in the step.

Syntax:

```
phaseunselecttablerow <row>
```

Operands:

<row>

The row number to unselect in the table.

picksize

This command allows the user to choose the size of the pick box that is used for picking atoms/bonds/residues etc.

Syntax:

picksize $\langle \text{size} \rangle$

Operands:

$\langle \text{size} \rangle$

This operand actually defines the size of the pick box. Allowed sizes are 7X7, 10X10 and 15X15 which are specified at command line as 7, 10 & 15 respectively.

place

Place the current fragment on screen at the $\langle x \rangle$, $\langle y \rangle$ and $\langle z \rangle$ positions given by the operands.

Syntax:

place $\langle x \rangle \langle y \rangle \langle z \rangle$

Operands:

$\langle x \rangle \langle y \rangle \langle z \rangle$

The operands are three real numbers which are the x, y and z coordinates where the new fragment is to be placed.

plotxyarrangecolumn

Puts all of the plots in a column.

Syntax:

plotxyarrangecolumn

plotxyarrangerow

Displays the plots in a row.

Syntax:

plotxyarrangerow

plotxyarrangetiled

Tiles all of the displayed plots.

Syntax:

plotxyarrangetiled

plotxyaspectratiolock

Displays with/without the aspect ratio maintained for the given plot.

Syntax:

plotxyaspectratiolock \langle plotname \rangle yes|no

Operands:

\langle plotname \rangle yes|no

The name of the plot to operate on. yes to display the plot with aspect ratio locked or no to display it as it is

plotxyaxis

Creates or modifies an XY axis.

Syntax:

plotxyaxis *axis*=x | y *maximum*= \langle text \rangle *minimum*= \langle text \rangle
nummarkers= \langle n \rangle *plot*= \langle text \rangle \langle title \rangle

Options:

axis This option indicates whether this axis is an X axis or a Y axis.

Valid values: x

y

Default value: **x**

maximum This option sets the maximum value for the axis.

Valid values: text strings

	Default value:
<i>minimum</i>	This option sets the minimum value for the axis. Valid values: text strings Default value:
<i>nummarkers</i>	This option sets the number of scale markers for the axis. Valid values: integers Default value: 10 Minimum: 0
<i>plot</i>	This option is the name of the plot containing the axis. Valid values: text strings Default value:

Operands:

$\langle \text{title} \rangle$

The name of the axis to create or modify.

plotxyaxisautorange

Toggles autorange for the given axis on or off.

Syntax:

```
plotxyaxisautorange autorange=yes | no axis=x | y
plot= $\langle \text{text} \rangle$   $\langle \text{axisname} \rangle$ 
```

Options:

<i>autorange</i>	This option indicates whether or not to enable autorange for this axis. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>axis</i>	This option indicates whether this axis is an X axis or a Y axis. Valid values: x y Default value: x
<i>plot</i>	This option is the name of the plot containing the axis. Valid values: text strings Default value:

Operands:

⟨axisname⟩

The name of the axis to set autorange for.

plotxyaxisdelete

Deletes an XY plot axis.

Syntax:

plotxyaxisdelete *axis*=x | y *plot*=⟨text⟩ ⟨axisname⟩

Options:

axis This option indicates whether this axis is an X axis or a Y axis.

Valid values: x

y

Default value: **x**

plot This option is the name of the plot containing the axis.

Valid values: text strings

Default value:

Operands:

⟨axisname⟩

The name of the axis to delete.

plotxyaxisdisplay

Toggles display of the given axis on or off.

Syntax:

plotxyaxisdisplay *autorange*=yes | no *axis*=x | y *plot*=⟨text⟩
⟨axisname⟩

Options:

autorange This option indicates whether to display or undisplay the given axis.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

<i>axis</i>	This option indicates whether this axis is an X axis or a Y axis. Valid values: x y Default value: x
<i>plot</i>	This option is the name of the plot containing the axis. Valid values: text strings Default value:

Operands:

\langle axisname \rangle

The name of the axis to display or undisplay.

plotxyaxisrename

Renames an existing axis.

Syntax:

```
plotxyaxisrename axis=x | y newname= $\langle$ text $\rangle$  plot= $\langle$ text $\rangle$ 
                      $\langle$ title $\rangle$ 
```

Options:

<i>axis</i>	This option indicates whether this axis is an X axis or a Y axis. Valid values: x y Default value: x
<i>newname</i>	This option is the new name for the axis. Valid values: text strings Default value:
<i>plot</i>	This option is the name of the plot containing the axis. Valid values: text strings Default value:

Operands:

\langle title \rangle

The name of the axis to rename.

plotxycaption

Sets the caption for the given plot.

Syntax:

plotxycaption \langle plotname \rangle \langle caption \rangle

Operands:

\langle plotname \rangle \langle caption \rangle

The name of the plot to operate on. The new caption for the plot.

plotxycaptionposition

Sets the position of the caption in the given plot.

Syntax:

plotxycaptionposition \langle plotname \rangle top|bottom

Operands:

\langle plotname \rangle top|bottom

The name of the plot to operate on. top to display the caption at the top of the plot or bottom to display it at the bottom.

plotxycopy

Copies the selected plots.

Syntax:

plotxycopy

plotxydelete

Deletes the selected plots.

Syntax:

plotxydelete

plotxydeleteall

Deletes all XY plots from the project.

Syntax:

plotxydeleteall

plotxydisplay

Displays the selected plots in addition to any currently displayed plots.

Syntax:

plotxydisplay

plotxydisplaycaption

Shows or hides the caption for the given plot.

Syntax:

plotxydisplaycaption \langle plotname \rangle yes|no

Operands:

\langle plotname \rangle yes|no

The name of the plot to operate on. yes to display the caption or no to not display the caption.

plotxydisplayincluded

Sets the Display Included Markers state for the given plot.

Syntax:

plotxydisplayincluded $\langle \text{plotname} \rangle$ yes|no

Operands:

$\langle \text{plotname} \rangle$ yes|no

The name of the plot to operate on. yes to display included markers or no to not display them.

plotxydisplaylegend

Displays or hides the legend for the given plot.

Syntax:

plotxydisplaylegend $\langle \text{plotname} \rangle$ yes|no [legend_psotion]

Operands:

$\langle \text{plotname} \rangle$ yes|no [legend_psotion]

The name of the plot to operate on. yes to display the legend or no to not display it. Positions are upper_right, upper_left, lower_left, lower_right, center_left, center_right, lower_center, upper_center.

plotxydisplayname

Shows or hides the name in the caption for the given plot.

Syntax:

plotxydisplayname $\langle \text{plotname} \rangle$ yes|no

Operands:

$\langle \text{plotname} \rangle$ yes|no

The name of the plot to operate on. yes to display the plot name in the caption or no to not display the name in the caption.

plotxydisplayonly

Displays only the selected plots.

Syntax:

plotxydisplayonly

plotxydisplaypointlabels

Displays or hides the point labels for the given plot.

Syntax:

plotxydisplaypointlabels \langle plotname \rangle yes|no

Operands:

\langle plotname \rangle yes|no

The name of the plot to operate on. yes to display the point labels or no to not display them.

plotxydisplayselected

Displays or hides the selected entry markers for the given plot.

Syntax:

plotxydisplayselected \langle plotname \rangle yes|no

Operands:

\langle plotname \rangle yes|no

The name of the plot to operate on. yes to display the selected entry markers or no to not display them.

plotxyequalizeaxisrange

Displays with/without the axis ranges equalized for the given plot.

Syntax:

plotxyequalizeaxisrange \langle plotname \rangle yes|no

Operands:

\langle plotname \rangle yes|no

The name of the plot to operate on. yes to display the plot with the axis range equalized to be the same.

plotxyhidesidebar

Hides the plotxy side bar.

Syntax:

plotxyhidesidebar

plotxyhidetoolbar

Hides the plotxy tool bar.

Syntax:

plotxyhidetoolbar

plotxylabel

Turns on labels for the given data point.

Syntax:

plotxylabel *entryname*=yes | no *plot*= \langle text \rangle *series*= \langle text \rangle
title=yes | no *xaxis*=yes | no *yaxis*=yes | no \langle entry \rangle

Options:

entryname

If this is set to true then the entry name is displayed as part of the label.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

<i>plot</i>	<p>This option is the name of the plot containing the series and entry.</p> <p>Valid values: text strings</p> <p>Default value:</p>
<i>series</i>	<p>This option is the name of the series containing the entry.</p> <p>Valid values: text strings</p> <p>Default value:</p>
<i>title</i>	<p>If this is set to true then the entry's title is displayed as part of the label.</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: true</p>
<i>xaxis</i>	<p>If this is set to true then the X-axis property is displayed as part of the label.</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: true</p>
<i>yaxis</i>	<p>If this is set to true then the Y-axis property is displayed as part of the label.</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: true</p>

Operands:

⟨entry⟩

The name of the entry to label.

plotxymove

Moves the selected plots.

Syntax:

plotxymove ⟨location⟩

Operands:

⟨location⟩

The location to move the selected plots to. If the location is 0, then the plots are moved to the beginning of the plot list. Otherwise, the plots are moved after the given plot number in the list.

plotxynew

Creates a new XY plot

Syntax:

plotxynew $\langle \text{plotname} \rangle$

Operands:

$\langle \text{plotname} \rangle$

The name of the plot to create.

plotxypan

Moves the selected plots.

Syntax:

plotxypan $\langle x \rangle \langle y \rangle$

Operands:

$\langle x \rangle \langle y \rangle$

This is a percentage of the full range to offset by. This is a percentage of the full range to offset by.

plotxypanplot

Moves the specified plot.

Syntax:

plotxypanplot $\langle \text{plot} \rangle \langle x \rangle \langle y \rangle$

Operands:

$\langle \text{plot} \rangle \langle x \rangle \langle y \rangle$

This is the name of the plot to pan. This is a percentage of the full range to offset by. This is a percentage of the full range to offset by.

plotxyrename

Renames the given plot to the new name.

Syntax:

```
plotxyrename <plotname> <newname>
```

Operands:

<plotname> <newname>

The name of the plot to operate on. The new name for the plot.

plotxyresetview

Resets the zoom and pan for the selected plots.

Syntax:

```
plotxyresetview
```

plotxysaveimage

Capture the current XY plotting window and save it to an image file.

Syntax:

```
plotxysaveimage format=tiff | jpeg <file_name>
```

Options:

<i>format</i>	Specifies the format of the saved image.
Valid values:	tiff jpeg
Default value:	tiff

Operands:

<file_name>

The file where the image will be saved.

plotxyselect

Selects the given plots.

Syntax:

```
plotxyselect <plotname>
```

Operands:

<plotname>

Names of the plots to select, or all, or displayed.

plotxyseries

Creates or modifies an XY data series.

Syntax:

```
plotxyseries color=black | red | green | blue | purple | orange |  
blue_green | light_green | red_purple | yellow | cyan  
color-by-property=<text> color-map=autumn | blues | grays |  
greens | oranges | purples | rainbow | red_blue | red | spring  
| summer | white_black | winter | yellow_green_blue  
line_color=black | red | green | blue | purple | orange |  
blue_green | light_green | red_purple | yellow | cyan  
marker=filled_square | square | filled_circle | circle | cross |  
point | diamond | filled_diamond | none marker-size=<n>  
plot=<text> rsquared=<x> show_color_bar=yes | no  
size-by-property=<text> style=none | solid | dash  
unityaspect=yes | no width=<n> xaxis=<text> xbuckets=<n>  
xproperty=<text> yaxis=<text> ybuckets=<n>  
yproperty=<text> <seriesname>
```

Options:

<i>color</i>	The color of the symbols and lines used on the plot for this data series.
--------------	---

Valid values: black
 red
 green
 blue
 purple
 orange
 blue_green
 light_green
 red_purple
 yellow
 cyan
 Default value: **black**

color_by_property

This option sets the name of the property the points are to be colored by.

Valid values: text strings
 Default value:

color_map The color maps to be used when coloring by properties.

Valid values: autumn
 blues
 grays
 greens
 oranges
 purples
 rainbow
 red_blue
 red
 spring
 summer
 white_black
 winter
 yellow_green_blue
 Default value: **rainbow**

line_color The color for lines to be used on the plot.

Chapter 5: Commands

	Valid values:	black red green blue purple orange blue_green light_green red_purple yellow cyan
	Default value:	black
<i>marker</i>	The marker for this data series.	
	Valid values:	filled_square square filled_circle circle cross point diamond filled_diamond none
	Default value:	filled_circle
<i>marker_size</i>	The marker size for this data series.	
	Valid values:	integers
	Default value:	3
	Minimum:	1
	Maximum:	8
<i>plot</i>	This option is the name of the plot containing the series.	
	Valid values:	text strings
	Default value:	
<i>rsquared</i>	This option sets the correlation coefficient or rsquared value for this data series. This option is for Strike BuildQSAR plots.	
	Valid values:	reals
	Default value:	0
<i>show_color_bar</i>	Settings this option to yes will display a color bar when coloring by a property.	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	true

size_by_property

This option sets the name of the property to be used to size the markers.

Valid values: text strings

Default value:

style

The line style for this data series.

Valid values: none
solid
dash

Default value: **none**

unityaspect

Settings this option to yes causes the plot to display with an equal increment on each axis mapping to the same distance on the screen.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

width

This option sets the width of the line for this data series.

Valid values: integers

Default value: **1**

Minimum: 1

xaxis

This option selects an existing X-axis.

Valid values: text strings

Default value:

xbuckets

This option sets the number of buckets for the X-axis.

Valid values: integers

Default value: **100**

Minimum: 0

xproperty

This option sets the name of the property for the X-axis.

Valid values: text strings

Default value:

yaxis

This option selects an existing Y-axis.

Valid values: text strings

Default value:

ybuckets

This option sets the number of buckets for the Y-axis.

Valid values: integers

Default value: **100**

Minimum: 0

yproperty

This option sets the name of the property for the Y-axis.

Valid values: text strings
Default value:

Operands:

⟨seriesname⟩

The name of the series to create or modify.

plotxyseriesdelete

Deletes the given XY data series.

Syntax:

plotxyseriesdelete *plot*=⟨text⟩ ⟨seriesname⟩

Options:

plot This option is the name of the plot containing the series.
Valid values: text strings
Default value:

Operands:

⟨seriesname⟩

The name of the series to delete.

plotxyseriesdisplay

Displays or undisplay the given series.

Syntax:

plotxyseriesdisplay ⟨plotname⟩ ⟨series⟩ yes|no

Operands:

⟨plotname⟩ ⟨series⟩ yes|no

The name of the plot to operate on. The name of the series. yes to display the given data series or no to not display it.

plotxyseriesdisplay45

Shows or hides a 45-degree reference line for the given series.

Syntax:

```
plotxyseriesdisplay45 <plotname> <series> yes|no [slope]
```

Operands:

```
<plotname> <series> yes|no [slope]
```

The name of the plot to operate on, the series to operate on, and a yes or no to display or hide the 45-degree reference line. The slope is optional but should be 1 or -1.

plotxyseriesdisplayregression

Shows or hides the regression line for the given series.

Syntax:

```
plotxyseriesdisplayregression <plotname> <series> yes|no
```

Operands:

```
<plotname> <series> yes|no
```

The name of the plot to operate on, the series to operate on, and a yes or no to display or hide the regression line.

plotxyseriesrename

Renames an existing series.

Syntax:

```
plotxyseriesrename newname=<text> plot=<text> <title>
```

Options:

newname This option is the new name for the series.

Valid values: text strings

Default value:

plot This option is the name of the plot containing the series.

Valid values: text strings
Default value:

Operands:

⟨title⟩

The name of the series to rename.

plotxyseriesselect

Selects only the entries corresponding to the given series.

Syntax:

plotxyseriesselect *plot*=⟨text⟩ ⟨seriesname⟩

Options:

plot This option is the name of the plot containing the series.

Valid values: text strings
Default value:

Operands:

⟨seriesname⟩

The name of the series to select.

plotxyseriesselectadd

Adds the entries corresponding to the given series to the current selection.

Syntax:

plotxyseriesselectadd *plot*=⟨text⟩ ⟨seriesname⟩

Options:

plot This option is the name of the plot containing the series.

Valid values: text strings
Default value:

Operands:

⟨seriesname⟩

The name of the series to select.

plotxyseriessetaxesequal

Sets the ranges of the X and Y axes to be equal.

Syntax:

plotxyseriessetaxesequal ⟨plotname⟩ ⟨series⟩

Operands:

⟨plotname⟩ ⟨series⟩

The name of the plot to operate on. The name of the series.

plotxyshowsidebar

Displays the plotxy side bar.

Syntax:

plotxyshowsidebar

plotxyshowtoolbar

Displays the plotxy tool bar.

Syntax:

plotxyshowtoolbar

plotxytoggledisplay

Toggles the display of the plot on or off.

Syntax:

`plotxytoggedisplay` \langle plotname \rangle

Operands:

\langle plotname \rangle

The name of the plot to toggle the display of.

plotxyundisplay

Undisplays the selected plots.

Syntax:

`plotxyundisplay`

plotxyunlabel

Turns off the label for the given data point.

Syntax:

`plotxyunlabel` *plot*= \langle text \rangle *series*= \langle text \rangle \langle entry \rangle

Options:

plot This option is the name of the plot containing the series and entry.

Valid values: text strings

Default value:

series This option is the name of the series containing the entry.

Valid values: text strings

Default value:

Operands:

\langle entry \rangle

The name of the entry to unlabel.

plotxyunselect

Unselects the given plots.

Syntax:

```
plotxyunselect <plotname>
```

Operands:

<plotname>

Names of the plots to unselect, or all, or displayed.

plotxyupdate

Updates the selected plots.

Syntax:

```
plotxyupdate
```

plotxyzzoom

Scales the selected plots.

Syntax:

```
plotxyzzoom <x> <y>
```

Operands:

<x> <y>

This is the amount of scaling to use in the horizontal direction. This is the amount of scaling to use in the vertical direction.

plotxyzzoompan

Scales and offsets the given plot.

Syntax:

plotxyzzoompan $\langle \text{plot} \rangle$ $\langle \text{zoomx} \rangle$ $\langle \text{zoomy} \rangle$ $\langle \text{panx} \rangle$ $\langle \text{pany} \rangle$

Operands:

$\langle \text{plot} \rangle$ $\langle \text{zoomx} \rangle$ $\langle \text{zoomy} \rangle$ $\langle \text{panx} \rangle$ $\langle \text{pany} \rangle$

This is the name of the plot to apply the zoom and pan factors to. This is the amount of scaling to use in the horizontal direction. This is the amount of scaling to use in the vertical direction. This is the amount to offset in the horizontal direction as a percentage of the full range. This is the amount to offset in the vertical direction as a percentage of the full range.

potential

Set various options associated with the definition of the potential energy to be used in a MacroModel job.

Syntax:

potential *cele*= $\langle x \rangle$ *charges*=force_field | structure_file *chnd*= $\langle x \rangle$
cutoff=normal | extended | user_defined | none *cvdw*= $\langle x \rangle$
debug= $\langle \text{text} \rangle$ *dielectric*= $\langle x \rangle$ *electrostatics*=field_field |
 constant | distance_dependant *field*=mm2* | mm3* | amber*
 | opls* | amber94 | mmff | mmffs | oplsa | oplsa2005 |
 oplsa2008 *mutualinteractions*=yes | no *solvent*=none | water |
 chcl3 | octanol *substructure*=yes | no

Options:

<i>cele</i>	This option determines what cutoff will be used for the electrostatic part of the energy calculation. Valid values: reals Default value: 12 Minimum: 0.0 Maximum: 99999.0
<i>charges</i>	This option determines where the charges to be used in the energy calculation will come from. Valid values: force_field structure_file Default value: force_field
<i>chnd</i>	This option determines what cutoff will be used for the hydrogen bond part of the energy calculation. Valid values: reals Default value: 4

	Minimum: 0.0
	Maximum: 99999.0
<i>cutoff</i>	This option determines what type of non-bonded cutoff will be used in the energy calculation.
	Valid values: normal extended user_defined none
	Default value: normal
<i>cvdw</i>	This option determines what VDW cutoff will be used in the energy calculation.
	Valid values: reals
	Default value: 7
	Minimum: 0.0
	Maximum: 99999.0
<i>debug</i>	Debug opcodes
	Valid values: text strings
	Default value:
<i>dielectric</i>	The dielectric constant to be used in the electrostatic part of the energy calculation.
	Valid values: reals
	Default value: 1
	Minimum: 0.9999999999
<i>electrostatics</i>	The electrostatic treatment to be used in the energy calculation.
	Valid values: field_field constant distance_dependant
	Default value: constant
<i>field</i>	The force field to be used for the energy calculation.
	Valid values: mm2* mm3* amber* opls* amber94 mmff mmffs oplsaa opls2005 opls2008
	Default value: opls2005

Chapter 5: Commands

mutualinteractions

Controls whether or not to use constrained-atom mutual interactions

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

solvent

The solvent model to be used for the energy calculation

Valid values: none
 water
 chcl3
 octanol

Default value: **water**

substructure

[NOTE: This option is no longer used.]

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

pprep

This keyword is used to set various options associated with running protein preparation.

Syntax:

```
pprep cavity=default | liaison | no_neutralization  
      fix_receptor=yes | no procedure=both | prepare | refine  
      rmsd=<x>
```

Options:

cavity Which preparation procedure to run. Both preparation and refinement by default.

Valid values: default
 liaison
 no_neutralization

Default value: **no_neutralization**

fix_receptor

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

procedure

Which preparation procedure to run. Both preparation and refinement by default.

	Valid values:	both prepare refine
	Default value:	both
<i>rmsd</i>	Minimum RMSD.	
	Valid values:	reals
	Default value:	0.3
	Minimum:	0.0

pprepreceptorligand

Defines a on-screen molecule or entry to be treated as the ligand for a Pprep calculation

Syntax:

pprepreceptorligand \langle molecule_number | entry_id \rangle

Operands:

\langle molecule_number | entry_id \rangle

The molecule number or entry ID to be included as the ligand.

pprepwrite

Write the files required for protein structure preparation.

Syntax:

pprepwrite

prefer

Set the global preferences for handling markers (derived graphical objects). Specifically, specify under which conditions markers will be deleted.

Syntax:


```
prefer 2dbondlinewidth=<n> 2dbondspacing=<n>  
      2delementlabelmargin=<x> 2dfontsize=<x>  
      2dhashspacing=<n> 2dlabelallcarbons=yes | no  
      2dmaxatoms=<n> 2dmaxscalefactor=<x> 2dscalefactor=<x>  
      2dshowimplicithydrogens=yes | no 2dusecolor=yes | no  
      adjustmentsfontbold=yes | no adjustmentsfontitalics=yes | no  
      adjustmentsfontname=<text> adjustmentsfontsize=<x>  
      annotationstoolbarshow=yes | no autofit=never | singleentry |  
always avoidhardlinks=never | afs | always | auto  
beep=yes | no buildbackbonesubjobs=<n>  
buildtoolbarshow=yes | no buildtoolbarstyle=icononly |  
textonly | textbeside | textunder changedatanames=yes | no  
checkproteinprepared=yes | no cleanupusingmmap=yes | no  
clipdistance=<x> clipincrement=<x> cmdhistory=unlimited |  
limited cmdhistorylimit=<n> cmdlogfile=<text>  
collapsed=<n> colorhypothesesbyentry=yes | no  
commandcompletion=yes | no commandinputshow=yes | no  
deletemarkers=mismatch | missing displayangleprecision=<n>  
displayatomstoolbarshow=yes | no  
displayatomstoolbarstyle=icononly | textonly | textbeside |  
textunder displaydihedralprecision=<n>  
displaydistanceprecision=<n> displayeditwarning=yes | no  
displayprecision=<n> displaywithinincludesnonpolarh=yes | no  
dockingpanels=yes | no docklocation=mainwindow | floating  
drawmeasurementborder=yes | no edittoolbarshow=yes | no  
edittoolbarstyle=icononly | textonly | textbeside | textunder  
enablejobdebugoutput=yes | no entryedittoolbarstyle=icononly |  
textonly | textbeside | textunder entryfeedbackshow=yes | no  
entryfeedbackshownames=yes | no entrytoolbarstyle=icononly |  
textonly | textbeside | textunder eplayertoolbarshow=yes | no  
eplayertoolbarstyle=icononly | textonly | textbeside | textunder  
feedbackproperties=<text> feedbackproperty=<text>  
feedbackshow=yes | no filerefreshinterval=<n>  
findtoolbarshow=yes | no fitenhance=yes | no  
fitenhancedensity=<x> fitenhancefar=<x> fitenhancenear=<x>  
fitgrow=yes | no fitoffset=<x> fluorine=yes | no  
fragmentstoolbarshow=yes | no fragmentstoolbarstyle=icononly  
| textonly | textbeside | textunder  
geometrycleanupmethod=builtin | uff  
graphicsmemorycachelimit=<n>  
hidewstransformationcmds=yes | no improvelighting=yes | no  
incorporatejobs>manual | prompt | auto initworkdir=<text>  
jobnamefilter=<text> jobnameretain=yes | no  
jobpanelupdateinterval=<n> jobstatusfeedbackshow=yes | no  
jobstatusupdateinterval=<n> jobstealth=yes | no  
kbrotationincrement=<n> kbtransincrement=<n>  
keepcombiglidejobfiles=yes | no  
keepcontinuationrefinementproperties=yes | no  
keepphasejobfiles=yes | no keepprimejobfiles=yes | no
```

Chapter 5: Commands

Options:

2dbondlinewidth

The width for the lines used to draw the bonds.

Valid values: integers
Default value: **1**
Minimum: 1
Maximum: 20

2dbondspacing

The spacing between multiple bonds in the 2D structure drawing.

Valid values: integers
Default value: **2**
Minimum: 1
Maximum: 20

2delementlabelmargin

The separation between the label and structure. Values less than 1.0 will decrease the space between the label and the structures, values greater than 1.0 will increase it. Note there may be clipping effects for values very different from 1.0.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 2.0

2dfontsize

The font size used when drawing the 2D structure

Valid values: reals
Default value: **30**
Minimum: 10.0
Maximum: 100.0

2dhashspacing

The spacing used for hashed lines in the 2D structures

Valid values: integers
Default value: **40**
Minimum: 1
Maximum: 100

2dlabelallcarbons

Whether all carbons are labeled when drawing the 2D structures.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

2dmaxatoms

This option is used as the maximum number of atoms we will ever attempt to draw a 2D structure for.

Valid values: integers
 Default value: **200**
 Minimum: 10
 Maximum: 5000

2dmaxscalefactor

A scaling factor representing the maximum amount small structures will be scaled up.

Valid values: reals
 Default value: **0.15**
 Minimum: 0.0
 Maximum: 500.0

2dscalefactor

A scaling factor used when drawing the 2D structure

Valid values: reals
 Default value: **10**
 Minimum: 1.0
 Maximum: 500.0

2dshowimplicithydrogens

Whether all hydrogens are shown during the structure drawing.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

2dusecolor

Whether to use color when performing the 2D structure drawing.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

adjustmentsfontbold

A bool which determines whether bold font is to be used for adjustments markers

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

adjustmentsfontitalics

A bool which determines whether italics font is to be used for adjustments markers

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

adjustmentsfontname

A string which determines the font to be used for adjustments markers

Valid values: text strings
 Default value:

Chapter 5: Commands

adjustmentsfontsize

A double which determines the font size to be used for adjustments markers

Valid values: reals

Default value: **14**

Minimum: 0.0

annotationstoolbarshow

Whether or not to show the Annotations toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

autofit

Specifies whether workspace entries fit to Workspace automatically or not. Valid values are “never”, “singleentry”, and “always”. These cause Maestro never fit to Workspace, fit when only one entry in workspace and always fit to Workspace respectively.

Valid values: never
singleentry
always

Default value: **singleentry**

avoidhardlinks

Specifies conditions for which files should be copied rather than hard-linked. Valid values are “never”, “afs”, “auto”, and “always”. For “afs”, copy the file if the original file appears to be on an AFS mount (having /afs/ in the file path), since hard links are not permitted between files in different directories under AFS (“afs” setting is deprecated, use “auto”). For “auto”, try to create a link, if link fails then try to copy the file.

Valid values: never
afs
always
auto

Default value: **auto**

beep

Whether or not the system beeps for picking feedback.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

buildbackbonesubjobs

This is the number of templates to run at one time in Build Backbone.

Valid values: integers

Default value: **1**

Minimum: 1

buildtoolbarshow

Whether or not to show the Build toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

buildtoolbarstyle

Determines the appearance of the Build toolbar

Valid values: icononly
 textonly
 textbeside
 textunder

Default value: **icononly**

changedatanames

If set to true, then when a non-standard property is renamed, change the internal data name for the property to match. Otherwise, the data name is not changed when the external (user-visible) property name is changed. A non-standard property is one that entries are not required to have, is added to project only after project is initialized and can be removed from the project.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

checkproteinprepared

Whether or not protein structure, to be used for job input, should be checked for having been prepared using Protein Preparation Wizard.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

cleanupusingmmap

On value uses memory mapped IO (mmap) when running project_cleanup. Off uses traditional read system calls. Performance can vary depending on your file system and on your nfs configuration (if you are using an nfs-mounted disk). Using mmap should give better performance on traditional file systems and where the maestro-project is stored on file systems conforming to POSIX.1b. In some cases such as GPFS you may get better performance by turning mmap off.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

clipdistance

Adjust clipping planes to within the distance when focusing on a substructure.

Chapter 5: Commands

Valid values: reals
Default value: **4**
Minimum: 0.0

clipincrement

Move the clipping planes by this percentage (of slab width) when moving forward, back, together or apart.

Valid values: reals
Default value: **3**
Minimum: 0.0

cmdhistory

Specifies whether to store limited or unlimited commands in command history. The default option is 'unlimited' and if user specifies 'limited' then 'cmdhistorylimit' option will be used.

Valid values: unlimited
 limited
Default value: **limited**

cmdhistorylimit

Specifies number of commands to be stored in command history. This option will be used when 'cmdhistory=limited'.

Valid values: integers
Default value: **500**
Minimum: 0

cmdlogfilename

If MM_OPTION_PREF_WRITE_CMDLOG_ON_EXIT is true, then commands will be written to this file.

Valid values: text strings
Default value: **maestrolog.cmd**

collapsed

Whether the status bar is collapsed or not.

Valid values: integers
Default value: **0**

colorhypothesesbyentry

If this option is set to true, then Maestro will color entry-based hypotheses by the entry coloring scheme.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

commandcompletion

If this option is set to true then command completion will be performed in the command input area

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

commandinputshow

Whether the command input area is show or hide.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

deletemarkers

This option sets the condition for marker deletion. Valid values are “mismatch” or “missing”.

Valid values: mismatch

missing

Default value: **mismatch**

displayangleprecision

Default display precision for angle measurement

Valid values: integers

Default value: **1**

Minimum: 1

Maximum: 6

displayatomstoolbarshow

Whether or not to show the Display Atoms toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

displayatomstoolbarstyle

Determines the appearance of the Display Atoms toolbar

Valid values: icononly

textonly

textbeside

textunder

Default value: **icononly**

displaydihedralprecision

Default display precision for dihedral measurement

Valid values: integers

Default value: **1**

Minimum: 1

Maximum: 6

displaydistanceprecision

Default display precision for distance measurement

Valid values: integers

Default value: **3**

Minimum: 1

Maximum: 6

Chapter 5: Commands

displayeditwarning

If set to true, then a warning message is displayed when user edits value of an entry property in PT.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

displayprecision

Default display precision for PT real type columns

Valid values: integers

Default value: **3**

Minimum: 0

Maximum: 15

displaywithinincludesnonpolarh

A bool which determines whether workspace sequence viewer has to be aligned by residue number or not

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

dockingpanels

If this option is set to true then most non-application panels will be dockable in the main window.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

docklocation

If set to 0 (mainwindow), dockable panels will be docked within the main window. If it is set to 1 (floating), they will be docked into a separate top-level window.

Valid values: mainwindow
floating

Default value: **mainwindow**

drawmeasurementborder

If this flag is on, a rectangle will be drawn in the background color around measurement text. If this is off, then the text will be drawn with no border.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

edittoolbarshow

Whether or not to show the Edit toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

edittoolbarstyle

Determines the appearance of the Edit toolbar

Valid values: icononly
 textonly
 textbeside
 textunder
 Default value: **textunder**

enablejobdebugoutput

This option if enabled would print job control diagnostic output to the job log file.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

entryedittoolbarstyle

Determines the appearance of the PT Entry Edit toolbar

Valid values: icononly
 textonly
 textbeside
 textunder
 Default value: **textunder**

entryfeedbackshow

Whether the workspace entry feedback is activated.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

entryfeedbackshownames

Whether the workspace entry feedback includes the name of the property along with its value.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

entrytoolbarstyle

Determines the appearance of the PT Entry toolbar

Valid values: icononly
 textonly
 textbeside
 textunder
 Default value: **textunder**

eplayertoolbarshow

Whether or not to show the EPlayer toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

eplayertoolbarstyle

Determines the appearance of the ePlayer toolbar

Chapter 5: Commands

Valid values: icononly
 textonly
 textbeside
 textunder
Default value: **icononly**

feedbackproperties

These are the entry properties displayed for a single entry in workspace. The property names should be m2io data names with no spaces, separated by spaces in the option string.

Valid values: text strings
Default value: **s_m_title s_pdb_PDB_ID**

feedbackproperty

This is the entry property displayed at the end of the atom feedback string in workspace.

Valid values: text strings
Default value: **Title**

feedbackshow

Whether the workspace feedback is activated.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

filerefreshinterval

Refresh interval (in secs) for file dialogs

Valid values: integers
Default value: **10**
Minimum: 0

findtoolbarshow

Whether or not to show the Find toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

fitenhance If set to true, automatically enhance depth cues when we fit to Workspace.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

fitenhancedensity

The extra density added to fog density for enhance depth cues when we fit to Workspace.

Valid values: reals
Default value: **0.5**
Minimum: 0.0

- fitenhancefar* The far position of structures for enhance depth cues when we fit to Workspace.
 Valid values: reals
 Default value: **1**
- fitenhancenear* The near position of structures for enhance depth cues when we fit to Workspace.
 Valid values: reals
 Default value: **-1**
- fitgrow* If this option is set to true, then Maestro will fit to Workspace after every grow operation in the builder.
 Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**
- fitoffset* Offset used when fitting to structures to screen. A non-zero offset will provide more of a buffer around the structure. This value is in Angstroms.
 Valid values: reals
 Default value: **2**
 Minimum: 0.0
- fluorine* This option determines whether Treat organofluorines as H-Bond acceptors.
 Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**
- fragmentstoolbarshow* Whether or not to show the Fragments toolbar.
 Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**
- fragmentstoolbarstyle* Determines the appearance of the Fragments toolbar
 Valid values: icononly
 textonly
 textbeside
 textunder
 Default value: **icononly**
- geometrycleanupmethod* Determines how geometry cleanup is performed - using built-in or UFF method.
 Valid values: builtin
 uff

Default value: **builtin**

graphicsmemorycachelimit

Specify number of atoms for which graphics library should cache the memory. If entries are included and later excluded, then graphics library would not release memory until this threshold is reached.

Valid values: integers

Default value: **50000**

hidewstransformationcmds

To hide or show translate, rotate, zoom and transformation commands in command history

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

improvelighting

If this flag is on, clipping planes settings and color palette will be set closer to PyMOL defaults.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

incorporatejobs

Specifies the conditions under which completed jobs that are incorporatable in the current project are incorporated. This setting does not affect the automatic incorporation of monitored jobs in Maestro. It affects the incorporation of jobs that are found to be incorporatable (incorporatable jobs) when the status of all jobs listed in the monitoring panel is updated. Valid values for this setting are “manual”, “prompt”, and “auto”. For the “manual” setting, incorporatable jobs are incorporated only if they are monitored. For the “prompt” setting, incorporatable jobs are incorporated if user approval is given. If approval is not given to incorporate these jobs, then they can only be incorporated in the current Maestro session if they are monitored. Otherwise, they can be incorporated from a different Maestro session. For the “auto” setting, all incorporatable jobs are incorporated automatically, except for those that the user chose not to incorporate under the “prompt” setting

Valid values: manual
prompt
auto

Default value: **prompt**

initworkdir

The Maestro input/output directory is used for starting jobs and for other file input and output. The i/o directory is always

displayed in the title bar of the Main Application Window. The i/o dir can be set to a number of values which affect Maestro's behavior when reading and writing files and when running jobs. You can specify "currentdir" ("startdir" is a synonym) to set the i/o directory to use Maestro's current working directory (cwd), "project" to change to the project directory, "projectparent" to change to the parent directory which contains the project directory, or "projectjobs" to change to the "jobs" directory within the project directory. Otherwise, specify a directory path to be used. Maestro has a current working directory (cwd), similar to what a Unix shell has. When "currentdir" is in effect job files will be placed in the cwd. This is located at the very top of the window next to the window's border. The cwd is changed whenever a cd (changedirectory) command

Valid values: text strings
 Default value: **currentdir**

jobnamefilter

The filter for identifying Maestro job names.

Valid values: text strings
 Default value: ***.inp *.in *.com**

jobnameretain

This option is used to record if last used job name should be the default job name for the next job.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

jobpanelupdateinterval

The time interval, in seconds, of how often status gets updated for unfinished jobs listed in the Monitor panel.

Valid values: integers
 Default value: **30**

jobstatusfeedbackshow

Whether the workspace job status feedback is activated.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

jobstatusupdateinterval

The time interval of how often job status gets updated in the Workspace status bar.

Valid values: integers
 Default value: **30**

jobstealth

Whether or not job monitoring and incorporation is expected to proceed with minimal disruption to the Workspace and entry selection.

Chapter 5: Commands

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

kbrotationincrement

Uses this value when performing workspace rotation with keyboard.

Valid values: integers
Default value: **5**

kbtransincrement

Uses this value when performing workspace translation with keyboard.

Valid values: integers
Default value: **5**

keepcombiglidejobfiles

If this is set to true, then CombiGlide job files will be preserved after a job completes. Otherwise, the job files will be removed when each wizard-based job completes.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

keepentrydisassociateproperties

If set to true, then copy properties to new entries created by entrydisassociate commands. Normally, individual molecules or chains would not have the same properties as the original structure, but this provides the option to preserve the original structure properties in the new entries.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

keepphasejobfiles

If this is set to true, then Phase job files will be preserved after a job completes. Otherwise, the job files will be removed when each wizard-based job completes.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

keepprimejobfiles

If this is set to true, then Prime job files will be preserved after a job completes. Otherwise, the job files will be removed when each wizard-based job completes.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

labelstoolbarstyle

Determines the appearance of the Labels toolbar

Valid values: icononly
 textonly
 textbeside
 textunder
 Default value: **textunder**

lastphasedb

This records the last Phase 3D Database which was opened.

Valid values: text strings
 Default value:

lastphasematchfile

This records the last match file which was created from Phase.

Valid values: text strings
 Default value:

lastproject This records the last project which was opened.

Valid values: text strings
 Default value:

limitanimationdistance

Whether or not view animation should be suppressed when the distance traversed would be too large.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

macromodelsuffix

The default suffix or extension to be used for exported Macro-Model files. The suffix will automatically be appended to a file name if it has no suffix.

Valid values: text strings
 Default value: **dat**

maestrocompressedsuffix

The default suffix or extension to be used for exported compressed Maestro files. The suffix will automatically be appended to a file name if it has no suffix.

Valid values: text strings
 Default value: **maegz**

maestrosuffix

The default suffix or extension to be used for exported Maestro files. The suffix will automatically be appended to a file name if it has no suffix.

Valid values: text strings
 Default value: **mae**

Chapter 5: Commands

mainmetatoolbarshow

Whether or not to show the Meta toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

mainprojecttoolbarshow

Whether or not to show the Project toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

maxnumcons

Default maximum number of constraints allowed

Valid values: integers

Default value: **5000**

Minimum: 1

Maximum: 20000

measurementsfontbold

A bool which determines whether bold font is to be used for measurements markers

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

measurementsfontitalics

An integer which determines whether italics font is to be used for measurements markers

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

measurementsfontname

A string which determines the font to be used for measurements markers

Valid values: text strings

Default value:

measurementsfontsize

A double which determines the font size to be used for measurements markers

Valid values: reals

Default value: **14**

Minimum: 0.0

mol2suffix

The default suffix or extension to be used for exported Mol2 files. The suffix will automatically be appended to a file name if it has no suffix.

Valid values: text strings

Default value: **mol2**

monitorjobs

Whether or not jobs launched from maestro are automatically monitored.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

mouseactionset

Sets the current mouse action set.

Valid values: 3buttonsandscrollwheel

3buttonsonly

2buttonsandscrollwheel

2buttonsonly

pymolmode

trackpad

Default value: **3buttonsandscrollwheel**

mruprojectlistlimit

This is the maximum number of most recently used project to be shown in the Project menu.

Valid values: integers

Default value: **5**

Minimum: 0

Maximum: 10

nonlinearanimation

Whether or not nonlinear animation is used for slow-in action near the starting pose and slow-out action near the ending pose.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

nonlineardegree

The degree of a NURBS curve which is used to control the non-linear animation. matching tolerance for the custom(Z) feature. The valid range is 0.0 to 100.0

Valid values: integers

Default value: **3**

Minimum: 2

Maximum: 4

nonlineareffect

The factor that controls one of slow-in and slow-out effects by changing the shape of a NURBS curve, which is used to control the nonlinear animation.

Valid values: reals

Default value: **10**

Minimum: 0.0

Maximum: 50.0

Chapter 5: Commands

obfuscatepaths

If this flag is set, then the postmortem utility will obfuscate paths automatically.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

openlastproject

If set to true, then when Maestro starts it will re-open the last project that was open when Maestro was closed.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

pdbsuffix

The default suffix or extension to be used for exported PDB files. The suffix will automatically be appended to a file name if it has no suffix.

Valid values: text strings

Default value: **pdb**

phasedefaultfeaturedefinitions

This allows the user to override the default `pharma_feature.ini` file in Phase application panels.

Valid values: text strings

Default value:

phasefeedbackproperties

These are the Phase properties displayed for a single Phase ligand in the workspace. The property names should be m2io data names with no spaces, separated by spaces in the option string.

Valid values: text strings

Default value:

phasetolerance_a

The feature matching tolerance for the hydrogen bond acceptor feature. The valid range is 0.0 to 100.0

Valid values: reals

Default value: **1**

Minimum: 0.0

Maximum: 100.0

phasetolerance_d

The feature matching tolerance for the hydrogen bond donor feature. The valid range is 0.0 to 100.0

Valid values: reals

Default value: **1**

Minimum: 0.0

Maximum: 100.0

phasetolerance_h

The feature matching tolerance for the hydrophobic feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1.5**
 Minimum: 0.0
 Maximum: 100.0

phasetolerance_n

The feature matching tolerance for the negative feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **0.75**
 Minimum: 0.0
 Maximum: 100.0

phasetolerance_p

The feature matching tolerance for the positive feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **0.75**
 Minimum: 0.0
 Maximum: 100.0

phasetolerance_r

The feature matching tolerance for the aromatic ring feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1.5**
 Minimum: 0.0
 Maximum: 100.0

phasetolerance_x

The feature matching tolerance for the custom(X) feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 100.0

phasetolerance_y

The feature matching tolerance for the custom(Y) feature. The valid range is 0.0 to 100.0

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 100.0

Chapter 5: Commands

phasetolerance_z

The feature matching tolerance for the custom(Z) feature. The valid range is 0.0 to 100.0

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 100.0

picktomovetiles

If this flag is on, then tile picking will be turned on when tile mode is entered. This allows tiles to be rearranged without changing the content of the Workspace. If this flag is off, tile picking will be left off when tile mode is entered.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

primextoolbarshow

Whether or not to show the PrimeX toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

primextoolbarstyle

Determines the appearance of the PrimeX toolbar

Valid values: icononly
 textonly
 textbeside
 textunder
Default value: **textunder**

prioritizedredraw

If this flag is on, then redraws will ignore pending events. If this flag is off, redraws will be skipped if there are pending events.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

projectbackuplimit

This is deprecated and calls: projectbackup maxbackups=<int>. Please use that and 'projectbackup instead.

Valid values: integers
Default value: **0**
Minimum: 0

projectdbsync

Holds the project database synchronous pragma flag value. This synchronous pragma flag will be set to the project database when we open an existing project or create a new project.

Valid values: off
 normal
 full
 Default value: **off**

projectmemorycachelimit

Specifies Maestro project memory cache limit in MBs. If memory usage exceeds from given cache limit, then Maestro dumps entry related data (structure information, volume, and surface) into the disk.

Valid values: integers
 Default value: **500**

projectsave

Specifies information to be saved to disk when project is closed. The “small” option saves only the compressed opening state. Because it takes longer to open and close projects when the current state must be expanded or deleted, this option is not recommended unless you are low on disk space. The “medium” option saves only the expanded current project state. This option requires more disk space than the “small” option, but is the fastest because it does not do compression, or save an opening state. There will generally not be an opening state saved for a project when the “medium” option is used, so the ability to revert to the opening state while the project is open is lost. The “large” option saves the compressed opening state and the expanded current project state. This option requires more disk space, but provides faster and easier access (compared to the “small” option) and also the safeguard of redundant project data. Because of the time required to save the opening state,

Valid values: small
 medium
 large
 Default value: **medium**

projectsuffix

The default suffix or extension to be used for project directories. When a project selector is used in Maestro to choose a project, the filter is automatically set to match this suffix. The suffix is automatically appended, if missing, to the returned project path when there is potential for creating a new project (projectnew, projectrename, or projectcopy). However, when a project command is issued without using a project selector (e.g. in a macro, script, or the command input area), the suffix is not automatically appended. Also, in cases where an existing project must be opened (projectopen, projectmerge) the suffix is not appended.

Valid values: text strings
Default value: **.prj**

projectsync

Specifies whether workspace changes to project entries are automatically saved or not. Valid values are “auto”, “prompt”, and “manual”. These cause Maestro to save changes automatically, prompt to save changes, or save changes only when explicitly directed by the user.

Valid values: auto
prompt
manual
Default value: **auto**

projecttabletoolbarstyle

Determines the appearance of the Project Table toolbar

Valid values: icononly
textonly
textbeside
textunder
Default value: **textunder**

projecttoolbarshow

Whether or not to show the Project Table toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

projecttoolbarstyle

Determines the appearance of the Project toolbar

Valid values: icononly
textonly
textbeside
textunder
Default value: **textunder**

propertyedittoolbarstyle

Determines the appearance of the PT Property Edit toolbar

Valid values: icononly
textonly
textbeside
textunder
Default value: **textunder**

propertyselectorsubset

When a new property selector instance is created, it should query this preference. If the subset is available, it should be used to initialize the Subset option menu for that instance.

Valid values: text strings
 Default value: **allprimaryproperties**

propertytoolbarstyle

Determines the appearance of the PT Property toolbar

Valid values: icononly
 textonly
 textbeside
 textunder
 Default value: **textunder**

pspinterface

Default user interface for Prime Structure Prediction. Can be user , msv or wizard .

Valid values: text strings
 Default value: **user**

pspsequenceviewerfontsize

Prime sequence viewer font size.

Valid values: integers
 Default value: **10**
 Minimum: 6
 Maximum: 48

pspsequenceviewerwrap

Sequence wrapping in Prime sequence viewer.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

ptcellcolor

A string which is a 6-hexadecimal-digit RGB string for the default color of cells in the Project Table.

Valid values: text strings
 Default value: **#FFFFFF**

ptcolorcolumnheader

If this flag is on, the column header cells in the Project Table will be colored with a different color.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

ptcolornonstandardfixedarea

If this flag is on, the cells in the nonstandard fixed area of the Project Table will be colored with a different color.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

ptcolorrownumber

If this flag is on, the row number cells in the Project Table will be colored with a different color.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

ptcolorstandardfixedarea

If this flag is on, the cells in the standard fixed area of the Project Table will be colored with a different color.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

ptcolumnheadercolor

A string which is a 6-hexadecimal-digit RGB string for the color of column header cells in the Project Table.

Valid values: text strings
Default value: **#87CEFA**

ptentryedittoolbarshow

Whether or not to show the Entry Edit toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

ptentryselectcolor

A string which is a 6-hexadecimal-digit RGB string for the color of cells for selected entries in the Project Table.

Valid values: text strings
Default value: **#FFFFC0**

ptentrytoolbarshow

Whether or not to show the Entry toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

ptfitcolumnbuffer

If set to 0, then no buffer will be added when fitting column width to data, else the given pixels will be added as buffer.

Valid values: integers
Default value: **0**

ptgroupfullselectcolor

A string which is a 6-hexadecimal-digit RGB string for the color of cells for groups with all entries selected in the Project Table.

Valid values: text strings
Default value: **#CDBA96**

ptgrouppartialselectcolor

A string which is a 6-hexadecimal-digit RGB string for the color of cells for groups with some entries selected in the Project Table.

Valid values: text strings
Default value: **#FFE7BA**

pthighlightalternatingrows

If this flag is on, the cells in the Project Table will be colored to highlight alternating rows.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

ptnonstandardfixedareacolor

A string which is a 6-hexadecimal-digit RGB string for the color of cells for nonstandard properties in the fixed area in the Project Table.

Valid values: text strings

Default value: **#FF8CFF**

ptnumberalignment

If set to 1, then numeric data in project table will be center aligned. Setting 0 will left align the numbers and 2 will right align.

Valid values: left
 center
 right

Default value: **right**

ptpropertyedittoolbarshow

Whether or not to show the Property Edit toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

ptpropertytoolbarshow

Whether or not to show the Property toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

ptrownumbercolor

A string which is a 6-hexadecimal-digit RGB string for the color of row number cells in the Project Table.

Valid values: text strings

Default value: **#87CEFA**

ptstandardfixedareacolor

A string which is a 6-hexadecimal-digit RGB string for the color of cells for standard properties in the fixed area in the Project Table.

Valid values: text strings

Default value: **#9DB3CE**

Chapter 5: Commands

pttextalignment

If set to 1, then string property data in project table will be center aligned. Setting 0 will left align the data and 2 will right align the data.

Valid values: left
 center
 right
Default value: **left**

ptvscrollbar

If set to 0, then vertical scrollbar will be displayed on the left side of project table, else it will be displayed on the right side of project table.

Valid values: left
 right
Default value: **left**

refinebackbonesubjobs

This is the number of refinement jobs to run at one time in Refine Backbone.

Valid values: integers
Default value: **1**
Minimum: 1

renamecrystalmates

If this option is set to true, Maestro will assign unique names to generates crystal mate chains.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

representationtoolbarshow

Whether or not to show the Representaion toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

representationtoolbarstyle

Determines the appearance of the Representation toolbar

Valid values: icononly
 textonly
 textbeside
 textunder
Default value: **icononly**

resetclipprefs

When true, resets clipping preferences (adjustclip and zoom_topview) to their default values, and sets itself to false.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

resetfitmarginprefs

When true, resets the fit margin to the default value.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

rightclicktimerperiod

The time delay for workspace menu to appear in milliseconds.

Valid values: integers
 Default value: **300**

rotamerfitselectedresidue

If this option is set to true then selected residue will be fit to Workspace.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

savedviewstoolbarstyle

Determines the appearance of the Saved Views toolbar

Valid values: icononly
 textonly
 textbeside
 textunder
 Default value: **icononly**

savelayoutonexit

If this option is set to true, the panel layout will be saved when you exit Maestro.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

savescratchproject

If set to true, then pop up 'Save Scratch Project' dialog to Save/Discard the scratch project. Else the scratch project will be discarded.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

savetempjobdir

This option if enabled would not remove temporary job dir when job completes.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

scriptupdatesws

If this option is set to true, force redraw of main window if redraw was requested by running a command in a script.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

sdcompressedsuffix

The default suffix or extension to be used for exported compressed SD files. The suffix will automatically be appended to a file name if it has no suffix.

Valid values: text strings
 Default value: **sdfgz**

sdsuffix

The default suffix or extension to be used for exported SD files. The suffix will automatically be appended to a file name if it has no suffix.

Valid values: text strings
 Default value: **sdf**

selectexcludeselected

If this option is set, then expanding the selection using the Workspace Selection menu will exclude the current selection.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

selectexpandresidues

If this option is set, then expanding the selection using the Workspace Selection menu will expand to residues rather than just atoms.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

sequenceviewerfontsize

The font size for the Workspace sequence viewer

Valid values: integers
 Default value: **10**

sequenceviewerproximity

Cutoff distance for proximity coloring in the sequence viewer

Valid values: reals
 Default value: **4**
 Minimum: 0.0

sequenceviewershow

Whether the sequence viewer is show or hide.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

show2dstructures

If set to true, then 2D column is shown in the PT for the new/scratch project otherwise its hidden by default.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

show2dtooltip

If set to true, then hovering the mouse over Project table 'Row' column shows the 2D structure image as the column tooltip otherwise no tooltip is shown for the 'Row' column.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

showdialogatcursor

If this option is set to true then the dialogs will be shown at mouse cursor location. If it is false then the dialogs will be displayed at their last shown position.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

showdynamictoolbar

Display the dynamic toolbar

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

showentryname

If set to true, then Entry Name property column is shown in the PT otherwise its hidden by default.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

showjobstatus

Whether or not job status is shown in the Workspace status bar.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

showmonitorpanel

If this is set to true, and monitorjobs is true, then Monitor panel will be displayed after a job is launched. Otherwise, the Monitor panel will not be displayed on job launch.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

shownewproperties

If set to 0, show all the new properties, for 1 shows only primary properties and for 2 hides new properties

Valid values: all
 primaryonly
 hidenew

Default value: **primaryonly**

Chapter 5: Commands

showobjectindex

Whether show the index of the object being picked in atom feedback.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

showpanelsontop

If this option is set to true then the panels will always be shown on top of the main window.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

showprojecttable

If set to true, then when a project is opened, the project table will automatically be displayed.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

showptpropertytree

If set to true, then property tree will be shown when showing the project table.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

showstars If set to true, then Starts property column is shown in the PT for new projects otherwise its hidden by default.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

spherequalityreduction

Controls how many steps we can reduce the quality by

Valid values: integers

Default value: **56**

Minimum: 0

statusbarshow

Whether the status bar is show or hide.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

styletoolbarstyle

Determines the appearance of the Style toolbar

Valid values: icononly
 textonly
 textbeside
 textunder

Default value: **textunder**

taskbasedmenu

If this option is set to true, the Applications and Workflows menus will be replaced by the Tasks menu.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

tbcollapsed

Whether the tool bar is collapsed or not.

Valid values: integers

Default value: **0**

tblocation

Whether the tool bar is located left or right.

Valid values: integers

Default value: **1**

tempprojectlocation

The path to create the Maestro temporary project

Valid values: text strings

Default value:

titlebardirectory

If this option is set to true, the main window's title bar will include the current working directory path.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

titlebarmmversion

If this option is set to true, the main window's title bar will include mmshare's version number.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

titlebarprofile

If this option is set to true, the main window's title bar will include current profile name.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

titlebarversion

If this option is set to true, the main window's title bar will include Maestro's version number.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

toolbarshow

Whether the tool bar is show or hide.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

Chapter 5: Commands

topviewshow

Whether the top view window is show or hide.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

trackpadrotsensitivity

Trackpad rotate sensitivity. The valid range is 1.0 to 100.0

Valid values: reals

Default value: **15**

Minimum: 1.0

Maximum: 100.0

trackpadzsensitivity

Trackpad Pinch gesture sensitivity for zoom. The valid range is 1.0 to 100.0. Larger is more sensitive

Valid values: reals

Default value: **25**

Minimum: 1.0

Maximum: 100.0

turnoffreapplystylewhenclosingproject

If this value is on then the Workspace Style reapply style setting will be turned off when a project is closed.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

useclipdistance

Set the flag of using clipping distance or not when focusing on a substructure.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

useremotedriver

If this flag is set, jobs launched on remote hosts will also have their driver script run on the remote host.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

viewanimation

Whether or not view animation is used on actions like fit to Workspace, spot centering, and view reset.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

viewanimationdistance

This value sets the maximum distance for view animation, if it is turned on.

Valid values: reals
 Default value: **10**
 Minimum: 0.0

viewanimationduration

This value sets the duration of each view animation, if it is turned on. The valid range is 0.0 to 5.0 seconds.

Valid values: reals
 Default value: **2**
 Minimum: 0.0
 Maximum: 5.0

viewname When saving a view, if set to 0, then 'Save View' dialog pops up asking for view name , else view will be stored with default name without prompting for name.

Valid values: prompt
 default
 Default value: **prompt**

viewtoolbarshow

Whether or not to show the Views toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

viewtoolbarshow

Whether or not to show the View toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

viewtoolbarstyle

Determines the appearance of the View toolbar

Valid values: icononly
 textonly
 textbeside
 textunder
 Default value: **icononly**

warnclosescratchproject

This option is used to record if warn users before closing scratch project

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

warnoverwritejobfiles

This option is used to record if warn users before existing job files being overwritten.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

Chapter 5: Commands

workspacemenushow

Whether the workspace menu should be displayed on right click.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

workspacetoolbarshow

Whether or not to show the Workspace toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

workspacetoolbarstyle

Determines the appearance of the Workspace toolbar

Valid values: icononly
 textonly
 textbeside
 textunder

Default value: **textunder**

writcmdlogonexit

If this option is set to true, cmd log will be saved when you exit Maestro.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

wsassistanttoolbarshow

Whether or not to show the Workspace Assistant toolbar.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

wsfeedbackfontbold

A bool which determines whether bold font is to be used for workspace feedback text

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

wsfeedbackfontitalics

A bool which determines whether italics font is to be used for workspace feedback text

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

wsfeedbackfontname

A string which determines the font to be used for workspace feedback text

Valid values: text strings

Default value:

wsfeedbackfontsize

A double which determines the font size to be used for workspace feedback text

Valid values: reals

Default value: **14**

Minimum: 0.0

wsselectpickstate

What kind of pick state the Workspace selection contains.

Valid values: atom
residue
molecule
chain
entry

Default value: **atom**

wsseqvieweralignbyresnum

A bool which determines whether workspace sequence viewer has to be aligned by residue number or not

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

wsseqviewerdisplaynonprotein

A bool which determines whether to display non-protein molecules in workspace sequence viewer.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

wsseqviewerdisplayssa

A bool which determines whether display SSA in workspace sequence viewer or not

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

wsseqviewerreorderresidues

Enables reordering residues by residue number.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

wsseqviewerwrap

A bool which determines whether workspace sequence viewer has to be wrapped or not

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

zeroorderbondtraversal

Whether Traverse zero order bonds in molecules.

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

previousresidue

If only a single residue is selected in the Workspace, then find the previous one.

Syntax:

previousresidue

primexaddwaters

Launches an add waters job.

Syntax:

primexaddwaters *bfactor_maximum*=⟨x⟩ *cutoff_maximum*=⟨x⟩
cutoff_minimum=⟨x⟩ *hbond_distance*=⟨x⟩ *peak_height*=⟨x⟩
⟨structure file⟩

Options:

bfactor_maximum

The maximum b-factor.

Valid values: reals

Default value: **100**

cutoff_maximum

The maximum distance cutoff, in angstroms.

Valid values: reals

Default value: **4**

Minimum: 0.0

cutoff_minimum

The minimum distance cutoff, in angstroms.

Valid values: reals

Default value: **2.6**

Minimum: 0.0

hbond_distance

Valid values: reals

Default value: **2.2**

peak_height

The minimum peak height.

Valid values: reals

Default value: **4**

Minimum: 0.0

Operands:

⟨structure file⟩

The name of a PDB structure file to use as the starting structure.

primexcalculatedensityfit

Calculates density fit for the given atoms.

Syntax:

primexcalculatedensityfit**primexcalculatedensitypeaks**

Calculates density peaks.

Syntax:

primexcalculatedensitypeaks *peak_threshold*=⟨x⟩

Options:

peak_threshold

Threshold for calculating density peaks. Only peaks that are greater than or equal to this number will be displayed.

Valid values: reals

Default value: **3**

Minimum: 1.0

Maximum: 10.0

primexcalculaterfactors

Calculate R-factors for the included structure

Syntax:

primexcalculatorfactors

primexcreatemap

Create map for Primex.

Syntax:

```
primexcreatemap map_2fofc=yes | no map_3fo2fc=yes | no
    map_box_xmax=<x> map_box_xmin=<x> map_box_ymax=<x>
    map_box_ymin=<x> map_box_zmax=<x> map_box_zmin=<x>
    map_extent=model | unitcell | asymmetric | box
    map_extent_box=fractional | orthogonal map_fo=yes | no
    map_fofc=yes | no map_grid_size=size033 | size025 | other
    map_grid_size_other=<x> map_weighting=sigma | unweighted
    maptype=regular | omit | composite max_displacement=<x>
    method=bhatcohen | kicked molecule_plus=<x>
    number_averaged=<n> omit=<x> scale_map=yes | no
```

Options:

map_2fofc Whether or not to create a 2Fo-Fc map
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

map_3fo2fc Whether or not to create a 3Fo-2Fc map
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

map_box_xmax The maximum X-coordinate for the map extent box.
Valid values: reals
Default value: **0**

map_box_xmin The minimum X-coordinate for the map extent box.
Valid values: reals
Default value: **0**

map_box_ymax The maximum Y-coordinate for the map extent box.
Valid values: reals
Default value: **0**

<i>map_box_ymin</i>	<p>The minimum Y-coordinate for the map extent box.</p> <p>Valid values: reals</p> <p>Default value: 0</p>
<i>map_box_zmax</i>	<p>The maximum Z-coordinate for the map extent box.</p> <p>Valid values: reals</p> <p>Default value: 0</p>
<i>map_box_zmin</i>	<p>The minimum Z-coordinate for the map extent box.</p> <p>Valid values: reals</p> <p>Default value: 0</p>
<i>map_extent</i>	<p>The extent of the map.</p> <p>Valid values: model unitcell asymmetric box</p> <p>Default value: model</p>
<i>map_extent_box</i>	<p>Which space the box extent is in.</p> <p>Valid values: fractional orthogonal</p> <p>Default value: fractional</p>
<i>map_fo</i>	<p>Whether or not to create an Fo map</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>
<i>map_fofc</i>	<p>Whether or not to create an Fo-Fc map</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>
<i>map_grid_size</i>	<p>The grid size for the maps.</p> <p>Valid values: size033 size025 other</p> <p>Default value: size033</p>
<i>map_grid_size_other</i>	<p>A user-specified grid size.</p> <p>Valid values: reals</p> <p>Default value: 0.1</p>

Chapter 5: Commands

map_weighting

The weighting for the maps.

Valid values: sigma
 unweighted

Default value: **sigma**

maptype

The map types of Primex create map.

Valid values: regular
 omit
 composite

Default value: **regular**

max_displacement

The maximum displacement for kicked method.

Valid values: reals

Default value: **0.1**

Minimum: 0.0

method

The methods of Primex create map.

Valid values: bhatcohen
 kicked

Default value: **bhatcohen**

molecule_plus

The amount in angstroms to increase the model by when calculating the extent for generating maps.

Valid values: reals

Default value: **5**

Minimum: 0.0

number_averaged

The number averaged for kicked method.

Valid values: integers

Default value: **10**

Minimum: 0

omit

The percentage of data to omit.

Valid values: reals

Default value: **5**

Minimum: 0.0

scale_map

Indicates whether or not to scale the map to sigma units

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

primexdecreaseisovalue

Decrease isovalue for currently selected map surface.

Syntax:

```
primexdecreaseisovalue entry=<text> numsigma=<x>
                        volume=<text> <surface>
```

Options:

<i>entry</i>	The entry name of the entry that the surface belongs to. Valid values: text strings Default value:
<i>numsigma</i>	The amount of decreasing isovalue, e.g. 0.1*sigma. Valid values: reals Default value: 0.1 Minimum: 0.01
<i>volume</i>	The name of the volume that the surface belongs to. Valid values: text strings Default value:

Operands:

<surface>

The name of the surface to decrease the isovalue for.

primexdeleteallrigidbodies

Deletes all rows in the rigid bodies table.

Syntax:

```
primexdeleteallrigidbodies
```

primexdeleteoccupancy

Deletes the partial occupancy data for the given atom

Syntax:

primexdeleteoccupancy \langle atom index \rangle

Operands:

\langle atom index \rangle

An atom index for the main CT.

primexdeleterigidbodies

Deletes the selected rows in the rigid bodies table.

Syntax:

primexdeleterigidbodies

primexentrywscreate

Creates a project entry from atoms in the workspace. The entry name and title are derived from the current WS entry, and all non-structure-specific PrimeX properties are copied to the new entry.

Syntax:

primexentrywscreate

primexexport

Export the included structure as a formatted PDB file

Syntax:

primexexport \langle file name \rangle

Operands:

\langle file name \rangle

The name of the file to export the PrimeX data to.

primexincreaseisovalue

Increase isovalue for currently selected map surfaces.

Syntax:

```
primexincreaseisovalue entry=⟨text⟩ numsigma=⟨x⟩
                        volume=⟨text⟩ ⟨surface⟩
```

Options:

entry The entry name of the entry that the surface belongs to.

Valid values: text strings
Default value:

numsigma The amount of increasing isovalue, e.g. 0.1*sigma.

Valid values: reals
Default value: **0.1**
Minimum: 0.01

volume The name of the volume that the surface belongs to.

Valid values: text strings
Default value:

Operands:

⟨surface⟩

The name of the surface to increase the isovalue for.

primexinputdata

Adds or replaces the PrimeX data with the given data.

Syntax:

```

primexinputdata assign_bond_orders=yes | no cella=⟨x⟩
                 cellalpha=⟨x⟩ cellb=⟨x⟩ cellbeta=⟨x⟩ cellc=⟨x⟩
                 cellgamma=⟨x⟩ constant_bfactor=⟨x⟩ foproperty=⟨text⟩
                 generate_testset=yes | no map_file=⟨text⟩ map_type=mapfo |
                 mapfofc | map2fofc | map3fo2fc reflectionfile=⟨text⟩
                 sequence_source=selected_entries | workspace | file
                 sequencefile=⟨text⟩ sigmaproperty=⟨text⟩ spacegroup=⟨n⟩
                 spacegroupname=⟨text⟩ structure_source=selected_entries |
                 workspace | file test_percentage=⟨x⟩ testproperty=⟨text⟩
                 use_bfactors=current | constant ⟨structure file⟩

```

Options:

assign_bond_orders

This option assigns correct bond orders to the input structure.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

cella

The unit cell a parameter.

Valid values: reals

Default value: **0**

cellalpha

The unit cell alpha parameter.

Valid values: reals

Default value: **90**

cellb

The unit cell b parameter.

Valid values: reals

Default value: **0**

cellbeta

The unit cell beta parameter.

Valid values: reals

Default value: **90**

cellc

The unit cell c parameter.

Valid values: reals

Default value: **0**

cellgamma

The unit cell gamma parameter.

Valid values: reals

Default value: **90**

constant_bfactor

The constant value to set the B-factors to.

Valid values: reals

Default value: **20**

foproperty The Fo property name for mtz reflection data file.

Valid values: text strings

Default value: **FP**

generate_testset

This option indicates whether or not a new random test set should be generated in the reflection file.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

map_file The map data file.

Valid values: text strings

Default value:

map_type What kind of map the map file contains.

Valid values: mapfo
mapfofc
map2fofc
map3fo2fc

Default value: **mapfo**

reflectionfile

The crystal reflection data file.

Valid values: text strings

Default value:

sequence_source

Whether to get the sequence from an external file or from the contents of the Workspace.

Valid values: selected_entries
workspace
file

Default value: **file**

sequencefile

The sequence file.

Valid values: text strings

Default value:

sigmaproperty

The Sigma(Fo) property name for mtz reflection file.

Valid values: text strings

Default value: **SIGFP**

spacegroup

The space group number.

Valid values: integers

Default value: **0**

Chapter 5: Commands

spacegroupname

The space group name.

Valid values: text strings

Default value:

structure_source

Whether to get the structure from an external file or from the contents of the Workspace.

Valid values: selected_entries
workspace
file

Default value: **file**

test_percentage

This is the percentatge of reflections to use as the test set.

Valid values: reals

Default value: **5**

Minimum: 0.01

Maximum: 100

testproperty

The test set property name for mtz reflection file.

Valid values: text strings

Default value: **FreeR_flag**

use_bfactors

Whether to keep the existing B-factor values or to replace them with a constant value.

Valid values: current
constant

Default value: **current**

Operands:

⟨ structure file ⟩

The name of a PDB structure file to use as the starting structure.

primexmarkersettings

Set graphical data of PrimeX peak markers.

Syntax:

```

primexmarkersettings ambient=⟨x⟩ diffuse=⟨x⟩ emission=⟨x⟩
                     num_slices=⟨n⟩ num_stacks=⟨n⟩ shininess=⟨x⟩
                     specular=⟨x⟩ transparency=⟨x⟩

```

Options:

<i>ambient</i>	Set the ambient material property for the peak markers. Valid values: reals Default value: 0.5 Minimum: 0.0 Maximum: 1.0
<i>diffuse</i>	Set the diffuse material property for the peak markers. Valid values: reals Default value: 0.4 Minimum: 0.0 Maximum: 1.0
<i>emission</i>	Set the emission material property for the peak markers. Valid values: reals Default value: 0.05 Minimum: 0.0 Maximum: 1.0
<i>num_slices</i>	Set the number of slices for the peak markers. Valid values: integers Default value: 18 Minimum: 2
<i>num_stacks</i>	Set the number of stacks for the peak markers. Valid values: integers Default value: 9 Minimum: 2
<i>shininess</i>	Set the shininess material property for the peak markers. Valid values: reals Default value: 80 Minimum: 0.0 Maximum: 128.0
<i>specular</i>	Set the specular material property for the peak markers. Valid values: reals Default value: 0 Minimum: 0.0 Maximum: 1.0

transparency

The transparency percentage of the peak markers.

Valid values: reals
Default value: **50**
Minimum: 0.0
Maximum: 100.0

primexmutatemodeltosequence

Mutate the model in the workspace (which must be a single entry) to match the reference sequence specified for PrimeX.

Syntax:

primexmutatemodeltosequence *delete_gap_residues*=yes | no

Options:

delete_gap_residues

This options controls whether residues that are aligned with gaps are deleted.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

primexncsaddrow

Adds a row with default settings to the NCS table.

Syntax:

primexncsaddrow

primexncsdeleterows

Deletes the selected rows from the NCS tab.

Syntax:

primexncsdeleterows

primexncseditcell

Edit a cell in the NCS table

Syntax:

```
primexncseditcell column=chains | residues | bfactor | coord
                    row=⟨n⟩ ⟨value⟩
```

Options:

column The column of the NCS cell to edit.

Valid values: chains
 residues
 bfactor
 coord

Default value: **chains**

row The row of the NCS cell to edit.

Valid values: integers
 Default value: **1**
 Minimum: 1

Operands:

⟨value⟩

The new value for the cell.

primexncsselectextendrow

Extends the selection to this row in the NCS table in PrimeX.

Syntax:

```
primexncsselectextendrow ⟨row⟩
```

Operands:

⟨row⟩

The row number to extend the select to.

primexncsselectonlyrow

Selects only this row in the NCS table in PrimeX.

Syntax:

```
primexncsselectonlyrow <row>
```

Operands:

<row>

The row number to select only in the table row.

primexncsselectrow

Selects the given row in the NCS table in PrimeX.

Syntax:

```
primexncsselectrow <row>
```

Operands:

<row>

The row number to select in the table.

primexncstogglebbonly

Toggles the BB Only value for the given row in the NCS table in PrimeX.

Syntax:

```
primexncstogglebbonly <row>
```

Operands:

<row>

The row number to toggle the BB Only value for.

primexncsunselectrow

Unselects the given row in the NCS table in PrimeX.

Syntax:

```
primexncsunselectrow <row>
```

Operands:

<row>

The row number to unselect in the table.

primexpartialoccupancy

Sets the partial occupancy for selected atoms

Syntax:

```
primexpartialoccupancy occupancy=<x> <ASL>
```

Options:

occupancy The partial occupancy.

Valid values: reals

Default value: **0.5**

Minimum: 0.0

Maximum: 1.0

Operands:

<ASL>

An ASL representing the atoms to set the partial occupancy for.

primexplaceligand

Launches a place ligand/solvent job.

Syntax:

```
primexplaceligand box_buffer=<x> energy_weight=<x>  
generate_mates=yes | no ligand=<text> ligandsfile=<text>  
ligandsource=file | entries | solvent refine_poses=yes | no  
scoring_weight=<x> use_ligprep=yes | no <job name>
```

Options:

Chapter 5: Commands

box_buffer Controls the size of the buffer added to the box which encloses the selected density peaks in order to create the grid box for ligand placement.

Valid values: reals

Default value: **16**

energy_weight

Energy density weight.

Valid values: reals

Default value: **120**

Minimum: 1.0

Maximum: 1000.1

generate_mates

Indicates whether or not to generate symmetry mates near the site.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

ligand

The entry name of the ligand to use.

Valid values: text strings

Default value:

ligandsfile

The file containing one or more ligands.

Valid values: text strings

Default value:

ligandsource

The source of the ligands to be used.

Valid values: file

entries

solvent

Default value: **file**

refine_poses

Set this option to refine the poses that Glide generates using the X-ray data.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

scoring_weight

Scoring density weight.

Valid values: reals

Default value: **50**

Minimum: 1.0

Maximum: 1000.1

use_ligprep

Indicates whether or not to run LigPrep on the ligand before placing it.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

Operands:

⟨job name⟩

The name of the job.

primexpublish

Exports the structure in the Workspace as the final structure for a crystallographic refinement.

Syntax:

primexpublish ⟨file⟩

Operands:

⟨file⟩

The name of a PDB file to save as the final structure.

primexrealspace

Launches a real-space refinement job.

Syntax:

```
primexrealspace bias_existing=yes | no grid_size=⟨x⟩
                  loop_bfactorsfit=⟨x⟩ loop_chain=⟨text⟩ loop_end=⟨text⟩
                  loop_setbfactorsfit=yes | no loop_start=⟨text⟩
                  minimize_bfactorsfit=⟨x⟩ minimize_grid_size=⟨x⟩
                  minimize_setbfactorsfit=yes | no sidechain_bfactorsfit=⟨x⟩
                  sidechain_grid_size=⟨x⟩ sidechain_setbfactorsfit=yes | no
                  ⟨structure file⟩
```

Options:

Chapter 5: Commands

bias_existing

If set to true, then PrimeX will bias loops towards the existing loop structure. If set to false, then PrimeX will bias loops towards the reflection data.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

grid_size

The grid size used for real-space loop refinement jobs.

Valid values: reals

Default value: **1**

Minimum: 0.5

Maximum: 1.0

loop_bfactorsfit

The radius of setting b_factors of atoms to be fit for a real-space loop refinement job.

Valid values: reals

Default value: **5**

Minimum: 0.0

loop_chain

This is the chain for loop refinement.

Valid values: text strings

Default value:

loop_end

This is the ending residue for loop refinement.

Valid values: text strings

Default value:

loop_setbfactorsfit

An option which determines if the real-space loop refinement job will set b-factors of atoms to be fit.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

loop_start

This is the starting residue for loop refinement.

Valid values: text strings

Default value:

minimize_bfactorsfit

The radius of setting b-factors of atoms to be fit for a real-space minimization job.

Valid values: reals

Default value: **20**

Minimum: 0.0

minimize_grid_size

The grid size used for real-space minimization refinement jobs.

Valid values:	reals
Default value:	0.5
Minimum:	0.5
Maximum:	1.0

minimize_setbfactorstfit

An option which determines if the real-space minimization job will set b-factors of atoms to be fit.

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

sidechain_bfactorstfit

The radius of setting b-factors of atoms to be fit for a real-space sidechain refinement job.

Valid values:	reals
Default value:	20
Minimum:	0.0

sidechain_grid_size

The grid size used for sidechain refinement jobs.

Valid values:	reals
Default value:	0.5
Minimum:	0.5
Maximum:	1.0

sidechain_setbfactorstfit

An option which determines if the real-space loop refinement job will set b-factors of atoms to be fit.

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

Operands:

⟨structure file⟩

The name of a PDB structure file to use as the starting structure.

primexreciprocal

Launches a reciprocal-space refinement job.

Syntax:

```

primexreciprocal bb_angles=⟨x⟩ bb_bonds=⟨x⟩
    bfactor=individual | anisotropic bfactor_high_limit=⟨x⟩
    bfactor_low_limit=⟨x⟩ bfactor_restraint_weighting=⟨x⟩
    constant_weight=⟨x⟩ constantvalue=⟨x⟩ cooling_steps=⟨n⟩
    energy_model=approximate | complete final_steps=⟨n⟩
    final_temp=⟨x⟩ heating_steps=⟨n⟩ high_temp=⟨x⟩
    initial_steps=⟨n⟩ initial_temp=⟨x⟩ map_2fofc=yes | no
    map_3fo2fc=yes | no map_fo=yes | no map_fofc=yes | no
    md_steps=⟨n⟩ md_time_step=⟨x⟩ method=minimization |
    rigidbodies | simulatedannealing
    min_minimizer=truncatednewton | conjugategradient |
    quasineutron | optimal minimization_max_cycles=⟨n⟩
    minimization_max_steps=⟨n⟩ minimize=coordinates | bfactors
    | occupancies | groupedbfactors nonbonded_cutoff=⟨x⟩
    rb_minimizer=truncatednewton | conjugategradient |
    quasineutron | optimal rfree=⟨x⟩ rigidmaxcycles=⟨n⟩
    rigidmaxsteps=⟨n⟩ side_angles=⟨x⟩ side_bonds=⟨x⟩
    target=likelihood | leastsquares ungroupedatoms=fix | group
    usevalues=current | constant weight_multiplication_factor=⟨x⟩
    xray_weight=multiply | constant ⟨structure file⟩

```

Options:

<i>bb_angles</i>	The target sigma value for B-factor restraints for the backbone angles. Valid values: reals Default value: 2
<i>bb_bonds</i>	The target sigma value for B-factor restraints for the backbone bonds. Valid values: reals Default value: 1.5
<i>bfactor</i>	The atomic B-factor for minimization Valid values: individual anisotropic Default value: individual
<i>bfactor_high_limit</i>	The high-resolution limit for B-factors. Valid values: reals Default value: 100 Minimum: 0.0
<i>bfactor_low_limit</i>	The low-resolution limit for B-factors.

Valid values: reals
 Default value: **2**
 Minimum: 0.0

bfactor_restraint_weighting

The B-factor restraint weighting factor

Valid values: reals
 Default value: **1**

constant_weight

A constant weight to replace the automatically calculated weight with.

Valid values: reals
 Default value: **0.5**

constantvalue

Size of the constant B-factor value to set to.

Valid values: reals
 Default value: **20**

cooling_steps

The number of cooling steps.

Valid values: integers
 Default value: **1250**
 Minimum: 1

energy_model

The energy model for simulated annealing.

Valid values: approximate
 complete
 Default value: **approximate**

final_steps The number of final minimization steps.

Valid values: integers
 Default value: **30**
 Minimum: 1

final_temp The final temperature, in Kelvin.

Valid values: reals
 Default value: **300**
 Minimum: 0

heating_steps

The number of heating steps.

Valid values: integers
 Default value: **125**
 Minimum: 1

Chapter 5: Commands

<i>high_temp</i>	The high temperature, in Kelvin. Valid values: reals Default value: 700 Minimum: 0
<i>initial_steps</i>	The number of initial minimization steps. Valid values: integers Default value: 30 Minimum: 1
<i>initial_temp</i>	The initial temperature, in Kelvin. Valid values: reals Default value: 50 Minimum: 0
<i>map_2fofc</i>	Whether or not to create a 2Fo-Fc map Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>map_3fo2fc</i>	Whether or not to create a 3Fo-2Fc map Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>map_fo</i>	Whether or not to create an Fo map Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>map_fofc</i>	Whether or not to create an Fo-Fc map Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>md_steps</i>	The number of MD energy scale estimation steps. Valid values: integers Default value: 100 Minimum: 1
<i>md_time_step</i>	The MD time step. Valid values: reals Default value: 0.008
<i>method</i>	The refinement method to use. Valid values: minimization rigidbodies simulatedannealing

	Default value: minimization
<i>min_minimizer</i>	Which minimizer to use for Minimization.
	Valid values: truncatednewton conjugategradient quasinewton optimal
	Default value: optimal
<i>minimization_max_cycles</i>	The maximum number of minimization cycles for minimization.
	Valid values: integers
	Default value: 3
	Minimum: 1
<i>minimization_max_steps</i>	The maximum number of minimization steps per cycle for minimization.
	Valid values: integers
	Default value: 8
	Minimum: 0
<i>minimize</i>	The minimize option for Minimization.
	Valid values: coordinates bfactors occupancies groupedbfactors
	Default value: coordinates
<i>nonbonded_cutoff</i>	The non-bonded interactions cutoff, in angstroms.
	Valid values: reals
	Default value: 9.5
<i>rb_minimizer</i>	Which minimizer to use for Rigid Bodies.
	Valid values: truncatednewton conjugategradient quasinewton optimal
	Default value: optimal
<i>rfree</i>	The R-free percentage.
	Valid values: reals
	Default value: 0

Chapter 5: Commands

rigidmaxcycles

Maximum number of minimization cycles for a rigid body refinement.

Valid values: integers

Default value: **3**

Minimum: 1

rigidmaxsteps

Maximum number of minimization steps per cycle for a rigid body refinement.

Valid values: integers

Default value: **8**

Minimum: 1

side_angles

The target sigma value for B-factor restraints for the side chain angles.

Valid values: reals

Default value: **2.5**

side_bonds

The target sigma value for B-factor restraints for the side chain bonds.

Valid values: reals

Default value: **2**

target

The refinement target.

Valid values: likelihood
least-squares

Default value: **likelihood**

ungroupedatoms

How should ungrouped atoms in a Rigid Body refinement be treated?

Valid values: fix
group

Default value: **fix**

usevalues

Whether to use the existing B-factor values or whether to replace them with a constant value.

Valid values: current
constant

Default value: **current**

weight_multiplication_factor

The factor to multiply the automatically calculated weight by.

Valid values: reals

Default value: **1**

xray_weight

Whether to multiply the weight by a factor or replace it with a constant term.

Valid values: multiply
constant

Default value: **constant**

Operands:

⟨structure file⟩

The name of a PDB structure file to use as the starting structure.

primexselectdensityblob

Selects a row in the density blobs table.

Syntax:

primexselectdensityblob

primexselectdensitypeak

Selects a row in the density peaks table.

Syntax:

primexselectdensitypeak

primexselectextenddensitypeak

Extends the selection in the density peaks table to the given row.

Syntax:

primexselectextenddensitypeak

primexselectextendligandsolvent

Extends the selection to this row in the ligand / solvents table.

Syntax:

```
primexselectextendligandsolvent
```

primexselectligandsolvent

Selects a row in the ligand / solvents table.

Syntax:

```
primexselectligandsolvent
```

primexselectonlydensitypeak

Selects only the given row in the density peaks table.

Syntax:

```
primexselectonlydensitypeak
```

primexselectonlyligandsolvent

Selects a row in the ligand / solvents table.

Syntax:

```
primexselectonlyligandsolvent
```

primexselectonlyrigidbody

Selects only the given rigid body row.

Syntax:

primexselectonlyrigidbody

primexselectrigidbody

Selects the given rigid body row.

Syntax:

primexselectrigidbody

primexsetoccupancy

Sets the partial occupancy data for the given atom

Syntax:

primexsetoccupancy *occupancy*=⟨x⟩ ⟨atom index⟩

Options:

occupancy The partial occupancy value to set.

Valid values:	reals
Default value:	0.5
Minimum:	0.0
Maximum:	1.0

Operands:

⟨atom index⟩

An atom index for the main CT.

primexsettings

Holds options for PrimeX jobs.

Syntax:

```

primexsettings covalent_radius=<x> filterf=<x> filterrms=<x>
                 filtersigma=<x> high_resolution=<x>
                 implicit_solvation=yes | no include_hatoms=yes | no
                 ion_radius=<x> low_resolution=<x> maprproperty=<text>
                 mapsource=model | file mapsourcefile=<text>
                 max_memory=<x> phiproperty=<text>
                 planargrouprestraints=low | normal | high postrefine=yes | no
                 rejectreflections=yes | no resolution_type=calculate | define
                 scaling=none | isotropic | anisotropic | both
                 shrink_factor=<x> solventmethod=mask | babinet | none
                 weightingfactor=<x> <none>

```

Options:

covalent_radius

The VdW radius for covalent atoms.

Valid values: reals

Default value: **1.4**

filterf

Reject reflections with $F < \text{this value}$.

Valid values: reals

Default value: **0**

filterrms

Reject reflections with $F / \text{rms}(F) > \text{this value}$.

Valid values: reals

Default value: **10000**

filtersigma

Reject reflections with $F / \text{sigma}(F) < \text{this value}$.

Valid values: reals

Default value: **0**

high_resolution

The high-resolution limit, in angstroms.

Valid values: reals

Default value: **1**

Minimum: 0.2

implicit_solvation

Whether or not to use implicit solvation.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

include_hatoms

Whether or not to include H atoms.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

- ion_radius* The radius for ions.
 Valid values: reals
 Default value: **0.8**
- low_resolution*
 The low-resolution limit, in angstroms.
 Valid values: reals
 Default value: **30**
 Minimum: 2.0
- maprproperty*
 The map F property name for mtz reflection data file.
 Valid values: text strings
 Default value:
- mapsource* The map source of Primex create map.
 Valid values: model
 file
 Default value: **model**
- mapsourcefile*
 The input map source file.
 Valid values: text strings
 Default value:
- max_memory*
 Maximum memory to use for calculations.
 Valid values: reals
 Default value: **400**
 Minimum: 10.0
- phiproperty*
 The phi property name for mtz reflection file.
 Valid values: text strings
 Default value:
- planargrouprestraints*
 The planar group restraints level (low, normal or high).
 Valid values: low
 normal
 high
 Default value: **normal**
- postrefine* Turns on running a reciprocal-space minimization job after a real-space refinement job.
 Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

Chapter 5: Commands

rejectreflections

Indicates whether or not to filter reflections.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

resolution_type

Specifies how the low and high resolution limit being defined: 1) calculate from data and cell constants, and 2) use data set by default or by users.

Valid values: calculate
 define

Default value: **calculate**

scaling

The overall B-factor scaling.

Valid values: none
 isotropic
 anisotropic
 both

Default value: **anisotropic**

shrink_factor

The shrink factor.

Valid values: reals

Default value: **1.4**

solventmethod

The solvent correction method.

Valid values: mask
 babinet
 none

Default value: **mask**

weightingfactor

The weighting factor.

Valid values: reals

Default value: **1**

Minimum: 0.0

Operands:

⟨ none ⟩

primexswitchaltpositions

Swaps the real and alternate positions (occupancies) for the given atoms.

Syntax:

`primexswitchaltpositions` \langle ASL \rangle

Operands:

\langle ASL \rangle

The atoms to swap positions.

primextoggleselectrigidbody

Toggles the selection of the given rigid body row.

Syntax:

`primextoggleselectrigidbody`

primexunitcell

This keyword is used to set options associated with the PrimeX unit cell.

Syntax:

`primexunitcell` *showmarker*=yes | no

Options:

showmarker

This determines whether the unit cell is shown or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

primexunselectdensitypeak

Unselects a row in the density peaks table.

Syntax:

`primexunselectdensitypeak`

primexunselectligandsolvent

Unselects a row in the ligand / solvents table.

Syntax:

```
primexunselectligandsolvent
```

primexwriteaddwaters

Write the inputs files for an add waters job.

Syntax:

```
primexwriteaddwaters
```

primexwritecalculaterfactors

Write a calculate R-factors job for the included structure

Syntax:

```
primexwritecalculaterfactors
```

primexwritecreatemap

Write a create map job for Primex.

Syntax:

```
primexwritecreatemap
```

primexwriteexport

Write an export job for the included structure

Syntax:

primexwriteexport

primexwriteplaceligand

Writes a place ligand/solvent job.

Syntax:

primexwriteplaceligand

primexwriterealspace

Writes the input files for a real-space refinement job.

Syntax:

primexwriterealspace

primexwritereciprocal

Write the input files for a reciprocal-space refinement job.

Syntax:

primexwritereciprocal

profile

Specifies a profile of Maestro

Syntax:

profile *current*=yes | no *default*=yes | no *directory*= \langle text \rangle
parent= \langle text \rangle *select*=yes | no *standard*=yes | no

Options:

<i>current</i>	State indicates whether the profile is a current profile or not.
Valid values:	boolean (true false; yes no; y n; on off)
Default value:	false

Chapter 5: Commands

<i>default</i>	State indicates whether the profile is a default profile or not. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>directory</i>	This is profile directory path. Valid values: text strings Default value: dir
<i>parent</i>	This is the parent of profile. Valid values: text strings Default value: parent
<i>select</i>	State indicates whether the profile is a currently elected or not.. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>standard</i>	State indicates whether the profile is a standard profile or not. Valid values: boolean (true false; yes no; y n; on off) Default value: false

projectbackup

Perform project backup. See projectbackupprefs for the preferences that control backups. Logging is done to schrodinger.autobackup.log in your home directory on Linux or to your %LOCALAPPDATA% on Windows.

Syntax:

projectbackup <comment_string>

Operands:

<comment_string>

The comment to include in the backup

projectbackuppreference

Set project backup preferences. These are persistent settings.

Syntax:

```
projectbackuppreference auto=yes | no dialog=never | always
                        | errors emaildebug=yes | no emailfrom=<text>
                        emailnotification=yes | no emailpassword=<text>
                        emailport=<n> emailsecurity=none | starttls | ssltls
                        emailsmtplibserver=<text> emailto=<text>
                        limitbackups=yes | no maxbackups=<n> time=<n>
```

Options:

- auto* When on perform an automatic backup once a day at the specified time. If off, no automatic backup is performed.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**
- dialog* When to show a notification dialog. This can be 'never' in which case no dialog is shown, 'always' in which case a dialog is always shown regardless of success or failure, and 'errors' meaning only show the dialog when there are errors. This setting applies to projectbackup. It applies to projectcheck when projectcheck is run from the automated backup. Otherwise, projectcheck always displays a dialog when an error is encountered.
Valid values: never
always
errors
Default value: **errors**
- emaildebug* If on, debug output is printed at each step of the EMail process. This value is not persistent and is reset to its default every time Maestro starts.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**
- emailfrom* Full e-mail From address
Valid values: text strings
Default value:
- emailnotification* Upon completion of the backup attempt whether or not to send an e-mail notification indicating the status. Only sending email via an smtp server email is supported.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**
- emailpassword* If sending email needs a password, specify it here.
Valid values: text strings

Chapter 5: Commands

	Default value:
<i>emailport</i>	Port to use for SMTP EMail Server
	Valid values: integers
	Default value: 587
<i>emailsecurity</i>	Type of security to use for EMail
	Valid values: none starttls ssltls
	Default value: none
<i>emailsmtpserver</i>	Full name of the smtp server
	Valid values: text strings
	Default value:
<i>emailto</i>	Full e-mail address to which EMail notifications will be sent
	Valid values: text strings
	Default value:
<i>limitbackups</i>	If false, the number of backups is unlimited and the maxbackups option is ignored. If enabled, then maxbackups is honored.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>maxbackups</i>	Specifies maximum number of backups allowed for projects. If limitbackups is enabled and this number is exceeded when a new backup is made, then the oldest backups will be deleted to stay within this limit. This option replaces the 'prefer' command's projectbackuplimit option.
	Valid values: integers
	Default value: 3
	Minimum: 1
<i>time</i>	When automatic backup is enabled start the backup some time at this time (may not trigger exactly at this time because it's on a timer but will happen within the hour of the specified time). Time is in 24-hour format.
	Valid values: integers
	Default value: 2
	Minimum: 0
	Maximum: 23

projectcheck

Check a project for consistency. Currently, only prjzip's are supported. Further, it has to be a project saved with the 'medium' (Fast close) projectsave setting, i.e. one which has a .mmpproj-admin directory. See the `prefer` command's `projectsave` option for further information.

Syntax:

```
projectcheck checklevel=none | exists | simple | requiredfiles |
               count | filesandoffsets <project path>
```

Options:

checklevel Specifies how much to check the project Used internally by Maestro to check a backup. Can be used by user, but use values of 3 and greater carefully. 0: no checking, 1: see if prjzip exists, 2: check prjzip to make sure key key files and directories are present, 3: Does 1 and 2 but also checks to make sure that the number of entries and properties matches the currently open project. This is fairly quick and does not require unzipping the full project. 4: Does 1-3 and also ensures that for each entry there is a corresponding structure file and that the starting offset can be reached (is within bounds). This requires unzipping the whole prjzip. So you may need a lot of disk and it may be slow. Note that 1, 2 and 4 can be used on any project that is not open in Maestro. 3, however, uses Maestro's currently opened project to compare against. So make sure that when using 3 the currently opened project and the one you name are correctly paired! This option is persistent across Maestro sessions.

Valid values: none
 exists
 simple
 requiredfiles
 count
 filesandoffsets

Default value: **simple**

Operands:

<project path>

The project which we are checking. Currently only prjzip files are supported.

projectclose

Close the current project and open a new scratch project.

Syntax:

`projectclose`

projectcopy

Make a copy of the current project.

Syntax:

`projectcopy combiglide=yes | no entry=yes | no job=yes | no
oldformat=yes | no phase=yes | no plot=yes | no
prime=yes | no prjzip=yes | no run=yes | no
scenes=yes | no snapshot=yes | no source=all | selected |
included table=yes | no user=yes | no <to_dir_path>`

Options:

<i>combiglide</i>	Enable/disable copying of project combiglide data. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>entry</i>	Enable/disable copying of project entry files. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>job</i>	Enable/disable copying of project job files. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>oldformat</i>	Enable/disable copying of project in old format. For any Maestro version, the old format will be the immediate previous version. Say for Suite2011 Maestro, if this option is set to true, then the project will be copied in Suite2011 format. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>phase</i>	Enable/disable copying of project phase data. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>plot</i>	Enable/disable copying of project plots. If this option is enabled, entries should also be copied.

	Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>prime</i>	Enable/disable copying of project prime data. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>prjzip</i>	If on, create a .prjzip. If off, create a .prj. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>run</i>	Enable/disable copying of project run files. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>scenes</i>	Enable/disable copying of scenes. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>snapshot</i>	Enable/disable copying of project saved state. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>source</i>	Selects the source of entries. Valid values: all selected included Default value: all
<i>table</i>	Enable/disable copying of project table files. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>user</i>	Enable/disable copying of project user files. Valid values: boolean (true false; yes no; y n; on off) Default value: false

Operands:

⟨to_dir_path⟩

The path (name and location) of the directory to be used to copy the current project. This can either be a new or existing project directory.

projectdelete

Delete the current project and open a new scratch project.

Syntax:

projectdelete

projectmerge

Merge data from another project into the current Maestro project.

Syntax:

projectmerge *entry*=yes | no *job*=yes | no *run*=yes | no
user=yes | no *<from_dir_path>*

Options:

<i>entry</i>	Enable/disable merging of project entry files. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>job</i>	Enable/disable merging of project job files. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>run</i>	Enable/disable merging of project run files. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>user</i>	Enable/disable merging of project user files. Valid values: boolean (true false; yes no; y n; on off) Default value: false

Operands:

<from_dir_path>

The path (name and location) of the project directory to be merged into the current project.

projectnew

Create a new project and open it in Maestro.

Syntax:

projectnew $\langle \text{dir_path} \rangle$

Operands:

$\langle \text{dir_path} \rangle$

The path (name and location) of the directory to be created as the new project directory.

projectopen

Open an existing project into Maestro.

Syntax:

projectopen $\langle \text{dir_path} \rangle$

Operands:

$\langle \text{dir_path} \rangle$

The path (name and location) of the project directory to be opened.

projectprefer

Set project-specific (rather than user-specific) preferences.

Syntax:

projectprefer *2dstructureheight*= $\langle n \rangle$ *2dstructurewidth*= $\langle n \rangle$
entryfeedbackshow=yes | no *entryfeedbackshownames*=yes | no
feedbackproperties= $\langle \text{text} \rangle$ *feedbackproperty*= $\langle \text{text} \rangle$
feedbackshow=yes | no *jobstatusfeedbackshow*=yes | no
phasefeedbackproperties= $\langle \text{text} \rangle$ *show2dstructures*=yes | no

Options:

2dstructureheight

The height of the 2D structure in the table.

Valid values: integers

Default value: **200**

Minimum: 16

Maximum: 1000

2dstructurewidth

The width of the 2D structure in the table.

Chapter 5: Commands

Valid values: integers
Default value: **300**
Minimum: 1
Maximum: 5000

entryfeedbackshow

Whether the workspace entry feedback is activated.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

entryfeedbackshownames

Whether the workspace entry feedback includes the name of the property along with its value.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

feedbackproperties

These are the entry properties displayed for a single entry in workspace. The property names should be m2io data names with no spaces, separated by spaces in the option string.

Valid values: text strings
Default value: **s_m_title s_pdb_PDB_ID**

feedbackproperty

This is the entry property displayed at the end of the atom feedback string in workspace.

Valid values: text strings
Default value: **Title**

feedbackshow

Whether the workspace feedback is activated.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

jobstatusfeedbackshow

Whether the workspace job status feedback is activated.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

phasefeedbackproperties

These are the Phase properties displayed for a single Phase ligand in the workspace. The property names should be m2io data names with no spaces, separated by spaces in the option string.

Valid values: text strings
Default value:

show2dstructures

Whether the 2D structures are shown in the project table.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

projectpublish

Publish this project to an html file with links to files in a corresponding subdirectory. The files include among other things a .prjzip version of the current project. The html and corresponding directory may then be used to post this project for use on the web.

Syntax:

```
projectpublish annotation=yes | no publishsummary=<text>  

               <html_file>
```

Options:

annotation

Enable/disable inclusion of the project annotation. This text appears in the Annotation section of the Project Summary panel.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

publishsummary

Summary to be included in the published project. The text to use as the project summary. If you save successive projects, you must update this string as the current value will be written the the project you are publishing.

Valid values: text strings
 Default value:

Operands:

<html_file>

The path (name and location) of the html file. A corresponding subdirectory using the base name of the html file will be crated.

projectrename

Rename the current project directory and/or move it to a new location.

Syntax:

projectrename $\langle \text{to_dir_path} \rangle$

Operands:

$\langle \text{to_dir_path} \rangle$

The new path (name and location) for the current project directory.

projectrestore

Replace specified project with a backup created earlier for that project. Does not open the specified project. If the specified project is the current project, it will be closed first. This command is not undoable. This command expects the project to be specified as a .prj directory, not as a .prjzip or .prj.zip file.

Syntax:

projectrestore $\langle \text{dir_path} \rangle$ $\langle \text{backup_name} \rangle$

Operands:

$\langle \text{dir_path} \rangle$ $\langle \text{backup_name} \rangle$

The path (name and location) of the project .prj directory to be replaced. Both absolute and relative paths are accepted. The name of the backup project file that will be used to replace the project.

projectrevertopen

Revert state of current project to that it had when opened in Maestro.

Syntax:

projectrevertopen

projectrevertsnapshot

Revert state of current project to that previously stored. Deprecated in favor of `projectrevertstate`.

Syntax:

projectrevertsnapshot

projectrevertstate

Revert state of current project to that previously stored.

Syntax:

projectrevertstate

projectsaveas

Save project and place user into that project

Syntax:

projectsaveas \langle to_dir_path \rangle

Operands:

\langle to_dir_path \rangle

Name to which the project will be saved. Saves all data from the current project (table, plots, etc.)

projectsceneclose

Close the open scene and return to the master project.

Syntax:

projectsceneclose

projectscenedelete

Deletes the selected scenes.

Syntax:

projectscenedelete

projectscenedescription

Sets the description for the given row in the scenes table.

Syntax:

projectscenedescription *description*=⟨text⟩ ⟨row⟩

Options:

description

This is the text to set as the description for the given row.

Valid values: text strings

Default value:

Operands:

⟨row⟩

The row number of the scene to set the description for.

projectscenedragrows

Drags the selected rows in the scene table to the target row.

Syntax:

projectscenedragrows ⟨target row⟩

Operands:

⟨target row⟩

The index (row number) to drop the selected rows.

projectsceneexportpresentation

Exports the selected scenes as a presentation (project).

Syntax:

projectsceneexportpresentation \langle project path \rangle

Operands:

\langle project path \rangle

Project path is the full path to the new presentation, which is stored as a project.

projectsceneexportproject

Exports the current scene as a project.

Syntax:

projectsceneexportproject \langle project path \rangle

Operands:

\langle project path \rangle

Project path is the full path to the new project.

projectsceneextendselect

Extends the selection to the given row in the Scenes table.

Syntax:

projectsceneextendselect \langle row \rangle

Operands:

\langle row \rangle

The row number of the scene to select.

projectscenenew

Create a new scene.

Chapter 5: Commands

Syntax:

```
projectscenenew description=⟨text⟩ source=all | selected |  
included
```

Options:

description

A description for the new scene

Valid values: text strings

Default value:

source

Selects the source of entries.

Valid values: all
selected
included

Default value: **included**

projectsceneoptions

Holds options for project scenes.

Syntax:

```
projectsceneoptions thumbnailheight=⟨n⟩ thumbnailwidth=⟨n⟩  
warnbeforeclosing=yes | no
```

Options:

thumbnailheight

The height for saving thumbnails.

Valid values: integers

Default value: **100**

thumbnailwidth

The width for saving thumbnails.

Valid values: integers

Default value: **100**

warnbeforeclosing

If set, Maestro will put up a warning message before closing a modified scene.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

projectscenereplace

Replaces the current scene with the current Workspace and project contents.

Syntax:

```
projectscenereplace
```

projectsceneselect

Selects only the given row in the Scenes table.

Syntax:

```
projectsceneselect <row>
```

Operands:

<row>

The row number of the scene to select.

projectscenetoggleselect

Toggle-selects the given row in the Scenes table.

Syntax:

```
projectscenetoggleselect <row>
```

Operands:

<row>

The row number of the scene to select.

projectsceneview

View a scene.

Syntax:

projectsceneview <scene index>

Operands:

<scene index>

The index (row number) of the scene to view.

projectsceneviewnamed

View a specific scene in the table.

Syntax:

projectsceneviewnamed <scene>

Operands:

<scene>

Scene can either be: first, last, previous, or next to view the appropriate scene.

projectstoresnapshot

Save copy of current project state for reversion. Deprecated in favor of projectstorestate.

Syntax:

projectstoresnapshot

projectstorestate

Save copy of current project state for reversion.

Syntax:

projectstorestate

projectsynchronize

Save changes in workspace to current project.

Syntax:

`projectsynchronize`

projecttablefind

Searches for a string in the project table and make that cell editable.

Syntax:

```
projecttablefind direction=up | down findtext=<text>
                matchcase=yes | no matchword=yes | no
                searchhiddenrows=yes | no selectmatchentries=yes | no
                selectproperty=all | selected selectrow=all | selected
                showtoolbar=yes | no toolbararea=top | bottom | none
                <find_string>
```

Options:

direction With this option user can specify the direction of search.

Valid values: up
down

Default value: **down**

findtext Text to find.

Valid values: text strings

Default value:

matchcase With this option user can specify for case sensitive search.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

matchword

With this option user can specify for matching word as a whole.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

searchhiddenrows

With this option user can specify whether to search in collapsed group rows or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Chapter 5: Commands

selectmatchentries

With this option user can select entries with matches.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

selectproperty

With this option user can select properties to search

Valid values: all
 selected

Default value: **all**

selectrow

With this option user can select rows to search

Valid values: all
 selected

Default value: **all**

showtoolbar

With this option user can specify for showing the find toolbar or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

toolbararea

With this option user can specify the dock area of the find toolbar.

Valid values: top
 bottom
 none

Default value: **bottom**

Operands:

<find_string>

Its the string for which the user want to search in the project table.

projecttablereplace

Replace the <findstring> in the project table with the <replacestring> and finds the next <findstring>

Syntax:

projecttablereplace <find_string> <replace_string>

Operands:

<find_string> <replace_string>

Its the string to be searched Its the string to be replaced

projecttablereplaceall

Replace all the instances of findstring in the project table with the <replaces-string>

Syntax:

projecttablereplaceall <find_string> <replace_string>

Operands:

<find_string> <replace_string>

It's the string to be searched It's the string to be replaced

projectupdatecoordinates

Update the coordinates for the included entries in the project. This applies any current transformations in the Workspace to the original coordinates and places the result back into the project, overwriting the original coordinates.

Syntax:

projectupdatecoordinates

projectupdateviews

Update any open views on project data.

Syntax:

projectupdateviews

propertycalculate

This command calculates the given property for the for the entries that are specified with the ESL.

Syntax:

```
propertycalculate buried=yes | no gridspacing=<x>  
hydrophobic=yes | no proberadius=<x> propertyname=<text>  
recalculate=yes | no receptor=<text> surfacetype=vdw |  
extended | molecular | pbrsp vdwscale=<x> <ESL>
```

Options:

buried This option determines whether to calculate buried surface area.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

gridspacing

This is the grid spacing for generating molecular surface.

Valid values: reals
Default value: **0.6**
Minimum: 0.05

hydrophobic

This option determines whether to calculate hydrophobic surface area.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

proberadius

This is the probe radius for generating extended radius or molecular surface.

Valid values: reals
Default value: **1.4**

propertyname

This option is the name of the property that need to be calculated. The valid values for this option are numatoms , numresidues , nummolecules , molweight , spin , molcharge , secstruct and numheavyatoms , for calculating the properties Number of atoms , Number of residues , Number of molecules , Molecular weight , Spin multiplicity , Molecular charge , Secondary Structure Content and Number of heavy atoms respectively. A new property will be added to the project this <propertyname> and the values will be assigned to those entries which match the ESL.

	Valid values: text strings
	Default value: numatoms
<i>recalculate</i>	The valid values are true and false. With this option the user can specify whether to re-calculate the property data, which has been calculated already. This option will not have any impact on the entries which does not have the property data. If this option value is true, then the property value will be calculated and the new value will be assigned. If the options value is false, then the already existing property values will not be calculated.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>receptor</i>	The name of the receptor which is used for buried surface area calculation.
	Valid values: text strings
	Default value:
<i>surfacetyp</i>	This is the surface type for generating surface.
	Valid values: vdw extended molecular pbrsp
	Default value: extended
<i>vdwscaling</i>	This is the scaling of VdW radius for generating extended radius or molecular surface.
	Valid values: reals
	Default value: 1

Operands:

⟨ ESL ⟩

A valid ESL specification. Calculates specified property for those entries which match the ESL description.

propertycalculatechirality

This command calculates the chiral centers for the entries that are specified with the ESL.

Syntax:

propertycalculatedchirality *recalculate*=yes | no \langle ESL \rangle

Options:

recalculate The valid values are true and false. With this option the user can specify whether to re-calculate the property data, which has been calculated already. This option will not have any impact on the entries which do not have the property data. If this option value is true, then the property value will be calculated and the new value will be assigned. If the options value is false, then the already existing property values will not be calculated.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

\langle ESL \rangle

A valid ESL specification. Calculates the chiral centers for those entries which match the ESL description.

propertycalculatedmolformula

This command calculates the molecular formula for the entries that are specified with the ESL.

Syntax:

propertycalculatedmolformula *maxatoms*= \langle n \rangle
recalculate=yes | no \langle ESL \rangle

Options:

maxatoms This option is used for calculating molecular formula for an entry which has only one molecule. If the entry contains more than the specified maxatom then the molecular formula will not be calculated for that entry.

Valid values: integers

Default value: **100**

recalculate The valid values are true and false. With this option the user can specify whether to re-calculate the property data, which has been calculated already. This option will not have any impact on the entries which does not have the property data. If this

option value is true, then the property value will be calculated and the new value will be assigned. If the options value is false, then the already existing property values will not be calculated.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

⟨ESL⟩

A valid ESL specification. Calculates molecular formula for those entries which match the ESL description.

propertycalculatepickpka

This command calculates the pickPka atom for the entries that are specified with the ESL.

Syntax:

```
propertycalculatepickpka atomname=⟨text⟩
                        recalculate=yes | no ⟨ESL⟩
```

Options:

atomname This option is the name of the atom that is to be set as the pKa atom for the Jaguar calculation.

Valid values: text strings

Default value:

recalculate The valid values are true and false. With this option the user can specify whether to re-calculate the property data, which has been calculated already. This option will not have any impact on the entries which does not have the property data. If this option value is true, then the property value will be calculated and the new value will be assigned. If the options value is false, then the already existing property values will not be calculated.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

⟨ESL⟩

A valid ESL specification. Selects the pickPka atom for those entries which match the ESL description.

propertycalculatesubstructs

This command calculates the number of sub-structures for the entries that are specified with the ESL.

Syntax:

```
propertycalculatesubstructs name=<text>  
recalculate=yes | no substructure=<text> <ESL>
```

Options:

name This option holds the name of substructure. A valid property name is the valid value for this option. A new property of the integer type will be created with this <name>, which holds the count of specified sub-structures in an entry.

Valid values: text strings

Default value:

recalculate The valid values are true and false. With this option the user can specify whether to re-calculate the property data, which has been calculated already. This option will not have any impact on the entries which do not have the property data. If this option value is true, then the property value will be calculated and the new value will be assigned. If the options value is false, then the already existing property values will not be calculated.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

substructure

This describes the sub-structure definition. This value is used for matching substructures in each entry. The count of these sub-structures in each entry is saved as a property in the table.

Valid values: text strings

Default value:

Operands:

<ESL>

A valid ESL specification. Calculates the number of sub-structures for those entries which match the ESL description.

propertyclearvalue

This command clear values of the given properties

Syntax:

propertyclearvalue *allentries*=false | true <propertynames>

Options:

allentries The property value to be cleared for selected entries or all the entries. This has two valid values, true for clearing all property values and false for clearing property values for the selected entries only.

Valid values: false
true

Default value: **false**

Operands:

<propertynames>

The names of properties for which values need to be cleared.

propertycreate

Create a new property for entries in the current project without assigning any entry property values.

Syntax:

propertycreate *author*=<text> *displayprecision*=<n> *type*=bool | int | double | string <property_name>

Options:

author The m2io signature of the owner (e.g. person, group, or software module), or authority, that defines the property (meaning and range of values) and generally assigns its values. The author for user-defined properties should normally be user, as these properties can be freely edited in the project table. Be careful if setting author to m, mmod, i, j, qp, sd, or other names reserved for properties generated by existing programs.

Valid values: text strings
Default value: **user**

displayprecision

The Display precision for the real data type variable for display use only.

Valid values: integers

Default value: **4**
Minimum: 0
Maximum: 15

type The manner in which the property values are to be stored and represented, either bool (for Boolean), int (for integer), double (for double precision floating point, real numbers), or string (for text character strings).

Valid values: bool
int
double
string
Default value: **int**

Operands:

⟨property_name⟩

The user-assigned name of the property, which `propertyrename` can alter. The property name must be unique within the current project. If it is, a unique m2io data name will be generated by combining the type, author, and name, possibly modified (such as replacing spaces with underscores and adding a number at the end) to make it valid and unique.

propertydelete

This command deletes the given properties.

Syntax:

propertydelete ⟨propertynames⟩

Operands:

⟨propertynames⟩

The names of properties to delete.

Aliases:

deleteproperty (see [\[deleteproperty\]](#), page 109)

propertygeneratecontacts

This command applies a Contacts measurement as a property to a number of selected entries in a project. It is valid at least `Contactset1` has been set.

If only Contactset1 is defined, then Contactset2 = Contactset1. Normally these two sets are required having been set before this command, otherwise no property will be generated.

Syntax:

```
propertygeneratecontacts
```

propertygeneratehbond

This command applies an HBond measurement as a property to a number of selected entries in a project. It is valid at least HBondset1 has been set. If only HBondset1 is defined, then HBondset2 = HBondset1. Normally these two sets are required having been set before this command, otherwise no property will be generated.

Syntax:

```
propertygeneratehbond
```

propertymeasurementsetting

This command sets whether a measurement is applied as a property to a number of selected entries in a project.

Syntax:

```
propertymeasurementsetting applytoselectedentries=yes | no
```

Options:

applytoselectedentries

This option determines whether to apply a measurement as a property to a number of selected entries in a project.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

propertymove

Moves the property in the first column to the second column.

Syntax:

propertymove $\langle \text{from} \rangle \langle \text{to} \rangle \langle \text{table} \rangle$

Operands:

$\langle \text{from} \rangle \langle \text{to} \rangle \langle \text{table} \rangle$

The source column number. The destination column number. The name of the table to use as the source for the column numbers. If the table operand is missing, the current or default table (1) will be used.

propertyprecision

This command sets the precision of the property.

Syntax:

propertyprecision $\langle \text{property_name} \rangle \langle \text{precision} \rangle$

Operands:

$\langle \text{property_name} \rangle \langle \text{precision} \rangle$

The name of the property. Precision for that property.

propertyrename

This command renames the given property to the given name.

Syntax:

propertyrename $\langle \text{from} \rangle \langle \text{to} \rangle$

Operands:

$\langle \text{from} \rangle \langle \text{to} \rangle$

The name of the property to rename. The name to rename the property to.

Aliases:

renameproperty (see [\[renameproperty\]](#), page 670)

propertyshowall

This command creates a property subset consisting of all the properties in the project.

Syntax:

```
propertyshowall
```

propertysuperimposesetting

This command sets whether a superimposition is applied as a property to any entry in the project; and if applied, to which entries it is applied i.e. to the selected entries or to the included entries in the project.

Syntax:

```
propertysuperimposesetting applytoentries=included | selected
                           applytoincludedentries=yes | no createproperty=yes | no
```

Options:

applytoentries

This option determines to which entries the superimposition is applied as a property. Valid values are “selected”, or “included”.

Valid values: included
 selected

Default value: **included**

applytoincludedentries

This option determines whether to apply a superimposition as a property to a number of included entries in a project. NOTE: This option is deprecated now onwards and is supported just for backward compatibility. Internally the option will be interpreted as: true => applytoentries = included createproperty = true
false => createproperty = false

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

createproperty

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

protassign

This keyword is used to set various options associated with running protassign jobs.

Syntax:

```
protassign input_file=⟨text⟩ structure_source=selected_entries |  
workspace | file
```

Options:

input_file The name of the structure input file.

Valid values: text strings

Default value:

structure_source

Whether to use the selected entries in the current project, or what is in the workspace, or a specified file with multiple structures as structure input for the job.

Valid values: selected_entries
workspace
file

Default value: **workspace**

protassignresidues

Defines a set of atoms for protein assignment.

Syntax:

```
protassignresidues ⟨ASL⟩
```

Operands:

⟨ASL⟩

The ASL expression which defines the atoms that will be used to define the residues for protein assignment.

protassignstart

Start a protein assignment job with the current settings.

Syntax:

protassignstart

protassignwrite

Write a protassign input file with the current settings.

Syntax:

protassignwrite

proteinsculpting

Define an atom for a protein sculpting operation and/or translate it by X, Y, Z angstroms

Syntax:

proteinsculpting $\langle \text{atom} \rangle$ [$\langle X \rangle$ $\langle Y \rangle$ $\langle Z \rangle$]

Operands:

$\langle \text{atom} \rangle$ [$\langle X \rangle$ $\langle Y \rangle$ $\langle Z \rangle$]

An atom number to be transformed and the X, Y and Z extents for transformation.

psp_helical

Set residues to define a Prime refinement helical constraint.

Syntax:

psp_helical *constrain*=yes | no *helixend*= $\langle \text{text} \rangle$
helixstart= $\langle \text{text} \rangle$ *loopend*= $\langle \text{text} \rangle$ *loopresidues*= $\langle \text{text} \rangle$
loopstart= $\langle \text{text} \rangle$

Options:

<i>constrain</i>	Whether or not to use the constraint in loop prediction jobs.
Valid values:	boolean (true false; yes no; y n; on off)
Default value:	false

Chapter 5: Commands

<i>helixend</i>	The residue where the helix ends. Valid values: text strings Default value:
<i>helixstart</i>	The residue where the helix starts. Valid values: text strings Default value:
<i>loopend</i>	The residue where the loop ends. Valid values: text strings Default value:
<i>loopresidues</i>	All residue of the loop. Valid values: text strings Default value:
<i>loopstart</i>	The residue where the loop starts. Valid values: text strings Default value:

psp_pairwise

Specifies a pair of atoms to define a Prime refinement pairwise constraint.

Syntax:

psp_pairwise *distance*= $\langle x \rangle$ *forceconstant*= $\langle n \rangle$ $\langle \text{atom1} \rangle$ $\langle \text{atom2} \rangle$

Options:

<i>distance</i>	The distance from atom1 to atom2 of the pairwise constraint. Valid values: reals Default value: 0 Minimum: 0.0
<i>forceconstant</i>	The force constant for the pairwise constraint. Valid values: integers Default value: 350 Minimum: 0

Operands:

$\langle \text{atom1} \rangle$ $\langle \text{atom2} \rangle$

The two atoms that define a Prime refinement pairwise constraint.

psp_spatial

Specifies a residue to define a Prime refinement spatial constraint.

Syntax:

```
psp_spatial caatom=<n> constrain=yes | no distance=<x>
           xcoord=<x> xcoordold=<x> ycoord=<x> ycoordold=<x>
           zcoord=<x> zcoordold=<x>
           <chain_name>:<residue_number><insertion_code>
```

Options:

<i>caatom</i>	The atom number of C-alpha atom of the residue. Valid values: integers Default value: 0 Minimum: 0
<i>constrain</i>	Whether or not to use the constraint in loop prediction jobs. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>distance</i>	The radius of sphere centered at the C-alpha atom of the residue. Valid values: reals Default value: 3 Minimum: 0.0
<i>xcoord</i>	The X coordinate of the C-alpha atom of the residue. Valid values: reals Default value: 0
<i>xcoordold</i>	The X coordoldinate of the C-alpha atom of the residue. Valid values: reals Default value: 0
<i>ycoord</i>	The Y coordinate of the C-alpha atom of the residue. Valid values: reals Default value: 0
<i>ycoordold</i>	The Y coordoldinate of the C-alpha atom of the residue. Valid values: reals Default value: 0

<i>zcoord</i>	The Z coordinate of the C-alpha atom of the residue. Valid values: reals Default value: 0
<i>zcoordold</i>	The Z coordoldinate of the C-alpha atom of the residue. Valid values: reals Default value: 0

Operands:

`<chain_name>:<residue_number><insertion_code>`

The residue that define a Prime refinement spatial constraint.

pspalignaddanchor

Adds an anchor at the given position for Edit Alignment.

Syntax:

pspalignaddanchor `<anchor position>`

Operands:

`<anchor position>`

The position at which to add an anchor.

pspaligndeleteanchor

Deletes the anchor at the given position for Edit Alignment.

Syntax:

pspaligndeleteanchor `<anchor position>`

Operands:

`<anchor position>`

The position at which to delete an anchor.

pspaligninsertgaps

Inserts gaps into the sequence at the given position.

Syntax:

pspaligninsertgaps *number*=⟨n⟩ *position*=⟨n⟩ ⟨sequence name⟩

Options:

<i>number</i>	This is the number of gaps to insert.
	Valid values: integers
	Default value: 1
	Minimum: 1
<i>position</i>	This is the position at which to insert a gap.
	Valid values: integers
	Default value: 1
	Minimum: 1

Operands:

⟨sequence name⟩

The name of the sequence to insert gaps into.

pspalignlockgaps

Locks gaps for all alignments in the Edit Alignment step.

Syntax:

pspalignlockgaps

pspalignmoveleft

Moves the alignment left freely starting at the given position. This will shift gaps from the left of the given position to the right of the given position.

Syntax:

pspalignmoveleft *number*=⟨n⟩ *position*=⟨n⟩ ⟨sequence name⟩

Options:

<i>number</i>	This is the number of spaces to move left.
	Valid values: integers
	Default value: 1
	Minimum: 1

position This is the position at which to move left.

Valid values: integers

Default value: **1**

Minimum: 1

Operands:

⟨sequence name⟩

The name of the sequence to move left.

pspalignmoveleftblock

Moves the alignment left as a block starting at the given position. This will close up gaps to the left of the given position while preserving them to the right.

Syntax:

```
pspalignmoveleftblock number=⟨n⟩ position=⟨n⟩  
⟨sequence name⟩
```

Options:

number This is the number of spaces to move left.

Valid values: integers

Default value: **1**

Minimum: 1

position This is the position at which to move left.

Valid values: integers

Default value: **1**

Minimum: 1

Operands:

⟨sequence name⟩

The name of the sequence to move left.

pspalignmoveright

Moves the alignment right freely starting at the given position. This will shift gaps from the right of the given position to the left of the given position. If necessary, this will open up gaps to the left of the given position.

Syntax:

pspalignmoveright *number*=⟨n⟩ *position*=⟨n⟩ ⟨sequence name⟩

Options:

number This is the number of spaces to move right.

Valid values: integers

Default value: **1**

Minimum: 1

position This is the position at which to move right.

Valid values: integers

Default value: **1**

Minimum: 1

Operands:

⟨sequence name⟩

The name of the sequence to move right.

pspalignstructures

Runs 'structalign' to align the selected template structures in the Find Homologs step.

Syntax:

pspalignstructures

pspalignunlockgaps

Unlocks gaps for all alignments in the Edit Alignment step.

Syntax:

pspalignunlockgaps

pspbstogglehetatom

Toggles the given hetatom between included and excluded.

Syntax:

pspbstogglehetatom \langle hetatom name \rangle

Operands:

\langle hetatom name \rangle

The name of the hetatom to toggle.

pspbuildbackbone

Runs the build backbone backend.

Syntax:

pspbuildbackbone

pspbuildconsensusmodel

Runs the consensus homology modeling backend.

Syntax:

pspbuildconsensusmodel

pspbuildheteromultimer

Runs the hetero-multimer building command.

Syntax:

pspbuildheteromultimer

pspbuildhomomultimer

Runs the build structure backend to build a homo-multimer model.

Syntax:

pspbuiidhomomultimer

pspbuiidstructure

Runs the build structure backend.

Syntax:

pspbuiidstructure

pspconstraintmarkersettings

Set graphical data of Phase excluded volume markers.

Syntax:

```
pspconstraintmarkersettings helicalblue=⟨ x ⟩ helicalgreen=⟨ x ⟩
helicalred=⟨ x ⟩ pairwiseblue=⟨ x ⟩ pairwisegreen=⟨ x ⟩
pairwisehighlightwidth=⟨ x ⟩ pairwisered=⟨ x ⟩
spatialambient=⟨ x ⟩ spatialblue=⟨ x ⟩ spatialdiffuse=⟨ x ⟩
spatialemission=⟨ x ⟩ spatialgreen=⟨ x ⟩ spatialhighlightblue=⟨ x ⟩
spatialhighlightgreen=⟨ x ⟩ spatialhighlightred=⟨ x ⟩
spatiallineblue=⟨ x ⟩ spatiallinegreen=⟨ x ⟩ spatiallinered=⟨ x ⟩
spatiallinewidth=⟨ x ⟩ spatialred=⟨ x ⟩ spatialshininess=⟨ x ⟩
spatialslices=⟨ n ⟩ spatialspecular=⟨ x ⟩ spatialstacks=⟨ n ⟩
spatialtransparency=⟨ x ⟩
```

Options:

helicalblue The blue color component of helical markers.

Valid values:	reals
Default value:	1
Minimum:	0.0
Maximum:	1.0

helicalgreen

The green color component of helical markers.

Valid values:	reals
Default value:	0
Minimum:	0.0
Maximum:	1.0

helicalred The red color component of helical markers.

Chapter 5: Commands

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

pairwiseblue

The blue color component of pairwise markers.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

pairwisegreen

The green color component of pairwise markers.

Valid values: reals
Default value: **0**
Minimum: 0.0
Maximum: 1.0

pairwisehighlightwidth

The highlight line width of pairwise markers.

Valid values: reals
Default value: **1**
Minimum: 0.0

pairwisered

The red color component of pairwise markers.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

spatialambient

Set material property - ambient, to its red, green, and blue components, for front face.

Valid values: reals
Default value: **0.4**
Minimum: 0.0
Maximum: 1.0

spatialblue

The blue color component of spatial markers.

Valid values: reals
Default value: **0**
Minimum: 0.0
Maximum: 1.0

spatialdiffuse

Set material property - diffuse, to its red, green, and blue components, for front face.

Valid values: reals
 Default value: **0.4**
 Minimum: 0.0
 Maximum: 1.0

spatialemission

Set material property - emission, to its red, green, and blue components, for front face.

Valid values: reals
 Default value: **0.1**
 Minimum: 0.0
 Maximum: 1.0

spatialgreen

The green color component of spatial markers.

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 1.0

spatialhighlightblue

The blue color component of spatial markers if the spatial is highlighted.

Valid values: reals
 Default value: **0**
 Minimum: 0.0
 Maximum: 1.0

spatialhighlightgreen

The green color component of spatial markers if the spatial is highlighted.

Valid values: reals
 Default value: **0**
 Minimum: 0.0
 Maximum: 1.0

spatialhighlightred

The red color component of spatial markers if the spatial is highlighted.

Valid values: reals
 Default value: **1**
 Minimum: 0.0
 Maximum: 1.0

Chapter 5: Commands

spatiallineblue

The blue color component of spatial line markers.

Valid values: reals
Default value: **0**
Minimum: 0.0
Maximum: 1.0

spatiallinegreen

The green color component of spatial line markers.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

spatiallinered

The red color component of spatial line markers.

Valid values: reals
Default value: **0**
Minimum: 0.0
Maximum: 1.0

spatiallinewidth

The red color component of spatial line markers.

Valid values: reals
Default value: **1**
Minimum: 0.0

spatialred

The red color component of spatial markers.

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

spatialshininess

Set material property - shininess, for front face.

Valid values: reals
Default value: **80**
Minimum: 0.0
Maximum: 128.0

spatialslices

Set the slices of drawing sphere.

Valid values: integers
Default value: **18**
Minimum: 2

spatialspecular

Set material property - specular, to its red, green, and blue components, for front face.

Valid values: reals
 Default value: **0.1**
 Minimum: 0.0
 Maximum: 1.0

spatialstacks

Set the stacks of drawing sphere.

Valid values: integers
 Default value: **9**
 Minimum: 2

spatialtransparency

The transparency of QSAR markers.

Valid values: reals
 Default value: **20**
 Minimum: 0.0
 Maximum: 100.0

pspeditquerysequence

Edits query sequence in Input Sequence step.

Syntax:

pspeditquerysequence

pspexcludetable1row

Excludes a composite structure from the table shown in the Build Backbone step. Excludes a structure, for the specified row in the input table shown in the current Structure Prediction step, from the Sequence Viewer (and Workspace). This applies to the Refine Backbone step.

Syntax:

pspexcludetable1row <row_number>

Operands:

<row_number>

The row number in the table which is to be excluded.

pspexcludetablerow

Excludes a composite structure from the table shown in the Build Backbone step. Excludes a structure, for the specified row in the table shown in the current Structure Prediction step, from the Sequence Viewer (and Workspace). This applies to the Fold Recognition, Build Backbone, Refine Backbone, and Refine Structure steps (which have structure or template tables, and the ability to specify rows in Workspace independently from selected rows).

Syntax:

pspexcludetablerow \langle row_number \rangle

Operands:

\langle row_number \rangle

The row number in the table which is to be excluded.

pspexportalignment

exports alignments selected into a file.

Syntax:

pspexportalignment \langle filename \rangle

Operands:

\langle filename \rangle

The name of the file to which the alignments are to be written.

pspfindfamily

Finds the family for the query sequence currently active in the sequence viewer.

Syntax:

pspfindfamily

pspfindhomologs

Runs Blast to find homologs for the query sequence currently active in the sequence viewer.

Syntax:

pspfindhomologs

pspfoldrecognitionoptions

Sets some options associated with the Fold Recognition step.

Syntax:

pspfoldrecognitionoptions *max_results*= $\langle n \rangle$

Options:

max_results

Specifies how many results will be shown in the fold recognition results table.

Valid values: integers

Default value: **100**

pspfoldrecognitionsearch

Runs the fold recognition search on the current query sequence.

Syntax:

pspfoldrecognitionsearch

pspgetquerysequence

Creates new query sequence from text in Input Sequence step.

Syntax:

`pspgetquerysequence`

pspimportalignment

imports alignments in a file and applies to the templates. since there can be more than one alignment stored in a file, the command gives messages on which templates have been altered.

Syntax:

`pspimportalignment` \langle filename \rangle

Operands:

\langle filename \rangle

The name of the file from which the alignments are to be read.

pspimporthomolog

Imports the sequences from a file as homologs to the current query sequence.

Syntax:

`pspimporthomolog` \langle filename \rangle

Operands:

\langle filename \rangle

The name of the file from which the sequences are to be read.

pspimportssp

Adds the secondary structure prediction from a file to those for the current query sequence.

Syntax:

pspimportssp \langle filename \rangle

Operands:

\langle filename \rangle

The name of the file from which the SSP is to be read.

pspincludetable1row

Includes a structure from the input table shown in the current Structure Prediction step into the Sequence Viewer (and Workspace), independent of selection. This applies to the Refine Backbone step.

Syntax:

pspincludetable1row \langle row_number \rangle

Operands:

\langle row_number \rangle

The row number in the table which is to be included.

pspincludetablerow

Includes a structure from the table shown in the current Structure Prediction step into the Sequence Viewer (and Workspace). This applies to the Fold Recognition, Build Backbone, Refine Backbone, and Refine Structure steps (which have structure or template tables, and the ability to specify rows in Workspace independently from selected rows).

Syntax:

pspincludetablerow \langle row_number \rangle

Operands:

\langle row_number \rangle

The row number in the table which is to be included.

pspminimizationresidues

Defines a set of atoms for minimization refinement.

Syntax:

pspminimizationresidues \langle ASL \rangle

Operands:

\langle ASL \rangle

The ASL expression which defines the atoms that will be used to define the residues for minimization refinement.

pspoptimizealignment

optimizes alignment of a template sequence.

Syntax:

pspoptimizealignment \langle rowindex \rangle

Operands:

\langle rowindex \rangle

The row index of the sequence which is aligned

pspquickbuild

Runs the build structure backend with fixed options.

Syntax:

pspquickbuild

pspquickhomomultimer

Runs the build homo-multimer backend with fixed options.

Syntax:

pspquickhomomultimer

psprefinebackbone

Runs the refine backbone backend.

Syntax:

```
psprefinebackbone
```

psprefinestructure

Runs the structure refinement program with the currently set refinement options on the currently selected structure.

Syntax:

```
psprefinestructure
```

psprsdefaulthelixloops

Rebuilds the Helices & Loops table for structure refinement with the default loop and helix refinement options for the currently selected structure.

Syntax:

```
psprsdefaulthelixloops
```

psprshelixlooptogglerefine

Toggle refinement of loop or helix.

Syntax:

```
psprshelixlooptogglerefine <row_index>
```

Operands:

<row_index>

The row index for the loop or helix in the Helices & Loops table that is to be toggled for refinement.

psprshelixlooptogglerefineonly

Toggle refinement of loop or helix.

Syntax:

```
psprshelixlooptogglerefineonly <row_index>
```

Operands:

<row_index>

The row index for the loop or helix in the Helices & Loops table that is the only one to be toggled for refinement.

psprshelixoptions

Specifies the refinement options for the currently selected helix in the Helices & Loops table.

Syntax:

```
psprshelixoptions roll_range=<x> roll_resolution=<x>  
                  set_roll=yes | no set_shift=yes | no set_tilt=yes | no  
                  set_trans=yes | no shift_range=<x> shift_resolution=<x>  
                  tilt_range=<x> tilt_resolution=<x> trans_range=<x>  
                  trans_resolution=<x>
```

Options:

roll_range Range, in degrees, for helix roll.

Valid values: reals

Default value: **20**

roll_resolution

Resolution, in degrees, for helix roll.

Valid values: reals

Default value: **20**

set_roll Set the helix roll range and resolution.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

set_shift Set the helix shift range and resolution.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

<i>set_tilt</i>	Set the helix tilt range and resolution. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>set_trans</i>	Set the helix translation range and resolution. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>shift_range</i>	Range, in angstroms, for helix shift. Valid values: reals Default value: 2
<i>shift_resolution</i>	Resolution, in angstroms, for helix shift. Valid values: reals Default value: 2
<i>tilt_range</i>	Range, in degrees, for helix tilt. Valid values: reals Default value: 5
<i>tilt_resolution</i>	Resolution, in degrees, for helix tilt. Valid values: reals Default value: 5
<i>trans_range</i>	Range, in angstroms, for helix translation. Valid values: reals Default value: 2
<i>trans_resolution</i>	Resolution, in angstroms, for helix translation. Valid values: reals Default value: 2

psprsloopoptions

Specifies the refinement options for the currently selected loop in the Helices & Loops table.

Syntax:

```
psprsl options ca_max=<x> membrane=none | inside |
               outside overlap_min=<x> restrict_ca=yes | no
               sphere_size=<x> use_sphere=yes | no
```

Options:

ca_max Maximum CA atom movement, in angstroms, from initial.

Valid values: reals

Default value: **3**

membrane Membrane constraint type for loop building

Valid values: none
inside
outside

Default value: **none**

overlap_min

Minimum overlap, in angstroms, from initial.

Valid values: reals

Default value: **0.7**

Minimum: 0.1

Maximum: 1.0

restrict_ca Restrict loop CA atom movement.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

sphere_size

Distance from loop for sidechain inclusion, in angstroms.

Valid values: reals

Default value: **7.5**

use_sphere

Include nearby sidechains in loop refinement.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

psprsnontemplatehelixloops

Rebuilds the Helices & Loops table for structure refinement with the loops that were not constructed from a template during homology modeling, for the currently selected structure.

Syntax:

psprsnontemplatehelixloops

psprunalign

Runs secondary structure prediction (if necessary) and then aligns the query to the active templates in Edit Alignment.

Syntax:

psprunalign *gpcr*=yes | no

Options:

gpcr Specifies whether or not to use GPCR alignment.
 Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

pspruncreate

Creates the run with the given name.

Syntax:

pspruncreate <run name>

Operands:

<run name>

The name of the new run to create.

psprundelele

Deletes the current run from the project.

Syntax:

psprundelele

psprunopen

Opens the run with the given name.

Syntax:

psprunopen ⟨run name⟩

Operands:

⟨run name⟩

The name of the run to open.

psprunquickalign

Runs quick multiple sequence alignment algorithm and to align the query to the active templates in Edit Alignment.

Syntax:

psprunquickalign

psprunrename

Changes the current run's name to the given name.

Syntax:

psprunrename ⟨run name⟩

Operands:

⟨run name⟩

The name to change the current run's name to.

psprunsaveas

Saves a copy of the current run under the given name.

Syntax:

psprunsaveas \langle run name \rangle

Operands:

\langle run name \rangle

The name of the run to save as.

pspsecstructprediction

Runs the available secondary structure prediction programs on the current query sequence.

Syntax:

pspsecstructprediction

pspselectextendtable1row

Extend the selection from the selected table row to joining up with an existing selection. This applies to the Refine Backbone step

Syntax:

pspselectextendtable1row \langle row_number \rangle

Operands:

\langle row_number \rangle

The row number in the table from which the selection is to begin.

pspselectextendablerow

Extend the selection from the selected table row to joining up with an existing selection. This applies to the Find Homologs, Fold Recognition, Build Backbone, Refine Backbone, Edit Alignment and Refine Structure steps.

Syntax:

pspselectextendablerow $\langle \text{row_number} \rangle$

Operands:

$\langle \text{row_number} \rangle$

The row number in the table from which the selection is to begin.

pspselecthelixlooprow

Select structure for refinement.

Syntax:

pspselecthelixlooprow $\langle \text{row_index} \rangle$

Operands:

$\langle \text{row_index} \rangle$

The row index for the loop or helix in the Helices & Loops table that is to be selected for refinement option editing.

pspselectonlytable1row

Select a row from the input structure or template table shown in the current Structure Prediction step, deselecting all other rows. This applies to the Refine Backbone step.

Syntax:

pspselectonlytable1row $\langle \text{row_number} \rangle$

Operands:

$\langle \text{row_number} \rangle$

The row number in the table which is to be selected.

pspselectonlyablerow

Select a row from the structure or template table shown in the current Structure Prediction step, deselecting all other rows. This applies to the Find Homologs, Fold Recognition, Build Backbone, Refine Backbone, Edit Alignment and Refine Structure steps.

Syntax:

pspselectonlytablerow \langle row_number \rangle

Operands:

\langle row_number \rangle

The row number in the table which is to be selected.

pspselectrscontext

Specify the context for structure refinement. This command may be executed instead of showing the psp or refinement panels. Be cautious about using this commands while either of these panels is shown, as setting the wrong context will interfere with the execution of the refinement commands issued from the panel.

Syntax:

pspselectrscontext pspstep|standalone

Operands:

pspstep|standalone

The context in which refine structure commands are executed.

pspselectrsrefinement

Selects the current psp structure refinement type.

Syntax:

pspselectrsrefinement \langle refine_type \rangle

Operands:

\langle refine_type \rangle

The type of structure refinement to be performed. Value should be sidechains , helixloops , or minimization .

pspselecttable1row

Selects a row from the input structure or template table shown in the current Structure Prediction step. This applies to the Refine Backbone step.

Syntax:

```
pspselecttable1row <row_number>
```

Operands:

<row_number>

The row number in the table which is to be selected.

pspselecttablerow

Selects a row from the structure or template table shown in the current Structure Prediction step. This applies to the Find Homologs, Fold Recognition, Build Backbone, Refine Backbone, Edit Alignment and Refine Structure steps.

Syntax:

```
pspselecttablerow <row_number>
```

Operands:

<row_number>

The row number in the table which is to be selected.

pspsequenceaddfile

Adds the sequences from the given file to the Select Sequence step in PSP.

Syntax:

```
pspsequenceaddfile <file>
```

Operands:

<file>

The file name of the sequence file to add.

pspsequenceaddworkspace

Adds the sequences from the Workspace to the Select Sequence step in PSP.

Syntax:

pspsequenceaddworkspace

pspsequencecrop

Crops the given sequence to the given residue index, towards the closer end of the sequence.

Syntax:

pspsequencecrop *res*=⟨n⟩ ⟨sequence name⟩ ⟨residue index⟩

Options:

<i>res</i>	This is the position to crop to
Valid values:	integers
Default value:	1
Minimum:	1

Operands:

⟨sequence name⟩ ⟨residue index⟩

The name of the sequence to crop.

pspsequenceselect

Selects the given sequence number as input for the next step in PSP.

Syntax:

pspsequenceselect ⟨sequence number⟩

Operands:

⟨sequence number⟩

The index of the sequence to select.

pspsequenceviewerexport

Exports all of the visible sequences in the sequence viewer to the given file.

Syntax:

```
pspsequenceviewerexport sequenceviewer=<text> <file name>
```

Options:

sequenceviewer

Specifies whether to save the contents of the Workspace or Prime sequence viewer.

Valid values: text strings

Default value: **prime**

Operands:

<file name>

The name of the file to export the sequences to.

pspsethelixloopresidues

Specify residues to define loop or helix feature for structure refinement.

Syntax:

```
pspsethelixloopresidues <row_index> <R1> <R2> <R3> <R4>
```

Operands:

<row_index> <R1> <R2> <R3> <R4>

The row index for the loop or helix in the Helices & Loops table, followed by the strings for the four residue column values, as they are to appear in the table. For a loop, only two end residues (<R1> <R2>) are required, so the <R3> and <R4> values can be missing or left blank. For a helix, <R1> and <R4> specify the outer residues for the loops at the two ends of the helix, while <R2> and <R3> specify the ends of the helix itself. The format for the residue strings is the residue number followed by the insertion code, if any. For example, 12C represent residue number 12, insertion code C.

pspsidechainresidues

Defines a set of atoms for the side chain refinement.

Syntax:

```
pspsidechainresidues prediction_algorithm=default |  
vectorsampling | backbonesampling <ASL>
```

Options:

prediction_algorithm

Determines which side chain sampling algorithm is used for Prime sidechain prediction jobs.

Valid values: default
 vectorsampling
 backbonesampling

Default value: **default**

Operands:

<ASL>

The ASL expression which defines the atoms that will be used to define the residues for side chain refinement.

pspsortbbtable

Resort the Build Backbone table based on the data in the specified column

Syntax:

```
pspsortbbtable <column_name>
```

Operands:

<column_name>

The name of the column to be sorted.

pspsortfindhomologstable

Resort the find homologs table based on the data in the specified column

Syntax:

pspsortfindhomologstable $\langle \text{column_name} \rangle$

Operands:

$\langle \text{column_name} \rangle$

The name of the column to be sorted.

pspsortfoldtable

Resort the fold recognition table based on the data in the specified column

Syntax:

pspsortfoldtable $\langle \text{column_name} \rangle$

Operands:

$\langle \text{column_name} \rangle$

The name of the column to be sorted.

pspsortrbtable

Resort the Refine Backbone output structure table based on the data in the specified column

Syntax:

pspsortrbtable $\langle \text{column_name} \rangle$

Operands:

$\langle \text{column_name} \rangle$

The name of the column to be sorted.

pspsortrbtable1

Resort the Refine Backbone composite structure table based on the data in the specified column

Syntax:

pspsortrbtable1 \langle column_name \rangle

Operands:

\langle column_name \rangle

The name of the column to be sorted.

pspsortrstable

Resort the Refine Structure table based on the data in the specified column

Syntax:

pspsortrstable \langle column_name \rangle

Operands:

\langle column_name \rangle

The name of the column to be sorted.

pspsspdelete

Deletes the given secondary structure prediction.

Syntax:

pspsspdelete \langle SSP name \rangle

Operands:

\langle SSP name \rangle

The name of the secondary structure prediction to delete.

pspsspexport

Exports the given secondary structure to the named file.

Syntax:

pspsspexport *ssp*=⟨text⟩ ⟨file name⟩

Options:

ssp This option is the name of the secondary structure prediction to export.

Valid values: text strings

Default value:

Operands:

⟨file name⟩

The name of the file to export to.

pspssprevert

Reverts a modified secondary structure prediction back to its original form.

Syntax:

pspssprevert ⟨SSP name⟩

Operands:

⟨SSP name⟩

The name of the secondary structure prediction to revert.

pspsspset

Sets the given range of positions in the secondary structure position to the given code, which must be H , E , or - .

Syntax:

pspsspset *code*=⟨text⟩ *from*=⟨n⟩ *to*=⟨n⟩ ⟨SSP name⟩

Options:

code This is the code for setting the secondary structure position to. It must be one of H , E , or - .

Valid values: text strings

Default value:

<i>from</i>	This is the starting position for setting codes.
	Valid values: integers
	Default value: 1
	Minimum: 1
<i>to</i>	This is the ending position for setting codes.
	Valid values: integers
	Default value: 1
	Minimum: 1

Operands:

⟨SSP name⟩

The name of the step to set data for.

pspstepforward

Moves forward to the next named step. Deletes any steps after the current step and copies data forward to create the specified step.

Syntax:

pspstepforward ⟨step name⟩

Operands:

⟨step name⟩

The name of the step to switch to.

pspstepgoto

Moves to an existing step in the current run.

Syntax:

pspstepgoto ⟨step name⟩

Operands:

⟨step name⟩

The name of the step to switch to.

pspstructureaddentry

Create project entry from selected Prime structure. This operation is only permitted for the Build Structure step.

Syntax:

pspstructureaddentry

psptemplatesetregion

Sets the region for the given template.

Syntax:

psptemplatesetregion *end*=⟨n⟩ *start*=⟨n⟩ ⟨template name⟩

Options:

end This is the position at which the region ends.

Valid values: integers

Default value: **1**

Minimum: 1

start This is the position at which to start the region.

Valid values: integers

Default value: **1**

Minimum: 1

Operands:

⟨template name⟩

The name of the template to set a region for.

psptmplsecstructprediction

Runs the available secondary structure prediction programs on the template sequence identified by row number

Syntax:

psptmplsecstructprediction $\langle \text{row_number} \rangle$

Operands:

$\langle \text{row_number} \rangle$

The row number in the table which is to be selected.

psptranslatespatial

Translate in Angstroms the selected Prime refinement spatial constraint.

Syntax:

psptranslatespatial $x=\langle x \rangle$ $y=\langle x \rangle$ $z=\langle x \rangle$

Options:

x Amount in Angstroms to translate in X

Valid values: reals

Default value: **0**

y Amount in Angstroms to translate in Y

Valid values: reals

Default value: **0**

z Amount in Angstroms to translate in Z

Valid values: reals

Default value: **0**

pspunselecttable1row

Unselects a row from the input table shown in the current Structure Prediction step. This applies to the Refine Backbone step.

Syntax:

pspunselecttable1row $\langle \text{row_number} \rangle$

Operands:

$\langle \text{row_number} \rangle$

The row number in the table which is to be unselected.

pspunselecttablerow

Unselects a row from the table shown in the current Structure Prediction step. This applies to the Find Homologs, Fold Recognition, Build Backbone, Refine Backbone, and Edit Alignment steps (which support multiple row selection).

Syntax:

```
pspunselecttablerow <row_number>
```

Operands:

<row_number>

The row number in the table which is to be unselected.

pspupdatehetatomlist

Updates list of ligands and cofactors in Build Structure step.

Syntax:

```
pspupdatehetatomlist showall=yes | no
```

Options:

showall Specifies whether or not to list all ligands and waters, or only ligands within 5Å from a chain.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

pspupdatescores

Updates Scores of all the alignments

Syntax:

```
pspupdatescores
```


pspwritebuildconsensusmodel

Writes input files for the build structure backend in consensus modeling mode.

Syntax:

```
pspwritebuildconsensusmodel
```

pspwritebuildhomomultimer

Writes input files for the build structure backend in homo-multimer mode.

Syntax:

```
pspwritebuildhomomultimer
```

pspwritebuildstructure

Writes input files for the build structure backend.

Syntax:

```
pspwritebuildstructure
```

pyeval

This is a standard alias for **pythoneval** (see [\[pythoneval\]](#), page 648).

pyimp

This is a standard alias for **pythonimport** (see [\[pythonimport\]](#), page 648).

pyrun

This is a standard alias for **pythonrun** (see [\[pythonrun\]](#), page 648).

pythoneval

Evaluate the python expression using the built-in Python interpreter.
KEY_OPERAND_SYNOPOSIS: <python expression>

Syntax:

pythoneval

Aliases:

pyeval (see [\[pyeval\]](#), page 647)

pythonimport

Imports the specified python module. If the module has already been imported then it will be reloaded. Under normal usage this command is not needed as pythonrun automatically imports the module specified as its operand. However during development of modules it is useful to be able to reload them using pythonimport
KEY_OPERAND_SYNOPOSIS: <module>

Syntax:

pythonimport

Aliases:

pyimp (see [\[pyimp\]](#), page 647)

pythonrun

Runs the function in <function name> in module <modulename>
KEY_OPERAND_SYNOPOSIS: <function name> [function parameters]

Syntax:

pythonrun

Aliases:

pyrun (see [\[pyrun\]](#), page 647)

pythonrunbuiltin

Runs the function in <function name> in module <modulename>
KEY_OPERAND_SYNOPOSIS: <function name> [function parameters]

Syntax:

`pythonrunbuiltin`

qikprop

This keyword is used to set various options associated with running QikProp jobs.

Syntax:

`qikprop` *fastmode*=yes | no *input_file*=⟨text⟩ *nsim*=⟨n⟩
sim=yes | no *structure_source*=selected_entries | workspace |
 file

Options:

<i>fastmode</i>	The option to identify the similar molecules Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>input_file</i>	The name of the structure input file. Valid values: text strings Default value:
<i>nsim</i>	 Valid values: integers Default value: 5 Minimum: 1 Maximum: 999
<i>sim</i>	The number of similar molecules to identify Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>structure_source</i>	Whether to use the selected entries in the current project or what is in the workspace, or a specified file with multiple structures as input for the job. Valid values: selected_entries workspace file Default value: selected_entries

qsarmarkerdump

Print out the current option values of the QSAR marker command.

Syntax:

```
qsarmarkerdump
```

qsarmarkersettings

Set graphical data of Phase QSAR markers.

Syntax:

```
qsarmarkersettings ambient=<x> coefficient_feature=<text>  
  combine_effects=yes | no diffuse=<x> emission=<x>  
  negative_blue=<x> negative_green=<x> negative_red=<x>  
  negativecoefficient=<x> numberpls=<n> positive_blue=<x>  
  positive_green=<x> positive_red=<x> positivecoefficient=<x>  
  roundingeffect=<x> selectedatomclass=<text> shininess=<x>  
  specular=<x> step=<n> transparency=<x>  
  volumeoccupied=workspacelig | qsarmodel
```

Options:

ambient Set material property - ambient, to its red, green, and blue components, for front face.

Valid values: reals
Default value: **0.4**
Minimum: 0.0
Maximum: 1.0

coefficient_feature

Sets which feature the coefficient options apply to. Can also be set to combined .

Valid values: text strings
Default value: **combined**

combine_effects

Whether or not to combine the effects from all classes.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

diffuse

Set material property - diffuse, to its red, green, and blue components, for front face.

	Valid values:	reals
	Default value:	0.4
	Minimum:	0.0
	Maximum:	1.0
<i>emission</i>	Set material property - emission, to its red, green, and blue components, for front face.	
	Valid values:	reals
	Default value:	0.1
	Minimum:	0.0
	Maximum:	1.0
<i>negative_blue</i>	Blue component of negative effects color.	
	Valid values:	reals
	Default value:	0
	Minimum:	0.0
	Maximum:	1.0
<i>negative_green</i>	Green component of negative effects color.	
	Valid values:	reals
	Default value:	0
	Minimum:	0.0
	Maximum:	1.0
<i>negative_red</i>	Red component negative effects color.	
	Valid values:	reals
	Default value:	1
	Minimum:	0.0
	Maximum:	1.0
<i>negativecoefficient</i>	Set the QSAR visualization option of negative coefficient threshold.	
	Valid values:	reals
	Default value:	-0.044
	Maximum:	0.0
<i>numberpls</i>	Set the QSAR visualization option of number of PLS factors.	
	Valid values:	integers
	Default value:	1
	Minimum:	1
<i>positive_blue</i>	Blue component of positive effects color.	

Chapter 5: Commands

Valid values: reals
Default value: **1**
Minimum: 0.0
Maximum: 1.0

positive_green

Green component of positive effects color.

Valid values: reals
Default value: **0**
Minimum: 0.0
Maximum: 1.0

positive_red

Red component of positive effects color.

Valid values: reals
Default value: **0**
Minimum: 0.0
Maximum: 1.0

positivecoefficient

Set the QSAR visualization option of positiv coefficient threshold.

Valid values: reals
Default value: **0.044**
Minimum: 0.0

roundingeffect

This determines the rounding effect of edges of a cell, for front face.

Valid values: reals
Default value: **5**
Minimum: 0.0

selectedatomclass

Set the QSAR visualization option of selected atom class.

Valid values: text strings
Default value: **D**

shininess

Set material property - shininess, for front face.

Valid values: reals
Default value: **80**
Minimum: 0.0
Maximum: 128.0

specular

Set material property - specular, to its red, green, and blue components, for front face.

Valid values: reals

	Default value:	0.1
	Minimum:	0.0
	Maximum:	1.0
<i>step</i>	The step domain tolerance of cells.	
	Valid values:	integers
	Default value:	3
	Minimum:	1
<i>transparency</i>	The transparency of QSAR markers.	
	Valid values:	reals
	Default value:	50
	Minimum:	0.0
	Maximum:	100.0
<i>volumeoccupied</i>	Set the QSAR visualization option of viewing volume occupied.	
	Valid values:	workspacelig qsarmodel
	Default value:	workspacelig

qsitebasis

Specifies an atom for QM basis atoms

Syntax:

```
qsitebasis basis=<text> diaplayname=<text> select=yes | no
          <atom>
```

Options:

<i>basis</i>	The basis set for the atom to be included in the QM basis of an Impact QSite simulation.	
	Valid values:	text strings
	Default value:	lacvp*
<i>diaplayname</i>	The display name for the atom to be included in the QM basis of an Impact QSite simulation.	
	Valid values:	text strings
	Default value:	
<i>select</i>	Whether or not the atom is selected.	

Chapter 5: Commands

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

Operands:

$\langle \text{atom} \rangle$

The atom that defines a basis set

qsitehcap

Specifies a pair of atoms to define a hydrogen cap

Syntax:

qsitehcap *basis*= $\langle \text{text} \rangle$ $\langle \text{atom1} \rangle$ $\langle \text{atom2} \rangle$

Options:

basis The basis set for the hydrogen cap to be included in the QM region of an Impact QSite simulation.

Valid values: text strings
Default value: **lacvp***

Operands:

$\langle \text{atom1} \rangle$ $\langle \text{atom2} \rangle$

The two atoms that define a hydrogen cap

qsitehcapmarkersettings

Set graphical data of QM hydrogen cap markers.

Syntax:


```

qsitehcapmarkersettings ambient=<x> blue=<x>
    cornradius=<x> cylinderheight=<x> cylinderradius=<x>
    diffuse=<x> drawstyle=solid | line emission=<x> green=<x>
    linewidth=<x> red=<x> selectblue=<x> selectgreen=<x>
    selectred=<x> shininess=<x> sliceline=<n> slicesolid=<n>
    specular=<x> stackline=<n> stacksolid=<n>
    transparency=<x>

```

Options:

ambient Set material property - ambient, to its red, green, and blue components, for front face.

Valid values: reals
 Default value: **0.5**
 Minimum: 0.0
 Maximum: 1.0

blue The blue component of attachment markers.

Valid values: reals
 Default value: **0.16**
 Minimum: 0.0
 Maximum: 1.0

cornradius

The radius of corn of attachment markers.

Valid values: reals
 Default value: **0.55**
 Minimum: 0.0

cylinderheight

The cylinder height ratio of attachment markers.

Valid values: reals
 Default value: **0.6**
 Minimum: 0.0
 Maximum: 1.0

cylinderradius

The radius of cylinder of attachment markers.

Valid values: reals
 Default value: **0.26**
 Minimum: 0.0

diffuse

Set material property - diffuse, to its red, green, and blue components, for front face.

Valid values: reals
 Default value: **0.4**
 Minimum: 0.0

Chapter 5: Commands

	Maximum: 1.0
<i>drawstyle</i>	The styles of rendering attachment markers, they are: 1 - solid, and 2 - lines. Default is solid. Valid values: solid line Default value: solid
<i>emission</i>	Set material property - emission, to its red, green, and blue components, for front face. Valid values: reals Default value: 0.05 Minimum: 0.0 Maximum: 1.0
<i>green</i>	The green component of attachment markers. Valid values: reals Default value: 0.16 Minimum: 0.0 Maximum: 1.0
<i>linewidth</i>	Set the width of lines in drawing attachment. Valid values: reals Default value: 1.5 Minimum: 0.0001
<i>red</i>	The red component of attachment markers. Valid values: reals Default value: 0.65 Minimum: 0.0 Maximum: 1.0
<i>selectblue</i>	The blue component of selected attachment markers. Valid values: reals Default value: 1 Minimum: 0.0 Maximum: 1.0
<i>selectgreen</i>	The green component of selected attachment markers. Valid values: reals Default value: 0.9 Minimum: 0.0 Maximum: 1.0
<i>selectred</i>	The red component of selected attachment markers. Valid values: reals

	Default value: 0.2
	Minimum: 0.0
	Maximum: 1.0
<i>shininess</i>	Set material property - shininess, for front face.
	Valid values: reals
	Default value: 80
	Minimum: 0.0
	Maximum: 128.0
<i>sliceline</i>	Set the slices of drawing line attachment.
	Valid values: integers
	Default value: 10
	Minimum: 2
<i>slicesolid</i>	Set the slices of drawing solid attachment.
	Valid values: integers
	Default value: 36
	Minimum: 2
<i>specular</i>	Set material property - specular, to its red, green, and blue components, for front face.
	Valid values: reals
	Default value: 0.2
	Minimum: 0.0
	Maximum: 1.0
<i>stackline</i>	Set the stacks of drawing line attachment.
	Valid values: integers
	Default value: 8
	Minimum: 2
<i>stacksolid</i>	Set the stacks of drawing solid attachment.
	Valid values: integers
	Default value: 18
	Minimum: 2
<i>transparency</i>	The transparency of rendering attachment markers.
	Valid values: reals
	Default value: 20
	Minimum: 0.0
	Maximum: 100.0

qsiteion

Defines an ion as part of the QM region for an Impact QSite simulation

Syntax:

```
qsiteion basis=<text> <molecule_num>
```

Options:

basis The basis set for the ligand to be included in the QM region of an Impact QSite simulation.

Valid values: text strings

Default value: **lacvp***

Operands:

<molecule_num>

The molecule number of an ion to be included in the QM region in a Impact QSite simulation.

qsiteresidue

Defines a residue as part of the QM region for an Impact QSite simulation

Syntax:

```
qsiteresidue backboneresidue=yes | no backbonestring=<text>  
             basis=<text> cuttype=<n>  
             <chain>:<molnum>:<resnum>:<inscode>
```

Options:

backboneresidue

The flag indicates if this residue is a lower or upper bound residue of Free ligand/ion residues.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

backbonestring

The string for all residues of a Free ligand/ion section.

Valid values: text strings

Default value:

basis

The basis set for the residue to be included in the QM region of an Impact QSite simulation.

Valid values: text strings
 Default value: **lacvp***

cuttype The cut type for the residue to be included in the QM region of an Impact QSite simulation.

Valid values: integers
 Default value: **1**
 Minimum: 0
 Maximum: 5

Operands:

$\langle \text{chain} \rangle : \langle \text{molnum} \rangle : \langle \text{resnum} \rangle : \langle \text{inscode} \rangle$

The number of a residue to be included in the QM region in a Impact QSite simulation.

qsiteset

Settings associated with QSite simulations in Impact.

Syntax:

Chapter 5: Commands

```
qsiteset alie=yes | no boxadjustment=⟨x⟩ charge=⟨n⟩  
electrondensity=yes | no enableip472=yes | no endalpha=⟨n⟩  
endalphabase=homominus | lumoplus endbeta=⟨n⟩  
endbetabase=homominus | lumoplus esp=yes | no  
espunits=kcalmolelectron | ktelectron | hartrees | kcalmol | kt  
| ev griddensity=⟨x⟩ hcap=none | point | gaussian |  
gaussgrid maxiter=⟨n⟩ method=dft | hf | lmp2 | pwb6k |  
m06 | m062x | m06l | m06hf | m05 | m052x | user | rm1 |  
am1 | pm3 | mndo | mndod multiplicity=⟨n⟩  
nddo_options=⟨text⟩ noncovgriddensity=⟨x⟩  
noncovinteraction=yes | no numberprocessors=⟨n⟩  
optimization=singlepoint | minimization | transitionstate  
optimize=yes | no options=⟨text⟩ orbitals=yes | no  
pathfraction=⟨x⟩ product_entry=⟨text⟩ reactant_entry=⟨text⟩  
spin_unrestricted=yes | no spindensity=yes | no  
startalpha=⟨n⟩ startalphabase=homominus | lumoplus  
startbeta=⟨n⟩ startbetabase=homominus | lumoplus  
tsguess_entry=⟨text⟩ tsmethod=standard | lst | qst  
use_guess=yes | no use_hessian=yes | no
```

Options:

<i>alie</i>	Whether to calculate Average local ionization energy Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>boxadjustment</i>	The box size adjustment, in angstroms, per side. Valid values: reals Default value: 0
<i>charge</i>	The total charge on the QM part in an Impact QSite simulation. Valid values: integers Default value: 0
<i>electrondensity</i>	Whether to calculate electron density surface. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>enableip472</i>	Whether enable ip472 keyword or not. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>endalpha</i>	The ending alpha orbital number or offset. Valid values: integers

	Default value: 0
	Minimum: 0
<i>endalphabase</i>	The base for the ending alpha orbital number.
	Valid values: homominus lumoplus
	Default value: lumoplus
<i>endbeta</i>	The ending beta orbital number or offset.
	Valid values: integers
	Default value: 0
<i>endbetabase</i>	The base for the ending beta orbital number.
	Valid values: homominus lumoplus
	Default value: lumoplus
<i>esp</i>	Whether to calculate electrostatic potential surface.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>espunits</i>	Which esp units to use for surfaces (kcal/mol/electron, kT/electron at 298.15K, hartrees, kcal/mol, kT at 298.15K, eV).
	Valid values: kcalmolelectron ktelectron hartrees kcalmol kt ev
	Default value: kcalmol
<i>griddensity</i>	The number of grid points per angstrom.
	Valid values: reals
	Default value: 5
<i>hcap</i>	The QM hydrogen cap electrostatics type in a QSite job in Maestro.
	Valid values: none point gaussian gaussgrid
	Default value: gaussgrid

Chapter 5: Commands

<i>maxiter</i>	The maximum number of iterations for the QM optimization. Valid values: integers Default value: 100
<i>method</i>	The QM Method used in a QSite job in Maestro Valid values: dft hf lmp2 pwb6k m06 m062x m06l m06hf m05 m052x user rm1 am1 pm3 mndo mndod Default value: dft
<i>multiplicity</i>	The multiplicity of the QM part in an Impact QSite simulation. Valid values: integers Default value: 1
<i>nddo_options</i>	Any additional Jaguar options which can be used during the Impact QSite job with NDDO method. Valid values: text strings Default value: mmok=1 noxyz=1 nointer=1 trust=0.5
<i>noncovgriddensity</i>	Noncovalent grid density for QSite. Valid values: reals Default value: 20
<i>noncovinteraction</i>	Whether to enable noncovalent interactions for QSite. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>numberprocessors</i>	The number of processors to be used for the Jaguar Part of the calculation.

	Valid values: integers Default value: 1
<i>optimization</i>	What type of calculation is to be performed in the QM part. Valid values: singlepoint minimization transitionstate Default value: minimization
<i>optimize</i>	[NOTE: This option is no longer used.] An option which determines if the QM part will be geometry optimized during a QSite job. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>options</i>	Any additional Jaguar options which can be used during the Impact QSite job. Valid values: text strings Default value: iacc=1 vshift=1.0 maxit=100 trust=0.5
<i>orbitals</i>	Whether to calculate molecular orbital surfaces. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>pathfraction</i>	The fraction of the path between the reactant and product which the TS is along. Valid values: reals Default value: 0.5 Minimum: -0.000000 Maximum: 1.000001
<i>product_entry</i>	The name of the entry which represents the product in a transition state calculation. Valid values: text strings Default value:
<i>reactant_entry</i>	The name of the entry which represents the entry in a transition state calculation. Valid values: text strings Default value:
<i>spin_unrestricted</i>	Whether spin is restricted or not.

Chapter 5: Commands

	Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>spindensity</i>	Whether to calculate spin density surface. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>startalpha</i>	The starting alpha orbital number or offset. Valid values: integers Default value: 0 Minimum: 0
<i>startalphabase</i>	The base for the starting alpha orbital number. Valid values: homominus lumoplus Default value: homominus
<i>startbeta</i>	The starting beta orbital number or offset. Valid values: integers Default value: 0
<i>startbetabase</i>	The base for the starting beta orbital number. Valid values: homominus lumoplus Default value: homominus
<i>tsguess_entry</i>	The name of the entry which represents the transtion state guess in a transition state calculation. Valid values: text strings Default value:
<i>tsmethod</i>	For a transition state calculation, how that is to be performed. Valid values: standard lst qst Default value: standard
<i>use_guess</i>	Whether use wave function in input file or not. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>use_hessian</i>	Whether use hessian in input file or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

quicktorsion

Adjust the dihedral specified by 4 atoms to the given value. Adjust the marker arrow direction via reverse operand.

Syntax:

```
quicktorsion contactdisplay=yes | no dihedral=⟨x⟩
             hbonddisplay=yes | no reverse=yes | no movelarger|⟨atom1⟩
             ⟨atom2⟩ ⟨atom3⟩ ⟨atom4⟩
```

Options:

contactdisplay

This option determines whether Contact markers will be displayed for quick torsion.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

dihedral

Value to set the torsion

Valid values: reals
 Default value: **0**
 Minimum: -180.1
 Maximum: 180.1

hbonddisplay

This option determines whether H-bond markers will be displayed for quick torsion.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

reverse

To set the torsion marker arrow direction

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

Operands:

```
movelarger|⟨atom1⟩ ⟨atom2⟩ ⟨atom3⟩ ⟨atom4⟩
```

Four atoms which are used to adjust a torsion. The *movelarger* operand will flip the direction of already defined dihedral marker arrow toward the larger atom structure.

quit

Quit the program. To quit issue just the quit command without any options.

Syntax:

`quit confirm=yes | no`

Options:

confirm If this option has been set to “false” then the program will exit without prompting the user in any way. So, for example, the user is not prompted to save any changed macros nor will the Quit panel be displayed. Note that invoking quit with any options only sets the option it does not also try to quit Maestro.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

raisepanels

Raise all visible panels

Syntax:

`raisepanels`

read

This is a standard alias for **fileread** (see [\[fileread\]](#), page 195).

readpotential

Read potential settings from a command file.

Syntax:

`readpotential <file_name>`

Operands:

`<file_name>`

The name of the file from which the potential settings will be read. If no name is specified, the default settings will be used.

reagentprep

Options for Reagent Preparation jobs.

Syntax:

```
reagentprep gen_conform=<n> gen_ionization=yes | no
            gen_stereo=<n> gen_tautomers=yes | no
            group_name_long=<text> group_name_short=<text>
            input_file=<text> ph=<x> ph_tolerance=<x>
            sd_title_property=<text> sd_title_source=molecule_name |
            property structure_source=selected_entries | workspace | file
            use_epik=yes | no
```

Options:

gen_conform

The percentage of generating low energy ring conformations.

Valid values: integers

Default value: **1**

Minimum: 1

gen_ionization

An option which allows generating ionization.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

gen_stereo The percentage of generating stereoisomers.

Valid values: integers

Default value: **10**

Minimum: 1

gen_tautomers

An option which allows generating tautomers.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

group_name_long

The long name of selected functional group for reagent preparation job.

Valid values: text strings

Default value:

Chapter 5: Commands

group_name_short

The short name of selected functional group for reagent preparation job.

Valid values: text strings

Default value:

input_file

The name of the structure input file.

Valid values: text strings

Default value:

ph

The reagent ionization pH value.

Valid values: reals

Default value: **7**

Minimum: 0.0

Maximum: 14.0

ph_tolerance

The reagent ionization pH tolerance.

Valid values: reals

Default value: **2**

Minimum: 0.0

Maximum: 7.0

sd_title_property

The property to be used to construct titles for reagents from a SD format file, if *sd_title_source* is 2 - property.

Valid values: text strings

Default value:

sd_title_source

The source of titles for reagents from a SD format file (1 - molecule_name or 2 - property).

Valid values: molecule_name

property

Default value: **molecule_name**

structure_source

Whether to use the selected entries in the current project, or what is in the workspace, or a specified file with multiple structures as structure input for the job.

Valid values: selected_entries

workspace

file

Default value: **file**

use_epik

An option which specifies whether or not Epik should be used for ionization and tautomerization.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

refinestart

Start a Refine input file with the current settings.

Syntax:

```
refinestart
```

refinewrite

Write a Refine input file with the current settings.

Syntax:

```
refinewrite
```

removeatomsfromtable

Removes the selected atoms from the atom-property table.

Syntax:

```
removeatomsfromtable
```

rename

Rename a named object. The object type is the same as the command which is used to create that type of object.

For example to rename a set named “alpha” use: `rename set alpha beta` .

Syntax:

rename <object_type> <current_object_name> <new_object_name>

Operands:

<object_type> <current_object_name> <new_object_name>

The first operand is the name of the existing object.

renameproperty

This is a standard alias for **propertyrename** (see [\[propertyrename\]](#), page 606).

repall

Set global representation properties

Syntax:

repall *ballhresolution*=<n> *balllresolution*=<n> *ballsize*=<x>
bondoutline=yes | no *border*=auto | on | true | yes | off |
false | no *borderscale*=<x> *bstyle*=split | blend
cpkhresolution=<n> *cpklresolution*=<n> *cpksize*=<x>
depthcutoff=<n> *depthfactor*=<x>
drawspheresaspoints=yes | no *enhanceddepth*=yes | no
maxwirewidth=<x> *minwirewidth*=<x> *resolution*=high | low
rstyle=multiple | thick *scalewirewidth*=yes | no
showaltpositions=yes | no *simplifymoving*=yes | no
smooth=yes | no *smooth_sgl*=yes | no *stickradius*=<x>
stickradiusmulti=<x> *thintuberadius*=<x>
thintuberadiusmulti=<x> *tubehresolution*=<n>
tubelresolution=<n> *tuberadius*=<x> *tuberadiusmulti*=<x>
tubestickborder=auto | on | off *usethicklines*=yes | no
wirebondsmooth=auto | on | off *wirethickness*=<n>
wirewidthpoint=<x>

Options:

ballhresolution

Set Ball high resolution

Valid values: integers

Default value: **16**

	Minimum: 4
	Maximum: 53
<i>ballresolution</i>	Set Ball resolution
	Valid values: integers
	Default value: 4
	Minimum: 4
	Maximum: 53
<i>ballsize</i>	Percentage to draw balls out at
	Valid values: reals
	Default value: 16
	Minimum: 4.0
	Maximum: 200.0
<i>bondoutline</i>	Enhance appearance of bond overlaps by drawing in the background color a thin border around each bond
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>border</i>	Set drawing of bond order to be on, off, or automatically determined by viewing scale.
	Valid values: auto
	on
	true
	yes
	off
	false
	no
	Default value: auto
<i>borderscale</i>	The drawing scale at which bond orders will be displayed
	Valid values: reals
	Default value: 30
<i>bstyle</i>	Set the bond style
	Valid values: split
	blend
	Default value: split
<i>cpkhresolution</i>	Set CPK high resolution
	Valid values: integers
	Default value: 20

Chapter 5: Commands

Minimum: 8
Maximum: 57

cpklresolution

Set CPK low resolution

Valid values: integers
Default value: **8**
Minimum: 8
Maximum: 57

cpksize

Percentage to draw CPK spheres out at

Valid values: reals
Default value: **85**
Minimum: 4.0
Maximum: 200.0

depthcutoff

Depth cutoff of ambient occlusion. Used when enhanceddepth is on.

Valid values: integers
Default value: **3**
Minimum: 2
Maximum: 10

depthfactor

Depth factor of ambient occlusion. Used when enhanceddepth is on.

Valid values: reals
Default value: **0.7**
Minimum: 0.0
Maximum: 1.0

drawspheresaspoints

Draw spheres using point shaders (true) or using triangle shaders (false)

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

enhanceddepth

Flag for whether or not enhance depth view.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

maxwirewidth

Maximum thickness of scaled wireframe lines

Valid values: reals
Default value: **3**

	Minimum:	2.0
	Maximum:	50.0
<i>minwirewidth</i>	Minimum thickness of scaled wireframe lines	
	Valid values:	reals
	Default value:	1
	Minimum:	0.1
	Maximum:	2.0
<i>resolution</i>	Set the overall resolution for drawing molecules. low or high.	
	Valid values:	high low
	Default value:	high
<i>rstyle</i>	Set the bond render style	
	Valid values:	multiple thick
	Default value:	multiple
<i>scalewirewidth</i>	Change wire width when zooming.	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	true
<i>showaltpositions</i>	Show the alternate occupancy positions	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	true
<i>simplifymoving</i>	Enable/disable use of simplified moving representation	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	false
<i>smooth</i>	Enable/disable line and polygon antialiasing	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	true
<i>smooth_sgl</i>	Enable/disable line and polygon antialiasing for workspace transformation operations, if Maestro was started with -SGL option.	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	true
<i>stickradius</i>	Radius of sticks in Angstroms	

Chapter 5: Commands

Valid values: reals
Default value: **0.12**
Minimum: 0.01
Maximum: 1.0

stickradiusmulti

Radius of multiple bonds sticks in Angstroms

Valid values: reals
Default value: **0.08**
Minimum: 0.01
Maximum: 0.5

thintuberadius

Radius of thin tubes in Angstroms

Valid values: reals
Default value: **0.08**
Minimum: 0.01
Maximum: 1.0

thintuberadiusmulti

Radius of multiple bonds thin tubes in Angstroms

Valid values: reals
Default value: **0.05**
Minimum: 0.01
Maximum: 0.5

tubehresolution

Set Tube high resolution

Valid values: integers
Default value: **16**
Minimum: 8
Maximum: 57

tubelresolution

Set Tube resolution

Valid values: integers
Default value: **8**
Minimum: 8
Maximum: 57

tuberadius

Radius of tubes in Angstroms

Valid values: reals
Default value: **0.16**
Minimum: 0.01
Maximum: 1.0

tuberadiusmulti

Radius of multiple bonds tubes in Angstroms

Valid values: reals
 Default value: **0.08**
 Minimum: 0.01
 Maximum: 0.5

tubestickborder

Set drawing of tube/stick bond order to be on, off, or automatically determined by viewing scale.

Valid values: auto
 on
 off
 Default value: **auto**

usethicklines

Enables the use of thick lines to replace sticks when rotating.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

wirebondsmooth

Set wire bond smoothing as automatically, on, or off.

Valid values: auto
 on
 off
 Default value: **on**

wirethickness

Thickness of wireframe lines

Valid values: integers
 Default value: **2**
 Minimum: 1
 Maximum: 20

wirewidthpoint

The linewidth stage. Points will be displayed at atoms places if current linewidth greater or equal to the stage.

Valid values: reals
 Default value: **3.5**
 Minimum: 2.0

repatom

Change the representation used to display a group of atoms in the main structure window.

Syntax:

repatom *rep*=none | circle | cpk | ballnstick <ASL>

Options:

rep Type of atom representation

Valid values: none
 circle
 cpk
 ballnstick

Default value: **none**

Operands:

<ASL>

A string in the atom specification language. All atoms which match this specification will have their representation changed.

repatombonds

Set representation of all of atoms' bonds. Uses bond rep specified in the repbond command.

Syntax:

repatombonds <ASL>

Operands:

<ASL>

A string in the atom specification language. All atoms which match this string will have the representation of all bonds to them changed.

repbond

Change the representation used to draw an on-screen bond.

Syntax:

repbond *rep*=none | wire | tube | thin_tube <atom1> <atom2>

Options:

<i>rep</i>	Type of bond representation
Valid values:	none wire tube thin_tube
Default value:	wire

Operands:

<atom1> <atom2>

The operands represent the numbers of the two atoms which define the bond to have its representation changed.

repdefault

Set representation of all atoms and bonds

Syntax:

repdefault *style*=wire | cpk | ballnstick | tube | thin_tube

Options:

<i>style</i>	Default representation
Valid values:	wire cpk ballnstick tube thin_tube
Default value:	wire

replacefromhold

This command replaces the on-screen structure with a copy of the the structure from the specified hold set.

Syntax:

replacefromhold <hold_name>

Operands:

<hold_name>

The name of the hold. This must be the name which was specified when the hold was created using the “hold” command.

repquick

Set representation of all atoms and bonds

Syntax:

repquick *default*=default | wire | cpk | ballnstick | tube |
 thin_tube *style*=default | wire | cpk | ballnstick | tube |
 thin_tube [update]

Options:

<i>default</i>	Default representation
	Valid values: default wire cpk ballnstick tube thin_tube
	Default value: wire
<i>style</i>	Default representation
	Valid values: default wire cpk ballnstick tube thin_tube
	Default value: default

Operands:

[update]

When update is present the currently set representation will be applied to the on-screen structure.

resetcsearch

Deletes all the variables used in a conformational search.

Syntax:

```
resetcsearch
```

residue name

Set the residue name for all atoms which match the ASL specification.

Syntax:

```
residue name <PDBNAME> <ASL>
```

Operands:

```
<PDBNAME> <ASL>
```

The first operand is the PDB residue name which will be used for all atoms which match the specification. The second operand is a valid ASL string which defines the set of atoms which are to have their residue names changed.

residue number

Set the residue number for all atoms which match the ASL specification.

Syntax:

```
residue number <res_num> <ASL>
```

Operands:

```
<res_num> <ASL>
```

The first operand is an integer, optionally followed by a single-character insertion code, that specifies the residue number for all the atoms which match the ASL specification. If the trailing alphabetic character is omitted, the residue insertion code will be set blank. The second operand must be a valid ASL string which specifies all the atoms to have their residue number changed.

residuerenumber

Renumber the residues, starting with the starting number, for all residues which match the ASL specification.

Syntax:

```
residuerenumber  $\langle$  starting_res_num  $\rangle$   $\langle$  ASL  $\rangle$ 
```

Operands:

```
 $\langle$  starting_res_num  $\rangle$   $\langle$  ASL  $\rangle$ 
```

The first operand is an integer which represents the starting residue number. The second operand must be a valid ASL string which specifies a set of residues to renumber.

restorepanels

Restores panel locations

Syntax:

```
restorepanels
```

retype

Change the atom type of the atom number specified by the operand to whatever type or element has previously been made current with the atom command.

Syntax:

```
retype  $\langle$  atom_num  $\rangle$ 
```

Operands:

```
 $\langle$  atom_num  $\rangle$ 
```

An atom number representing the atom which is to have its type changed to the current type.

ribbon

Creates a new Ribbon.

Syntax:

```
ribbon ambient=<x> ambientback=<x> backcoloroffset=<x>
      blend=yes | no calphalinelwidth=<n> calphatubesteps=<n>
      calphatubewidth=<x> color=black | gray | darkblue | blue |
      lightblue | aquamarine | turquoise | springgreen | darkgreen |
      green | limegreen | yellowgreen | yellow | orange | maroon |
      red | pink | plum | purple | bluepurple | white
      curvedlinesteps=<n> curvedlinewidth=<n> diffuse=<x>
      diffuseback=<x> display=ribbononly | atomonly | both
      emission=<x> emissionbackblue=<x> emissionbackgreen=<x>
      emissionbackred=<x> flatstrand=yes | no helixcolor=oncolor
      | twocolors hide=none | allatoms | backboneatoms
      includecalpha=yes | no ladderuknot=<x> ladderwidth=<x>
      lowerresidue=<n> minmolsize=<n> perpendicular=yes | no
      resolution=high | low ribbonendweight=<x>
      ribbonhasthick=yes | no ribbonsteps=<n> ribbonthick=<x>
      ribbonweight=<x> ribbonwidth=<x> scheme=constant |
      secondarystructure | chain | calphaatom | residuecharge |
      residueproperty | residuetype | residueposition | entry
      setlowerresidue=yes | no setupperresidue=yes | no
      shininess=<x> shininessback=<x> showladder=yes | no
      shownormals=yes | no simplifymoving=yes | no specular=<x>
      specularback=<x> sphereslices=<n> spherestacks=<n>
      strandarrowweight=<x> strandarrowwidth=<x>
      strandendweight=<x> strandendweight1=<x> strandsteps=<n>
      strandthick=<x> strandweight=<x> strandwidth=<x>
      style=none | cartoon | ribbon | tube | thintube | curvedline |
      calphaline | calphatube thintubesteps=<n>
      thintubeweight=<x> thintubewidth=<x> tubesteps=<n>
      tubeweight=<x> tubewidth=<x> upperresidue=<n>
      useshader=yes | no <ASL-definition>
```

Options:

<i>ambient</i>	Set material property - ambient, to its red, green, and blue components, for front face.
	Valid values: reals
	Default value: 0.5
	Minimum: 0.0
	Maximum: 1.0

Chapter 5: Commands

ambientback

Set material property - ambient, to its red, green, and blue components, for back face.

Valid values: reals
Default value: **0.05**
Minimum: 0.0
Maximum: 1.0

backcoloroffset

The offset for drawing back ribbons with gray color.

Valid values: reals
Default value: **0.5**
Minimum: 0.0
Maximum: 1.0

blend

The flag indicates whether get smooth color changes along the ribbon.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

calphalinelinewidth

Set linewidth for drawing CA Trace ribbons.

Valid values: integers
Default value: **2**
Minimum: 1
Maximum: 40

calphatubesteps

Set the steps of drawing ribbon CA Tube Trace from one node point (mapped from one CA atom) to another. In the both U and V directions. It is used when the resolution option is HIGH.

Valid values: integers
Default value: **3**
Minimum: 1
Maximum: 10

calphatubewidth

Set CA tube width for drawing CA Trace Tube ribbons.

Valid values: reals
Default value: **0.4**
Minimum: 0.001

color

An option which controls the ribbon constant color.

Valid values: black
 gray
 darkblue
 blue
 lightblue
 aquamarine
 turquoise
 springgreen
 darkgreen
 green
 limegreen
 yellowgreen
 yellow
 orange
 maroon
 red
 pink
 plum
 purple
 bluepurple
 white

Default value: **green**

curvedlinesteps

Set the steps of drawing ribbon Curved Line from one node point (mapped from one CA atom) to another. It is used when the resolution option is HIGH.

Valid values: integers
 Default value: **9**
 Minimum: 1
 Maximum: 100

curvedlinewidth

Set curve width for drawing Curved Line ribbons.

Valid values: integers
 Default value: **2**
 Minimum: 1
 Maximum: 40

diffuse

Set material property - diffuse, to its red, green, and blue components, for front face.

Valid values: reals
 Default value: **0.4**
 Minimum: 0.0
 Maximum: 1.0

diffuseback

Set material property - diffuse, to its red, green, and blue components, for back face.

Valid values: reals

Default value: **0.85**

Minimum: 0.0

Maximum: 1.0

display

An option which controls whether the atoms which define the ribbon are to be shown. The three options are: Ribbon Only, Atoms Only, and Both.

Valid values: ribbononly
 atomsonly
 both

Default value: **ribbononly**

emission

Set material property - emission, to its red, green, and blue components, for front face.

Valid values: reals

Default value: **0.05**

Minimum: 0.0

Maximum: 1.0

emissionbackblue

Set material property - emission, to its blue component, for back face. Emission R, G, B values are used to control the color of back face. The great the values, the lighter the back face.

Valid values: reals

Default value: **0.1**

Minimum: 0.0

Maximum: 1.0

emissionbackgreen

Set material property - emission, to its green component, for back face. Emission R, G, B values are used to control the color of back face. The great the values, the lighter the back face.

Valid values: reals

Default value: **0.1**

Minimum: 0.0

Maximum: 1.0

emissionbackred

Set material property - emission, to its red component, for back face. Emission R, G, B values are used to control the color of back face. The great the values, the lighter the back face.

Valid values: reals

	Default value: 0.1
	Minimum: 0.0
	Maximum: 1.0
<i>flatstrand</i>	The flag indicates whether use flat strand arrows.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>helixcolor</i>	An option which controls whether ribbons in the helical part are to be colored with a single color or with two colors and the inside of the helix has a contrasting color.
	Valid values: onecolor twocolors
	Default value: twocolors
<i>hide</i>	An option which controls whether the atoms which define the ribbon are to be hide when creating ribbons. The three options are: None, All associated atoms, and associated backbone atoms.
	Valid values: none allatoms backboneatoms
	Default value: none
<i>includecalpha</i>	The flag indicates whether ribbons pass through the alpha carbons.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>ladderuknot</i>	The value of parameter in the U direction that determines the ladder end of the residue along DNA/RNA backbone ribbon: [0.0, 1.0]. The great the value, the small the radius.
	Valid values: reals
	Default value: 0.75
	Minimum: 0.0
	Maximum: 1.0
<i>ladderwidth</i>	Set tube width for drawing ladders.
	Valid values: reals
	Default value: 0.4
	Minimum: 0.05
	Maximum: 5.0
<i>lowerresidue</i>	The lower residue limit for residue position color scheme.

Chapter 5: Commands

Valid values: integers
Default value: **1**
Minimum: 1

minmolsize

The minimum size of molecule (number of molecule atoms), which determines if a molecule can be displayed as ribbons.

Valid values: integers
Default value: **20**
Minimum: 1

perpendicular

Set DNA ribbon as perpendicular to base.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

resolution

Set the resolution for drawing ribbons.

Valid values: high
low
Default value: **high**

ribbonendweight

This weight value is used to control the radius of the Ribbon ends. The greater the value, the smaller the radius.

Valid values: reals
Default value: **0.6**
Minimum: 0.001

ribbonhasthick

This bool value is used to determine which kind of ribbon will be drawn, a flat sheet or a solid ribbon has thick.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

ribbonsteps

Set the steps of drawing Ribbon ribbon from one node point (mapped from one CA atom) to another. In the both U and V directions. It is used when the resolution option is HIGH.

Valid values: integers
Default value: **2**
Minimum: 1
Maximum: 10

ribbonthick

This value is used to control the thickness of ribbon, if the ribbon is displayed in solid ribbon.

Valid values: reals

Default value: **0.1**
 Minimum: 0.1
 Maximum: 2.0

ribbonweight

This weight value is used to control the radius of the rectangle cross-section of a strand ribbon. The great the value, the small the radius.

Valid values: reals
 Default value: **3**
 Minimum: 0.001

ribbonwidth

Set ribbon width for drawing Ribbon ribbons.

Valid values: reals
 Default value: **2**
 Minimum: 0.05
 Maximum: 5.0

scheme

An option which controls the ribbon coloring scheme.

Valid values: constant
 secondarystructure
 chain
 calphaatom
 residuecharge
 residueproperty
 residuetype
 residueposition
 entry
 Default value: **secondarystructure**

setlowerresidue

The flag indicates whether use lower residue limit for residue position color scheme.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

setupperresidue

The flag indicates whether use upper residue limit for residue position color scheme.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

shininess

Set material property - shininess, for front face.

Valid values: reals
 Default value: **80**
 Minimum: 0.0

Chapter 5: Commands

Maximum: 128.0

shininessback

Set material property - shininess, for back face.

Valid values: reals

Default value: **80**

Minimum: 0.0

Maximum: 128.0

showladder

Set the flag for showing DNA ladders.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

shownormals

Set the flag for showing normals.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

simplifymoving

Enable/disable use of simplified moving representation.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

specular

Set material property - specular, to its red, green, and blue components, for front face.

Valid values: reals

Default value: **0.7**

Minimum: 0.0

Maximum: 1.0

specularback

Set material property - specular, to its red, green, and blue components, for back face.

Valid values: reals

Default value: **0.5**

Minimum: 0.0

Maximum: 1.0

sphereslices

Set the slices of drawing sphere in the tube style. We use `gluSphere()` to draw a sphere. The higher the slices, the finer the sphere.

Valid values: integers

Default value: **18**

Minimum: 2

spherestacks

Set the stacks of drawing sphere in the tube style. We use `gluSphere()` to draw a sphere. The higher the stacks, the fine the sphere.

Valid values: integers
 Default value: **18**
 Minimum: 2

strandarrowweight

This weight value is used to control the radius of the strand arrow. The great the value, the small the radius.

Valid values: reals
 Default value: **10**
 Minimum: 0.001

strandarrowwidth

This ratio value is used to control the width of the strand arrow.

Valid values: reals
 Default value: **2**
 Minimum: 1.0

strandendweight

This weight value is used to control the radius of the strand ends. The great the value, the small the radius.

Valid values: reals
 Default value: **0.6**
 Minimum: 0.001

strandendweight1

This weight value is used to control the radius of the strand ends in Cartoon. The great the value, the small the radius.

Valid values: reals
 Default value: **50**
 Minimum: 0.001

strandsteps

Set the steps of drawing ribbon Strand from one node point (mapped from one CA atom) to another. In the both U and V directions. It is used when the resolution option is HIGH.

Valid values: integers
 Default value: **3**
 Minimum: 1
 Maximum: 10

strandthick

Set strand thick for drawing Strand ribbons.

Valid values: reals

Chapter 5: Commands

Default value: **0.2**
Minimum: 0.1
Maximum: 2.0

strandweight

Ribbons are defined with NURBS curve or surfaces that are defined by 4D homogeneous coordinate (x, y, z, w) control points array or mesh. The weight w can push/pull away/towards the curve/surface part near the control point, by decreasing/increasing weight value. This weight value is used to control the radius of the rectangle cross-section of a strand ribbon. The great the value, the small the radius.

Valid values: reals
Default value: **5**
Minimum: 0.001

strandwidth

Set strandwidth for drawing Strand ribbons.

Valid values: reals
Default value: **2**
Minimum: 0.05
Maximum: 5.0

style

An option which controls the ribbon style representation.

Valid values: none
cartoon
ribbon
tube
thintube
curvedline
calphaline
calphatube
Default value: **cartoon**

thintubesteps

Set the steps of drawing ribbon Thin Tube from one node point (mapped from one CA atom) to another. In the both U and V directions. It is used when the resolution option is HIGH.

Valid values: integers
Default value: **3**
Minimum: 1
Maximum: 10

thintubeweight

This weight value is used to control the radius of the Thin Tube ends. The great the value, the small the radius.

Valid values: reals
 Default value: **1**
 Minimum: 0.001

thintubewidth

Set thin tube width for drawing Thin Tube ribbons.

Valid values: reals
 Default value: **0.4**
 Minimum: 0.05
 Maximum: 5.0

tubesteps

Set the steps of drawing ribbon Tube from one node point (mapped from one CA atom) to another. In the both U and V directions. It is used when the resolution option is HIGH.

Valid values: integers
 Default value: **3**
 Minimum: 1
 Maximum: 10

tubeweight

This weight value is used to control the radius of the Tube ends. The great the value, the small the radius.

Valid values: reals
 Default value: **1**
 Minimum: 0.001

tubewidth

Set tube width for drawing Tube ribbons.

Valid values: reals
 Default value: **0.9**
 Minimum: 0.05
 Maximum: 5.0

upperresidue

The upper residue limit for residue position color scheme.

Valid values: integers
 Default value: **1**
 Minimum: 1

usesshader

The flag indicates whether using ribbon specified shaders.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

Operands:

⟨ ASL-definition ⟩

The operand must be a valid string in the atom specification language. It will define which atoms are to have a ribbon drawn.

ribbondump

Print out the current option values of the ribbon command.

Syntax:

ribbondump

ringclosure

A command which defines a ring closure to be used during a conformational search.

Syntax:

ringclosure *maximum*=⟨x⟩ *minimum*=⟨x⟩ ⟨atom1⟩ ⟨atom2⟩
 ⟨atom3⟩ ⟨atom4⟩

Options:

maximum The maximum distance between the ends of the ring which will be accepted as a candidate for “closure”.

Valid values: reals
Default value: **2.5**
Minimum: 0.0

minimum The minimum distance between the ends of the ring which will be accepted as a candidate for “closure”.

Valid values: reals
Default value: **0.5**
Minimum: 0.0

Operands:

⟨atom1⟩ ⟨atom2⟩ ⟨atom3⟩ ⟨atom4⟩

The four atom numbers which define a point in a ring which is to be opened while new structures are generated in a conformational search. The actually opening takes place between the second and third atoms specified. Note that specifying a-b-c-d is the same as specifying d-c-b-a.

ringflip

Conversion of aliphatic ring substituents between axial and equatorial

Syntax:

ringflip $\langle \text{atom_num} \rangle$

Operands:

$\langle \text{atom_num} \rangle$

An atom number representing the atom which is a ring atom to have the substituents flipped.

rotate

Rotate in degrees whatever is specified in the transform set. This is either global (all atoms) or a local grouping defined via an ASL in the transform command.

Syntax:

rotate $x=\langle x \rangle$ $y=\langle x \rangle$ $z=\langle x \rangle$ [reset]

Options:

x	Amount in degrees to rotate in X
	Valid values: reals
	Default value: 0
y	Amount in degrees to rotate in Y
	Valid values: reals
	Default value: 0
z	Amount in degrees to rotate in Z
	Valid values: reals
	Default value: 0

Operands:

[reset]

If reset is present the global rotation matrix will be reset.

run

This is a standard alias for **scriptrun** (see [\[scriptrun\]](#), page 698).

saveimage

Capture the current main structure window and save to an image file.

Syntax:

```
saveimage format=tiff | jpeg | png jpeg-quality=<n>  
          png-compression=<n> png-gamma=<x> showoptions=yes | no  
          smooth=yes | no transparentbg=yes | no <file_name>
```

Options:

format Specifies the format of the saved image.

Valid values: tiff
 jpeg
 png

Default value: **png**

jpeg-quality

Quality of the JPEG file (1-100)

Valid values: integers

Default value: **95**

png-compression

Compression ratio of png format

Valid values: integers

Default value: **6**

png-gamma

Gamma of the image in png format

Valid values: reals

Default value: **0**

Minimum: 0.0

Maximum: 1.0

showoptions

If enabled, show expanded options when dialog is posted.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

smooth

Image saved will have smooth curves in it if enabled.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

transparentbg

Image saved will have transparent background color if enabled.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

⟨file_name⟩

The file where the image will be saved.

saveimageheight

Sets the height of the image saved due to saveimage command

Syntax:

saveimageheight ⟨height⟩

Operands:

⟨height⟩

Height of the image to be saved.

saveimagewidth

Sets the width of the image saved due to saveimage command

Syntax:

saveimagewidth ⟨width⟩

Operands:

⟨width⟩

Width of the image to be saved.

savelayout

Save the size and position of all currently visible panels.

Syntax:

savelayout

savemovie

Syntax:

savemovie <filename>

Operands:

<filename>

The file where the movie will be saved.

savemovieoptions

Sets the options for the save movie while eplaying entries

Syntax:

savemovieoptions *addframes*=<n> *duration*=*frameduration* |
movieduration *frameduration*=<x> *movieduration*=<n>
quality=low | medium | high *resolution*=low | medium | high
smooth=yes | no *speed*=<n>

Options:

addframes This option controls the number of frames to add when 'smooth' is enabled.

Valid values: integers
Default value: **20**

duration This option sets the duration option for the movie

Valid values: *frameduration*
movieduration
Default value: **frameduration**

frameduration

This option sets the number of seconds that each frame to be displayed for.

Valid values: reals
Default value: **0.04**
Minimum: 0.04
Maximum: 5.0

movieduration

This option sets the number of seconds that the recorded movie should be played for.

	Valid values: integers
	Default value: 1
	Minimum: 1
	Maximum: 300
<i>quality</i>	This option sets the quality of movie file.
	Valid values: low medium high
	Default value: medium
<i>resolution</i>	This option sets the resolution of movie file.
	Valid values: low medium high
	Default value: medium
<i>smooth</i>	This option controls whether additional interpolated frames are added to the movie or not. Entries should be conformers. This is useful for creating a smoother animation in the movie for things like molecular dynamics (MD) results. Keep in mind this is strictly a linear interpolation of coordinates.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>speed</i>	This option sets the number of structures to be displayed per second.
	Valid values: integers
	Default value: 1

scanmode

Set the Coordinate Scan mode.

Syntax:

scanmode *modetype*=distance | angle | dihedral

Options:

<i>modetype</i>	Determines the mode for coordiante scan softlimit=<n> <input_file_name>
	Valid values: distance angle dihedral

Default value: **distance**

scriptlogfile

Commands are by default logged to a temporary file which is deleted when the program terminates. If the “logfile” command is used to name the logfile then the commands will be logged to that file and that file will not be deleted when the program ends.

Syntax:

scriptlogfile *logging*=yes | no <logfile_name>

Options:

logging If this option is set to “false” then no logging will be done to a file or in memory.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

Operands:

<logfile_name>

The name of the file to which the commands will be logged. The full name of the file, including any suffix must be included.

Aliases:

logfile (see [\[logfile\]](#), page 338)

scriptrun

Run the script file whose name is given as the operands.

Syntax:

scriptrun <script_name>

Operands:

<script_name>

The name of the script which is to be executed. The full name, including any suffix must be given. If the file is not in the local directory then a full pathname must be given.

Aliases:

run (see [\[run\]](#), page 693)

searchdbconfgen

Defines settings for Find Matches to Hypothesis Generate Conformers job.

Syntax:

```
searchdbconfgen amidebonds=vary | retain | trans
                  eliminate=atom_deviation | rmsd field=mmffs | opls2005
                  incorporate=append | replace | ignore | appendungrouped |
workspace max_rmsd=⟨x⟩ maxdist=⟨x⟩ method=default |
mixed minimizationsteps=⟨n⟩ numrotatablesteps=⟨n⟩
numsteps=⟨n⟩ postmaxiter=⟨n⟩ postprocessing=yes | no
postprocessingmethod=mini | filter | rce premaxiter=⟨n⟩
preprocessing=yes | no sampling=standard | rapid | complete
| thorough skipconformergeneration=⟨n⟩ solvation=gbsa |
distance_dependent window=⟨x⟩
```

Options:

amidebonds

This determines whether to vary amide bond conformation, retain original amide bond conformation, or set amide bond conformation to trans.

Valid values: vary
 retain
 trans

Default value: **vary**

eliminate

The method to use for eliminating redundant conformers: maximum atom deviation or RMSD.

Valid values: atom_deviation
 rmsd

Default value: **rmsd**

field

This determines which force field mmffs|mmff|opls2001 is used. Currently we always use mmffs, so it will have only one option value.

Valid values: mmffs
 opls2005

Default value: **opls2005**

Chapter 5: Commands

incorporate

This option controls the incorporation of the results (replace or append).

Valid values: append
 replace
 ignore
 appendungrouped
 workspace

Default value: **replace**

max_rmsd Maximum RMSD for considering two structures equal.

Valid values: reals

Default value: **1**

Minimum: 0.0

maxdist Maximum distance between atoms in equal structures.

Valid values: reals

Default value: **2**

Minimum: 0.0

method This determines whether MacroModel uses the ligand torsion search method (default) or the mixed MCMM/LMOD search method (mixed) to generate conformers. Currently Find Matches always uses the default method, so it will have only one option value.

Valid values: default
 mixed

Default value: **default**

minimizationsteps

This option determines the maximum number of minimization steps for Mixed MCMM/LMOD conformer generation.

Valid values: integers

Default value: **100**

Minimum: 1

numrotatablesteps

An option which sets the number of steps which will be performed during the ConfGen conformational search.

Valid values: integers

Default value: **10**

Minimum: 1

numsteps An option which sets the number of steps which will be performed during the conformational search. This also limits number of conformations generated.

Valid values: integers
 Default value: **100**
 Minimum: 0

postmaxiter

This option determines the maximum number of iterations for post-minimization of generated structures.

Valid values: integers
 Default value: **50**
 Minimum: 0
 Maximum: 9999999

postprocessing

Indicates whether or not to perform MacroModel postprocessing.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

postprocessingmethod

This determines which type of postprocessing method to use (minimization, filtering and redundant conformer elimination, or redundant conformer elimination only).

Valid values: mini
 filter
 rce
 Default value: **rce**

premaxiter

This option determines the maximum number of iterations for pre-minimization of input structures.

Valid values: integers
 Default value: **100**
 Minimum: 0
 Maximum: 9999999

preprocessing

Indicates whether or not to perform MacroModel preprocessing.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

sampling

This determines whether rapid (standard) or thorough (complete) sampling will be used.

Valid values: standard
 rapid
 complete
 thorough
 Default value: **standard**

skipconformergeneration

Skip conformer generation for structures with more than the given number of rotatable bonds.

Valid values: integers

Default value: **15**

Minimum: 0

solvation

This determines whether GB/SA Water (gbsa) or Distance Dependent Dielectric (distance_dependent) solvation treatment is used.

Valid values: gbsa
distance_dependent

Default value: **distance_dependent**

window

The energy window (in kcal/mol) within which structures will be saved.

Valid values: reals

Default value: **10**

Minimum: 0.0

selectatomsforediting

Selects atoms for editing their properties.

Syntax:

selectatomsforediting <ASL>

Operands:

<ASL>

The operand is a valid ASL string which specifies which atoms are to be selected.

selecteditatom

Selects a row in the atom-level property table.

Syntax:

selecteditatom \langle row number \rangle

Operands:

\langle row number \rangle

The operand is the row number to select.

selecteditproperty

Selects an atom-level property from the edit atom-level property list.

Syntax:

selecteditproperty \langle property \rangle

Operands:

\langle property \rangle

The operand is a property from the edit property list.

sequenceviewersaveimage

Capture the sequence viewer contents and save it to an image file.

Syntax:

sequenceviewersaveimage *format*=tiff | jpeg | png
sequenceviewer= \langle text \rangle \langle file_name \rangle

Options:

format Specifies the format of the saved image.

Valid values: tiff
 jpeg
 png
 Default value: **tiff**

sequenceviewer

Specifies whether to save the contents of the Workspace or Prime sequence viewer.

Valid values: text strings
 Default value: **workspace**

Operands:

⟨file_name⟩

The file where the image will be saved.

set

Creates a new named set. The set name must be a single token (or “quoted” if multiple tokens). A set can be redefined by specifying a new definition.

Syntax:

set ⟨set_name⟩ ⟨ASL-definition⟩

Operands:

⟨set_name⟩ ⟨ASL-definition⟩

The name which will be applied to the set. If the name contains embedded spaces then it must be enclosed in double quotation marks.

setread

Read set definitions from the file whose name is given as the operand. The file name usually has a “.set” suffix.

Syntax:

setread ⟨set_file_name⟩

Operands:

⟨set_file_name⟩

The name of the file from which the set definitions are to be read. The full name of the file (including any .set suffix) must be specified.

setwrite

Write the currently defined sets to the file whose name is given as the operand. The file name usually has a “.set” suffix

Syntax:

setwrite <set_file_name>

Operands:

<set_file_name>

The name of the file to which the current set definitions are to be written .
The full name of the file (including any .set suffix) must be specified.

showdockablepanels

Show visible dockable panels previously hidden with
MM_ID_HIDE_DOCKABLE_PANELS.

Syntax:

showdockablepanels

showfirstselectedentry

Scroll to make first selected entry visible in Project Table.

Syntax:

showfirstselectedentry

showhwstereosetup

Write out the hardware stereo setup.

Syntax:

showhwstereosetup

showmarkers

This command displays the given marker.

Syntax:

showmarkers <marker_name>

Operands:

<marker_name>

The name of the marker.

showpanel

Show the panel whose name is given by the operands.

Syntax:

showpanel <panel_name> [:<tab_name>]

Operands:

<panel_name> [:<tab_name>]

The first operand is the name of the panel which is to be displayed. The name must match to all characters. The names of the panels to be used in the “showpanel” command are displayed in parentheses after each item in the main menu bar. The optional second argument is the name of a tab folder within that panel which is to be made the current tab folder. The name of a tab folder to be used in the second optional operand is displayed in the associated panel. The name must match to all characters, but it is not case sensitive.

showpanels

Make visible panels previously hidden with HIDEPANELS

Syntax:

showpanels

showpropertiesselectedentries

Displays only those properties for which at least one of the selected entries has a value.

Syntax:

showpropertyselectedentries

showproperty

This command shows the given property by creating a subset of all the properties in show state. This function also switches to the property subset view.

Syntax:

showproperty \langle propertyname \rangle

Operands:

\langle propertyname \rangle

The name of the property to show.

showtoolbar

Show the toolbar of given id displayed under given panel.

Syntax:

showtoolbar \langle panel_name:toolbar_id \rangle

Operands:

\langle panel_name:toolbar_id \rangle

The first part is the name of the panel under which toolbar would be displayed. The name must match to all characters and it is a case sensitive. The second part is the name of a toolbar id which needs to be displayed. Id must match to all characters and it is a case sensitive. Any toolbar can not be displayed under any panel. There is a fix set of toolbars which can be displayed under given panel.

sleep

This is a standard alias for **energysleep** (see [\[energysleep\]](#), page 134).

sorteditatomcolumn

Sorts the given column in the atom-property table.

Syntax:

sorteditatomcolumn \langle column \rangle

Operands:

\langle column \rangle

The column to sort.

specifiedname

Set the specified name to that specified for all atoms which match the ASL specification.

Syntax:

specifiedname \langle specified_name \rangle \langle ASL \rangle

Operands:

\langle specified_name \rangle \langle ASL \rangle

The first operand is the specified name which is to be applied to the atom. Only the first 20 characters of the specified name will be used. The second operand is the ASL specification for all the atoms which are to have the specified name applied.

spotcenter

Center the given atom in the Workspace. Make it the center of global rotation.

Syntax:

spotcenter \langle atom_num \rangle

Operands:

\langle atom_num \rangle

The number of the atom which is to be centered in the Workspace and to become the center of rotation.

spotcenterpoint

Center the given point in the Workspace. Make it the center of global rotation.

Syntax:

```
spotcenterpoint <x y z>
```

Operands:

```
<x y z>
```

The point is to be centered in the Workspace and to become the center of rotation.

startdiagnosis

This commands collects diagnostic information such as total system memory, available system memory, running maestro process peak virtual memory usage, current virtual memory usage, peak resident memory usage (only Linux), current resident memory usage (Only Linux), total swap disk space, free swap disk space, page fault etc. It also starts timer to collect user, sys, and real time information when explicitly called to stop it. All information are written into the file residing LINUX : \$HOME/.schrodinger/maestro-diagnosis-[unique id](#).txt Winows : App-Data/Schrodinger//maestro-diagnosis-[unique id](#).txt Unique file name is generated using second, minute, hour, day, month, and year. This file is created for each session.

Syntax:

```
startdiagnosis
```

stop

This is a standard alias for **energystop** (see [\[energystop\]](#), page 135).

stopdiagnosis

This commands collects diagnostic information such as total system memory, available system memory, running maestro process peak virtual mem-

ory usage, current virtual memory usage, peak resident memory usage (only Linux), current resident memory usage (Only Linux), total swap disk space, free swap disk space, page fault etc. It also stops timer which is started by MM_ID_START_DIAGNOSIS command and write user, system, and real time information for the tasks performed between start and stop operation. All information are written into the file residing LINUX : \$HOME/.schrodinger/maestro-diagnosis-*<unique id>*.txt Winows : AppData/Schrodinger//maestro-diagnosis-*<unique id>*.txt Unique file name is generated using second, minute, hour, day, month, and year. This file is created for each session.

Syntax:

stopdiagnosis

strikebuildqsar

Runs a Build QSAR job.

Syntax:

```
strikebuildqsar activity_property=<text> max_pls_factors=<n>
                    method=pls | pca | mlr remove_outliers=yes | no
                    supyintercept=yes | no which_descriptors=all | subset
                    <job name>
```

Options:

activity_property

The name of the activity property.

Valid values: text strings

Default value:

max_pls_factors

How many factors for Partial Least Squares.

Valid values: integers

Default value: **1**

Minimum: 0

method

The regression method to use.

Valid values: pls
pca
mlr

Default value: **pls**

remove_outliers

Set to true to automatically remove outliers.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

supyintercept

Set to true to force the y-intercept in MLR fitting (MLRO & MLRS) to be zero.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

which_descriptors

Which descriptors to use.

Valid values: all
 subset

Default value: **all**

Operands:

⟨job name⟩

The name of the job to run.

strikedeletemodel

Deletes the given model.

Syntax:

strikedeletemodel ⟨model name⟩

Operands:

⟨model name⟩

The name of the model to delete.

strikeexportmodel

Exports a model from Strike.

Syntax:

strikeexportmodel $\langle \text{file} \rangle$

Operands:

$\langle \text{file} \rangle$

The file to export the Strike model to.

strikeextendselectdescriptor

Extends the selected descriptors via the given descriptor.

Syntax:

strikeextendselectdescriptor $\langle \text{M2IO descriptor name} \rangle$

Operands:

$\langle \text{M2IO descriptor name} \rangle$

The M2IO data name of the descriptor to select.

strikeimportmodel

Imports a model to the Strike panels.

Syntax:

strikeimportmodel $\langle \text{file} \rangle$

Operands:

$\langle \text{file} \rangle$

The file to import the Strike model from.

strikeplotmodel

Plots the given model.

Syntax:

strikeplotmodel \langle model name \rangle

Operands:

\langle model name \rangle

The name of the model to plot.

strikepredict

Runs a prediction job.

Syntax:

strikepredict \langle job name \rangle

Operands:

\langle job name \rangle

The name of the job to run.

strikeselectdescriptor

Selects the given descriptor

Syntax:

strikeselectdescriptor \langle M2IO descriptor name \rangle

Operands:

\langle M2IO descriptor name \rangle

The M2IO data name of the descriptor to select.

strikeselectmodel

Selects only the given model.

Syntax:

strikeselectmodel \langle model name \rangle

Operands:

\langle model name \rangle

The name of the model to select.

strikesimilarity

Runs a similarity job.

Syntax:

strikesimilarity *job_type*=atompairs | descriptors \langle job name \rangle

Options:

job_type Whether to calculate atom pair or descriptor similarities.

Valid values: atompairs
 descriptors

Default value: **atompairs**

Operands:

\langle job name \rangle

The name of the job to run.

striketoggleselectdescriptor

Toggles the selection of the given descriptor on or off.

Syntax:

striketoggleselectdescriptor \langle M2IO descriptor name \rangle

Operands:

\langle M2IO descriptor name \rangle

The M2IO data name of the descriptor to select.

structalignatoms

Sets the ASL that the next struct align job will operate on.

Syntax:

```
structalignatoms <ASL>
```

Operands:

<ASL>

The residues to align.

structalignatomsalt

Sets the separate ASL that the next struct align job will operate on. If the separate asl is set, the original ASL will be used as reference residues.

Syntax:

```
structalignatomsalt <ASL>
```

Operands:

<ASL>

The residues to align.

structalignstart

Start a struct align job with the current settings.

Syntax:

```
structalignstart
```

substructure

Specifies a set of atoms to be used as the “substructure” during a substructure minimization or dynamics simulation.

Syntax:

substructure *fillres*=yes | no *radius*=⟨x⟩ ⟨ASL⟩

Options:

fillres A boolean option that determines if the substructure definition will be expended to complete residue boundaries.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

radius The radius from the basic substructure definition within which atoms will be included in the substructure.

Valid values: reals

Default value: **0**

Minimum: 0.0

Operands:

⟨ASL⟩

A string in the atom specification language. All atoms which match this will be treated as part of the substructure in a substructure minimization.

substructurefile

This is a standard alias for **substructurefilewrite** (see [\[substructurefilewrite\]](#), [page 717](#)).

substructurefileread

Will read a .sbc file and replace the current substructure and constrained/fixed atom shells.

Syntax:

substructurefileread ⟨sbc_file_name⟩

Operands:

⟨sbc_file_name⟩

The complete name (including the *.sbc suffix) of the file from which the substructure information will be read.

substructurefilewrite

Will write a .sbc file with the current substructure and constrained/fixed atom shells.

Syntax:

```
substructurefilewrite absolutecoords=yes | no  
                    aslformat=yes | no <sbcs_file_name>
```

Options:

absolutecoords

A boolean option which determines whether the atoms are considered to be at absolute coordinates (fixed at their current workspace locations) or at relative coordinates (at their locations in the input structures).

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

aslformat

A boolean option which determines whether the atoms specification to be written in ASL format. If user specifies 'aslformat=yes' then absolutecoords option will be ignored.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

<sbcs_file_name>

The complete name (including the *.sbc suffix) of the file to which the substructure information will be written.

Aliases:

substructurefile (see [\[substructurefile\]](#), page 716)

substructureshell

Specifies a set of atoms forming a “shell” around the substructure defined by the “subs” command.

Syntax:

```
substructureshell addatoms=⟨text⟩ constant=⟨x⟩
                    fillres=yes | no frozen=yes | no radius=⟨x⟩ ⟨shell_number⟩
                    ⟨ASL-definition⟩
```

Options:

<i>addatoms</i>	A string in the atom specification language. All atoms which match this will be treated as part of the shell. Valid values: text strings Default value:
<i>constant</i>	The harmonic force constant to be applied to the constrained shell. Valid values: reals Default value: 200 Minimum: 0.0
<i>fillres</i>	A boolean option which determines whether the shell is to be made up of complete residues. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>frozen</i>	A boolean option which determines whether the atoms in the shell are to be “frozen” during the substructure calculation. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>radius</i>	The radius of a “shell” of atoms around the “substructure” during a substructure energy procedure. Valid values: reals Default value: 0 Minimum: 0.0

Operands:

⟨shell_number⟩ ⟨ASL-definition⟩

The number of the shell. Shells are usually numbered from 1, but any sequence will work. The shell with the lowest number is defined relative to the atom in the substructure. Subsequent shells are defined relative to next lowest numbered shell.

superimpose

Perform a superposition using previously defined atoms or all corresponding atoms in all onscreen entries if all is used.

Syntax:

superimpose *inplace*=yes | no [all]

Options:

inplace Enable superposition to be performed in place (effectively just calculates the RMS and doesn't actually move any entries.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

Operands:

[all]

If all is present then an attempt will be made to superimpose all onscreen entries. Otherwise only the atom pairs defined by superimposeatom commands will be used to perform a superposition.

superimposeatom

Define an atom pair for which will be superimposed in a subsequent superposition operation.

Syntax:

superimposeatom \langle atom1 \rangle \langle atom2 \rangle

Operands:

\langle atom1 \rangle \langle atom2 \rangle

Two atom numbers which represent an atom pair to be superimposed by a subsequent superimpose command. The two atoms must be from different entries and all superimposeatom commands must specify entries in the same order.

superimposeset

Use an ASL expression to define superposition atom pairs. The ASL set must define exactly the same number of atoms in each on-screen entry.

Syntax:

superimposeset $\langle \text{ASL} \rangle$

Operands:

$\langle \text{ASL} \rangle$

A string in the atom specification language. This set must define exactly the same number of atoms in each on-screen entry and those atoms will become the basis for superposition.

superimposesmarts

Use a SMARTS expression to define superposition atom pairs. The SMARTS expression must define exactly the same number of atoms in each on-screen entry.

Syntax:

superimposesmarts $\langle \text{SMARTS} \rangle$

Operands:

$\langle \text{SMARTS} \rangle$

A SMARTS expression. This must define exactly the same number of atoms in each on-screen entry and those atoms will become the basis for superposition.

surfaceactivegrid

Changes the settings of active grid for the given surface.

Syntax:

surfaceactivegrid *centerx*= $\langle x \rangle$ *centery*= $\langle x \rangle$ *centerz*= $\langle x \rangle$
entry= $\langle \text{text} \rangle$ *mode*=center | size | use | center_point
size= $\langle x \rangle$ *surface*= $\langle \text{text} \rangle$ *use*=yes | no *volume*= $\langle \text{text} \rangle$
 $\langle \text{ASL} \rangle$

Options:

centerx The X coordinate of active grid center for the surface.

Valid values: reals

Default value: **0**

centery The Y coordinate of active grid center for the surface.

	Valid values: reals Default value: 0
<i>centerz</i>	The Z coordinate of active grid center for the surface. Valid values: reals Default value: 0
<i>entry</i>	The entry name of the entry that the surface belongs to. Valid values: text strings Default value:
<i>mode</i>	Sets the mode of setting active grid for the given surface, mode is center, size, or use. Valid values: center size use center_point Default value: center
<i>size</i>	The size of active grid for the surface. Valid values: reals Default value: 10 Minimum: 0.001
<i>surface</i>	The name of the surface. Valid values: text strings Default value:
<i>use</i>	Whether or not use active grid for eletron density map for the surface . Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>volume</i>	The name of the volume that the surface belongs to. Valid values: text strings Default value:

Operands:

⟨ ASL ⟩

The ASL defines the center of active grid.

surfacecontours

Control surface contour settings Edges are contoured when the distance along z to the nearest surface point behind it is at least the distance specified in the contour 'depth' value. The contour color is controlled by the red,

green, and blue option values. The contour line's thickness is controlled via the 'thickness' option. The color rendered for the contour is additionally affected by 'intensity' which affects the final color by blending the contour color and the surface color according to 'intensity'.

Syntax:

```
surfacecontours blue=<x> contours=yes | no depth=<x>
               green=<x> intensity=<x> red=<x> thickness=<x>
```

Options:

<i>blue</i>	Sets blue component of contour color Valid values: reals Default value: 0 Minimum: 0.0 Maximum: 1.0
<i>contours</i>	Display or hide contours on opaque surfaces. For transparent surfaces these are never shown regardless of what value this setting has. Valid values: boolean (true false; yes no; y n; on off) Default value: false
<i>depth</i>	Display contours where the distance to the closest surface point is at least this distance in z. Valid values: reals Default value: 1.2 Minimum: 0.0 Maximum: 5.0
<i>green</i>	Sets green component of contour color Valid values: reals Default value: 0 Minimum: 0.0 Maximum: 1.0
<i>intensity</i>	Intensity of the contour. A value of 1.0 is the most intense showing only the contour color. A value of zero is the least intense showing none of the contour color and only the surface color. Valid values: reals Default value: 0.2 Minimum: 0.0 Maximum: 1.0
<i>red</i>	Sets red component of contour color

Valid values:	reals
Default value:	0
Minimum:	0.0
Maximum:	1.0

thickness Sets thickness of contour

Valid values:	reals
Default value:	0.1
Minimum:	0.0
Maximum:	1.0

surfacedarkencolor

Darkens the color by cavity depth for the given surface.

Syntax:

```
surfacedarkencolor bycavitydepth=yes | no entry=⟨text⟩
                  ⟨surface⟩
```

Options:

bycavitydepth

Darkens the color by surface cavity depth.

Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

entry The entry that the surface belongs to.

Valid values:	text strings
Default value:	

Operands:

⟨surface⟩

The name of the surface to set whether to darken colors by cavity depth when drawing the surface.

surfacedelete

Deletes the given surface.

Chapter 5: Commands

Syntax:

surfacedelete *entry*=⟨text⟩ ⟨surface⟩

Options:

entry The entry that the surface belongs to.
Valid values: text strings
Default value:

Operands:

⟨surface⟩

The name of the surface to delete.

surfacedisplay

Displays or undisplay the given surface

Syntax:

surfacedisplay *display*=yes | no *entry*=⟨text⟩ ⟨surface⟩

Options:

display Sets whether or not to display the surface.
Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

entry The entry that the surface belongs to.
Valid values: text strings
Default value:

Operands:

⟨surface⟩

The name of the surface to display.

surfaceduplicate

Duplicates the given surface.

Syntax:

surfaceduplicate *entry*=⟨text⟩ ⟨surface⟩

Options:

entry The entry that the surface belongs to.
 Valid values: text strings
 Default value:

Operands:

⟨surface⟩

The name of the surface to duplicate.

surfaceextended

Creates a new extended radius surface for the current workspace.

Syntax:

surfaceextended *atom_radius_scaling*=⟨x⟩ *context*=entry |
 molecule | workspace | none | asl *grid_spacing*=⟨x⟩
probe_radius=⟨x⟩ *transparency*=⟨x⟩ name of the surface,
 followed by ASL defining the atoms to be surfaced

Options:

atom_radius_scaling
 This is scaling of VdW radii.
 Valid values: reals
 Default value: **1**
 Minimum: 0.5
 Maximum: 10.0

context This option sets the context for generating a surface.
 Valid values: entry
 molecule
 workspace
 none
 asl
 Default value: **entry**

grid_spacing
 This option sets the grid spacing for the surface.

Valid values:	reals
Default value:	0.8
Minimum:	0.01

probe_radius

This is the probe radius.

Valid values:	reals
Default value:	1.4
Minimum:	1.0

transparency

The transparency option controls how transparent the surface appears. A value of 100 means 100 percent transparent. 0 means completely opaque.

Valid values:	reals
Default value:	0
Minimum:	0.0
Maximum:	100.0

Operands:

name of the surface, followed by ASL defining the atoms to be surfaced
Name of the extended radius surface to be created

surfaceextendedradiuscontext

Defines a set of atoms for which a surface can be clipped against with the `surfaceextendedradius` command.

Syntax:

surfaceextendedradiuscontext <ASL>

Operands:

<ASL>

The ASL expression which defines the atoms the surface will be clipped against.

surfaceextendedradiusset

Defines a set of atoms for which a surface can be created for with the `surfaceextendedradius` command.

Syntax:

surfaceextendedradiusset \langle ASL \rangle

Operands:

\langle ASL \rangle

The ASL expression which defines the atoms the surface will be created for.

surfacemolecular

Creates a new molecular surface for the current workspace.

Syntax:

surfacemolecular *atom_radius_scaling*= \langle x \rangle
boolean_operation=none | union | intersection | difference
context=entry | molecule | workspace | none | asl
edgesmoothing=yes | no *grid_spacing*= \langle x \rangle *probe_radius*= \langle x \rangle
transparency= \langle x \rangle \langle surface name \rangle

Options:

atom_radius_scaling

This is scaling of VdW radii.

Valid values: reals
 Default value: **1**
 Minimum: 0.5
 Maximum: 10.0

boolean_operation

This option sets the context for generating a surface.

Valid values: none
 union
 intersection
 difference
 Default value: **none**

context

This option sets the context for generating a surface.

Valid values: entry
 molecule
 workspace
 none
 asl
 Default value: **none**

Chapter 5: Commands

edgesmoothing

Enables edge smoothing of truncated surfaces.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

grid_spacing

This option sets the grid spacing for the surface.

Valid values: reals

Default value: **0.8**

Minimum: 0.01

probe_radius

This is the probe radius.

Valid values: reals

Default value: **1.4**

Minimum: 1.0

transparency

The transparency option controls how transparent the surface appears. A value of 100 means 100 percent transparent. 0 means completely opaque.

Valid values: reals

Default value: **0**

Minimum: 0.0

Maximum: 100.0

Operands:

⟨surface name⟩

Name of the molecular surface to be created

surfacemolecularcontext

Defines a set of atoms for which a surface will be clipped against with the `surfacemolecular` command.

Syntax:

surfacemolecularcontext ⟨ASL⟩

Operands:

⟨ASL⟩

The ASL expression which defines the atoms the surface will be clipped against.

surfacemolecularset

Defines a set of atoms for which a surface can be created for with the `surfacemolecular` command.

Syntax:

```
surfacemolecularset <ASL>
```

Operands:

<ASL>

The ASL expression which defines the atoms the surface will be created for.

surfacemolecularset2

Defines a set of atoms for surface boolean operation.

Syntax:

```
surfacemolecularset2 <ASL>
```

Operands:

<ASL>

The ASL expression which defines the atoms for surface boolean operation.

surfaceoption

Set surface options.

Syntax:

```

surfaceoption cavitydepth=yes | no cavityiterations=⟨n⟩
                defaultttransparency=⟨n⟩ defaultttransparencyback=⟨n⟩
                defaultttransparencytogether=yes | no extentstep=⟨x⟩
                fastrender=yes | no keepcenter=yes | no
                lowqualitytransparency=rotating | always | never
                maxmeshwidth=⟨x⟩ meshwidth=⟨n⟩ minmeshwidth=⟨x⟩
                pairedcolor1=black | gray | dark_blue | blue | light_blue |
                aquamarine | turquoise | spring_green | dark_green | green |
                lime_green | yellow_green | yellow | orange | maroon | red |
                pink | plum | magenta | blue_violet | white
                pairedcolor2=black | gray | dark_blue | blue | light_blue |
                aquamarine | turquoise | spring_green | dark_green | green |
                lime_green | yellow_green | yellow | orange | maroon | red |
                pink | plum | magenta | blue_violet | white
                randomstipple=interlaced | always | never
                scalemeshwidth=yes | no showtable=yes | no
                singlecolor=black | gray | dark_blue | blue | light_blue |
                aquamarine | turquoise | spring_green | dark_green | green |
                lime_green | yellow_green | yellow | orange | maroon | red |
                pink | plum | magenta | blue_violet | white
                singlelayertransparency=yes | no smoothing=⟨n⟩ style=solid
                | mesh | dots usecheesy=yes | no

```

Options:

cavitydepth

Enables weighting surface colors by cavity depth.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

cavityiterations

Sets number of iterations of cavity depth coloring algorithm.

Valid values: integers

Default value: **50**

Minimum: 5

Maximum: 200

defaultttransparency

Sets default front surface transparency.

Valid values: integers

Default value: **0**

Minimum: 0

Maximum: 100

defaultttransparencyback

Sets default back surface transparency.

	Valid values: integers
	Default value: 0
	Minimum: 0
	Maximum: 100
<i>defaulttransparencytogether</i>	Set the flag of adjusting default front and back surfaces together.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>extentstep</i>	Step size of increasing or decreasing density map extent.
	Valid values: reals
	Default value: 1
	Minimum: 0.1
<i>fastrender</i>	Enables optimized rendering of transparent surfaces.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>keepcenter</i>	Set the flag of keeping electron density map center while translation is performed.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>lowqualitytransparency</i>	Specifies drawing quality for transparent surfaces when rotating.
	Valid values are ‘When rotating’, ‘Always’, and ‘Never’.
	Valid values: rotating
	always
	never
	Default value: never
<i>maxmeshwidth</i>	Maximum thickness of scaled mesh lines.
	Valid values: reals
	Default value: 1
	Minimum: 0.5
	Maximum: 10.0
<i>meshwidth</i>	Thickness of mesh lines.
	Valid values: integers
	Default value: 2
	Minimum: 1
	Maximum: 20
<i>minmeshwidth</i>	Minimum thickness of scaled mesh lines.

Chapter 5: Commands

Valid values:	reals
Default value:	0.5
Minimum:	0.1
Maximum:	0.5

pairedcolor1

The default first color for a paired surface

Valid values:	black gray dark_blue blue light_blue aquamarine turquoise spring_green dark_green green lime_green yellow_green yellow orange maroon red pink plum magenta blue_violet white
Default value:	blue

pairedcolor2

The default second color for a paired surface

Valid values: black
 gray
 dark_blue
 blue
 light_blue
 aquamarine
 turquoise
 spring_green
 dark_green
 green
 lime_green
 yellow_green
 yellow
 orange
 maroon
 red
 pink
 plum
 magenta
 blue_violet
 white

Default value: **red**

randomstipple

Specifies when to use random stipple mask for lower-quality transparency. The default is to only use the random stipple pattern for interlaced stereo. Valid values are ‘interlaced’, ‘always’, and ‘never’.

Valid values: interlaced
 always
 never

Default value: **interlaced**

scalemeshwidth

Change mesh width when zooming.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

showtable Set the flag of showing surface manager (table).

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

singlecolor The default color for a standard surface.

Valid values: black
 gray
 dark_blue
 blue
 light_blue
 aquamarine
 turquoise
 spring_green
 dark_green
 green
 lime_green
 yellow_green
 yellow
 orange
 maroon
 red
 pink
 plum
 magenta
 blue_violet
 white
 Default value: **gray**

singlelayertransparency

Enables single layer transparency.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

smoothing

Sets number of surface color smoothing iterations.

Valid values: integers
 Default value: **3**
 Minimum: 0
 Maximum: 20

style

Sets the drawing style.

Valid values: solid
 mesh
 dots
 Default value: **solid**

usecheesy

Off (the default) uses the high-quality but possibly slower-to-draw transparency. On uses the lower-quality cheesy translucency, but in some cases this will be significantly faster to draw. Cheesy translucency mimics transparency by mapping level of transparency to amount of holes in the surface.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

surfacepbrsp

Creates a new Electrostatic Potential Surface for the current workspace.

Syntax:

```
surfacepbrsp coarsedepth=⟨n⟩ coarsespacing=⟨x⟩
               espmayscale=⟨x⟩ focusdepth=⟨n⟩ focusspacing=⟨x⟩
               grid_spacing=⟨x⟩ isosurfaces=yes | no isovalue=isovalue1 |
isovalue2 | isovalue3 | isovalue4 mgextension=⟨x⟩
               pbrspmap=yes | no soluteconstant=⟨x⟩ solventconstant=⟨x⟩
               solventradius=⟨x⟩ temperature=⟨x⟩ transparency=⟨x⟩ Name
of the surface
```

Options:

coarsedepth

This is the MG depth (nlev) for PBE coarse grid.

Valid values: integers
 Default value: **2**
 Minimum: 1

coarsespacing

This option sets the grid spacing for the PBE coarse grid.

Valid values: reals
 Default value: **1.2**
 Minimum: 0.01

espmayscale

This option sets the scale of ESP mapping range.

Valid values: reals
 Default value: **0.25**
 Minimum: 0.01
 Maximum: 1.0

focusdepth This is the MG depth (nlev) in PBE. for PBE focused grid.

Valid values: integers
 Default value: **4**
 Minimum: 1

focusspacing

This option sets the grid spacing for the PBE focused grid.

Valid values: reals

Chapter 5: Commands

Default value: **0.8**
Minimum: 0.01

grid_spacing

This option sets the grid spacing for the surface.

Valid values: reals
Default value: **0.8**
Minimum: 0.01

isosurfaces

This option indicates whether a paired (+ve and -ve isovalue) surfaces will be generated from the PB calculation.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

isovalue

This is the initial isovalue in PBE.

Valid values: isovalue1
isovalue2
isovalue3
isovalue4
Default value: **isovalue1**

mgextension

This option sets the grid extension for the PBE.

Valid values: reals
Default value: **5**
Minimum: 0.01

pbrspmap

This option indicates whether the PBRSP from the PB calculation will be mapped onto a molecular surface.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

soluteconstant

This is the dielectric constant for solute.

Valid values: reals
Default value: **1**
Minimum: 1.0

solventconstant

This is the dielectric constant for solvent.

Valid values: reals
Default value: **80**
Minimum: 1.0

solventradius

This is the solvent radius.

Valid values:	reals
Default value:	1.4
Minimum:	0.01

temperature

This is the temperature in PBE.

Valid values:	reals
Default value:	298
Minimum:	0.01

transparency

The transparency option controls how transparent the surface appears. A value of 100 means 100 percent transparent. 0 means completely opaque.

Valid values:	reals
Default value:	0
Minimum:	0.0
Maximum:	100.0

Operands:

Name of the surface

Name of the Electrostatic Potential Surface to be created

surfacebrspset

Defines a set of atoms for which a surface can be created for with the surfacebrspset command.

Syntax:

surfacebrspset \langle ASL \rangle

Operands:

\langle ASL \rangle

The ASL expression which defines the atoms the surface will be created for.

surfacereaname

Renames the given surface.

Syntax:

surfacerename *entry*=⟨text⟩ *newname*=⟨text⟩ ⟨surface⟩

Options:

entry The entry that the surface belongs to.

Valid values: text strings

Default value:

newname The new name for the surface.

Valid values: text strings

Default value:

Operands:

⟨surface⟩

The name of the surface to rename.

surfaceresolution

Set surface displaying resolution options.

Syntax:

surfaceresolution *lowresstyle*=solid | mesh | dot
resolution=high | low

Options:

lowresstyle

Specifies drawing style for solid surfaces when rotating. Valid values are ‘solid’, ‘mesh’, and ‘dot’.

Valid values: solid
 mesh
 dot

Default value: **solid**

resolution Set the resolution for drawing solid surfaces to low or high .

Valid values: high
 low

Default value: **high**

surfacescheme

Sets the color scheme for the given surface.

Syntax:

```
surfacescheme color=<text> colorramp=<text>
               datarange=on_surface | entire_volume defaultcolor=<n>
               dotradius=<x> entry=<text> espcolorramp=<text>
               espmax=<x> espmin=<x> iterations=<n> linewidth=<x>
               mapmax=<x> mapmin=<x> negativecolor=<text>
               scheme=<text> schemevolume=<text> showlegend=yes | no
               smooth=yes | no <surface>
```

Options:

<i>color</i>	Sets the color for a constant color. Valid values: text strings Default value: gray
<i>colorramp</i>	Sets the color ramp for the surface when color scheme is Map Values From Volume. Valid values: text strings Default value: redwhiteblue
<i>datarange</i>	Sets the data range to which we apply color ramp. Valid values: on_surface entire_volume Default value: on_surface
<i>defaultcolor</i>	Sets default constant color for a new surface. Valid values: integers Default value: 2 Minimum: 1 Maximum: 21
<i>dotradius</i>	Sets radius for surface dots, unit is Angstroms. Valid values: reals Default value: 0.1 Minimum: 0.0001
<i>entry</i>	The entry that the surface belongs to. Valid values: text strings Default value:
<i>espcolorramp</i>	Sets the color ramp for the surface when color scheme is electro-static potential (ESP).

Chapter 5: Commands

	Valid values: text strings Default value: redwhiteblue
<i>espmax</i>	Sets maximum ESP value for ESP color scheme. Valid values: reals Default value: 0.3 Minimum: 0.001
<i>espmin</i>	Sets minimum ESP value for ESP color scheme. Valid values: reals Default value: -0.3 Maximum: -0.001
<i>iterations</i>	Specify the number of iterations for Laplacian smoothing. Valid values: integers Default value: 35 Minimum: 0 Maximum: 100
<i>linewidth</i>	Sets line width for surface mesh, unit is pixel. Valid values: reals Default value: 0.025 Minimum: 0.0001
<i>mapmax</i>	Sets maximum mapping value. Valid values: reals Default value: 0
<i>mapmin</i>	Sets minimum mapping value. Valid values: reals Default value: 0
<i>negativecolor</i>	Sets the color for the negative surface. Valid values: text strings Default value: white
<i>scheme</i>	Sets the color scheme. Valid values: text strings Default value:
<i>schemevolume</i>	Sets the color scheme volume name. Valid values: text strings Default value:
<i>showlegend</i>	Sets whether or not show legends in workspace.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

smooth Turn on/off surface smoothing processing.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

Operands:

⟨ surface ⟩

The name of the surface to set the color scheme for.

surfacesetcomment

Changes the comments of a given surface.

Syntax:

surfacesetcomment *comment*=⟨ text ⟩ *entry*=⟨ text ⟩ ⟨ surface ⟩

Options:

comment The new comments for the surface.

Valid values: text strings
 Default value:

entry The entry that the surface belongs to.

Valid values: text strings
 Default value:

Operands:

⟨ surface ⟩

The name of the surface to change the comments for.

surfacesetisovalue

Changes the isovalue of a given surface.

Syntax:

surface**setiso****value** *entry*=⟨text⟩ *iso***value**=⟨x⟩ ⟨surface⟩

Options:

entry The entry that the surface belongs to.

Valid values: text strings

Default value:

*iso***value** The new iso**value** for the surface.

Valid values: reals

Default value: **0.3**

Operands:

⟨surface⟩

The name of the surface to change the iso**value** for.

surface**settings**

This keyword is used to set various options associated with surface types.

Syntax:

surface**settings** *default_color*=black | gray | dark_blue | blue |
light_blue | aquamarine | turquoise | spring_green | dark_green
| green | lime_green | yellow_green | yellow | orange | maroon
| red | pink | plum | magenta | blue_violet | white
default_paired_color=black | gray | dark_blue | blue |
light_blue | aquamarine | turquoise | spring_green | dark_green
| green | lime_green | yellow_green | yellow | orange | maroon
| red | pink | plum | magenta | blue_violet | white
⟨model_name⟩

Options:

default_color

The default color for the surface.

Valid values:	black gray dark_blue blue light_blue aquamarine turquoise spring_green dark_green green lime_green yellow_green yellow orange maroon red pink plum magenta blue_violet white
Default value:	yellow_green

default_paired_color

The default color for the paired surface.

Valid values:	black gray dark_blue blue light_blue aquamarine turquoise spring_green dark_green green lime_green yellow_green yellow orange maroon red pink plum magenta blue_violet white
Default value:	red

Operands:

`<model.name>`

The name of a surface type which has settings like default color and default paired color.

surfaceasetviewasl

This command sets a surface to only display the portions of the surface within a given distance of a given ASL.

Syntax:

```
surfaceasetviewasl distance=<x> entry=<text> surface=<text>  
use=yes | no ASL
```

Options:

<i>distance</i>	The distance to the ASL.
Valid values:	reals
Default value:	5

<i>entry</i>	The entry that the surface belongs to. Valid values: text strings Default value:
<i>surface</i>	The name of the surface to operate on. Valid values: text strings Default value:
<i>use</i>	Determines whether or not to use the view by ASL. Valid values: boolean (true false; yes no; y n; on off) Default value: true

Operands:

ASL

The asl which controls the visible surface.

surfacesplitpair

Split paired surfaces into two separate surfaces.

Syntax:

surfacesplitpair *entry*=⟨text⟩ *volume*=⟨text⟩ ⟨surface⟩

Options:

<i>entry</i>	The entry name of the entry that the surface belongs to. Valid values: text strings Default value:
<i>volume</i>	The name of the volume that the surface belongs to. Valid values: text strings Default value:

Operands:

⟨surface⟩

The name of the surface to be splitted.

surfacestyle

Sets the drawing style for the given surface.

Syntax:

surfacestyle *entry*=⟨text⟩ *style*=solid | mesh | dots ⟨surface⟩

Options:

entry The entry that the surface belongs to.

Valid values: text strings

Default value:

style Sets the drawing style.

Valid values: solid
 mesh
 dots

Default value: **solid**

Operands:

⟨surface⟩

The name of the surface to set the style for.

surfacetransparency

Sets the transparency for the given surface.

Syntax:

surfacetransparency *entry*=⟨text⟩ *transparency*=⟨x⟩ ⟨surface⟩

Options:

entry The entry that the surface belongs to.

Valid values: text strings

Default value:

transparency Sets the transparency percentage.

Valid values: reals

Default value: **0**

Minimum: 0.0

Maximum: 100.0

Operands:

⟨surface⟩

The name of the surface to set transparency on.

surfacetransparencyback

Sets the transparency for the given surface (back).

Syntax:

```
surfacetransparencyback entry=⟨text⟩ transparency=⟨x⟩
                        ⟨surface⟩
```

Options:

entry The entry that the surface belongs to.
 Valid values: text strings
 Default value:

transparency
 Sets the transparency percentage.
 Valid values: reals
 Default value: **0**
 Minimum: 0.0
 Maximum: 100.0

Operands:

⟨surface⟩

The name of the surface to set transparency on.

surfacevdw

Creates a new vdW surface for the current workspace.

Syntax:

```
surfacevdw asl=⟨text⟩ atom_radius_scaling=⟨x⟩ context=entry |
molecule | workspace | none | asl grid_spacing=⟨x⟩
transparency=⟨x⟩ Name of the surface
```

Options:

asl This option sets the ASL specification for the atoms which define the surface.

Chapter 5: Commands

Valid values: text strings
Default value:

atom_radius_scaling

This is scaling of VdW radii.

Valid values: reals
Default value: **1**
Minimum: 0.5
Maximum: 10.0

context

This option sets the context for generating a surface.

Valid values: entry
molecule
workspace
none
asl
Default value: **entry**

grid_spacing

This option sets the grid spacing for the surface.

Valid values: reals
Default value: **0.8**
Minimum: 0.01

transparency

The transparency option controls how transparent the surface appears. A value of 100 means 100 percent transparent. 0 means completely opaque.

Valid values: reals
Default value: **0**
Minimum: 0.0
Maximum: 100.0

Operands:

Name of the surface

Name of the vdW surface to be created

surfacevdwcontext

Defines a set of atoms for which a surface will be clipped against with the surfacevdw command.

Syntax:

surfacevdwcontext \langle ASL \rangle

Operands:

\langle ASL \rangle

The ASL expression which defines the atoms the surface will be clipped against.

surfacevdwset

Defines a set of atoms for which a surface can be created for with the surfacevdw command.

Syntax:

surfacevdwset \langle ASL \rangle

Operands:

\langle ASL \rangle

The ASL expression which defines the atoms the surface will be created for.

surfaceviewaslset

Defines a set of atoms for the View by ASL property which will be used by the surfacesetviewasl command.

Syntax:

surfaceviewaslset \langle ASL \rangle

Operands:

\langle ASL \rangle

The ASL expression which defines the atoms that will be used for the View by ASL property.

swapptworkspace

Swaps the project table and the Workspace.

Syntax:

`swapptworkspace`

symmetrizeworkspace

Symmetrizes the workspace finding the point groups.

Syntax:

`symmetrizeworkspace tolerance=<x> analyze|update`

Options:

<i>tolerance</i>	Tolerance for finding the pointing groups and symmetrizing the workspace
Valid values:	reals
Default value:	0.04
Minimum:	0.04
Maximum:	1.00

Operands:

`analyze|update`

If `analyze`, finds the point groups for the current tolerance. If `update`, symmetrize the Workspace with the current tolerance.

system

Execute a command from the system.

Syntax:

`system <command>`

Operands:

`<command>`

A command which is to be executed from the current shell.

tablealigncolumn

Set the alignment of the specified column in the current table.

Syntax:

```
tablealigncolumn alignment=left | center | right <columnname>
                  <alignment>
```

Options:

alignment The alignment to be set for the column. Valid values are “left”, “center”, or “right”.

Valid values:	left
	center
	right
Default value:	left

Operands:

<columnname> <alignment>

The name of the column to align. This is the name displayed in the column header.

tablecopyprop

Set property values for selected project entries from a value associated with entry row in the current Project Table, and indirectly associated with the entry. This is similar to `entrycopyprop` command, but the from property value comes from a table row or from an entry group associated with the table row, not from the entry itself. The to property value will be copied from the from property for the entries selected by the ESL expression. If the to property does not already exist, it will be created. When the data type for the two properties is not the same, a reasonable attempt is made to convert between them. If a selected entry has no value for the from property, or the conversion from a string value fails, the value is generally cleared for the to value for that entry. Entry names are never cleared.

Syntax:

```
tablecopyprop from=<text> to=<text> <ESL>
```

Options:

<i>from</i>	<p>The user name of the table property which provides the values to be copied. Currently supported table properties are Group Title and Row .</p> <p>Valid values: text strings</p> <p>Default value:</p>
<i>to</i>	<p>The name of the property to be modified. This can be either the user name or the m2io data name for the property, if the to property exists. If the to property does not exist, it will be created using the data type of the from property, with user as the author. If the property is the entry name (e.g. Entry Name or s_m_entry_name), an entryrename with replace=no will be done. If the property is the 'included in Workspace' property (In or b_m_entry_in_workspace), then entrywsinclude or entrywsexclude will be done, if needed.</p> <p>Valid values: text strings</p> <p>Default value:</p>

Operands:

⟨ ESL ⟩

The ESL expression specifies for which entries the property values are to be copied.

tablegroupsort

Reorder the groups in the specified table.

Syntax:

```
tablegroupsort field=⟨text⟩ order=ascending | descending  
                  source=group | firstentry ⟨table⟩ all | ⟨group name list⟩
```

Options:

<i>field</i>	<p>The name of the group/entry property to be used as sort-key attribute. Use grouptitle to reorder groups based on group's title. Use property name to reorder groups based on first entry's property value.</p> <p>Valid values: text strings</p> <p>Default value: grouptitle</p>
<i>order</i>	<p>This option sets order of sorted values to be either ascending or descending.</p>

	Valid values:	ascending descending
	Default value:	ascending
<i>source</i>	This option determines whether the 'field' option takes 'group property name' or 'entry property name'. Valid values are group or firstentry .	
	Valid values:	group firstentry
	Default value:	group

Operands:

⟨table⟩ all|⟨group name list⟩

The name of the table to sort. Comma separated list of group names. Or 'all' to sort all groups.

tableresizecolumn

Set the width of the specified column in the current table.

Syntax:

tableresizecolumn ⟨columnname⟩ ⟨width⟩

Operands:

⟨columnname⟩ ⟨width⟩

The name of the column to resize. This is the name displayed in the column header. The width to set the column to. This is the new width for the column.

tablesort

Sort current entry selection in the specified table.

Syntax:

tablesort *field*=⟨text⟩ *order*=ascending | descending ⟨table⟩

Options:

<i>field</i>	The name of the property to be sorted on. This option has been deprecated in favor of the <code>tablesortfields</code> command. This option does nothing now. Valid values: text strings Default value: s_m_entry_name
<i>order</i>	This option sets order of sorted values to be either ascending or descending. This option has been deprecated in favor of the <code>tablesortfields</code> command. This option does nothing now. Valid values: ascending descending Default value: ascending

Operands:

⟨table⟩

The name of the table to sort.

tablesortall

Sort all rows in the specified table.

Syntax:

tablesortall ⟨table⟩

Operands:

⟨table⟩

The name of the table to sort.

tablesortfields

Sets multiple fields for sorting.

Syntax:

```
tablesortfields groupbypropertyname=⟨text⟩
    groupentriesbeforesort=yes | no selectoption=nentries |
    npercentofentries | nentriesineachgroup |
    npercentofentriesineachgroup selecttop=yes | no
    sortentries=yes | no sortgroups=yes | no topn=⟨n⟩ ⟨field⟩
    ⟨ascending|descending⟩
```

Options:

groupbypropertyname

If *groupentriesbeforesort* is true then entries will be grouped based on this property.

Valid values: text strings

Default value: **Title**

groupentriesbeforesort

This option sets if we should group the entries before sorting.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

selectoption

An option which controls the selection of entries after sorting.

Valid values: nentries
npercentofentries
nentriesineachgroup
npercentofentriesineachgroup

Default value: **nentries**

selecttop

This option sets if we should select the given number of entries from the top, in the project table, after sorting the rows.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

sortentries

Sorts the entries when *tablesort* or *tablesortall* command is executed.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

sortgroups

Sorts the groups when *tablesort* or *tablesortall* command is executed.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

topn

This option along with *selecttop* option sets number of entries to be selected from the top, in the project table, after sorting the rows.

Valid values: integers
Default value: **5**

Operands:

⟨field⟩ ⟨ascending|descending⟩

The fields to sort together with the sort order: ascending or descending.

tablesortgroupfields

Sets multiple fields for sorting groups.

Syntax:

```
tablesortgroupfields sortgroupoption=grouptitle |  
firstentrypropertysortfield | firstentrypropertygroupsortfield
```

Options:

sortgroupoption

An option which controls the sorting of groups.

Valid values: group**title**
firstentryproperty**sortfield**
firstentryproperty**groupsortfield**
Default value: **group**title****

tableunselectnonsubset

Unselects all entries which are not in the current table's subset.

Syntax:

```
tableunselectnonsubset
```

targetentry

Designate the target entry for Workspace entry feedback, subject to the entryfeedbackshow and feedbackproperties preferences.

Syntax:

targetentry *target*=⟨text⟩

Options:

target This option sets the name (ID string) of the target entry. If there is no target entry, this string should be empty. Generally, the target entry should be set by `eplayer` commands and by `entrywsincludeonly` commands.

Valid values: text strings

Default value:

tile

Spread the on-screen entries out in a tile-pattern.

Syntax:

tile

tileapplyview

Set the view for all tiles from that of the specified Workspace tile.

Syntax:

tileapplyview ⟨tile_name⟩

Operands:

⟨tile_name⟩

The name which identifies the contents of the tile. This name will be the index of an entry group (greater than 0), or 0: followed by the entry ID for an ungrouped entry (or Scratch entry). A tile named 0: (with no entry ID) does not belong to any entry or entry group. In tile by surface mode, the entry ID will be followed by the name of the special surface for that tile.

tileapplyzoom

Set the zoom for all tiles from that of the specified Workspace tile.

Syntax:

tileapplyzoom <tile_name>

Operands:

<tile_name>

The name which identifies the contents of the tile. This name will be the index of an entry group (greater than 0), or 0: followed by the entry ID for an ungrouped entry (or Scratch entry). A tile named 0: (with no entry ID) does not belong to any entry or entry group. In tile by surface mode, the entry ID will be followed by the name of the special surface for that tile.

tilemode

Control over mode which effectively creates separate workspaces for different entries.

Syntax:

tilemode *max_columns*=<n> *mode*=entry | group | surface
tile=yes | no *transformall*=yes | no

Options:

max_columns

Specifies the maximum number of tile columns.

Valid values: integers

Default value: **4**

Minimum: 1

Maximum: 10

mode

Specifies method used to assign entries to tiles

Valid values: entry
group
surface

Default value: **entry**

tile

A boolean option which determines whether to display entries in tile mode.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

transformall

A boolean option which determines whether to apply transforms globally (on) or to tiles individually (off).

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

tilemove

Move the tile identified by the first operand to the location of the tile identified by the second operand.

Syntax:

```
tilemove <move_tile_name> <to_tile_name>
```

Operands:

```
<move_tile_name> <to_tile_name>
```

The names of two Workspace tiles, the first of which is to be moved to the location currently occupied by the second. This name will be the index of an entry group (greater than 0), or 0: followed by the entry ID for an ungrouped entry (or Scratch entry). A tile named 0: (with no entry ID) does not belong to any entry or entry group. In tile by surface mode, the entry ID will be followed by the name of the special surface for that tile.

tilereset

Reset transform for the specified Workspace tile to that which was in effect the last time the current tile mode was entered.

Syntax:

```
tilereset <tile_name>
```

Operands:

```
<tile_name>
```

The name which identifies the contents of the tile. This name will be the index of an entry group (greater than 0), or 0: followed by the entry ID for an ungrouped entry (or Scratch entry). A tile named 0: (with no entry ID) does not belong to any entry or entry group. In tile by surface mode, the entry ID will be followed by the name of the special surface for that tile.

tileresetall

Reset transform for all tiles in Workspace to that which was in effect the last time the current tile mode was entered. This also sets the global transform.

Syntax:

tileresetall

timingsetup

Set timing experiment variables.

Syntax:

timingsetup *duration*=⟨x⟩ *file*=⟨text⟩ *period*=⟨x⟩

Options:

duration This option sets the maximum duration (in seconds) for the timing experiment. A value of 0.0 means that there is no limit on the time, and that the timing experiment will continue until a **timingstop** command is done, **maestro** quits, or a new **timingstart** command is done.

Valid values: reals
Default value: **0**

file The timing data file.

Valid values: text strings
Default value:

period This option sets the sampling period (in seconds) for accumulating and reporting timing results. A value of 0.0 means that there is no periodic sampling, so that results are accumulated for the entire duration of the timing experiment and reported at the end.

Valid values: reals
Default value: **0**
Minimum: 0.0

timingstart

Start a timing experiment with the settings from the last `timingsetup` command. This first stops any timing experiment currently in progress (from a previous `timingstart` command).

Syntax:

```
timingstart
```

timingstop

Stop any timing experiment currently in progress (from a previous `timingstart` command). Timing results for a partially completed sampling period are not reported.

Syntax:

```
timingstop
```

toggleseeditatomboolean

Toggles the boolean property for the given row in the atom-property table.

Syntax:

```
toggleseeditatomboolean property=⟨text⟩ ⟨row number⟩
```

Options:

property The M2IO data name of the property to toggle
 Valid values: text strings
 Default value:

Operands:

⟨row number⟩

The operand is the row number to toggle.

toggleselecteditatom

Toggles the selection of the given row in the edit atom-property table.

Syntax:

toggleselecteditatom \langle row number \rangle

Operands:

\langle row number \rangle

The operand is the row number to toggle-select.

toggleselecteditproperty

Toggles the selection of the given property in the edit property list.

Syntax:

toggleselecteditproperty \langle property \rangle

Operands:

\langle property \rangle

The operand is a property from the edit property list.

torsioncheck

Specifies four atoms which define a torsion to be checked during a conformational search.

Syntax:

torsioncheck *maximum*= \langle x \rangle *minimum*= \langle x \rangle \langle atom1 \rangle \langle atom2 \rangle
 \langle atom3 \rangle \langle atom4 \rangle

Options:

maximum The maximum allowed value for the torsion. during the torsion check.

Valid values: reals
Default value: **180**
Minimum: 0.0

minimum The minimum allowed value for the torsion check (Degrees)

Valid values: reals
Default value: **0**

Operands:

`<atom1> <atom2> <atom3> <atom4>`

The numbers of four atoms which define a torsion angle to be checked during the conformational search. Note that specifying a-b-c-d is the same as specifying d-c-b-a.

torsiongroup

Sets the current torsion group and the conformation within that group

Syntax:

`torsiongroup <group_name> <conformation_name>`

Operands:

`<group_name> <conformation_name>`

The first operand must be the name of a torsion group within the current fragment mode. The second operand must be the name of a conformation within that torsion group.

trajectoryexport

Export selected frames from the current project to a file or to PT

Syntax:

`trajectoryexport`

trajectoryexportatoms

Specify the set of atoms to include in each exported frame

Syntax:

`trajectoryexportatoms <ASL>`

Operands:

`<ASL>`

A string in the atom specification language.

trajectoryexportmovieoptions

Sets the options for the export movie for trajectory player.

Syntax:

```
trajectoryexportmovieoptions format=avi | mpeg | quicktime  
                             frames=selected | all
```

Options:

<i>format</i>	This option sets the format of movie file. Valid values: avi mpeg quicktime Default value: mpeg
<i>frames</i>	This option sets the number of frames to be exported. The two options are: export selected frames or all frames Valid values: selected all Default value: all

trajectoryexportoptions

Sets the options for the export structure for trajectory player.

Syntax:

```
trajectoryexportoptions exportto=pt | file filename=⟨text⟩  
                       frames=selected | currentsingle | currentmultiple
```

Options:

<i>exportto</i>	This option sets the destination of the frame to be written. Valid values are 'pt' or 'file'. Valid values: pt file Default value: file
<i>filename</i>	If the frames are to be exported to a file, specifies filename to be used for the file to be exported Valid values: text strings Default value:

frames This option sets the number of frames to be exported. The three options are: export selected frames or current frame only as single entry or current frame only as multiple entries.

Valid values: selected
 currentsingle
 currentmultiple

Default value: **currentsingle**

trajectoryplayerdisplayatoms

Specify the set of atoms to be displayed in each frame

Syntax:

trajectoryplayerdisplayatoms \langle ASL \rangle

Operands:

\langle ASL \rangle

A string in the atom specification language.

trajectoryplayergoto

Go to the specified frame in the ordered sequence of trajectory frames.

Syntax:

trajectoryplayergoto \langle frame_number \rangle

Operands:

\langle frame_number \rangle

The frame number within the ordered sequence of trajectory frames, which is to be included in the workspace. This frame is recorded as an option to the **trajectoryplayersettings** command.

trajectoryplayergotofirst

Go to the start frame in the ordered sequence of trajectory frames.

Syntax:

trajectoryplayergotofirst

trajectoryplayergotolast

Go to the end frame in the ordered sequence of trajectory frames.

Syntax:

trajectoryplayergotolast

trajectoryplayersettings

Set trajectoryplayer state variables.

Syntax:

trajectoryplayersettings *arep*=⟨n⟩ *brep*=⟨n⟩ *crep*=⟨n⟩
drawfaster=yes | no *end*=⟨n⟩ *entry*=⟨text⟩ *frame*=⟨n⟩
frameduration=⟨x⟩ *hideclipping*=yes | no *playmode*=loop |
reverse | once *position*=none | superimpose | center |
workspace *referframe*=⟨n⟩ *showaxes*=yes | no
showonlyspecified=yes | no *showsimbox*=yes | no
showtimes=yes | no *smoothing*=⟨n⟩ *start*=⟨n⟩ *step*=⟨n⟩
update_ssa=yes | no

Options:

<i>arep</i>	This option sets the number of replications in the a direction. Valid values: integers Default value: 1 Minimum: 1
<i>brep</i>	This option sets the number of replications in the b direction. Valid values: integers Default value: 1 Minimum: 1
<i>crep</i>	This option sets the number of replications in the c direction. Valid values: integers Default value: 1 Minimum: 1

<i>drawfaster</i>	<p>If true, user lower quality drawing to speed up continuous trajectory play.</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: false</p>
<i>end</i>	<p>This option sets the number of the trajectory end frame.</p> <p>Valid values: integers</p> <p>Default value: 0</p> <p>Minimum: 0</p>
<i>entry</i>	<p>This option sets the name (entry ID) of the entry whose trajectory is displayed.</p> <p>Valid values: text strings</p> <p>Default value:</p>
<i>frame</i>	<p>This option sets the number of the trajectory frame being displayed.</p> <p>Valid values: integers</p> <p>Default value: 0</p> <p>Minimum: 0</p>
<i>frameduration</i>	<p>This option determines the minimum duration, in seconds, of each displayed (entry) trajectory frame during continuous play. The actual frame duration may be longer than the specified value, due to time required for drawing and screen update.</p> <p>Valid values: reals</p> <p>Default value: 0</p> <p>Minimum: 0.0</p> <p>Maximum: 10.0</p>
<i>hideclipping</i>	<p>If true, hide clipping planes during continuous trajectory play.</p> <p>TRAJECTORYEXPORTOPTIONS</p> <p>Valid values: boolean (true false; yes no; y n; on off)</p> <p>Default value: true</p>
<i>playmode</i>	<p>This option sets the mode for continuous play. Valid values are “loop”, “reverse”, or “once”. These cause play to continue, change direction, or stop, respectively, when reaching either end of the current trajectory frame range.</p> <p>Valid values: loop reverse once</p> <p>Default value: once</p>

- position* This option sets the mode for repositioning frame coordinates. Valid values are “none”, “superimpose”, “workspace”, or “center”. For the value “none”, do not adjust positions. For the value “superimpose”, superimpose current trajectory frame onto specified reference frame, aligning atoms that match the current superimpose ASL expression. For the value “workspace”, superimpose current trajectory frame onto specified workspace atoms, aligning atoms that match the current superimpose ASL expression. For the value “center”, recenter selected molecules back to simulation box, using the same superimpose ASL expression to select atoms in the molecules of interest.
- Valid values: none
superimpose
center
workspace
- Default value: **center**
- referframe* This option sets the number of the trajectory reference frame for superimpose.
- Valid values: integers
- Default value: **0**
- Minimum: 0
- showaxes* If true, display x, y, z axes for trajectory.
- Valid values: boolean (true|false; yes|no; y|n; on|off)
- Default value: **false**
- showonlyspecified*
- If true, display only the atoms specified by trajectoryplayerdisplayatoms ASL (and optionally those within a specified distance) when a trajectory frame is shown. If false, display whichever atoms have currently been made visible in Maestro.
- Valid values: boolean (true|false; yes|no; y|n; on|off)
- Default value: **false**
- showsimbox*
- If true, display simulation box for trajectory.
- Valid values: boolean (true|false; yes|no; y|n; on|off)
- Default value: **true**
- showtimes* If true, perform timing tests while showing trajectory frames, and display timing results.
- Valid values: boolean (true|false; yes|no; y|n; on|off)
- Default value: **false**
- smoothing* This option sets the trajectory smoothing.

	Valid values: integers
	Default value: 1
	Minimum: 1
<i>start</i>	This option sets the number of the trajectory start frame.
	Valid values: integers
	Default value: 0
	Minimum: 0
<i>step</i>	This option sets the number of the trajectory frames per step (skipping 1 less than this).
	Valid values: integers
	Default value: 1
	Minimum: 1
<i>update_ssa</i>	If true, recalculate secondary structure assignment for each frame.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false

trajectoryplayerstepahead

Go to the next frame in the ordered sequence of trajectory frames, if there is one after the frame specified in the trajectoryplayersettings command. The step value in the trajectoryplayersettings command determines how many frames to go forward.

Syntax:

```
trajectoryplayerstepahead
```

trajectoryplayerstepback

Go to the previous frame in the ordered sequence of trajectory frames, if there is one before the frame specified in the trajectoryplayersettings command. The step value in the trajectoryplayersettings command determines how many frames to go back.

Syntax:

trajectoryplayerstepback

trajectoryplayersuperimposeatoms

Specify the set of atoms used in superimposition

Syntax:

trajectoryplayersuperimposeatoms \langle ASL \rangle

Operands:

\langle ASL \rangle

A string in the atom specification language.

trajectoryresolvemeasurements

Sets the options for the resolving trajectory measurements

Syntax:

trajectoryresolvemeasurements *use*=wsonly | trajonly | both

Options:

use This option sets the preference of the user whether to use trajectory measurements or ws measurements or both

Valid values: wsonly
 trajonly
 both

Default value: **wsonly**

transform

Specify what is to be transformed.

Syntax:

```

transform centerlocal=centroid | atom centerrotation=yes | no
contactdisplay=yes | no gui=none | rotate | translate
hbonddisplay=yes | no includealternate=yes | no
rsensitivity=⟨x⟩ scope=global | local
showrotationcenter=yes | no smootherrotation=yes | no
threshold=⟨n⟩ tsensitivity=⟨x⟩ [reset]

```

Options:

centerlocal

Whether local center is an atom or centroid of group of atoms

Valid values: centroid
 atom

Default value: **centroid**

centerrotation

Set the flag of center rotation.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

contactdisplay

This option determines whether Contact markers will be displayed for local transformations.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

gui

Whether rotation, translation or nothing is performed

Valid values: none
 rotate
 translate

Default value: **rotate**

hbonddisplay

This option determines whether H-bond markers will be displayed for local transformations.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

includealternate

Whether to include alternate positions in local transformations to atom coordinates.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

rsensitivity

Mouse rotation sensitivity. Larger is more sensitive.

Valid values: reals

Chapter 5: Commands

	Default value: 75
	Minimum: 1.0
	Maximum: 500.0
<i>scope</i>	Whether global or local transformations are done
	Valid values: global local
	Default value: global
<i>showrotationcenter</i>	Set the flag of displaying rotate center.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: false
<i>smootherrotation</i>	Set the flag of smoother local rotation.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>threshold</i>	How many pixels a mouse must move from the initial mouse down position before any transformation will occur
	Valid values: integers
	Default value: 5
	Minimum: 1
	Maximum: 30
<i>tsensitivity</i>	[NOTE: This option is no longer used.] Mouse translation sensitivity. Larger is more sensitive.
	Valid values: reals
	Default value: 3.33
	Minimum: 1.0
	Maximum: 500.0

Operands:

[reset]

Specifies that the transformations are to be reset - this must be reset .

translate

Translate in Angstroms whatever is specified in the transform set. This is either global (all atoms) or a local grouping defined using an ASL expression in the transform command.

Syntax:

translate $x=\langle x \rangle$ $y=\langle x \rangle$ $z=\langle x \rangle$

Options:

x	Amount in Angstroms to translate in X
	Valid values: reals
	Default value: 0
y	Amount in Angstroms to translate in Y
	Valid values: reals
	Default value: 0
z	Amount in Angstroms to translate in Z
	Valid values: reals
	Default value: 0

undisplayatom

Undisplay atoms in the set described by the ASL.

Syntax:

undisplayatom $\langle \text{ASL} \rangle$

Operands:

$\langle \text{ASL} \rangle$

A string in the atom specification language which describes the set of atoms which are to be undisplayed.

undo

Undo the effect of the last change on the on-screen structure.

Syntax:

undo

ungroupentries

Ungroups the entries that match the given ESL expression.

Syntax:

ungroupentries \langle ESL \rangle

Operands:

\langle ESL \rangle

\langle ESL \rangle A valid ESL expression to specify which entries are to be ungrouped. Ungroups the entries that match the given ESL expression, and move them all to the end of ungrouped section. In project table, all the ungrouped entries (if any) will be present at the top i.e. before all the groups.

unhookimport

Restores normal Import behavior.

Syntax:

unhookimport

uniquename

With the set of specified atoms, make the atom name unique. This is done by adding a “~N” to each duplicate name where the “N” is a digit which represents how often this name is repeated in the set.

Syntax:

uniquename

uniquepdb

Set unique PDB atom names (by residue) for all atoms that match the ASL specification.

Syntax:

uniquepdb \langle ASL \rangle

Operands:

\langle ASL \rangle

The operand is a valid ASL string that defines the set of atoms that are to have their PDB atom names changed.

update

This is a standard alias for **energyupdate** (see [\[energyupdate\]](#), page 136).

updateribbons

Update existing ribbons with current ribbon style and color scheme.

Syntax:

updateribbons

varymolecule

A command which defines a rotatable/translatable molecule during a conformational search. The molecule is defined by specifying any atom which belongs to it.

Syntax:

varymolecule $rmax=\langle x \rangle$ $rmin=\langle x \rangle$ $tmax=\langle x \rangle$ $tmin=\langle x \rangle$
 \langle atom_number \rangle

Options:

<i>rmax</i>	The maximum value for the molecule rotation.
	Valid values: reals
	Default value: 180
	Minimum: 0.0
<i>rmin</i>	The minimum value for the molecule rotation.
	Valid values: reals
	Default value: 0
	Minimum: 0.0

Chapter 5: Commands

<i>tmax</i>	The maximum value for the molecule translation
	Valid values: reals
	Default value: 1
	Minimum: 0.0
<i>tmin</i>	The minimum value for the molecule translation
	Valid values: reals
	Default value: 0
	Minimum: 0.0

Operands:

$\langle \text{atom_number} \rangle$

The number of an atom which is a member of the molecule which is to be translated or rotated.

varytorsion

Defines a rotatable bond to be used in a conformational search.

Syntax:

varytorsion *maximum*= $\langle x \rangle$ *minimum*= $\langle x \rangle$ $\langle \text{atom1} \rangle$ $\langle \text{atom2} \rangle$

Options:

<i>maximum</i>	The maximum value for the torsional rotation.
	Valid values: reals
	Default value: 180
	Minimum: 0.0
<i>minimum</i>	The minimum value for the torsional rotation.
	Valid values: reals
	Default value: 0
	Minimum: 0.0

Operands:

$\langle \text{atom1} \rangle$ $\langle \text{atom2} \rangle$

The two atom numbers which define a bond to be rotated in a conformational search. These two atoms must have a bond (usually single) between them. Note that specifying a-b is the same as specifying b-a.

vcsaddattachment

Adds an attachment to the core molecule using the given atoms.

Syntax:

vcsaddattachment *atom1*=⟨n⟩ *atom2*=⟨n⟩ ⟨attachment name⟩

Options:

atom1 The atom number of the atom in the original core to set as an attachment point. This is the atom which will be kept.

Valid values: integers

Default value: **1**

Minimum: 1

atom2 The atom number of the atom in the original core to set as an attachment point. This is the atom which will be removed.

Valid values: integers

Default value: **1**

Minimum: 1

Operands:

⟨attachment name⟩

The name of the attachment.

vcsaddcorefromproject

Adds the given entry as a core pose in CombiGlide.

Syntax:

vcsaddcorefromproject ⟨entry name⟩

Operands:

⟨entry name⟩

The entry name.

vcsaddmincapcore

Adds the minimally capped core to the poses for Define Core Poses in CombiGlide.

Syntax:

vcsaddmincapcore

vcsaddoriginalcore

Adds the original core to the poses for Define Core Poses in CombiGlide.

Syntax:

vcsaddoriginalcore

vcscanceldockjob

Cancels the currently running Dock Library job associated results.

Syntax:

vcscanceldockjob

vcsclearreagentfile

Clears the reagent file for the selected rows.

Syntax:

vcsclearreagentfile

vcscombiexportdockingfile

Exports the combinatorial docking results from CombiGlide to the given file.

Syntax:

vcscombiexportdockingfile <file name>

Operands:

<file name>

The file name.

vcscombiexportdockingproject

Exports the combinatorial docking results from CombiGlide.

Syntax:

```
vcscombiexportdockingproject
```

vcscombiexportoptions

Holds the options for exporting combinatorial docking results from CombiGlide.

Syntax:

```
vcscombiexportoptions includereceptor=yes | no  
                      numreagents=⟨ n ⟩
```

Options:

includereceptor

Indicates whether or not to include the receptor in the exported results

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

numreagents

The number of reagents to export

Valid values: integers

Default value: **100**

Minimum: 1

vcsconfiguredocking

Allows the settings of some values that determine how the overall VCS job runs.

Syntax:

```

vcconfiguredocking bondrotation=⟨text⟩ gridfilename=⟨text⟩
gridfiletype=⟨text⟩ inputring=yes | no lig_ccut=⟨x⟩
lig_vscale=⟨x⟩ ligandwithmetal=charged | either | neutral
maxatom=⟨n⟩ maxrotbonds=⟨n⟩ ninvert=yes | no
numreqgroup1=⟨n⟩ numreqgroup2=⟨n⟩ numreqgroup3=⟨n⟩
numreqgroup4=⟨n⟩ penalizeamidebondrotations=yes | no
reqmodegroup1=all | at least reqmodegroup2=all | at least
reqmodegroup3=all | at least reqmodegroup4=all | at least
ringconf=yes | no

```

Options:

bondrotation

The option of amide bond rotation.

Valid values: text strings

Default value: **penal**

gridfilename

The base name for the file which the receptor grid is to be written to or read from.

Valid values: text strings

Default value:

gridfiletype

The type for the file which the receptor grid is to be written to or read from. The type is either zip or grid .

Valid values: text strings

Default value: **zip**

inputring

An option which allows input ring conformation

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

lig_ccut

The partial atomic charge below which ligand atoms are considered to be non-polar and will have their VDW radii scaled.

Valid values: reals

Default value: **0.15**

Minimum: 0.00000001

lig_vscale

The scaling factor for the VDW radii of non-polar ligand atoms.

Valid values: reals

Default value: **0.8**

Minimum: 0.00000001

ligandwithmetal

Controls which ligand atoms can interact with metal sites.

	Valid values:	charged either neutral
	Default value:	charged
<i>maxatom</i>	Any ligands in the input with more than this number of atoms will be skipped.	
	Valid values:	integers
	Default value:	300
	Minimum:	1
	Maximum:	300
<i>maxrotbonds</i>	Any ligands in the input with more than this number of rotatable bonds will be skipped.	
	Valid values:	integers
	Default value:	50
	Minimum:	1
<i>ninvert</i>	An option which allows nitrogen inversions	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	true
<i>numreggroup1</i>	Number of constraints to be required for group 1 in docking.	
	Valid values:	integers
	Default value:	1
	Minimum:	0
	Maximum:	4
<i>numreggroup2</i>	Number of constraints to be required for group 2 in docking.	
	Valid values:	integers
	Default value:	1
	Minimum:	0
	Maximum:	4
<i>numreggroup3</i>	Number of constraints to be required for group 3 in docking.	
	Valid values:	integers
	Default value:	1
	Minimum:	0
	Maximum:	4
<i>numreggroup4</i>	Number of constraints to be required for group 4 in docking.	
	Valid values:	integers

Chapter 5: Commands

Default value: **1**
Minimum: 0
Maximum: 4

penalizeamidebondrotations

An option that penalizes twisted (non-planar) amide bonds.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

reqmodegroup1

The mode determines how to set the number of required constraints for group 1 in docking.

Valid values: all
atleast
Default value: **atleast**

reqmodegroup2

The mode determines how to set the number of required constraints for group 2 in docking.

Valid values: all
atleast
Default value: **atleast**

reqmodegroup3

The mode determines how to set the number of required constraints for group 3 in docking.

Valid values: all
atleast
Default value: **atleast**

reqmodegroup4

The mode determines how to set the number of required constraints for group 4 in docking.

Valid values: all
atleast
Default value: **atleast**

ringconf

An option which allows ring flips

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

vcscoreoptions

Sets the options for Define Core Poses.

Syntax:

```
vcscoreoptions centroidx=⟨x⟩ centroidy=⟨x⟩ centroidz=⟨x⟩
                constrainradius=⟨x⟩ maxglidermsd=⟨x⟩ maxrmsd=⟨x⟩
                poseconstraint=box | sphere | glidecore
```

Options:

centroidx The X-coordinate of the location to constrain the core center of mass to.

Valid values: reals

Default value: **0**

centroidy The Y-coordinate of the location to constrain the core center of mass to.

Valid values: reals

Default value: **0**

centroidz The Z-coordinate of the location to constrain the core center of mass to.

Valid values: reals

Default value: **0**

constrainradius

The radius of the constraint for the core's center of mass.

Valid values: reals

Default value: **5**

Minimum: 1.0

maxglidermsd

The maximum RMSD that the Glide core can move.

Valid values: reals

Default value: **1**

Minimum: 0.0

maxrmsd The maximum RMSD that the core can move.

Valid values: reals

Default value: **2**

Minimum: 0.0

poseconstraint

Controls the allowed placement for core structures in CombiGlide docking.

Valid values: box
 sphere
 glidecore

Default value: **box**

vcscreatedockedlibrary

Start the creation and docking of the library in the Analyze Library step of CombiGlide. This is similar to `vcstrunenumerateddocking`, except that it uses only the selected reagents at each position and does not remove or replace contents of the Combinatorial Screening run.

Syntax:

```
vcscreatedockedlibrary
```

vcsdeleteattachment

Deletes all of the selected attachments.

Syntax:

```
vcsdeleteattachment
```

vcsdeletecore

Deletes the selected cores in the Define Core Poses step in CombiGlide.

Syntax:

```
vcsdeletecore
```

vcsdeleteresults

Deletes the given results.

Syntax:

```
vcsdeleteresults <name>
```

Operands:

<name>

The name of the results file.

vcsdisplayreceptor

This function displays the receptor for the current CombiGlide run in the Workspace.

Syntax:

vcsdisplayreceptor

vcsdocking

Sets docking job options for CombiGlide.

Syntax:

vcsdocking *applystrain*=yes | no *postdock*=yes | no
task=combination | single | combinatorial | enumerated

Options:

applystrain

An option which determines if the CombiGlide docking job will apply strain correction terms.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

postdock

An option which determines if the CombiGlide docking job will do post-docking minimization.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

task

Determines what kind of docking task will be used in CombiGlide docking, 1) user-capped core and 2) minimally capped core.

Valid values: combination
 single
 combinatorial
 enumerated

Default value: **combination**

vcsenumeratedockoptions

Holds the options for enumerate and dock for CombiGlide.

Syntax:

vcseenumeratedockoptions *mode*=⟨text⟩

Options:

mode An option controlling what type of docking will be done after combinatorial enumeration. The allowed values are *combi* , *xp* , *sp* , and *htvs* .

Valid values: text strings

Default value: **combi**

vcsexcludetablerow

Excludes the given row in the structure table in the step from the Workspace.

Syntax:

vcsexcludetablerow ⟨row⟩

Operands:

⟨row⟩

The row number to exclude in the Workspace.

vcsexportdefinition

Stores the current core molecule and attachments in a file.

Syntax:

vcsexportdefinition ⟨file name⟩

Operands:

⟨file name⟩

The name of the file to store the core definition in.

vcsexportresults

Stores the settings from the Filter and Select dialog, and all the results, in a human-readable text file.

Syntax:

vcsexportresults *<file name>*

Operands:

<file name>

The name of the file to store the settings and results in.

vcsexportspreadsheet

Export actives found by libselector to a file for use in a spreadsheet. This can be comma-separated value (.csv) format or tab-delimited format.

Syntax:

vcsexportspreadsheet *delimiter=<text>* *<filename>*

Options:

delimiter This option sets the delimiter to use to separate columns.

Valid values: text strings

Default value: ,

Operands:

<filename>

The name of the file to which actives will be written. If no name is specified, then no export will be done.

vcimportdefinition

Reads a core molecule and attachments from the given file.

Syntax:

vcimportdefinition *<file name>*

Operands:

<file name>

The name of the file to read the core definition from.

vcsimportrun

Incorporates the combiglide results into the Maestro project for the specified run

Syntax:

```
vcsimportrun <directory name> <run name>
```

Operands:

<directory name> <run name>

The first operand is the name of the directory to read the combiglide output files. The second operand is the run name to be imported into maestro

vcsincludeextendtablerow

Extends the rows included in the workspace to include this one.

Syntax:

```
vcsincludeextendtablerow <row>
```

Operands:

<row>

The row number to include in the Workspace.

vcsincludeonlytablerow

Includes only the given row in the structure table in the step into the Workspace.

Syntax:

```
vcsincludeonlytablerow <row>
```

Operands:

<row>

The row number to include in the Workspace.

vcsincludetablerow

Includes the given row in the structure table in the step into the Workspace.

Syntax:

```
vcsincludetablerow <row>
```

Operands:

<row>

The row number to include in the Workspace.

vcsinverttableselection

Inverts the row selection in the first table in the step.

Syntax:

```
vcsinverttableselection
```

vcsoptions

This command holds general options for CombiGlide.

Syntax:

```
vcsoptions mode=<text> sidechainnode=<n> untangle=yes | no
```

Options:

<i>mode</i>	An option controlling what type of job CombiGlide will run.
Valid values:	text strings
Default value:	sidechain

<i>sidechainnode</i>	Which sidechain node we are working on.
Valid values:	integers
Default value:	1
Minimum:	1

<i>untangle</i>	An option which allows post-combgen minimization (for Create Library, enumeration only).
Valid values:	boolean (true false; yes no; y n; on off)
Default value:	true

vcsrefreshstructure

This function refreshes the structure in the Workspace from the current core structure in CombiGlide.

Syntax:

vcsrefreshstructure *viewcappedcore*=yes | no

Options:

viewcappedcore

An option which allows viewing of the minimally capped core, rather than the original core, in the Workspace.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

vcsrenameattachment

Renames the attachment in CombiGlide to the new name.

Syntax:

vcsrenameattachment *row*=⟨n⟩ ⟨new name⟩

Options:

row The row to rename.

Valid values: integers

Default value: **1**

Minimum: 1

Operands:

⟨new name⟩

The new name for the attachment.

vcsrestorerresults

Restores the filter and selection settings, and all associated results.

Syntax:

```
vcsrestorereresults <name>
```

Operands:

<name>

The name of the results file.

vcsruncombinationdocking

Runs a both single-position and combinatorial docking job for CombiGlide.

Syntax:

```
vcsruncombinationdocking
```

vcsruncombinatorialdocking

Runs a combinatorial docking job for CombiGlide.

Syntax:

```
vcsruncombinatorialdocking maxresults=<n> mode=sp | xp  
                        <job name>
```

Options:

maxresults

The number of combinatorial results to return.

Valid values: integers

Default value: **1000**

Minimum: 0

mode

Determines what kind of docking will be taken, 1) SP docking (faster), and 2) CombiGlide XP docking.

Valid values: sp

xp

Default value: **sp**

Operands:

<job name>

The job name.

vcsruncbinatorialselection

Runs a combinatorial selection job for CombiGlide.

Syntax:

```
vcsruncbinatorialselection cmdargs=⟨text⟩
```

Options:

cmdargs The arguments for the reagent selection command.

Valid values: text strings

Default value:

vcsruncreate

Creates the run with the given name.

Syntax:

```
vcsruncreate ⟨run name⟩
```

Operands:

⟨run name⟩

The name of the new run to create.

vcsrundelele

Deletes the current run from the project.

Syntax:

```
vcsrundelele
```

vcsruneenumerateddocking

Runs an enumerated docking job in the Dock Library step of CombiGlide.

Syntax:

vcsruneenumerateddocking

vcsruneopen

Opens the run with the given name.

Syntax:

vcsruneopen \langle run name \rangle

Operands:

\langle run name \rangle

The name of the run to open.

vcsrunerename

Changes the current run's name to the given name.

Syntax:

vcsrunerename \langle run name \rangle

Operands:

\langle run name \rangle

The name to change the current run's name to.

vcsruneasaveas

Saves a copy of the current run under the given name.

Syntax:

vcstrunsaveas \langle run name \rangle

Operands:

\langle run name \rangle

The name of the run to save as.

vcstrunsingledocking

Runs a single position docking job for CombiGlide.

Syntax:

vcstrunsingledocking *usecore*=usercapped | minimallycapped

Options:

usecore Determines what kind of core will be used in single-position docking, 1) user-capped core, 2) minimally capped core.

Valid values: usercapped
 minimallycapped

Default value: **usercapped**

vcstrunsingleselection

Runs a single position selection job for CombiGlide.

Syntax:

vcstrunsingleselection *cmdargs*= \langle text \rangle

Options:

cmdargs The arguments for the reagent selection command.

Valid values: text strings
Default value:

vcssaveresults

Saves the current filter and selection settings, and all associated results.

Syntax:

vcssaveresults $\langle \text{name} \rangle$

Operands:

$\langle \text{name} \rangle$

The name of the results file.

vcselectalltablerows

Selects all rows in the first table in the step.

Syntax:

vcselectalltablerows

vcselectextendablerow

Extends the selection to this row in the table.

Syntax:

vcselectextendablerow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to extend the select to.

vcselectonlyablerow

Selects only this row in the table.

Syntax:

vcselectonlyablerow $\langle \text{row} \rangle$

Operands:

$\langle \text{row} \rangle$

The row number to select only in the table row.

vcsselecttablerow

Selects the given row in the first table in the step.

Syntax:

vcsselecttablerow *<row>*

Operands:

<row>

The row number to select in the table.

vcssetattachmentfile

Sets the reagent file for the given attachment.

Syntax:

vcssetattachmentfile *file*=*<text>* *<attachment name>*

Options:

file The file name of the reagent file to add to the given attachment.
Valid values: text strings
Default value:

Operands:

<attachment name>

The name of the attachment.

vcssetmolecule

Sets the core molecule for the current CombiGlide run to the molecule containing the given atom.

Syntax:

vcssetmolecule *title*=*<text>* *<atom number>*

Options:

title This option sets the title for the core molecule.

Valid values: text strings
Default value: **core**

Operands:

⟨ atom number ⟩

The atom number of the molecule.

vcsetreagentfile

Sets the reagent file for the selected rows.

Syntax:

vcsetreagentfile ⟨ reagent name ⟩

Operands:

⟨ reagent name ⟩

The name of the reagent file.

vcssingleexportdockingfile

Exports the single-position docking results from CombiGlide to the given file.

Syntax:

vcssingleexportdockingfile ⟨ file name ⟩

Operands:

⟨ file name ⟩

The file name.

vcssingleexportdockingproject

Exports the single-position docking results from CombiGlide.

Syntax:

`vcssingleexportdockingproject`

vcssingleexportoptions

Holds the options for exporting single-position docking results from CombiGlide.

Syntax:

`vcssingleexportoptions` *includereceptor*=yes | no
numreagents= $\langle n \rangle$

Options:

includereceptor

Indicates whether or not to include the receptor in the exported results

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

numreagents

The number of reagents to export

Valid values: integers

Default value: **100**

Minimum: 1

vcssorttable

Resort the given CombiGlide table based on the data in the specified column

Syntax:

`vcssorttable` *table*= $\langle n \rangle$ \langle column.name \rangle

Options:

table The table to set.

Valid values: integers

Default value: **10**

Operands:

⟨column_name⟩

The name of the column to be sorted.

vcsstepforward

Moves forward to the next VCS step. Deletes any steps after the current step, then creates the next step, using the data from previous steps.

Syntax:

vcsstepforward

vcsstepgoto

Moves to an existing step in the current project.

Syntax:

vcsstepgoto ⟨step name⟩

Operands:

⟨step name⟩

The name of the step to switch to.

vcsundisplayreceptor

This function undisplays the receptor for the current CombiGlide run in the Workspace.

Syntax:

vcsundisplayreceptor

vcsunselecttablerow

Unselects the given row in the first table in the step.

Syntax:

```
vcsunselecttablerow <row>
```

Operands:

<row>

The row number to unselect in the table.

vcswritecombinationdocking

Writes a both single-position and combinatorial docking job input file for CombiGlide.

Syntax:

```
vcswritecombinationdocking
```

vcswritecombinatorialdocking

Writes a combinatorial docking job input file for CombiGlide.

Syntax:

```
vcswritecombinatorialdocking
```

vcswriteenumerateddocking

Writes an enumerated docking job input file in the Dock Library step of CombiGlide.

Syntax:

```
vcswriteenumerateddocking
```

vcswritesingledocking

Writes a single position docking job input file for CombiGlide.

Syntax:

`vcswritesingledocking`

viewcrystalmates

Set the options for removing and adding crystal mates in Workspace

Syntax:

viewcrystalmates *groupradius*=⟨x⟩ *renamechains*=yes | no
viewallcrystalmates=yes | no *viewcrystalmateswithin*=⟨x⟩
viewmates=yes | no

Options:

groupradius

Controls the group radius used in the crystal mates calculation. This controls the radius within which a symmetric unit is considered to be in contact with the AsU. This is different to the *viewcrystalmateswithin* value which is only control the display and not the generation of the mates.

Valid values: reals
 Default value: **14**
 Minimum: 0.0

renamechains

Enabling this option will automatically rename crystal mates as they are generated.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

viewallcrystalmates

Indicates whether or not to view all crystal mates.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

viewcrystalmateswithin

This feature displays crystal mates within the given distance. This only takes affect if *viewallcrystalmates* is set to false.

Valid values: reals
 Default value: **10**
 Minimum: 0.0

viewmates Indicates whether or not to include crystal mates in Workspace for included project entries that have the required crystal settings.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

viewdelete

Deletes transform data and view data corresponding to the given view name.

Syntax:

viewdelete <view_name>

Operands:

<view_name>

A valid view name existing in the current views table.

viewdeletesected

Delete selected views in the current views table.

Syntax:

viewdeletesected

viewdragselection

Move the current view selection to the specified table row.

Syntax:

viewdragselection <row>

Operands:

<row>

A valid destination row number in the views table.

viewexport

Exports views from the current views table to given file. If non of the views are selected then exports all the views, else exports only selected views.

Syntax:

viewexport <file>

Operands:

<file>

The views file name to export.

viewextendselect

Extend current view selection to encompass the specified view in the table.

Syntax:

viewextendselect <view_name>

Operands:

<view_name>

A valid view name existing in the current views table. A range selection will be done for the views between the currently selected view and the specified view name.

viewextendselectrow

Extend current view selection to encompass the specified view table row.

Syntax:

viewextendselectrow <row>

Operands:

<row>

A row number between 1 and the total number of rows in the views table. A range selection will be done for the views between currently selected view rows and the specified view row number.

viewimport

Imports views from the given file.

Syntax:

```
viewimport <file>
```

Operands:

<file>

The views file name to import.

viewmatrix

Sets workspace rotation matrix to the matrix supplied. When 'inverse' is supplied, it sets the inverse rotation matrix and when 'nocenter', no center rotation matrix. And in the absence of both the keywords, it sets the rotation matrix. The transpose of this matrix is used in a call to `glMultMatrixd` to set the user rotation. Points in molecule coordinate space (x, y, z, 1) are pre-multiplied by this matrix to generate rotated coordinates for display, subject to additional view transforms for perspective or stereo rotation, and mapping into normalized view space and to window or screen pixel space. The `viewmatrix` command takes, as operands, 16 floating point numbers, to define a 4-by-4 viewing matrix in row major order, plus an optional matrix identifier string (inverse or nocenter). Maestro reads each line from the `wsvview` file and processes it as a separate command, silently ignoring any errors. The `viewmatrix` settings are all inter-dependent. If one is set, all three should be set.

Syntax:

```
viewmatrix
```

viewplayergotofirst

Goto the first view in the current views table, for which play state is set.

Syntax:

viewplayergotofirst

viewplayergotolast

Goto the last view in the current views table, for which play state is set.

Syntax:

viewplayergotolast

viewplayerplaybackward

Start playing views backward from the recently restored view in the ordered sequence of play toggled on views in the current views table.

Syntax:

viewplayerplaybackward

viewplayerplayforward

Start playing views forward from the recently restored view in the ordered sequence of play toggled on views in the current views table.

Syntax:

viewplayerplayforward

viewplayersettings

Set viewplayer state variables.

Syntax:

viewplayersettings *playmode=loop | reverse | once*

Options:

playmode This option sets the mode for continuous play. Valid values are “loop”, “reverse”, or “once”. These cause play to continue, change direction, or stop, respectively, when reaching either end of the current play toggle on views.

Valid values: loop
 reverse
 once

Default value: **once**

viewplayerstepahead

Goto the next view in the ordered sequence of play toggled on view in the current views table.

Syntax:

viewplayerstepahead

viewplayerstepback

Goto the previous view in the ordered sequence of play toggled on views in the current views table.

Syntax:

viewplayerstepback

viewplayerstop

Stop playing views/saving views movie.

Syntax:

viewplayerstop

viewrecordmovie

Save a movie with views for which Play is set.

Syntax:

viewrecordmovie $\langle \text{file} \rangle$

Operands:

$\langle \text{file} \rangle$

The file where the views movie will be saved.

viewrename

Renames the given view name.

Syntax:

viewrename $\langle \text{curr_name} \rangle \langle \text{new_name} \rangle$

Operands:

$\langle \text{curr_name} \rangle \langle \text{new_name} \rangle$

A valid view name existing in the current views table. A valid new view name.

viewreset

Resets the viewing transform to original.

Syntax:

viewreset

viewresetnofit

Resets the viewing transform to original but not fit to window.

Syntax:

viewresetnofit

viewrestore

Restores the viewing transform corresponding to the view name.

Syntax:

viewrestore \langle view_name \rangle

Operands:

\langle view_name \rangle

A valid view name existing in the current views table.

viewrestorerow

Restore the transform of the view corresponding to the given row.

Syntax:

viewrestorerow \langle row \rangle

Operands:

\langle row \rangle

A valid row in the views table.

viewsave

Saves the current viewing transform.

Syntax:

viewsave *pause*= \langle n \rangle *play*=yes | no *select*=yes | no \langle view_name \rangle

Options:

<i>pause</i>	Pause value of the view.
Valid values:	integers
Default value:	1
Minimum:	0

	Maximum:	5
<i>play</i>	Play state of the view.	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	true
<i>select</i>	Selection state of the view.	
	Valid values:	boolean (true false; yes no; y n; on off)
	Default value:	true

Operands:

⟨view_name⟩

A valid view name with which the current transform will be saved.

viewselect

Select the given view in the current views table.

Syntax:

viewselect ⟨view_name⟩

Operands:

⟨view_name⟩

A valid view name existing in the current views table.

viewselectonly

Select only the given view in the current views table.

Syntax:

viewselectonly ⟨view_name⟩

Operands:

⟨view_name⟩

A valid view name existing in the current views table.

viewselectonlyrow

Select only the view corresponding to the row in the current views table and unselect remaining views.

Syntax:

```
viewselectonlyrow <row>
```

Operands:

<row>

A valid row in the views table.

viewselectrow

Select the view corresponding to the given row in the current views table.

Syntax:

```
viewselectrow <row>
```

Operands:

<row>

A valid row in the views table.

viewsetpause

Sets the pause value to the given view.

Syntax:

```
viewsetpause <view_name> <pause_value>
```

Operands:

<view_name> <pause_value>

A valid view name existing in the current views table. A valid pause value greater than -1.

viewsetplay

Set the play state of the given view.

Syntax:

```
viewsetplay <view_name>
```

Operands:

<view_name>

A valid view name existing in the current views table.

viewsetrowpause

Set the given pause value to the view corresponding to the row in the current views table.

Syntax:

```
viewsetrowpause <row> <pause>
```

Operands:

<row> <pause>

A row number between 1 and the total number of rows in the views table.

A valid pause value greater than -1.

viewtrajectory

Attempt to view the Chorus trajectory associated with the specified entry.

Syntax:

```
viewtrajectory <entry_name>
```

Operands:

<entry_name>

The name (entry ID) of the entry, whose trajectory data is to be viewed.

viewunselect

Unselect the given view in the current views table.

Syntax:

```
viewunselect <view_name>
```

Operands:

<view_name>

A valid view name existing in the current views table.

viewunselectrow

Unselect the view corresponding to the given row in the current views table.

Syntax:

```
viewunselectrow <row>
```

Operands:

<row>

A valid row in the views table.

viewunsetplay

Unset the play state of the given view.

Syntax:

```
viewunsetplay <view_name>
```

Operands:

<view_name>

A valid view name existing in the current views table.

viewvolume

Sets workspace view bounding box to the supplied values

Syntax:

```
viewvolume bottom=⟨x⟩ far=⟨x⟩ left=⟨x⟩ near=⟨x⟩ right=⟨x⟩
               top=⟨x⟩
```

Options:

<i>bottom</i>	Bottom co-ordinate of the bounding box Valid values: reals Default value: -5
<i>far</i>	Far co-ordinate of the bounding box Valid values: reals Default value: -10
<i>left</i>	Left co-ordinate of the bounding box Valid values: reals Default value: -5
<i>near</i>	Near co-ordinate of the bounding box Valid values: reals Default value: 10
<i>right</i>	Right co-ordinate of the bounding box Valid values: reals Default value: 5
<i>top</i>	Top co-ordinate of the bounding box Valid values: reals Default value: 5

visexport

Export a volume to the given file.

Syntax:

```
visexport entry=⟨text⟩ format=cnsunformatted | cnsformatted |
           cnsccp4 surface=⟨text⟩ surfacetype=⟨text⟩
           surfacetypecomment=⟨text⟩ volume=⟨text⟩ ⟨file name⟩
```

Options:

<i>entry</i>	The entry name of the entry that the surface belongs to. Valid values: text strings Default value:
--------------	---

Chapter 5: Commands

<i>format</i>	The format of the file that will be exported. Valid values: cnsunformatted cnsformatted cnscpp4 Default value: cnsformatted
<i>surface</i>	The name of the surface. Valid values: text strings Default value:
<i>surfacetype</i>	The surface type of the surface that will be exported. Valid values: text strings Default value:
<i>surfacetypecomment</i>	The surface type comment of the surface that will be exported. Valid values: text strings Default value:
<i>volume</i>	The name of the volume that the surface belongs to. Valid values: text strings Default value:

Operands:

⟨ file name ⟩

The name of the visio file to export.

visimport

Creates all volumes and surfaces in the given file.

Syntax:

```
visimport entry=⟨text⟩ isovalue=⟨x⟩ transparency=⟨x⟩  
          ⟨file.vis⟩ [ :⟨surface name⟩ ]
```

Options:

<i>entry</i>	This is the entry to associate the objects with. Valid values: text strings Default value:
--------------	---

isovalue If any of the volumes do not have associated surfaces or suggested isovalues, then this isovalue will be used to isosurface the volumes.

Valid values: reals

Default value: **0.03**

transparency

This option sets the transparency of any surfaces.

Valid values: reals

Default value: **0**

Minimum: 0.0

Maximum: 100.0

Operands:

$\langle \text{file.vis} \rangle$ [$\langle \text{surface name} \rangle$]

The first operand is the name of the visio or Jaguar plot file to load. The optional second operand is the specified surface name.

volumedelete

Deletes the given volume.

Syntax:

volumedelete *entry*= $\langle \text{text} \rangle$ $\langle \text{volume} \rangle$

Options:

entry The entry that the volume belongs to.

Valid values: text strings

Default value:

Operands:

$\langle \text{volume} \rangle$

The name of the volume to delete.

volumedisplay

Displays or undisplay the given volume

Syntax:

```
volumedisplay display=yes | no entry=⟨text⟩ ⟨volume⟩
```

Options:

<i>display</i>	Sets whether or not to display the volume. Valid values: boolean (true false; yes no; y n; on off) Default value: true
<i>entry</i>	The entry that the volume belongs to. Valid values: text strings Default value:

Operands:

⟨volume⟩

The name of the volume to display.

volumeisosurface

Isosurfaces the given volume.

Syntax:

```
volumeisosurface entry=⟨text⟩ isovalue=⟨x⟩ surface=⟨text⟩  
                  ⟨volume⟩
```

Options:

<i>entry</i>	The entry that the volume belongs to. Valid values: text strings Default value:
<i>isovalue</i>	The isovalue to use for the isosurfacing. Valid values: reals Default value: 0.1
<i>surface</i>	The name for the new surface. Valid values: text strings Default value:

Operands:

⟨volume⟩

The name of the volume to isosurface.

volumerename

Renames the given volume.

Syntax:

volumerename *entry*=⟨text⟩ *newname*=⟨text⟩ ⟨volume⟩

Options:

entry The entry that the volume belongs to.

Valid values: text strings

Default value:

newname The new name for the volume.

Valid values: text strings

Default value:

Operands:

⟨volume⟩

The name of the volume to rename.

volumesetcolors

Changes color map for the specified volume. The color points has to be specified after the volume name as five floating point numbers. The numbers correspond to data value, alpha value (opacity), R, G, B color components.

Syntax:

volumesetcolors *entry*=⟨text⟩ ⟨volume⟩

Options:

entry The entry that the volume belongs to.

Valid values: text strings

Default value:

Operands:

⟨ volume ⟩

The name of the volume to change colors.

wake

This is a standard alias for **energywake** (see [\[energywake\]](#), page 137).

workspacecopy

Currently this is only used by Edit->Copy to determine if it generates in the current working directory a number of files corresponding the data it puts on the clipboard. It does not actually cause the copy to execute. This has to be done graphically.

Syntax:

workspacecopy *writefiles*=yes | no

Options:

writefiles If true, write out files corresponding to the clipboard data. If false, do not write out the files.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

workspacefontdecrease

Decrease the size of the Workspace fonts

Syntax:

workspacefontdecrease

workspacefontincrease

Increase the size of the Workspace fonts

Syntax:

`workspacefontincrease`

workspacepaste

If there is a mime type of application/x-qt-windows-mime;value=SMILES on the clipboard, a plain text SMILES string in the global mouse selection, a plain text SMILES string on the clipboard, or an ISIS/Draw MDLCT, then use that as the basis for pasting into center of the Workspace a 2D structure. The structure is created as a project entry and included into the contents of the Workspace. If atoms were cut or copied within Maestro, then the pasted atoms are placed in the Workspace but not added as an entry to the Project Table.

Syntax:

`workspacepaste`

workspacepastesettings

Settings that affect workspacepaste

Syntax:

`workspacepastesettings enableselect=yes | no`

Options:

enableselect

Enable this to allow pasting of selected text on Linux. This only applies to systems running X Windows (such as Linux). If true, allows user to do text selection in windows (such as an xterm) and does not require copying the text to the clipboard (via a text editor's Copy, for example). There have been cases where text selection in python-based tkinter panels can be extremely slow when this feature is enabled. As of Suite 2011 most, if not all, GUI panels done in python use PyQt. So this should not be an issue.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

workspacescopeoff

Turn off the Workspace scope

Syntax:

workspacescopeoff

workspacescopeon

Turn on the Workspace scope

Syntax:

workspacescopeon

workspacescopesettings

Settings for the Workspace scope

Syntax:

workspacescopesettings *alpha*=⟨x⟩ *blendcolor*=background |
user *blendfactor*=⟨x⟩ *blue*=⟨x⟩ *green*=⟨x⟩ *radius*=⟨x⟩
red=⟨x⟩ *track*=yes | no *x*=⟨x⟩ *y*=⟨x⟩

Options:

alpha Alpha component. 0 is invisible. 1 is opaque.

Valid values: reals

Default value: **1**

Minimum: 0.0

Maximum: 1.0

blendcolor For pixels outside of the scope the color source to use.

Valid values: background

user

Default value: **background**

blendfactor

How much the pixels outside of the scope should be blended with the background color. 1.0 means 100%.

Valid values: reals

	Default value: 0.75
	Minimum: 0.0
	Maximum: 1.0
<i>blue</i>	Blue component when blendcolor is 'user'.
	Valid values: reals
	Default value: 0.5
	Minimum: 0.0
	Maximum: 1.0
<i>green</i>	Green component when blendcolor is 'user'.
	Valid values: reals
	Default value: 0.5
	Minimum: 0.0
	Maximum: 1.0
<i>radius</i>	Radius in pixels
	Valid values: reals
	Default value: 100
	Minimum: 0.0
	Maximum: 1000.0
<i>red</i>	Red component when blendcolor is 'user'.
	Valid values: reals
	Default value: 0.5
	Minimum: 0.0
	Maximum: 1.0
<i>track</i>	Enable or disable tracking of the cursor. This is only use if the workspace scope is on.
	Valid values: boolean (true false; yes no; y n; on off)
	Default value: true
<i>x</i>	X position
	Valid values: reals
	Default value: 0
<i>y</i>	Y position
	Valid values: reals
	Default value: 0

workspacescopeview

Determine which things are visible inside and outside of the scope. Only used when the scope is enabled.

Syntax:

```
workspacescopeview molecules=nowhere | inside | outside |  
everywhere ribbons=nowhere | inside | outside | everywhere  
surfaces=nowhere | inside | outside | everywhere  
volumes=nowhere | inside | outside | everywhere
```

Options:

<i>molecules</i>	<p>View molecules nowhere, inside the scope, outside the scope or everywhere (inside and outside).</p> <p>Valid values: nowhere inside outside everywhere</p> <p>Default value: everywhere</p>
<i>ribbons</i>	<p>View ribbons nowhere, inside the scope, outside the scope or everywhere (inside and outside).</p> <p>Valid values: nowhere inside outside everywhere</p> <p>Default value: everywhere</p>
<i>surfaces</i>	<p>View surfaces nowhere, inside the scope, outside the scope or everywhere (inside and outside).</p> <p>Valid values: nowhere inside outside everywhere</p> <p>Default value: everywhere</p>
<i>volumes</i>	<p>View volumes nowhere, inside the scope, outside the scope or everywhere (inside and outside).</p> <p>Valid values: nowhere inside outside everywhere</p> <p>Default value: everywhere</p>

workspaceselectionadd

Extends workspace selection with the atoms confirming to the supplied ASL.

Syntax:

workspaceselectionadd $\langle \text{ASL} \rangle$

Operands:

$\langle \text{ASL} \rangle$

ASL representing the set of atoms

workspaceselectionclear

Resets current workspace selection to none.

Syntax:

workspaceselectionclear

workspaceselectioninvert

Flips the state of an atom in the workspace selection. i.e., removes if the atom already exists or adds if it doesn't.

Syntax:

workspaceselectioninvert $\langle \text{ASL} \rangle$

Operands:

$\langle \text{ASL} \rangle$

ASL representing the set of atoms

workspaceselectionreplace

Overwrites the workspace selection with the supplied set of atoms.

Syntax:

workspaceselectionreplace $\langle \text{ASL} \rangle$

Operands:

$\langle \text{ASL} \rangle$

ASL representing the set of atoms

workspaceselectionsubtract

Removes atoms confirming to the supplied ASL from the workspace selection.

Syntax:

workspaceselectionsubtract \langle ASL \rangle

Operands:

\langle ASL \rangle

ASL representing the set of atoms

workspacetile

The workspacetile command specifies contents for a single tile in workspace tile mode.

Syntax:

workspacetile \langle tile_name \rangle

Operands:

\langle tile_name \rangle

The name which identifies the contents of the tile. This name will be the index of an entry group (greater than 0), or 0: followed by the entry ID for an ungrouped entry (or Scratch entry). A tile named 0: (with no entry ID) does not belong to any entry or entry group. In tile by surface mode, the entry ID will be followed by the name of the special surface for that tile.

write

This is a standard alias for **filewrite** (see [\[filewrite\]](#), page 197).

writeangle

Write the angles in the angle table to a file.

Syntax:

```
writeangle atomnumber=entry | molecule | workspace
           delimiter=comma | tab | userdefined userdelimiter=⟨text⟩
           ⟨file_name⟩
```

Options:

atomnumber

This option is deprecated.

Valid values: entry
 molecule
 workspace

Default value: **entry**

delimiter This option sets what delimiter will be used in the output file:
 comma=1, tab=2, and userdefined=3.

Valid values: comma
 tab
 userdefined

Default value: **comma**

userdelimiter

The delimiter string defined by users. It is valid only when option *delimiter* is userdefined.

Valid values: text strings

Default value:

Operands:

⟨file_name⟩

The name of the file to which the angles will be written.

writecontact

Write Contacts information to a file.

Syntax:

```
writecontact atomnumber=entry | molecule | workspace  
             delimiter=comma | tab | userdefined userdelimiter=⟨text⟩  
             ⟨file_name⟩
```

Options:

atomnumber

This option is deprecated.

Valid values: entry
 molecule
 workspace

Default value: **entry**

delimiter This option sets what delimiter will be used in the output file:
comma=1, tab=2, and userdefined=3.

Valid values: comma
 tab
 userdefined

Default value: **comma**

userdelimiter

The delimiter string defined by users. It is valid only when
option *delimiter* is userdefined.

Valid values: text strings
Default value:

Operands:

⟨file_name⟩

The name of the file to which Contacts information will be written.

writecoupling

Write the coupling constants in the coupling table to a file.

Syntax:

```
writecoupling atomnumber=entry | molecule | workspace  
             delimiter=comma | tab | userdefined userdelimiter=⟨text⟩  
             ⟨file_name⟩
```

Options:

atomnumber

This option is deprecated.

	Valid values:	entry molecule workspace
	Default value:	entry
<i>delimiter</i>	This option sets what delimiter will be used in the output file: comma=1, tab=2, and userdefined=3.	
	Valid values:	comma tab userdefined
	Default value:	comma
<i>userdelimiter</i>	The delimiter string defined by users. It is valid only when option delimiter is userdefined.	
	Valid values:	text strings
	Default value:	

Operands:

⟨file_name⟩

The name of the file to which the coupling constants will be written.

writedihtedral

Write the dihedrals in the dihedral table to a file.

Syntax:

```
writedihtedral atomnumber=entry | molecule | workspace
               delimiter=comma | tab | userdefined userdelimiter=⟨text⟩
               ⟨file_name⟩
```

Options:

atomnumber

This option has been deprecated.

Valid values: entry
molecule
workspace

Default value: **entry**

delimiter This option sets what delimiter will be used in the output file:
comma=1, tab=2, and userdefined=3.

Valid values: comma
 tab
 userdefined
Default value: **comma**

userdelimiter

The delimiter string defined by users. It is valid only when option *delimiter* is *userdefined*.

Valid values: text strings
Default value:

Operands:

$\langle \text{file_name} \rangle$

The name of the file to which the dihedrals will be written.

writedistance

Write the distances in the distance table to a file.

Syntax:

writedistance *atomnumber*=entry | molecule | workspace
 delimiter=comma | tab | userdefined *userdelimiter*= $\langle \text{text} \rangle$
 $\langle \text{file_name} \rangle$

Options:

atomnumber

This option is deprecated.

Valid values: entry
 molecule
 workspace

Default value: **entry**

delimiter

This option sets what delimiter will be used in the output file: comma=1, tab=2, and userdefined=3.

Valid values: comma
 tab
 userdefined

Default value: **comma**

userdelimiter

The delimiter string defined by users. It is valid only when option *delimiter* is *userdefined*.

Valid values: text strings
 Default value:

Operands:

`< file_name >`

The name of the file to which the distances will be written.

writehbond

Write H-Bonds information to a file.

Syntax:

```
writehbond atomnumber=entry | molecule | workspace
           delimiter=comma | tab | userdefined userdelimiter=< text >
           < file_name >
```

Options:

atomnumber

This option is deprecated.

Valid values: entry
 molecule
 workspace

Default value: **entry**

delimiter This option sets what delimiter will be used in the output file:
 comma=1, tab=2, and userdefined=3.

Valid values: comma
 tab
 userdefined

Default value: **comma**

userdelimiter

The delimiter string defined by users. It is valid only when option *delimiter* is userdefined.

Valid values: text strings
 Default value:

Operands:

`< file_name >`

The name of the file to which H-Bonds information will be written.

wsassistantapply

Command that applies workspace style settings.

Syntax:

```
wsassistantapply
```

wsassistantapplycontacthbond

Command that applies contacts and hbonds options of workspace assistant.

Syntax:

```
wsassistantapplycontacthbond
```

wsassistantapplyfocus

Command that applies workspace style settings and sets focus on ligand.

Syntax:

```
wsassistantapplyfocus
```

wsassistantbindingsurface

Command that handles workspace assistant binding site surface options.

Syntax:

```
wsassistantbindingsurface colorscheme=constant |  
    partialcharge | atomtype | chainname | element | molnum |  
    molnumcarbon | residuecharge | residueproperty |  
    residueposition | residuetype | entry | atomcolor | elepotential  
    | pbelepotential createligandsurfaces=yes | no  
    createreceptorsurfaces=yes | no radius=<x> spacing=<x>  
    style=solid | mesh | dot transparency=<n>
```

Options:

colorscheme

What kind of color scheme for the binding surface.

Valid values: constant
 partialcharge
 atomtype
 chainname
 element
 molnum
 molnumcarbon
 residuecharge
 residueproperty
 residueposition
 residuetype
 entry
 atomcolor
 elepotential
 pbelepotential

Default value: **constant**

createligandsurfaces

Create ligand surfaces or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

createreceptorsurfaces

Create receptor surfaces or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

radius

The radius used to truncate receptor surface.

Valid values: reals
 Default value: **5**
 Minimum: 0.01
 Maximum: 20.0

spacing

The surface grid spacing value.

Valid values: reals
 Default value: **0.6**
 Minimum: 0.01
 Maximum: 10.0

style

What kind of style for the binding surface: solid (0), Mesh (1), or Dot (2).

Valid values: solid
 mesh
 dot
 Default value: **solid**

transparency

The percentage of surface transparency.

Valid values: integers

Default value: **70**

Minimum: 0

Maximum: 100

wsassistantcontacthbond

Command that handles workspace assistant contacts and hbonds options.

Syntax:

```
wsassistantcontacthbond contacttype=badandugly | ugly | bad |  
good intraligandcontacts=yes | no  
receptorligandcontacts=yes | no receptorligandhbonds=yes | no
```

Options:

contacttype

What kind of contact type the workspace assistant contains.

Valid values: badandugly
ugly
bad
good

Default value: **badandugly**

intraligandcontacts

Show intra-ligand contacts or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

receptorligandcontacts

Show receptor-ligand contacts or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

receptorligandhbonds

Show reseptor-ligand hbonds or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

wsassistantcreatesurfaces

Command that creates binding site surfaces.

Syntax:

```
wsassistantcreatesurfaces
```

wsassistantdefaultstyle

Command that sets default workspace style.

Syntax:

```
wsassistantdefaultstyle
```

wsassistantfixreceptor

Command to fix the receptor in the Workspace.

Syntax:

```
wsassistantfixreceptor
```

wsassistanthideallsurfaces

Command that Hides all surfaces.

Syntax:

```
wsassistanthideallsurfaces
```

wsassistantliganddetection

Command that handles workspace assistant ligand detection options.

Syntax:

```
wsassistantliganddetection allowaminoacids=yes | no  
                        allowsmallions=yes | no includedresiduecode= $\langle$ text $\rangle$   
                        includedresiduecodes= $\langle$ text $\rangle$  maximumatomcount= $\langle$ n $\rangle$   
                        minimumatomcount= $\langle$ n $\rangle$  residuecode= $\langle$ text $\rangle$   
                        residuecodes= $\langle$ text $\rangle$ 
```

Options:

allowaminoacids

Allow molecules containing only amino acids or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

allowsmallions

Allow molecules that are small ions or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

includedresiduecode

The 3-letter residue code for included residues.

Valid values: text strings

Default value:

includedresiduecodes

All 3-letter residue codes for included residues.

Valid values: text strings

Default value:

maximumatomcount

The maximum atom count of ligand molecule.

Valid values: integers

Default value: **130**

Minimum: 2

minimumatomcount

The minimum atom count of ligand molecule.

Valid values: integers

Default value: **5**

Minimum: 1

residuecode

The 3-letter residue code for excluded residues.

Valid values: text strings

Default value:

residuecodes

All 3-letter residue codes for excluded residues.

Valid values: text strings
 Default value:

wsassistantmaterial

Command to set the the appearance to one of a set of preset definitions. These control the lighting and material properties used in the Workspace.

Syntax:

wsassistantmaterial

wsassistantprimarystyle

Command that handles workspace assistant primary style options.

Syntax:

wsassistantprimarystyle *colorscheme*=element | elementgreen | elementlight | elemententry | bfactor | densityfit | moleculesize | residuetype | entry *distalribbonscheme*=residueposition | secondarystructure | chain | residuetype | residueproperty | entry | calphaatom *distalribbonstyle*=none | cartoon | ribbon | tube | thintube | curvedline | calphaline | calphatube *fit*=yes | no *pickstate*=atom | residue | molecule | chain | entry *proximalribbonscheme*=residueposition | secondarystructure | chain | residuetype | residueproperty | entry | calphaatom *proximalribbonstyle*=none | cartoon | ribbon | tube | thintube | curvedline | calphaline | calphatube *representation*=default | wire | cpk | ballbstick | tube

Options:

colorscheme

What kind of color scheme the primary style contains.

Chapter 5: Commands

Valid values: element
elementgreen
elementlight
elemententry
bfactor
densityfit
moleculesize
residuetype
entry
Default value: **element**

distalribbonscheme

What kind of distal ribbon scheme the primary style contains.

Valid values: residueposition
secondarystructure
chain
residuetype
residueproperty
entry
calphaatom
Default value: **residueposition**

distalribbonstyle

What kind of distal ribbon style the primary style contains.

Valid values: none
cartoon
ribbon
tube
thintube
curvedline
calphaline
calphatube
Default value: **calphatube**

fit

Fit Workspace to picked entity or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **true**

pickstate

What kind of pick state the primary style contains.

Valid values: atom
residue
molecule
chain
entry
Default value: **residue**

proximalribbonscheme

What kind of proximal ribbon scheme the primary style contains.

Valid values: residueposition
 secondarystructure
 chain
 residuetype
 residueproperty
 entry
 calphaatom

Default value: **residueposition**

proximalribbonstyle

What kind of proxiaml ribbon style the primary style contains.

Valid values: none
 cartoon
 ribbon
 tube
 thintube
 curvedline
 calphaline
 calphatube

Default value: **cartoon**

representation

What kind of representation the primary style contains.

Valid values: default
 wire
 cpk
 ballbstick
 tube

Default value: **default**

wsassistantrendering

Command to set the quality of workspace rendering. This allows a tradeoff between quality and speed.

Syntax:

wsassistantrendering

wsassistantshowallsurfaces

Command that shows all surfaces.

Syntax:

wsassistantwsstyle

Command that handles workspace assistant workspace style options.

wsassistantshowallsurfaces

Syntax:


```

wsassistantwsstyle applycolorschemes=yes | no
    applyrepresentations=yes | no applysettings=yes | no
    atomradius=<x> ballpercentage=<n>
    createligandsurfaces=yes | no createreceptorsurfaces=yes | no
    defaultrepresentation=wire | cpk | ballnstick | tube |
thin_tube defaultscheme=element | elementgreen |
elementlight | elemententry | bfactor | densityfit | moleculesize
| residuetype | entry displayperspective=yes | no
displayribbons=yes | no displaywithin=yes | no
fittoligand=yes | no focusligands=yes | no
ionrepresentation=wire | cpk | ballnstick | tube | thin_tube
ionscheme=element | elementgreen | elementlight |
elemententry | bfactor | densityfit | moleculesize | residuetype
| entry ligandhydrogen=none | polaronly | all
ligandrepresentation=wire | cpk | ballnstick | tube | thin_tube
ligandscheme=element | elementgreen | elementlight |
elemententry | bfactor | densityfit | moleculesize | residuetype
| entry materialeffect=default | aluminum | sundown |
dramatic | gloss multibondradius=<x> reapply=yes | no
reapplyfocus=yes | no reapplyligand=yes | no
reapplyto=workspace | ligand receptorhydrogen=none |
polaronly | all ribbonlimit=yes | no ribbonradius=<x>
ribbonscheme=residueposition | secondarystructure | chain |
residuetype | residueproperty | entry | calphaatom
ribbonstyle=none | cartoon | ribbon | tube | thintube |
curvedline | calphaline | calphatube singlebondradius=<x>
stickradius=<x> suppressdisplay=yes | no waterdisplay=none
| all | beyond3a | beyond5a | beyond10a
waterrepresentation=wire | cpk | ballnstick | tube | thin_tube
waterscheme=element | elementgreen | elementlight |
elemententry | bfactor | densityfit | moleculesize | residuetype
| entry within=<x>

```

Options:

applycolorschemes

Apply color schemes or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

applyrepresentations

Apply representations or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

Chapter 5: Commands

applysettings

When Workspace changes applying Workspace style settings or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

atomradius

The radius of atom in representation.

Valid values: reals

Default value: **5**

Minimum: 0.0

Maximum: 1000.0

ballpercentage

The percentage of ball radius in representation.

Valid values: integers

Default value: **16**

Minimum: 4

Maximum: 200

createligandsurfaces

Create ligand surface or not when reapply.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

createreceptorsurfaces

Create receptor surface or not when reapply.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

defaultrepresentation

What kind of default representation the workspace style contains.

Valid values: wire
 cpk
 ballnstick
 tube
 thin_tube

Default value: **wire**

defaultscheme

What kind of default color scheme the workspace style contains.

Valid values: element
 elementgreen
 elementlight
 elemententry
 bfactor
 densityfit
 moleculesize
 residuetype
 entry
 Default value: **element**

displayperspective

Display in perspective or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

displayribbons

Display ribbons or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

displaywithin

Display residues within a distance or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

fittoligand When Workspace changes Fitting to ligand or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

focusligands

Focus display on ligands or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **false**

ionrepresentation

What kind of ion representation the workspace style contains.

Valid values: wire
 cpk
 ballnstick
 tube
 thin_tube
 Default value: **cpk**

ionscheme What kind of ion color scheme the workspace style contains.

Chapter 5: Commands

Valid values: element
elementgreen
elementlight
elemententry
bfactor
densityfit
moleculesize
residuetype
entry
Default value: **element**

ligandhydrogen

Options for displaying hydrogens of ligand: None (0), Polar only (1), or All (2).

Valid values: none
polaronly
all
Default value: **all**

ligandrepresentation

What kind of ligand representation the workspace style contains.

Valid values: wire
cpk
ballnstick
tube
thin_tube
Default value: **ballnstick**

ligandscheme

What kind of ligand color scheme the workspace style contains.

Valid values: element
elementgreen
elementlight
elemententry
bfactor
densityfit
moleculesize
residuetype
entry
Default value: **elemententry**

materialeffect

What kind of lighting/material effect the workspace style contains.

Valid values: default
aluminum
sundown
dramatic
gloss
Default value: **default**

multibondradius

The radius of multiple tube bond in representation.

Valid values: reals
Default value: **0.05**
Minimum: 0.01
Maximum: 1.0

reapply

Reapply workspace style settings for receptor or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

reapplyfocus

Reapply workspace style focus or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

reapplyligand

Reapply workspace style settings for ligand or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

reapplyto

Reapply workspace style to which structures in the Workspace, entire Workspace (0) or ligand only (1).

Valid values: workspace
ligand
Default value: **workspace**

receptorhydrogen

Options for displaying hydrogens of receptor: None (0), Polar only (1), or All (2).

Valid values: none
polaronly
all
Default value: **polaronly**

ribbonlimit

Display ribbons within a range or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
Default value: **false**

Chapter 5: Commands

ribbonradius

The additional radius of ribbon in representation.

Valid values: reals
Default value: **5**
Minimum: 0.0
Maximum: 1000.0

ribbonscheme

What kind of ribbon scheme the workspace style contains.

Valid values: residueposition
 secondarystructure
 chain
 residuetype
 residueproperty
 entry
 calphaatom
Default value: **residueposition**

ribbonstyle

What kind of ribbon style the primary style contains.

Valid values: none
 cartoon
 ribbon
 tube
 thintube
 curvedline
 calphaline
 calphatube
Default value: **cartoon**

singlebondradius

The radius of single tube bond in representation.

Valid values: reals
Default value: **0.1**
Minimum: 0.01
Maximum: 1.0

stickradius

The radius of stick bond in representation.

Valid values: reals
Default value: **0.12**
Minimum: 0.01
Maximum: 1.0

suppressdisplay

Suppress nonpolar-hydrogen display or not.

Valid values: boolean (true|false; yes|no; y|n; on|off)
 Default value: **true**

waterdisplay

What kind of water display suppression the workspace style contains.

Valid values: none
 all
 beyond3a
 beyond5a
 beyond10a
 Default value: **none**

waterrepresentation

What kind of water representation the workspace style contains.

Valid values: wire
 cpk
 ballnstick
 tube
 thin_tube
 Default value: **ballnstick**

waterscheme

What kind of water color scheme the workspace style contains.

Valid values: element
 elementgreen
 elementlight
 elemententry
 bfactor
 densityfit
 moleculesize
 residuetype
 entry
 Default value: **element**

within

The radius within which structures will be displayed.

Valid values: reals
 Default value: **8**
 Minimum: 0.0
 Maximum: 999.0

xcluster

This keyword is used to set various options associated with starting XCluster jobs from Maestro.

Syntax:

```
xcluster cluster_by=atom | torsion compare_enantiomers=yes | no  
          input_file=⟨text⟩ job=⟨text⟩ rms_in_place=yes | no  
          structure_source=selected_entries | workspace | file
```

Options:

cluster_by Whether to use atomic RMS or torsional RMS differences for clustering analysis.

Valid values: atom
torsion

Default value: **atom**

compare_enantiomers

A boolean which controls whether XCluster will perform superposition in such a way so as to compare enantiomers.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **true**

input_file The name of the structure input file.

Valid values: text strings

Default value:

job The name for the XCluster job.

Valid values: text strings

Default value: **xclustmp**

rms_in_place

A boolean which controls whether XCluster will calculate RMS differences without first doing a superposition.

Valid values: boolean (true|false; yes|no; y|n; on|off)

Default value: **false**

structure_source

Whether to use the selected entries in the current project, or a specified file with multiple structures as structure input for the job.

Valid values: selected_entries
workspace
file

Default value: **selected_entries**

xclusterstart

Start a XCluster job with the current settings.

Syntax:

```
xclusterstart
```

xclusterwrite

Write a XCluster input file with the current settings.

Syntax:

```
xclusterwrite
```

zoom

Zoom in or out (i.e. scale up or down)

Syntax:

```
zoom factor=⟨x⟩ in | out
```

Options:

factor A zooming (scaling) factor which determines by how much the on-screen structures is scaled. 1.0 is 1 percent, 99.0 is 99 percent, etc. Values are used relative to current size

Valid values: reals

Default value: **10**

Minimum: 1.0

Maximum: 99.0

Operands:

in | out

Must specify whether to zoom in and make larger or out and make smaller.

Command Index

1

1ddataset.....	29
1dplot.....	31
1drescale.....	33
1dtable.....	33

2

2ddataset.....	33
2dplot.....	35
2drescale.....	36

A

addatomproperty.....	36
addfromhold.....	37
adjustangle.....	37
adjustchirality.....	38
adjustdihedral.....	38
adjustdistance.....	39
alias.....	40
alignaxis.....	40
alignplane.....	41
angle.....	42
appendribbons.....	42
application.....	43
assigncharges.....	43
atom.....	43
atomname.....	44
atomproplabel.....	44
attach.....	45
attachmentmarkerdump.....	45
attachmentmarkersettings.....	45
autosetup.....	49

B

beginundoblock.....	49
bmincomfile.....	50
bond.....	50
bondorder.....	50
bondoutline.....	51
bondtonew.....	53
buildoptions.....	53

C

calcenergy.....	54
canonicalname.....	54
caption.....	55
captionbackground.....	55
captioncolor.....	55
captiondelete.....	56
captionfont.....	56
captionhide.....	57
captionposition.....	57
captionselect.....	58
captionselectadd.....	58
captionshow.....	58
captionsoff.....	58
captionson.....	58
captiontext.....	59
captionunselect.....	59
captionunselectall.....	59
cascadepanels.....	59
cellsmarkerdump.....	60
cellsmarkersettings.....	60
centeratom.....	63
centerbond.....	64
centercoordinates.....	64
centroid.....	65
centroidatom.....	65
centroidatomset.....	65
cglidedockconstraintposition.....	66
cglidedockconstraintregion.....	67
chainname.....	68
changedirectory.....	69
chiralatom.....	69
clip.....	70
clipatomset.....	71
clipsurfaceset.....	71
cliptoallatoms.....	72
cliptoligand.....	72
cliptoselectedatoms.....	72
clusteratom.....	72
clusterheavyatoms.....	73
clusterset.....	73
clustertorsion.....	73
color.....	74
coloratom.....	75
coloratombonds.....	75
coloratomrgb.....	76
colorbond.....	76

Command Index

colorscheme	77
combilibenum	78
combilibenumaddattachment	78
combilibenumclearreagentfile	79
combilibenumdeleteattachment	79
combilibenumexportdefinition	79
combilibenumimportdefinition	79
combilibenumoptions	80
combilibenumrefreshstructure	80
combilibenumrenameattachment	80
combilibenumselectextendtablerow ..	81
combilibenumselectonlytablerow	81
combilibenumselecttablerow	82
combilibenumsetmolecule	82
combilibenumsetreagentfile	82
combilibenumunselecttablerow	83
compareatom	83
compareset	83
confelim	84
confelimstart	85
confelimwrite	85
confgenltsearch	85
confgenmini	88
confgenpotential	89
confgenreadpotential	91
confgenstart	91
confgenstd	91
confgenstdsettings	92
confgenstdwrite	93
confgenwrite	93
confsearch	94
connect	97
connectfuseatom	97
constrainedangle	97
constrainedatom	98
constraineddist	99
constrainedset	100
constrainedtorsion	101
contactcriteria	102
contactset1	103
contactset2	103
coordinatescan	104
coupling	105
createcustompropertygroup	106
createlibrary	106
createpropsubset	107
createsubset	107

D

defaultfc	107
-----------------	-----

delete	108
deletcustompropertygroup	108
deleteemptyentrygroups	109
deleteslectedproperties	109
deletezeroorderbonds	109
densityblobupdate	110
dialogoptions	110
dihedral	111
dihedraldrive	111
dipolemoment	112
displayatom	118
displayonlyatom	118
displayopt	118
distance	123
distancecheck	124
dockgridbox	125
doplanealignment	126
dynamics	126

E

editatomproperty	128
ejob	129
elementlabels	130
embrace	130
embracecsearch	132
emptyviewreset	133
endundoblock	134
energykill	134
energymonitor	134
energysleep	135
energystart	135
energystop	136
energytask	136
energyupdate	136
energyupdatejobstatus	137
energywake	137
enhance3d	138
enterposeviewermode	138
entryadddisplay	138
entrycolorbyproperty	139
entrycombine	140
entrycopyprop	141
entrycreatesmoothed	141
entrydelete	142
entrydisassociate	142
entrydisassociatebychain	143
entrydisassociatebyligand	143
entrydisplaysurfaces	144
entrydragselection	144
entryduplicate	144

entryexport	145	entryshowall	177
entryexportcanvas	149	entrystars	177
entryexportspreadsheet	149	entrytable	178
entryextendsselection	150	entryundisplaysurfaces	178
entryextendselectrow	151	entryunselect	178
entryextendwsinclude	151	entryunselectall	179
entrygroupbyproperty	151	entryunselectrow	179
entrygroupcollapse	152	entrywatermapexamineresults	179
entrygroupcountcontacts	153	entrywscreate	180
entrygroupcounthbonds	153	entrywsexclude	180
entrygroupcreate	154	entrywsexcludefixed	180
entrygroupdelete	154	entrywsinclude	181
entrygroupdragselection	155	entrywsincludelock	181
entrygroupduplicate	155	entrywsincludeonly	181
entrygroupexpand	156	entrywsincludeunlock	182
entrygroupexpandonly	156	epik	182
entrygroupextendsselection	157	epikreadinput	184
entrygroupinvertselection	157	epikstart	184
entrygroupmove	157	epikwriteinput	185
entrygrouprename	158	eplayergoto	185
entrygroupselect	158	eplayergotofirst	185
entrygroupselectonly	159	eplayergotolast	186
entrygroupsettitle	159	eplayernextgroup	186
entrygroupungroup	159	eplayerpreviousgroup	186
entrygroupunselect	160	eplayersettings	187
entrygroupviewposes	160	eplayerstepahead	189
entrygroupwsexclude	161	eplayerstepback	189
entrygroupwsinclude	162	errorcheck	189
entrygroupwsincludeonly	162	excludedvolumesmarkerssettings	190
entryimport	163	exportrotatablegroups	191
entryimportspreloadsheet	167	extendsselecteditatom	192
entryimportvibration	168		
entryinvertselection	168	F	
entrymergeprop	169	fieldqsarmarkersettings	192
entrymovetogroup	170	fileread	195
entryremovedisplay	170	filewrite	198
entryrename	170	filter	200
entryresetcolor	171	find	201
entryselect	171	findnext	204
entryselectall	172	findprevious	204
entryselectonly	172	fit	205
entryselectonlyrow	172	fitbox	205
entryselectrandom	172	fitligand	205
entryselectrow	173	fitnextligand	205
entrysetcolor	174	fitoptions	206
entrysetdeletable	174	fitpreviousligand	206
entrysetnondeletable	175	flexibleligandalignment	207
entrysetprop	175	forcefield	207
entrysetreadonly	176	forcefieldbend	207
entrysettitle	176	forcefielddele	208
entrysetwritable	177		

Command Index

forcefieldgbolv 209
forcefieldimproper 209
forcefieldsasolv 210
forcefieldstretch 211
forcefieldtorsion 211
forcefieldvdw 212
forcefieldview 213
forcefieldwilson 213
formalcharge 214
fragment 214
fragmentfuse 215
fragmenttype 215
frozenatom 215
frozenet 216
fuse 216

G

geometrycleanup 216
glideactivesiteres 217
glideactivesiteset 217
glidecalcbboxfromligand 217
glideconstraintatomlabel 218
glideconstraintposition 218
glideconstraintregion 219
glideconstraintregionactivecell 219
glideconstraintregiongrow 219
glideconstraintregioninvisible 220
glideconstraintregionshrink 220
glideconstraintregionvisible 220
glidecoremolecule 221
glidecorermssubset 221
glidedisplayreceptor 221
glidedockconstraintposition 222
glidedockconstraintregion 223
glidedockexcludedvolume 224
glidedockingreset 225
glidedockregionnumatoms 225
glidedockregionselect 226
glidedockregionunselect 226
glideexcludedvolume 227
glidegridhydrophobic 227
glidegridreset 228
glideligand 228
glideligandfile 231
glideligandgrid 231
glidelocatehydrophobiccells 233
glideoutput 233
glidereceptorligand 235
glidereferenceligand 236
glideringsampling 237

gliderotatablegroups 237
gliderotatablegroupsextenduserow
..... 238
gliderotatablegroupstoggleatom... 238
gliderotatablegroupstogglerow... 238
gliderotatablegroupsuseexisting.. 239
glidescoring 239
glidesettings 241
glidesite 246
glidetorsionalconstraint 248
glidetorsionalconstraintallbonds
..... 249
glidetorsionalconstraintpatternadd
..... 249
glidetorsionalconstraintpatterndelete
..... 250
glidetorsionalconstraintpatterndeleteall
..... 250
glidetorsionalconstraintpatterndrag
..... 250
glidetorsionalconstraintpatternreplace
..... 250
glidetorsionalconstraintpatternselect
..... 251
glidetorsionalconstrainttorsionadd
..... 251
glidetorsionalconstrainttorsiondelete
..... 251
glidetorsionalconstrainttorsiondeleteall
..... 252
glidetorsionalconstrainttorsiondeleteinvalid
..... 252
glidetorsionalconstrainttorsionreplace
..... 252
glidetorsionalconstrainttorsionselect
..... 252
glidetorsionalconstrainttorsionsetangle
..... 253
glideundisplayreceptor 253
googleprofiling 253
graphicsdebugdisable 254
graphicsdebugenable 254
grow 255
growbond 255
growbond2 255
growdirection 256
growname 256

H

hbondcriteria 257

hbondset1	257
hbondset2	258
helpauto	258
helpcategory	259
helpsearch	259
helptopic	260
hidedockablepanels	260
hideentries	261
hidemarkers	261
hidepanel	261
hidepanels	261
hideproperty	262
hidetoolbar	262
highlight	262
highlightarrowatoms	263
highlightarrowcolor	263
highlightarrowsettings	264
highlightatoms	264
highlightbackground	264
highlightcolor	265
highlightdelete	265
highlightdirection	266
highlightduplicate	266
highlighthide	266
highlightmethod	267
highlightpreference	267
highlightramp	268
highlightrename	268
highlightselect	268
highlightshow	269
highlightsoff	269
highlightson	269
highlighttext	269
highlighttextcolor	269
highlighttextfont	270
highlighttextposition	271
highlighttextselect	272
highlighttextselectadd	272
highlighttextsettings	272
highlighttextunselect	272
highlighttextunselectall	273
highlightwidth	273
historyvisible	274
hold	280
hppmap	280
hppmapbox	281
hppmapset	282
hppmapstart	282
hppmapwrite	282
hydrogenapply	283
hydrogentreat	283

I

impactbufferedatom	284
impactbufferedset	284
impactconstraints	284
impactcontinuumsolvent	285
impactdynamics	287
impactfastmultipole	289
impactfrozenatom	289
impactfrozenset	290
impacthybridmc	290
impactjob	291
impactmdparams	292
impactminimization	294
impactperiodicboundary	296
impactpotential	297
impactread	299
impactrepexch	299
impactselectextendablerow	300
impactselectonlytablerow	301
impactstart	301
impacttask	302
impacttruncation	302
impactwrite	303
invert	303
invertset	303

J

jaguarassignatomnames	304
jaguarimportgeometry	304
jaguarinputfilesjob	305
jaguarjob	306
jaguarselectextendablerow	306
jaguarselectonlytablerow	307
jobcleanup	307
jobsettings	308

L

labelatom	311
labelatomfieldsnone	321
labelatomoffset	321
labelatomoffsetreset	321
labelclear	322
labeldisplayusertext	322
labelformat	322
labelupdate	323
labelusertextatomupdate	323
liaisonanalysis	324
liaisonparameters	326

Command Index

liaisonselectlig 329
liaisonsettings 329
liaisonssystem 330
ligandbond 331
light 332
lightambient 333
lightdiffuse 333
lightposition 333
lightspecular 334
ligprep 334
ligprepread 337
ligprepstart 337
ligprepwrite 337
localatoms 337
localbitset 338
localcenter 338
logp 338

M

macrodefine 339
macrorun 339
makedirectory 340
markers 340
material 341
materialdefault 342
materialmolecular 342
mcscd 343
minienergy 344
minimize 345
minta 345
monitorangle 346
monitordistance 346
monitorrhbond 347
monitorsetsurf 348
monitorsettings 348
monitorsurf 349
monitortorsion 349
mouse 350
move 350
multiplemini 351
mutate 352

N

nextresidue 352

O

optimizefog 353

P

partialcharge 353
pause 353
pausecommands 354
peptidecistrans 354
peptiderotate 355
perresiduedisplay 355
perresiduesubset 356
perresiduetype 356
phaseaddconstraint 356
phaseaddcustomfeature 357
phaseaddhypotheses 357
phaseaddligands 358
phaseaddligandsfromfile 358
phaseaddligandsfromrun 359
phaseaddsite 360
phasealignhypotheses 360
phasebuildqsar 361
phasechangesitetype 364
phasechooseactiveset 364
phasecleanupstructures 365
phaseclearcentroidatoms 365
phaseclusterhypotheses 365
phaseconfgen 366
phasecopyhypothesisfrombuilder 369
phasecreateevactives 369
phasecreateevactivessettings 369
phasecreateevreceptor 370
phasecreateevreceptoratoms 370
phasecreateevreceptorsettings 371
phasecreateevreference 372
phasecreateevreferenceatoms 372
phasecreateevreferencesettings 373
phasecreateexcludedvolume 373
phasecreatefreestylehypothesis 374
phasedbaddconfset 374
phasedbaddligands 374
phasedbaddligandsfromdb 375
phasedbconfgen 376
phasedbcreatesites 379
phasedbcreatesubsetfromhitfile 379
phasedbcreatesubsetfromselected 380
phasedbcreatesubsetfromtextfile 380
phasedbdeleteligands 380
phasedbdeletesubset 381
phasedbdisplayrange 381
phasedbexportligands 381
phasedbfilter 382
phasedbfilteradd 382
phasedbgenerateconformers 383
phasedbnewfromfile 383

phasedbopen	383	phasefieldqsarselectcustom	400
phasedbremoveconformers	384	phasefieldqsarselectligandrow....	401
phasedbselectextendtablerow	384	phasefieldqsarsetactivity	401
phasedbselectonlytablerow	384	phasefieldqsarsettrainingrows....	402
phasedbselectsuubsetrow	385	phasefieldqsarsort	402
phasedbselecttablerow	385	phasefieldqsartest	402
phasedbunselecttablerow	385	phasefieldqsartoggleinclude	403
phasedeletecustomfeature	386	phasefieldqsartoggleligandrow....	403
phasedeleteexcludedvolumes	386	phasefieldqsartogglerow	403
phasedeletehypothesis	386	phasefieldqsartoggleselectcustom	
phasedeletehypothesisfrombuilder		404
.....	387	phasefieldqsarvisselectcontour...	404
phasedeleteselectdconstraints...	387	phasefieldqsarvisselectintensity	
phasedeleteselectdligands	387	404
phasedeletesites	387	phasefieldqsarvistoggleselectcontour	
phasedisplayproperty	388	405
phaseeoptions	388	phasefindmatches	405
phaseexcludetablerow	388	phasefindpharmacophores	411
phaseexportalignmentstofile	389	phasegenerateconformers	411
phaseexportconformerstofile	389	phasegroupligands	411
phaseexportfeature	389	phasegroupligandsbytitle	412
phaseexporthypothesis	390	phasehookimport	412
phaseexporthypothesisfrombuilder		phasehypothesiscreateexcludedvolume	
.....	390	412
phaseexportselectedalignments....	390	phasehypothesisdeleteexcludedvolumes	
phaseexportselectedhypotheses....	391	412
phaseexportselectedhypothesesfrombuilder		phasehypothesisselectevrow	413
.....	391	phasehypothesisselectextendevrow	
phaseexportselectedligands	391	413
phaseextendsite	391	phasehypothesisselectonlyevrow...	413
phasefieldqsar	392	phasehypothesisselectrow	414
phasefieldqsaraddligands	395	phasehypothesissetexcludedvolumes	
phasefieldqsaraddligandsfromfile		414
.....	396	phasehypothesissetid	414
phasefieldqsaraddtoproject	396	phasehypothesissetsitemask	415
phasefieldqsarbuild	396	phasehypothesissettolerances	415
phasefieldqsardelete	396	phasehypothesistoggleinclude	417
phasefieldqsardeleteall	397	phasehypothesisunselectevrow....	418
phasefieldqsarexport	397	phaseimportfeature	418
phasefieldqsarexportligandstofile		phaseimporthypothesis	418
.....	397	phaseincludeextendtablerow	419
phasefieldqsarextendinclude	398	phaseincludeonlytablerow	419
phasefieldqsarextendselectligandrow		phaseinclugetablerow	419
.....	398	phaseinverttableselection	420
phasefieldqsarimport	398	phaseligprep	420
phasefieldqsarinclude	399	phaselocatesites	421
phasefieldqsarinvertselection....	399	phasemarkerdump	421
phasefieldqsarplotactivity	399	phasemarkersettings	422
phasefieldqsarpredict	400	phasemarkfeature	433
phasefieldqsarrandomtraining....	400	phasemergestereoisomers	433
phasefieldqsarselectall	400	phasemergetitles	433

Command Index

phaseoptions	434	phaserunsetseed	457
phasepatterndelete	436	phaserunsetstereo	457
phasepatternedit	436	phasesavematching	457
phasepatternmovedown	437	phasescorehypotheses	458
phasepatternmoveup	438	phasescoreinactives	462
phasepatternnew	438	phasescoreselecthypothesis	462
phasepatternsetoptions	439	phasesearchformatches	463
phaseplotactivity	440	phaseselectalltablerows	463
phasepreview	440	phaseselectconstraintrow	464
phaseqsar	441	phaseselectevrow	464
phaseqsaraddligands	442	phaseselectexcludedfeature	464
phaseqsaraddligandsfromfile	443	phaseselectextendconstraintrow	465
phaseqsaraddtohypothesis	443	phaseselectextendevrow	465
phaseqsaraddtoproject	444	phaseselectextendhypothesisrow	465
phaseqsarbuild	444	phaseselectextendtablerow	466
phaseqsardelete	444	phaseselecthypothesisrow	466
phaseqsardeleteall	444	phaseselectonlyevrow	466
phaseqsarexport	445	phaseselectonlyhypothesisrow	467
phaseqsarexportligandstofile	445	phaseselectonlytablerow	467
phaseqsarextendinclude	445	phaseselectsite	467
phaseqsarextendselectligandrow	446	phaseselecttablerow	468
phaseqsarimport	446	phaseseparatestereoisomers	468
phaseqsarinclude	446	phaseseparatetitles	468
phaseqsarinvertselection	447	phasesetactiverows	468
phaseqsarplotactivity	447	phasesetactivity	469
phaseqsarpredict	447	phasesetactivityproperty	469
phaseqsarrandomtraining	448	phasesetactivitythresholds	470
phaseqsarsearch	448	phasesetalignmentoptions	470
phaseqsarselectall	448	phasesetconstrainttolerance	471
phaseqsarselecthypothesis	448	phasesetexcludedvolumes	471
phaseqsarselectligandrow	449	phasesetfrequency	472
phaseqsarsetactivity	449	phasesettrainingrows	472
phaseqsarsettrainingrows	450	phaseshowclusters	473
phaseqsarsort	450	phasesiteoptions	473
phaseqsartest	450	phasesortable	474
phaseqsartoggleinclude	451	phasetstepforward	475
phaseqsartoggleligandrow	451	phasetstepgoto	475
phaseqsartogglerow	451	phasetoggleactivetablerow	475
phaserandomtraining	452	phasetogglecentroidatom	476
phaserandomtrainingoptions	452	phasetoggleconstraintmark	476
phaserrefreshfrequencytable	452	phasetoggleexcludedfeature	476
phaseremoveligands	453	phasetogglerequiredmatch	477
phaserenameligand	453	phasetoggleselectconstraintrow	477
phaserescorehypotheses	453	phasetoggleselectsite	477
phaserescoreweighting	454	phasetoggletrainingrow	478
phaseresetfeatures	455	phasetranslateexcludedvolumes	478
phaseruncreate	455	phasetranslatesite	478
phaserundelete	456	phaseundisplayproperty	479
phaserunopen	456	phaseungroupactives	479
phaserunrename	456	phaseunmarkfeature	479
phaserunsaveas	456	phaseunselectevrow	480

phaseunselecthypothesisrow	480	plotxyunlabel	502
phaseunselecttablerow	480	plotxyunselect	503
picksize	481	plotxyupdate	503
place	481	plotxyzoom	503
plotxyarrangecolumn	481	plotxyzzoompan	504
plotxyarrangerow	482	potential	504
plotxyarrangetiled	482	pprep	506
plotxyaspectratiolock	482	pprepreceptorligand	507
plotxyaxis	482	pprepwrite	507
plotxyaxisautorange	483	prefer	508
plotxyaxisdelete	484	previousresidue	544
plotxyaxisdisplay	484	primexaddwaters	544
plotxyaxisrename	485	primexcalculateddensityfit	545
plotxycaption	486	primexcalculateddensitypeaks	545
plotxycaptionposition	486	primexcalculaterfactors	546
plotxycopy	486	primexcreatemap	546
plotxydelete	487	primexdecreaseisovalue	549
plotxydeleteall	487	primexdeleteallrigidbodies	549
plotxydisplay	487	primexdeleteoccupancy	550
plotxydisplaycaption	487	primexdeleterigidbodies	550
plotxydisplayincluded	488	primexentrywscreate	550
plotxydisplaylegend	488	primexexport	550
plotxydisplayname	488	primexincreaseisovalue	551
plotxydisplayonly	489	primexinputdata	552
plotxydisplaypointlabels	489	primexmarkersettings	555
plotxydisplayselected	489	primexmutatemodeltosequence	556
plotxyequalizeaxisrange	490	primexncsaddrow	556
plotxyhidesidebar	490	primexncsdeleterows	557
plotxyhidetoolbar	490	primexncseditcell	557
plotxylabel	490	primexncsselecttextendrow	557
plotxymove	491	primexncsselectonlyrow	558
plotxynew	492	primexncsselectrow	558
plotxypan	492	primexncstogglobbonly	558
plotxypanplot	492	primexncsunselectrow	559
plotxyrename	493	primexpartialoccupancy	559
plotxyresetview	493	primexplaceligand	559
plotxysaveimage	493	primexpublish	561
plotxyselect	494	primexrealspace	561
plotxyseries	494	primexreciprocal	564
plotxyseriesdelete	498	primexselectdensityblob	569
plotxyseriesdisplay	498	primexselectdensitypeak	569
plotxyseriesdisplay45	499	primexselecttextenddensitypeak	569
plotxyseriesdisplayregression	499	primexselecttextendligandsolvent	570
plotxyseriesrename	499	primexselectligandsolvent	570
plotxyseriesselect	500	primexselectonlydensitypeak	570
plotxyseriesselectadd	500	primexselectonlyligandsolvent	570
plotxyseriesselectaxesequal	501	primexselectonlyrigidbody	571
plotxyshowsidebar	501	primexselectrigidbody	571
plotxyshowtoolbar	501	primexsetoccupancy	571
plotxytoggledisplay	502	primexsettings	572
plotxyundisplay	502	primexswitchaltpositions	575

Command Index

primextoggleselectrigidbody	575	propertycalculate	598
primexunitcell	575	propertycalculatetchirality	600
primexunselectdensitypeak	575	propertycalculatemolformula	600
primexunselectligandsolvent	576	propertycalculatpickpka	601
primexwriteaddwaters	576	propertycalculatesubstructs	602
primexwritecalculaterfactors	576	propertyclearvalue	603
primexwritecreatemap	576	propertycreate	603
primexwriteexport	577	propertydelete	604
primexwriteplaceligand	577	propertygeneratecontacts	605
primexwriterealspace	577	propertygeneratehbond	605
primexwriterecipocal	577	propertymeasurementsetting	605
profile	577	propertymove	606
projectbackup	578	propertyprecision	606
projectbackuppreference	579	propertyrename	606
projectcheck	581	propertyshowall	607
projectclose	582	propertysuperimposesetting	607
projectcopy	582	protassign	608
projectdelete	584	protassignresidues	608
projectmerge	584	protassignstart	609
projectnew	585	protassignwrite	609
projectopen	585	proteinsculpting	609
projectprefer	585	psp_helical	609
projectpublish	587	psp_pairwise	610
projectrename	588	psp_spatial	611
projectrestore	588	pspalignaddanchor	612
projectrevertopen	588	pspaligndeleteanchor	612
projectrevertsnapshot	589	pspaligninsertgaps	613
projectrevertstate	589	pspalignlockgaps	613
projectsaveas	589	pspalignmoveleft	613
projectsceneclose	589	pspalignmoveleftblock	614
projectscenedelete	590	pspalignmoveright	615
projectscenedescription	590	pspalignstructures	615
projectscenedragrows	590	pspalignunlockgaps	615
projectsceneexportpresentation	591	pspbstogglehetatom	616
projectsceneexportproject	591	pspbuiltbackbone	616
projectsceneextendselect	591	pspbuiltconsensusmodel	616
projectscenenew	592	pspbuiltheteromultimer	616
projectsceneoptions	592	pspbuiltathomomultimer	617
projectscenereplace	593	pspbuiltstructure	617
projectsceneseselect	593	pspconstraintmarkersettings	617
projectscenetoggleselect	593	pspeditquerysequence	621
projectsceneview	594	pspexcludetablelrow	621
projectsceneviewnamed	594	pspexcludetablelrow	622
projectstoresnapshot	594	pspexportalignment	622
projectstorestate	594	pspfindfamily	623
projectstynchronize	595	pspfindhomologs	623
projecttablefind	595	pspfoldrecognitionoptions	623
projecttablereplace	597	pspfoldrecognitionsearch	623
projecttablereplaceall	597	pspgetquerysequence	624
projectupdatecoordinates	597	pspimportalignment	624
projectupdateviews	597	pspimporthomolog	624

[illegible]

Command Index

retype 680
 ribbon 681
 ribbondump 692
 ringclosure 692
 ringflip 693
 rotate 693

S

saveimage 694
 saveimageheight 695
 saveimagewidth 695
 savelayout 695
 savemovie 696
 savemovieoptions 696
 scanmode 697
 scriptlogfile 698
 scriptrun 698
 searchdbconfgn 699
 selectatomsforediting 702
 selecteditatom 703
 selecteditproperty 703
 sequenceviewersaveimage 703
 set 704
 setread 704
 setwrite 705
 showdockablepanels 705
 showfirstselectedentry 705
 showhwstereosetup 705
 showmarkers 706
 showpanel 706
 showpanels 706
 showpropertiesselectedentries 707
 showproperty 707
 showtoolbar 707
 sorteditatomcolumn 708
 specifiedname 708
 spotcenter 708
 spotcenterpoint 709
 startdiagnosis 709
 stopdiagnosis 710
 strikebuildqsar 710
 strikedeletemodel 711
 strikeexportmodel 712
 strikeextendselectdescriptor 712
 strikeimportmodel 712
 strikeplotmodel 713
 strikepredict 713
 strikeselectdescriptor 713
 strikeselectmodel 714
 strikesimilarity 714

striketoggleselectdescriptor 714
 structalignatoms 715
 structalignatomsalt 715
 structalignstart 715
 substructure 716
 substructurefileread 716
 substructurefilewrite 717
 substructureshell 718
 superimpose 719
 superimposeatom 719
 superimposeset 720
 superimposesmarts 720
 surfaceactivegrid 720
 surfacecontours 722
 surfacedarkencolor 723
 surfacedelete 724
 surfacedisplay 724
 surfaceduplicate 725
 surfaceextended 725
 surfaceextendedradiuscontext 726
 surfaceextendedradiusset 727
 surfacemolecular 727
 surfacemolecularcontext 728
 surfacemolecularset 729
 surfacemolecularset2 729
 surfaceoption 730
 surfacepbrsp 735
 surfacepbrspset 737
 surfacereaname 738
 surfaceresolution 738
 surfacescheme 739
 surfacesetcomment 741
 surfacesetisovalue 742
 surfacesettings 742
 surfacesetviewas1 744
 surfacesplitpair 745
 surfacestyle 746
 surfacetransparency 746
 surfacetransparencyback 747
 surfacevdw 747
 surfacevdwcontext 749
 surfacevdwset 749
 surfaceviewas1set 749
 swapptworkspace 750
 symmetrizeworkspace 750
 system 750

T

tablealigncolumn 751
 tablecopyprop 751

tablegroupsort	752
tableresizecolumn	753
tablesort	753
tablesortall	754
tablesortfields	755
tablesortgroupfields	756
tableunselectnonsubset	756
targetentry	757
tile	757
tileapplyview	757
tileapplyzoom	758
tilemode	758
tilemove	759
tilereset	759
tileresetall	760
timingsetup	760
timingstart	761
timingstop	761
toggleseditatomboolean	761
togglesselecteditatom	762
togglesselecteditproperty	762
torsioncheck	762
torsiongroup	763
trajectoryexport	763
trajectoryexportatoms	763
trajectoryexportmovieoptions	764
trajectoryexportoptions	764
trajectoryplayerdisplayatoms	765
trajectoryplayergoto	765
trajectoryplayergotofirst	766
trajectoryplayergotolast	766
trajectoryplayersettings	766
trajectoryplayerstepahead	769
trajectoryplayerstepback	770
trajectoryplayersuperimposeatoms	770
trajectoryresolvemeasurements	770
transform	771
translate	773

U

undisplayatom	773
undo	773
ungroupentries	774
unhookimport	774
uniquepdb	774
uniquepdb	775
updateribbons	775

V

varymolecule	775
varytorsion	776
vcsaddattachment	777
vcsaddcorefromproject	777
vcsaddmincapcore	778
vcsaddoriginalcore	778
vcscanceldockjob	778
vcsclearreagentfile	778
vcscombiexportdockingfile	778
vcscombiexportdockingproject	779
vcscombiexportoptions	779
vcsconfiguredocking	780
vcscoreoptions	783
vcscreatedockedlibrary	784
vcsdeleteattachment	784
vcsdeletecore	784
vcsdeleterresults	784
vcsdisplayreceptor	785
vcsdocking	785
vcsenumeratedockoptions	786
vcsexcludetablerow	786
vcsexportdefinition	786
vcsexportresults	787
vcsexportspreadsheet	787
vcsimportdefinition	787
vcsimportrun	788
vcsincludeextendtablerow	788
vcsincludeonlytablerow	788
vcsincludetablerow	789
vcsinverttableselection	789
vcsoptions	789
vcsrefreshstructure	790
vcsrenameattachment	790
vcsrestorerresults	791
vcsruncombinationdocking	791
vcsruncombinatorialdocking	791
vcsruncombinatorialeselection	792
vcsruncreate	792
vcsrundelete	792
vcsrunenumerateddocking	793
vcsrunopen	793
vcsrunrename	793
vcsrunsas	794
vcsrunsingledocking	794
vcsrunsingleselection	794
vcssaverresults	795
vcselectalltablerows	795
vcselectextendtablerow	795
vcselectonlytablerow	795
vcselecttablerow	796

Command Index

vcsetattachmentfile 796
vcsetmolecule 796
vcsetreagentfile 797
vcssingleexportdockingfile 797
vcssingleexportdockingproject 798
vcssingleexportoptions 798
vcssorttable 798
vcstepforward 799
vcstepgoto 799
vcundisplayreceptor 799
vcunselecttable 800
vcswritecombinationdocking 800
vcswritecombinatorialdocking 800
vcswriteenumerateddocking 800
vcwritesingledocking 801
viewcrystal 801
viewdelete 802
viewdeleteselect 802
viewdragselection 802
viewexport 803
viewextendselect 803
viewextendselectrow 803
viewimport 804
viewmatrix 804
viewplayergotofirst 805
viewplayergotolast 805
viewplayerplaybackward 805
viewplayerplayforward 805
viewplayersettings 806
viewplayersteppahead 806
viewplayerstepback 806
viewplayerstop 806
viewrecordmovie 807
viewrename 807
viewreset 807
viewresetnofit 808
viewrestore 808
viewrestorerow 808
viewsave 808
viewselect 809
viewselectonly 809
viewselectonlyrow 810
viewselectrow 810
viewsetpause 810
viewsetplay 811
viewsetrowpause 811
viewtrajectory 811
viewunselect 812
viewunselectrow 812
viewunsetplay 812
viewvolume 813
visexport 813

visimport 814
volumedelete 815
volumedisplay 816
volumeisosurface 816
volumerename 817
volumesetcolors 817

W

workspacecopy 818
workspacefontdecrease 818
workspacefontincrease 819
workspacepaste 819
workspacepastesettings 819
workspacescopeoff 820
workspacescopeon 820
workspacescopesettings 820
workspacescopeview 822
workspaceselectionadd 823
workspaceselectionclear 823
workspaceselectioninvert 823
workspaceselectionreplace 823
workspaceselectionsubtract 824
workspacetile 824
writeangle 825
writecontact 826
writecoupling 826
writedihehdral 827
writedistance 828
writebond 829
wsassistantapply 830
wsassistantapplycontactbond 830
wsassistantapplyfocus 830
wsassistantbindingsurface 830
wsassistantcontactbond 832
wsassistantcreatesurfaces 833
wsassistantdefaultstyle 833
wsassistantfixreceptor 833
wsassistanthideallsurfaces 833
wsassistantligandddetection 834
wsassistantmaterial 835
wsassistantprimarystyle 835
wsassistantrendering 837
wsassistantshowallsurfaces 838
wsassistantwsstyle 839

X

xcluster 846
xclusterstart 847
xclusterwrite 847

Table of Contents

1	The Maestro Command Language	1
2	Conventions Used in this Manual	3
3	The Atom Specification Language	5
3.1	Why an Atom Specification Language?	5
3.2	The ASL Hierarchy	5
3.3	Atom Specification	5
3.4	Operators	13
3.4.1	The Boolean <code>and</code> Operator	13
3.4.2	The Boolean <code>or</code> Operator	14
3.4.3	The Boolean <code>not</code> Operator	14
3.4.4	The <code>fillres</code> and <code>fillmol</code> Operators	14
3.4.5	The <code>within</code> and <code>beyond</code> Operators	14
3.4.6	The <code>withinbonds</code> Operator	15
3.4.7	The <code>beyondbonds</code> Operator	15
3.5	Operator Priority	16
3.6	Implicit Operators	16
3.7	Creating New Sets from Existing Ones	17
3.8	Special Specifications	17
3.9	Matching PDB Atom Names	18
3.10	Miscellaneous	18
3.10.1	Atoms not yet present in a structure	18
3.10.2	Aliasing	18
3.11	Useful Hints when using ASL with the Project Facility	19
3.12	ASL Examples	20
4	The Entry Specification Language	25
4.1	Why an Entry Selection Language?	25
4.2	Entry Properties	25
4.3	Logical Operators	26
4.4	Entry Property Comparisons	26
4.5	Examples	27
5	Commands	29
	<code>lddataset</code>	29
	<code>ldplot</code>	30
	<code>ldrescale</code>	32
	<code>ldtable</code>	33

2ddataset	33
2dplot	34
2drescale	36
addatomproperty	36
addfromhold	37
adjustangle	37
adjustchirality	38
adjustdihedral	38
adjustdistance	39
alias	40
alignaxis	40
alignplane	41
angle	41
appendribbons	42
application	42
assigncharges	43
atom	43
atomname	44
atomproplabel	44
attach	45
attachmentmarkerdump	45
attachmentmarkersettings	45
autosetup	48
beginundoblock	49
bmincomfile	49
bond	50
bondorder	50
bondoutline	50
bondtonew	52
buildoptions	53
calcenergy	53
canonicalname	54
caption	54
captionbackground	55
captioncolor	55
captiondelete	56
captionfont	56
captionhide	57
captionposition	57
captionselect	57
captionselectadd	58
captionshow	58
captionsoff	58
captionson	58
captiontext	58

captionunselect	59
captionunselectall	59
cascadepanels	59
cd	59
cellsmarkerdump	59
cellsmarkersettings	60
centeratom	63
centerbond	63
centercoordinates	64
centroid	64
centroidatom	65
centroidatomset	65
cglidedockconstraintposition	65
cglidedockconstraintregion	67
chainname	68
changedirectory	68
chdir	69
chiralatom	69
clip	69
clipatomset	71
clipsurfaceset	71
cliptoallatoms	72
cliptoligand	72
cliptoselectedatoms	72
clusteratom	72
clusterheavyatoms	72
clusterset	73
clustertorsion	73
color	74
coloratom	74
coloratombonds	75
coloratomrgb	75
colorbond	76
colorscheme	77
combilibenum	77
combilibenumaddattachment	78
combilibenumclearreagentfile	78
combilibenumdeleteattachment	79
combilibenumexportdefinition	79
combilibenumimportdefinition	79
combilibenumoptions	79
combilibenumrefreshstructure	80
combilibenumrenameattachment	80
combilibenumselectextendtable	81
combilibenumselectonlytable	81

combibibenumselecttablerow	81
combibibenumsetmolecule	82
combibibenumsetreagentfile	82
combibibenumunselecttablerow	82
compareatom	83
compareset	83
confelim	84
confelimstart	85
confelimwrite	85
confgenltsearch	85
confgenmini	88
confgenpotential	89
confgenreadpotential	91
confgenstart	91
confgenstd	91
confgenstdsettings	91
confgenstdwrite	93
confgenwrite	93
confsearch	93
connect	96
connectfuseatom	97
constrainedangle	97
constrainedatom	98
constraineddist	99
constrainedset	100
constrainedtorsion	100
contactcriteria	101
contactset1	103
contactset2	103
coordinatescan	104
coupling	105
createcustompropertygroup	106
createlibrary	106
createpropsubset	106
createsubset	107
defaultfc	107
delete	107
deletcustompropertygroup	108
deleteemptyentrygroups	109
deleteproperty	109
deleteselectedproperties	109
deletezeroorderbonds	109
densityblobupdate	109
dialogoptions	110
dihedral	110

dihedraldrive	111
dipolemoment	112
displayatom	117
displayonlyatom	118
displayopt	118
distance	122
distancecheck	124
dockgridbox	125
doplanealignment	126
dynamics	126
ecalc	127
editatomproperty	128
ejob	128
elementlabels	130
embrace	130
embracecsearch	131
emptyviewreset	133
endundoblock	133
energykill	134
energymonitor	134
energysleep	134
energystart	135
energystop	135
energytask	136
energyupdate	136
energyupdatejobstatus	137
energywake	137
enhance3d	137
enterposeviewermode	138
entryadddisplay	138
entrycolorbyproperty	138
entrycombine	140
entrycopyprop	140
entrycreatesmoothed	141
entrydelete	142
entrydisassociate	142
entrydisassociatebychain	143
entrydisassociatebyligand	143
entrydisplaysurfaces	143
entrydragselection	144
entryduplicate	144
entryexport	145
entryexportcanvas	149
entryexportspreadsheets	149
entryextendsselection	150

entryextendselectrow	150
entryextendwsinclude	151
entrygroupbyproperty	151
entrygroupcollapse	152
entrygroupcountcontacts	153
entrygroupcounthbonds	153
entrygroupcreate	153
entrygroupdelete	154
entrygroupdragselection	154
entrygroupduplicate	155
entrygroupexpand	155
entrygroupexpandonly	156
entrygroupextendselection	156
entrygroupinvertselection	157
entrygroupmove	157
entrygrouprename	158
entrygroupselect	158
entrygroupselectonly	159
entrygroupsettitle	159
entrygroupungroup	159
entrygroupunselect	160
entrygroupviewposes	160
entrygroupwsexclude	161
entrygroupwsinclude	161
entrygroupwsincludeonly	162
entryimport	162
entryimportspreadsheet	167
entryimportvibration	168
entryinvertselection	168
entrymergeprop	168
entrymovetogroup	169
entryremovedisplay	170
entryrename	170
entryresetcolor	171
entryselect	171
entryselectall	171
entryselectonly	172
entryselectonlyrow	172
entryselectrandom	172
entryselectrow	173
entrysetcolor	174
entrysetdeletable	174
entrysetnondeletable	175
entrysetprop	175
entrysetreadonly	176

entrysettitle	176
entrysetwritable	176
entryshowall	177
entrystars	177
entrytable	177
entryundisplaysurfaces	178
entryunselect	178
entryunselectall	178
entryunselectrow	179
entrywatermapexaminerresults	179
entrywscreate	179
entrywsexclude	180
entrywsexcludeunlocked	180
entrywsinclude	180
entrywsincludeunlock	181
entrywsincludeonly	181
entrywsincludeunlock	181
epik	182
epikreadinput	184
epikstart	184
epikwriteinput	184
eplayergoto	185
eplayergotofirst	185
eplayergotolast	185
eplayernextgroup	186
eplayerpreviousgroup	186
eplayersettings	186
eplayersteppahead	189
eplayerstepback	189
errorcheck	189
excludedvolumesmarkerssettings	190
exportrotatablegroups	191
extendselecteditatom	192
fieldqsarmarkerssettings	192
fileread	195
filewrite	197
filter	200
find	200
findnext	204
findprevious	204
fit	204
fitbox	205
fitligand	205
fitnextligand	205
fitoptions	205

fitpreviousligand	206
flexibleligandalignment	206
forcefield	207
forcefieldbend	207
forcefieldele	208
forcefieldgbsolv	209
forcefieldimproper	209
forcefieldsasolv	210
forcefieldstretch	211
forcefieldtorsion	211
forcefieldvdw	212
forcefieldview	213
forcefieldwilson	213
formalcharge	214
fragment	214
fragmentfuse	214
fragmenttype	215
frozenatom	215
frozenset	216
fuse	216
geometrycleanup	216
glideactivesiteres	216
glideactivesiteset	217
glidecalcboxfromligand	217
glideconstraintatomlabel	217
glideconstraintposition	218
glideconstraintregion	218
glideconstraintregionactivecell	219
glideconstraintregiongrow	219
glideconstraintregioninvisible	219
glideconstraintregionshrink	220
glideconstraintregionvisible	220
glidecoremolecule	220
glidecorermssubset	221
glidedisplayreceptor	221
glidedockconstraintposition	221
glidedockconstraintregion	223
glidedockexcludedvolume	224
glidedockingreset	225
glidedockregionnumatoms	225
glidedockregionselect	225
glidedockregionunselect	226
glideexcludedvolume	226
glidegridhydrophobic	227
glidegridreset	228

glideligand	228
glideligandfile	231
glideligandgrid	231
glidelocatehydrophobiccells	232
glideoutput	233
glidereceptorligand	235
glidereferenceligand	236
glideringsampling	237
gliderotatablegroups	237
gliderotatablegroupsextenduserow	238
gliderotatablegroupstoggleatom	238
gliderotatablegroupstogglerow	238
gliderotatablegroupsuseexisting	239
glidescoring	239
glidesettings	240
glidesite	245
glidetorsionalconstraint	248
glidetorsionalconstraintallbonds	249
glidetorsionalconstraintpatternadd	249
glidetorsionalconstraintpatterndelete	249
glidetorsionalconstraintpatterndeleteall	250
glidetorsionalconstraintpatterndrag	250
glidetorsionalconstraintpatternreplace	250
glidetorsionalconstraintpatternselect	250
glidetorsionalconstrainttorsionadd	251
glidetorsionalconstrainttorsiondelete	251
glidetorsionalconstrainttorsiondeleteall	251
glidetorsionalconstrainttorsiondeleteinvalid	252
glidetorsionalconstrainttorsionreplace	252
glidetorsionalconstrainttorsionselect	252
glidetorsionalconstrainttorsionsetangle	253
glideundisplayreceptor	253
googleprofiling	253
graphicsdebugdisable	254
graphicsdebugenable	254
grow	254
growbond	255
growbond2	255
growdirection	255
growname	256
happly	256
hbondcriteria	256
hbondset1	257
hbondset2	258
help	258

helpauto	258
helpcategory	259
helpsearch	259
helptopic	260
hidedockablepanels	260
hideentries	260
hidemarkers	261
hidepanel	261
hidepanels	261
hideproperty	261
hidetoolbar	262
highlight	262
highlightarrowatoms	262
highlightarrowcolor	263
highlightarrowsettings	263
highlightatoms	264
highlightbackground	264
highlightcolor	265
highlightdelete	265
highlightdirection	265
highlightduplicate	266
highlighthide	266
highlightmethod	266
highlightpreference	267
highlightramp	267
highlightrename	268
highlightselect	268
highlightshow	268
highlightsoff	269
highlightson	269
highlighttext	269
highlighttextcolor	269
highlighttextfont	270
highlighttextposition	271
highlighttextselect	271
highlighttextselectadd	272
highlighttextsettings	272
highlighttextunselect	272
highlighttextunselectall	272
highlightwidth	273
historyvisible	273
hold	280
hppmap	280
hppmapbox	281
hppmapset	282

hppmapstart	282
hppmapwrite	282
htreat	282
hydrogenapply	283
hydrogentreat	283
impactbufferedatom	283
impactbufferedset	284
impactconstraints	284
impactcontinuumsolvent	285
impactdynamics	286
impactfastmultipole	288
impactfrozenatom	289
impactfrozenet	289
impacthybridmc	290
impactjob	290
impactmdparams	291
impactminimization	293
impactperiodicboundary	296
impactpotential	297
impactread	299
impactrepexch	299
impactselectextendablerow	300
impactselectonlyablerow	301
impactstart	301
impacttask	301
impacttruncation	302
impactwrite	303
invert	303
invertset	303
jaguarassignatomnames	303
jaguarimportgeometry	304
jaguarinputfilesjob	305
jaguarjob	305
jaguarselectextendablerow	306
jaguarselectonlyablerow	307
jobcleanup	307
jobsettings	308
kill	310
labelatom	310
labelatomfieldsnone	320
labelatomoffset	321
labelatomoffsetreset	321
labelclear	321
labeldisplayusertext	322
labelformat	322

labelupdate	323
labelusertextatomupdate	323
liaisonanalysis	323
liaisonparameters	326
liaisonselectlig	329
liaisonsettings	329
liaisonsystem	330
ligandbond	331
light	332
lightambient	332
lightdiffuse	333
lightposition	333
lightspecular	334
ligprep	334
ligprepread	336
ligprepstart	337
ligprepwrite	337
localatoms	337
localbitset	337
localcenter	338
logfile	338
logp	338
macrodefine	339
macrorun	339
makedirectory	340
markers	340
material	341
materialdefault	342
materialmolecular	342
mcsd	342
mini	343
minienergy	343
minimize	344
minta	345
monitor	346
monitorangle	346
monitordistance	346
monitorhbond	347
monitorsetsurf	347
monitorsettings	348
monitorsurf	349
monitortorsion	349
mouse	349
move	350
multiplemini	351

mutate	352
nextresidue	352
optimizefog	352
partialcharge	353
pause	353
pausecommands	353
peptidecistrans	354
peptiderotate	354
perresiduedisplay	355
perresiduesubset	355
perresiduetype	356
phaseaddconstraint	356
phaseaddcustomfeature	357
phaseaddhypotheses	357
phaseaddligands	357
phaseaddligandsfromfile	358
phaseaddligandsfromrun	359
phaseaddsite	359
phasealignhypotheses	360
phasebuildqsar	361
phasechangesitetype	364
phasechooseactiveset	364
phasecleanupstructures	364
phaseclearcentroidatoms	365
phaseclusterhypotheses	365
phaseconfgn	365
phasecopyhypothesisfrombuilder	369
phasecreateevactives	369
phasecreateevactivesettings	369
phasecreateevreceptor	370
phasecreateevreceptoratoms	370
phasecreateevreceptorsettings	370
phasecreateevreference	372
phasecreateevreferenceatoms	372
phasecreateevreferencesettings	373
phasecreateexcludedvolume	373
phasecreatefreestylehypothesis	373
phasedbaddconfset	374
phasedbaddligands	374
phasedbaddligandsfromdb	375
phasedbconfgn	375
phasedbcreatesites	379
phasedbcreatesubsetfromhitfile	379
phasedbcreatesubsetfromselected	380
phasedbcreatesubsetfromtextfile	380

phasedbdeleteligands	380
phasedbdeletesubset	381
phasedbdisplayrange	381
phasedbexportligands	381
phasedbfilter	382
phasedbfilteradd	382
phasedbgenerateconformers	382
phasedbnewfromfile	383
phasedbopen	383
phasedbremoveconformers	383
phasedbselectextendablerow	384
phasedbselectonlyablerow	384
phasedbselectsuubsetrow	385
phasedbselectablerow	385
phasedbunselectablerow	385
phasedeletecustomfeature	385
phasedeleteexcludedvolumes	386
phasedeletehypothesis	386
phasedeletehypothesisfrombuilder	386
phasedeletesectedconstraints	387
phasedeletesectedligands	387
phasedeletesites	387
phasedisplayproperty	387
phaseevoptions	388
phaseexcludablerow	388
phaseexportalignmentstofile	388
phaseexportconformerstofile	389
phaseexportfeature	389
phaseexporthypothesis	389
phaseexporthypothesisfrombuilder	390
phaseexportselectedalignments	390
phaseexportselectedhypotheses	390
phaseexportselectedhypothesesfrombuilder	391
phaseexportselectedligands	391
phaseextendsite	391
phasefieldqsar	392
phasefieldqsaraddligands	394
phasefieldqsaraddligandsfromfile	395
phasefieldqsaraddtoproject	396
phasefieldqsarbuild	396
phasefieldqsardelete	396
phasefieldqsardeleteall	396
phasefieldqsarexport	397
phasefieldqsarexportligandstofile	397
phasefieldqsarexextendinclude	397

phasefieldqsarextendselectligandrow	398
phasefieldqsarimport	398
phasefieldqsarinclude	398
phasefieldqsarinvertselection	399
phasefieldqsarplotactivity	399
phasefieldqsarpredict	399
phasefieldqsarrandomtraining	400
phasefieldqsarselectall	400
phasefieldqsarselectcustom	400
phasefieldqsarselectligandrow	401
phasefieldqsarsetactivity	401
phasefieldqsarsettrainingrows	401
phasefieldqsarsort	402
phasefieldqsartest	402
phasefieldqsartoggleinclude	402
phasefieldqsartoggleligandrow	403
phasefieldqsartogglerow	403
phasefieldqsartoggleselectcustom	403
phasefieldqsarvisselectcontour	404
phasefieldqsarvisselectintensity	404
phasefieldqsarvistoggleselectcontour	404
phasefindmatches	405
phasefindpharmacophores	410
phasegenerateconformers	411
phasesgroupligands	411
phasesgroupligandsbytitle	411
phasehookimport	412
phasehypothesiscreateexcludedvolume	412
phasehypothesisdeleteexcludedvolumes	412
phasehypothesisselectevrow	412
phasehypothesisselectextendevrow	413
phasehypothesisselectonlyevrow	413
phasehypothesisselectrow	413
phasehypothesissetexcludedvolumes	414
phasehypothesissetid	414
phasehypothesissetsitemask	415
phasehypothesissettolerances	415
phasehypothesistoggleinclude	417
phasehypothesisunselectevrow	417
phaseimportfeature	418
phaseimportthypothesis	418
phaseincludeextendablerow	418
phaseincludeonlyablerow	419
phaseincludetablerow	419
phaseinverttableselection	419

phaseligprep	420
phaselocatesites	420
phasemarkerdump	421
phasemarkersettings	421
phasemarkfeature	432
phasemergestereoisomers	433
phasemergetitles	433
phaseoptions	433
phasepatterndelete	436
phasepatternedit	436
phasepatternmovedown	437
phasepatternmoveup	438
phasepatternnew	438
phasepatternsetoptions	439
phaseplotactivity	440
phasepreview	440
phaseqsar	440
phaseqsaraddligands	442
phaseqsaraddligandsfromfile	443
phaseqsaraddtohypothesis	443
phaseqsaraddtoproject	443
phaseqsarbuild	444
phaseqsardelete	444
phaseqsardeleteall	444
phaseqsarexport	444
phaseqsarexportligandstofile	445
phaseqsarextendinclude	445
phaseqsarextendselectligandrow	445
phaseqsarimport	446
phaseqsarinclude	446
phaseqsarinvertselection	446
phaseqsarplotactivity	447
phaseqsarpredict	447
phaseqsarrandomtraining	448
phaseqsarsearch	448
phaseqsarselectall	448
phaseqsarselecthypothesis	448
phaseqsarselectligandrow	449
phaseqsarsetactivity	449
phaseqsarsettrainingrows	449
phaseqsarsort	450
phaseqsartest	450
phaseqsartoggleinclude	450
phaseqsartogglegligandrow	451
phaseqsartogglrow	451

phaserandomtraining	451
phaserandomtrainingoptions	452
phaserrefreshfrequencytable	452
phaserremoveligands	452
phaserrenameligand	453
phaserescorehypotheses	453
phaserescoreweighting	453
phaserresetfeatures	455
phaseruncreate	455
phaserundelete	455
phaserunopen	456
phaserunrename	456
phaserunsaveas	456
phaserunsetseed	457
phaserunsetstereo	457
phasesavematching	457
phasescorehypotheses	457
phasescoreinactives	462
phasescoreselecthypothesis	462
phasesearchformatches	463
phaseselectalltablerows	463
phaseselectconstraintrow	463
phaseselectevrow	464
phaseselectexcludedfeature	464
phaseselectextendconstraintrow	464
phaseselectextendevrow	465
phaseselectextendhypothesisrow	465
phaseselectextendtablerow	465
phaseselecthypothesisrow	466
phaseselectonlyevrow	466
phaseselectonlyhypothesisrow	466
phaseselectonlytablerow	467
phaseselectsite	467
phaseselecttablerow	467
phaseseparatestereoisomers	468
phaseseparatetitles	468
phasesetactiverows	468
phasesetactivity	469
phasesetactivityproperty	469
phasesetactivitythresholds	469
phasesetalignmentoptions	470
phasesetconstrainttolerance	471
phasesetexcludedvolumes	471
phasesetfrequency	472
phasesettrainingrows	472

phaseshowclusters	473
phasesiteoptions	473
phasesorttable	474
phasetepforward	474
phasetepgoto	475
phasetoggleactivetablerow	475
phasetogglecentroidatom	475
phasetoggleconstraintmark	476
phasetoggleexcludedfeature	476
phasetogglerequiredmatch	476
phasetoggleselectconstraintrow	477
phasetoggleselectsite	477
phasetoggletrainingrow	477
phasetranslateexcludedvolumes	478
phasetranslatesite	478
phaseundisplayproperty	479
phaseungroupactives	479
phaseunmarkfeature	479
phaseunselectevrow	480
phaseunselecthypothesisrow	480
phaseunselecttablerow	480
picksize	480
place	481
plotxyarrangecolumn	481
plotxyarrangerow	481
plotxyarrangetiled	482
plotxyaspectratiolock	482
plotxyaxis	482
plotxyaxisautorange	483
plotxyaxisdelete	484
plotxyaxisdisplay	484
plotxyaxisrename	485
plotxycaption	485
plotxycaptionposition	486
plotxycopy	486
plotxydelete	486
plotxydeleteall	487
plotxydisplay	487
plotxydisplaycaption	487
plotxydisplayincluded	487
plotxydisplaylegend	488
plotxydisplayname	488
plotxydisplayonly	488
plotxydisplaypointlabels	489
plotxydisplayselected	489

plotxyequalizeaxisrange	489
plotxyhidesidebar	490
plotxyhidetoolbar	490
plotxylabel	490
plotxymove	491
plotxynew	491
plotxypan	492
plotxypanplot	492
plotxyrename	492
plotxyresetview	493
plotxysaveimage	493
plotxyselect	493
plotxyseries	494
plotxyseriesdelete	498
plotxyseriesdisplay	498
plotxyseriesdisplay45	498
plotxyseriesdisplayregression	499
plotxyseriesrename	499
plotxyseriesselect	500
plotxyseriesselectadd	500
plotxyseriessetaxesequal	501
plotxyshowsidebar	501
plotxyshowtoolbar	501
plotxytoggledisplay	501
plotxyundisplay	502
plotxyunlabel	502
plotxyunselect	502
plotxyupdate	503
plotxyzoom	503
plotxyzoompan	503
potential	504
pprep	506
pprepreceptorligand	507
pprepwrite	507
prefer	507
previousresidue	544
primexaddwaters	544
primexcalculateddensityfit	545
primexcalculateddensitypeaks	545
primexcalculaterfactors	545
primexcreatemap	546
primexdecreaseisovalue	548
primexdeleteallrigidbodies	549
primexdeleteoccupancy	549
primexdeleterigidbodies	550

primexentrywscrate	550
primexexport	550
primexincreaseisovalue	550
primexinputdata	551
primexmarkersettings	554
primexmutatemodeltosequence	556
primexncsaddrow	556
primexncsdeleterows	556
primexncseditcell	557
primexncsselecttextendrow	557
primexncsselectonlyrow	558
primexncsselectrow	558
primexncstogglebbonly	558
primexncsunselectrow	558
primexpartialoccupancy	559
primexplaceligand	559
primexpublish	561
primexrealspace	561
primexreciprocal	563
primexselectdensityblob	569
primexselectdensitypeak	569
primexselecttextenddensitypeak	569
primexselecttextendligandsolvent	569
primexselectligandsolvent	570
primexselectonlydensitypeak	570
primexselectonlyligandsolvent	570
primexselectonlyrigidbody	570
primexselectrigidbody	571
primexsetoccupancy	571
primexsettings	571
primexswitchaltpositions	574
primextoggleselectrigidbody	575
primexunitcell	575
primexunselectdensitypeak	575
primexunselectligandsolvent	575
primexwriteaddwaters	576
primexwritecalculatorfactors	576
primexwritecreatemap	576
primexwriteexport	576
primexwriteplaceligand	577
primexwriterealspace	577
primexwritereciprocal	577
profile	577
projectbackup	578
projectbackuppreference	578

projectcheck	581
projectclose	581
projectcopy	582
projectdelete	583
projectmerge	584
projectnew	584
projectopen	585
projectprefer	585
projectpublish	587
projectrename	587
projectrestore	588
projectrevertopen	588
projectrevertsnapshot	588
projectrevertstate	589
projectsaveas	589
projectsceneclose	589
projectscenedelete	589
projectscenedescription	590
projectscenedragrows	590
projectsceneexportpresentation	590
projectsceneexportproject	591
projectsceneextendselect	591
projectscenenew	591
projectsceneoptions	592
projectscenereplace	593
projectsceneselect	593
projectscenetoggleselect	593
projectsceneview	593
projectsceneviewnamed	594
projectstoresnapshot	594
projectstorestate	594
projectsynchronize	594
projecttablefind	595
projecttablereplace	596
projecttablereplaceall	597
projectupdatecoordinates	597
projectupdateviews	597
propertycalculate	597
propertycalculatechirality	599
propertycalculatemolformula	600
propertycalculatepickpka	601
propertycalculatesubstructs	602
propertyclearvalue	602
propertycreate	603
propertydelete	604

propertygeneratecontacts	604
propertygeneratehbond	605
propertymeasurementsetting	605
propertymove	605
propertyprecision	606
propertyrename	606
propertyshowall	606
propertysuperimposesetting	607
protassign	607
protassignresidues	608
protassignstart	608
protassignwrite	609
proteinsculpting	609
psp_helical	609
psp_pairwise	610
psp_spatial	611
pspalignaddanchor	612
pspaligndeleteanchor	612
pspaligninsertgaps	612
pspalignlockgaps	613
pspalignmoveleft	613
pspalignmoveleftblock	614
pspalignmoveright	614
pspalignstructures	615
pspalignunlockgaps	615
pspbstogglehetatom	615
pspbuiltbackbone	616
pspbuiltconsensusmodel	616
pspbuiltheteromultimer	616
pspbuilthomomultimer	616
pspbuiltstructure	617
pspconstraintmarkersettings	617
pspeditquerysequence	621
pspexcludetable1row	621
pspexcludetablerow	622
pspexportalignment	622
pspfindfamily	622
pspfindhomologs	623
pspfoldrecognitionoptions	623
pspfoldrecognitionsearch	623
pspgetquerysequence	623
pspimportalignment	624
pspimporthomolog	624
pspimportssp	624
pspincludetable1row	625

pspincludetable	625
pspminimizationresidues	625
pspoptimizealignment	626
pspquickbuild	626
pspquickhomomultimer	626
psprefinebackbone	626
psprefinestructure	627
psprsdefaulthelixloops	627
psprshelixlooptogglerefine	627
psprshelixlooptogglerefineonly	627
psprshelixoptions	628
psprslloopoptions	629
psprsnontemplatehelixloops	630
psprunalign	631
pspruncreate	631
psprundelete	631
psprunopen	631
psprunquickalign	632
psprunrename	632
psprunsaveas	632
pspsecstructprediction	633
pspselectextendtable1row	633
pspselectextendtable	633
pspselecthelixlooprow	634
pspselectonlytable1row	634
pspselectonlytable	634
pspselectrscontext	635
pspselectrsrefinement	635
pspselecttable1row	635
pspselecttable	636
pspsequenceaddfile	636
pspsequenceaddworkspace	637
pspsequencecrop	637
pspsequenceselect	637
pspsequenceviewerexport	637
pspsethelixloopresidues	638
pspsidechainresidues	638
pspsortbbtable	639
pspsortfindhomologstable	639
pspsortfoldtable	640
pspsortrbtable	640
pspsortrbtable1	640
pspsortrstable	641
pspsppdelete	641
pspsppexport	641

pspssprevert	642
pspsspset	642
pspstepforward	643
pspstepgoto	643
pspstructureaddentry	644
psptemplatesetregion	644
psptmplsecstructprediction	644
psptranslatespatial	645
pspunselecttablelrow	645
pspunselecttablerow	646
pspupdatehetatomlist	646
pspupdatescores	646
pspwritebuildconsensusmodel	646
pspwritebuildhomomultimer	647
pspwritebuildstructure	647
pyeval	647
pyimp	647
pyrun	647
pythoneval	648
pythonimport	648
pythonrun	648
pythonrunbuiltin	648
qikprop	649
qsarmarkerdump	649
qsarmarkersettings	650
qsitebasis	653
qsitehcap	654
qsitehcapmarkersettings	654
qsiteion	657
qsiteresidue	658
qsiteset	659
quicktorsion	665
quit	666
raisepanels	666
read	666
readpotential	666
reagentprep	667
refinestart	669
refinewrite	669
removeatomsfromtable	669
rename	669
renameproperty	670
repall	670
repatom	675
repatombonds	676

repbond	676
repdefault	677
replacefromhold	677
repquick	678
resetcsearch	679
residuenam	679
residuenumber	679
residuerenumber	679
restorepanels	680
retype	680
ribbon	680
ribbondump	692
ringclosure	692
ringflip	692
rotate	693
run	693
saveimage	694
saveimageheight	695
saveimagewidth	695
savelayout	695
savemovie	695
savemovieoptions	696
scanmode	697
scriptlogfile	698
scriptrun	698
searchdbconfg	699
selectatomsforediting	702
selecteditatom	702
selecteditproperty	703
sequenceviewersaveimage	703
set	704
setread	704
setwrite	704
showdockablepanels	705
showfirstselectedentry	705
showhwstereosetup	705
showmarkers	705
showpanel	706
showpanels	706
showpropertiesselectedentries	706
showproperty	707
showtoolbar	707
sleep	707
sorteditatomcolumn	708
specifiedname	708

spotcenter	708
spotcenterpoint	709
startdiagnosis	709
stop	709
stopdiagnosis	709
strikebuildqsar	710
strikedeletemodel	711
strikeexportmodel	711
strikeextendselectdescriptor	712
strikeimportmodel	712
strikeplotmodel	712
strikepredict	713
strikeselectdescriptor	713
strikeselectmodel	713
strikesimilarity	714
striketoggleselectdescriptor	714
structalignatoms	714
structalignatomsalt	715
structalignstart	715
substructure	715
substructurefile	716
substructurefileread	716
substructurefilewrite	717
substructureshell	717
superimpose	718
superimposeatom	719
superimposeset	719
superimposesmarts	720
surfaceactivegrid	720
surfacecontours	721
surfacedarkencolor	723
surfacedelete	723
surfacedisplay	724
surfaceduplicate	724
surfaceextended	725
surfaceextendedradiuscontext	726
surfaceextendedradiusset	726
surfacemolecular	727
surfacemolecularcontext	728
surfacemolecularset	729
surfacemolecularset2	729
surfaceoption	729
surfacepbrsp	735
surfacepbrspset	737
surfacereaname	737

surfaceresolution	738
surfacescheme	738
surfacecomment	741
surfaceisovalue	741
surfacesettings	742
surfaceviewasl	744
surfacesplitpair	745
surfacestyle	745
surfaceopacity	746
surfaceopacityback	747
surfacevdw	747
surfacevdwcontext	748
surfacevdwset	749
surfaceviewaslset	749
swapptworkspace	749
symmetrizeworkspace	750
system	750
tablealigncolumn	750
tablecopyprop	751
tablegroupsort	752
tablesizecolumn	753
tablesort	753
tablesortall	754
tablesortfields	754
tablesortgroupfields	756
tableunselectnonsubset	756
targetentry	756
tile	757
tileapplyview	757
tileapplyzoom	757
tilemode	758
tilemove	759
tilereset	759
tileresetall	759
timingsetup	760
timingstart	760
timingstop	761
toggleseditatomboolean	761
toggleseditatom	761
toggleseditproperty	762
torsioncheck	762
torsiongroup	763
trajectoryexport	763
trajectoryexportatoms	763
trajectoryexportmovieoptions	764

trajectoryexportoptions	764
trajectoryplayerdisplayatoms	765
trajectoryplayergoto	765
trajectoryplayergotofirst	765
trajectoryplayergotolast	766
trajectoryplayersettings	766
trajectoryplayerstepahead	769
trajectoryplayerstepback	769
trajectoryplayersuperimposeatoms	770
trajectoryresolvemeasurements	770
transform	770
translate	772
undisplayatom	773
undo	773
ungroupentries	773
unhookimport	774
uniquename	774
uniquepdb	774
update	775
updateribbons	775
varymolecule	775
varytorsion	776
vcsaddattachment	776
vcsaddcorefromproject	777
vcsaddmincapcore	777
vcsaddoriginalcore	778
vcscanceldockjob	778
vcsclearreagentfile	778
vcscombiexportdockingfile	778
vcscombiexportdockingproject	779
vcscombiexportoptions	779
vcsconfiguredocking	779
vcscoreoptions	782
vcscreatedockedlibrary	784
vcsdeleteattachment	784
vcsdeletecore	784
vcsdeleteresults	784
vcdisplayreceptor	784
vcsdocking	785
vcsenumeratedockoptions	785
vcsexcludetablerow	786
vcsexportdefinition	786
vcsexportresults	786
vcsexportspreadsheet	787
vcimportdefinition	787

vcsimportrun	788
vcsincludeextendtablerow	788
vcsincludeonlytablerow	788
vcsincludetablerow	788
vcsinverttableselection	789
vcsoptions	789
vcsrefreshstructure	790
vcsrenameattachment	790
vcsrestorerresults	790
vcsruncombinationdocking	791
vcsruncombinatorialdocking	791
vcsruncombinatorialselection	792
vcsruncreate	792
vcsrundelete	792
vcsrunenumerateddocking	792
vcsrunopen	793
vcsrunrename	793
vcsrunsas	793
vcsrunsingledocking	794
vcsrunsinglselection	794
vcssaveresults	794
vcssselectalltablerows	795
vcssselectextendtablerow	795
vcssselectonlytablerow	795
vcssselecttablerow	796
vcsssetattachmentfile	796
vcsssetmolecule	796
vcssetreagentfile	797
vcssingleexportdockingfile	797
vcssingleexportdockingproject	797
vcssingleexportoptions	798
vcssorttable	798
vcsstepforward	799
vcsstepgoto	799
vcsundisplayreceptor	799
vcsunselecttablerow	799
vcswritecombinationdocking	800
vcswritecombinatorialdocking	800
vcswriteenumerateddocking	800
vcswritesingledocking	800
viewcrystalmates	801
viewdelete	802
viewdeleteselectd	802
viewdragselection	802
viewexport	803

viewextendselect	803
viewextendselectrow	803
viewimport	804
viewmatrix	804
viewplayergotofirst	804
viewplayergotolast	805
viewplayerplaybackward	805
viewplayerplayforward	805
viewplayersettings	805
viewplayerstepahead	806
viewplayerstepback	806
viewplayerstop	806
viewrecordmovie	806
viewrename	807
viewreset	807
viewresetnofit	807
viewrestore	808
viewrestorerow	808
viewsave	808
viewselect	809
viewselectonly	809
viewselectonlyrow	809
viewselectrow	810
viewsetpause	810
viewsetplay	810
viewsetrowpause	811
viewtrajectory	811
viewunselect	811
viewunselectrow	812
viewunsetplay	812
viewvolume	812
visexport	813
visimport	814
volumedelete	815
volumedisplay	815
volumeisosurface	816
volumerename	817
volumesetcolors	817
wake	818
workspacecopy	818
workspacefontdecrease	818
workspacefontincrease	818
workspacepaste	819
workspacepastesettings	819
workspacescopeoff	820

workspacescopeon	820
workspacescopesettings	820
workspacescopeview	821
workspaceselectionadd	822
workspaceselectionclear	823
workspaceselectioninvert	823
workspaceselectionreplace	823
workspaceselectionsubtract	824
workspacetile	824
write	824
writeangle	824
writecontact	825
writecoupling	826
writedihehdral	827
writedistance	828
writehbond	829
wsassistantapply	830
wsassistantapplycontacthbond	830
wsassistantapplyfocus	830
wsassistantbindingsurface	830
wsassistantcontacthbond	832
wsassistantcreatesurfaces	832
wsassistantdefaultstyle	833
wsassistantfixreceptor	833
wsassistanthideallsurfaces	833
wsassistantliganddetection	833
wsassistantmaterial	835
wsassistantprimarystyle	835
wsassistantrendering	837
wsassistantshowallsurfaces	837
wsassistantwsstyle	838
xcluster	845
xclusterstart	847
xclusterwrite	847
zoom	847

Command Index 849

